



## Technique Spotlight

### Intra-operative Pancreatoscopy for Patient Undergoing Whipple Procedure

Presented by:  
Dr. Steven Trocha  
Surgical Oncologist  
Prisma Health-Greenville Memorial Hospital  
Greenville, NC



#### Patient Presentation:

---

70 year old caucasian female presented with jaundice. She had a gastric bypass less than a year old. Patient had an EDGE procedure (EUS-directed transgastric ERCP) where cytology brushings were positive for Pancreatic Adenocarcinoma. PreOperative imaging included a pancreas protocol CT which demonstrated a resectable pancreas cancer, double duct sign and no evidence of regional or metastatic disease.

#### Procedure:

---




This case was unusual as the patient had a Roux En Y gastric bypass which necessitated a subtotal gastrectomy with a standard whipple. In addition, the patient had a replaced right hepatic artery. After standard mobilization, the neck of the pancreas was divided with Bovie cautery.



The dilated duct proximating distally was visualized, the sterile SpyGlass Discover Digital Catheter was used to cannulate the proximal and distal ducts. (See video). No intraluminal lesions were identified. The standard anastomoses were then completed.

## Video

---



## Conclusion:

---

Pancreatotomy added less than 5 minutes to the 5 hour procedure and provided confidence there was no intraluminal extension of tumor. For patients undergoing a whipple procedure, confirming no intraluminal lesions can be of great value to any surgeon. For patients presenting with Roux En Y gastric bypass with positive stones on IOC, a transcystic approach to complete LCBDE may be used to eliminate the need of an ERCP. This was a useful adjunct to ductal evaluation in this whipple case.

**IMPORTANT INFORMATION:** These materials are intended to describe common clinical considerations and procedural steps for the use of referenced technologies but may not be appropriate for every patient or case. Decisions surrounding patient care depend on the physician's professional judgment in consideration of all available information for the individual case. Boston Scientific (BSC) does not promote or encourage the use of its devices outside their approved labeling. Case studies are not necessarily representative of clinical outcomes in all cases as individual results may vary.

This educational video was produced in cooperation with Dr. Trocha. Results from case studies are not predictive of results in other cases. Results in other cases may vary.

**CAUTION:** The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labelling supplied with each device. Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. This material not intended for use in France. Rx Only.

All trademarks are the property of their respective owners.

**Boston  
Scientific**  
Advancing science for life™

Boston Scientific Corporation  
300 Boston Scientific Way  
Marlborough, MA 01752-1234  
www.bostonscientific.com

©2020 Boston Scientific Corporation  
or its affiliates. All rights reserved.

ENDO-888407-AA