



## Ampullectomy with EXALT™ Model D Single-Use Duodenoscope

**Bhavesh Shah, M.D., FASGE**  
Chief of Endoscopy,  
Director of Advanced Endoscopy,  
MetroHealth Medical Center, Cleveland, Ohio



### Patient History & Assessment


A 78-year-old woman underwent upper endoscopy for reflux symptoms. A polypoid lesion was visualized at the ampulla. This was biopsied with a cold forceps. Pathology revealed tubular adenoma. The patient was referred for endoscopic ultrasound and consideration for subsequent ampullectomy.

### Procedure


Endoscopy was first performed with a duodenoscope demonstrating the previously visualized ampullary lesion adjacent to a duodenal diverticulum. Endoscopic ultrasound revealed an 11.7mm x 9.2mm ampullary lesion consistent with the patient's adenoma. No other abnormalities were identified.

ERCP with endoscopic ampullectomy utilizing EXALT™ Model D was then performed. The pancreatic duct was cannulated and an 0.025 inch guidewire was advanced into the main pancreatic duct. An endoscopic snare was placed over the guidewire and resection of the ampulla was performed utilizing cautery. A 5Fr x 5cm pancreatic stent was placed 4cm into the pancreatic duct. The remainder of the ampulla was removed with a second snare resection.


The ampullary sample was retrieved with a RescueNet™ Retrieval Net and sent for pathologic analysis. The common bile duct was then cannulated and a 10Fr x 7cm plastic stent was placed 6cm into the bile duct. Pathology revealed tubular adenoma of the duodenum with high-grade dysplasia. The margins were free of adenomatous change/dysplasia. The patient was discharged home without complication.




*Figure 1*  
Initial Image of the Ampulla  
Prior to Instrumentation




*Figure 2*  
Passing the Snare Over the Wire  
and Closing Down on the Ampulla



*Figure 3*  
Initial Resection of the Ampulla  
over the Wire



*Figure 4*  
Immediate Post Resection with  
Ampulla Inside Snare



*Figure 5*  
Placement of Pancreatic Stent  
Prior to Resection

## Case Outcome/Discussion

The patient tolerated the procedure well without complication and was discharged the same day. EXALT™ Model D was used for multiple different interventions during this procedure, each intervention with different positions and maneuvers and all completed successfully.

The scope reacted well to adjustments in positioning and torque and was stable during the procedure. EXALT Model D provided additional stiffness which assisted in this procedure where multiple interventions were performed sequentially.

Other devices used during the procedure included: Autotome™ 39 Cannulating Sphincterotome, Jagwire™ Revolution High Performance Guidewire, Advanix™ Pancreatic Stent, Captivator™ II Single-Use Snare.

**IMPORTANT INFORMATION:** These materials are intended to describe common clinical considerations and procedural steps for the use of referenced technologies but may not be appropriate for every patient or case. Decisions surrounding patient care depend on the physician's professional judgment in consideration of all available information for the individual case. Boston Scientific (BSC) does not promote or encourage the use of its devices outside their approved labeling. Case studies are not necessarily representative of clinical outcomes in all cases as individual results may vary.

Images provided courtesy of Bhavesh Shah, M.D.

**CAUTION:** The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings, and instructions for use can be found in the product labeling supplied with each device or at <http://www.ifu-bsci.com>. Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. This material not intended for use in France. Rx Only.

2022 Copyright © Boston Scientific Corporation or its affiliates. All rights reserved.

ENDO-1431505-AA