

□ □ ■ □ □ □ □ □ □ □ □ □ □ □ □ □ □ ■

# WIRE-GUIDED CANNULATION TECHNIQUE GUIDE

□ □ □ □ □ □ □ □ □ □ □ ■ □ □ □ □ □ □ □ □

*As described by*

**Steven A. Edmundowicz,  
MD, FASGE**

**Thomas E. Kowalski, MD**

**Peter D. Stevens, MD**

Dear Colleagues:

Cannulation is the initial step in all diagnostic and therapeutic Endoscopic Retrograde Cholangiography and Pancreatography (ERCP). It can be a major hurdle for individuals performing this procedure. While reported success rates for cannulation are very high, even experts encounter difficult and sometimes unsuccessful cannulation attempts.

Technologic advancements have modified the tools available to achieve cannulation. Recent developments include systems that allow for a physician controlled cannulation technique. With this technique, Physician Controlled Wire-Guided Cannulation (PCWGC), the physician controls the guidewire and is provided with direct tactile feedback from the wire during cannulation. This affords the operator additional control over the cannulation process that was not possible with previous cannulation techniques, during which the assistant manipulated the guidewire. Additionally, developments in wire technology have allowed the production of coated nitinol core wires with hydrophilic tips that are designed to provide excellent push-ability, precise control, and a less atraumatic “lumen seeking” tip. It is our belief that these features limit ampullary trauma and inadvertent pancreatic duct injections, factors that in theory could reduce the incidence of post procedure pancreatitis.<sup>2,3,7</sup>

We feel PCWGC is an important skill that may assist the physician in obtaining cannulation of the desired duct at ERCP. During discussions of this technique at meetings and live courses we have learned there are different PCWGC approaches. This guide is an effort to communicate and share with you the PCWGC approaches we have found most helpful in our practices. It is our hope that by carefully describing the process with step-by-step instructions and tips, you will find it easier to develop your personalised technique for PCWGC to adopt in your practice. Whether you elect to use Physician Controlled Wire-Guided Cannulation in every case or as an option if your initial cannulation technique fails, we have created this guide to assist you and hopefully benefit both you and your patients.

Sincerely,

Steven A. Edmundowicz, MD, FASGE

Thomas E. Kowalski, MD

Peter D. Stevens, MD

---



## TABLE OF CONTENTS

---

Device Selection .....	2
Procedural Technique .....	3
Alternative Procedural Methods .....	9
Indications & Contraindications.....	14

**Steven A. Edmundowicz, MD, FASGE**

*Professor of Medicine  
Chief of Endoscopy  
Washington University School of Medicine  
St. Louis, Missouri*

**Thomas E. Kowalski, MD**

*Associate Professor of Medicine  
Director, Gastrointestinal Endoscopy  
Thomas Jefferson University  
Philadelphia, Pennsylvania*

**Peter D. Stevens, MD**

*Assistant Professor of Clinical Medicine  
Columbia University College of Physicians  
and Surgeons  
Director of Endoscopy  
Director of Interventional Endoscopy  
Columbia University Medical Centre  
New York Presbyterian Hospital  
New York, New York*



# PROCEDURAL TECHNIQUE

1. A sphincterotome pre-loaded with a guidewire is prepared and passed through the endoscope and below the major papilla.<sup>A</sup> (Figures 1, 2)
2. The guidewire is then “peeled” (removed from inside the cannula) from the insertion point, to a point on the sphincterotome just above the operating channel port (to allow manipulation by the physician.) (Figure 3)




Figure 1 – Sphincterotome, preloaded with guidewire, is inserted into scope.




Figure 2 – Sphincterotome and guidewire are passed into the major papilla.




Figure 3 – Guidewire is “peeled” with forefinger from sphincterotome (other hand holds sphincterotome and wire at device head).

**A** ORIENTATION – Orientation to the major papilla is a key component of all cannulation attempts. Wire-guided cannulation requires the operator to place the sphincterotome in the optimum cannulation position for cannulation but uses the wire instead of contrast or the tip of the sphincterotome to actually engage the bile or pancreatic duct.



### Technique 1: Sphincterotome Outside the Ampullary Orifice

The guidewire is advanced gently into the papillary orifice and, under fluoroscopic guidance, manipulated into the bile duct. Once the wire advances beyond the papillary orifice, fluoroscopy can be activated and its position in the desired duct confirmed. If prompt and resistance-free passage into the bile duct does not occur, the guidewire is withdrawn and the angle of approach is changed—by manipulating the flex on the sphincterotome or changing the elevator or scope position.<sup>C, D</sup> (Figures 6, 7, 8)

*Selective cannulation with the wire.*



Figure 6 – Cannulation of CBD with sphincterotome outside ampullary orifice.



Figure 7 – Side view




Figure 8 – En face view

**C** APPROACH TO THE AMPULLARY ORIFICE (Technique 1) – When using this technique, the guidewire is advanced out of the sphincterotome so that it extends 2-5 mm beyond the sphincterotome tip toward the papillary orifice. In small motions, the wire is gently moved back and forth to help the operator further realise the vector of the wire in relationship to the major papilla and underlying duct systems. Small adjustments to the scope position and the sphincterotome tip position help to determine the desired angle of entry. The operator analyses the surface of the papilla to the point of entry. At times an orifice is clearly visible, but in other situations it is more difficult to determine. In this case, the surface of the papilla may provide clues for several putative openings. The tip of the wire is used to locate the orifice by gently touching the papilla, feeling and watching for the wire to slip into an orifice. Multiple gentle attempts to enter putative openings are made in succession as the operator changes the point of wire contact on the papilla to precisely target a specific point or to change to a new target. If prompt and resistance free passage into the bile duct does not occur, the guidewire is withdrawn and the angle of approach is changed by manipulating the flex on the sphincterotome or changing the elevator and/or scope position.

**D** WIRE AS CANNULATION TOOL (Technique 1) – Commonly, the wire tip slips into the common channel at an angle that is not suited for deep cannulation of the common bile duct. Multiple, subtle changes in direction of the wire are required to establish deep cannulation. Using the wire as the cannulation tool in this way allows for additional methods of altering the direction of the wire's tip. The first technique is to use the ampullary orifice as a fulcrum, and to “see-saw” the tip of the wire into a new direction. A second useful technique is to use the wall of the duct to carom the tip of the wire into a new vector. If a gentle bow forms in the wire, the tip of the cannula can be used to rotate the bow and new angles are discovered. At times, more stability is required to for deep cannulation of the CBD, so the endoscopist may “follow” the wire with the sphincterotome into the papilla. Once the tip of the sphincterotome engages the papilla approximately 1-2 mm, the sphincterotome is used as a steer-able conduit for the wire and the cannulation process is then identical to technique 2.

When using these wire cannulation techniques, it is essential to be gentle to avoid trauma to the papilla. One welcomed feature when keeping a short length of wire outside of the papilla, with technique 1, is that the wire itself can absorb excess forward pressure by bowing. Once no “free wire” is visible between the tome and the papilla, bowing is constrained by the duct walls and more force may be applied to the tip of the wire. CAUTION: If only 1 mm or less of the wire is out of the catheter the wire may act like a needle and pierce the epithelium, causing a false tract. This may be prevented by carefully monitoring the force on the wire using the excellent tactile feedback provided by the wire itself.



## Wire Positioning Tips

### USING THE SCOPE FOR POSITIONING




Figure 12 – Scope is positioned for cannulation just below the papilla. Scope is put in “shortened position” and may be pushed inward and outward to facilitate sphincterotome positioning.




Figure 13 – Big wheel is rolled in “up” direction (counter-clockwise) to move scope toward ampulla or “down” (clockwise) to move scope away from ampulla. Small wheel facilitates movement to the left or right of papilla.

### USING THE SPHINCTEROTOME FOR POSITIONING

*Sphincterotome flex may be used to adjust the position of the sphincterotome.*



Figure 14 – Side view (technique 1)




Figure 15 – En face view (technique 1)



Figure 16 – Assistant flexes sphincterotome tip.



## ALTERNATIVE METHODS

*In some cases, technical and anatomical challenges will present or the physician may prefer to have an assistant manipulate the wire. If either of these situations arises, there are two alternative wire-guided cannulation methods to consider.*

## DOUBLE WIRE METHOD

### Alternative method for technical or anatomical challenges

If repeated attempts at cannulation result in continued wire placement in the pancreatic duct, a double wire technique may be used to facilitate bile duct cannulation. With this method, the guidewire is left in the pancreatic duct while a second guidewire is passed to cannulate the CBD. This technique may be especially useful before proceeding with more aggressive cannulation measures such as pancreatic duct stent placement or pre-cut sphincterotomy.

1. If this method is used, the pancreatic wire is advanced deep into the main pancreatic duct, if possible, so that the stiff portion of the wire is within the ampullary segment. Once the wire is advanced to the tail of the pancreas, it is locked in place and the sphincterotome is removed from the duodenoscope. <sup>F</sup> (Figures 19, 20)



Figure 19 – Guidewire is placed in Pancreatic Duct (PD).




Figure 20 – Guidewire placed in PD is locked into place with locking device.




Figure 21 – Alpha loop in PD to facilitate passage.

**F** PANCREATIC WIRE ADVANCEMENT – When advancing the wire in the pancreatic duct care must be taken not to advance a wire in a straight conformation against resistance so as to avoid potential side branch perforation. If there is no resistance, the wire may be advanced to the tail of the pancreas with limited concerns of pancreatic injury. If the tip of the wire enters a side branch, continued gentle advancement may allow the wire to find the path of least resistance by folding on itself and adopting an alpha (or “safety”) loop conformation. In the alpha loop conformation the wire is then advanced to the tail, typically without engaging any additional side branches. (Figure 21)



4. The tip of the sphincterotome is manipulated to direct the guidewire slightly above and coming from the right to the left of the pancreatic wire. Using one of the two wire-guided cannulation techniques described in the first section of this guide, the second wire is used to gain access to the bile duct. <sup>G</sup> (Figures 24, 25, 26, 27)
  
5. Upon achievement of CBD cannulation, the second wire is locked into place using the external locking device, and the procedure is continued. Depending on the clinical situation and physician preference, the first wire may be removed, used for additional therapy, or locked into place until the end of the procedure for prophylactic pancreatic stent placement.




Figure 24 – Side view (technique 2)

*Guidewire is positioned to perceived CBD orientation. Initial guidewire is still placed in PD.*




Figure 25 – En face view (technique 2)




Figure 26 – Guidewires placed in CBD and PD. (technique 2)



Figure 27 – Fluoroscopy image of guidewire placement in CBD and PD. (technique 2)

**G PD WIRE INTENDED BENEFIT** – The wire in the pancreas lengthens and stabilises the papilla facilitating biliary cannulation. If the intraduodenal biliary segment has a “sigmoid” morphology the pancreatic wire straightens the “S-turn” again facilitating biliary cannulation.



## Resource List

---

To learn more about wire-guided cannulation, below is a list of published literature on this topic for your reference.

1. Freeman ML, DiSario JA, Nelson DB et al. Risk factors for post-ERCP pancreatitis: a prospective, multicenter study. *Gastrointestinal Endoscopy* 2001; 54:425-434.
2. Johnson KG, Geenen JE, Johanson JF et al. Evaluation of post-ERCP pancreatitis: potential causes noted during a controlled study of differing contrast media. *Gastrointestinal Endoscopy* 1997; 46:217-222.
3. Freeman ML and Guda NM. Prevention of post-ERCP pancreatitis: a comprehensive review. *Gastrointestinal Endoscopy* 2004; 59:845-864
4. Maede S, Hayashi H, Hosokawa O, et al., Prospective randomised pilot trial of selective biliary cannulation using pancreatic guide-wire placement. *Endoscopy* 2003; 35:721-724.
5. Lella F, Bagnolo F, Columbo E, and Bonassi U. A simple way of avoiding post ERCP pancreatitis. *GI Endoscopy* 2004; 59:830-834.
6. Gyokeres T, Duhl J, Varsanyi M, et al., Double guidewire placement for endoscopic pancreaticobiliary procedures. *Endoscopy* 2003; 35:95-96.
7. Freeman ML, Nalini M, and Guda NM. ERCP Cannulation: A Review of Reported Techniques. *Gastrointestinal Endoscopy* 2005; 61:112-125

This educational technique guide was developed and produced in cooperation with Steven A. Edmundowicz, MD, FASGE, Thomas E. Kowalski, MD, and Peter D. Stevens, MD. The opinions, recommendations and techniques reflected in this guide are those of these physicians and do not necessarily reflect the opinions and recommendations of Boston Scientific Corporation, its employees or its affiliates.

\* Clinical images courtesy of Peter D. Stevens, MD.

**European Headquarters – Paris**

T: +33 1 57 66 80 00 F: +33 1 57 66 84 99

**Argentina**

T: +5411 4896 8556 F: +5411 4896 8550

**Australia/New Zealand**

T: +61 2 8336 5555 F: +61 2 8335 0404

**Austria**

T: +43 1 726 30 05 1022 F: +43 1 726 30 05 1030

**Balkans**

T: +30 210 95 37 890 F: +30 210 95 79 836

**Belgium (Freephone)**

T: 080094 494 F: 080093 343

**Brazil**

T: +55 11 5502 8500 F: +55 11 5103 2212

**Bulgaria**

T: +359 2 986 50 48 F: +359 2 986 57 09

**Canada**

T: +1 888 359 9691 F: +1 888 575 7396

**Chile**

T: +562 445 4904 F: +562 445 4915

**China – Beijing**

T: +86 10 8525 1588 F: +86 10 8525 1566

**China – Guangzhou**

T: +86 20 8767 9791 F: +86 20 8767 9789

**China – Shanghai**

T: +86 21 6141 5959 F: +86 21 6141 5900

**Colombia**

T: +57 1 629 5045 F: +57 1 612 4761

**Czech Republic**

T: +420 2 3536 2911 F: +420 2 3536 4334

**Denmark (Freephone)**

T: 80 30 80 02 F: 80 30 80 05

**Finland (Freephone)**

T: 09 622 42 85 F: 09 622 42 86

**France**

T: +33 1 39 30 49 30 F: +33 1 39 30 49 31

**Germany (Freephone)**

T: 0800 072 3303 F: 0800 072 3319

**Greece**

T: +30 210 9542 300 F: +30 210 9542 310

**Hong Kong**

T: +852 2960 7100 F: +852 2563 5276

**Hungary**

T: +36 1 456 30 40 F: +36 1 456 30 41

**India – Bangalore**

T: +91 80 5112 1104/5 F: +91 80 5112 1106

**India – Chennai**

T: +91 44 2648 0318 F: +91 44 2641 4695

**India – Delhi**

T: +91 11 2618 0445/6 F: +91 11 2618 1024

**India – Mumbai**

T: +91 22 5677 8844 F: +91 22 2617 2783

**Italy**

T: +39 010 60 60 1 F: +39 010 60 60 200

**Korea**

T: +82 2 3476 2121 F: +82 2 3476 1776

**Mexico**

T: +52 55 5687 63 90 F: +52 55 5687 62 28

**Middle East / Gulf / North Africa**

T: +961 1 805 282 F: +961 1 805 445

**The Netherlands**

T: +31 43 356 8260 F: +31 43 356 8265

**Norway (Freephone)**

T: 800 104 04 F: 800 101 90

**Poland**

T: +48 22 435 1414 F: +48 22 435 1410

**Portugal**

T: +351 1 381 25 40 F: +351 1 381 28 29

**South Africa**

T: +27 11 840 8600 F: +27 11 463 6077

**South East Asia – Malaysia**

T: +60 3 2284 1499 F: +60 3 2284 1507

**South East Asia – Philippines**

T: +63 2 687 3239 F: +63 2 687 3047

**South East Asia – Singapore**

T: +65 6418 8888 F: +65 6418 8899

**South East Asia – Thailand**

T: +66 2 2654 3810 F: +66 2 2654 3818

**Spain**

T: +34 901 11 12 15 F: +34 93 405 9045

**Sweden (Freephone)**

T: 020 65 25 30 F: 020 55 25 35

**Switzerland (Freephone)**

T: 0800 826 786 F: 0800 826 787

**Taiwan**

T: +886 2 2747 7278 F: +886 2 2747 7270

**Turkey – Istanbul**

T: +90 216 464 36 66 F: +90 216 464 36 67

**Uruguay**

T: +59 82 900 6212 F: +59 82 900 6212

**UK & Eire**

T: +44 1727 866633 F: +44 1727 830371

**Venezuela**

T: +58 212 959 8106 F: +58 212 959 5328

**Boston  
Scientific**

*Delivering what's next.™*

[www.bostonscientific.com/international](http://www.bostonscientific.com/international)

[www.bostonscientific.ca](http://www.bostonscientific.ca)

[www.bostonscientific.co.uk](http://www.bostonscientific.co.uk)

[www.bostonscientific.de](http://www.bostonscientific.de)

[www.bostonscientific.es](http://www.bostonscientific.es)

[www.bostonscientific.fr](http://www.bostonscientific.fr)

[www.bostonscientific.ie](http://www.bostonscientific.ie)

[www.bostonscientific.it](http://www.bostonscientific.it)

[www.bostonscientific.nl](http://www.bostonscientific.nl)

RCS Nanterre B420 668 420

© 2007 Boston Scientific Corporation

or its affiliates. All rights reserved.

DINEND2178EA

All cited trademarks are the property of their respective owners. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labelling supplied with each device. Information for the use only in countries with applicable health authority product registrations.

PSST 4318 Printed in the UK by Gosling.