



# TheraSphere™ Y-90 Glass Microspheres | RASER STUDY

## Radiation segmentectomy for curative intent of unresectable very early to early stage hepatocellular carcinoma (RASER): a single-center, single-arm study

Kim E, Sher A, Abboud G, et al. Radiation segmentectomy for curative intent of unresectable very early to early stage hepatocellular carcinoma (RASER): a single-centre, single-arm study [published online ahead of print, 2022 May 23]. *Lancet Gastroenterol Hepatol.* 2022;S2468-1253(22)00091-7. doi:10.1016/S2468-1253(22)00091-7

First prospective study to assess outcomes after radiation segmentectomy (RS) in patients with very early or early stage hepatocellular carcinoma which showed 100% initial objective response with TheraSphere Y-90 glass microspheres.

### OVERVIEW

The aim of this study was to assess the safety and efficacy of RS in patients with unresectable hepatocellular carcinoma (HCC) deemed unfavorable for ablation. The study provides a strong rationale for new randomized trials comparing RS to ablation and supports inclusion of RS in BCLC guidelines.

### OBJECTIVE

**Primary endpoint:** Target tumor response measured by modified RECIST (mRECIST).

**Secondary endpoints:** Time to progression (TTP) of the target lesion and overall disease, and adverse events using CTCAE version 5.0.

CTCAE = National Cancer Institute Common Terminology Criteria for Adverse Events

### STUDY DESIGN/METHODS

The study included adults (>18 years) with solitary HCC with unfavorable location for ablation, without metastasis or macrovascular invasion. Eligibility criteria included measurable disease ≤ 3 cm in diameter, Child-Pugh score A–B7, an ECOG score of 0, and adequate hematological and organ function. Of the 44 individuals assessed for eligibility, 29 patients were included in the study. Patients were followed up with imaging and office visits for up to 24 months.

## KEY RESULTS

#### PARTICIPANT DEMOGRAPHICS

ECOG	0
Sex	Male: 23/29 (79%) Female: 6/29 (21%)
Child-Pugh	A5: 14/29 (48%) A6: 12/29 (41%) B7: 3/29 (10%)
Mean Perfused Liver Volume	153.6 mL (mean SD 99.2)
Median Tumoral Dose Delivered	1004.6 Gy (95% CI [190.8–3730.0])


#### TARGET LESION RESPONSE; % OF PARTICIPANTS (n=29)

Initial Objective Response	Initial Complete Response	Partial Response	Sustained Complete Response
100% (29)	83% (24)	17% (5)	90% (26) Median Time: 43 days Median Duration: 635 days

#### TARGET LESION & OVERALL DISEASE PROGRESSION; % OF PARTICIPANTS (n=29)

Target Lesion Progression	Cumulative Incidence of Target Lesion Progression	Actuarial Overall Survival	Overall Progression	Cumulative Incidence of Overall Progression
10% (3)	Year 1: 4% Year 2: 12%	Year 1: 96% Year 2: 96%	31% (9)	Year 1: 14% Year 2: 27%

#### LIVER TRANSPLANT NECROSIS



#### LOW PROPORTION OF HIGH-GRADE ADVERSE EVENTS\* (AEs); % OF PARTICIPANTS (n=29)

Grade 3 Leukopenia	14% (4)
Grade 3 Thrombocytopenia	7% (2)
Grade 3 Non-laboratory-related adverse events	7% (2)

\*A complete list of adverse events can be found in Table 2 in the publication

Eight patients received a liver transplant. Pathology results show all eight target lesions had 100% necrosis.

### CONCLUSION

In this study, radiation segmentectomy with TheraSphere Y-90 glass microspheres was shown to be effective, with a low proportion of high-grade adverse events in patients with unresectable very early to early-stage hepatocellular carcinoma with suboptimal location for ablation.

Sustained complete response rates and local progression of the target lesion were similar to the previously reported rates after thermal ablation.\*\*

Given complete pathological necrosis of the explanted tumors, larger investigative studies on the curative potential of radiation segmentectomy are warranted.

\*\*References results cited in the publication.

# TheraSphere™ Y-90 Glass Microspheres | RASER STUDY

## TheraSphere™ Yttrium-90 Glass Microspheres

**INDICATION FOR USE:** TheraSphere is indicated for use as selective internal radiation therapy (SIRT) for local tumor control of solitary tumors (1-8 cm in diameter), in patients with unresectable hepatocellular carcinoma (HCC), Child-Pugh Score A cirrhosis, well-compensated liver function, no macrovascular invasion, and good performance status. **CONTRAINDICATIONS:** TheraSphere is contraindicated in patients: whose Tc-99m macroaggregated albumin (MAA) hepatic arterial perfusion scintigraphy shows any deposition to the gastrointestinal tract that may not be corrected by angiographic techniques • who show shunting of blood to the lungs that could result in delivery of greater than 16.5 mCi (0.61 GBq) of Y-90 to the lungs. Radiation pneumonitis has been seen rarely in patients receiving doses to the lungs greater than 30 Gy in a single treatment • in whom hepatic artery catheterization is contraindicated, such as patients with vascular abnormalities or bleeding diathesis • who have pulmonary insufficiency (conventionally defined by an arterial oxygen pressure (PaO<sub>2</sub>) of < 60 mmHg or oxygen saturation (SaO<sub>2</sub>) of < 90%) or severe liver dysfunction, including hepatic encephalopathy, clinically evident ascites or treatment with diuretics for ascites • with portal vein thrombosis (PVT) Type 4 involvement and lack of Tc-99m MAA deposition on the PVT seen on the Tc-99m MAA imaging with >70% tumor replacement in the liver • with comorbidities or poor overall health (e.g., ECOG performance status rating ≥ 2) which may make the patient a poor candidate for locoregional radiation treatment. • who are pregnant. **WARNINGS:** The following pre-treatment, high-risk factors (disease characteristics) have been associated with serious adverse events deemed possibly related to use of the device: infiltrative tumor type • tumor nodules too numerous to count • AST or ALT > 5 times ULN • bilirubin > 2 mg/dL • tumor volume > 50% combined with albumin < 3 g/dL. Keep the TheraSphere dose vial upright and stored in its lead pot before and during patient treatment, except as required for radiation measurement. Do not open the dose vial acrylic shield prior to patient treatment. Post-treatment, waste materials require caution to prevent contamination and beta shielding due to residual glass microspheres. **PRECAUTIONS: GENERAL PRECAUTIONS:** As in any intra-arterial procedure, aseptic technique should be practiced, and care should be taken to ensure minimum patient anesthesia exposure extraneous to therapeutic objective. • Consideration of patient comorbidities should be used when determining the type and volume of fluid to infuse via catheter to avoid electrolyte imbalance, fluid shift, and hyperglycemia. • It is important to avoid any aggressive arterial procedure that may lead to arterial spasm that impairs TheraSphere distribution into the perfused liver target volume which may lead to underdosing or non-target deposition of TheraSphere. **PRECAUTION IN PATIENTS WITH IMPAIRED LIVER FUNCTION:** No efficacy or safety data from the LEGACY study are available to support the use of the device in patients with Child-Pugh score B or C cirrhosis. **PRECAUTION IN VULNERABLE PATIENTS:** No effectiveness or safety data are available to support the use of the device in children or breast-feeding women. **ENDOCRINE DISRUPTION, CARINOGENICITY, MUTAGENICITY, TOXICITY TO REPRODUCTION:** Ideally the use of this radioactive device in women of childbearing capability should be performed during the first few (approximately 10) days following the onset of menses. **RADIATION SAFETY:** Radioactive products should be used only by healthcare professionals who are qualified by specific training in the safe use and handling of radionuclides and whose experience and training have been approved by the appropriate government agency authorized to license the use of radionuclides. • As in the use of any radioactive material, ensure minimum radiation exposure to the patient extraneous to the therapeutic objective, and to minimize radiation exposure to workers and others in contact with the patient. **RELEASE AND POST-TREATMENT PRECAUTIONS:** Post-treatment patient care: use universal precautions for body fluid contact. Tare Y-90 may be detectable in blood and urine; handle with gloves and dispose as normal body fluids. The radiation field is expected to be less than 1 mrem/h (10 µSv/h) at 3 ft (1 m) from the patient's abdomen. Supplemental shielding and segregation of the patient are not required to maintain exposure to others below regulated limits. • Release instructions: The patient should follow good hygiene (e.g., proper hand washing). Caregivers, family, and others do not require restrictions on patient contact; however, they can minimize their radiation exposure by avoiding prolonged time (>12 hours per day) within 1 ft (0.3 m) of the patient's abdomen for the first week post therapy. Patients should be advised that radiation emitted from the patient may be detectable at security screening (e.g., international travel). • Special precautions post-administration: If the patient requires hospitalization, surgery, medical assessment or treatment regarding any part of their thorax or abdomen within first 2 weeks of treatment, the patient should advise the hospital and treating physician of the Y-90 TheraSphere implant. The physician should consult their radiation safety staff for handling and disposal of liver tissue. • Special liver tissue handling: Special liver tissue handling may be required for post-treatment surgery, explant, or transplant since the glass microspheres remain permanently implanted in the liver tissue. Disclosure of the treatment will be required if cremation is considered. **POTENTIAL ADVERSE EVENTS:** The use of this product leads to irradiation of both tumorous and normal liver tissue. As a result, patients with compromised liver function may be at greater risk of liver function impairment and hence could experience complications. Clinical side effects usually occur within the first 4 to 6 weeks after treatment. Based on clinical trial data, literature reviews and post-market surveillance, adverse events potentially associated with treatment using Y-90 microspheres, including TheraSphere, may include the following: Allergic reaction • Altered liver function, acute or chronic • Anorexia • Anxiety • Ascites • Bile Duct Injury • Bleeding/hemorrhage • Chills • Rigors • Cholecystitis (inflammatory or infectious) • Colitis • Death • Dehydration • Diarrhea • Dizziness • Dyspnea • Edema (any location) • Electrolyte abnormalities • Elevated BUN/creatinine • Fall • Fatigue • Fever • Gastrointestinal bleeding / hemorrhage • Gastrointestinal ulcer or ulceration • Hepatic encephalopathy • Hepatorenal failure • Hiccups • Hypertension • Hypotension • Infection (any location) • Liver failure, acute or chronic • Lymphopenia • Malaise • Mood alteration • Muscle weakness • Nausea • Neutropenia • Pain (any location) • Pancreatitis • Platelet count abnormalities • Pleural effusion • Portal hypertension • Pre-existing chronic liver disease decompensation • Pulmonary edema • Pulmonary fibrosis • Radiation hepatitis • Radiation induced disease, acute • Radio Embolization Induced Liver Disease (REILD) • Sepsis • Supraventricular arrhythmia • Thrombosis (arterial or venous) • Tumor inflammation (including tumor edema) • Tumor-lysis syndrome • Vomiting • Weight loss. Complications related to the administration procedure itself may include: Allergic reaction: Arterial injury including vessel dissection • Aspiration pneumonia • Bruising/bleeding/hematoma at site • Constipation/abdominal distension • Fatigue • Flushing • Infection • Nausea • Nerve damage. **CAUTION:** Federal (USA) law restricts this device to sale by or on order of a physician. PI-992004-AA. **Note:** Dose to the liver does not exceed 150 Gy. TheraSphere is a registered trademark of Theragnostics Corporation used under license by Biocompatibles UK Ltd. All other trademarks are property of their respective owners.

**Boston Scientific**  
Advancing science for life™

**Peripheral Interventions**  
300 Boston Scientific Way  
Marlborough, MA 01752-1234  
[www.bostonscientific.com](http://www.bostonscientific.com)

To order product or for more information  
contact customer service at 1.888.272.1001.

© 2022 Boston Scientific Corporation  
or its affiliates. All rights reserved.

PI-1373002-AA