



AMS 700™

Inflatable Penile Prosthesis with TENACIO™ Pump

It all starts with a squeeze.

Powering the AMS 700™ Inflatable Penile Prosthesis (IPP), the new TENACIO™ Pump is purpose-built from the inside out with advanced engineering, and a bold design that aims to impact outcomes and streamline the patient experience.¹



Engineered to redefine reliability.


The TENACIO Pump introduces independent valve systems, thoughtfully engineered to maximize reliability across a wider user profile and to address pump stiction.



- **Independent valving** enables pump activation regardless of fast or slow squeeze rates.¹
- **Direct fluid pathways** designed to enable efficient inflation. In bench testing conducted by BSC (n=5), the TENACIO pump was found to have a 27% faster refill rate as compared to the MS Pump™.*¹
- **Independent valve functions** designed to limit mechanical failure.¹
- **Valve design** intended to reduce oscillation and increase tactile feedback.¹
- **Dedicated anti-auto inflate valve** designed to prevent unexpected inflation.¹

Designed to unlock the feel of confidence.

The TENACIO Pump design is intended to help shorten the patient learning curve and increase satisfaction.¹



IDENTIFICATION

The prominent pump bulb grips are intended to facilitate identification and inflation squeeze efficiency.¹




INFLATION

Engineered for reliable and efficient inflation, regardless of pumping style.¹



DEFLATION

High-profile, one-press release button provides tactile feedback designed to help users confirm that deflation has started.¹



Designed for flexibility.

AMS 700 line has three design options – the most variations of any penile implant on the US market. All are available with InhibiZone, clinically proven to reduce revision surgery due to infection, and offer the same innovative silicone construction to fit ED patients' needs.¹

Ordering Information

Product: AMS 700 LGX™ Inflatable Penile Prosthesis	InhibiZone item number
LGX PRECONNECT TENACIO 12CM, IP, IZ	72404455
LGX PRECONNECT TENACIO 12CM, PS, IZ	72404450
LGX PRECONNECT TENACIO 15CM, IP, IZ	72404456
LGX PRECONNECT TENACIO 15CM, PS, IZ	72404451
LGX PRECONNECT TENACIO 18CM, IP, IZ	72404457
LGX PRECONNECT TENACIO 18CM, PS, IZ	72404452
LGX PRECONNECT TENACIO 21CM, IP, IZ	72404458
LGX PRECONNECT TENACIO 21CM, PS, IZ	72404453

Product: AMS 700 CXR Inflatable Penile Prosthesis	InhibiZone item number
CXR PRECONNECT TENACIO 10CM, IP, IZ	72404465
CXR PRECONNECT TENACIO 10CM, PS, IZ	72404460
CXR PRECONNECT TENACIO 12CM, IP, IZ	72404466
CXR PRECONNECT TENACIO 12CM, PS, IZ	72404461
CXR PRECONNECT TENACIO 14CM, IP, IZ	72404467
CXR PRECONNECT TENACIO 14CM, PS, IZ	72404462
CXR PRECONNECT TENACIO 16CM, IP, IZ	72404468
CXR PRECONNECT TENACIO 16CM, PS, IZ	72404463
CXR PRECONNECT TENACIO 18CM, IP, IZ	72404469
CXR PRECONNECT TENACIO 18CM, PS, IZ	72404464

Product: AMS 700 CX Inflatable Penile Prosthesis	InhibiZone item number
CX PRECONNECT TENACIO 12CM, IP, IZ	72404435
CX PRECONNECT TENACIO 12CM, PS, IZ	72404430
CX PRECONNECT TENACIO 15CM, IP, IZ	72404436
CX PRECONNECT TENACIO 15CM, PS, IZ	72404431
CX PRECONNECT TENACIO 18CM, IP, IZ	72404437
CX PRECONNECT TENACIO 18CM, PS, IZ	72404432
CX PRECONNECT TENACIO 21CM, IP, IZ	72404438
CX PRECONNECT TENACIO 21CM, PS, IZ	72404433
CX PRECONNECT TENACIO 24CM, IP, IZ	72404439
CX PRECONNECT TENACIO 24CM, PS, IZ	72404434

Product: Standalone Pump	InhibiZone item number
PUMP, TENACIO, IZ	72404429

References

1. Data on file with Boston Scientific

Caution: U.S. Federal law restricts this device to sale by or on the order of a physician.

Prior to use, please review the Instructions for Use for a complete listing of indications, contraindications, warnings, precautions and potential adverse events.

Indications for Use: The AMS 700™ Series Inflatable Penile Prosthesis product line is intended for use in the treatment of chronic, organic, male erectile dysfunction (impotence).

Contraindications: The AMS 700 Series Inflatable Penile Prostheses are contraindicated in the patients that have active urogenital infections or active skin infections in the region of surgery or (for the AMS 700 prosthesis with InhibiZone™ Antibiotic Surface Treatment) have a known sensitivity or allergy to rifampin, minocycline or other tetracyclines, or patients with lupus erythematosus because minocycline has been reported to aggravate this condition. **Warnings:** Implantation of the device will make latent natural or spontaneous erections, as well as other interventional treatment options, impossible. Men with diabetes, spinal cord injuries or open sores may have an increased risk of infection associated with the implantation of a prosthesis. Failure to evaluate and promptly treat erosion may result in a substantial worsening of the condition leading to infection and loss of tissue. Implantation may result in penile curvature, or scarring. Pre-existing abdominal or penile scarring or contracture may make surgical implantation more complicated or impractical. If a hypersensitivity reaction develops to a device coated with InhibiZone, the penile prosthesis should be removed and the patient treated appropriately. **Precautions:** Migration of the device components can occur if the cylinders are improperly sized, if the pump or the reservoir is not positioned properly, or if the tubing lengths are incorrect.

Potential Adverse Events: May include device malfunction/failure leading to additional surgery, device migration potentially leading to exposure through the tissue, device/tissue erosion, infection, unintended-inflation of the device and pain/soreness. MH-545408-AB

All images are the property of Boston Scientific. All trademarks are the property of their respective owners.

Boston Scientific
Advancing science for life™

Boston Scientific Corporation
300 Boston Scientific Way
Marlborough, MA 01752
www.bostonscientific.com

All images are property of Boston Scientific. All trademarks are the property of their respective owners.

© 2023 Boston Scientific Corporation or its affiliates. All rights reserved.
URO-1680704-AA NOV 2023