



Orca™

Single Use Air/Water and Suction Valves

OrcaPod™

*Single Use Air/Water, Suction and
Biopsy Valve Set with Water Jet Connector*

**Boston
Scientific**
Advancing science for life™

Refer to the device directions for use for complete instructions on device use.

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

Warning

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

Orca Air/Water and Suction Valves

Indications For Use

The Orca Air/Water Valve is intended to be used to control the air/water function on an endoscope during a GI (gastrointestinal) endoscopic procedure.

The Orca Suction Valve is intended to be used to control the suction function on an endoscope during a GI endoscopic procedure.

Contraindications

Contraindications include those specific to any endoscopic procedure.

Adverse Events

None known.



Orca™

Single Use Air/Water and Suction Valves

OrcaPod™

Single Use Air/Water, Suction and
Biopsy Valve Set with Water Jet Connector

**Boston
Scientific**
Advancing science for life™

Seal Single Use Biopsy Valve

Indications For Use

The Seal Single Use Biopsy Valve provides access for endoscopic device passage and exchange, helps maintain insufflation, minimizes leakage of biomaterial from the biopsy port throughout the endoscopic procedure, and provides access for irrigation.

Contraindications

Contraindications include those specific to any endoscopic procedure.

Adverse Events

None known.

Hydra™ Disposable Auxiliary Water Jet Connector

Indications For Use

The Hydra Disposable Auxiliary Water Jet Connector is used in conjunction with Hydra irrigation tubing (tubing and accessories to accommodate various endoscopes and irrigation pumps) and is intended to provide irrigation via sterile water supply during GI endoscopic procedures when used in conjunction with an irrigation pump.

Contraindications

Contraindications include those specific to any endoscopic procedure.

All trademarks are the property of their respective owners.

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.

Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device.

2023 Copyright © Boston Scientific Corporation or its affiliates. All rights reserved.

ENDO-1578313-AA