


Increasing HeartLogic Alert Level Leads to Proactive Check with Intervention and Rapid Resolution of Heart Failure Symptoms

Summary

HeartLogic Heart Failure Index



Age and Gender: 85 y/o male

Patient History: Chronic HFrEF stable at NYHA class II symptoms, paroxysmal atrial flutter/Afib, previous CRT-D implantation, CAD, ischemic cardiomyopathy, LVEF reported at 20% in 2017, CHF symptoms upon hospital admittance and pacemaker dependency

Day 1: CRT-D change out due to ERI indicator, implanted MOMENTUM CRT-D with HeartLogic diagnostic. Atrial flutter at implant with successful pace termination.

- Latitude enabled with HeartLogic tailored for device reporting to EP/device clinicians. HeartLogic reporting tailored to HF cardiologist in practice. Weight scale and BP cuff ordered.
- Post-op/wound check with recurrent atrial flutter and successful pace termination. Sleep incline calibrated. Patient symptomatic.


(A) Day 80: Pro-active device check at remote clinic due to patient symptoms with increasing HeartLogic score (56). Patient back in atrial flutter with successful pace termination.

- Patient referred back to EP for arrhythmia management and possible ablation.
- HeartLogic data evaluated by HF physician with addition of Torsemide 20mg and Metolazone 2.5mg and d/c Furosemide 40mg.




(B) Day 98: HeartLogic index returns to 0.

Trend Graphs



Trend Graphs




Trend Graphs

—●— Daily Value — 3 Day Average

Day 8 Day 38 Day 69 Day 100 Day 130


V Therapy

Day 155




RV Rate during AT/AF

Day 155




Mean Heart Rate

Day 155



% LV Paced

Day 155



Heart Rate Variability (SDANN)

Data Table & Timeline

LATITUDE™ NXT Patient Management – Heart Failure Management Report

Trend Data									
Date	HeartLogic Heart Failure Index	S3 (mG)	S1 (mG)	Thoracic Impedance (Ω)	Respiratory Rate (rpm)	Night Heart Rate (bpm)	Sleep Incline (degrees)	Weight (lbs)	Blood Pressure (mmHg)
Day 98	0	1.61	1.16	52.3	17.8	71	10	177	114/67
Day 97	2	1.76	1.25	52.3	16.8	70	5	178	125/78
Day 96	5	1.63	1.22	52.2	17.8	70	4	179	136/79
Day 95	6	1.45	1.27	54.8	18.3	71	5	179	126/81
Day 94	10	1.72	1.24	50.2	18.0	70	11		127/81
Day 93	13	1.75	1.07	48.4	17.3	70	11		125/78
Day 92	16	1.88	1.14	48.0	18.0	71	7		143/91
Day 91	18	1.77	1.13	46.8	18.8	71	5		134/83
Day 90	18	1.54	1.09	45.3	18.0	70	8		133/78
Day 89	20	1.60	1.10	44.6	18.0	70	8		117/74
Day 88	20	1.61	1.12	44.7	18.0	70	8		126/77
Day 87	20	1.49	1.13	43.4	19.3	71	11	145	125/72
Day 86	20	1.54	1.33	40.8	N/R	71	7		
Day 85	20	1.85	1.66	42.9	N/R	70	25		147/91
Day 84	31	1.94	1.36	40.6	18.5	73	22		
Day 83	46	1.28	1.75	38.3	19.0	80	4		
Day 82	61	0.93	1.55	37.5	19.3	80	7		116/71
Day 81	61	N/R	N/R	N/R	19.3	96	11		120/82
Day 80	56	2.25	1.17	38.9	20.8	93	14		120/77
Day 79	56	2.23	1.20	38.3	20.3	94	8		105/66
Day 78	51	2.57	1.30	40.2	19.3	95	9		111/65
Day 77	50	3.07	1.65	43.8	19.0	94	9		
Day 76	47	2.13	1.26	41.4	19.0	92	18		110/85
Day 75	47	2.54	1.44	42.8	18.3	93	16		
Day 74	46	2.54	1.48	43.2	18.5	96	15		112/87
Day 73	46	2.36	1.33	41.4	18.0	95	8	197	
Day 72	45	2.32	1.54	42.8	19.0	88	17		104/73
Day 71	45	2.15	1.15	42.2	19.0	71	19		
Day 70	43	2.30	1.17	42.8	18.3	91	12		108/78
Day 69	41	2.66	1.49	44.4	18.8	93	13		
Day 68	33	2.46	1.24	43.6	18.0	91	12		114/79
Day 67	28	2.33	1.52	45.2	19.0	91	15	192	122/78
Day 66	21	2.04	1.52	44.7	19.0	93	11		
Day 65	16	2.83	1.34	44.3	18.3	92	17		
Day 66					19.0				
Day 65	16	2.83	1.34	44.3	18.3	92	17		
Day 64	15	2.07	1.70	45.0	19.0	93	26		
Day 63	30	1.90	1.83	44.4	18.5	95	8		
Day 62	42	1.89	1.74	44.6	18.5	95	10		
Day 61	45	1.40	1.48	44.8	18.3	94	12		
Day 60	45	2.11	1.43	45.0	18.5	91	23		
Day 59	44	2.23	1.44	45.8	18.3	88	27		
Day 58	43	2.47	1.45	46.7	18.3	91	24		
Day 57	40	2.66	1.33	45.6	18.3	89	20		
Day 56	40	2.28	1.56	46.2	18.0	88	11		
Day 55	38	2.46	1.44	45.7	18.0	89	33		
Day 54	38	2.23	1.47	45.7	17.3	93	16		
Day 53	35	2.10	1.54	47.9	18.0	93	14		
Day 52	33	2.57	1.52	45.0	19.3	91	27		
Day 51	24	2.23	1.44	43.7	19.3	93	23		
Day 50	14	2.29	1.21	42.3	18.5	91	13		
Day 49	12	2.52	1.45	43.8	17.8	92	16		
Day 48	12	1.59	1.52	46.7	18.8	93	23		
Day 47	12	2.53	1.63	47.6	18.5	92	19		
Day 46	11	2.38	1.53	46.8	17.8	92	14		
Day 45	11	1.76	1.70	46.2	18.8	91	30		

Data Table & Timeline

LATITUDE™ NXT Patient Management – Heart Failure Management Report

Trend Data									
Date	HeartLogic Heart Failure Index	S3 (mG)	S1 (mG)	Thoracic Impedance (Ω)	Respiratory Rate (rpm)	Night Heart Rate (bpm)	Sleep Incline (degrees)	Weight (lbs)	Blood Pressure (mmHg)
Day 45	11	1.76	1.70	46.2	18.8	91	30		
Day 44	10	2.38	1.51	46.3	18.0	92	27		
Day 43	21	1.78	1.69	46.5	17.8	92	23		
Day 42	25	1.86	1.61	45.9	18.8	92	22		
Day 41	25	1.68	1.40	45.1	17.8	90	20		
Day 40	21	2.43	1.37	46.8	18.3	88	6		
Day 39	7	2.56	1.38	45.2	18.0	90	18		
Day 38	5	2.62	1.42	45.5	18.5	88	17		
Day 37	5	1.54	1.46	46.1	18.0	90	19		

Date	Activity Level (hour(s))	AT/AF Burden (hour(s), events)	RV Rate during AT/AF (max, mean) (bpm)	Mean Heart Rate (bpm)	% LV Paced (%)	Heart Rate Variability (SDANN) (ms)
Day 98	0.3	24.0, 0	145, 71	71	97	N/R
Day 97	0.2	24.0, 0	144, 71	71	96	N/R
Day 96	0.2	24.0, 0	150, 71	71	96	N/R
Day 95	0.2	24.0, 0	150, 71	71	96	N/R
Day 94	0.3	24.0, 0	138, 71	71	96	N/R
Day 93	0.5	24.0, 0	156, 71	71	96	N/R
Day 92	0.6	24.0, 0	133, 71	71	96	N/R
Day 91	0.3	24.0, 0	120, 71	71	96	N/R
Day 90	0.5	24.0, 0	129, 71	71	95	N/R
Day 89	0.4	24.0, 0	158, 71	71	96	N/R
Day 88	0.2	24.0, 0	129, 71	71	94	N/R
Day 87	0.3	24.0, 0	123, 71	71	94	N/R
Day 86	0.3	24.0, 0	148, 72	72	91	N/R
Day 85	0.1	24.0, 0	153, 71	71	94	N/R
Day 84	0.5	1.4, 0	113, 72	81	96	N/R
Day 83	0.4	0.0, 0	N/R, N/R	81	96	N/R
Day 82	0.4	4.2, 31	132, 72	84	96	N/R
Day 81	0.2	3.9, 14	153, 73	90	97	N/R
Day 80	0.3	5.3, 24	132, 72	89	97	N/R
Day 79	0.4	4.5, 25	166, 72	90	97	N/R
Day 78	0.4	4.1, 19	126, 72	90	97	N/R
Day 77	0.5	8.1, 20	140, 72	86	97	N/R
Day 76	0.2	5.9, 30	148, 72	88	97	N/R
Day 75	0.4	5.1, 19	135, 73	89	96	N/R
Day 74	0.3	5.5, 20	121, 72	89	97	N/R
Day 73	0.6	9.3, 33	134, 71	84	95	N/R
Day 72	0.4	19.1, 25	120, 71	75	96	N/R
Day 71	0.6	17.5, 2	142, 71	77	96	N/R
Day 70	0.6	7.4, 26	155, 72	87	96	N/R
Day 69	0.7	10.5, 16	119, 72	83	96	N/R
Day 68	0.7	11.2, 26	116, 71	82	97	N/R
Day 67	0.8	8.3, 33	145, 72	85	96	N/R
Day 66	0.5	8.1, 16	115, 73	86	96	N/R
Day 65	0.4	6.0, 28	119, 72	87	97	N/R

Data Table & Timeline

LATITUDE™ NXT Patient Management – Heart Failure Management Report

Date	Activity Level (hour(s))	AT/AF Burden (hour(s), events)	RV Rate during AT/AF (max, mean) (bpm)	Mean Heart Rate (bpm)	% LV Paced (%)	Heart Rate Variability (SDANN) (ms)
Day 64	1.1	5.4, 22	150, 72	89	97	N/R
Day 63	0.7	3.8, 30	124, 72	90	97	N/R
Day 62	0.7	3.5, 44	153, 72	90	96	N/R
Day 61	0.6	8.1, 30	118, 72	85	96	N/R
Day 60	0.5	9.5, 29	123, 71	83	97	N/R
Day 59	0.4	11.5, 14	130, 71	82	97	N/R
Day 58	0.4	9.2, 18	119, 71	84	96	N/R
Day 57	0.5	9.1, 20	130, 71	84	97	N/R
Day 56	0.7	15.3, 28	144, 71	79	94	N/R
Day 55	0.9	10.2, 32	141, 71	84	96	N/R
Day 54	0.7	6.5, 36	120, 72	87	97	N/R
Day 53	0.8	7.0, 33	134, 72	86	96	N/R
Day 52	0.5	5.9, 30	115, 72	87	97	N/R
Day 51	0.5	6.6, 28	138, 71	86	98	N/R
Day 50	0.6	8.9, 62	135, 72	85	98	N/R
Day 49	0.5	8.6, 62	120, 71	85	98	N/R
Day 48	0.8	5.2, 27	116, 72	88	97	N/R
Day 47	0.8	6.6, 27	123, 72	87	97	N/R
Day 46	0.4	5.0, 30	116, 72	89	98	N/R
Day 45	0.8	5.4, 27	113, 71	87	98	N/R
Day 44	0.8	4.2, 36	138, 71	89	98	N/R
Day 43	0.5	3.5, 26	156, 72	90	98	N/R
Day 42	1.0	6.5, 34	109, 71	86	98	N/R
Day 41	0.6	10.1, 48	115, 71	83	98	N/R
Day 40	0.7	10.0, 26	118, 71	83	98	N/R
Day 39	0.6	5.9, 22	124, 71	86	98	N/R
Day 38	0.4	7.3, 39	112, 71	85	98	N/R
Day 37	1.2	11.3, 15	153, 71	81	97	N/R

Brief Summary Statement

RESONATE™ HF, RESONATE™ EL, PERCIVA™ HF, PERCIVA™, VIGILANT™ EL, MOMENTUM™ EL ICD – Manual 360199-003

CAUTION

Federal law (USA) restricts this device to sale by or on the order of a physician. Rx only. Prior to use, please see the complete “Directions for Use” for more information on Indications, Contraindications, Warnings, Precautions, Adverse Events, and Operator’s Instructions.

INDICATIONS AND USAGE

Boston Scientific implantable cardioverter defibrillators (ICDs) are intended to provide ventricular antitachycardia pacing (ATP) and ventricular defibrillation for automated treatment of life-threatening ventricular arrhythmias.

CONTRAINDICATIONS

Use of these Boston Scientific pulse generators are contraindicated for the following: patients whose ventricular tachyarrhythmias may have reversible cause, such as: digitalis intoxication, electrolyte imbalance, hypoxia, sepsis; or patients whose ventricular tachyarrhythmias have a transient cause, such as: acute myocardial infarction (MI), electrocution, drowning; or patients who have a unipolar pacemaker.

WARNINGS

Read this manual thoroughly before implantation to avoid damage to the pulse generator and/or lead. For single patient use only. Do not reuse, reprocess, or resterilize. Always have external defibrillation equipment available during implant and electrophysiologic testing. Ensure that an external defibrillator and medical personnel skilled in CPR are present during post-implant device testing should the patient require external rescue. Do not use this pulse generator with another pulse generator. Program the pulse generator Tachy Mode(s) to Off during implant, explant, or postmortem procedures to avoid inadvertent high voltage shocks. Do not kink, twist, or braid the lead with other leads as doing so could cause lead insulation abrasion damage or conductor damage. For leads that require the use of a Connector Tool, use caution handling the lead terminal when the Connector Tool is not present on the lead. Do not directly contact the lead terminal with any surgical instruments or electrical connections such as PSA (alligator) clips, ECG connections, forceps, hemostats, and clamps. Do not contact any other portion of the DF4–LLHH or DF4–LLHO lead terminal, other than the terminal pin, even when the lead cap is in place. Do not use atrial tracking modes in patients with chronic refractory atrial tachyarrhythmias. Tracking of atrial arrhythmias could result in ventricular tachyarrhythmias. Advise patients to seek medical guidance before entering environments that could adversely affect the operation of the active implantable medical device, including areas protected by a warning notice that prevents entry by patients who have a pulse generator. RESONATE HF, RESONATE, PERCIVA HF, PERCIVA, VIGILANT and MOMENTUM devices are considered MR Conditional. For these devices, unless all of the MRI Conditions of Use are met, MRI scanning of the patient does not meet MR Conditional requirements for the implanted system, and significant harm to or death of the patient and/or damage to the implanted system may result. Do not expose patients with non-MR conditional devices to MRI scanning. For potential adverse events applicable when the Conditions of Use are met or not met, refer to the MRI Technical Guide. Do not subject a patient with an implanted pulse generator and/or lead to diathermy. If desired, ensure that Patient Triggered Monitor is enabled prior to sending the patient home. Once the Patient Triggered Monitor feature has been triggered by the magnet and an EGM has been stored, or after 60 days have elapsed from the day that Store EGM was enabled, the patient should not apply the magnet.

PRECAUTIONS

For specific information on precautions, refer to the following sections of the product labeling: clinical considerations, sterilization and storage, implantation, device programming, environmental and medical therapy hazards, hospital and medical environments, home and occupational environments, follow up testing, explant and disposal, supplemental precautionary information.

POTENTIAL ADVERSE EVENTS

Based on the literature and on pulse generator and/or lead implant experience, the following alphabetical list includes the possible adverse events associated with the included devices: Air embolism; Allergic reaction; Bleeding; Bradycardia; Cardiac tamponade; Chronic nerve damage; Component failure; Conductor coil fracture; Death; Elevated thresholds; Erosion; Excessive fibrotic tissue growth; Extracardiac stimulation (muscle/nerve stimulation); Failure to convert an induced arrhythmia; Fluid accumulation; Foreign body rejection phenomena; Formation of hematomas or seromas; Heart block; Heart failure following chronic RV apical pacing; Inability to defibrillate or pace; Inappropriate therapy (e.g., shocks and antitachycardia pacing (ATP) where applicable, pacing; Incisional pain; Incomplete lead connection with pulse generator; Infection including endocarditis; Insulating myocardium during defibrillation with internal or external paddles; Lead dislodgement; Lead fracture; Lead insulation breakage or abrasion; Lead perforation; Lead tip deformation and/or breakage; Local tissue reaction; Loss of capture; Myocardial infarction (MI); Myocardial necrosis; Myocardial trauma (e.g., tissue damage, valve damage); Myopotential sensing; Oversensing/undersensing; Pacemaker-mediated tachycardia (PMT); Pericardial rub, effusion; Pneumothorax; Pulse generator migration; Shunting current during defibrillation with internal or external paddles; Syncope; Tachyarrhythmias, which include acceleration of arrhythmias and early, recurrent atrial fibrillation; Thrombosis/thromboemboli; Valve damage; Vasovagal response; Venous occlusion; Venous trauma (e.g., perforation, dissection, erosion); Worsening heart failure.

For a list of potential adverse events associated with MRI scanning, refer to the MRI Technical Guide

Patients may develop psychological intolerance to a pulse generator system and may experience the following: Dependency; Depression; Fear of premature battery depletion; Fear of a device malfunction.

Indications, Safety and Warnings

CRT-D Systems – RESONATE™HF, RESONATE™, RESONATE™X4, VIGILANT™, VIGILANT™X4, MOMENTUM™, MOMENTUM™ X4

INDICATIONS AND USAGE

These Boston Scientific Cardiac Resynchronization Therapy Defibrillators (CRT-Ds) are indicated for patients with heart failure who receive stable optimal pharmacologic therapy (OPT) for heart failure and who meet any one of the following classifications: Moderate to severe heart failure (NYHA Class III-IV) with EF ≤ 35% and QRS duration ≥ 120 ms; or left bundle branch block (LBBB) with QRS duration ≥ 130 ms, EF ≤ 30%, and mild (NYHA Class II) ischemic or nonischemic heart failure or asymptomatic (NYHA Class I) ischemic heart failure

CONTRAINDICATIONS

There are no contraindications for this device.

WARNINGS

Read this manual thoroughly before implantation to avoid damage to the pulse generator and/or lead. For single patient use only. Do not reuse, reprocess, or resterilize. Always have external defibrillation equipment available during implant and electrophysiologic testing. Ensure that an external defibrillator and medical personnel skilled in CPR are present during post-implant device testing should the patient require external rescue. Do not use this pulse generator with another pulse generator. Program the pulse generator Tachy Mode(s) to Off during implant, explant, or postmortem procedures. Do not kink, twist, or braid the lead with other leads. For leads that require the use of a Connector Tool, use caution handling the lead terminal when the Connector Tool is not present on the lead. Do not directly contact the lead terminal with any surgical instruments or electrical connections such as PSA (alligator) clips, ECG connections, forceps, hemostats, and clamps. Do not contact any other portion of the DF4–LLHH or DF4–LLHO lead terminal, other than the terminal pin, even when the lead cap is in place. Do not contact any other portion of the IS4–LLLL lead terminal, other than the terminal pin, even when the lead cap is in place. When implant a system that uses both a DF4–LLHH or DF4–LLHO and IS4–LLLL lead, ensure that the leads are inserted and secured in the appropriate ports. Do not use atrial tracking modes in patients with chronic refractory atrial tachyarrhythmias. Do not use atrial-only modes in patients with heart failure. Left ventricular

lead dislodgement to a position near the atria can result in atrial oversensing and left ventricular pacing inhibition. Physicians should use medical discretion when implanting this device in patients who present with slow VT. Advise patients to seek medical guidance before entering environments that could adversely affect the operation of the active implantable medical device, including areas protected by a warning notice that prevents entry by patients who have a pulse generator. RESONATE HF, RESONATE, and VIGILANT devices with an IS-1/DF4/IS4 lead connection are considered MR Conditional. For these devices, unless all of the MRI Conditions of Use are met, MRI scanning of the patient does not meet MR Conditional requirements for the implanted system, and significant harm to or death of the patient and/or damage to the implanted system may result. For potential adverse events applicable when the Conditions of Use are met or not met, refer to the MRI Technical Guide. Do not subject a patient with an implanted pulse generator and/or lead to diathermy. If desired, ensure that Patient Triggered Monitor (PTM) is enabled prior to sending the patient home by confirming the magnet response is programmed to Store EGM. Once the PTM feature has been triggered and the magnet response set to Inhibit therapy the patient should not reapply the magnet.

PRECAUTIONS

For specific information on precautions, refer to the following sections of the product labeling: clinical considerations, sterilization and storage, implantation, device programming, environmental and medical therapy hazards, hospital and medical environments, follow up testing, explant and disposal, supplemental precautionary information. Advise patients to avoid sources of EMI because EMI may cause the pulse generator to deliver inappropriate therapy or inhibit appropriate therapy.

POTENTIAL ADVERSE EVENTS

Based on the literature and on pulse generator and/or lead implant experience, the following alphabetical list includes the possible adverse events associated with the included devices: Air embolism; Allergic reaction; Bleeding; Bradycardia; Cardiac tamponade; Chronic nerve damage; Component failure; Conductor coil fracture; Death; Electrolyte imbalance/dehydration; Elevated thresholds; Erosion; Excessive fibrotic tissue growth; Extracardiac stimulation (muscle/nerve stimulation); Failure to convert an induced arrhythmia; Fluid accumulation; Foreign body rejection phenomena; Formation of hematomas or seromas; Heart block; Inability to defibrillate or pace; Inappropriate therapy (e.g., shocks and antitachycardia pacing (ATP) where applicable, pacing); Incisional pain; Incomplete lead connection with pulse generator; Infection including endocarditis; Insulating myocardium during defibrillation with internal or external paddles; Lead dislodgement; Lead fracture; Lead insulation breakage or abrasion; Lead perforation; Lead tip deformation and/or breakage; Local tissue reaction; Loss of capture; Myocardial infarction (MI); Myocardial necrosis; Myocardial trauma (e.g., tissue damage, valve damage); Myopotential sensing; Oversensing/undersensing; Pacemaker-mediated tachycardia (PMT); Pericardial rub, effusion; Pneumothorax; Pulse generator migration; Shunting current during defibrillation with internal or external paddles; Syncope; Tachyarrhythmias, which include acceleration of arrhythmias and early, recurrent atrial fibrillation; Thrombosis/thromboemboli; Valve damage; Vasovagal response; Venous occlusion; Venous trauma (e.g., perforation, dissection, erosion); Worsening heart failure.

For a list of potential adverse events associated with MRI scanning, refer to the MRI Technical Guide

Patients may develop psychological intolerance to a pulse generator system and may experience the following: Dependency; Depression; Fear of premature battery depletion; Fear of a device malfunction.

Additionally, potential adverse events associated with the implantation of a coronary venous lead system include: Allergic reaction to contrast media; Breakage/failure of implant instruments; Prolonged exposure to fluoroscopic radiation; Renal failure from contrast media used to visualize coronary veins.

Refer to the product labeling for specific indications, contraindications, warnings/precautions and adverse events. Rx only.

92436222 (Rev A)

Results from the case studies are not necessarily predictive of results in other cases. Results in other cases may vary.

All trademarks are the property of their respective owners.

Prior to use, review DFU for indications, contraindications, warnings, precautions, adverse events and operating instructions.

Call us: 1-800-227-3422

**Boston
Scientific**
Advancing science for life™

Rhythm Management
300 Boston Scientific Way
Marlborough, MA 01752-1234
www.bostonscientific.com

Medical Professionals:
1.800.CARDIAC (227.3422)
Patients and Families:
1.866.484.3268

© 2019 Boston Scientific Corporation
or its affiliates. All rights reserved.

CRM-701305-AA