# **SYNERGY™**

## **MONORAIL™**

### Everolimus-Eluting Platinum Chromium Coronary Stent System

<b>Directions for Use</b>	<b>2</b>
<b>Instrucciones de uso</b>	<b>11</b>
<b>Mode d'emploi</b>	<b>21</b>
<b>Gebrauchsanweisung</b>	<b>31</b>
<b>Istruzioni per l'uso</b>	<b>41</b>
<b>Gebruiksaanwijzing</b>	<b>51</b>
<b>Instruções de Utilização</b>	<b>61</b>

## TABLE OF CONTENTS

<b>WARNING</b> .....	3
<b>DEVICE DESCRIPTION</b> .....	3
Table 1. SYNERGY™ Stent System Product Description .....	3
Device Component Description .....	3
Contents.....	3
Drug Eluting Coating Description .....	3
Everolimus.....	3
Figure 1. The Chemical Structure of Everolimus .....	3
Polymer Carrier.....	4
Figure 2. The Chemical Structure of PLGA.....	4
Product Matrix.....	4
Table 2. SYNERGY Stent System Product Matrix and Everolimus Content .....	4
<b>INTENDED USE/INDICATIONS FOR USE</b> .....	4
<b>CONTRAINDICATIONS</b> .....	4
<b>WARNINGS</b> .....	4
<b>PRECAUTIONS</b> .....	5
General Precautions.....	5
Stent System Handling (also see, Operational Instructions).....	5
Stent Placement.....	5
Preparation.....	5
Placement.....	5
Stent System Removal - Pre-deployment.....	5
Stent System Removal - Post-deployment .....	5
Table 3. Representative System Deflation Times (seconds) .....	6
Post-Procedure .....	6
Brachytherapy.....	6
Magnetic Resonance Imaging (MRI).....	6
3.0 Tesla Temperature Information.....	6
1.5 Tesla Temperature Information.....	6
Image Artifact Information .....	6
Pre- and Post-Procedure Antiplatelet Regimen .....	6
Drug Interactions .....	6
Use in Special Populations:.....	7
Pregnancy .....	7
Use of Multiple Stents .....	7
Drug Information: .....	7
Mechanism of Action .....	7
Drug Interactions .....	7
Carcinogenicity, Genotoxicity, and Reproductive Toxicology .....	7
<b>ADVERSE EVENTS</b> .....	7
<b>HOW SUPPLIED</b> .....	8
Handling and Storage.....	8
<b>OPERATIONAL INSTRUCTIONS</b> .....	8
Device Selection .....	8
Inspection Prior to Use .....	8
Material Required (not included in stent system package).....	8
Preparation.....	8
Package Removal .....	8

Guidewire Lumen Flush.....	8
Balloon Preparation.....	9
Delivery Procedure .....	9
Deployment Procedure .....	9
Removal Procedure and Completion .....	9
<b>IN VITRO INFORMATION</b> .....	10
SYNERGY Stent System Compliance .....	10
Table 4. SYNERGY Compliance Chart.....	10
<b>WARRANTY</b> .....	10

# SYNERGY™

MONORAIL™

## Everolimus-Eluting Platinum Chromium Coronary Stent System

### Rx ONLY

**Caution:** Federal Law (USA) restricts this device to sale by or on the order of a physician.

#### WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

#### DEVICE DESCRIPTION

The SYNERGY Everolimus-Eluting Platinum Chromium Coronary Stent System (SYNERGY Stent System) is a device/drug combination product comprised of two regulated components: a device (coronary stent system) and a drug product (a formulation of everolimus contained in a biodegradable polymer coating). SYNERGY is uniquely designed with a low initial polymer load, abluminal coating and bioabsorbable polymer which may reduce the risk of thrombosis and the need for prolonged dual antiplatelet therapy. The characteristics of the SYNERGY Stent System are described in Table 1. SYNERGY Stent System Product Description.

**Table 1. SYNERGY Stent System Product Description**

Characteristic	SYNERGY Stent System
Available stent lengths (mm)	8, 12, 16, 20, 24, 28, 32, 38
Available stent diameters (mm)	2.25, 2.50, 2.75, 3.00, 3.50, 4.00
Stent material	Platinum Chromium (PtCr) Alloy
Drug product	An abluminal (outer surface of the stent) coating of a polymer carrier with approximately 1 µg of everolimus per mm <sup>2</sup> of total stent surface area and a maximum nominal polymer content of 351 µg and drug content of 287 µg on the largest stent (4.00 x 38 mm)
Delivery system effective length	144 cm
Delivery system Y-adapter ports	Single access port to inflation lumen. Guidewire exit port is located approximately 25 cm from tip. Designed for guidewire ≤ 0.014 in (0.36 mm).
Average stent length change upon deployment at nominal diameter	Small Vessel (SV) average: -0.15 mm Workhorse (WH) average: -0.20 mm Large Vessel (LV) average: 0.15 mm
Stent delivery balloon	A balloon with two radiopaque markers nominally 0.4 mm longer than the stent at each end.

Characteristic	SYNERGY Stent System
Balloon inflation pressure	Nominal inflation pressure: 11 atm (1117 kPa) Rated burst pressure: 2.25 - 2.75 mm: 18 atm (1827 kPa) 3.00 - 4.00 mm: 16 atm (1620 kPa)
Guide catheter inner diameter	≥ 0.056 in (1.42 mm)
Catheter shaft outer diameter	Proximal: 2.1F (0.70 mm) Distal: 2.25 - 2.75 mm: 2.6F (0.90 mm) 3.00 mm: • 8 - 28 mm: 2.6F (0.90 mm) • 32 - 38 mm: 2.7F (0.95 mm) 3.50 mm: • 8 - 20 mm: 2.6F (0.90 mm) • 24 - 38 mm: 2.7F (0.95 mm) 4.00 mm: 2.7F (0.95 mm)
Stent strut thickness	2.25 - 2.75 mm: 0.074 mm 3.00 - 3.50 mm: 0.079 mm 4.00 mm: 0.081 mm

#### Device Component Description

The SYNERGY Stent System consists of a platinum chromium stent platform with an abluminal drug/polymer coating mounted onto a Monorail Delivery System.

The SYNERGY stent is available in 3 stent models each designed for specific diameters as follows:

- Small Vessel (SV): 2.25, 2.50 and 2.75 mm
- Workhorse (WH): 3.00 and 3.50 mm
- Large Vessel (LV): 4.00 mm

#### Contents

Qty	Material
One (1)	SYNERGY Stent System
One (1)	Flushing needle with Luer fitting

#### Drug Eluting Coating Description


The SYNERGY stent is a stent with a drug/polymer coating. The coating is comprised of a polymer matrix that contains an active pharmaceutical ingredient.

See **Everolimus** and **Polymer Carrier** sections for descriptions of drug and polymer, respectively.

#### Everolimus

Everolimus is the active pharmaceutical ingredient in the SYNERGY stent.

The everolimus chemical name is 42-O-(2-hydroxyethyl)-rapamycin, and its chemical structure is shown below in Figure 1.

**Figure 1. The Chemical Structure of Everolimus**

## Polymer Carrier

The SYNERGY™ stent is coated on the abluminal stent surface with a biodegradable drug matrix. The biodegradable drug matrix is composed of PLGA [poly(DL-lactide-co-glycolide)] mixed with everolimus. The chemical structure of PLGA is shown below in Figure 2.


Figure 2. The Chemical Structure of PLGA

## Product Matrix

Table 2. SYNERGY Stent System Product Matrix and Everolimus Content

UPN Number	Nominal Expanded Inner Diameter (mm)	Nominal Unexpanded Stent Length (mm)	Nominal Everolimus Content (µg)
H7493926208220	2.25	8	38.9
H7493926208250	2.50	8	38.9
H7493926208270	2.75	8	38.9
H7493926208300	3.00	8	46.5
H7493926208350	3.50	8	46.5
H7493926208400	4.00	8	67.5
H7493926212220	2.25	12	58.3
H7493926212250	2.50	12	58.3
H7493926212270	2.75	12	58.3
H7493926212300	3.00	12	66.3
H7493926212350	3.50	12	66.3
H7493926212400	4.00	12	96.2
H7493926216220	2.25	16	77.6
H7493926216250	2.50	16	77.6
H7493926216270	2.75	16	77.6
H7493926216300	3.00	16	92.7
H7493926216350	3.50	16	92.7
H7493926216400	4.00	16	124.8
H7493926220220	2.25	20	96.9
H7493926220250	2.50	20	96.9
H7493926220270	2.75	20	96.9
H7493926220300	3.00	20	112.5
H7493926220350	3.50	20	112.5
H7493926220400	4.00	20	153.5
H7493926224220	2.25	24	121.1
H7493926224250	2.50	24	121.1
H7493926224270	2.75	24	121.1
H7493926224300	3.00	24	132.3
H7493926224350	3.50	24	132.3

UPN Number	Nominal Expanded Inner Diameter (mm)	Nominal Unexpanded Stent Length (mm)	Nominal Everolimus Content (µg)
H7493926224400	4.00	24	182.2
H7493926228220	2.25	28	140.5
H7493926228250	2.50	28	140.5
H7493926228270	2.75	28	140.5
H7493926228300	3.00	28	158.7
H7493926228350	3.50	28	158.7
H7493926228400	4.00	28	210.8
H7493926232220	2.25	32	159.8
H7493926232250	2.50	32	159.8
H7493926232270	2.75	32	159.8
H7493926232300	3.00	32	178.5
H7493926232350	3.50	32	178.5
H7493926232400	4.00	32	239.5
H7493926238220	2.25	38	188.9
H7493926238250	2.50	38	188.9
H7493926238270	2.75	38	188.9
H7493926238300	3.00	38	211.6
H7493926238350	3.50	38	211.6
H7493926238400	4.00	38	287.2

## INTENDED USE/INDICATIONS FOR USE

The SYNERGY Stent System is indicated for improving coronary luminal diameter in patients with symptomatic ischemic heart disease, including patients with acute myocardial infarction (AMI), due to discrete de novo native coronary artery lesions. The treated lesion length should be less than the nominal stent length (8 mm, 12 mm, 16 mm, 20 mm, 24 mm, 28 mm, 32 mm and 38 mm) with a reference vessel diameter of 2.25 mm - 4.0 mm.

## CONTRAINDICATIONS

Use of the SYNERGY Stent System is contraindicated in patients with the following:

- Known hypersensitivity to platinum, the platinum chromium alloy, or similar alloy types such as stainless steel.
- Known hypersensitivity to everolimus or structurally related compounds.
- Known hypersensitivity to the polymer or its individual components (see details in **Polymer Carrier**).
- Known severe reaction to contrast agents that cannot be adequately pre-medicated prior to the SYNERGY stent placement procedure.

Coronary artery stenting is contraindicated for use in the following:

- Patients who cannot receive recommended anti-platelet and/or anticoagulant therapy.
- Patients judged to have a lesion that prevents complete inflation of an angioplasty balloon or proper placement of the stent or delivery device.

## WARNINGS

- This product should not be used in patients who are not likely to comply with recommended anti-platelet therapy.
- To maintain sterility, the package should not be opened or damaged prior to use. The package should be opened as described in **Operational Instructions**.
- The use of this product carries the risks associated with coronary artery stenting, including stent thrombosis, vascular complications, and/or bleeding events.

- Patients with known hypersensitivity to platinum or stainless steel may suffer an allergic reaction to this implant.

## PRECAUTIONS

### General Precautions

- Only physicians who have received adequate training should perform implantation of the stent.
- Stent placement should only be performed at hospitals where emergency coronary artery bypass graft (CABG) surgery is readily available.
- Potential interactions of the SYNERGY™ stent with other drug-eluting or coated stents have not been evaluated.
- Subsequent restenosis may require repeat dilatation of the arterial segment containing the stent. The long-term outcome following repeat dilatation of coronary stents is unknown at present.
- Consideration should be given to the risks and benefits of use in patients with a history of severe reaction to contrast agents.
- Do not expose the stent delivery system to organic solvents such as alcohol or detergents.
- Care should be taken to control the position of the guide catheter tip during stent delivery, deployment, and balloon withdrawal. Before withdrawing the stent delivery system (SDS), visually confirm complete balloon deflation by fluoroscopy (see Table 3, Representative System Deflation Times). Failure to do so may cause increased SDS withdrawal forces and result in guide catheter movement into the vessel and subsequent arterial damage.
- Orally administered everolimus combined with cyclosporine is associated with increased serum cholesterol and triglyceride levels.

### Stent System Handling (also see, Operational Instructions)

- For single use only. Do not resterilize or reuse this product. Note product “Use By” date and do not use after the “Use By” date.
- The SYNERGY stent and its delivery system are designed for use as a unit. The stent is not to be removed from its delivery balloon. The stent is not designed to be crimped onto another balloon. Removing the stent from its delivery balloon may damage the stent and coating and/or lead to stent embolization.
- Prior to angioplasty, carefully examine all equipment to be used during the procedure including the dilatation catheter to verify proper function.
- Special care must be taken not to handle or in any way disrupt the stent position on the delivery balloon. This is most important during catheter removal from packaging, placement over the guidewire, and advancement through the hemostatic valve and guide catheter hub.
- Excessive manipulation or handling may cause coating damage, contamination, or dislodgment of the stent from the delivery balloon.
- Use only the appropriate balloon inflation media (see **Operational Instructions, Balloon Preparation**). DO NOT use air or any gas medium to inflate the balloon.
- In the event that the SYNERGY stent is not deployed, follow product return procedures and avoid handling of the stent with bare hands.
- Stent contact with any fluid prior to placement is not recommended as there is a possibility of drug release. However, if it is absolutely necessary to flush or soak the stent with sterile/isotonic saline, contact time should be limited (1 minute maximum).

### Stent Placement

#### Preparation

- DO NOT PREPARE OR PRE-INFLATE BALLOON PRIOR TO STENT DEPLOYMENT OTHER THAN AS DIRECTED. Use the balloon purging technique described in **Operational Instructions, Balloon Preparation**.
- If unusual resistance is felt at any time during lesion access before stent implantation, the stent system and guide catheter should be removed as a single unit (see **Stent System Removal - Pre-deployment** and **Stent System Removal - Post-deployment**).
- An unexpanded stent should be introduced into the coronary arteries **one time only**. An unexpanded stent should not be subsequently moved in and out through the distal end of the guide catheter as stent or coating damage or stent dislodgment from the balloon may occur.

### Placement

- The vessel should be pre-dilated with an appropriately sized balloon. Failure to do so may increase the risk of placement difficulty and procedural complications.
- Do not expand the stent if it is not properly positioned in the vessel (see **Precautions, Stent System Removal - Pre-deployment**).
- Balloon pressures should be monitored during inflation. Do not exceed the rated burst pressure as indicated on the product label (see Table 4, SYNERGY Compliance Chart). Use of pressures higher than specified on the product label may result in a ruptured balloon or shaft. This may result in potential intimal damage, dissection or vessel rupture.
- The stent inner diameter should approximate 1.1 times the distal reference vessel diameter.
- Placement of the stent has the potential to compromise side branch patency if stenting near a side branch.
- Implanting a stent may lead to dissection of the vessel distal and/or proximal to the stented portion and may cause acute closure of the vessel requiring additional intervention (e.g., CABG, further dilation, placement of additional stents, or other).
- When treating multiple lesions, the distal lesion should generally be stented first, followed by stenting of the more proximal lesion(s). Stenting in this order avoids the requirement to cross the proximal stent when placing the distal stent and reduces the chances of the stent dislodgment.

### Stent System Removal - Pre-deployment

- If unusual resistance is felt at any time during lesion access before stent implantation, the stent system and the guide catheter should be removed as a single unit (see note below).
- Do not attempt to pull an unexpanded stent back into the guide catheter, as stent or coating damage or stent dislodgment from the balloon may occur.
- Stent retrieval methods (use of additional wires, snares, and/or forceps) may result in additional trauma to the vascular site. Complications can include bleeding, hematoma, or pseudoaneurysm.

**Note:** When removing the entire stent system and guide catheter as a single unit the following steps should be executed in the order indicated under direct visualization using fluoroscopy.

- If greater than usual resistance is felt during delivery system withdrawal, pay particular attention to guide catheter position. In some cases it may be necessary to pull back slightly on the guide catheter in order to prevent deep seating (unplanned advancement) of the guide catheter and subsequent vessel damage. In cases where unplanned guide catheter movement has occurred, angiographic assessment of the coronary tree should be undertaken to ensure that there is no damage to the coronary vasculature.
- Maintain guidewire placement across the lesion during the entire removal process. Carefully pull back the stent system until the proximal balloon marker of the stent system is just distal to the guide catheter distal tip.
- The stent system and the guide catheter should be pulled back as a single unit until the tip of the guide catheter is just distal to the arterial sheath, allowing the guide catheter to straighten. Carefully retract the un-deployed stent into the tip of the guide catheter and remove the stent system and the guide catheter from the patient again as a single unit while leaving the guidewire across the lesion.

### Stent System Removal - Post-deployment

- Following stent placement, confirm complete balloon deflation (Table 3, Representative System Deflation Times).
- If greater than usual resistance is felt during delivery system withdrawal, pay particular attention to guide catheter position. In some cases it may be necessary to pull back slightly on the guide catheter in order to prevent deep seating (unplanned advancement) of the guide catheter and subsequent vessel damage. In cases where unplanned guide catheter movement has occurred, angiographic assessment of the coronary tree should be undertaken to ensure that there is no damage to the coronary vasculature.

- If greater than usual resistance is felt during delivery system withdrawal into the guide catheter, the stent system and guide catheter should be removed as a single unit (see note in above section).

**Table 3. Representative System Deflation Times (seconds)**

Balloon Length / Diameter	8 mm	12 mm	16 mm	20 mm	24 mm	28 mm	32 mm	38 mm
2.25 mm	≤ 16	≤ 16	≤ 16	≤ 16	≤ 16	≤ 16	≤ 16	≤ 21
2.50 mm								
2.75 mm								
3.00 mm				≤ 21	≤ 21	≤ 21	≤ 30	
3.50 mm								
4.00 mm								

#### Post-Procedure

- Care must be exercised when crossing a newly deployed stent with ancillary devices to avoid disrupting the stent placement, apposition, geometry, and/or coating.

If the patient requires Magnetic Resonance Imaging (MRI), see **Magnetic Resonance Imaging**.

#### Brachytherapy

The safety and effectiveness of the SYNERGY™ stent in patients with prior brachytherapy of the target lesion have not been established.

The safety and effectiveness of the use of brachytherapy to treat in-stent restenosis in the SYNERGY stent has not been established.

Both vascular brachytherapy and the SYNERGY stent alter arterial remodeling. The interaction, if any, between these two treatments has not been determined.

#### Magnetic Resonance Imaging (MRI)

Through non-clinical testing, the SYNERGY stent has been shown to be MR Conditional (poses no known hazards under specified conditions). The conditions are as follows:

- Field strengths of 3.0 and 1.5 Tesla with
  - Static magnetic field gradient < 11 T/m (extrapolated).
  - Product of static magnetic field and static magnetic field gradient < 25 T<sup>2</sup>/m (extrapolated).
- A calculated rate of change of magnetic field (dB/dt) of 60 T/s or less.
- A maximum whole body averaged specific absorption rate (SAR) of lower than 2.0 W/kg for a total active MR scan time (with RF exposure) of 15 minutes or less. The SYNERGY stent should not migrate in this MRI environment. MR imaging within these conditions may be performed immediately following the implantation of the stent. This stent has not been evaluated to determine if it is MR Conditional beyond these conditions.

#### 3.0 Tesla Temperature Information

Non-clinical testing of RF-induced heating was performed at 123 MHz in a 3.0 Tesla Magnetom Trio™, Siemens Medical Solutions MR system, software version Numaris/4, syngo™ MR A30A. The stents were in a location and orientation in the phantom that produced the worst case Radio Frequency (RF) heating. RF power was applied for 15 minutes and the measured conductivity of the phantom material was about 0.49 S/m. The phantom average SAR calculated using calorimetry was 2.3 W/kg. The maximal in-vitro temperature rise was calculated as 2.6°C when the local SAR was scaled to 2.0 W/kg for a measured stent length up to 74 mm. Predicted in-vivo heating based on these non-clinical tests and computer simulation of the patient exposure to the electromagnetic fields in MRI yielded to the following maximal in vivo rises: for landmarks at the chest level, the calculated temperature rise was 2.6°C with a calculated uncertainty upper bound temperature of 4.7°C for a whole body average SAR value of 2.0 W/kg and a continuous scan time of 15 minutes.

The actual in vivo rise is expected to be less than these values as the calculations did not include the cooling effects due to blood flow in the lumen of the stent and blood perfusion in the tissue outside the stent.

#### 1.5 Tesla Temperature Information

Non-clinical testing of RF-induced heating was performed at 64 MHz in a 1.5 Tesla Intera™ Philips Medical Systems, software version Release 12.6.1.3, 2010-12-02 whole body coil MR scanner. The stents were in a location and orientation in the phantom that produced the worst case RF heating. RF power was applied for 15 minutes and the measured conductivity of the phantom material was about 0.50 S/m. The phantom average SAR calculated using calorimetry was 2.3 W/kg. The maximal in-vitro temperature rise was calculated as 2.6°C when the local SAR was scaled to 2.0 W/kg for a measured stent length up to 74 mm. Predicted in-vivo heating based on these non-clinical tests and computer simulation of the patient exposure to the electromagnetic fields in MRI yielded to the following maximal in vivo rises: for landmarks at the chest level, the calculated temperature rise was 2.6°C with an uncertainty upper bound temperature of 4.8°C for a whole body average SAR value of 2.0 W/kg and a continuous scan time of 15 minutes.

The actual in vivo rise is expected to be less than these values as the calculations did not include the cooling effects due to blood flow in the lumen of the stent and blood perfusion in the tissue outside the stent.

In vivo, local SAR depends on MR Field strength and may be different than the estimated whole body averaged SAR, due to body composition, stent position within the imaging field, and scanner used, thereby affecting the actual temperature rise. No tests have been performed on possible nerve or other tissue stimulation possible to be activated by strong gradient magnetic fields and resulting induced voltages.

#### Image Artifact Information

The calculated image artifact extends approximately 7 mm from the perimeter of the device diameter and 5 mm beyond each end of the length of the stent when scanned in non-clinical testing using a Spin Echo sequence. With a Gradient Echo sequence the calculated image artifact extends 7 mm beyond the perimeter of the diameter and 6 mm beyond each end of the length with both sequences partially shielding the lumen in a 3.0 Tesla Intera (Achieva Upgrade), Philips Medical Solutions, software version Release 2.6.3.5 2009-10-12 MR system with a transmit/receive head coil. This testing was completed using ASTM F2119-07 test method.

#### Pre- and Post-Procedure Antiplatelet Regimen

The device carries an associated risk of acute, subacute, or late thrombosis, vascular complications, and/or bleeding events. Therefore, the patient should be carefully selected, and a P2Y<sub>12</sub> inhibitor (i.e., clopidogrel, ticlopidine, prasugrel, or ticagrelor) must be prescribed post procedure to reduce risk of thrombosis. Aspirin must be administered concomitantly with P2Y<sub>12</sub> inhibitor, and then continued indefinitely to reduce the risk of thrombosis. SYNERGY is designed with a low initial polymer load, abluminal coating and bioabsorbable polymer which may reduce the risk of thrombosis and the need for prolonged dual antiplatelet therapy. It is strongly advised that the treating physician consider the European Society of Cardiology recommendations (or other applicable country guidelines) for antiplatelet therapy pre- and post-procedure to reduce the risk of thrombosis. In selected patients, it may be reasonable to interrupt or discontinue P2Y<sub>12</sub> inhibitor therapy after 3 months.

It is very important that the patient be compliant with the post-procedural antiplatelet recommendations. Premature discontinuation of prescribed antiplatelet medication could result in a higher risk of thrombosis, myocardial infarction, or death. This should be carefully considered by the treating physicians prior to Percutaneous Coronary Intervention (PCI) for patients who may require premature cessation of antiplatelet therapy, e.g., for surgical or dental procedures. Patients who require premature discontinuation of antiplatelet therapy due to significant active bleeding or the expectation of significant active bleeding should be monitored carefully for cardiac events and once stabilised have their antiplatelet therapy restarted without unnecessary delay.

#### Drug Interactions

When taken orally, everolimus is extensively metabolized by the cytochrome P4503A4 (CYP3A4) in the gut wall and liver and is a substrate for the countertransporter P-glycoprotein. Therefore, absorption and subsequent elimination of everolimus may be influenced by drugs that affect these pathways. Concurrent treatment with strong 3A4 inhibitors and inducers is not recommended unless the benefits outweigh the risk. Inhibitors of P-glycoprotein may decrease the efflux of everolimus from intestinal cells and increase everolimus blood concentrations. In vitro, everolimus was a competitive inhibitor of CYP3A4 and of CYP2D6, potentially increasing the concentrations of drugs eliminated by these enzymes. Thus, caution should be exercised when coadministering everolimus with 3A4 and 2D6 substrates with a

narrow therapeutic index. Everolimus has also been shown to reduce the clearance of some prescription medications when it was administered orally along with cyclosporine (CsA).

Everolimus, when prescribed as an oral medication, may interact with the following drugs or substances.

**Note:** The list below describes interactions for orally administered everolimus at significantly higher doses than are present on the SYNERGY™ Stent System. Interactions observed at these higher, oral doses may not be relevant to the SYNERGY Stent System.

- CYP3A4 isozyme inhibitors (ketoconazole, itraconazole, ritonavir, erythromycin, clarithromycin, fluconazole, calcium channel blockers)
- Inducers of CYP3A4 isozyme (rifampicin, rifabutin, carbamazepin, phenobarbital, phenytoin)
- Antibiotics (Ciprofloxacin, ofloxacin)
- Glucocorticoids
- HMGCoA reductase inhibitors (simvastatin, lovastatin)
- Digoxin
- Cisapride (theoretical potential interaction)
- Sildenafil (Viagra™) (theoretical potential interaction)
- Antihistaminics (terfenadine, astemizole)
- Grapefruit juice

Because systemic everolimus levels are below the lower limit of detection in pre-clinical studies after two days, formal drug interaction studies have not been performed with SYNERGY Stent System. Therefore, due consideration should be given to the potential for both systemic and local drug interactions in the vessel wall when deciding to place the SYNERGY stent in a subject taking a drug with known interaction with everolimus.

#### Use in Special Populations:

##### Pregnancy

This product has not been tested in pregnant women or in men intending to father children; effects on the developing foetus have not been studied. While there is no contraindication, the risks and reproductive effects are unknown. It is not recommended that the SYNERGY Stent System be used in women attempting to conceive, or who are pregnant.

##### Use of Multiple Stents

Potential interactions of the SYNERGY stent with other drug-eluting or coated stents have not been evaluated in vivo. Patients should be treated with no more than 2 planned SYNERGY stents. Additional stents may be placed if bailout stenting is required. The use of multiple drug-eluting stents will expose the patient to larger amounts of drug and polymer.

When more than one stent is required and results in stent-to-stent contact, stent materials should be of similar composition to avoid the possibility of corrosion due to the presence of dissimilar metals in a conducting medium. Placing multiple stents of different metals in contact with each other may increase the potential for corrosion, though in vitro tests to assess stent-to-stent contact using a platinum chromium alloy stent in combination with a 316L stainless steel or cobalt-chromium alloy stent suggest there is no increased risk of corrosion with this pair. If more than one SYNERGY stent is needed to cover the lesion, it is suggested that, to avoid the potential for gap restenosis, the stents be adequately overlapped (with a minimum of 2 mm overlap).

##### Drug Information:

##### Mechanism of Action

The mechanism by which the SYNERGY stent inhibits neointimal growth has not been established. At the cellular level, everolimus inhibits growth factor-stimulated cell proliferation. At the molecular level, everolimus forms a complex with the cytoplasmic protein FKBP-12 (FK 506 Binding Protein). This complex binds to and interferes with FRAP (FKBP-12 Rapamycin Associated Protein), also known as mTOR (mammalian Target of Rapamycin), leading to inhibition of cell metabolism, growth and proliferation by arresting the cell cycle at the late G1 stage.

##### Drug Interactions

See **Precautions, Drug Interactions.**

#### Carcinogenicity, Genotoxicity, and Reproductive Toxicology

Carcinogenicity, genotoxicity, and reproductive toxicology of SYNERGY stent were not evaluated. However, testing has been completed on PROMUS (Xience V™). PROMUS (Xience V) and SYNERGY use the same drug (everolimus) and release profile. A 26-week carcinogenicity study was conducted to evaluate the carcinogenic potential of PROMUS (Xience V) everolimus-eluting stents following subcutaneous implantation in transgenic mice. During the course of the study, there were no abnormal clinical observations that suggested a carcinogenic effect of the test group PROMUS (Xience V). The test group did not demonstrate an increased incidence of neoplastic lesions when compared to the negative control group. However, the positive control and the experimental positive control groups demonstrated notable increases in the incidence of neoplastic lesions compared to either the test or the negative control group. Based on the results of this study, the PROMUS (Xience V) stent does not appear to be carcinogenic when implanted in transgenic mice for 26 weeks.

In addition, a reproductive toxicity (teratology) study was conducted to demonstrate that implantation of PROMUS (Xience V) stents in female Sprague-Dawley rats does not affect their fertility or reproductive capability and shows a lack of any reproductive toxicity on their offspring. The PROMUS (Xience V) stent did not affect the fertility or reproductive capability of female Sprague-Dawley rats. There was no statistical difference between the test article PROMUS (Xience V) stent and the control system in terms of any of the evaluated parameters. The test article had no effect on litter size and caused no increase of in utero mortality. Additionally, the PROMUS (Xience V) stent did not cause any reproductive toxicity in the offspring in this study.

#### ADVERSE EVENTS

Potential adverse events (in alphabetical order) which may be associated with the implantation of a coronary stent in a native coronary artery include those risks associated with percutaneous transluminal coronary angioplasty as well as additional risks related to the use of the stent as listed below.

- Abrupt stent closure
- Acute myocardial infarction
- Allergic reaction to anti-coagulant and/or antiplatelet therapy, contrast medium, or stent materials
- Angina
- Arrhythmias, including ventricular fibrillation and ventricular tachycardia
- Arteriovenous fistula
- Bleeding
- Cardiac tamponade
- Cardiogenic shock/pulmonary edema
- Coronary aneurysm
- Death
- Dissection
- Emboli, distal (air, tissue, or thrombotic material or material from device(s) used in the procedure)
- Heart failure
- Hematoma
- Hemorrhage, which may require transfusion
- Hypotension/hypertension
- Infection, local or systemic
- Ischemia, myocardial
- Pain, access site
- Perforation or rupture of coronary artery
- Pericardial effusion
- Pseudoaneurysm, femoral
- Renal insufficiency or failure
- Respiratory failure

- Restenosis of stented segment
- Stent deformation
- Stent embolization or migration
- Stent fracture
- Stent thrombosis/vessel occlusion
- Stroke/cerebrovascular accident/transient ischemic attack
- Total occlusion of coronary artery
- Vessel spasm
- Vessel trauma requiring surgical repair or reintervention

Adverse events associated with daily oral administration of everolimus (or potential adverse events not captured above, that may be unique to the everolimus drug coating):

- Abdominal pain
- Acne
- Allergic/immunologic reaction to drug (everolimus or structurally related compounds) or the polymer stent coating or its individual components (see **Drug Component Description**)
- Anemia
- Coagulopathy
- Diarrhea
- Edema
- Hemolysis
- Hypercholesterolemia
- Hyperlipidemia
- Hypertension
- Hypertriglyceridemia
- Hypogonadism male
- Leukopenia
- Liver function test abnormal
- Lymphocele
- Myalgia/Arthralgia
- Nausea
- Pain
- Pneumonia
- Pyelonephritis
- Rash
- Renal tubular necrosis
- Sepsis
- Surgical wound complication
- Thrombocytopenia
- Urinary tract infection
- Venous thromboembolism
- Viral, bacterial, and fungal infections
- Vomiting
- Wound infection

There may be other potential adverse events that are unforeseen at this time.

#### HOW SUPPLIED

Non-pyrogenic.

Do not use if package is opened or damaged.

Do not use if labeling is incomplete or illegible.

#### Handling and Storage

Keep dry and protect from light.

Store at 25°C (77°F); excursions permitted to 15-30°C (59-86°F).

Store product in outer carton until ready for use.

DO NOT REMOVE FROM FOIL POUCH UNTIL READY TO USE.

THE FOIL POUCH IS NOT A STERILE BARRIER.

Do not store devices where they are directly exposed to organic solvents or ionizing radiation.

The foil pouch contains nitrogen gas (N<sub>2</sub>) and desiccant as a storage medium.

#### OPERATIONAL INSTRUCTIONS

##### Device Selection

Select device(s) with nominal stent length(s) and diameter(s) appropriate for the lesion.

##### Inspection Prior to Use

Check foil pouch for “Use By” date. Carefully inspect the foil pouch and the sterile package before opening. Do not use the product after the “Use By” date. If the integrity of the foil pouch or the sterile package has been compromised prior to the product “Use By” date (e.g., damage of the package), contact your local Boston Scientific representative for return information. Do not use if any defects are noted.

##### Material Required (not included in stent system package)

Quantity	Material
1	Guide catheter with ≥ 0.056 in (1.42 mm) inner diameter
2-3	20-mL (cc) syringe
1000 u / 500 cc	Normal heparinised saline
1	≤ 0.014 in (0.36 mm) guidewire
1	Rotating hemostatic valve
1	Diluted contrast medium 1:1 with normal heparinised saline
1	Inflation device (with luer fitting)
1	Torque device (optional)
1	Pre-deployment dilatation catheter
1	Three-way stopcock
1	Appropriate arterial sheath

##### Preparation

##### Package Removal

##### Step Action

1. Open the outer box to reveal the foil pouch and carefully inspect the foil pouch for damage.
2. Carefully open the foil pouch by tearing along the tear strip as indicated on the foil pouch to access the sterile barrier package containing the stent delivery system.
3. Carefully inspect the sterile barrier package for damage.
4. Carefully peel open the sterile barrier using aseptic techniques and extract the stent delivery system.
5. Carefully remove the stent delivery system from its protective tubing for preparation of the delivery system. Do not bend or kink the device during removal.
6. Remove the product mandrel and stent protector by grasping the catheter just proximal to the stent (at the proximal balloon bond site) and, with the other hand, grasp the stent protector and gently remove distally.
7. Examine the device for any damage. If it is suspected that the sterility or performance of the device has been compromised, the device should not be used.

##### Guidewire Lumen Flush

##### Step Action

1. Flush stent system guidewire lumen with normal heparinised saline using the flushing needle supplied for the Monorail™ Delivery System at the distal end.

- Verify that the stent is positioned between the proximal and distal balloon markers. Check for bends, kinks and other damage. Do not use if any defects are noted.

---

**Note:** Use caution while flushing guidewire lumen with flushing needle to avoid damage to catheter tip.

---

**Note:** Avoid manipulation of the stent during flushing of the guidewire lumen, as this may disrupt the placement of the stent on the balloon.

---

**Note:** Stent contact with any fluid is not recommended as there is a possibility of initiating drug release. However, if it is absolutely necessary to flush the stent with saline, contact time should be limited (1 minute maximum).

---

#### Balloon Preparation

##### Step Action

- Prepare inflation device/syringe with diluted contrast medium.
- Attach inflation device/syringe to stopcock; attach to inflation port. Do not bend the hypotube when connecting to inflation device/syringe.
- With tip down, orient stent system vertically.
- Open stopcock to stent system; pull negative pressure for 15 seconds; release to neutral for contrast fill.
- Close stopcock to stent system; purge inflation device/syringe of all air.
- Repeat steps 4 through 6 until all air is expelled. If bubbles persist, do not use product.
- If a syringe was used, attach a prepared inflation device to stopcock.
- Open stopcock to stent system.
- Leave on neutral.

#### Delivery Procedure

##### Step Action

- Prepare the vascular access site according to standard practice.
- Predilate the lesion/vessel with appropriate diameter balloon.
- Maintain neutral pressure on inflation device attached to stent system.
- Backload stent system onto proximal portion of guidewire while maintaining guidewire position across target lesion.
- Fully open rotating hemostatic valve to allow for easy passage of the stent and prevent damage to the stent.
- Carefully advance the stent system into the hub of the guide catheter. Be sure to keep the hypotube straight. Ensure guide catheter stability before advancing the stent system into the coronary artery.

---

**Note:** If unusual resistance is felt before the stent exits the guide catheter, do not force passage. Resistance may indicate a problem. Use of excessive force may result in stent damage or stent dislodgment from the balloon. Maintain guidewire placement across the lesion, and remove the stent system and guide catheter as a single unit.

---

- Advance the stent system over the guidewire to target lesion under direct fluoroscopic visualisation. Use the proximal and distal radiopaque balloon markers as a reference point. Fully cover the entire lesion and balloon treated area. The stent should adequately cover healthy vessel proximal and distal to the lesion. If the position of the stent is not optimal, it should be carefully repositioned or removed (see also **Precautions, Stent System Removal - Pre-deployment**). The inside edges of the marker bands indicate both the stent edges and balloon shoulders. Expansion of the stent should not be undertaken if the stent is not properly positioned in the target lesion segment of the vessel.

---

**Note:** If unusual resistance is felt at any time during lesion access before stent implantation, the stent system and the guide catheter should be removed as a single unit (see also **Precautions, Stent System Removal - Pre-deployment**). Once the stent delivery system has been removed, do not reuse.

---

- Sufficiently tighten the rotating hemostatic valve. The stent is now ready to be deployed.

#### Deployment Procedure

##### Step Action

- Inflate the delivery system expanding the stent to a minimum pressure of 11 atm - 1117 kPa. Higher pressure may be necessary to optimise stent apposition to the arterial wall. Accepted practice generally targets an initial deployment pressure that would achieve a stent inner diameter of about 1.1 times the distal reference vessel diameter (see Table 4.). Balloon pressure must not exceed rated burst pressure of 18 atm - 1827 kPa for the 2.25–2.75 mm sizes and 16 atm - 1620 kPa for the 3.00-4.00 mm sizes (see Table 4.).
- Maintain inflation pressure for 15-30 seconds for full expansion of the stent.
- Deflate balloon by pulling negative pressure on inflation device until balloon is fully deflated.
- Confirm stent position and deployment using standard angiographic techniques. For optimal results, the entire stenosed arterial segment should be covered by the stent. Fluoroscopic visualisation during stent expansion should be used in order to properly judge the optimum expanded stent diameter as compared to the proximal and distal coronary artery diameter(s). Optimal expansion requires that the stent be in full contact with the artery wall.
- If stent sizing/apposition requires optimisation, readvance the stent system balloon, or another high-pressure balloon catheter of the appropriate size, to the stented area using standard angioplasty techniques.
- Inflate the balloon to the desired pressure while observing under fluoroscopy. Deflate the balloon (refer to product labeling and/or see Balloon Compliance Chart supplied with device).
- If more than one SYNERGY™ stent is needed to cover the lesion and balloon treated area, it is suggested that, to avoid the potential for gap restenosis, the stents be adequately overlapped. To ensure that there are no gaps between stents, the balloon marker bands of the second SYNERGY stent should be positioned inside of the deployed stent prior to expansion.
- Reconfirm stent position and angiographic result. Repeat inflations until optimal stent deployment is achieved, or remove stent delivery system for larger post-dilatation balloon catheter.

#### Removal Procedure and Completion

##### Step Action

- Ensure balloon is fully deflated before delivery system withdrawal.
- Fully open rotating hemostatic valve.
- While maintaining guidewire position and negative pressure on inflation device, withdraw delivery system.
- Repeat angiography to assess the stented area.
- If an adequate expansion has not been obtained, exchange back to the original stent delivery catheter or exchange to another balloon catheter of appropriate balloon diameter to achieve proper stent apposition to the vessel wall.
  - The stent delivery balloon may be used for post-dilatation up to stent diameters indicated on the compliance chart (Table 4).
  - A post-dilatation balloon catheter may be used to expand the stent up to the post-dilatation limits indicated in the following table.

#### Post-Deployment Dilatation of Stented Segments

---

**Precaution:** Do not dilate the stent beyond the limit tabulated below.

---

Nominal Stent Diameter (ID)	Post-Dilatation Limits (ID)*
2.25 mm, 2.50 mm, 2.75 mm	3.50 mm
3.00 mm, 3.50 mm	4.25 mm
4.00 mm	5.75 mm

\*Max Stent Inner Diameter

**Note:** All efforts should be taken to assure that the stent is not underdilated. If the deployed stent size is still inadequate with respect to vessel diameter, or if full contact with the vessel wall is not achieved, a larger post-dilatation balloon catheter may be used to expand the stent further. The balloon should be centered within the stent and should not extend outside of the stented region.

**Note:** Care must be exercised when crossing a newly deployed stent with an intravascular ultrasound (IVUS) catheter, a coronary guidewire, or a balloon catheter to avoid disrupting the stent placement, apposition, geometry, and/or coating. If recrossing with a guidewire, the stented segment should be recrossed carefully with a prolapsed tip to avoid dislodging the stent.

- Complete angiographic confirmation, remove devices, and close vascular access site according to standard practice.

## IN VITRO INFORMATION

### SYNERGY™ Stent System Compliance

Table 4. SYNERGY Compliance Chart

Pressure atm - kPa	Stent I.D. (mm)					
	2.25	2.50	2.75	3.00	3.50	4.00
8 - 814	---	2.35	2.57	2.89	3.30	3.81
9 - 910	2.13	2.42	2.65	2.96	3.40	3.91
10 - 1014	2.19	2.48	2.72	3.02	3.48	3.98
11 - 1117	2.24	2.54	2.79	3.08	3.55	4.06
12 - 1213	2.28	2.59	2.85	3.13	3.61	4.12
13 - 1317	2.31	2.63	2.89	3.17	3.66	4.17
14 - 1420	2.35	2.67	2.93	3.20	3.70	4.22
15 - 1517	2.37	2.70	2.96	3.24	3.74	4.26
16 - 1620	2.40	2.73	3.00	3.27*	3.79*	4.30*
17 - 1724	2.43	2.76	3.03	3.32	3.83	4.36
18 - 1827	2.45*	2.79*	3.06*	3.37	3.87	4.42
19 - 1924	2.48	2.82	3.10	3.43	3.93	4.52
20 - 2027	2.51	2.85	3.13	3.49	3.99	---
21 - 2130	2.54	2.90	3.19	---	---	---
22 - 2227	2.58	2.95	3.23	---	---	---

\*Rated Burst Pressure. DO NOT EXCEED  
Nominal Pressure = 11.0 atm - 1117 kPa

Pressure atm - kPa	Stent O.D. (mm)					
	2.25	2.50	2.75	3.00	3.50	4.00
8 - 814	---	2.51	2.73	3.05	3.46	3.99
9 - 910	2.29	2.58	2.81	3.12	3.56	4.09
10 - 1014	2.35	2.64	2.88	3.18	3.64	4.16
11 - 1117	2.40	2.70	2.95	3.24	3.71	4.24
12 - 1213	2.44	2.75	3.01	3.29	3.77	4.30
13 - 1317	2.47	2.79	3.05	3.33	3.82	4.35
14 - 1420	2.51	2.83	3.09	3.36	3.86	4.40
15 - 1517	2.53	2.86	3.12	3.40	3.90	4.44
16 - 1620	2.56	2.89	3.16	3.43*	3.95*	4.48*
17 - 1724	2.59	2.92	3.19	3.48	3.99	4.54
18 - 1827	2.61*	2.95*	3.22*	3.53	4.03	4.60
19 - 1924	2.64	2.98	3.26	3.59	4.09	4.70
20 - 2027	2.67	3.01	3.29	3.65	4.15	---
21 - 2130	2.70	3.06	3.35	---	---	---
22 - 2227	2.74	3.11	3.39	---	---	---

\*Rated Burst Pressure. DO NOT EXCEED  
Nominal Pressure = 11.0 atm - 1117 kPa

## WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

Viagra is a trademark of Pfizer, Inc.

Xience V is a trademark of Abbott Cardiovascular Systems, Inc.

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Magnetom Trio and syngo are trademarks of Siemens Aktiengesellschaft.
Intera is a trademark of Koninklijke Philips Electronics NV.
© 2014 Boston Scientific Corporation or its affiliates. All rights reserved.
Boston Scientific International S.A.
Boston Scientific (Australia) Pty Ltd