

**Boston
Scientific**

CLINICAL SUMMARY

COMPANION

CAUTION: Federal law restricts this device to sale by or on the order of a physician trained or experienced in device implant and follow-up procedures.

Boston Scientific Corporation acquired Guidant Corporation in April 2006. During our transition period, you may see both the Boston Scientific and Guidant names on product and patient materials. As we work through the transition, we will continue to offer doctors and their patients technologically advanced and high quality medical devices and therapies.

TABLE OF CONTENTS

Clinical Study Populations	1
Summary	1
Observed Adverse Events	2
Prior History	2
Study Background	2
Adverse Event Definitions	3
Deaths	9
Study Design	10
Inclusion/Exclusion Criteria	10
Results	16

CLINICAL STUDY - COMPANION

CLINICAL STUDY POPULATIONS

Guidant CRT-Ds, when compared to OPT alone, have been demonstrated with reasonable assurance, to be safe and effective in significantly reducing: the risk of a composite of all-cause mortality or first hospitalization by 20%, the risk of all-cause mortality by 36%, and heart failure symptoms in patients who have moderate to severe heart failure (NYHA III/IV) including left ventricular dysfunction ($EF \leq 35\%$) and QRS duration ≥ 120 ms and remain symptomatic despite stable, optimal heart failure drug therapy, based on the Guidant sponsored COMPANION clinical study. (Guidant devices were the only devices studied in the COMPANION clinical trial.)

SUMMARY

The COMPANION clinical study was designed to determine whether combined all-cause mortality or first hospitalization in heart failure patients receiving optimal pharmacologic therapy (OPT) can be reduced by combining OPT and either of the following:

- Biventricular pacing therapy alone (CRT-P)
- Biventricular pacing with defibrillation (CRT-D)

All-cause mortality or first hospitalization (time to first event) analyzed from the time of randomization, was the primary endpoint of the study.

Guidant conducted the COMPANION study in part to demonstrate the safety and effectiveness of Guidant CRT-D and CRT-P devices in the COMPANION patient population. Trial objectives included establishing that OPT combined with biventricular pacing with defibrillation [CONTAK CD] is superior to OPT alone in improving exercise performance (Sub-study), reducing combined all-cause mortality or first hospitalization (Primary endpoint), reducing cardiac morbidity (Secondary endpoint) and reducing all-cause mortality (Secondary endpoint).

The COMPANION trial utilized a Steering Committee, Data Safety Monitoring Board (DSMB), and Morbidity and Mortality Committee for study conduct, safety, and event adjudication respectively.

The clinical study began January 20, 2000 and was conducted at 128 centers within the United States.

The COMPANION clinical study was monitored using a sequential design and on November 18, 2002, after review of the data by the Data Safety and Monitoring Board, enrollment in the study was stopped. The CRT-D arm of the trial had reached the target number of events for the combined primary all-cause mortality or first hospitalization endpoint, as well as the secondary all-cause mortality endpoint. All effectiveness follow-ups ended by December 1, 2002.

OBSERVED ADVERSE EVENTS

Prior History

The CONTAK RENEWAL family of CRT-Ds and CONTAK CD devices provide the same defibrillation and cardiac resynchronization therapy (biventricular pacing) and have the same Indications for Use.

Therefore, the Comparison of Medical, Pacing, and Defibrillation Therapies in Heart Failure (COMPANION) clinical trial data (based on CONTAK CD devices) used to support expanding Guidant CRT-D indications to the COMPANION patient population, are applicable to all Guidant/Boston Scientific CRT-Ds including LIVIAN devices.

The primary difference between CONTAK CD devices and CONTAK RENEWAL family of CRT-Ds is that CONTAK CD utilizes an electrically common RV and LV sensing/pacing circuit whereas CONTAK RENEWAL family of CRT-Ds incorporate an independent RV and LV sensing/pacing circuit. Additional clinical analysis was conducted with CONTAK RENEWAL in a European study, to provide confirmation that the independent sensing and pacing capability did not adversely affect the ability of the device to detect ventricular tachyarrhythmias or provide continuous biventricular pacing therapy.

Study Background

The Comparison of Medical Therapy, Pacing, and Defibrillation in Heart Failure (COMPANION) Study was a prospective, open-label, randomized, controlled, multi-center, unblinded study conducted at 128 sites and enrolled a total of 1638 patients, of which 1520 were randomized. Patients were randomly assigned 1:2:2 to receive optimal pharmacological therapy (OPT, 308 patients) or a cardiac resynchronization therapy pacemaker (CRT-P, 617 patients) or a cardiac resynchronization therapy pacemaker with defibrillator (CRT-D, 595 patients). Of the 1520 patients randomized, 903 were randomized to OPT and CRT-D. This summary focuses on

data and analyses for the CRT-D and OPT groups, only, with the exception of the Exercise Performance results, which are based on pooled CRT-D and CRT-P data.

The CRT-D devices (CONTAK CD) in this trial, were approved for commercial distribution via the CONTAK CD study, which provided a reasonable assurance of safety. A similar safety analysis was applied to the COMPANION patient population. The results were consistent with safety measurements obtained in the CONTAK CD trial. See “Data Analysis and Results: CRT-D System Safety” on page 30.

Adverse Event Definitions

Adverse events were defined as any undesirable clinical occurrence, whether it was related to the device or not. Table 1 includes adverse events occurring in the first six months related to the device (pulse generator and leads) and implant procedures (including attempts). Table 2 includes adverse events occurring in the first six months related to patient condition (i.e., worsening heart failure). Adverse events are listed in descending order by total number of patients experiencing the event.

Adverse events related to the device were further reported using two sub-categories based on the nature of the intervention. These events were defined as a *complication* if the event resulted in invasive intervention, loss of significant device function, and death or permanent disability. An *observation* was a device-related adverse event that was resolved non-invasively. Forty-nine percent of CRT-D patients reported a device and/or procedure-related adverse event.

Table 1. Device- and Procedure-Related Adverse Events [Occurring During the First Six Months Post Randomization^a (N = 588)]

	Total Events (Patients)	% Comps (Patients)	% Obs (Patients)
Total Adverse Events	498 (290)	13.1 (77)	43.4 (255)
Post surgical wound discomfort	68 (62)	0.0 (0)	10.5 (62)
Phrenic nerve/diaphragm stimulation	77 (59)	1.4 (8)	9.0 (53)
Brady capture - LV	38 (36)	4.3 (25)	2.2 (13)
Hematoma	37 (34)	0.3 (2)	5.4 (32)
Inappropriate shock above rate cutoff	26 (24)	0.0 (0)	4.1 (24)
Multiple counting - tachy	22 (17)	0.3 (2)	2.9 (17)
Pocket infection	19 (17)	0.5 (3)	2.6 (15)
Dissection, coronary sinus	15 (15)	0.0 (0)	2.6 (15)
Brady capture - atrium	14 (12)	1.5 (9)	0.5 (3)

Table 1. Device- and Procedure-Related Adverse Events [Occurring During the First Six Months Post Randomization^a (N = 588)]

	Total Events (Patients)	% Comps (Patients)	% Obs (Patients)
Inappropriate shock due to oversensing	11 (11)	0.0 (0)	1.9 (11)
Pneumothorax	10 (10)	1.0 (6)	0.7 (4)
Hypotension	10 (9)	0.2 (1)	1.4 (8)
Brady capture - RV	8 (8)	0.9 (5)	0.5 (3)
Physical trauma	8 (8)	0.2 (1)	1.2 (7)
AV Block - heart block, complete	7 (7)	0.2 (1)	1.0 (6)
Pacemaker-mediated tachycardia (PMT)	7 (6)	0.0 (0)	1.0 (6)
Physiological reaction ^b	6 (6)	0.0 (0)	1.0 (6)
Arrhythmia - atrial fibrillation	5 (5)	0.0 (0)	0.9 (5)
Bleeding/fluid accumulation	5 (5)	0.0 (0)	0.9 (5)
Perforation, coronary venous	5 (5)	0.5 (3)	0.3 (2)
Renal failure	5 (5)	0.0 (0)	0.9 (5)
Thrombosis	5 (5)	0.0 (0)	0.9 (5)
Vascular related	5 (5)	0.0 (0)	0.9 (5)
Oversensing - atrium pace sense	4 (4)	0.3 (2)	0.3 (2)
Allergic reaction	3 (3)	0.0 (0)	0.5 (3)
Congestive heart failure	3 (3)	0.0 (0)	0.5 (3)
Nausea (2), Constipation (1)	3 (3)	0.0 (0)	0.5 (3)
High DFTs - tachy	3 (3)	0.2 (1)	0.3 (2)
Oversensing - ventricle rate - tachy	3 (3)	0.2 (1)	0.3 (2)
Respiratory related	3 (3)	0.2 (1)	0.3 (2)
Ventricular tachycardia	3 (3)	0.2 (1)	0.3 (2)
Cardiac tamponade	2 (2)	0.3 (2)	0.0 (0)
Dyspnea (shortness of breath)	2 (2)	0.0 (0)	0.3 (2)
Electrolyte/lab	2 (2)	0.0 (0)	0.3 (2)
Hemorrhage	2 (2)	0.2 (1)	0.2 (1)
Insulation breach suspected	2 (2)	0.3 (2)	0.0 (0)
Migration of device	2 (2)	0.0 (0)	0.3 (2)
Muscle stimulation	2 (2)	0.0 (0)	0.3 (2)
Myocardial infarction	2 (2)	0.0 (0)	0.3 (2)
Numbness	2 (2)	0.0 (0)	0.3 (2)
Perforation, venous	2 (2)	0.0 (0)	0.3 (2)

Table 1. Device- and Procedure-Related Adverse Events [Occurring During the First Six Months Post Randomization^a (N = 588)]

	Total Events (Patients)	% Comps (Patients)	% Obs (Patients)
Phantom shock	2 (2)	0.0 (0)	0.3 (2)
Undersensing - atrium pace sense - brady	2 (2)	0.2 (1)	0.2 (1)
Altered hemodynamic status	1 (1)	0.0 (0)	0.2 (1)
Arrhythmia	1 (1)	0.0 (0)	0.2 (1)
Arrhythmia - sinus tachycardia	1 (1)	0.0 (0)	0.2 (1)
Bruise	1 (1)	0.0 (0)	0.2 (1)
Cardiac arrest	1 (1)	0.2 (1)	0.0 (0)
Change in arrhythmia - SVT	1 (1)	0.0 (0)	0.2 (1)
Change in arrhythmia - brady	1 (1)	0.0 (0)	0.2 (1)
Change in arrhythmia - junctional	1 (1)	0.0 (0)	0.2 (1)
Change in physical status	1 (1)	0.0 (0)	0.2 (1)
Chest pain	1 (1)	0.0 (0)	0.2 (1)
Dizziness, cause undetermined	1 (1)	0.0 (0)	0.2 (1)
Edema	1 (1)	0.0 (0)	0.2 (1)
Fatigue	1 (1)	0.0 (0)	0.2 (1)
Febrile	1 (1)	0.0 (0)	0.2 (1)
Unable to urinate	1 (1)	0.0 (0)	0.2 (1)
Helix related (screw tip), broken or stretched	1 (1)	0.2 (1)	0.0 (0)
Hemoglobin drop	1 (1)	0.2 (1)	0.0 (0)
Hypertension	1 (1)	0.0 (0)	0.2 (1)
Infection	1 (1)	0.2 (1)	0.0 (0)
Insulation breach observed	2 (1)	0.2 (1)	0.0 (0)
Malfunction, memory problem	1 (1)	0.2 (1)	0.0 (0)
Materials unretrieved in body	1 (1)	0.2 (1)	0.0 (0)
Pacemaker mediated tachycardia (PMT)	1 (1)	0.0 (0)	0.2 (1)
Pacemaker syndrome	1 (1)	0.0 (0)	0.2 (1)
Pericardial effusion	1 (1)	0.2 (1)	0.0 (0)
Pericarditis	2 (1)	0.0 (0)	0.2 (1)
Placement difficulty, stylet related	1 (1)	0.2 (1)	0.0 (0)
Pleural effusion	1 (1)	0.2 (1)	0.0 (0)

Table 1. Device- and Procedure-Related Adverse Events [Occurring During the First Six Months Post Randomization^a (N = 588)]

	Total Events (Patients)	% Comps (Patients)	% Obs (Patients)
Pleurisy	2 (1)	0.0 (0)	0.2 (1)
Pocket erosion/extrusion	1 (1)	0.2 (1)	0.0 (0)
Anxiety	1 (1)	0.0 (0)	0.2 (1)
Respiratory arrest	1 (1)	0.2 (1)	0.0 (0)
Ventricular fibrillation	1 (1)	0.0 (0)	0.2 (1)

- a. Observations and complications may not sum to total because some patient may have events in both categories.
- b. Physiological reaction includes: swelling, rash, and/or drainage.

Table 2. Patient-Related Adverse Events [Occurring During the First Six Months Post Randomization^a (N = 588)]

	Total Events (Patients)		% of Patients with Events		Events/Patient Year	
	CRT-D	OPT	CRT-D N = 595 Patients	OPT N = 308 Patients	CRT-D 281 Years	OPT 134 Years
Total Patient Related Adverse Events	1437 (443)	625 (207)	74.5	67.2	5.11 (1437)	4.66 (625)
Cardiovascular Related Events	814 (351)	399 (176)	59.0	57.1	2.90 (814)	2.98 (399)
Congestive heart failure ^b	269 (166)	185 (111)	27.9	36.0	0.96 (269)	1.38 (185)
Chest pain	83 (65)	50 (37)	10.9	12.0	0.30 (83)	0.37 (50)
Supraventricular tachyarrhythmia	69 (56)	11 (11)	9.4	3.6	0.25 (69)	0.08 (11)
Ventricular tachyarrhythmia	66 (51)	16 (15)	8.6	4.9	0.23 (66)	0.12 (16)
Electrolyte/lab	51 (42)	17 (16)	7.1	5.2	0.18 (51)	0.13 (17)
Hypotension	42 (40)	16 (15)	6.7	4.9	0.15 (42)	0.12 (16)
Dizziness, cause undetermined	33 (30)	26 (23)	5.0	7.5	0.12 (33)	0.19 (26)
Renal failure	40 (29)	16 (14)	4.9	4.5	0.14 (40)	0.12 (16)
Fatigue	27 (25)	12 (12)	4.2	3.9	0.10 (27)	0.09 (12)
Bradycardia	32 (30)	5 (5)	5.0	1.6	0.11 (32)	0.04 (5)
Vascular	14 (11)	11 (10)	1.8	3.2	0.05 (14)	0.08 (11)
Syncope	12 (12)	7 (7)	2.0	2.3	0.04 (12)	0.05 (7)
GI bleed	14 (13)	4 (4)	2.2	1.3	0.05 (14)	0.03 (4)
Arrhythmia	12 (10)	6 (6)	1.7	1.9	0.04 (12)	0.04 (6)
Hypertension	12 (9)	6 (5)	1.5	1.6	0.04 (12)	0.04 (6)
Palpitations	9 (9)	3 (3)	1.5	1.0	0.03 (9)	0.02 (3)
Myocardial infarction	7 (7)	4 (4)	1.2	1.3	0.02 (7)	0.03 (4)
Stroke syndrome or CVA	7 (7)	2 (2)	1.2	0.6	0.02 (7)	0.01 (2)
Deep vein thrombosis	4 (4)	0 (0)	0.7	0.0	0.01 (4)	0.00 (0)

Table 2. Patient-Related Adverse Events [Occurring During the First Six Months Post Randomization^a (N = 588)]

	Total Events (Patients)		% of Patients with Events		Events/Patient Year	
	CRT-D	OPT	CRT-D N = 595 Patients	OPT N = 308 Patients	CRT-D 281 Years	OPT 134 Years
Transient ischemic attack (TIA)	3 (3)	1 (1)	0.5	0.3	0.01 (3)	0.01 (1)
Hematuria	3 (3)	0 (0)	0.5	0.0	0.01 (3)	0.00 (0)
Ischemia	2 (2)	1 (1)	0.3	0.3	0.01 (2)	0.01 (1)
Coagulopathy	2 (2)	0 (0)	0.3	0.0	0.01 (2)	0.00 (0)
Bleeding/fluid accumulation	1 (1)	0 (0)	0.2	0.0	0.00 (1)	0.00 (0)
Non-cardiovascular Related Events	623 (293)	226 (119)	49.2	38.6	2.22 (623)	1.69 (226)
Respiratory related ^c	130 (108)	53 (41)	18.2	13.3	0.46 (130)	0.40 (53)
GI ^d	124 (95)	30 (24)	16.0	7.8	0.44 (124)	0.22 (30)
Pain	82 (66)	40 (32)	11.1	10.4	0.29 (82)	0.30 (40)
Physiological reaction ^e	76 (61)	20 (18)	10.3	5.8	0.27 (76)	0.15 (20)
Infection	54 (37)	18 (15)	6.2	4.9	0.19 (54)	0.13 (18)
Endocrine	41 (35)	16 (14)	5.9	4.5	0.15 (41)	0.12 (16)
Psychological effects	24 (19)	13 (12)	3.2	3.9	0.09 (24)	0.10 (13)
Change in physical status	20 (18)	9 (9)	3.0	2.9	0.07 (20)	0.07 (9)
Physical trauma	26 (22)	4 (4)	3.7	1.3	0.09 (26)	0.03 (4)
Neurologic	14 (14)	6 (6)	2.4	1.9	0.05 (14)	0.04 (6)
Genitourinary	9 (7)	5 (4)	1.2	1.3	0.03 (9)	0.04 (5)
Cancer, other	5 (5)	6 (5)	0.8	1.6	0.02 (5)	0.04 (6)
Febrile	7 (7)	0 (0)	1.2	0.0	0.02 (7)	0.00 (0)
Respiratory failure	4 (4)	1 (1)	0.7	0.3	0.01 (4)	0.01 (1)
Tumors, growths	1 (1)	2 (2)	0.2	0.6	0.00 (1)	0.01 (2)
Ulceration	2 (1)	2 (2)	0.2	0.6	0.01 (2)	0.01 (2)

Table 2. Patient-Related Adverse Events [Occurring During the First Six Months Post Randomization^a (N = 588)]

	Total Events (Patients)		% of Patients with Events		Events/Patient Year	
	CRT-D	OPT	CRT-D N = 595 Patients	OPT N = 308 Patients	CRT-D 281 Years	OPT 134 Years
Diabetes complications	2 (2)	0 (0)	0.3	0.0	0.01 (2)	0.00 (0)
Pulmonary embolism	1 (1)	1 (1)	0.2	0.3	0.00 (1)	0.01 (1)
Pneumonia (respiratory infection)	1 (1)	0 (0)	0.2	0.0	0.00 (1)	0.00 (0)

- a. Observations and complications may not sum to total because some patients may have events in both categories.
- b. Congestive heart failure includes: congestive heart failure, dyspnea, volume overload, edema, pulmonary edema, change in drug therapy.
- c. The most frequent three events in this category were: upper respiratory infection, bronchitis, and influenza.
- d. The most frequent three events in this category were: nausea, diarrhea, and abdominal pain.
- e. The most frequent three events in this category were: swelling, rash, and weight gain.

Deaths

There were a total of 182 deaths (77 OPT, 105 CRT-D) that occurred during the trial and recorded through November 30, 2002. Table 3 presents cause of death stratified by treatment group.

Table 3. CRT-D and OPT Cause of Death

Cause of Death	OPT Arm (N = 308)	CRT-D Arm (N = 595)	Total (N = 903)
Cardiac	58 (18.8%)	76 (12.8%)	134 (14.8%)
Vascular	0 (0.0%)	3 (0.5%)	3 (0.3%)
Non-Cardiac	11 (3.6%)	21 (3.5%)	32 (3.5%)
Unknown/ Unclassified	8 (2.6%)	5 (0.8%)	13 (1.4%)
Total Deaths	77 (25.0%)	105 (17.6%)	182 (20.2%)

NOTE: After the study was stopped in November 2002, follow-up for safety continued for approximately one more year on 151 OPT and 449 CRT-D patients with the final data cut-off on November 26, 2003. During this post-trial follow-up period, an additional 54 deaths were reported, consisting of 14/151 (9.3%) OPT

patients and 40/449 (8.9%) CRT-D patients.

The mortality rates are approximately equal during this post-trial follow-up period. This may be because CRT devices were made available to OPT patients. Thus, most patients were receiving the same therapy during this interval.

STUDY DESIGN

The COMPANION study design and study results have been previously described in the medical literature.^{1,2}

The COMPANION study was a prospective, randomized [1:2:2 to OPT, CRT-P (delivered by the CONTAK TR device), or CRT-D (delivered by the CONTAK CD device)], controlled, multi-center study. Both of these devices became commercially available during the course of the study.

Randomization was stratified by centers and by beta-blocker use to assure proper balance between the treatment groups within each center. Each randomized patient remained counted as a member of the original randomization assignment (intention-to-treat) regardless of subsequent crossover or protocol adherence.

Eligible patients were also enrolled in a sub-study designed to measure improvement in exercise performance in patients randomized to CRT (CRT-P and CRT-D pooled data) therapy compared to OPT.

Inclusion/Exclusion Criteria

The study population consisted of patients with moderate to severe heart failure, New York Heart Association Classification III or IV, left ventricular ejection fraction \leq 35%, and QRS width \geq 120 ms due to ischemic or non-ischemic cardiomyopathy.

All patients were required to have been treated with a stable dose of beta-blocker, ACE inhibitor or ARB, diuretic, and aldosterone antagonist. A stable dose was defined as 30 days for all drugs except beta-blocker, which required 90 days stabilization from last up titration prior to randomization. Diuretic dosage could be adjusted any time by the investigator using medical discretion.

1. Bristow MR, Feldman AM, Saxon LA. Heart failure management using implantable devices for ventricular resynchronization: Comparison of Medical Therapy, Pacing, and Defibrillation in Chronic Heart Failure (COMPANION) trial. *J Card Fail.* 2000;6(3):276-285.
2. Bristow MR, Saxon LA, Boehmer J, et al. Cardiac resynchronization therapy with or without an implantable defibrillator in advanced chronic heart failure. *N Engl J Med.* 2004;350:2140-2150.

Patients enrolled in the study were required to meet the following inclusion criteria:

- Moderate or severe heart failure, defined as symptomatic heart failure for at least six months with NYHA class III or IV symptoms at the time of enrollment, and at least one of the following events in the previous 12 months:
 - Hospitalization for heart failure management
 - Outpatient visit in which intravenous (IV) inotropes or vasoactive infusion were administered continuously for at least 4 hours
 - Emergency room visit of at least twelve hours duration in which IV heart failure medications were administered (including diuretics)
- QRS \geq 120 ms and PR interval $>$ 150 ms from any two leads of a 12-lead ECG
- Left ventricular ejection fraction \leq 35%
- Left ventricular end diastolic dimension \geq 60 mm (required only if LVEF measured by echo) or $>$ 3.0 cm/m² (The cm/m² is calculated by LVEDD [in cm] divided by BSA [body surface area])
- Age \geq 18 years
- Optimal pharmacologic therapy for heart failure (beta blocker, ACE inhibitor, diuretics, and spironolactone)

Additional eligibility criteria for the Exercise Performance sub-study:

- Understand the nature of the sub-study and provide informed consent
- Have been enrolled at a participating sub-study investigational center
- Have no neuromuscular or vascular disability that prevents normal walking (e.g., intermittent claudication, arthritis, residual stroke weakness)
- Have no history of angina during previous exercise testing
- FEV₁/FVC \geq 50%
- 150m \leq Six-minute walk distance \leq 425 m

- Baseline Peak $\text{VO}_2 < 22 \text{ ml/kg/min}$
- Have no cardiac disabilities that would ordinarily contraindicate exercise testing:
 - Changing pattern on the ECG
 - Changing pattern of chest discomfort
 - Decompensated heart failure
 - Uncontrolled arrhythmias

Patients were excluded from the investigation if they met any of the following criteria:

- Unable or unwilling to undergo device implant and follow-up testing
- Patients with a hypersensitivity to 0.7 mg nominal dose of dexamethasone acetate
- Meet the general indications for an implantable cardioverter defibrillator
- Surgically uncorrected primary valvular disease
- Meet the general indications for antibradycardia pacing
- Coronary artery disease (CAD) in which surgical or percutaneous correction is recent (within 60 days of randomization)
- Expected to receive a heart transplant in the next six months
- Women who are pregnant or not using medically acceptable birth control
- Chronic, medically refractory atrial tachyarrhythmias
- Hypertrophic obstructive cardiomyopathy
- Unexplained syncope
- Amyloid disease
- Myocardial infarction within 60 days of randomization

- Hospitalization for heart failure or IV inotropic or vasoactive therapy in excess of 4 hours in the 30 days prior to enrollment
- History of non-compliance with oral heart failure therapy
- Involved in any other investigational studies
- Progressive or unstable angina
- Life expectancy < 6 months due to any other medical conditions
- Uncontrolled blood pressure: Systolic BP > 160 mm Hg or < 85 mm Hg or diastolic BP > 90 mm Hg

Endpoints

This summary focuses on the CRT-D vs. OPT contrast, providing evidence of safety and effectiveness for Guidant CRT-Ds in the COMPANION patient population¹. The clinical data and analyses herein address the following study endpoints for all patients randomized to CRT-D and OPT only, unless otherwise stated as a sub-study measurement (where CRT-D and CRT-P data were pooled for exercise performance):

Primary Endpoint

The primary endpoint was a composite consisting of all-cause mortality or first hospitalization (time to first event) as analyzed from the date of randomization on an intention-to-treat basis. The study was designed to demonstrate a 25% relative reduction with CRT-D when compared to an estimated 40% annual rate in the OPT cohort. All-cause mortality was defined as death from any cause. Hospitalization is defined below:

Qualifying Duration for Hospitalization The intent behind hospitalization was to capture hospitalizations that were of sufficient duration to enter into a composite with all-cause mortality. Thus, hospitalization was defined as care provided at a hospital in which hospital admission and discharge occurred on separate dates. Patients excluded from this definition were those who received care at a hospital, but were discharged on the same day as admission. In addition to hospitalizations, the use of intravenous inotropes or vasoactive agents for a duration of greater than

1. Guidant CRT-Ps are already approved for use in the COMPANION patient population, P030005, approved 01/26/04.

four hours was also considered to be of significant importance to be treated as an instance of hospitalization.

Hospitalizations Related to the Implant Procedure Hospitalizations associated with device implant (initial and reattempted for unsuccessful initial implant) were not considered to be an event for evaluating the primary endpoint. Similarly, hospitalizations associated with elective implant of devices (i.e., absence of an electrophysiological indication or an ongoing hospitalization requiring intravenous therapy) in the OPT cohort also were not considered to be a primary endpoint event. Surgical revisions of a previous implanted system did count as a primary endpoint event if the revision was of a sufficient duration to result in different admission and discharge dates. Table 4 summarizes the criteria for determining which hospitalization events were considered as a primary endpoint event.

Table 4. Hospitalizations Contributing to Primary Endpoint

Event Description	CRT-D	OPT
Initial implant/reattempts	No	No
Surgical revisions of system	Yes ^a	Yes ^a
Hospitalization with no calendar date change	No	No
Hospitalization with a calendar date change	Yes	Yes
IV inotrope and/or vasoactive drug use > 4 hours	Yes	Yes

a. If calendar date change.

Secondary Endpoints

All-cause mortality The all-cause mortality (death from any cause) endpoint was designed to show a 25% reduction in mortality in the CRT-D arm from an OPT annual mortality rate of 24%. Difference in mortality was determined by contrasting patients randomized to CRT-D in addition to OPT versus patients randomized to OPT alone.

Cardiac morbidity The hospitalization component of the primary endpoint included non-cardiac events that may not be impacted by CRT-D. The cardiac morbidity endpoint was unique to the COMPANION study. It is a more specific outcome measure intended to determine whether CRT-D when compared to OPT would reduce the type of events that are pertinent to a hospitalization for heart failure.

Cardiac morbidity was defined as the occurrence of one or more of the following events:

- Worsening heart failure resulting in use of intravenous vasoactive or inotropic therapy exceeding four hours
- Mechanical respiratory or cardiac support
- Any cardiac surgery, including heart transplant
- Resuscitated cardiac arrest or sustained ventricular tachycardia requiring intervention (e.g., chest thump, external cardioversion, or external defibrillation)
- Hospitalization for acute decompensation of heart failure
- Hospitalization that results in death from cardiac causes
- Significant device-related events resulting in:
 - Permanent disability
 - Hospitalization for pending death or permanent disability

Safety

CRT-D system-related complication-free rate is determined by measuring complications related to any of the implanted components or their associated implant procedure in those patients who were successfully implanted with the CRT-D system.

NOTE: *During the course of the COMPANION clinical study, the EASYTRAK Coronary Venous pace/sense lead was established as safe and effective in a separate clinical study and was approved for commercial distribution (P010012, 05/02/02). Refer to the commercially available EASYTRAK Coronary Venous pace/sense lead labeling for clinical safety and performance characteristics.*

Sub-study Primary Endpoint and Additional Tertiary Endpoints

Exercise performance The co-primary endpoint, which consists of Peak VO₂ and Six-Minute Walk, is designed to demonstrate improvement in exercise performance with CRT (CONTAK TR and CONTAK CD pooled data) compared to OPT at six months post-baseline.

Additional tertiary endpoints included Quality of Life as measured by the Minnesota Living with Heart Failure Questionnaire[®] and NYHA Class.

Follow-up Schedule

Enrollment	Initial assessment of patient eligibility; taking of patient history.
Baseline Screening	Special testing ^a
Randomization	Randomization status (OPT, CRT-P, or CRT-D) was assigned.
Implant (CRT-P or CRT-D arm)	Implant of investigational devices and acute device testing for those randomized to a CRT therapy arm.
Routine Follow-up	Routine evaluation of device function and patient condition at pre-discharge, one week, and one month.
Three- and six-month Visits	Evaluation of randomized therapy with Special Testing ^a and device function at three and six months after the Post-Recovery Visit.
Quarterly Visits	After the six-month visit, patients were seen for routine evaluation of device function and patient condition.

a. Special Testing included a Symptom-Limited Treadmill Test with measurement of oxygen uptake (Peak VO₂), a Six-Minute Walk, Quality of Life (QOL) questionnaire and New York Heart Association Classification.

RESULTS**Demographic Data**

All baseline patient characteristics are presented in Table 5.

Table 5. Characteristics of Patient Population for COMPANION (OPT and CRT-D)

Characteristic		OPT (N = 308)	CRT-D (N = 595)	P-value
Age (years)	Mean ± SD	66.7 ± 10.7	65.6 ± 11.2	0.14
Gender [N (%)]	Female	97 (31.4)	194 (32.6)	0.73
	Male	211 (68.5)	401 (67.3)	
NYHA Classification [N (%)]	Class III	253 (82.1)	512 (86.1)	0.12
	Class IV	55 (17.8)	83 (13.9)	

Table 5. Characteristics of Patient Population for COMPANION (OPT and CRT-D)

Characteristic		OPT (N = 308)	CRT-D (N = 595)	P-value
Ischemic Etiology (%)	Ischemic	58.7	54.6	0.13
	Non-ischemic	41.3	45.4	
LVEF (%)	Mean ± SD	22.8 ± 7.2	22.5 ± 6.8	0.47
Resting Heart Rate (bpm)	Mean ± SD	72 ± 12	73 ± 13	0.37
QRS Width (ms)	Mean ± SD	156 ± 24	159 ± 24	0.09
Conduction Abnormality (%)	LBBB	69.8	72.9	0.21
	Non-specific	21.4	16.8	
	RBBB	8.77	10.2	
Duration of Heart Failure (years)	Mean ± SD	4.86 ± 4.41	4.44 ± 3.83	0.43
Heart Failure Medications [(%)]	Diuretic	94.4	96.6	0.12
	ACE inhibitor or ARB	88.6	89.6	0.66
	Beta Blockers	66.2	67.6	0.69
	Aldosterone Antagonist	54.8	55.1	0.94
	Digoxin	67.2	70.9	0.25

Patient Accountability and Follow-up Duration

The COMPANION study enrolled 1638 patients, with 1520 patients randomized to a treatment group and one hundred eighteen patients (118) not randomized due to changes in patient condition or consent between time of enrollment and time of randomization, such that the inclusion criteria were no longer satisfied. Of the 1520 patients, 595 were randomized to CRT-D with a mean follow-up of 1.3 years and 308 were randomized to OPT with a mean follow-up of 1.1 years. Figure 1 provides an overview of patient enrollment.

Table 6 gives a summary (by treatment group) of patient disposition over time through 12 months after randomization. This does not account for patients that had a hospitalization or death event that contributed to the primary endpoint or secondary endpoint of all-cause mortality. For events contributing to the primary

endpoint or the secondary endpoint of all-cause mortality, please refer to Figure 2 on page 20 and Figure 3 on page 21.

Table 6. Patient Follow-up Disposition 12 Months Post Randomization

	CRT-D				OPT			
	# of With-drawn Patients	# of Deceased Patients	(N = 595) # Reached end of study (Nov. 30, 2002)	# of Active Patients at end of time interval	# of With-drawn Patients	# of Deceased Patients	(N = 308) # Reached end of study (Nov. 30, 2002)	# of Active Patients at end of time interval
1 Day - 7 Days	4	3	0	588	6	0	0	302
7 Days - 1 Month	4	3	5	576	10	3	1	288
1 Month - 3 Months	4	15	6	551	11	11	1	265
3 Months - 9 Months	12	28	49	462	26	22	29	188
9 Months - 12 Months	1	12	35	414	11	11	19	147




Figure 1. Study Patient Enrollment and Randomization for CRT-D and OPT.

Event Contributing to Primary Endpoint and Secondary Endpoint of All-cause Mortality

A total of 903 COMPANION patients in the CRT-D (595) and OPT (308) groups were eligible for the primary endpoint. Figure 2 provides patient randomization and status for the primary endpoint and Figure 3 provides patient randomization and status for the secondary mortality endpoint.




Figure 2. CRT-D and OPT Patient Randomization for Primary Endpoint.




Figure 3. CRT-D and OPT Patient Randomization for Mortality Endpoint.

Data Analysis and Results for Primary Endpoint and Secondary All-Cause Mortality Endpoint

Sequential Monitoring

The COMPANION DSMB met approximately every six months to review the trial's progress and to review the safety and effectiveness data collected. An "O'Brien-Fleming" type boundary as implemented by Lan and DeMets was used in monitoring the trial. The Group sequential procedure ensured that the total alpha spent across repeated analyses did not exceed the total type I error, in this case $\alpha = 0.03$.

On November 18, 2002 the DSMB reviewed the study progress for the final time. The CRT-D arm of the Study had reached the target number of events for both the combined mortality and hospitalization endpoint as well as the all-cause mortality endpoint prompting the DSMB to recommend to the Steering Committee that enrollment be stopped. All effectiveness follow-ups ended on December 1, 2002.

Primary Endpoint: All-cause Mortality or First Hospitalization

The Kaplan-Meier curves illustrating the time to all-cause mortality or first hospitalization are shown in Figure 4. There were 216 primary endpoint events observed in the OPT arm and 390 in the CRT-D arm ($p = 0.010$; $p = 0.011$ after adjustment for interim analyses). The median time to first event was 209 days in the OPT group and 269 days in the CRT-D group. The annual event rates for OPT and CRT-D, respectively, were 68.0% and 55.9%, with a hazard ratio of 0.80; 95% CI (0.68, 0.95). This result demonstrated that CRT-D significantly reduced the relative risk of all-cause mortality or first hospitalization by 20% when compared to OPT alone.




Figure 4. Primary Endpoint: All-cause Mortality or First Hospitalization.

In addition to the hazard ratio, point estimates of risk reduction were also calculated (Table 7). These estimates will vary with time from the true treatment effect, and thus should be interpreted with caution.

Table 7. Primary Endpoint [Risk Reduction Point Estimates] (Overall Hazard Ratio = 0.80; p = 0.010)

	% Failure		Absolute Risk Reduction	Relative Risk Reduction
	OPT	CRT-D		
6 months	44.9% (38.9%, 50.3%)	42.9%(38.7%, 46.7%)	2.0%	4.5%
12 months	68.0% (61.7%, 73.2%)	55.9% (51.6%, 59.8%)	12.1%	17.8%
18 months	77.8% (71.6%, 82.7%)	69.0% (64.5%, 73.1%)	8.8%	11.3%

Secondary Endpoints

All-cause Mortality Deaths from any cause were reported in 77 patients randomized to OPT and 105 patients randomized to CRT-D (p = 0.003, p = 0.004 after adjusting for interim analyses). The Kaplan-Meier curves are depicted in Figure 5. These numbers correspond to an annual mortality rate of 19% in the OPT arm and 12% in the CRT-D arm, with a hazard ratio of 0.64, 95% CI (0.48, 0.86). These results demonstrated that CRT-D was associated with a 36% reduction in the risk of all-cause mortality when compared to OPT alone.




Figure 5. Secondary Endpoint: All-cause Mortality.

In addition to the hazard ratio, point estimates of risk reduction were also calculated (Table 8). These estimates will vary with time from the true treatment effect, and thus should be interpreted with caution.

Table 8. Mortality Endpoint Risk Reduction Point Estimates (Overall Hazard Ratio = 0.64; $p = 0.003$)

	% Failure		Absolute Risk Reduction	Relative Risk Reduction
	OPT	CRT-D		
6 months	9.0% (5.7%, 12.2%)	7.3% (5.1%, 9.3%)	1.7%	18.9%
12 months	18.9% (14.1%, 23.5%)	12.1% (9.3%, 14.8%)	6.8%	36.0%
18 months	28.4% (22.3%, 34.1%)	18.0% (14.4%, 21.5%)	10.4%	36.6%

Results for Secondary Cardiac Morbidity Endpoint

As previously mentioned in the Cardiac Morbidity section on page 14, cardiac morbid events were reported during hospitalizations.

During a hospitalization more than one of the pre-specified cardiac morbid events could occur. The Anderson-Gill extension to the Cox proportional hazard model was used to analyze time to multiple cardiac morbid events. Caution must be used in interpreting p-values in this analysis because this analysis does not account for the competing risk of death.

In Figure 6, the frequency and duration of cardiac morbid events are illustrated. CRT-D was associated with a 36% reduction ($p < 0.0001$) in the proportion of patients with at least one event, a 52% reduction ($p < 0.0001$) in events on an annual basis, and a 41% reduction ($p < 0.0001$) in the hospital duration on an annual basis. These reductions are primarily due to the reduction of hospitalizations for acute decompensation of heart failure, worsening heart failure resulting in IV inotrope or vasoactive therapy > 4 hours (during an inpatient hospitalization) and cardiac surgery (including percutaneous intervention), as shown in Figure 7.




Figure 6. Secondary Endpoint of Cardiac Morbidity.

Caution must be used in interpreting p-values; analysis does not account for competing risk of death.




Figure 7. Cardiac Morbidity by Major Component.

For a given cardiac hospitalization, patients may have events in more than one category, and if there are multiple occurrences in a single category, then only the first occurrence was counted.

Other Analyses: Implant Disposition

Table 9 identifies the number of initial and subsequent implant procedures attempted in patients randomized to CRT-D and the rate of success for each additional implant procedure. There were 81 CRT-D patients that had an unsuccessful initial implant for the CRT-D system. Fifty (50) of these patients had a second implant procedure, of which 33 were successful and 17 were unsuccessful. Three patients had a third implant procedure, of which one was successful. Therefore, there were 541 patients implanted with the CRT-D system.

Table 9. CRT-D System Implant Disposition

		Attempt successful	Failed implant	Reattempt not done after this procedure
Initial implants	588 (98.8%)	507 (85.0%)	81 (14.0%)	31 (5.2%)
First reattempt	50 (8.4%)	33 (5.5%)	17 (2.9%)	14 (2.3%)
Second reattempt	3 (0.5%)	1 (0.2%)	2 (0.3%)	2 (0.34%)

Additional Outcome Measures

First Heart Failure Hospitalizations An additional outcome that was not pre-specified in the protocol provides further insight into the results observed in the composite primary endpoint. This post-hoc analysis was conducted using cause-specific hospitalizations as adjudicated by the morbidity and mortality committee and therefore should be interpreted with caution.

The outcome of all-cause mortality or first heart failure hospitalization was analyzed on an intention-to-treat basis and time to first event.

Hospitalizations were defined per the following:

- Care provided at a hospital for any reason in which the duration is associated with a date change, or
- Use of intravenous inotropes and/or vasoactive drugs for a duration > 4 hours (inpatient or outpatient).

NOTE: Hospitalizations associated with a device implant attempt or re-attempt are excluded.

Those contributing to the heart failure hospitalization outcome were required by the Morbidity and Mortality committee to meet at least one of the following additional criteria:

- IV diuretics
- IV inotrope/vasoactive therapy
- Other parenteral therapy for the treatment of heart failure
- Significant alterations in oral therapy for the treatment of heart failure

All-cause Mortality or First Heart Failure Hospitalization The Kaplan-Meier curves for all-cause mortality or first heart failure hospitalization is shown in Figure 8. OPT and CRT-D had annual event rates of 45% and 29%, respectively with a hazard ratio of 0.60, 95% CI (0.49-0.75), $p < 0.001$. Therefore, CRT-D was associated with a 40% relative reduction in the risk of all-cause mortality or first heart-failure hospitalization when compared to OPT alone.




Figure 8. Additional Outcome: All-cause Mortality or First Heart Failure Hospitalization.

Disposition of Hospitalization Implantation of the CRT-D system generally requires hospitalization. To differentiate between the hospitalization required to implant the system and those hospitalizations that occurred after the system was implanted, the following terms are used:

- **Implant hospitalization:** The elective hospitalization associated with either the implant procedure or a repeat implant procedure if the initial procedure was unsuccessful.
- **All other hospitalizations:** Patients who required a revision for an implanted system (e.g., lead dislodgment or infection) were included in this category as were hospitalizations for non-elective device related implants.

The hospitalizations analysis illustrated in Figure 9 and hospitalization days analysis in Figure 10 depicts hospitalization data stratified by *implant* and *non-elective hospitalizations*. This analysis was on an intention-to-treat basis and includes patients who underwent an attempted implant procedure. Patients randomized to CRT-D had a follow-up duration approximately 30% longer than OPT patients. Thus, hospitalization data are normalized per patient-year of follow-up. An additional comparison of hospitalization days for heart failure hospitalizations is shown in Figure 11.

NOTE: CRT-D was associated with a reduction in all-cause mortality and therefore there is a competing risk for hospitalizations. This data should be interpreted with caution.




Figure 9. Hospitalizations/Patient year.




Figure 10. Hospitalization Days/Patient-Year: Hospitalizations.




Figure 11. Hospitalization Days/Patient-Year: Heart Failure Hospitalizations.

Data Analysis and Results: CRT-D System Safety

The system-related complication-free rate analysis was not a predefined endpoint in the protocol. The intent of this analysis is to provide reasonable assurance of safety of the CONTAK CD system in this patient population. The system-related complication-free rate was defined over a six-month follow-up period as the proportion of patients who are free of complications attributed to:

- Any implanted component (e.g, pulse generator, coronary venous lead, right atrial pace/sense lead, cardioversion/defibrillation lead)
- The surgical procedure required to implant the CRT-D system

In the COMPANION study, this analysis was performed on an intention-to-treat basis and also extends to those patients who underwent an implant procedure but did not ultimately receive a device. Of the 595 patients analyzed, 522 (87.7%) were free of system-related complications.

Of the 73 (12.3%) patients who experienced a system-related complication, the most common were loss of left ventricular capture (25 patients, 4.2%), loss of right atrial capture (9 patients, 1.5%), and phrenic nerve/diaphragmatic stimulation (8 patients, 1.3%).

When analyzed on a time-to-event basis, the system-related complication-free rate was 87.7%. The safety performance of the CONTAK CD system compares favorably with the safety performance observed in the prior CONTAK CD study (P010012, May 2, 2002).

Data Analysis and Results for COMPANION Sub-study

The Exercise Performance Sub-study consisted of:

CRT Effectiveness: Primary Co-primary endpoint consisting of Peak VO_2 derived from a symptom-limited exercise test and Six-Minute Walk, with CRT results pooled from the CONTAK TR and CONTAK CD arms.

Effectiveness was determined by assessing both Peak VO_2 and Six-Minute Walk distance improvements with CRT compared to OPT. Prospectively, success was defined as occurring if:

- Peak VO_2 improved ≥ 0.7 ml/kg/min ($p < 0.05$) and 6 MWD improvement resulted in $p < 0.10$, or
- Peak VO_2 improved ≥ 0.5 ml/kg/min ($p < 0.10$) and 6 MWD improvement resulted in $p < 0.05$.

Additional: Quality of Life as measured by the Minnesota Living with Heart Failure Questionnaire® and NYHA Class.

Patient Accountability




Figure 12. Enrollment and Follow-up of Randomized Patients.

Baseline Characteristics

Table 10. Characteristics of Patient Population

Characteristic		CRT (N = 318)	OPT (N = 87)	P-value ^a
Age (years)	Mean ± SD	62.1 ± 11.8	63.1 ± 10.6	0.48
	Range	32.0–86.0	27.0–85.0	
Gender [N (%)]	Female	109 (34.3)	24 (27.6)	0.24
	Male	209 (65.7)	63 (72.4)	
NYHA Classification [N (%)]	III	294 (92.5)	79 (90.8)	0.61
	IV	24 (7.5)	8 (9.2)	
Ischemic Etiology	Ischemic	141 (44.3)	42 (48.3)	0.51
	Non-ischemic	177 (55.7)	45 (51.7)	
LVEF (%)	Mean ± SD	22.5 ± 6.9	22.2 ± 8.0	0.79
	Range	5.0–35.0	5.0–35.0	
Resting Heart Rate (bpm)	Mean ± SD	73.1 ± 12.8	73.5 ± 11.5	0.78
	Range	46.0–122.0	54.0–103.0	
QRS Width (ms)	Mean ± SD	159.2 ± 25.0	155.7 ± 25.8	0.26
	Range	120.0–276.0	120.0–224.0	
LBBB/NSIVCD (%)	LBBB	230 (72.3)	62 (71.3)	0.60
	Nonspecific	54 (17.0)	18 (20.7)	
	RBBB	34 (10.7)	7 (8.0)	
Peak VO ₂ (ml/kg/min)	Mean ± SD	12.7 ± 3.3	12.4 ± 3.3	0.42
	Range	3.0–21.2	4.8–21.5	
Six-Minute Walk Distance (m)	Mean ± SD	292.4 ± 65.5	291.6 ± 70.5	0.92
	Range	152.0–411.5	162.4–414.0	
Quality of Life Score (points)	Mean ± SD	59.8 ± 23.1	55.4 ± 23.3	0.12
	Range	0.0–105.0	0.0–97.0	
Heart Failure Medications [N (%)]	Diuretic	300 (94.3)	82 (94.3)	0.98
	ACE Inhibitor or ARB	286 (89.9)	82 (94.3)	0.22
	Beta Blockers	240 (75.5)	60 (69.0)	0.22
	Aldosterone Antagonist	178 (56.0)	51 (58.6)	0.66
	Digoxin	239 (75.2)	65 (74.7)	0.93

a. Continuous data were analyzed using a two-tailed t-test procedure, and categorical data were analyzed using a two-tailed chi-square procedure. A p-value < 0.05 is considered significant.

CRT Effectiveness

Peak VO₂ Peak VO₂ was determined from a standardized protocol for exercise testing as a means of measuring a patient's capacity for performing physical activity.




Figure 13. Maximal Oxygen Consumption Results.

Table 11. Maximal Oxygen Consumption Results

Peak VO ₂ (ml/kg/min)	CRT		OPT		P-value ^a
	N	Mean \pm S.E.	N	Mean \pm S.E.	
Δ at 3 months	247	1.1 \pm 0.2	52	0.3 \pm 0.4	0.026
Δ at 6 months	230	1.2 \pm 0.2	46	0.6 \pm 0.4	0.074

a. P-values obtained using one-tailed longitudinal analysis methods.

The longitudinal analysis was performed on all available data. The percentages of missing data at the six-month endpoints for Peak VO₂ and Six-Minute Walk were 36 percent and 28 percent for the CRT arm and 47 percent and 34 percent for the OPT arm. The longitudinal analysis performed is most appropriate when missing data occurs at the percentages found in this trial.

Six-Minute Walk The Six-Minute Walk test is a measure of a patient’s ability to sustain exercise during an activity similar to that which a patient may typically perform on a daily basis. For this test, patients are instructed to walk as far as possible in 6 minutes in a level corridor.




Figure 14. Change in Six-Minute Walk.

Table 12. Change in Six-Minute Walk

Six-Minute Walk (m)	CRT		OPT		P-value ^a
	N	Mean ± S.E.	N	Mean ± S.E.	
Δ at 3 months	274	37 ± 5	63	13 ± 10	0.016
Δ at 6 months	260	41 ± 5	57	17 ± 10	0.017

a. P-values obtained using one-tailed longitudinal analysis methods.

NYHA Class The determination for New York Heart Association (NYHA) Class is based on mutual assessment, by the patient and physician, of the patient's heart failure symptoms both at rest and while performing ordinary physical activity.




Figure 15. Change in NYHA.

Table 13. Change in NYHA

NYHA Classification	Change	CRT		OPT		P-value ^a
		N	Patients	N	Patients	
3 months	Improve 2 Classes	294	22 (7.5%)	69	3 (4.4%)	< 0.01
	Improve 1 Class		142 (48.3%)		13 (18.8%)	
	No Change		122 (41.5%)		48 (69.6%)	
	Worsen 1 Class		8 (2.7%)		5 (7.3%)	
6 months	Improve 2 Classes	291	20 (6.9%)	65	2 (3.1%)	0.032
	Improve 1 Class		149 (51.2%)		28 (43.1%)	
	No Change		118 (40.6%)		34 (52.3%)	
	Worsen 1 Class		4 (1.4%)		1 (1.5%)	

a. P-values are not adjusted for multiplicity and were obtained using a one-tailed Mantel-Haenszel chi-square method.

Quality of Life Quality of Life (QOL) was assessed using the 21-question Minnesota Living with Heart Failure questionnaire. Each question, answered by the patient, is ranked on a scale ranging from 0 to 5. A lower total score indicates an improved quality of life.




Figure 16. Quality of Life Score.

Table 14. Quality of Life Score

Quality of Life (points)	CRT		OPT		P-value ^a
	N	Mean ± S.E. (95% CI)	N	Mean ± S.E. (95% CI)	
Δ at 3 months	289	23 ± 1 (20.1, 25.7)	72	6 ± 3 (0.6, 11.3)	< 0.001
Δ at 6 months	279	23 ± 1 (19.7, 25.4)	66	10 ± 3 (4.2, 15.2)	< 0.001

a. P-values are not adjusted for multiplicity and were obtained using one-tailed longitudinal analysis methods.

Additional Functional Capacity Data

In addition to the Exercise Performance sub-study, functional capacity was evaluated by means of NYHA Class, six-minute walk distance, and Minnesota Living with Heart Failure Questionnaire[®] QOL for the all patients randomized to OPT and CRT-D through 6-months of follow up.

As shown in Table 15, NYHA Class, six-minute walk distance, and QOL scores were significantly improved in the CRT-D group compared to the OPT group at 3 and

6 months. These findings are similar to those presented in the exercise performance sub-study and previous cardiac resynchronization therapy trials.

Table 15. Changes in Six-Minute Walk, QOL and NYHA

Six Minute Walk Distance	CRT-D		OPT		P-value ^a
	N	Mean ± SD	N	Mean ± SD	
Δ at 3 months	420	42 ± 98	172	8 ± 82	< 0.0001
Δ at 6 months	377	45 ± 98	141	2 ± 92	< 0.0001
QOL	N	Mean ± SD	N	Mean ± SD	
Δ at 3 months	514	-24 ± 28	243	-8 ± 21	< 0.0001
Δ at 6 months	479	-23 ± 28	207	-12 ± 23	< 0.0001
NYHA	N	% Improved	N	% Improved	
Δ at 3 months	543	55	242	24	< 0.0001
Δ at 6 months	498	57	199	38	< 0.0001

a. P-values are not adjusted for multiplicity and were obtained using t-tests for continuous data and chi-square for categorical data.

Boston Scientific



Manufacturer

Boston Scientific

4100 Hamline Avenue North
St. Paul, MN 55112-5798 USA

www.bostonscientific.com

1.800.CARDIAC (227.3422)
+1.866.484.3268

© 2011 Boston Scientific or its affiliates
All Rights Reserved.

358487-001 EN US 01/11

