



50736123-01

2018-10  
< en >

# endoglide<sup>+</sup>™

## Sterile Bacteriostatic Lubricating Gel

### WARNING

Contents supplied STERILE using a radiation process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

### DEVICE DESCRIPTION

Sterile, transparent, greaseless, water-soluble jelly containing glycerine, PEG, preservatives, methylparaben, propylparaben.

### INTENDED USE/INDICATIONS FOR USE

Intended for medical purposes, to lubricate body orifices to facilitate entry of diagnostic or therapeutic devices when a sterile field is required. Contents sterile in unopened, undamaged package.

### USER INFORMATION

For use by a professional in an endoscopy suite.

### CONTENTS

Provided in either a 56.7 g (2 oz.) tube or 32 g (1.1 oz.) pouch.



### WARNING

May cause eye irritation. If irritation or discomfort occurs, discontinue use and consult a doctor. Keep out of the reach of children. Keep out of eyes and ears.

### CONTRAINDICATIONS

- Sensitivity to ingredients listed in Device Description.

### ADVERSE EVENTS

- Irritation or allergic reaction to ingredients.

### HOW SUPPLIED

Do not use if package is opened or damaged.

Do not use if labeling is incomplete or illegible.

### Handling and Storage

Store in a cool, dry, dark place.

### OPERATIONAL INSTRUCTIONS

**Pouch:** Tear and remove tab and squeeze lubricating gel out of the pouch as needed.

**Tube:** Before first use, remove or puncture seal under cap.

Squeeze lubricating gel out of the tube as needed.



Catalog Number



Consult instructions for use.



Contents



EU Authorized Representative



Legal Manufacturer



Lot



Use By



For single use only. Do not reuse.



Do Not Resterilize



Sterilized using irradiation.



### EU Authorized Representative

Shanghai International Holding Corp. GmbH (Europe)  
Eiffestraße 80, 20537 Hamburg, Germany



### Australian Sponsor Address

Boston Scientific (Australia) Pty Ltd  
PO Box 332  
BOTANY  
NSW 1455  
Australia  
Free Phone 1800 676 133  
Free Fax 1800 836 666



### Argentina Local Contact

Para obtener información de contacto de Boston Scientific Argentina SA, por favor, acceda al link [www.bostonscientific.com/arg](http://www.bostonscientific.com/arg)



### Legal Manufacturer

Jiangsu Province JianErKang  
Medical Dressing Co., Ltd.  
Zhixi Town, 213251 Jintan/Jiangsu, China

### DISTRIBUTED BY:

Boston Scientific Corporation  
300 Boston Scientific Way  
Marlborough MA 01752 USA  
USA Customer Service 888-272-1001



Do not use if package is damaged.



Recyclable Package

© 2018 Boston Scientific Corporation or its affiliates.  
All rights reserved.