AMS Ambicor™

Penile Prosthesis

Instructions for Use

English	AMS Ambicor™ <i>Penile Prosthesis</i> Instructions For Use	3
Français	AMS Ambicor™ <i>Prothèse pénienne</i> Mode d'emploi	21
Deutsch	AMS Ambicor™ <i>Schwellkörperimplantat</i> Gebrauchsanweisung	39
Italiano	AMS Ambicor™ <i>Protesi peniena</i> Istruzioni per l'uso	57
Español	AMS Ambicor™ <i>Prótesis de pene</i> Instrucciones de uso	75
Português	AMS Ambicor™ <i>Prótese Peniana</i> Instruções de Utilização	93

Rx ONLY

AMS™

en English de Deutsch es Español fr Français it Italiano
pt Português

	<p>en Caution: Federal law (U.S.) restricts this device to sale by or on the order of a physician. de Vorsicht: Nach US-amerikanischem Recht darf dieses Gerät nur von einem Arzt oder auf ärztliche Anordnung verkauft werden. es Precaución: las leyes federales de los EE. UU. limitan la venta de este dispositivo a médicos o bajo prescripción facultativa. fr Mise en Garde : Selon la législation fédérale américaine, ce dispositif ne peut être vendu que sur prescription d'un médecin. it Attenzione: le leggi federali degli Stati Uniti limitano la vendita di questo dispositivo ai medici o dietro prescrizione medica. pt Atenção: a lei federal dos Estados Unidos só permite a venda deste dispositivo mediante receita médica.</p>
	<p>en Manufacturer de Hersteller es Fabricante fr Fabricant it Produttore pt Fabricante</p>
	<p>en Authorized Representative in the European Community de Autorisierte Vertretung in der EU es Representante autorizado en la Comunidad Europea fr Représentant autorisé dans la Communauté Européenne it Rappresentante autorizzato della comunità europea pt Representante autorizado na comunidade europeia</p>
	<p>en Date of manufacture de Herstellungsdatum es Fecha de fabricación fr Date de fabrication it Data di fabbricazione pt Data de fabrico</p>
	<p>en Use by YYYY-MM-DD de Verwendbar bis JJJJ-MM-TT es Usar antes de AAAA-MM-DD fr Utiliser avant le AAAA-MM-JJ it Usare entro AAAA-MM-GG pt Utilizar até AAAA-MM-DD</p>
	<p>en Lot number de Chargennummer es N.º de lote fr Numéro de lot it Numero di lotto pt Número de lote</p>
	<p>en Catalog number de Katalognummer es Número de catálogo fr Référence catalogue it Numero di catalogo pt Número de catálogo</p>
	<p>en Sterilized using ethylene oxide de Mit Ethylenoxid sterilisiert es Esterilizado con óxido de etileno fr Stérilisé à l'oxyde d'éthylène it Sterilizzato con ossido di etilene pt Esterilizado por óxido de etileno</p>
	<p>en Sterilized using steam de Mit Dampf sterilisiert es Esterilizar con vapor fr Stérilisé à la vapeur it Sterilizzato a vapore pt Esterilizado por vapor</p>
	<p>en Do not resterilize de Nicht reesterilisieren es No volver a esterilizar fr Ne pas résteriliser it Non risterilizzare pt Não reesterilizar</p>
	<p>en Do not reuse de Nicht wiederverwenden es No reutilizar fr Ne pas réutiliser it Non riutilizzare pt Não reutilizar</p>
	<p>en Consult Instructions for Use de Gebrauchsanweisung beachten es Consultar las instrucciones de uso fr Consulter le mode d'emploi it Consultare le istruzioni per l'uso pt Consultar as instruções de utilização</p>
	<p>en Recyclable packaging de Wiederverwertbare Verpackung es Embalaje reciclable fr Emballage recyclable it Confezione riciclabile pt Embalagem Reciclável</p>
	<p>en MR Conditional de Bedingt MR-tauglich es RM condicional fr Compatible IRM sous conditions it RM con riserva pt Condicional para RM</p>

AMS AMBICOR™

Penile Prosthesis

INSTRUCTIONS FOR USE

NOTE: Refer to the Operating Room Manual for further information on the AMS Ambicor Device and its implantation.

Brief Device Description

The AMS Ambicor Penile Prosthesis is a closed fluid-filled system consisting of two implanted components: a pair of cylinders implanted in the corpora cavernosa, and a pump implanted in the scrotum. Kink-resistant tubing connects the components. The device is delivered prefilled with normal saline and preconnected. The cylinders are inflated as fluid is pumped from the reservoirs, which are located at the proximal end of the cylinders, to the main portion of the cylinder, creating an erection. They are deflated as fluid is transferred back to the reservoirs, making the penis flaccid once again. This device contains solid silicone elastomer. These devices are for men who, after appropriate patient history and diagnostic evaluations as well as discussions with the urologist about other alternative treatment methods, are determined to be suitable candidates for implantation surgery.

Indications for Use

The AMS Ambicor Penile Prosthesis is intended for use in the treatment of chronic, organic, male erectile dysfunction (impotence).

Contraindication

The implantation of this device is contraindicated in patients who have active urogenital infections or active skin infections in the region of surgery.

Warnings

1. **Implantation of the device will make latent natural or spontaneous erections as well as other interventional treatment options, impossible.**
2. Men with diabetes, spinal cord injuries, or open sores may have an increased risk of infection associated with a prosthesis.
3. Failure to evaluate and promptly treat erosion may result in a substantial worsening of the condition leading to infection and loss of tissue.
4. Implantation of a penile prosthesis may result in penile shortening, curvature or scarring.
5. This device contains solid silicone elastomer. The risks and benefits of implanting this device in patients with documented sensitivity to silicone should be carefully considered.
6. Pre-existing abdominal or penile scarring or contracture may make surgical implantation more complicated or impractical.

Precautions

Surgery Related

1. Spontaneous unintended inflation or deflation of the cylinders may occur and may result in unintended partial or full erections.
2. Migration of the device components can occur if the cylinders are improperly sized, or if the pump is not positioned properly.
3. Removal of an implanted prosthesis without timely reimplantation of a new prosthesis may complicate subsequent reimplantation or may make it impossible.
4. Improper measurement technique, positioning or sizing may reduce cylinder life.
5. Unsuccessful outcomes have been reported due to improper surgical technique, anatomical misplacement of components, improper sizing of components, or tubing kinks.

Device Related

1. Some of the materials used in the construction of this device have been shown to cause minor irritation when implanted in animals. Therefore, implantation of this device may cause minor irritation or discomfort in some patients.
2. Do not use product that has damaged or open packaging, as sterility may be compromised.

Patient Related

1. A thorough preoperative consultation should include a discussion between patient and physician of all available treatment options and their risks and benefits.
2. Adequate patient manual dexterity and strength are required for proper device inflation and deflation.
3. Mental or psychological conditions, such as senile dementia, may inhibit the patient's successful operation of the prosthesis.
4. Trauma to the pelvic area, such as impact injuries associated with sports (e.g. bicycle riding), can result in damage of the implanted device and/or surrounding tissues. This damage may result in the malfunction of the device and may necessitate surgical correction, including replacement of the device.
5. The implantation of this device should only be considered in patients whom the physician determines are adequate surgical candidates.

Adverse Events

A clinical trial was conducted to determine the safety and effectiveness of the AMS Ambicor inflatable penile prostheses. A total of 179 patients were implanted with follow-up out to 18 months or longer for 141 patients. The Adverse Device Effects, detailed in the table below, were noted during the duration of this clinical trial for all enrolled patients.

AMS Ambicor Clinical Trial Device Related Adverse Events*

ADE	# Patients (%*)	# Events	Mean Onset time in Days (Range in Days)
Urogenital Pain (Typically Associated with Healing Process)	38 (21.2%)	41	83 (0 – 548)
Patient Dissatisfaction (With Length, Ability to Use and Nonspecific Reasons)	22 (12.3%)	25	269 (0 – 1105)
Mechanical Malfunction	11 (6.1%)	11	394 (0 – 1508)
Penile Sensation Change	11 (6.1%)	11	266 (10 – 573)
Urogenital Hematoma	9 (5.0%)	9	9 (2 – 34)
Urogenital Edema	8 (4.5%)	9	9 (0 – 29)
Infection/Device Infection	6 (3.4%)	8	207 (32 – 606)
Auto-inflation	4 (2.2%)	4	110 (27 – 224)
Penile Curvature	4 (2.2%)	4	323 (0 – 999)

ADE	# Patients (%*)	# Events	Mean Onset time in Days (Range in Days)
Abnormal Ejaculation	3 (1.7%)	3	256 (40 – 487)
Application Site Reaction (Wound Separation, Delay in Cutaneous Closure)	3 (1.7%)	3	26 (0 – 46)
Phimosis/Paraphimosis	3 (1.7%)	3	21 (3 – 38)
Post-Operative Bleeding	3 (1.7%)	3	4 (2 – 6)
Seroma	3 (1.7%)	3	242 (34 – 596)
Urinary retention	3 (1.7%)	3	121 (0 – 361)
Urogenital Ecchymosis	3 (1.7%)	3	123 (6 – 351)
Adhesion	2 (1.1%)	2	42 (34 – 50)
Auto-deflation	2 (1.1%)	2	486 (212 – 760)
Nocturia	2 (1.1%)	2	245 (105 – 385)
Penile Rash	2 (1.1%)	2	272 (43 – 501)

ADE	# Patients (%*)	# Events	Mean Onset time in Days (Range in Days)
Renal Failure	2 (1.1%)	2	476 (396 – 555)
Urethral Laceration	2 (1.1%)	2	1 (0 – 2)
Urinary Tract Infection	2 (1.1%)	2	283 (4 – 561)
Urination Impaired	2 (1.1%)	2	225 (64 – 385)
Urogenital Erythema	2 (1.1%)	2	193 (13 – 373)
Other urogenital	23 (12.8%)	24	97 (0 – 727)
Other non-urogenital	10 (5.6%)	14	176 (0 – 556)

* Percentages based on total number of patients implanted (n=179).

The following “other” urogenital adverse events each occurred in one patient: abnormal sexual function, allergic reaction, BPH, burning urination, crossover from left to right corpora, device migration, dribbling, dysuria, epididymitis, extrusion, incontinence, induration in penile area, lower urinary tract symptoms, malposition, meatal excoriation, micro-hematuria, penile scar tissue, prerenal azotemia, scab on corpora, soft glans, tunica perforation, uresepsis, wound dehiscence and wound infection.

Resolution of Device Related Adverse Events

ADE	# Patients (%)	# Events	Method of Resolution					Continuing
			Surgical	Medication	Other**	No Medical Intervention		
Urogenital Pain	38 (21.2 %)	41	3	9	1	29	8	
Patient Dissatisfaction	22 (12.3 %)	25	3	3	3	16	14	
Mechanical Malfunction	11 (6.1 %)	11	7	0	1	3	1	
Penile Sensation Change	11 (6.1 %)	11	0	2	0	10	4	
Urogenital Hematoma	9 (5.0 %)	9	0	1	3	6	0	
Urogenital Edema	8 (4.5 %)	9	0	2	2	5	0	
Infection/Device Infection	6 (3.4 %)	8	6	1	0	1	1	
Auto-inflation	4 (2.2 %)	4	0	0	0	4	2	

ADE	# Patients (%)	# Events	Method of Resolution					Continuing
			Surgical	Medication	Other**	No Medical Intervention		
Penile Curvature	4 (2.2 %)	4	2	0	0	2	2	
Abnormal Ejaculation	3 (1.7 %)	3	0	1	0	2	1	
Application Site Reaction	3 (1.7 %)	3	0	1	0	2	0	
Phimosis/Paraphimosis	3 (1.7 %)	3	1	2	0	0	0	
Post-operative Bleeding	3 (1.7 %)	3	0	0	0	3	0	
Seroma	3 (1.7 %)	3	0	0	1	2	0	
Urinary Retention	3 (1.7 %)	3	1	1	2	0	1	
Urogenital Ecchymosis	3 (1.7 %)	3	0	0	0	3	0	
Adhesion	2 (1.1 %)	2	1	0	0	1	0	

ADE	# Patients (%)	# Events	Method of Resolution					Continuing
			Surgical	Medication	Other**	No Medical Intervention		
Auto-deflation	2 (1.1 %)	2	0	0	1	1	0	
Nocturia	2 (1.1 %)	2	0	1	0	1	2	
Penile Rash	2 (1.1 %)	2	0	1	1	0	1	
Renal Failure	2 (1.1 %)	2	0	0	1	1	1	
Urethral Laceration	2 (1.1 %)	2	0	0	1	1	0	
Urinary Tract Infection	2 (1.1 %)	2	0	2	0	0	0	
Urination Impaired	2 (1.1 %)	2	0	1	1	0	2	
Urogenital Erythema	2 (1.1 %)	2	0	2	0	0	0	

** Other methods of resolution consist of standard medical treatments.

The following are adverse events associated with this product: Abscess, Bleeding, Erosion, Excessive Fibrous Capsular Growth, Exposure to Biohazardous Material, Genital Changes, Granulomas, Hemorrhage, Immune-Related Connective Tissue Disorders, Improper Size, Inguinal Hernia, Ischemia, Necrosis, Pain (which may be prolonged or severe), Perforation or Injury of Bladder, Corpus Cavernosum, Nerve, Tunica, or Urethra, Prolonged Procedure, Super Sonic Transporter (SST) (as a result of improper sizing), Ulceration, Unretrieved Device Fragment, Vascular Compromise, Ventral Chordee.

The following events occurred during this clinical study but were classified by the Investigators as not related to the device or the implantation surgery: urethral stricture (9 cases), prostate cancer (4 cases), BPH (3 cases), incontinence (3 cases), urinary retention (3 cases), urinary tract infection (3 cases), abnormal prostate exam/PSA (2 cases), abnormal sexual function (2 cases), balanitis (2 cases), hematuria (2 cases), pain (2 cases), renal calculus (2 cases), abnormal ejaculation (1 case), bladder instability (1 case), bladder spasms (1 case), blood in urine (1 case), dysuria (1 case), edema (1 case), epididymal cyst (1 case), eroded artificial sphincter (1 case), erythema (1 case), infection (1 cases), inguinal hernia (1 case), nocturia (1 case), prostate nodule (1 case), prostatitis (1 case), pyuria (1 case), renal colic (1 case), renal insufficiency (1 case), testicular cyst (1 case), testicular hydrocele (1 case), trichomonas (1 case), urethral mucosal atrophy (1 case), urination impaired (1 case), urgency (1 case), and urosepsis (1 case).

There were four patient deaths during the course of the trial. No deaths that occurred during the duration of the clinical study were attributed to the device implantation or use.

A total of 14 patients underwent revision surgeries in the study. Information on device revisions is described in the “Clinical Studies” section.

Clinical Studies

A clinical trial was undertaken to demonstrate that the AMS Ambicor Inflatable Penile Prosthesis provides an erection that is suitable for intercourse and has acceptable rates of surgical revision and of significant clinical events associated with the implantation and use of these devices. This trial was also designed to demonstrate the implantation of these devices does not negatively impact the sexual satisfaction, psychological well being, self-esteem or quality of life of patients who receive these devices. It was a prospective, multi-center cohort trial in which the patients served as their own control.

One hundred and seventy-nine (179) male patients, over 21 years of age, were enrolled and implanted in this study. All patients with diagnosed organic erectile dysfunction were eligible for enrollment, if they did not present with a history of allergy/sensitivity to silicone, preexisting autoimmune or connective tissue diseases or active urogenital infection.

All safety-related data, diagnoses and health status evaluations were captured on detailed case report forms. The Investigators' professional evaluation of the erections provided by the IPPs after implantation and their suitability for intercourse was the primary efficacy endpoint. The number of surgical revisions performed and reported by the Investigators was the primary safety endpoint. Patient self-evaluations on four validated outcome instruments were the secondary efficacy endpoints (concerning quality of life, self-esteem, and sexual satisfaction and functioning).

This clinical trial provided the following results through the 18-month evaluation for the first 141 patients to reach this post-surgical follow-up. Some patients were studied for longer than 18 months.

Physician Assessment of Device Function:

The function of 141 implanted devices was evaluated on or after the 18-month follow-up exam. Of these 141 devices, 140 could be inflated (99.3%). The device that could not be inflated was surgically replaced. All of the 140 devices that were successfully

inflated were determined to provide an erection suitable for intercourse. This assessment of device function does not include the following information regarding device malfunctions: Of the 140 devices found to be functioning properly on or after 18 months, 2 required revision surgery prior to this evaluation to correct a mechanical malfunction. These and other cases of revision surgery are discussed further in the next section.

Surgical Revisions:

The incidence of revisions was evaluated in the 141 patients with follow-up out to at least 18-months, and an additional 6 patients who were explanted prior to 18 months and subsequently discontinued from the study. (A revision is considered any urogenital surgical intervention that is related to the function, placement or site reaction to the implanted device.) Of these 147 patients, 14 (9.5%) experienced a total of 15 revision surgeries, and 133 (90.5%) have not been revised.

The average time to the first revision surgery was 12 months (ranging from 0 to 50 months). Of the fifteen (15) revision surgeries, there were six (6) revisions due to "Mechanical Malfunction"; five (5) due to "Infection"; one (1) due to "Migration/Extrusion" (pump migration and cylinder aneurysm/bulging/extrusion); one (1) due to "Malposition"; one (1) due to "Phimosis"; and one (1) due to "Patient Dissatisfaction".

Patient Evaluation of Quality of Life, Self-esteem, Psychological and Sexual Well-Being

No decline in overall patient quality of life, self-esteem, psychological well being, sexual functioning or sexual satisfaction was demonstrated during this clinical trial. The psychological well being of these patients, as judged by the Brief Symptom Inventory (BSI) was determined to be equivalent to the pre-implant state. The self-esteem of the patients, as evaluated by the Rosenberg Self-Esteem Scale, was also equivalent to the pre-implant state. The sexual functioning and sexual satisfaction did improve from the pre-implant state, as measured by the Sexual History Form (SHF).

Patient Counseling Information

Patients should be counseled in order to have a realistic expectation of the physical, psychological and functional outcome of the implantation. The risks, benefits and potential adverse events of all available treatment options should be discussed with the patient and considered by the physician and patient when choosing a treatment option.

An appropriate patient history, including history of personality disorders, and diagnostic work-up should be a part of the patient decision making process.

Some patients may become dissatisfied by the presence of the prosthetic device in their body. This issue should be discussed with the patient prior to the surgery. Patient dissatisfaction may lead to device removal.

Implantation of a penile prosthesis may result in penile shortening, curvature or scarring. The prosthetic erection may differ from the patient's original, natural erection in that it may be shorter, less firm, have less girth, and reduced sensations. Realistic cosmetic expectations should be communicated to the patient and should include the potential for skin scarring, scrotal deformity, pump bulge in the scrotum, lack of concealability and other possible adverse events. Patients should also be aware that penile prostheses are not considered to be lifetime implants.

Improper implantation of a penile prosthesis may not provide rigidity to the glans, which may result in a floppy glans and may result in a lack of rigidity of the corpus spongiosum. Penile flaccidity may be less than it was before implantation.

Patients who undergo revision surgery may notice a change in the character of their erection compared to their previous implant, which may include differences in sensation, length, girth, rigidity, and/or flaccidity. It is also important that the physician discusses with the patient the possibility of an allergic reaction to the materials in the device (See Silicone Information).

Silicone Information

This device is composed of a number of materials, including solid silicone elastomers and a fluorosilicone lubricant. Silicone gel is not a component in the materials of this device.

Silicone elastomers have been commonly used in a variety of biomedical devices for over 40 years and are used as a biocompatibility reference against which new materials are tested. Silicone fluids have an extensive history of use in medical devices.

Scientific literature has included reports of adverse events and other observations in patients with implantable silicone devices. As reported, these events/observations indicate “allergic-like” symptoms and in other cases a symptom complex associated with immunological disorders. No casual relationship has been established between these events and silicone elastomer or fluorosilicone lubricant.

There are reports of malignant tumor formation in laboratory animals only, not humans, associated with implants of relatively large size. Many different materials are associated with this effect in animals, silicone elastomers among them. No such effect has been described in humans.

Extensive testing has been conducted on all materials that comprise the prostheses in the AMS Ambicor. This testing has indicated no toxicological response attributable to the materials. However, some of the materials caused minor irritation when implanted in animals.

Silicone elastomer particulate shedding and particulate migrations to regional lymph nodes have been reported in the literature on penile implants. There are no known clinical sequelae to this phenomenon.

Magnetic Resonance Imaging (MRI) Important Safety Information

Non-clinical testing has demonstrated the

Ambicor product line is MR Conditional. The device can be scanned safely under the following conditions:

Static Magnetic Field	1.5 Tesla ^a	3.0 Tesla ^b
Spatial Gradient Field	450 Gauss/cm or less	720 Gauss/cm or less
Maximum whole body averaged Specific Absorption Rate (SAR)	1.5 W/kg for 15 minutes of scanning as assessed by calorimetry	2.9 W/kg for 15 minutes of scanning as assessed by calorimetry

MRI-Related Heating

Non-clinical testing has demonstrated the Ambicor product line produced the temperature rises during MRI performed for 15 minutes of scanning in the respective MR systems which would not pose a hazard to the human subject.

Static Magnetic Field	1.5 Tesla ^a	3.0 Tesla ^b
Highest Temperature Change	≤ + 0.4° C	≤ + 1.9 ° C
(a) 1.5T - 64 MHz MR System (General Electric Healthcare, Milwaukee, WI) (b) 3.0T MR Excite, General Electric Healthcare, software version 14X.M5		

Artifact Information

Non-clinical testing has demonstrated that the Ambicor product line may compromise the MR image quality if the area of interest is relatively close to the position of the implant. The maximum image artifact produced by a MR gradient echo pulse sequence was a “moderate” localized signal void in size and shape of the implant. Optimization of MR imaging parameters to compensate for the presence of the device may be necessary.

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	6,244 mm ²	1,589 mm ²	10,295 mm ²	2,779 mm ²
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

How Supplied and Storage

WARNING: Contents supplied STERILE. Do not use if sterile barrier is damaged. If damage is found, call your AMS representative.

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

Store device in a clean, dry, dark area at room temperature.

Inventory Returns and Product Replacement Information

In the United States

Before returning any components, whether explanted or unused (sterile or nonsterile), customers must fill out the Return Goods Form located on the last page of the Patient Information Form. Follow all of the instructions on the form carefully, and be sure that the components have been thoroughly cleaned before returning them to AMS.

In all cases, obtaining credit or percentage of credit for a returned component is subject to approval under the terms of the AMS Return Goods Policy and the AMS Product Replacement Policy. For complete information regarding these policies, contact the AMS Customer Service Department.

Outside the United States

Customers outside of the United States should contact their local AMS Representative prior to returning any product.

This document is written for professional medical audiences. Contact American Medical Systems for lay publications.

American Medical Systems periodically updates product literature. If you have questions about the currency of this information, contact American Medical Systems.

This page intentionally left blank.

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All rights reserved. All trademarks are property of the
American Medical Systems, Inc.
Minnetonka, MN 55343 U.S.A.
American Medical Systems Europe B.V.
©2017 Boston Scientific Corporation or its affiliates.