


- Three-year Results from the EVOLVE II Randomized Trial presented at ACC
- Synergy™ BP Stent Launches 48mm Length
- Boston Scientific announces agreement to acquire Symetis



- EuroPCR LIVE case: Rotablation of LAD
- InCathlab Live Case: Antegrade Recanalization

Learning with Clinical Studies:

- Guidezilla™: back-up support for calcified, complex and tortuous anatomy



- Early Discontinuation of Dual-Antiplatelet Therapy with Synergy™ Stents in High-Risk Patients Undergoing Complex PCI
- Rotational atherectomy: You will never regret using it!
- Synergy™ Stent: ST Rates from the Academic Research Consortium



- Join the Complex PCI community
- Download New App: Stroke-Bleed Risks Calculator!




Featured
Content

Congresses and Events

Three-year Results from the EVOLVE II Randomized Trial presented at ACC

Late Clinical Outcomes with Bioresorbable compared to Permanent Polymer Everolimus-Eluting Stents


Dean J. Kereiakes, Ian T. Meredith, Stephan Windecker, Josep Rodes-Cabau, Tommy Lee, Arthur Reitman, Andrejs Erglis, Mark Dorogy, Barry Bertolet, Louis Cannon, Juhani Airaksinen, Craig Siegel, Thomas Christen, Dominic J. Allocco, Keith D. Dawkins



Three-year follow-up from EVOLVE II trial presented by Dr. Dean Kereiakes supports longer-term safety and efficacy of the novel abluminal bioabsorbable polymer Synergy everolimus-eluting stent in a broad range of patients undergoing PCI.

Key data points:

- Definite/Probable ST was 0.8% in the Promus Element™ Plus Stent arm and 0.5% in the Synergy™ BP Stent arm
- Synergy™ BP Stent had 0.2% Definite/Probable ST after 24 hours to 3-years, while Promus Element™ Plus Stent had 0.7% (landmark analysis)
- Synergy™ BP Stent shows a 0.06% per year ST rate beyond 1 year (a ST rate 10x lower than contemporary permanent polymer stents)
- TLF was similar in both arms: 11% for Synergy™ BP Stent and 10% for Promus Element™ Plus Stent




[WATCH VIDEO »](#)




Congresses and Events

Excellent Safety Performance with Synergy™ BP Stent:



Numerically Lower ARC Definite/Probable ST with the Synergy™ BP-DES compared to the Promus™ Element™ PP-DES:





Featured
Content

Congresses and Events

Synergy™ BP Stent Launches 48mm Length

Long lesions represent a complex lesion subset that is associated with increased safety related events, reduced efficacy and increased procedure time. Synergy™ BP Stent was designed for quality healing and safety over time and with the addition of a 48mm length it has a complete size matrix to address these complex lesions. Synergy™ BP Stent, with its abluminal, bioabsorbable polymer, shows more complete stent coverage compared to permanent polymer stents.

Heal long lesions with confidence



[WATCH VIDEO »](#)

euro
PCR Symposium webcast
May 16th to 19th Paris

REGISTER NOW >

The Synergy Stent™: Complete revascularisation solutions in complex PCI
Live Streaming May 16th - 17h15-18h45 With an educational grant from Boston Scientific

[LEARN MORE »](#)



Featured
Content

Congresses and Events

Boston Scientific announces Agreement to acquire Symetis

Boston Scientific announced a definitive agreement to acquire Symetis SA, a privately-held Swiss structural heart company focused on minimally-invasive, for \$435 million in up-front cash.

The Symetis portfolio includes the ACURATE TA™ and ACURATE neo/TF valve* systems for use in the treatment of high-risk patients suffering from severe and symptomatic aortic valve stenosis, which are sold in Europe and in other geographies outside of the United States. Symetis is also developing the ACURATE neo/AS** next generation valve system, currently in a clinical trial intended to serve as the basis for a future CE mark application.



Lotus™ Valve System, Boston Scientific

This agreement follows the recent acquisition by Boston Scientific of certain Neovasc, Inc. manufacturing assets, and demonstrates the company's continued investment in structural heart through intellectual property, research and development, and manufacturing capabilities.

The acquisition is projected to close during the second quarter of 2017, subject to customary closing conditions.



Ian Meredith, M.D.,
Executive Vice President
and Global Chief Medical
Officer of Boston Scientific.

*„The steps we are taking reflect our commitment to being a leader in TAVI and structural heart technologies now and over the long-term, as we broaden our portfolio and pipeline to address the needs of our global health care providers and their patients. The ACURATE family of valve products is strongly complementary to our cornerstone Lotus™ valve*** platform, and this compelling combination of technologies will allow us to provide interventional cardiologists and cardiac surgeons with multiple TAVI offerings for varying patient pathologies and anatomy.”*

euro
PCR Symposium webcast
May 16th to 19th Paris

REGISTER NOW >

The Lotus™ Valve: Expanding treatment options in real-world practice
Live Streaming May 17th - 13h35-14h35 With an educational grant from Boston Scientific

*The ACURATE TA™ and ACURATE neo/TF valve systems are not available for use or sale in the US.

**The ACURATE neo/AS valve system is not available for use or sale.

***The Lotus™ valve is currently not available for use or sale

LEARN MORE >>



Educational Corner

EuroPCR LIVE case: Rotablation of LAD

Invasive physiological assessment can be used to guide revascularization of residual disease after PCI. From St Thomas' Hospital (London) a live demonstration of an optimal complete revascularization after primary PCI in patients with multivessel disease, to learn more about rotablation of LAD.



Simon Redwood, MD
St Thomas' Hospital,
London, UK



Bernard Prendergast, MD
St Thomas' Hospital,
London, UK

Educational objectives:

- Understand the evidence for full revascularization after primary PCI in patients with multivessel disease
- Know how to assess residual ischaemia using invasive techniques (FFR)
- Appreciate the thresholds for considering rotablation in calcific coronary artery disease

Patient History:

- 80-year old male
- Minimal strenuous activity, carries out light work and an independent lifestyle
- Hypertension, smoker, hypercholesterolemia, inferior STEMI
- 24/03: Inferior STEMI with PPCI to RCA 3 x drug-eluting stents
- 23/04: Admitted to local hospital for reassessment of residual left coronary artery disease; LAD more heavily calcified than originally appreciated; referred for physiological assessment +/- rotablation

Procedure:

1. Check right coronary angiography
2. Update left coronary angiography with physiological assessment of the LAD & LCx
3. PCI of functionally significant lesion(s)
4. Rotablation available if needed for LAD disease



Featured
Content

Educational Corner

Rotablation of LAD



[WATCH VIDEO »](#)

[LEARN MORE »](#)



Featured
Content

Educational Corner

InCathlab Live Case: Right Coronary Artery CTO. Antegrade Recanalization

Educational objectives:

An innovative operation, the Antegrade Recanalization with Hybrid Approach, conducted with a live demonstration from London Chest Hospital.

Patient history:

- 69 years old man, taxi driver
- PCI to RCA 2000
- 8/12 worsening angina
- Failed relicensing ETT
- Coronary angiogram June 2013
- RCA CTO (JCTO score 1)



Elliot Smith, MD
Chest Hospital, London



Simon Walsh, MD
Belfast Trust, Belfast



[WATCH VIDEO »](#)

[LEARN MORE »](#)




Featured
Content

Educational Corner

Learning with Case Studies: Guidezilla™: back-up support for calcified, complex and tortuous anatomy

This case illustrates the need of Guidezilla™ Guide Extension Catheter as a transradial back-up support for calcified, complex and tortuous anatomy.



Brian D. Arcement MD
Gulf Coast Hospital
Fort Myers, Florida, USA

Key Learnings:

- Complex, calcified and tortuous distal lesions and transradial cases are perfect cases to use the Guidezilla™ to increase the probability of success and reduce the time spent delivering devices to lesions
- Balloon anchor technique makes delivery of Guidezilla™ through complex, transradial anatomy easy.

MASTER THE COMPLEX

Optimizing revascularization through innovation, training, and education.

CASE STUDY
GUIDEZILLA™
TRANSRADIAL BACK-UP SUPPORT
FOR CALCIFIED, COMPLEX
AND TORTUOUS ANATOMY

Author: Dr. Brian D. Arcement MD
Gulf Coast Hospital
Fort Myers, Florida, USA

This case illustrates the need of GUIDEZILLA™ Guide Extension Catheter as a transradial back-up support for calcified, complex and tortuous anatomy.

Patient History

- 65-year-old male
- Prior stroke, severe depression
- Severe prior angina, CABG in April 2010

Presenting Problem (PC)

- Presented to PC with severe angina, inferior to anterior angiology
- Urgent PCI - could not cross the lesion with the catheter after using guide wires, a 13 mm dilation to pre-dilate lesion, and several different guiders
- Patient was on aspirin and clopidogrel

Diagnosis

- Patient presented with severe anginal symptoms
- Collateral flow present, no significant stenosis
- Anterior to posterior angiology, severe PC
- Severe angina, severe depression and severe tortuous anatomy (distal to proximal)
- Patient had PCI in RFA for the PC which resulted in a 13 mm dilation to pre-dilate lesion, and several different guiders
- Severe stenosis of the proximal to the PC
- Severe stenosis of the PC, distal to the stenosis to the PC
- Severe stenosis of the PC, distal to the stenosis to the PC
- Severe stenosis of the PC, distal to the stenosis to the PC
- Severe stenosis of the PC, distal to the stenosis to the PC

KEY LEARNINGS

This case highlights:

- Complex, calcified and tortuous distal lesions and transradial cases are perfect cases to use the GUIDEZILLA™ to increase the probability of success and reduce the time spent delivering devices to lesions
- Balloon anchor technique makes delivery of GUIDEZILLA™ through complex, transradial anatomy easy

[Download Case Study »](#)

[LEARN MORE »](#)




Clinical Studies

Early Discontinuation of Dual-Antiplatelet Therapy with Synergy™ Stents in High-Risk Patients Undergoing Complex PCI

Rebecca L. Noad, MB, PhD; Colm G. Hanratty, MD; Simon J. Walsh, MD. The Journal of Invasive Cardiology, Epub 2016 December 15

As more elderly and co-morbid patients require percutaneous revascularization, 1 year of dual-antiplatelet therapy (DAPT) becomes concerning. Synergy™ stents allow for early cessation of DAPT. This study assessed those in our unit who underwent percutaneous coronary intervention (PCI) with a Synergy™ stent to examine a minimum of 6 months of clinical outcomes after early discontinuation of DAPT.




Conclusion:

In this small cohort, the use of Synergy stents allows for early discontinuation of DAPT, reducing the risk of bleeding complications and facilitating non-cardiac procedures, without an increase in the incidence of ST. The results for TLF and clinical outcomes are excellent for a group of patients with significant co-morbidities and complex coronary lesions.



[Download Clinical Article »](#)

[LEARN MORE »](#)



Clinical Articles

Rotational atherectomy: You will never regret using it!

EuroIntervention 2016;12:1441-1442

Teresa Strisciuglio¹, MD; Emanuele Barbato^{1,2*}, MD, PhD, 1. Department of Advanced Biomedical Sciences, University of Naples Federico II, Naples, Italy; 2. Cardiovascular Research Center Aalst, Aalst, Belgium

Dr. Teresa Strisciuglio (Naples, Italy) and Dr. Emanuele Barbato (Aalst, Belgium) illustrate the advantages deriving from the contemporary rotational atherectomy technique.

Rotational Atherectomy has recently gone through a surge of interest. The contemporary RA technique in fact aims to smoothen the lumen and disconnect the calcified coronary ring, thus leaving the way clear for further balloon dilatation and stent implantation. This has implied smaller burr-to-artery ratio, lower rotational speed, and burr manipulation, aiming to reduce the friction and temperature increase within the ablated coronary segment.

In addition, the procedural outcome of RA has significantly improved thanks to the evolution from the original extensive “plaque debulking” to the current “plaque modification” technique (Table 1).



	Contemporary
Arterial access	Radial (6-7.5 Fr) or femoral (6-8 Fr), depending upon burr size requirement and operator experience
Guiding catheter	Single curve with strong support. Operator preference but stable catheter position required
Guidewire	Rotawire placement not always straightforward. Use of regular wire placement, with exchange using microcatheter placement often required
Burr size	Plaque modification with small burrs (1.25 mm to 1.5 mm) as initial strategy is default position. A step-up approach is encouraged to limit debris size and complications
Ablation speed	Plaque modification usually achieved at low speeds (135,000 to 180,000 rpm) to reduce risk of complications
Temporary pacemaker	Smaller burrs at lower speeds have led to lower incidence of transient heart block. Many operators use atropine to treat, avoiding any complications of temporary pacemaker placement
Rotablation flush	Rotablation cocktail with verapamil, nitrates and heparin in saline recommended

Table 1

EDITORIAL

Rotational atherectomy: you will never regret using it but you often regret not having used it!

Teresa Strisciuglio¹, MD, Emanuele Barbato^{1,2*}, MD, PhD
1. Department of Advanced Biomedical Sciences, University of Naples Federico II, Naples, Italy;
2. Cardiovascular Research Center Aalst, Aalst, Belgium

Rotational atherectomy (RA) was introduced nearly 30 years ago in order to remove coronary atheroma by “debulking” as opposed to conventional plain old balloon angioplasty (POBA) that was meant to enlarge the lumen by displacing the plaque¹. The initial enthusiasm around RA was quickly replaced by skepticism following the short-term and mid-term results which were far below expectations². Percutaneous coronary intervention (PCI) with RA, in fact, was not infrequently characterized by procedural complications (such as large dissections, no-reflow, etc.) and a high restenosis rate (up to 50% in some instances)³. The suboptimal procedural results observed in these early experiences were mostly related to the aggressive debulking RA strategy and to the less efficient antiproliferative therapy regimens used. In addition, many of these RA procedures were stand-alone, followed by POBA or by bare metal stents. More recently, RA has gone through a second stage of interest, mainly for the following three reasons: 1) PCI is increasingly being performed in heavily calcified coronary arteries as a consequence of the aging of the patients being referred to the cathlab; 2) excellent results achieved with current-generation drug-eluting stents (DES) have encouraged the performance of PCI in more challenging settings; 3) novel technologies such as bioresorbable vascular scaffolds require more extensive lesion preparation. In addition, the procedural outcome of RA has significantly improved thanks to the evolution from the original extensive “plaque debulking” to the current “plaque modification” technique (Table 1). The contemporary RA technique in fact aims to smoothen the lumen and disconnect the calcified coronary ring, thus leaving the way clear for further balloon dilatation (often with non-compliant or cutting balloons) and stent implantation. This has implied smaller burr-to-artery ratio, lower rotational speed, and burr manipulation, aiming to reduce the friction and temperature increase within the ablated coronary segment (pecking motion technique). In addition to these technical modifications, novel and more potent oral antiplatelet agents have allowed this complex procedure to be performed more safely (less no-reflow, fewer bleedings). In this issue of EuroIntervention, two reports from the ROTATE registry provide procedural data and clinical endpoints of the RA technique adopted in patients treated from 2002 to 2013^{4,5}. As one

Table 1. Contemporary rotational atherectomy technique (modified from Barbato et al⁶).

	Contemporary
Arterial access	Radial (6-7.5 Fr) or femoral (6-8 Fr), depending upon burr size requirement and operator experience
Guiding catheter	Single curve with strong support. Operator preference but stable catheter position required
Guidewire	Rotawire placement not always straightforward. Use of regular wire placement, with exchange using microcatheter placement often required
Burr size	Plaque modification with small burrs (1.25 mm to 1.5 mm) as initial strategy is default position. A step-up approach is encouraged to limit debris size and complications
Ablation speed	Plaque modification usually achieved at low speeds (135,000 to 180,000 rpm) to reduce risk of complications
Temporary pacemaker	Smaller burrs at lower speeds have led to lower incidence of transient heart block. Many operators use atropine to treat, avoiding any complications of temporary pacemaker placement
Rotablation flush	Rotablation cocktail with verapamil, nitrates and heparin in saline recommended

would expect, patients treated with RA were elderly and at high risk considering that one third had diabetes, one quarter had renal failure (15% treated with dialysis), 17% presented with multivessel disease, and 20% presented with acute coronary syndromes (ACS). The RA technique applied aimed at plaque modification (1.5 mm was the most frequent burr size used), and in 80% of the cases one burr was sufficient. Interestingly, the arterial access used was femoral in 70% of the cases. MACE occurred in-hospital in 8.6% of the patients and was mostly driven by periprocedural myocardial infarction. All patients were treated with classic dual antiplatelet therapy (aspirin plus ticagrelor or clopidogrel). We have no data on the use of bivalirudin in these patients, something which was previously associated with lower periprocedural myocardial infarction in patients undergoing RA⁷. Whether these results could be improved by the use of more potent P2Y₁₂ inhibitors still remains to be addressed. Unfortunately, the authors did not provide a temporal trend in the technique applied that would have enabled a better understanding

*Corresponding author: Cardiovascular Research Center Aalst, Moorslaan 164, B-9300 Aalst, Belgium.
E-mail: emanuele.barbato@ghc-aalst.be

© European Digital & Publishing 2016. All rights reserved.
The opinions expressed in this article are not necessarily those of the Editors of EuroIntervention or of the European Association of Percutaneous Cardiovascular Interventions.

Download Clinical Article »

LEARN MORE »



Featured
Content

Clinical Studies

Synergy™ Stent: ST Rates from the Academic Research Consortium

According to the Academic Research Consortium, Synergy™ has demonstrated very low stent thrombosis rates across multiple clinical trials.

Data in the table below show that Synergy™ Stent, with its abluminal and bioabsorbable polymer has reported consistently low sub-acute, late & very late ST in 18,000 patients across 9 studies.

	SWEET Registry	Fribourg Experience	Belfast Experience	EVOLVE II Trial	EVOLVE Trial	EVOLVE China	EVOLVE II QCA Study	SCAAR Registry	BIO-RESORT Trial
N:	820	671	185	846	94	205	100	14,979	1172
Acute	1.5%	0.3%	0%	0.2%	0%	0%	0%	0.08%*	0.1%
Sub-acute	0.1%	0.3%	0%	0%	0%	0%	0%	0.02%*	0.1%
Late	0.1%	0.1%	0%	0%	0%	0%	0%	0.2%*	0.2%
Very Late			0%	0.1%	0%			0.1%*	


SWEET Registry: Cook TCT 2015., Fribourg Experience: Arroyo CRT 2016. Belfast Experience: Noad TCT 2015., EVOLVE II: Kereiakes, et al. Circ Cardiovasc Interv. 2015;8:e002372.DOI:10.1161/CIRCINTERVENTIONS.114.002372., EVOLVE FHU: Meredith et al. J Am Coll Cardiol. 2012; 59 (15):1362., EVOLVE II QCA: Meredith ACC 2015., SCAAR Registry: James TCT 2016 BIO-RESORT presented by Clemens von Birgelen, MD PhD, TCT 2016.

The SWEET Registry, conducted at three centers in Switzerland, demonstrated a 1.5% ST rate in the first 24 hours. This rate is higher than we have seen in other clinical trials that have evaluated the performance of Synergy. Since SWEET is an unselected registry, it included by default most patients who were excluded from other clinical trials.

Acute: ≤ 1 day
Subacute: 2 – 30 days
Late: 30 days – 1 year
Very Late: Beyond 1 year


*Cumulative adjusted ARC def ST estimated from Kaplan Meier Curve

LEARN MORE »



Featured
Content

Featured Content



Join the Complex PCI Community

Stay up to date with the latest advances for optimizing revascularization.

Join the Complex PCI Community and download the EuroIntervention Clinical Article: *"Rationale and design of the SYNTAX II trial evaluating the short to long-term outcomes of state-of-the-art percutaneous coronary revascularisation in patients with de novo three-vessel disease"* by Prof. Escaned et al.

Register Now »

LEARN MORE »



NEW Stroke-Bleed Risks Calculator App

This App will help you easily balance the risk of stroke against the risk of bleeding prior to decide the best treatment option for your AF patients and consider the LAAC therapy.

Watch Video »

download new App



LEARN MORE »



**Boston
Scientific**
Advancing science for life™

www.bostonscientific.eu

* Lund experience with the™ Valve in an all-comers population, Matthias Götzberg, EuroPCR 2016

IC-465107-AA Printed in Germany by medicalvision.

© 2017 Boston Scientific Corporation
or its affiliates. All rights reserved.