


**AGENT™**  
Drug-Coated Balloon



**Leave The Right  
AGENT Behind**

AGENT makes Reliable Drug Transfer More Deliverable.

**Compliance Chart**

Pressure atm (kPa)	Balloon Size (mm)						
	2.00	2.25	2.50	2.75	3.00	3.50	4.00
3.0 (304)	1.86	2.06	2.28	2.53	2.76	3.19	3.66
4.0 (405)	1.93	2.14	2.37	2.61	2.85	3.30	3.80
5.0 (507)	1.99	2.20	2.44	2.68	2.93	3.39	3.88
6.0 (608)	2.03	2.26	2.50	2.75	3.00	3.46	3.96
7.0 (709)	2.07	2.31	2.55	2.81	3.06	3.52	4.04
8.0 (811)	2.10	2.34	2.59	2.85	3.11	3.57	4.14
9.0 (912)	2.13	2.38	2.62	2.88	3.15	3.61	4.00
10.0 (1013)	2.15	2.40	2.65	2.91	3.18	3.64	4.18
11.0 (1115)	2.18	2.42	2.67	2.94	3.21	3.68	4.22
12.0 (1216)	2.19	2.44	2.69	2.96	3.23	3.72	4.25
13.0 (1317)	2.21	2.46	2.72	2.99	3.26	-	-
14.0 (1419)	2.23	2.48	2.74	3.02	3.28	-	-


Nominal Pressure  Rated Burst Pressure. Do Not Exceed.

\*Drug coating designed to deliver therapeutic drug dose in association with balloon angioplasty treatment.

† Bench testing as compared to Apex™ and Maverick™ balloons. Data on file at Boston Scientific.

AGENT Drug-Coated Balloon vis distributed by Boston Scientific Corporation. The legal manufacturer is Hemoteq AG, Adenauerstrasse 15, 52146 Würselen Germany.

# AGENT™ Drug-Coated Balloon



## Product Details

Indications for Use	The AGENT DCB Balloon Catheter is indicated for Percutaneous Transluminal Coronary Angioplasty (PTCA) in the coronary arteries to treat in-stent restenosis (ISR) and de novo / small vessel disease.
Drug Coating	The TransPax drug coating is a proprietary formulation of Paclitaxel and a citrate ester excipient (acetyl tributyl citrate - ATBC). The AGENT balloon catheter is designed to prevent hyperplasia of smooth muscle cells at the treated lesion site.
Drug Dose	2.0 µg per mm <sup>2</sup> of the balloon surface
Recommended Inflation Time	At least 30 seconds
Antiplatelet Regimen	Recommended for 3 months after treatment
Available Balloon Lengths	In mm: 12, 15, 20, 30
Available Balloon Diameters	In mm: 2.00, 2.25, 2.50, 2.75, 3.00, 3.50, 4.00
Lesion Entry Profile	0.43 mm or 0.017"
Delivery System Effective Length	144 cm
Delivery System Port	Designed for 0.36 mm or 0.014" guidewire
Balloon Material	OptiLEAP™
Shelf Life	24 months for Use
Sterilization	Ethylene oxide (EO)
Single/Multiple Use	Single Use
Guide Catheter Inner Diameter	1.42 mm or 0.056"
Guide Catheter Compatibility	5 F or 1.67 mm
Catheter Shaft Outer Diameter	1.67 mm or 0.0656"
Marker Band Material	Platinum-Iridium
Marker Band Placement	The outside edges of the marker bands indicate the balloon's shoulders and coated region

## Ordering Information

ø (mm)	Balloon Length (mm)			
	12mm	15mm	20mm	30mm
2.00	H749 3922220121 0	H749 3922220151 0	H749 3922220201 0	H749 3922220301 0
2.25	H749 3922222121 0	H749 3922222151 0	H749 3922222201 0	H749 3922222301 0
2.50	H749 3922225121 0	H749 3922225151 0	H749 3922225201 0	H749 3922225301 0
2.75	H749 3922227121 0	H749 3922227151 0	H749 3922227201 0	H749 3922227301 0
3.00	H749 3922230121 0	H749 3922230151 0	H749 3922230201 0	H749 3922230301 0
3.50	H749 3922235121 0	H749 3922235151 0	H749 3922235201 0	H749 3922235301 0
4.00	H749 3922240121 0	H749 3922240151 0	H749 3922240201 0	H749 3922240301 0

**Boston  
Scientific**  
Advancing science for life™

[bostonscientific-international.com](http://bostonscientific-international.com)

All trademarks are the property of their respective owners.

CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings, and instructions for use can be found in the product labelling supplied with each device or at [www.IFU-BSCI.com](http://www.IFU-BSCI.com). Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. This material not intended for use in France.

2023 Copyright © Boston Scientific Corporation or its affiliates. All rights reserved.

Copyright © 2018 by  
Boston Scientific Corporation  
or its affiliates. All rights reserved.

IC-236507-AD