Expect™ Slimline (SL)

Endoscopic Ultrasound Aspiration Needle

Expect™ Slimline (SL)

F L E X I B L E

Endoscopic Ultrasound Aspiration Needle

Directions for Use	2
Instrucciones de uso	8
Mode d'emploi	15
Gebrauchsanweisung	22
Istruzioni per l'uso	29
Gebruiksaanwijzing	36
Instruções de Utilização	43

TABLE OF CONTENTS

WARNING3

DEVICE DESCRIPTION3

NEEDLE SPECIFICATIONS3

INTENDED USE/INDICATIONS FOR USE3

CONTRAINDICATIONS.....3

WARNINGS.....3

PRECAUTIONS.....3

ADVERSE EVENTS4

HOW SUPPLIED.....4

OPERATIONAL INSTRUCTIONS.....4

WARRANTY7

Expect™ Slimline (SL)

Endoscopic Ultrasound Aspiration Needle

Expect™ Slimline (SL)

FLEXIBLE

Endoscopic Ultrasound Aspiration Needle

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death.

Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

The needle tip is sharp. Take precautions to ensure that the needle is handled properly. The needle should be treated as infectious material and could create an infection risk.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

DEVICE DESCRIPTION

The Expect Slimline (SL) Endoscopic Ultrasound Aspiration Needle is an endoscopic ultrasound aspiration needle that can be coupled to the biopsy channel of a Curvilinear Array (CLA) Echoendoscope with a standard luer connection and delivered into the digestive tract. The device sheath length can be adjusted to accommodate different model echoendoscopes.

The needle is used to acquire aspiration samples from lesions and provide therapy within and adjacent to the digestive system's major lumens that can be identified and targeted using the echoendoscope. Both the sheath and needle length can be adjusted based on distance to the target lesion. The sheath and needle length adjustments are set and locked by the physician by using the locking knob mechanisms on the handle of the device.

An aspiration sample is obtained by penetrating the lesion with the needle while applying suction and manipulating the needle in a back and forth motion to aspirate cytology samples into the needle. The sample can be prepared per normal institutional protocol.

The Expect Slimline (SL) Needle has echogenic (visible under ultrasound) features at the distal end to facilitate real time visualization of the device under ultrasound.

NEEDLE SPECIFICATIONS

Working Length: 137.5 cm to 141.5 cm adjustable

Needle Length: 0 cm to 8 cm adjustable

Diameter: Three different needle sizes - 19 gauge, 22 gauge, and 25 gauge

M00555500	19 gauge OD: 1.10 mm	Minimum 2.8 mm Working Channel
M00555510	22 gauge OD: 0.72 mm	Minimum 2.4 mm Working Channel
M00555520	25 gauge OD: 0.52 mm	Minimum 2.4 mm Working Channel
M00555530	19 gauge flexible OD: 1.14 mm	Minimum 2.8 mm Working Channel

INTENDED USE/INDICATIONS FOR USE

The Expect Slimline (SL) Needle is designed to sample targeted submucosal and extramural gastrointestinal lesions through the accessory channel of a curvilinear echoendoscope. It can also be used for delivery of injectable materials (fluids) or fiducials into tissue or for passage of accessory devices.

CONTRAINDICATIONS

Contraindications for this device are those specific to the primary endoscopic procedure to be performed in gaining access to the desired site. Relative contraindications to submucosal and extramural aspiration include, but are not limited to: coagulopathy.

WARNINGS

The Expect Slimline (SL) Endoscopic Ultrasound Aspiration Needle should only be used to sample tissue where possible hemorrhage will not present a danger to patients. The device should be used with caution and only after careful consideration in patients with elevated bleeding times or coagulopathies.

Not for use in the heart or vascular system.

PRECAUTIONS

Read the entire directions for use before using the Expect Slimline (SL) Needle as well as the directions for use supplied by the manufacturers of injectable materials (fluids), fiducials, and accessory devices, if used. The Expect Slimline (SL) Needle should only be used by or under the supervision of physicians trained in Endoscopic Ultrasound procedures and Fine Needle Aspiration. A thorough understanding of the technical principles, clinical applications and risks associated with these procedures are necessary before using this device.

The packaging and the device should be inspected prior to use. Do not use the device if the product or packaging is damaged.

When targeting multiple sites, replace device for each site. Use of the device for multiple sites has not been established.

Do not use the same device to perform multiple indications. Use of the device for multiple indications has not been established.

ADVERSE EVENTS

Complications associated with the use of the Expect™ Slimline (SL) Needle may include:

- Bleeding
- Perforation
- Pancreatitis
- Infection
- Peritonitis
- Inflammation
- Aspiration
- Fever
- Allergic Reaction
- Hypotension
- Respiratory Depression or Arrest
- Cardiac Arrhythmia or Arrest
- Tumor Seeding
- Pain
- Fiducial Migration

HOW SUPPLIED

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if package is opened or damaged. Do not use if labeling is incomplete or illegible. If damage is found, call your Boston Scientific representative. Package contents include:

- One (1) Expect Slimline (SL) Needle
- One (1) Vacuum Syringe
- One (1) One-Way Stopcock

Storage

Store in a cool, dry, dark place.

OPERATIONAL INSTRUCTIONS

Device Selection

1. Choose appropriate size of device for the procedure.

Device Selection	19 Ga Flexible	19 Ga Standard	22 Ga	25 Ga
Sampling of submucosal and extramural gastrointestinal lesions (fine needle aspiration)	X	X	X	X
Delivery of injectable materials (fluids)	X	X	X	
Delivery of fiducials (luer loading)	X	X	X	
Delivery of fiducials (loading via needle tip)	X	X	X	
Passage of accessory devices*	X			

*19 Ga Flexible has demonstrated to be compatible for the delivery of the Mauna Kea Cellvizio AQ-Flex™ 19 Confocal Miniprobe™ accessory device.

See below table for needle size by fiducial dimensions:

Needle Size	Maximum Fiducial Outer Diameter (mm)	Maximum Fiducial Length (mm)
19 Ga Flexible and 19 Ga	0.80	10
22 Ga	0.35	10

Fiducials are not recommended for use with the 25 ga needle.

2. Inspect Boston Scientific device to ensure it is within acceptable shelf life constraints.

Product Preparation

1. Open device package and remove plastic insert containing FNA device and syringe.
2. Remove the device and syringe from the package.

Caution: Visually inspect the stopcock when removed from the package to confirm it is in the open position; otherwise, DO NOT USE IT. Call Boston Scientific Customer Service and return the product.

Caution: Visually inspect the device for loose, bent or broken parts, cracks, or other abnormalities. Inspect the catheter for any kinks or other damage. If an abnormality is detected, DO NOT USE. Broken parts, cracks, or kinks will hinder the mechanical operation of the Expect Slimline (SL) Needle. If the device fails to operate properly in any way or shows any signs of damage, DO NOT USE IT. Call Boston Scientific Customer Service and return the product.

Syringe Preparation

[For following indications: Sampling of submucosal and extramural gastrointestinal lesions (fine needle aspiration), delivery of injectable materials (fluids)]


Figure 1.

- 1(a): Normal Sliding Use 1(c): Locking Fins
1(b): Locks to Hold Vacuum 1(d): Stop Pin

1. Examine the stopcock (Figure 1(a)). The stopcock has two luer connections to attach to the needle and syringe. Air can be exchanged with the stopcock in the open position. The stopcock is open when it is aligned parallel to the syringe; it is closed when perpendicular to the syringe (as shown).

- Examine the syringe (Figure 1). The syringe barrel has one stop pin (Figure 1 (d)) and the syringe plunger has four locking fins (Figure 1(c)). The syringe plunger can be maneuvered within the syringe barrel to lock and unlock the syringe. To lock the syringe, pull back on the plunger until it aligns with the desired suction volume (Figure 1(a)). Turn the plunger clockwise so that the locking pin engages with the locking fins on the plunger (Figure 1(b)). Turn the plunger counterclockwise to release.

Note: The stopcock is required in order to maintain suction during the procedure.

Caution: If the stopcock is not set properly, adequate suction may not be achieved.

Product Use


Figure 2.

- | | |
|---------------------------|--------------------------------------|
| 2(a): Luer Attachment | 2(d): Needle Adjustment Lock |
| 2(b): Sheath Adjust Gauge | 2(e): Needle Depth Gauge |
| 2(c): Sheath Adjust Lock | 2(f): Aspiration Port and Stylet Cap |

Directions for Use

For additional instructions per indication, refer to specific instructions below.

Preparation Instructions:

- Remove the needle from the package and examine it for damage.
- Confirm that the needle is fully retracted and that the needle adjustment lock (Figure 2(d)) is secure in the zero position (Figure 2(e)).
- Determine the appropriate sheath length relative to the length of the echoendoscope. Use the sheath adjustment lock (Figure 2(c)) to set the desired length of the sheath and lock it in place. Turn the sheath adjustment lock clockwise to lock the sheath in place. The distal end of the sheath should be visible on the endoscopic image.

Note: The reference numbers and markings on the sheath adjust gauge (Figure 2(b)) are intended as reference only.

- Turn the endoscope elevator control knob to lower the elevator.

Caution: Failure to lower the elevator prior to insertion may cause damage to the device.

- Introduce the catheter into the echoendoscope working channel. Slowly advance the device through the working channel on the echoendoscope.

Caution: If resistance is encountered, stop pushing and reposition the catheter or the echoendoscope. Pushing too hard could cause damage to the echoendoscope.

- Continue to advance the device until it exits the echoendoscope and the sheath is visible in the endoscopic image. Tighten the luer attachment (Figure 2(a)) by turning clockwise to attach the device to the working channel port of the echoendoscope.

Note: If it is necessary to adjust the sheath length, loosen the sheath adjustment lock (Figure 2(c)) and reposition it to the appropriate reference point.

Caution: Do not tighten the luer connection too tightly to the echoendoscope as it may cause damage to the echoendoscope.

Caution: Prior to advancing the needle, ensure that the device is securely fastened to the echoendoscope and both the needle adjustment and sheath adjustment locks are secure. Failure to do so could result in damage to the echoendoscope.

Caution: Insufflation performance may be diminished when device is attached to the echoendoscope.

- Verify the distance from the distal end of the sheath to the target site using the ultrasound image.
- Adjust the needle penetration depth to the desired position using the needle adjustment lock (Figure 2(d)). To control the depth of needle penetration to the target site, loosen the needle adjustment lock by turning the knob counterclockwise. Align the needle adjustment lock with the appropriate reference number on the device handle. Lock in place by turning the needle adjustment lock clockwise.

Note: The reference numbers on the needle depth gauge (Figure 2(e)) are intended as reference only. The reference numbers represent needle extension (in centimeters) when the device is in a straight position.

- Advance the needle by sliding the handle toward the echoendoscope in a slow and controlled motion to penetrate the target site while observing the ultrasound image.

Note: If it is necessary to adjust the amount of needle penetration, loosen the needle adjustment lock, reposition, and lock it to the appropriate reference point.

- Remove the stylet (Figure 2(f)) from the aspiration port of the device by gently loosening and pulling it from the device handle.

Note: For delivery of fiducial via needle tip, stylet should be positioned per fiducial manufacturer's instructions. Continue to instructions below for delivery of fiducials (loading via needle tip).

Note: The stylet can be coiled and clipped together using the clip mechanism on the stylet cap (Figure 2(f)). The stylet should be retained for additional needle passes during the same procedure.

Caution: Failure to properly handle the stylet could result in damage to the device.

Caution: The stylet tip is sharp. Take precautions to ensure that the stylet is handled properly. After it has been removed from the needle, the stylet should be treated as infectious material and could create an infection risk.

After each procedure has been performed, follow instructions below for removal of needle.

Needle Removal:

1. Retract the needle fully into the sheath using the device handle by sliding the handle away from the echoendoscope until it stops moving. Secure the needle using the needle adjustment lock prior to withdrawing the device from the echoendoscope.

Caution: Ensure that the needle is fully retracted into the sheath. Failure to secure the needle could result in damage to the echoendoscope or injury to user.

2. Lower the elevator on the echoendoscope.

Caution: Failure to lower the elevator prior to withdrawal may cause damage to the device.

3. Detach the device from the echoendoscope by turning the luer connection counterclockwise. Slowly and carefully remove the device from the echoendoscope.

After device has been removed from endoscope, follow instructions below for disposal of needle.

Needle Disposal:

1. After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

Caution: Do not use the same device to perform multiple indications. Use of the device for multiple indications has not been established.

Specific Instructions per Indication:

Sampling of submucosal and extramural gastrointestinal lesions (fine needle aspiration)

1. Prepare device, supplied syringe, and stopcock as previously noted. Rotate the stopcock to the closed position, perpendicular to the syringe. Pull the syringe plunger to the desired volume and use the locking fins and stop pin on the syringe to secure the plunger in place. (Figures 1(c) and 1(d)).

Caution: Methods of providing suction other than the supplied syringe are not recommended with this device.

2. Connect the supplied syringe to the aspiration port (Figure 2 (f)) on the device handle.
3. Turn the stopcock to the open position (parallel to the device handle) to apply suction.
4. Maneuver the needle within the target site to maximize aspiration sample collection while observing needle penetration on the ultrasound image.
5. After an adequate number of passes have been made with the needle, close the stopcock by rotating until it is perpendicular to the syringe. This will stop suction.
6. Remove needle from echoendoscope (Refer to steps above).
7. After the device has been removed from the echoendoscope, release the needle adjustment lock and advance the device handle to extend the needle out of the sheath.

8. Remove the syringe and stopcock from the aspiration port.
9. Open the stopcock on the syringe by rotating it to the position that is parallel to the syringe. Pull the syringe plunger back to pull air into the syringe.
10. Reconnect the syringe to the aspiration port.
11. Push the syringe plunger forward to expel the aspiration sample from the needle.

Note: The locking fins and stop pin (Figures 1(c) and 1(d)) will have to be disengaged in order to pass air from the syringe.

Caution: Take precautions to ensure that the aspiration sample does not spray when it is expelled from the needle. The aspiration sample should be treated as infectious material and could create an infection risk.

12. Prepare the aspiration sample per institutional protocol.

Note: Multiple passes may be required to obtain an adequate aspiration sample.

13. If additional passes to the same target site are required, prepare the device by flushing the needle and wiping the stylet with sterile water or saline. Reinsert the stylet into the needle, examine the needle for damage, and repeat instructions.

Caution: Failure to flush the needle and wipe the stylet prior to reinserting the stylet into the needle could make the stylet difficult to pass or result in damage to the device.

14. Dispose of needle (Refer to steps above).

Delivery of injectable materials (fluids)

1. Prepare device and load the supplied syringe with fluid.
2. Connect the supplied syringe to the aspiration port (Figure 2(f)) on the device handle.
3. Use the supplied syringe to inject desired quantity of fluid into the targeted injection site.
4. Remove needle from echoendoscope (Refer to steps above).
5. Determine if additional fluid needs to be injected. If so, prepare the device by flushing the needle with sterile water or saline, repeat instructions for delivery of injectable materials (fluids).
6. Dispose of needle (Refer to steps above).

Delivery of fiducials (luer loading)

1. While preparing device per initial instructions and prior to turning the endoscope elevator control knob to lower the elevator, follow manufacturer's instructions to prepare the fiducial transfer device to load the fiducial.
2. Per the manufacturer's instructions, deliver the fiducial.

Caution: Ensure the stylet is not inadvertently advanced prior to or during needle insertion into the lesion. Advancing the stylet prior to or during needle insertion into the lesion may result in premature deployment of the fiducial.

Caution: Severe flexion of the endoscope may result in inability to deliver the fiducial.

3. Determine if additional fiducial needs to be placed. If so, repeat instructions for delivery of fiducials (luer loading).
4. Remove needle from echoendoscope (Refer to steps above).
5. Dispose of needle (Refer to steps above).

Delivery of fiducials (loading via needle tip)

Note: Stylet should be positioned per fiducial manufacturer's instructions.

1. Prepare device and prior to turning the endoscope elevator control knob to lower the elevator, follow manufacturer's instructions to load and deliver the fiducial.

Caution: Ensure the stylet is not inadvertently advanced prior to or during needle insertion into the lesion. Advancing the stylet prior to or during needle insertion into the lesion may result in premature deployment of the fiducial.

Caution: Severe flexion of the endoscope may result in inability to deliver the fiducial.

2. Remove needle from echoendoscope (Refer to steps above).

Warning: The needle tip is sharp. Take precautions to ensure that the needle is handled properly. The needle should be treated as infectious material and could create an infection risk.

3. Determine if additional fiducial needs to be placed. If so, prepare the device by flushing the needle with sterile water or saline, repeat instructions for delivery of fiducials (loading via needle tip).

Caution: Failure to properly clean and flush the needle could lead to difficulty in completing the procedure as intended.

4. Dispose of needle (Refer to steps above).

Passage of accessory devices

Note: The needle has demonstrated to be compatible for the delivery of the Mauna Kea Cellvizio AQ-Flex™ 19 Confocal Miniprob™ accessory device.

1. Prepare device.
2. Follow manufacturer's instructions to select, prepare and preload the accessory device into the needle prior to turning the endoscope elevator control knob to lower the elevator.
3. Follow manufacturer's instructions for operating and removal of accessory device.
4. Remove needle from echoendoscope and dispose of needle (Refer to steps above).

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

Mauna Kea Cellvizio AQ-Flex™ 19 Confocal Miniprob™ is a trademark of Mauna Kea Technologies.

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© 2016 Boston Scientific Corporation or its affiliates. All rights reserved.
Boston Scientific (Australia) Pty Ltd