

FRONTLINE GUIDEWIRE

Savion FLX™ Guidewires

UPN Order Code	Diameter (in)	Length (cm)	Tip Load (grams)	Tip Shape
H7493942318510	0.014	185	2.9	Straight
H749394231851J0	0.014	185	2.9	Angled
H7493942230010	0.014	300	2.9	Straight
H749394223001J0	0.014	300	2.9	Angled

DELIVERY SUPPORT GUIDEWIRE

Savion DLVR™ Guidewires

UPN Order Code	Diameter (in)	Length (cm)	Tip Load (grams)	Tip Shape
H7493942118210	0.014	182	0.8	Straight
H749394211821J0	0.014	182	0.8	Angled
H7493942030010	0.014	300	0.8	Straight
H749394203001J0	0.014	300	0.8	Angled

SAVION DLVR SAVION FLX Guidewires

**INTENDED USE/INDICATIONS FOR USE:** Boston Scientific SAVION DLVR™ and SAVION FLX™ Guidewires are intended to facilitate the placement of balloon dilatation catheters or other interventional therapeutic devices during PTCA, PTA, or other intravascular interventional procedures. They are not intended for use in the cerebral vasculature.

**CONTRAINDICATIONS:** None known. **WARNINGS:** Use extreme caution and careful judgment in patients for whom anticoagulation is not indicated. Severe reaction may occur in response to contrast agents that cannot be adequately premedicated. When the guidewire is in the body, it should be manipulated only under fluoroscopy. The hydrophilic coating of these guidewires increases the possibility of vessel wall perforation compared to non-hydrophilic coatings. **PRECAUTIONS:** Carefully check and match therapeutic device compatibility to the wire prior to use. **ADVERSE EVENTS:** Potential adverse events which may result from the use of the device include but are not limited to: Abrupt closure • Allergic reaction (to contrast, device or other) • Angina or unstable angina • Arrhythmias • Bleeding/Hemorrhage • Cardiac tamponade/pericardial effusion • Death • Embolization (plaque, thrombus, device, tissue, or other) • Hematoma • Infection/Sepsis • Myocardial infarction or ischemia • Pain at the access site • Pseudoaneurysm • Renal insufficiency or renal failure • Stroke/cerebral vascular accident (CVA)/transient ischemic attack (TIA) • Thrombosis/Thrombus • Vasospasm • Vessel trauma (dissection, perforation, rupture or injury). Some of the above potential adverse events may require additional urgent intervention or surgery.

**CAUTION:** Federal law (USA) restricts this device to sale by or on the order of a physician. Rx only. Prior to use, please see the complete "Directions for Use" for more information on Indications, Contraindications, Warnings, Precautions, Adverse Events, and Operator's Instructions.



**Peripheral Intervention**  
300 Boston Scientific Way  
Marlborough, MA 01752-1234  
**www.bostonscientific.com**

*To order product or for more information  
contact customer service at 1.888.272.1001.*

© 2018 Boston Scientific Corporation  
or its affiliates. All rights reserved.



SAVION FLX™ & SAVION DLVR™ GUIDEWIRES

TAKE  
CONTROL OF  
COMPLEXITY

Building on trusted technology

Savion Guidewires exhibit attributes ideally matched for BTK interventions to help you treat complex lesions successfully and restore blood flow to the leg.



# Savion FLX™

The **Savion FLX** guidewire allows you to take control—quickly, safely, and reliably—in challenging below-the-knee (BTK) cases. Featuring a resilient nitinol tip, this frontline wire enables you to confidently navigate your preferred path to reach the most distal and complex lesions without hesitation.


## Moderate Rail Support



# Savion DLVR™


The strength and structure of the **Savion DLVR** guidewire provides the confidence you need to deliver targeted therapies distally. Expect maximum support with easier exchanges for your toughest cases, including pedal and retrograde access.

## Maximum Rail Support



**Exceptional tip durability**  
Nitinol tip provides superior durability and resilience for the duration of most cases

Images of wire distal tips are captured before and after being advanced 20 mm into a cylindrical tube.




**Superior shape retention**  
Stainless steel shaping ribbon helps enable tip shape retention throughout the procedure

**Convenient crossing**  
Polymer sleeve tip with ICE™ hydrophilic coating facilitates lesion crossing

**Precise control**  
1:1 torque response for control\*

\* Bench Test results may not necessarily be indicative of clinical performance. Data on file.

**Improved performance**  
Maximum rail support with gradual distal transitions for improved performance



**Accurate feedback**  
Uncoated 3 cm atraumatic distal tip designed for precise tactile feedback

**Smooth tracking**  
ICE™ hydrophilic coating for excellent lubricity and smooth device tracking

\* excluding distal 3 cm (uncoated)