**Boston  
Scientific**

**Contour™**  
Embolization Particles

90495391-01

2010-03

<b>Directions for Use</b>	<b>2</b>
<b>Instrucciones de uso</b>	<b>13</b>
<b>Mode d'emploi</b>	<b>24</b>
<b>Gebrauchsanweisung</b>	<b>35</b>
<b>Istruzioni per l'uso</b>	<b>47</b>
<b>Gebruiksaanwijzing</b>	<b>58</b>
<b>Instruções de Utilização</b>	<b>69</b>

**TABLE OF CONTENTS**

<b>WARNING .....</b>	<b>3</b>
<b>DEVICE DESCRIPTION .....</b>	<b>3</b>
<b>INTENDED USE/INDICATIONS FOR USE .....</b>	<b>3</b>
<b>CONTRAINDICATIONS .....</b>	<b>3</b>
<b>CONTRAINDICATIONS SPECIFIC TO UTERINE FIBROID EMBOLIZATION (UFE) .....</b>	<b>3</b>
<b>POTENTIAL COMPLICATIONS .....</b>	<b>4</b>
<b>POTENTIAL COMPLICATIONS SPECIFIC TO UFE .....</b>	<b>4</b>
<b>WARNINGS AND PRECAUTIONS .....</b>	<b>4</b>
<b>PRECAUTIONS APPLICABLE TO ALL INDICATIONS .....</b>	<b>5</b>
<b>INSTRUCTIONS FOR USE .....</b>	<b>5</b>
<b>HOW SUPPLIED .....</b>	<b>6</b>
<b>CLINICAL STUDIES .....</b>	<b>6</b>
<b>SECONDARY ENDPOINTS INCLUDING SAFETY SUMMARY .....</b>	<b>9</b>
<b>ADVERSE EVENTS .....</b>	<b>9</b>
<b>LONG-TERM FOLLOW-UP .....</b>	<b>10</b>
<b>REFERENCES .....</b>	<b>11</b>
<b>WARRANTY .....</b>	<b>12</b>

# Contour™

## Embolization Particles

**Rx ONLY**

**Caution:** Federal Law (USA) restricts this device to sale by or on the order of a physician.

### WARNING

Contents supplied STERILE using a Radiation process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

### DEVICE DESCRIPTION

The Contour Embolization Particles are irregularly-shaped polyvinyl alcohol embolization devices. These devices are intended to provide vascular occlusion or reduction of blood flow within target vessels upon selective placement through a variety of infusion catheters.

The Contour Embolization Particles are packaged sterile, 1 cm<sup>3</sup> (1 cc) dry volume per vial, in a peel-away pouch.

Each vial is intended for single use only. Do not resterilize. Discard any unused material.

The Contour Embolization Particles size ranges are listed below and are available either as a box of two vials or a box of five vials.

### Contents

- Vials of Contour Embolization Particles

Size	Minimum Compatible Catheter ID
45-150 µm 150-250 µm 250-355 µm 355-500 µm	0.53 mm (0.021 in) (e.g., Renegade™ 18 Microcatheter, Renegade STC Microcatheter)
500-710 µm	0.69 mm (0.027 in) (e.g., Renegade HI-FLO™ Microcatheter)
710-1000 µm 1000-1180 µm	1.12 mm (0.044 in) (e.g., Imager™ II Catheter)

Contour Embolization Particles Size Range - Minimum Catheter I.D. Required.

### INTENDED USE/INDICATIONS FOR USE

The Contour Embolization Particles are used for the embolization of peripheral hypervascular tumors, including leiomyoma uteri and peripheral arteriovenous malformations (AVMs).

Do not use particles smaller than 355 microns for the treatment of leiomyoma uteri.

### CONTRAINDICATIONS

Contraindications Specific to All Peripheral Indications

1. Vascular anatomy or blood flow precludes stable, selective Contour Embolization Particles or catheter placement
2. Presence of vasospasm
3. Presence of hemorrhage
4. Presence of severe atheromatous disease
5. Presence of feeding arteries smaller than distal branches from which they emerge
6. Presence of collateral vessel pathways potentially endangering normal territories during embolization
7. Presence of arteries supplying the lesion not large enough to accept Contour Embolization Particles
8. Vascular resistance peripheral to the feeding arteries precluding passage of Contour Embolization Particles into the lesion
9. In large diameter arteriovenous shunts
10. In the pulmonary vasculature
11. Patient intolerance to occlusion procedures

### CONTRAINDICATIONS SPECIFIC TO UTERINE FIBROID EMBOLIZATION (UFE)

1. Pregnant women
2. Suspected pelvic inflammatory disease or any other active pelvic infection
3. Any malignancy of the pelvic region
4. Endometrial neoplasia or hyperplasia

5. Presence of one or more submucosal fibroid(s) with more than 50% growth into the uterine cavity
6. Presence of pedunculated serosal fibroid as the dominant fibroid(s)
7. Fibroids with significant collateral feeding by vessels other than the uterine arteries

#### POTENTIAL COMPLICATIONS

Complications specific to embolization include, but may not be limited to:

1. Foreign body reactions (i.e. pain, rash) necessitating medical intervention
2. Allergic reaction to contrast media
3. Infection necessitating medical intervention
4. Complications related to catheterization (e.g., hematoma at the site of entry, clot formation at the tip of the catheter and subsequent dislodgment, vasospasm and nerve and/or circulatory injuries, which may result in leg injury)
5. Undesirable reflux or passage of Contour™ Embolization Particles into arteries adjacent to the targeted lesion or through the lesion into other arteries or arterial beds
6. Ischemia at an undesirable location
7. Incomplete occlusion of vascular beds or territories may give rise to the possibility of post procedural hemorrhage, development of alternative vascular pathways, recanalization or recurrence of symptoms
8. Vessel or lesion rupture and hemorrhage
9. Recurrent hemorrhage
10. Ischemic stroke or myocardial infarction
11. Death

#### POTENTIAL COMPLICATIONS SPECIFIC TO UFE

1. Postembolization syndrome
2. Vaginal Discharge
3. Tissue passage, fibroid sloughing or fibroid expulsion post-UFE
4. Temporary or permanent stopping of menstrual bleeding
5. Infection of the pelvic region
6. Endometrial atrophy with amenorrhea despite normal ovarian function
7. Complications to pregnancy
8. Premature Ovarian Failure (i.e., menopause)
9. Uterine/Ovarian necrosis
10. Uterine Rupture
11. Post-UFE Intervention to remove necrotic fibroid tissue
12. Hysterectomy

---

#### WARNINGS AND PRECAUTIONS

Warnings Applicable to All Peripheral Indications

- **PRIOR TO EMBOLIZATION, PROSPECTIVE PATIENTS OR THEIR REPRESENTATIVES MUST BE PROVIDED AN INFORMED CONSENT DESCRIBING THE POSSIBLE COMPLICATIONS ASSOCIATED WITH THE USE OF THIS DEVICE. WRITTEN ACKNOWLEDGMENT IS WARRANTED.**
- The safety and effectiveness of Contour Embolization Particles for neurovascular use have not been established.
- As with any embolization device, patient injury, permanent disability or death may occur as a result of its use.
- Vascular occlusion should only be performed by physicians possessing skilled interventional occlusion experience in the territory intended to be embolized.
- A thorough evaluation of a patients medical condition, vascular pathways and the desired embolization goal is necessary to achieve successful occlusion. This evaluation should include baseline angiography to determine the presence of potentially dangerous collateral pathways. Do not proceed with embolization unless these pathways can be protected.
- Do not use if protective sterile packaging is opened or damaged.
- Smaller Contour Embolization Particles may be more likely to result in ischemic infarction because of the potential to block vessels at the precapillary level; however, Contour Embolization Particles of all sizes share this potential.
- Typically the artery will accept less Contour Embolization Particles as the treatment progresses. Proximal slowing or termination of Contour Embolization Particles passage may occur when the vessel or malformation is occluded by prior Contour Embolization Particles, or in the presence of severe atheromatous disease. Continued infusion may result in inadvertent reflux into critical arteries, creating the potential for undesirable ischemic infarction.
- Exercise conservative judgment in determining embolization endpoint. Terminate infusion before complete vessel occlusion has occurred.
- Avoid the use of a tapered tip delivery catheter. The Contour Embolization Particles may lodge at the catheter tip creating an obstruction.
- “Clumping” of Contour Embolization Particles or catheter obstruction may be a function of contrast dilution volume; ensure that sufficient volume of the appropriate contrast/saline mix is utilized such that Contour Embolization Particles are free floating and not observed as aggregates.
- Should catheter obstruction occur, remove the catheter from the patient. Do not use forceful injection, guidewires or other instruments to dislodge the blockage. Do not continue using a catheter which has been obstructed as damage to the device may have occurred.
- Incomplete occlusion of vascular beds or territories may give rise to the possibility of post procedural hemorrhage, development of alternative vascular pathways, or recurrence of symptoms.
- Postprocedural patient follow-up to assess the continued level of vascular occlusion is necessary. Angiography may be indicated.

- Serious radiation induced skin injury may occur to the patient due to long periods of fluoroscopic exposure, large patient diameter, angled x-ray projections, and multiple image recording runs or radiographs. Refer to your facility's institutional protocol to ensure the proper radiation dose is applied for each specific type of procedure performed. Physicians should monitor patients that may be at risk.
- Onset of radiation-induced injury to the patient may be delayed. Patients should be counseled on potential radiation side effects and whom they should contact if they show symptoms.

---

#### PRECAUTIONS APPLICABLE TO ALL INDICATIONS

- Patients with known allergy to contrast medium may require pre-medication prior to embolization.
- Additional evaluations or precautions may be necessary in managing periprocedural care for patients with the following conditions:
  - A. Bleeding diathesis or hypercoagulable state;
  - B. Immunocompromised.
- Sterile and single use product. Never reuse a vial that has been opened. Do not use if the vial, screw cap, or tray package appears damaged.
- Ensure cleanliness and attention to technique during preparation of the device to avoid introducing contaminants.
- The appropriate size particles must be chosen based upon the lesion to be treated and the measurements taken from the baseline angiography.
- The use of sophisticated imaging equipment is necessary for successful embolization therapy.
- Appropriate facilities should be available to treat potential complications of the procedure.
- While it is anticipated that long-term embolization of vascular structures with Contour™ Embolization Particles will be achieved, no guarantee of permanence, cure or benefit can be made.

#### UFE Specific Warnings for Pregnancy (Specific for Treatment of Leiomyoma Uteri)

- UFE is not intended for women who desire future pregnancy. The effects of UFE on the ability to become pregnant and carry a fetus to term, and on the development of the fetus, have not been determined. Therefore, this procedure should only be performed on women who do not intend future pregnancy.
- Women who become pregnant following UFE, should be aware that they may be at increased risk for preterm delivery, cesarean delivery, malpresentation (incorrect positioning of the baby), and postpartum hemorrhage (post-delivery bleeding).
- Devascularization of uterine myometrium resulting from UFE may put women who become pregnant following UFE at increased risk of uterine rupture.

#### Other UFE Specific Warnings

- Do not use particles smaller than 355 microns.
- An appropriate gynecologic work-up should be performed on all patients presenting for embolization of uterine fibroids (e.g., endometrial sampling to rule out carcinoma for patients with abnormal bleeding).

- Devascularization of uterine myometrium resulting from UFE may put women at increased risk of uterine rupture.
- The diagnosis of uterine sarcoma could be delayed by taking a non-surgical approach (such as UFE), to treat uterine fibroids. Conduct a more thorough work-up for patients with warning signs for sarcoma (e.g., prior pelvic radiation, MRI findings, rapid tumor growth, postmenopausal with new uterine enlargement). Recurrent or continued tumor growth following UFE should be considered a potential warning sign for sarcoma and surgery should be considered.

#### UFE Specific Precautions

- It is recommended that patients undergoing embolization of leiomyoma uteri be provided a clear understanding of who will provide post-procedure care prior to the embolization procedure.
- UFE should only be performed by physicians who have received appropriate training for treatment of uterine leiomyomata (fibroids).
- There is an increased chance of retro-migration of Contour Embolization Particles into unintended blood vessels as uterine artery flow diminishes. Embolization should be stopped when the vasculature surrounding the fibroid can no longer be visualized but before complete stasis in the uterine artery.

#### INSTRUCTIONS FOR USE

1. Inspect packaging prior to use to ensure seal integrity for maintenance of sterility.
2. Access the femoral artery using the standard Seldinger technique. Insert the valved sheath of choice to aid in control of arterial backflow.
3. Perform baseline angiographic evaluation of vascular network associated with the target lesion prior to beginning embolization procedure.
4. Choose appropriate Contour Embolization Particle size based on clinical presentation including vascular target/vessel size and desired clinical outcome.
5. Choose appropriate delivery catheter per standard protocol.
6. Under fluoroscopy, advance the delivery catheter to the desired lesion according to standard technique.
  - Inject contrast agent through the delivery catheter to confirm position and flow.
7. Open Contour Embolization Particles vial and pour contents into sterile container (e.g., bowl or barrel of a syringe).
  - Mix contents of one vial of Contour Embolization Particles with contrast agent/ saline to form a dilute solution.
  - Agitate by inverting syringe several times or gently use back-and-forth agitation through a stopcock or stir mixture in bowl.

If using syringe method, continue to next step. If using mixing bowl method: Draw the solution into delivery syringe. If a 3-way stop-cock delivery method is preferred, draw solution into "reservoir" syringe first. Ensure that all air is purged from the system prior to injection. Connect the reservoir syringe and/or delivery syringe as appropriate to the stop-cock and/or catheter. Observe contents of delivery syringe to ensure that Contour Embolization Particles are well distributed and suspended.

8. For all indications: Carefully inject Contour™ Embolization Particles under fluoroscopic visualization, observing contrast flow rate. Appreciate a reduction in flow acceptance, and decrease infusion accordingly. Exercise conservative judgment in determining embolization endpoint. Terminate infusion before complete vessel occlusion has occurred.

**The UFE specific end point** is generally described as complete stasis or near stasis, with the main uterine artery remaining patent, but with negligible residual flow. This endpoint generally corresponds to an angiographic image of a patent uterine artery with all its distal branches occluded. As with any embolic particle, in order to avoid a false endpoint with early recanalization, the embolization endpoint should be confirmed by leaving the catheter in the uterine artery for approximately 5 minutes after the apparent conclusion of the procedure. The endpoint should then be confirmed with an injection of contrast and observation with fluoroscopy. Additional particles can then be administered to reach the stated endpoint if flow restoration due to redistribution is identified on this contrast injection.

---

**Caution:** “Clumping” of Contour Embolization Particles or catheter obstruction may be a function of contrast dilution volume; ensure that enough contrast is utilized such that Contour Embolization Particles are free floating and not observed as aggregates.

---

9. Should catheter obstruction occur, remove the catheter from the patient. Do not use forceful injection, guidewires or other instruments to dislodge the blockage. Do not continue using a catheter which has been obstructed as damage to the device may have occurred.
10. Upon completion of the procedure, remove catheter from the patient with gentle suction to ensure any remaining particles are maintained within the catheter lumen as it is withdrawn.

#### HOW SUPPLIED

Boston Scientific products are sterile in unopened, undamaged packaging. Packaging is designed to maintain sterility unless the primary product pouch has been opened or damaged. Do not use if package is opened or damaged. Do not use if labeling is incomplete or illegible. Store in a cool, dry, dark place.

#### CLINICAL STUDIES

##### THE Contour Embolization Particles Clinical Study Summary

##### UFE Clinical Study Summary

##### Study Design

A prospective, multi-center trial was conducted to study uterine fibroid embolization (UFE) using Contour Embolization Particles for the treatment of symptomatic uterine fibroids. A total of 209 women were administered anesthesia and were defined as the intent-to-treat population (ITT), 149 in the UFE arm and 60 in the myomectomy control arm. Of the 149 UFE subjects, 2 subjects were never treated with embolic, making the number of subjects treated for their symptomatic uterine fibroids 147 for Contour Embolization Particles and 60 women were treated for their symptomatic uterine fibroids in a concurrent, non-randomized control arm with either abdominal or laparoscopic myomectomy. 16 sites participated in the study, 11 of which performed UFEs and 10 of which performed myomectomy.

The study was designed to evaluate the safety and effectiveness of Contour Embolization Particles when used in the uterine fibroid embolization procedure for the management of symptomatic uterine fibroids.

##### Primary Study Endpoint

The primary endpoint evaluated in this study was the improvement of fibroid related symptoms from baseline to 6-months post-procedure using the Overall Fibroid Symptom component of a validated Uterine Fibroid Quality of Life (UFQoL) Questionnaire. Subjects in the UFE cohort acted as their own control and were considered a success at 6-months post-procedure if they had an increased score from baseline to 6-months post-procedure of  $\geq 5$  points on the Overall Fibroid Symptom component of the UFQoL which includes assessment of:

- Heavy bleeding
- Abdominal swelling
- Pelvic pain
- Backaches
- Bowel problems
- Frequent urination
- Pain with urination
- Poor appetite
- Hot flashes
- Vaginal dryness
- Painful period

##### Secondary endpoints included:

- Overall adverse event rates for both the myomectomy and UFE cohorts
- Success for Overall Fibroid Symptom scores between myomectomy and UFE cohorts at 6-months post-procedure
- Quality of Life measures for UFE and myomectomy cohorts at 3 and 6-months post-procedure
- Time to return to normal activities between UFE and myomectomy cohorts
- Change in dominant fibroid size at 6-months post-procedure as compared to pre-procedure (UFE cohort)
- Relationship between change in dominant fibroid size and fibroid specific quality of life measures (UFE cohort)

##### Subject eligibility criteria for the study included:

##### Inclusion Criteria

- Diagnosis of symptomatic uterine myoma
- Record of a normal pap smear within the last 12 months
- Menstrual cycles within every 22 to 35 days (pre-menopausal)
- Age  $\geq 30$
- If subject was currently not taking hormonal contraceptives, subject had been off hormonal contraceptive for a minimum of 3 months prior to study enrollment and

agrees to continue without change in the regimen for the study duration (through 6-month follow-up)

- Subject was willing and able to complete the follow-up requirements outlined in the study design section of the protocol
- Subject was willing to sign a consent form

**Subjects in the UFE arm had to meet all of the following additional criteria to be eligible to be enrolled into the study:**

- Subject did not desire pregnancy
- Creatinine levels were within normal ranges
- If subject is currently on hormonal contraceptives, subject had taken hormonal contraceptive for  $\geq 3$  months prior to study enrollment and agreed to continue without change in the regimen for the study duration (through 6-month follow-up)

**Subjects in the myomectomy arm could also meet the following additional criteria and still be eligible to be enrolled into the study:**

- Subjects whose symptom was infertility related to fibroids
- Subjects to be treated GnRH agonists e.g., Lupron Depot™ Pharmaceutical (leuprolide acetate for depot suspension) may be enrolled if their therapy is started prior to the procedure but after taking the UFQoL questionnaire and remaining within the 60-day window
- Fibroids could be removed via laparoscopic myomectomy

**Exclusion Criteria**

**Subjects who met any of the following criteria were not enrolled into the study:**

- Fibroid(s) could only be removed by hysteroscopic myomectomy
- Active pelvic infection
- Known or suspected gynecologic malignancy
- Undiagnosed pelvic mass outside of the uterus (not expected to be functional in nature).
- Unexplained abnormal menstrual bleeding
- Acute or chronic infection (e.g., kidney, bladder, lung, etc.)
- Subject with coagulopathy
- Prior history of pelvic irradiation
- American Society of Anesthesiologists (ASA) score greater than or equal to IV
- Random FSH value  $\geq 40$  IU/L
- Participating in any other investigational device or drug study

**Subjects in the UFE arm who met any of the following additional criteria were not enrolled into the study:**

- Uterine arteriovenous fistula.
- Known severe contrast allergy.
- Pedunculated subserosal fibroid (attachment to uterus is less than 30% of the greatest diameter of the fibroid).

**Additional Inclusion Criteria to Continue Participation:**

- Subjects were required to complete the baseline UFQoL in order to participate in the study.
- Subjects were required to have a baseline score on the Overall Fibroid Symptom component of the UFQoL of  $\leq 90$  points with the exception of myomectomy subjects whose fibroids were a cause of infertility and who may not present with the symptoms of bulk/pressure or abnormal menstrual bleeding.
- Subjects from all cohorts could not be pregnant and this had to be confirmed by a pregnancy test on the day of, but prior to the procedure.

For subjects that met all of the inclusion and none of the exclusion criteria, subject demographics (Table 1), medical history, gynecological history, menstrual history, symptom assessment, imaging evaluation, and pre-procedure labs were collected prior to treatment.

All subjects had to have a baseline magnetic resonance imaging (MRI) prior to the UFE or myomectomy procedure. All MRI films were sent to a core laboratory for evaluation in order to provide consistent and non-biased measurement analysis.

Subjects were evaluated at 3 and 6 months during the study. The subjects in the UFE arm were then followed annually for up to three years to assess their long-term outcomes.

**Table 1. Subject Demographics**

	UFE cohort (N=149)	Myomectomy Cohort (N=60)
<b>Age</b>		
Mean	43.9	38.2
Std Dev	3.58	5.16
<b>Race/Ethnicity</b>		
Caucasian	71 (47.7%)	22 (36.7%)
African American	65 (43.6%)	29 (48.3%)
Hispanic	8 (5.4%)	4 (6.7%)
Asian	3 (2.0%)	3 (5.0%)
Native American	0 (0.0%)	1 (1.7%)
Other	2 (1.3%)	1 (1.7%)
<b>Height (inches)</b>		
Mean	64.9	64.4
Std Dev	2.78	3.23
<b>Weight (lbs)</b>		
Mean	159.7	165.9
Std Dev	32.10	41.45
<b>Prior Fibroid Treatment</b>		
Previous Myomectomy	24 (16.1%)	7 (11.7%)

	UFE cohort (N=149)	Myomectomy Cohort (N=60)
Myolysis	2 (1.3%)	0 (0.0%)
Dilatation & Curettage	44 (29.5%)	14 (23.3%)
Ovarian Cystectomy	6 (4.0%)	2 (3.3%)
Endometrial Ablation or Fulguration	6 (4.0%)	1 (1.7%)
Tubal Ligation	28 (18.8%)	1 (1.7%)
Salpingo-Oophorectomy	1 (0.7%)	2 (3.3%)
Other	30 (20.1%)	17 (28.3%)
<b>Length of Menstrual Bleeding (days)</b>		
Mean	7.0	6.2
Std Dev	2.67	2.47
<b>Dominant Symptom</b>		
Abnormal Bleeding	77 (51.7%)	20 (33.3%)
Bulk/Pressure	38 (25.5%)	16 (26.7%)
Pelvic Pain	29 (19.5%)	18 (30.0%)
Infertility	0 (0.0%)	2 (3.3%)
Other	5 (3.4%)	4 (6.7%)
<b>Duration of Dominant Symptom (Months)</b>		
Mean	27.7	36.5
Std Dev	37.64	43.73
Median	12.0	24.0
Minimum	0.03	0.07
Maximum	204	252
<b>Number of Fibroids Present</b>		
None	2 (1.3%)*	1 (1.7%)**
One	9 (6.0%)	5 (8.3%)
Two	10 (6.7%)	4 (6.7%)
Three	10 (6.7%)	8 (13.3%)
Four	10 (6.7%)	7 (11.7%)
Five	6 (4.0%)	2 (3.3%)
Six-Ten	27 (18.1%)	14 (23.3%)
> Ten	75 (50.3%)	13 (21.7%)
Missing	0 (0.0%)	6 (10.0%)
<b>Location of First Dominant Fibroid</b>		
Intramural	88 (59.1%)	26 (43.3%)

	UFE cohort (N=149)	Myomectomy Cohort (N=60)
Submucosal	1 (0.7%)	3 (5.0%)
Submucosal Pedunculated	17 (11.4%)	2 (3.3%)
Subserosal	8 (5.4%)	8 (13.3%)
Subserosal Pedunculated	31 (20.8%)	13 (21.7%)
Other	0 (0.0%)	1 (1.7%)
Cannot Determine	2 (1.3%)	0 (0.0%)
Missing	2 (1.3%)	7 (11.7%)
<b>Size of Dominant Fibroid (cm<sup>3</sup>)</b>		
N	147	53
Mean	182.12	226.92
Std Dev	208.978	196.394
Median	114.72	185.75
Minimum	0.63	6.45
Maximum	1116.91	869.95

Note: \*One subject had an ultrasound at baseline that showed evidence of the presence of fibroids, prior to her enrollment and UFE procedure. Another subject had a MRI at baseline that showed evidence of the presence of fibroids, prior to her enrollment and UFE procedure. Subsequent to enrollment and treatment, the Core Lab was unable to measure uterine fibroids for two subjects after reviewing each subject's baseline MRI. The Core Lab noted no measurable fibroids and the presence of adenomyosis. \*\*The Core Lab was unable to measure uterine fibroids for one myomectomy subject because the subject's baseline MRI was considered corrupt and could not be reviewed.

#### Study Results

One hundred forty nine (149) UFE subjects and 60 myomectomy subjects were considered intent-to-treat (ITT) and included in the FDA submission.

#### Procedure, Discharge, and Recovery Information

In the majority (96.6%) of the UFE procedures, the investigators used a right common femoral artery approach, while 3.4% used a bilateral common femoral artery approach. The majority (63.1%) of the Contour™ Embolization Particles used during the study was 355-500 microns, while 24.1% was 500-710 microns, and 12.7% was 250-355 microns. Many of the UFE procedures used more than one size of Contour Embolization Particles.

A majority of the UFE subjects (96.0%) were under conscious sedation during the procedure, 38.9% were under local anesthesia, 11.4% received epidural anesthesia, 0.7% received general anesthesia, and 6.0% used another form of anesthesia. The majority (98.3%) of the myomectomy subjects were under general anesthesia during the procedure, 1.7% were under local anesthesia, and 3.3% received epidural anesthesia. In some of the UFE and myomectomy procedures, more than one form of anesthesia was administered.

The mean total procedure time for the UFE subjects was 83.9 minutes (SD 47.44, range 13 to 240 minutes) which was less than the mean total procedure time in the myomectomy cohort of 141.6 minutes (SD 69.64, range 45 to 402 minutes).

The mean hospital stay for the UFE subjects was 23.8 hours (SD 11.02, range 5 to 60 hours) which was less than the mean hospital stay for the myomectomy cohort of 61.6 hours (SD 32.89, range 6 to 196 hours).

#### Primary Effectiveness Endpoint

Based on the 149 UFE ITT subjects, 130/149 (87.2%) UFE subjects completed up to 6-months of follow-up. Of the 130 subjects, 121/149 (81.2%, 95% CI (74.00%, 87.13%)) achieved success based on a  $\geq 5$  point improvement on the Overall Fibroid Symptom component of the UFQoL. There were 28 (18.8%) subjects not considered a success at 6-months; 8/149 (5.4%) did not achieve a  $\geq 5$  point increase on the Overall Fibroid Symptom component of a validated uterine fibroid quality of life questionnaire, 2/149 (1.3%) were never treated with embolic on the day of the UFE procedure and were considered technical failures, 10/149 (6.7%) were lost to follow-up, 4/149 (2.7%) withdrew consent, 1/149 (0.7%) had a post procedure myomectomy, 2/149 (1.3%) had a post procedure hysterectomy, and 1/149 (0.7%) were excluded from the analysis due to a significant inclusion criteria protocol deviation that may have confounded her data. Table # 2 presents the improvement in Overall Fibroid Symptom Score Success.

**Table 2. Primary Effectiveness Success Endpoint: Improvement in Overall Fibroid Symptom Score Success Rates**

Change in Overall Fibroid Symptom Score	6 Month Success Rates n/149 (%)
<b>Considered Success</b>	
$\geq 5$ point increase	121 (81.2%)
<b>Not Considered Success</b>	
$< 5$ point increase	8 (5.4%)
Technical Procedure Failure (never treated with embolic)	2 (1.3%)
Lost to Follow-up	10 (6.7%)
Withdrew Consent	4 (2.7%)
Post-procedure Myomectomy	1 (0.7%)
Post-procedure Hysterectomy	2 (1.3%)
Significant protocol deviation that may have confounded outcomes	1 (0.7%)

#### SECONDARY ENDPOINTS INCLUDING SAFETY SUMMARY

81.2% (121/149) UFE subjects through 6-months experienced an increase  $\geq 5$  points from baseline up to 6-months post procedure on the Overall Fibroid Symptom component of the UFQoL, while 75.0% (45/60) myomectomy subjects through 6-months experienced an increase  $\geq 5$  points from baseline up to 6-months post procedure on the Overall Fibroid Symptom component of the UFQoL. ( $p=0.315$ ).

Overall, the UFE cohort returned to normal activities after their procedures in a significantly shorter amount of time when compared to myomectomy cohort. The median number of days for the UFE cohort to return to normal activities was 9.5 days; however, the median number of days for the myomectomy cohort to return to normal activities was 39.5 days. However, when restricting this comparison to ages 36-40, to

remove the confounding with age, there was no significant difference between the UFE (median = 17.5) and myomectomy (median 26.5) cohorts.

Due to the significant non-Normality of dominant fibroid sizes at each of baseline, 3 months, and 6 months ( $p < 0.05$ ), the median change in fibroid size from baseline to 3 or 6 months was assessed using Wilcoxon Signed – Rank Test. The median size of the dominant fibroid in UFE subjects significantly reduced from baseline to 3 and 6-months post-procedure ( $p < 0.05$ ). The median fibroid size of the dominant fibroid was 116.9 cm<sup>3</sup> (at baseline) compared to 61.5 cm<sup>3</sup> (at 3-months) and 49.1 cm<sup>3</sup> (at 6-months). This represents a 39.2% decrease in fibroid size from baseline to 3-months and a 52.4% decrease in fibroid size from baseline to 6-months.

The total uterine volume for the UFE patients significantly reduced from baseline to 3 and 6-months post-procedure with  $p$ -values  $< 0.001$ . In the UFE cohort, there was a mean 30.0% reduction in uterine volume from baseline to 3-months and a mean 38.6% reduction from baseline to 6-months.

The menorrhagia bleeding scores for both the UFE and myomectomy patients significantly reduced from baseline to 3 and 6 months post-procedure with  $p$ -values  $< 0.001$ . In the UFE cohort, there was a 49.2% reduction in the menorrhagia bleeding score at 3-months and a 55.2% reduction in the menorrhagia bleeding score from baseline up to 6-months.

The UFE cohort missed a significantly shorter amount of days of work when compared to myomectomy cohort. The mean number of work days for the UFE cohort missed was 9.9 days; however, the mean number of missed work days for the myomectomy cohort was 37.0 days ( $p < 0.001$ ).

#### ADVERSE EVENTS

There were no deaths reported during the course of the study.

There were no reports of unanticipated adverse device effects occurring during the study.

There was a significantly higher adverse event rate in the myomectomy cohort when compared to the UFE cohort for number of subjects with at least one adverse event ( $p < 0.05$ ). Thirty seven out of 149 UFE subjects (24.8%) experienced at least one adverse event, whereas 24/60 (40.0%) myomectomy subjects experienced at least one adverse event.

Twenty four out of the 53 adverse events (45.3%) of the in the UFE cohort were related to the procedure, while 22/43 (51.2%) of the adverse events in the myomectomy cohort were related to the procedure. There were no reports of device-related adverse events in the UFE cohort occurring during the study.

There were 6/149 (4.0%) subjects with major adverse events for the UFE cohort and 1/60 (1.7%) for the myomectomy cohort ( $p > 0.05$ ). Of the 6 major adverse events in the UFE cohort, 3/6 were classified as related to the procedure, 1/6 was classified as not related to the procedure, and 2/6 were classified as prior condition/new onset. The one major adverse event in the myomectomy cohort was classified as procedure-related.

There were no UFE subjects enrolled with normal FSH values and regular menstrual cycles who post-procedure experienced both signs and symptoms of the ovarian failure (e.g., elevated FHS [sign], amenorrhea or hot flashes [symptoms]).

**Table 3. Timing and Type of Adverse Events Classified as UFE Procedure-Related**

Event Description	# Events
Hematoma	2
Post Embolization Syndrome	4
Prolapsed Fibroid	3
Urinary Tract Infection	4
Vaginal Discharge	2
Other*	9

\*Other includes events with only one occurrence; amenorrhea, abdominal cramping, arteriospasm left uterine artery, hot flashes, intractable pain, irregular period, ovarian failure, stinging pain on left side of pelvic area during intercourse, superficial venous thrombosis – L) calf pain.

**LONG-TERM FOLLOW-UP****Long-term Follow-up Subject Accountability**

Long-term follow-up information for the UFE cohort was collected at 1-year, 2-year and 3-year time intervals. A total of 120/149 (80.5%) UFE subjects returned for the 1-year follow-up visit, 102/149 (68.5%) returned for 2-year follow-up visit, and 90/149 (60.4%) UFE subjects returned for the 3-year follow-up visit. Table # 4 summarizes subject accountability throughout the 3-year time period.

**Table 4. Subject accountability**

	Contour™ Embolization Particles N
<b>Intent To Treat (ITT) Population</b>	149
6-Month Long-term Follow-up	130
Patient did not return for visit	-1
Lost to Follow-up	-9
12-Month Long-term Follow-up	120
Patient did not return for visit	-5
Patients that classified as LTF earlier - returned for 24-Month	2
Patients that voluntarily withdrew	-16
Missed 12-Month Follow-up - returned for 24-month	1
24-Month Long-term Follow-up	102
Lost to Follow-up	-17
Subjects who missed 24-month & came back for 36-month	5
36-Month Long-term Follow-up	90

**Long-term Effectiveness**

As the primary endpoint of the study was at 6 months, the long-term follow-up (1, 2, and 3 years) data was collected only on complications, menstrual bleeding questionnaire (Ruta), overall fibroid symptom component of the UFQOL questionnaire, new fibroid/recurrence symptoms, additional treatments and pregnancy information.

Change from baseline in all of the specific UFQOL measures (except hot flashes) were found to be statistically significant ( $p < 0.05$ ) at 1 year, 2 year, and 3 year follow-up in the UFE cohort of patients. This indicates the continuing trend in significant improvement in UFQOL measures over baseline through 3 year follow-up. Table # 5 summarizes the statistics on the Post-Market Change from Baseline Fibroid Specific Quality of Life Assessments for the UFE cohort.

**Table 5. Post-market Change From Baseline Fibroid Specific Quality of Life Assessments**

Overall Symptoms*	N	Mean	Std Dev	Min	Max
3 Months	135	27.9**	17.6	-24.2	81.5
6 Months	125	31.3**	18.1	-35.4	84.8
1 Year	115	35.6**	17.1	2.3	79.0
2 Years	100	36.4**	19.5	-7.2	84.8
3 Years	87	35.0**	18.7	-5.7	79.9

Note: \*Higher scores = Higher functioning/Less interference/Less symptoms/Less distress \*\*Statistically significant with  $p$ -value < 0.05.

Change from baseline in menorrhagia bleeding score was found to be statistically significant ( $p < 0.05$ ) at 1 year, 2 year, and 3 year follow-up in the UFE cohort of patients, indicating a continued trend of significant improvement over baseline through 3 year follow-up. Table # 6 summarizes statistics on the Post-Market Menorrhagia Questionnaire Bleeding Score Assessments for the UFE cohort.

**Table 6. Post-Market Menorrhagia Questionnaire Bleeding Score**

Change From Baseline Menorrhagia Score	N	Mean	% Change	Std Dev	Min	Max	p-value
3 Months	139	-24.5**	-49.3%	16.7	-63.7	20.3	< 0.001
6 Months	129	-26.6**	-55.3%	15.9	-65.8	11.5	< 0.001
1 Year	113	-28.6**	-61.1%	15.1	-65.0	3.3	< 0.001
2 Years	94	-28.6**	-61.0%	17.1	-85.7	9.0	< 0.001
3 Years	79	-28.9**	-60.9%	17.3	-75.2	11.9	< 0.001

Note: \*\*Statistically significant with  $p$ -value < 0.05.

### Symptom Recurrence, Complications, and Pregnancy Post-Market

Study subjects were evaluated for symptom recurrence, complication and pregnancy at 1, 2, and 3-years post-procedure. Table # 7 summarizes the results obtained from the clinician completed clinical assessment:

**Table 7. Symptom Recurrence Post-Market**

Recurrence/New Fibroid Symptoms	Contour™ Embolization Particles UFE arm (N=149)		
	Year 1 (n=120)	Year 2 (n=102)	Year 3 (n=88)
Yes	15 (12.5%)	14 (13.7%)	13 (14.8%)
No	105 (87.5%)	88 (86.3%)	75 (85.2%)
Missing	0 (0.0%)	0 (0.0%)	0 (0.0%)
Other Symptoms	Year 1	Year 2	Year 3
Yes	13 (10.8%)	16 (15.7%)	15 (17.0%)
No	107 (89.2%)	86 (84.3%)	73 (83.0%)
Missing	0 (0.0%)	0 (0.0%)	0 (0.0%)
New Treatments	Year 1	Year 2	Year 3
Yes	2 (1.7%)	6 (5.9%)	3 (3.4%)
No/Missing	118 (98.3%)	96 (94.1%)	85 (96.6%)
Specific Treatment	Year 1	Year 2	Year 3
Repeat UFE	0 (0.0%)	1 (1.0%)	0 (0.0%)
Myomectomy	0 (0.0%)	1 (1.0%)	1 (1.1%)
Hysteroscopy	0 (0.0%)	0 (0.0%)	0 (0.0%)
Hysterectomy	0 (0.0%)	2 (2.0%)	2 (2.3%)
Drug Treatment	2 (1.7%)	2 (2.0%)	0 (0.0%)
Other	0 (0.0%)	0 (0.0%)	0 (0.0%)

**Table 8. Complications and Pregnancy Post-Market**

Complications Related to UFE Procedure	Contour Embolization Particles UFE arm (N=149)		
	Year 1 (n=120)	Year 2 (n=102)	Year 3 (n=88)
Yes	*3 (2.5%)	0 (0.0%)	0 (0.0%)
No	117 (97.5%)	102 (100.0%)	88 (100.0%)
Missing	0 (0.0%)	0 (0.0%)	0 (0.0%)
Pregnant	Year 1	Year 2	Year 3
Yes	0 (0.0%)	0 (0.0%)	***1 (1.1%)
No	119 (99.2%)	102 (100.0%)	87 (98.9%)
Unknown	1 (0.8%)	0 (0.0%)	0 (0.0%)
Missing	0 (0.0%)	0 (0.0%)	0 (0.0%)
Specific Outcome	Year 1	Year 2	Year 3
Complicated	0 (0.0%)	0 (0.0%)	0 (0.0%)
Uncomplicated	0 (0.0%)	0 (0.0%)	1 (1.1%)
Therapeutic Abortion	0 (0.0%)	0 (0.0%)	0 (0.0%)
Spontaneous Abortion	0 (0.0%)	0 (0.0%)	0 (0.0%)

Note: \*03U06-Metallic taste in mouth last month. 09U04-"Something is there" in right upper leg below right groin area. 10U11-"Onset menopause symptoms". \*\*\*Subject 03U53-pregnancy reported at the 3-year follow-up visit.

### REFERENCES

#### UFE Indications

- Walker WJ and Pelage JP. Uterine Artery Embolization for Symptomatic Fibroids: Clinical Results in 400 Women with Imaging Follow-up. *BJOG* 2002; 109:1262-1272.
- Goodwin S, McLucas B, Lee M, et al. Uterine Artery Embolization for the Treatment of Uterine Leiomyomata Midterm Results. *JVIR* 1999; 10:1159-1165.
- Spies J, Spector A, Roth A, et al. Complications after Uterine Embolization for Leiomyomas. *Obstetrics & Gynecology* 2002; 100(5): 873-880.
- Pron G, Cohen M, Soucie J, Garvin G, Vanderburgh L, et al. The Ontario Uterine Fibroid Embolization Trial. Part 1. Baseline patient characteristics, fibroid burden, and impact of life. *Fertility and Sterility* 2003; 79(1): 112-119.
- Pron G, Bennett J, Common A, Wall J, et al. The Ontario Uterine Fibroid Embolization Trial. Part 2. Uterine fibroid reduction and symptom relief after uterine artery embolization for fibroids. *Fertility and Sterility* 2003; 79(1): 120-127.
- Siskin GP, Englander M, Stainken, et al. Embolic Agents Used for Uterine Artery Embolization. *AJR* 2000; 175: 767-773.
- Spies JB, Scialli AR, Jha LC, et al. Initial Results for Uterine Fibroid Embolization for Symptomatic Leiomyomata. *JVIR* 1999; 10: 1149-1157.

#### Other Indications

1. Solomon B, Soulen MC, Baum RA, et al. Hepatocellular carcinoma treated with CAM/Ethiodol/PVA™ chemoembolization: evaluation of survival and biologic and morphologic response in a U.S. population. *JVIR* 1999;10:793-798.
2. Tadavarthy SM, Knight L, Oviatt TW, et al. Therapeutic Transcatheter Arterial Embolization. *Radiology* 1974; 111: 13-16.
3. Brown KT, Koh BY, Brody LA, et al. Particle Embolization of Hepatic Neuroendocrine Metastases for Control of Pain and Hormonal Symptoms. *JVIR* 1999; 10: 397-403.
4. Sun S and Lang EV. Bone Metastases for Renal Cell Carcinoma: Preoperative Embolization. *JVIR* 1998; 9: 263-269.
5. Aina R, Oliva VL, Therasse E, et al. Arterial Embolotherapy for Upper Gastrointestinal Hemorrhage: Outcome Assessment. *JVIR* 2001; 12: 195-200.
6. Endovascular Treatment of Acute and Subacute Hemorrhage in the Head and Neck. *Arch Otolaryngol Head Neck Surg* 2000; 126: 1255-1262.
7. Leppanen M, Seppanen S, et al. Microcatheter embolization of intractable idiopathic epistaxis. *Cardiovasc Intervent Radiol*. 1999 Nov-Dec; 22(6):499-503.
8. Yu-Tang Goh P, Lin M, et al. Embolization for hemoptysis: a six-year review. *Cardiovasc Intervent Radiol*. 2002; 25(1):17-25.
9. Madoff DC, Hicks ME, Abdalla EK, et al. Portal vein embolization with polyvinyl alcohol particles and coils in preparation for major liver resection for hepatobiliary malignancy: safety and effectiveness--study in 26 patients. *Radiology* 2003; 227(1):251-60.

#### WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

Lupron Depot is a trademark of Tap Pharmaceuticals, Inc.

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## REFERENCIAS

### Indicaciones de la UFE

- Walker WJ and Pelage JP. Uterine Artery Embolization for Symptomatic Fibroids: Clinical Results in 400 Women with Imaging Follow-up. *BJOG* 2002; 109:1262-1272.
- Gowdin S, McLucas B, Lee M, et al. Uterine Artery Embolization for the Treatment of Uterine Leiomyomata Midterm Results. *JVIR* 1999; 10:1159-1165.
- Spies J, Spector A, Roth A, et al. Complications after Uterine Embolization for Leiomyomas. *Obstetrics & Gynecology* 2002; 100(5): 873-880.
- Pron G, Cohen M, Soucie J, Garvin G, Vanderburgh L, et al. The Ontario Uterine Fibroid Embolization Trial. Part 1. Baseline patient characteristics, fibroid burden, and impact of life. *Fertility and Sterility* 2003; 79(1): 112-119.
- Pron G, Bennett J, Common A, Wall J, et al. The Ontario Uterine Fibroid Embolization Trial. Part 2. Uterine fibroid reduction and symptom relief after uterine artery embolization for fibroids. *Fertility and Sterility* 2003; 79(1): 120-127.
- Siskin GP, Englander M, Stainken, et al. Embolic Agents Used for Uterine Artery Embolization. *AJR* 2000; 175: 767-773.
- Spies JB, Scialli AR, Jha LC, et al. Initial Results for Uterine Fibroid Embolization for Symptomatic Leiomyomata. *JVIR* 1999; 10: 1149-1157.

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- Solomon B, Soulen MC, Baum RA, et al. Hepatocellular carcinoma treated with CAM/Ethiodol/PVA™ chemoembolization: evaluation of survival and biologic and morphologic response in a U.S. population. *JVIR* 1999;10:793-798.
- Tadavarthi SM, Knight L, Oviatt TW, et al. Therapeutic Transcatheter Arterial Embolization. *Radiology* 1974; 111: 13-16.
- Brown KT, Koh BY, Brody LA, et al. Particle Embolization of Hepatic Neuroendocrine Metastases for Control of Pain and Hormonal Symptoms. *JVIR* 1999; 10: 397-403.
- Sun S and Lang EV. Bone Metastases for Renal Cell Carcinoma: Preoperative Embolization. *JVIR* 1998; 9: 263-269.
- Aina R, Oliva VL, Therasse E, et al. Arterial Embolotherapy for Upper Gastrointestinal Hemorrhage: Outcome Assessment. *JVIR* 2001; 12: 195-200.
- Endovascular Treatment of Acute and Subacute Hemorrhage in the Head and Neck. *Arch Otolaryngol Head Neck Surg* 2000; 126: 1255-1262.
- Leppanen M, Seppanen S, et al. Microcatheter embolization of intractable idiopathic epistaxis. *Cardiovasc Intervent Radiol.* 1999 Nov-Dec; 22(6):499-503.
- Yu-Tang Goh P, Lin M, et al. Embolization for hemoptysis: a six-year review. *Cardiovasc Intervent Radiol.* 2002; 25(1):17-25.
- Maddoff DC, Hicks ME, Abdalla EK, et al. Portal vein embolization with polyvinyl alcohol particles and coils in preparation for major liver resection for hepatobiliary malignancy: safety and effectiveness—study in 26 patients. *Radiology* 2003; 227(1):251-60.

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- Walker WJ and Pelage JP. Uterine Artery Embolization for Symptomatic Fibroids: Clinical Results in 400 Women with Imaging Follow-up. *BJOG* 2002; 109:1262-1272.
- Goodwin S, McLucas B, Lee M, et al. Uterine Artery Embolization for the Treatment of Uterine Leiomyomata Midterm Results. *JVIR* 1999; 10:1159-1165.
- Spies J, Spector A, Roth A, et al. Complications after Uterine Embolization for Leiomyomas. *Obstetrics & Gynecology* 2002; 100(5): 873-880.
- Pron G, Cohen M, Soucie J, Garvin G, Vanderburgh L, et al. The Ontario Uterine Fibroid Embolization Trial. Part 1. Baseline patient characteristics, fibroid burden, and impact of life. *Fertility and Sterility* 2003; 79(1): 112-119.
- Pron G, Bennett J, Common A, Wall J, et al. The Ontario Uterine Fibroid Embolization Trial. Part 2. Uterine fibroid reduction and symptom relief after uterine artery embolization for fibroids. *Fertility and Sterility* 2003; 79(1): 120-127.
- Siskin GP, Englander M, Stainken, et al. Embolic Agents Used for Uterine Artery Embolization. *AJR* 2000; 175: 767 773.
- Spies JB, Scialli AR, Jha LC, et al. Initial Results for Uterine Fibroid Embolization for Symptomatic Leiomyomata. *JVIR* 1999; 10: 1149-1157.

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- Solomon B, Soulen MC, Baum RA, et al. Hepatocellular carcinoma treated with CAM/Ethiodol/PVA™ thermoembolization: evaluation of survival and biologic and morphologic response in a U.S. population. *JVIR* 1999;10:793-798.
- Tadavarthi SM, Knight L, Oviatt TW, et al. Therapeutic Transcatheter Arterial Embolization. *Radiology* 1974; 111: 13-16.
- Brown KT, Koh BY, Brody LA, et al. Particle Embolization of Hepatic Neuroendocrine Metastases for Control of Pain and Hormonal Symptoms. *JVIR* 1999; 10: 397-403.
- Sun S and Lang EV. Bone Metastases for Renal Cell Carcinoma: Preoperative Embolization. *JVIR* 1998; 9: 263-269.
- Aina R, Oliva VL, Therasse E, et al. Arterial Embolotherapy for Upper Gastrointestinal Hemorrhage: Outcome Assessment. *JVIR* 2001; 12: 195-200.
- Endovascular Treatment of Acute and Subacute Hemorrhage in the Head and Neck. *Arch Otolaryngol Head Neck Surg* 2000; 126: 1255-1262.
- Leppanen M, Seppanen S, et al. Microcatheter embolization of intractable idiopathic epistaxis. *Cardiovasc Intervent Radiol.* 1999 Nov-Dec; 22(6):499-503.
- Yu-Tang Goh P, Lin M, et al. Embolization for hemoptysis: a six-year review. *Cardiovasc Intervent Radiol.* 2002; 25(1):17-25.
- Maddoff DC, Hicks ME, Abdalla EK, et al. Portal vein embolization with polyvinyl alcohol particles and coils in preparation for major liver resection for hepatobiliary malignancy: safety and effectiveness--study in 26 patients. *Radiology* 2003; 227(1):251-60.

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**LITERATURHINWEISE**

**UFE-Indikationen**

- Walker WJ und Pelage JP. Uterine Artery Embolization for Symptomatic Fibroids: Clinical Results in 400 Women with Imaging Follow-up. BJOG 2002; 109:1262-1272.
- Goodwin S, McLucas B, Lee M, et al. Uterine Artery Embolization for the Treatment of Uterine Leiomyomata Midterm Results. JIVR 1999; 10:1159-1165.
- Spies J, Spector A, Roth A, et al. Complications after Uterine Embolization for Leiomyomas. Obstetrics & Gynecology 2002; 100(5): 873-880.
- Pron G, Cohen M, Soucie J, Garvin G, Vanderburgh L, et al. The Ontario Uterine Fibroid Embolization Trial. Part 1. Baseline patient characteristics, fibroid burden, and impact of life. Fertility and Sterility 2003; 79(1): 112-119.

5. Pron G, Bennett J, Common A, Wall J, et al. The Ontario Uterine Fibroid Embolization Trial. Part 2. Uterine fibroid reduction and symptom relief after uterine artery embolization for fibroids. *Fertility and Sterility* 2003; 79(1): 120-127.
6. Siskin GP, Englander M, Stainken, et al. Embolic Agents Used for Uterine Artery Embolization. *AJR* 2000; 175: 767-773.
7. Spies JB, Scialli AR, Jha LC, et al. Initial Results for Uterine Fibroid Embolization for Symptomatic Leiomyomata. *JVIR* 1999; 10: 1149-1157.

#### Andere Indikationen

1. Solomon B, Soulen MC, Baum RA, et al. Hepatocellular carcinoma treated with CAM/Ethiodiol/PVA™ chemoembolization: evaluation of survival and biologic and morphologic response in a U.S. population. *JVIR* 1999; 10:793-798.
2. Tadavarthy SM, Knight L, Oviatt TW, et al. Therapeutic Transcatheter Arterial Embolization. *Radiology* 1974; 111: 13-16.
3. Brown KT, Koh BY, Brody LA, et al. Particle Embolization of Hepatic Neuroendocrine Metastases for Control of Pain and Hormonal Symptoms. *JVIR* 1999; 10: 397-403.
4. Sun S and Lang EV. Bone Metastases for Renal Cell Carcinoma: Preoperative Embolization. *JVIR* 1998; 9: 263-269.
5. Aina R, Oliva VL, Therasse E, et al. Arterial Embolotherapy for Upper Gastrointestinal Hemorrhage: Outcome Assessment. *JVIR* 2001; 12: 195-200.
6. Endovascular Treatment of Acute and Subacute Hemorrhage in the Head and Neck. *Arch Otolaryngol Head Neck Surg* 2000; 126: 1255-1262.
7. Leppanen M, Seppanen S, et al. Microcatheter embolization of intractable idiopathic epistaxis. *Cardiovasc Intervent Radiol.* 1999 Nov-Dec; 22(6):499-503.
8. Yu-Tang Goh P, Lin M, et al. Embolization for hemoptysis: a six-year review. *Cardiovasc Intervent Radiol.* 2002; 25(1):17-25.
9. Madoff DC, Hicks ME, Abdalla EK, et al. Portal vein embolization with polyvinyl alcohol particles and coils in preparation for major liver resection for hepatobiliary malignancy: safety and effectiveness--study in 26 patients. *Radiology* 2003; 227(1):251-60.

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## BIBLIOGRAFIA

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- Walker WJ and Pelage JP. Uterine Artery Embolization for Symptomatic Fibroids: Clinical Results in 400 Women with Imaging Follow-up. BJOG 2002; 109:1262-1272.
- Goodwin S, McLucas B, Lee M, et al. Uterine Artery Embolization for the Treatment of Uterine Leiomyomata Midterm Results. JVIR 1999; 10:1159-1165.
- Spies J, Spector A, Roth A, et al. Complications after Uterine Embolization for Leiomyomas. Obstetrics & Gynecology 2002; 100(5): 873-880.
- Pron G, Cohen M, Soucie J, Garvin G, Vanderburgh L, et al. The Ontario Uterine Fibroid Embolization Trial. Part 1. Baseline patient characteristics, fibroid burden, and impact of life. Fertility and Sterility 2003; 79(1): 112-119.
- Pron G, Bennett J, Common A, Wall J, et al. The Ontario Uterine Fibroid Embolization Trial. Part 2. Uterine fibroid reduction and symptom relief after uterine artery embolization for fibroids. Fertility and Sterility 2003; 79(1): 120-127.
- Siskin GP, Englander M, Stainken, et al. Embolic Agents Used for Uterine Artery Embolization. AJR 2000; 175: 767-773.
- Spies JB, Scialli AR, Jha LC, et al. Initial Results for Uterine Fibroid Embolization for Symptomatic Leiomyomata. JVIR 1999; 10: 1149-1157.

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- Solomon B, Soulen MC, Baum RA, et al. Hepatocellular carcinoma treated with CAM/Ethiodol/PVA™ chemoembolization: evaluation of survival and biologic and morphologic response in a U.S. population. JVIR 1999;10:793-798.
- Tadavarthy SM, Knight L, Oviatt TW, et al. Therapeutic Transcatheter Arterial Embolization. Radiology 1974; 111: 13-16.
- Brown KT, Koh BY, Brody LA, et al. Particle Embolization of Hepatic Neuroendocrine Metastases for Control of Pain and Hormonal Symptoms. JVIR 1999; 10: 397-403.
- Sun S and Lang EV. Bone Metastases for Renal Cell Carcinoma: Preoperative Embolization. JVIR 1998; 9: 263-269.
- Aina R, Oliva VL, Therasse E, et al. Arterial Embolotherapy for Upper Gastrointestinal Hemorrhage: Outcome Assessment. JVIR 2001; 12: 195-200.
- Endovascular Treatment of Acute and Subacute Hemorrhage in the Head and Neck. Arch Otolaryngol Head Neck Surg 2000; 126: 1255-1262.
- Leppanen M, Seppanen S, et al. Microcatheter embolization of intractable idiopathic epistaxis. Cardiovasc Intervent Radiol. 1999 Nov-Dec; 22(6):499-503.
- Yu-Tang Goh P, Lin M, et al. Embolization for hemoptysis: a six-year review. Cardiovasc Intervent Radiol. 2002; 25(1):17-25.
- Maddoff DC, Hicks ME, Abdalla EK, et al. Portal vein embolization with polyvinyl alcohol particles and coils in preparation for major liver resection for hepatobiliary malignancy: safety and effectiveness--study in 26 patients. Radiology 2003; 227(1):251-60.

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## LITERATUUR

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- Walker WJ and Pelage JP. Uterine Artery Embolization for Symptomatic Fibroids: Clinical Results in 400 Women with Imaging Follow-up. BJOG 2002; 109:1262-1272.
- Gowdin S, McLucas B, Lee M, et al. Uterine Artery Embolization for the Treatment of Uterine Leiomyomata Midterm Results. JVIR 1999; 10:1159-1165.
- Spies J, Spector A, Roth A, et al. Complications after Uterine Embolization for Leiomyomas. Obstetrics & Gynecology 2002; 100(5): 873-880.
- Pron G, Cohen M, Soucie J, Garvin G, Vanderburgh L, et al. The Ontario Uterine Fibroid Embolization Trial. Part 1. Baseline patient characteristics, fibroid burden, and impact of life. Fertility and Sterility 2003; 79(1): 112-119.
- Pron G, Bennett J, Common A, Wall J, et al. The Ontario Uterine Fibroid Embolization Trial. Part 2. Uterine fibroid reduction and symptom relief after uterine artery embolization for fibroids. Fertility and Sterility 2003; 79(1): 120-127.
- Siskin GP, Englander M, Stainken, et al. Embolic Agents Used for Uterine Artery Embolization. AJR 2000; 175: 767-773.
- Spies JB, Scialli AR, Jha LC, et al. Initial Results for Uterine Fibroid Embolization for Symptomatic Leiomyomata. JVIR 1999; 10: 1149-1157.

### Overige indicaties

- Solomon B, Soulen MC, Baum RA, et al. Hepatocellular carcinoma treated with CAM/Ethiodol/PVA™ chemoembolization: evaluation of survival and biologic and morphologic response in a U.S. population. JVIR 1999;10:793-798.
- Tadavarthi SM, Knight L, Oviatt TW, et al. Therapeutic Transcatheter Arterial Embolization. Radiology 1974; 111: 13-16.
- Brown KT, Koh BY, Brody LA, et al. Particle Embolization of Hepatic Neuroendocrine Metastases for Control of Pain and Hormonal Symptoms. JVIR 1999; 10: 397-403.
- Sun S and Lang EV. Bone Metastases for Renal Cell Carcinoma: Preoperative Embolization. JVIR 1998; 9: 263-269.
- Aina R, Oliva VL, Therasse E, et al. Arterial Embolotherapy for Upper Gastrointestinal Hemorrhage: Outcome Assessment. JVIR 2001; 12: 195-200.
- Endovascular Treatment of Acute and Subacute Hemorrhage in the Head and Neck. Arch Otolaryngol Head Neck Surg 2000; 126: 1255-1262.
- Leppanen M, Seppanen S, et al. Microcatheter embolization of intractable idiopathic epistaxis. Cardiovasc Intervent Radiol. 1999 Nov-Dec; 22(6):499-503.
- Yu-Tang Goh P, Lin M, et al. Embolization for hemoptysis: a six-year review. Cardiovasc Intervent Radiol. 2002; 25(1):17-25.
- Maddoff DC, Hicks ME, Abdalla EK, et al. Portal vein embolization with polyvinyl alcohol particles and coils in preparation for major liver resection for hepatobiliary malignancy: safety and effectiveness--study in 26 patients. Radiology 2003; 227(1):251-60.

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**REFERÊNCIAS**

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- Walker WJ and Pelage JP. Uterine Artery Embolization for Symptomatic Fibroids: Clinical Results in 400 Women with Imaging Follow-up. BJOG 2002; 109:1262-1272.
- Goodwin S, McLucas B, Lee M, et al. Uterine Artery Embolization for the Treatment of Uterine Leiomyomata Midterm Results. JIVIR 1999; 10:1159-1165.
- Spies J, Spector A, Roth A, et al. Complications after Uterine Embolization for Leiomyomas. Obstetrics & Gynecology 2002; 100(5): 873-880.
- Pron G, Cohen M, Soucie J, Garvin G, Vanderburgh L, et al. The Ontario Uterine Fibroid Embolization Trial. Part 1. Baseline patient characteristics, fibroid burden, and impact of life. Fertility and Sterility 2003; 79(1): 112-119.
- Pron G, Bennett J, Common A, Wall J, et al. The Ontario Uterine Fibroid Embolization Trial. Part 2. Uterine fibroid reduction and symptom relief after uterine artery embolization for fibroids. Fertility and Sterility 2003; 79(1): 120-127.
- Siskin GP, Englander M, Stainken, et al. Embolic Agents Used for Uterine Artery Embolization. AJR 2000; 175: 767-773.
- Spies JB, Scialli AR, Jha LC, et al. Initial Results for Uterine Fibroid Embolization for Symptomatic Leiomyomata. JIVIR 1999; 10: 1149-1157.

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- Solomon B, Soulen MC, Baum RA, et al. Hepatocellular carcinoma treated with CAM/ethiodiol/PVA™ chemoembolization: evaluation of survival and biologic and morphologic response in a U.S. population. JIVIR 1999;10:793-798.
- Tadavarthy SM, Knight L, Ovitt TW, et al. Therapeutic Transcatheter Arterial Embolization. Radiology 1974; 111: 13-16.
- Brown KT, Koh BY, Brody LA, et al. Particle Embolization of Hepatic Neuroendocrine Metastases for Control of Pain and Hormonal Symptoms. JIVIR 1999; 10: 397-403.
- Sun S and Lang EV. Bone Metastases for Renal Cell Carcinoma: Preoperative Embolization. JIVIR 1998; 9: 263-269.
- Aina R, Oliva VL, Therasse E, et al. Arterial Embolotherapy for Upper Gastrointestinal Hemorrhage: Outcome Assessment. JIVIR 2001; 12: 195-200.
- Endovascular Treatment of Acute and Subacute Hemorrhage in the Head and Neck. Arch Otolaryngol Head Neck Surg 2000; 126: 1255-1262.

7. Leppanen M, Seppanen S, et al. Microcatheter embolization of intractable idiopathic epistaxis. Cardiovasc Intervent Radiol. 1999 Nov-Dec; 22(6):499-503.
8. Yu-Tang Goh P, Lin M, et al. Embolization for hemoptysis: a six-year review. Cardiovasc Intervent Radiol. 2002; 25(1):17-25.
9. Madoff DC, Hicks ME, Abdalla EK, et al. Portal vein embolization with polyvinyl alcohol particles and coils in preparation for major liver resection for hepatobiliary malignancy: safety and effectiveness--study in 26 patients. Radiology 2003; 227(1):251-60.

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This page left intentionally blank.

This page left intentionally blank.

This page left intentionally blank.

This page left intentionally blank.

This page left intentionally blank.

This page left intentionally blank.

This page left intentionally blank.

This page left intentionally blank.

This page left intentionally blank.

This page left intentionally blank.

This page left intentionally blank.

This page left intentionally blank.

This page left intentionally blank.

This page left intentionally blank.

This page left intentionally blank.

This page left intentionally blank.

This page left intentionally blank.

This page left intentionally blank.

This page left intentionally blank.

This page left intentionally blank.

This page left intentionally blank.

This page left intentionally blank.

This page left intentionally blank.

This page left intentionally blank.

This page left intentionally blank.

This page left intentionally blank.

This page left intentionally blank.

This page left intentionally blank.

This page left intentionally blank.

This page left intentionally blank.

This page left intentionally blank.

This page left intentionally blank.

This page left intentionally blank.

This page left intentionally blank.

This page left intentionally blank.

This page left intentionally blank.

This page left intentionally blank.

This page left intentionally blank.

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**LOT**

Lot  
Lote  
Lot  
Charge  
Lotto  
Partij  
Lote

**UPN**

Product Number  
Número del producto  
Référence  
Produktnummer  
Codice prodotto  
Productnummer  
Número do Produto

Recyclable Package  
Envase reciclable  
Emballage recyclable  
Wiederverwertbare Verpackung  
Confezione riciclabile  
Recyclebare verpakking  
Embalagem Reciclável

Use By  
Fecha de caducidad  
Date limite d'utilisation  
Verwendbar bis  
Usare entro  
Uiterste gebruiksdatum  
Validade

**AUS**

Australian Sponsor Address  
Dirección del patrocinador australiano  
Adresse du promoteur australien  
Adresse des australischen Sponsors  
Indirizzo sponsor australiano  
Adres Australische sponsor  
Endereço do Patrocinador Australiano

Do Not Resterilize  
No reesterilizar  
Ne pas restériliser  
Nicht erneut sterilisieren  
Non risterilizzare  
Niet opnieuw steriliseren  
Não reesterilize

For single use only. Do not reuse.  
Para un solo uso. No reutilizar.  
À usage unique. Ne pas réutiliser.  
Für den einmaligen Gebrauch. Nicht wieder verwenden.  
Esclusivamente monouso. Non riutilizzare.  
Uitsluitend bestemd voor eenmalig gebruik. Niet opnieuw gebruiken.  
Apenas para uma única utilização. Não reutilize.

Do not use if package is damaged.  
No usar si el envase está dañado.  
Ne pas utiliser si l'emballage est endommagé.

Bei beschädigter Verpackung nicht verwenden.  
Non usare il prodotto se la confezione è danneggiata.  
Niet gebruiken als de verpakking is beschadigd.  
Não utilize se a embalagem estiver danificada.

**STERILE R**

Sterilized using irradiation.  
Esterilizado por radiación.  
Stérilisé par irradiation.  
Durch Bestrahlung sterilisiert.  
Sterilizzato mediante radiazioni.  
Gesteriliseerd met bestraling.  
Esterilizado por irradiação.

Open Here  
Abrir aquí  
Ouvrir ici  
Hier öffnen  
Aprire qui  
Hier openen  
Abra Aqui

Peripheral Use Only  
Uso periférico exclusivamente  
Usage périphérique uniquement  
Nur zur peripheren Verwendung  
Eclusivamente per uso periférico  
Uitsluitend voor gebruik in perifere vaten  
Apenas para uso periférico

**EC REP**

**EU Authorized Representative**

**Boston Scientific International S.A.**  
55 avenue des Champs Pierreux  
TSA 51101  
92729 NANTERRE CEDEX  
FRANCE

**AUS**

**Australian Sponsor Address**

**Boston Scientific (Australia) Pty Ltd**  
PO Box 332  
BOTANY  
NSW 1455  
Australia  
Free Phone 1800 676 133  
Free Fax 1800 836 666

**Legal Manufacturer**

**Boston Scientific Corporation**  
One Boston Scientific Place  
Natick, MA 01760-1537  
USA  
USA Customer Service 888-272-1001

**Do not use if package is damaged.**

**Recyclable Package**

**CE 0197**

© 2010 Boston Scientific Corporation or its affiliates. All rights reserved.