


WATCHMAN FLX™
LEFT ATRIAL APPENDAGE CLOSURE DEVICE

**Boston
Scientific**
Advancing science for life™





—
THE LEADER IN LAAC THERAPY




—
BUILT ON THE MOST STUDIED AND IMPLANTED LAAC
DEVICE IN THE WORLD — WATCHMAN FLX IS DESIGNED
TO **ADVANCE PROCEDURAL PERFORMANCE** AND **SAFETY**
WHILE **EXPANDING THE TREATABLE PATIENT POPULATION.**

WATCHMAN FLX DEVICE

Full recapture, reposition and redeploy capabilities for precise placement



77% reduced metal exposure



80% more contact points for sealing

Dual-row precision anchors designed to provide optimal device engagement with LAA tissue for long-term stability



WATCHMAN FLX ball – **fully rounded** designed to safely advance and maneuver within the LAA



—
ADVANCE SAFETY

ADVANCE PROCEDURAL PERFORMANCE


EXPAND THE TREATABLE PATIENT POPULATION

ADVANCE SAFETY




The PINNACLE FLX clinical trial demonstrated the procedural safety and closure efficacy of the WATCHMAN FLX device.


Primary Safety Endpoint*



*All-cause death, ischemic stroke, systemic embolism, or device- or procedure-related adverse events requiring surgery or major endovascular intervention within 7 days following the procedure or by hospital discharge, whichever is later.



*Based on the combined rate observed in PREVAIL¹ and CAP2², plus a clinically acceptable delta.



¹ Holmes, DR., et al. (2014). J Am Coll Cardiol 64(1): 1-12


² Holmes, DR., et al. JACC 2019

ADVANCE PROCEDURAL PERFORMANCE

Primary Effectiveness Endpoint

100%

of Subjects
Demonstrated
Effective LAA
Closure at
12 Months*



*LAA closure at 12 months is defined as any peri-device flow with jet size ≤ 5 mm per core laboratory-assessed TEE

**Performance goal based on the rates observed in PREVAIL¹ and CAP2², minus a clinically relevant delta

Procedure Performance

Procedure/Implant Success

99% Implant Success

Implant success defined as successful delivery and release of a WATCHMAN FLX device into the LAA

NOAC Discontinuation

96.2% of Patients Discontinued NOAC at 45-day Follow-up

Study/OAC	% Discontinuation
PINNACLE FLX/NOAC	96.2%
PREVAIL/warfarin ¹	92%
CAP2/warfarin ²	93%

¹ Holmes, DR., et al. (2014). J Am Coll Cardiol 64(1): 1-12

² Holmes DR et al, JACC 2019


EXPAND THE TREATABLE PATIENT POPULATION

WATCHMAN FLX is designed to treat the widest range of patient anatomies, with five device sizes treating ostia from 14 mm to 31.5 mm.



**Devices not shown to scale*

Greater Device Sizing Overlap



WATCHMAN FLX SYSTEM


WATCHMAN FLX Device



Delivery Catheter



WATCHMAN FXD Curve Access System



WATCHMAN FLX Device

Nitinol frame with Polyethylene Terephthalate (PET) fabric cover

WATCHMAN FLX Delivery Catheter

SHEATH MATERIAL

Braided Pebax® with PTFE liner and platinum/iridium marker band

WATCHMAN FXD Curve Access System

HUB MATERIAL

Pebax® with polycarbonate cap

SHEATH MATERIAL

Braided Pebax® with PTFE liner and platinum/iridium marker band

DILATOR



High density polyethylene (HDPE)/low density polyethylene (LDPE) 50:50 blend

ORDERING INFORMATION

WATCHMAN FLX LAAC DEVICE ORDERING INFORMATION

Reference Catalog No.	Description	Size	Order Number (GTIN)	ID	OD	Barcode
M635WU50200	WATCHMAN FLX LAAC Device and Delivery Catheter	20 mm	08714729860488	-	12F (4.0 mm)	
M635WU50240	WATCHMAN FLX LAAC Device and Delivery Catheter	24 mm	08714729860495	-	12F (4.0 mm)	
M635WU50270	WATCHMAN FLX LAAC Device and Delivery Catheter	27 mm	08714729860501	-	12F (4.0 mm)	
M635WU50310	WATCHMAN FLX LAAC Device and Delivery Catheter	31 mm	08714729860518	-	12F (4.0 mm)	
M635WU50350	WATCHMAN FLX LAAC Device and Delivery Catheter	35 mm	08714729860471	-	12F (4.0 mm)	

WATCHMAN FXD CURVE ACCESS SYSTEM ORDERING INFORMATION

Reference Catalog No.	Description	Curve	Order Number (GTIN)	ID	OD	Barcode
M635TU80010	WATCHMAN FXD Access System SGL US	Single	00191506013806	12F (4.2 mm)	15F (5.0 mm)	
M635TU80020	WATCHMAN FXD Access System DBL US	Double	00191506013813	12F (4.2 mm)	15F (5.0 mm)	

Please contact your Boston Scientific sales representative for ordering information.

WATCHMAN FLX is preloaded into the delivery catheter thus reducing the preparation time.

BRIEF SUMMARY

CAUTION: *Federal law (USA) restricts this device to sale by or on the order of a physician. Rx only. Prior to use, please see the complete "Instructions for Use" for more information on Indications, Contraindications, Warnings, Precautions, Adverse Events, and Operator's Instructions.*

INTENDED USE/INDICATIONS FOR USE

The WATCHMAN FLX™ Device is indicated to reduce the risk of thromboembolism from the left atrial appendage in patients with non-valvular atrial fibrillation who:

- Are at increased risk for stroke and systemic embolism based on CHA₂DS₂-VASc scores and are recommended for anticoagulation therapy;
- Are deemed by their physicians to be suitable for anticoagulation therapy;
- Have an appropriate rationale to seek a non-pharmacologic alternative to anticoagulation therapy, taking into account the safety and effectiveness of the device compared to anticoagulation therapy.

CONTRAINDICATIONS

Do not use the WATCHMAN FLX Device if:

- Intracardiac thrombus is present.
- An atrial septal defect repair or closure device or a patent foramen ovale repair or closure device is present.
- The LAA anatomy will not accommodate a Closure Device (see Table 62 of the eIFU).
- The patient has a known hypersensitivity to any portion of the device material or the individual components (see Device Description section of the eIFU) such that the use of the WATCHMAN FLX Device is contraindicated.
- Any of the customary contraindications for other percutaneous catheterization procedure (e.g.,

patient size too small to accommodate TEE probe or required catheters) or conditions (e.g., active infection, bleeding disorder) are present.

- There are contraindications to the use of anti-coagulation therapy, aspirin, or P2Y₁₂ inhibitor.

WARNINGS

Implantation of the WATCHMAN FLX Device should only be performed by interventional cardiologists and/or electrophysiologists who are proficient in percutaneous procedures, transseptal procedures, the imaging modality utilized and who have completed the WATCHMAN FLX Physician Training program.

- This device has not been studied in pregnant or breastfeeding women. Careful consideration should be given to use of the Closure Device in pregnant and/or breastfeeding women due to the risk of significant exposure to x-rays and the use of anticoagulation medication.
- Device selection should be based on accurate LAA measurements obtained using transesophageal or intracardiac echocardiographic imaging guidance in multiple views to avoid improper Closure Device sizing. For TEE recommended in multiple angles [e.g., 0°, 45°, 90°, 135°]; For ICE imaging, visualization of the LAA is recommended with the following anatomical structures: aortic valve (short-axis), mitral valve (long-axis), and pulmonary artery (short-axis), to assess the minimum and maximum diameter of the LAA ostium.
- Do not release (i.e., unscrew) the WATCHMAN FLX Device from the core wire unless all release criteria are satisfied to avoid suboptimal results.
- Potential for Closure Device embolization exists with cardioversion < 30 days following Closure Device implantation; verify Closure Device position after cardioversion during this period.
- Appropriate post-procedure drug therapy should

be followed. See Post-Procedure Information section (of the eIFU) for further detail.

PRECAUTIONS

- The safety and effectiveness (and benefit-risk profile) of the WATCHMAN FLX Device has not been established in patients for whom long-term anticoagulation is determined to be contraindicated.
- The LAA is a thin-walled structure. Use caution when accessing the LAA, and deploying, recapturing, and repositioning the Closure Device.
- Use caution when introducing a WATCHMAN Access System to prevent damage to cardiac structures.
- Use caution when introducing the Delivery System to prevent damage to cardiac structures.
- To prevent damage to the Delivery Catheter or Closure Device, do not allow the WATCHMAN FLX Device to protrude beyond the distal tip of the Delivery Catheter when inserting the Delivery System into the Access Sheath.
- If using a power injector, the maximum pressure should not exceed 690 kPa (100 psi).

PATIENT SELECTION FOR TREATMENT

In considering the use of the WATCHMAN FLX Device, the rationale for seeking an alternative to long-term anticoagulation therapy and the safety and effectiveness of the device compared to anticoagulation should be taken into account.

- The presence of indication(s) for long-term anticoagulation therapy, other than non-valvular atrial fibrillation (e.g., mechanical heart valve, hypercoagulable states, recurrent deep venous thrombosis). Details regarding the indications, contraindications, warnings, and precautions for oral anticoagulants approved for patients with non-valvular atrial fibrillation are provided in their respective Instructions for Use.

Of note:

- The safety and effectiveness (and benefit-risk profile) of the WATCHMAN FLX Device has not been established in patients for whom long-term anti-coagulation is determined to be contraindicated. Factors that need to be considered for the WATCHMAN FLX Device and implantation procedure include the following:
- Overall medical status, including conditions which might preclude the safety of a percutaneous, transcatheter procedure.
- Suitability for percutaneous, transseptal procedures, including considerations of:
 - Cardiac anatomy relating to the LAA size and shape.
 - Vascular access anatomy (e.g., femoral vein size, thrombus, or tortuosity).
 - Ability of the patient to tolerate general or local anesthesia.
 - Ability of the patient to undergo required imaging.
- Ability to comply with the recommended post-WATCHMAN FLX Device implant pharmacologic regimen (see Post-Procedure Information section) especially for patients at high risk for bleeding.

ADVERSE EVENTS

Potential adverse events (in alphabetical order) which may be associated with the use of a left atrial appendage closure device or implantation procedure include but are not limited to: Air embolism, Airway trauma, Allergic reaction to the contrast media, anesthetic, WATCHMAN Implant material, or medication, Altered mental status, Anemia requiring transfusion, Anesthesia risks, Angina, Anoxic encephalopathy, Arrhythmias, Atrial septal defect, Bruising, hematoma, or seroma near the catheter insertion site, Cardiac perforation,

Chest pain/discomfort, Confusion post procedure, Congestive heart failure, Contrast related nephropathy, Cranial bleed, Death, Decreased hemoglobin, Deep vein thrombosis, Device embolism, Device fracture, Device thrombosis, Edema, Embolism, Excessive bleeding, Fever, Fistula, Groin pain, Groin puncture bleed, Hematuria, Hemoptysis, Hypotension, Hypoxia, Improper wound healing, Inability to reposition, recapture, or retrieve the device, Infection/pneumonia, Interatrial septum thrombus, Intratracheal bleeding, Major bleeding requiring transfusion, Misplacement of the device/improper seal of the appendage/movement of device from appendage wall, Myocardial erosion, Myocardial infarction, Nausea, Oral bleeding, Pericardial effusion/tamponade, Pleural effusion, Prolonged bleeding from a laceration, Pseudoaneurysm, Pulmonary edema, Radiation injury, Renal failure, Respiratory insufficiency/failure, Stroke – Hemorrhagic, Stroke – Ischemic, Surgical removal of the device, TEE complications (e.g., throat pain, bleeding, esophageal trauma), Thrombocytopenia, Thrombosis, Transient ischemic attack (TIA), Valvular or vascular damage, Vasovagal reactions.

There may be other potential adverse events that are unforeseen at this time.

92574167 B.5



Interventional Cardiology

300 Boston Scientific Way
Marlborough, MA 01752-1234

www.bostonscientific.com

Medical Professionals:

1.800.CARDIAC (227.3422)

Patients and Families:

1.866.484.3268

©2020 Boston Scientific Corporation
or its affiliates. All rights reserved.

SH-1152208-AA

