FAQ

Device Replacement
Cambio de Dispositivo


Boston
Scientific

Defining tomorrow, today.™

Cardiology, Rhythm and Vascular
One Boston Scientific Place
Natick, MA 01760-1537 USA
www.bostonscientific.com

Medical Professionals:
1.800.CARDIAC (227.3422)
Patients and Families:
1.866.484.3268

© 2011 by Boston Scientific Corporation
or its affiliates. All rights reserved.
CRM9-1722-0411

Why Choose a Boston Scientific Device?

Boston Scientific Cardiac Rhythm Management (CRM) has dedicated our business to provide you with products and therapies that help you live a full and active life.

All Boston Scientific devices have enhanced therapy capabilities and information to help your doctor manage your medical condition.

- Pacemakers provide advanced minute ventilation sensor technology designed to do a better job of adjusting your heart rate while going uphill or doing other activities by sensing the associated increase in breathing.
- ICD devices can treat fast, abnormal and lethal ventricular rhythms with proven therapies. Additional features control the pacing response to atrial fibrillation, reduce unnecessary ventricular pacing, and evaluate heart rhythms to better determine which ventricular rhythms should be treated with shock therapy.
- The three-lead, CRT devices use cardiac resynchronization therapy (CRT) to deliver pacing pulses to the left and right ventricles to help them contract at the same time. This will help your heart pump more efficiently.
- If you are concerned about body image, Boston Scientific's newest generation of ICD and CRT-D devices are the smallest, thinnest, high-energy devices in the world.
- A complete line of leads are available to meet different patient needs.
- The LATITUDE® Patient Management system allows certain implanted devices to send specific device information and check your heart health from the comfort of your home as scheduled by your doctor.
- A dedicated patient support line is available to answer questions you have about your implanted device. Call 1.866.484.3268 for information about your LATITUDE system, connect with Technical Services or talk to Boston Scientific Medical Records.

For More Information

Contact Boston Scientific Patient Services at 1.866.484.3268.

Boston Scientific wants to help you lead a full, active life. We know that you may have questions as you adjust to living with your device. While we always recommend that you talk with your doctor first about your specific therapy, we also encourage you to:

Read the patient handbook that came with your device.

For additional copies of the handbook, call Boston Scientific Patient Services toll free at 1.866.484.3268 and select the option for technical questions.

Visit www.lifebeatonline.com

Visit www.lifebeatonline.com for stories of other heart patients, explanations of common heart conditions, resources, EMI, travel, and much more.

Our website offers:

- Animations that show how the heart works
- Clear explanations of common heart and blood vessel conditions
- Stories about other heart patients
- Pictures and descriptions of implantable devices
- Concise, printable summaries specific to your Boston Scientific device
- Information on how to be safe around appliances and equipment
- Worldwide clinic locator that identifies where devices can be checked

Call Boston Scientific Patient Services toll free at 1.866.484.3268.

You or your family can call to speak personally with a Patient Services representative.

Frequently Asked Questions About Device Replacement

All implanted cardiac devices [pacemakers, implantable cardioverter defibrillators (ICDs), and cardiac resynchronization therapy devices (CRT-P/CRT-D)] have batteries inside that supply the energy the devices need to do their jobs. When the energy in your device battery begins to get low, your doctor will make an appointment to discuss replacing your device.

Below are typical questions patients ask as they prepare for device replacement.

Q. How will I know when my device battery is getting low and needs to be changed?

A. Your doctor will regularly check the battery of your implanted device to determine how long it will last.

These checks will be done on a schedule set by your doctor and may include a combination of checks done in the office and from your home.

If your device is an ICD or CRT-D, it can be programmed to beep at replacement time.

Q. What conditions cause my device's battery to lose energy?

A. Like any battery, the energy in your implanted device will naturally decrease over time. The battery's life is affected both by how much energy is programmed and how often it is required to pace or deliver shock therapy to your heart. This varies patient-by-patient depending on the individual's medical condition.


Q. When and where will my device be replaced?

A. Implanted device batteries typically last between 4 and 8 years. Your doctor will schedule the replacement procedure based on the battery status information of the device. The device replacement is a surgical procedure and therefore is done in a sterile environment and requires you to be in a hospital setting.

Q. Will my entire device be replaced or just the battery?

A. Pacemakers, ICDs and cardiac resynchronization therapy (CRT) devices are "all-in-one" units. During manufacturing, the battery and other parts of your device are sealed inside a metal case. Nothing can get in or out of the metal case.

The implantable device cannot be opened and "re-sealed." When the battery runs down, a new device will be implanted. In most cases, your original lead wires will not need to be replaced.

Inside an ICD: The battery and other parts of the device are sealed inside the metal case. These include digital hybrid circuitry, a single-cell battery, and capacitors (for energy storage).

Q. Will I get the same type of device? Will it be from the same manufacturer?

A. Your doctor will select the type of device you need based on your medical condition. You may get the same type of device that you previously had, but if your medical condition has changed, the type of device you need may change. For example, you may go from a pacemaker to an ICD or CRT-D device if your heart condition changes.

- If your doctor chooses the same manufacturer, you will continue to have the same technical support and company representatives as before.
- If your doctor chooses a device from a different company, device features and how they are used may change. In addition, the company representatives and support systems will change.

Q. What happens to the old device that is taken out?


A. When the device is removed from your body, it is commonly sent back to the manufacturing company. Some hospitals require that it be sent to the hospital pathology lab for identification before being sent back to the manufacturing company.

Q. Will I need new leads as well?

A. Devices and leads made today use a common design to connect to each other, regardless of who made them. It is very common to use previously implanted leads with new devices.

- If your new device is from the same manufacturer, it is unlikely that you will need any additional equipment as the leads and devices are made to fit together.
- Older leads that may be of a different size may be adapted to fit a new device.

Your doctor will determine how your existing leads will be used.

Anatomy of a three-lead implant (CRT). There is a lead in the right atrium, a lead in the right ventricle, and a lead on the outside wall of the left ventricle.

Questions to ask your doctor:

We encourage you to talk with your healthcare provider about your condition and treatment options. To help with that conversation, we've listed some questions you may wish to discuss in addition to your own.

Procedure

- How will I prepare for device replacement?
- Will I have any discomfort afterwards?
- What are the risks of the procedure?
- What is the replacement procedure like?
- Will I need to be given medication to make me sleep?
- Do I need to stay overnight?

Medications

- Do I need to alter my medication schedule?

Device

- Will I get the same type of device?
- Will the device be from the same manufacturer?
- Will I need any new leads or adapters?

Follow-up

- Will my new device be monitored from home?
- How soon can I return to my regular activities?

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Important Safety Information

Cardiac Resynchronization Therapy Devices

Cardiac resynchronization therapy pacemakers (CRT-P) and defibrillators (CRT-D) are designed to treat heart failure patients who may or may not have symptoms (CRT-D), or who may have symptoms despite the best available drug therapy (CRT-P and CRT-D). They are also designed to help your heart pump more effectively and meet your body's need for blood flow. CRT-P and CRT-D systems are not for everyone, including people with separate implantable cardioverter-defibrillators (CRT-P only) or certain steroid allergies. Procedure risks include infection, tissue damage, and kidney failure. Patients who should not receive this device include: patients who have additional medical conditions that may not allow the pacemaker to function appropriately (CRT-P only) and patients whose ventricular rhythm disturbances or heart failure have a reversible or temporary cause. In some cases, the device may be unable to detect or appropriately respond to your heart rhythm (CRT-P and CRT-D) or may deliver inappropriate shocks (CRT-D only). In rare cases severe complications or device failures can occur. Electrical or magnetic fields can affect the device. Only your doctor knows what is right for you. These devices are available by prescription only. Individual results may vary. Your physician should discuss all potential benefits and risks with you. For further safety information, refer to the Physicians Instructions for Use or call 1/866-484-3268.

Implantable Cardioverter Defibrillators

An implantable cardioverter defibrillator is designed to monitor and treat heart rhythm problems, greatly reducing the risks associated with them. But it is not for everyone, including people with certain steroid allergies. Procedure risks include infection, tissue damage, and kidney failure. Patients who should not receive this device include: patients whose ventricular rhythm disturbances have a reversible or temporary cause and patients with certain types of atrial rhythm disturbances. Procedure risks include infection and tissue damage. In some cases, the device may not respond to irregular heartbeats or may deliver inappropriate shocks. In rare cases severe complications or device failures can occur. Electrical or magnetic fields can affect the device. Only your doctor knows what is right for you. This device is available by prescription only. Individual results may vary. Your physician should discuss all potential benefits and risks with you. For further safety information, refer to the Physicians Instructions for Use or call 1/866-484-3268.

Pacemakers

A pacemaker system is designed to monitor and treat your heart rhythm problems, greatly reducing the risks associated with them. But it is not for everyone, including patients with certain steroid allergies. Procedure risks include infection, tissue damage, and kidney failure. Patients who have additional medical conditions that may not allow the pacemaker to function appropriately should not receive a device. In rare cases severe complications or device failures can occur. Electrical or magnetic fields can affect the device. Only your doctor knows what is right for you. This device is available by prescription only. Individual results may vary. Your physician should discuss all potential benefits and risks with you. For further safety information, refer to the Physicians Instructions for Use or call 1/866-484-3268.

Device Quality and Reliability

It is Boston Scientific's intent to provide implantable devices of high quality and reliability. However, these devices may exhibit malfunctions that may result in lost or compromised ability to deliver therapy. Refer to Boston Scientific's CRM product performance report on www.bostonscientific.com for more information about device performance, including the types and rates of malfunctions that these devices have experienced historically. While historical data may not be predictive of future device performance, such data can provide important context for understanding the overall reliability of these types of products. Also, it is important that you talk with your doctor about the risks and benefits associated with the implantation of a device.

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