

# CONTAK RENEWAL<sup>®</sup> TR

Models H120 , H125

## Details About Your Boston Scientific CONTAK RENEWAL<sup>®</sup> TR Cardiac Resynchronization Therapy Pacemaker (CRT-P)

### What is a CRT-P?

CRT-P stands for *cardiac resynchronization therapy pacemaker*. A CRT-P is a small, battery-powered device that holds a tiny computer. It is implanted in the body to provide special pacing therapy to treat heart failure. Your doctor programs your CRT-P to respond appropriately based on your individual heart rhythm needs.

### Treating dyssynchrony

In heart failure, the chambers of the heart do not contract quite as they should. They are *dyssynchronous*, or *out of synch*. When that happens, your heart cannot provide your body with the blood flow it needs to function properly. Cardiac resynchronization therapy (CRT) is designed to help your heart pump better. It delivers specially timed electrical impulses to coordinate, or *resynchronize*, the timing of your heart's contractions. This can help restore your heart's pumping function or ability to fill with blood and work more efficiently.

### Treating slow heart rhythms

If an abnormally slow heart rhythm occurs, your CONTAK RENEWAL TR CRT-P can deliver small electrical signals that help your heart beat more normally. Your CRT-P can also help increase your heart rate to meet your body's needs when you are active.

### How small is my CRT-P?

Height	Width	Thickness	Weight	Volume
2.13 in. 5.40 cm	1.77 in. 4.50 cm	0.33 in. 0.85 cm	0.92 oz. 26 g	14 cc

### How long will my CRT-P last?

Like anything that operates on a battery, the life of your CRT-P will depend on how much the battery is used. The battery is used every time your device delivers therapy. So the amount of therapy you receive, and the amount of energy your doctor programs into those therapies, affects the life of the battery. A CONTAK RENEWAL TR CRT-P typically lasts 5–7 years.



ACTUAL SIZE  
(when printed at 100%)

**Boston  
Scientific**

## What happens after my implant procedure?

### *Care and recovery at home*

Your doctor or nurse will give you specific instructions to follow when returning home. Full recovery from surgery can take from several days to a few months, depending on your condition. It is important that you are actively involved in your own recovery. Follow your doctor's instructions about caring for the implant site and letting your other doctors, dentists, and emergency personnel know that you have an implanted CRT-P system.

### *Follow-up visits*

Once you have a CRT-P, it is important to follow your doctor's guidelines for follow-up visits. Your doctor establishes a periodic follow-up schedule based on his/her preference. Your doctor uses the follow-up visits to check your device and monitor your heart condition, which can change over time. Your doctor may also adjust your device settings.

### **Important Safety Information**

#### **Cardiac Resynchronization Therapy Pacemakers from Boston Scientific CRM**

Cardiac resynchronization therapy pacemakers (CRT-Ps) are used to treat heart failure patients who have symptoms despite the best available drug therapy. These patients also have an electrical condition in which the lower chambers of the heart contract in an uncoordinated way and a mechanical condition in which the heart pumps less blood than normal. CRT-Ps are not for everyone, including people with certain steroid allergies. Procedure risks include infection, tissue damage, and kidney failure. In some cases, the device may be unable to respond to irregular heartbeats or may deliver inappropriate shocks. In rare cases severe complications or device failures can occur. Electrical or magnetic fields can affect the device. Only your doctor knows what is right for you. Find further important safety information at [www.lifebeatonline.com](http://www.lifebeatonline.com). These devices are available by prescription only. Individual results may vary.

(Rev. C)

**Boston  
Scientific**

*Delivering what's next.™*

#### **Cardiac Rhythm Management**

Boston Scientific Corporation  
4100 Hamline Avenue North  
St. Paul, MN 55112-5798 USA  
Tel: 651.582.4000 Fax: 651.582.4166

Medical Professionals:  
1.800.CARDIAC (227.3422)  
Patients and Families:  
1.866.484.3268

**[www.bostonscientific.com](http://www.bostonscientific.com)**

Copyright © 2009 by  
Boston Scientific Corporation  
or its affiliates. All rights reserved.

C9-133-0409