

# SpyBite<sup>®</sup> Biopsy Forceps

*Compatible with the SpyGlass<sup>®</sup> Direct Visualization System*

- SpyBite Biopsy Forceps are compatible with the SpyScope<sup>®</sup> Access and Delivery Catheter and SpyGlass Probe, to provide an intraluminal view for directing therapeutic devices within the biliary duct system.
- Designed to collect tissue samples endoscopically for histologic examination
- Site-specific biopsy specimens can be retrieved for examination from the common bile duct under direct visualization
- Central spike in the specimen cup may aid in securing small tissue samples
- Delivered into the biliary anatomy through a port in the SpyScope Access and Delivery Catheter
- Compatible with delivery devices that have a minimum working channel inner diameter(ID) of 1.2mm

**For additional information on the SpyGlass System, please contact your Boston Scientific sales representative.**

## **INSTRUCTIONS FOR USE**

Refer to the operator's manual for complete instructions for use.

## **INDICATIONS**

The SpyBite Biopsy Forceps are indicated for tissue acquisition in the biliary system.

## **CONTRAINDICATIONS**

Contraindications for this device are those specific to endoscopic retrograde cholangiopancreatography (ERCP) procedures.

## **WARNINGS**

The SpyBite Biopsy Forceps should only be used to biopsy tissue where possible hemorrhage will not present a danger to patients.

The device should be used with caution and only after careful consideration in patients with elevated bleeding times or coagulopathies.

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found call your Boston Scientific representative.

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

## **POTENTIAL ADVERSE EFFECTS**

Complications associated with the use of the SpyBite Biopsy Forceps may include: bleeding, perforation, infection, peritonitis, and inflammation.

Please be aware that potential adverse effects may arise even with the proper use of medical devices. Accordingly, this device should only be used by persons qualified in the procedures for which it is indicated.

## **CAUTIONS**

Cautions can be found in the product labeling supplied with each device. CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

**Boston  
Scientific**

Boston Scientific Endoscopy  
Tel 508.650.8000  
www.bostonscientific.com

Ordering Information  
1.800.225.3226

©2007 Boston Scientific Corporation or its affiliates.  
All rights reserved.