WATCHMAN™

1 2 F

Left Atrial Appendage Closure Device with Delivery System

Directions for Use	3
Instrucciones de uso	8
Mode d'emploi	14
Gebrauchsanweisung	20
Istruzioni per l'uso	26
Gebruiksaanwijzing	32
Instruções de Utilização	38

TABLE OF CONTENTS

WARNING3

DEVICE DESCRIPTION3

 Contents3

INTENDED USE/ INDICATIONS FOR USE3

CONTRAINDICATIONS.....3

WARNINGS.....3

MAGNETIC RESONANCE IMAGING3

 3.0 Tesla Temperature Information3

 1.5 Tesla Temperature Information3

 Image Artifact Information4

PRECAUTIONS.....4

ADVERSE EVENTS4

HOW SUPPLIED.....4

HANDLING AND STORAGE.....4

OPERATIONAL INSTRUCTIONS.....4

 Pre-Procedural Instructions4

PROCEDURAL INSTRUCTIONS.....4

 Equipment Needed for Implantation Procedure.....4

 Figure 1. WATCHMAN™ Delivery System (Delivery Catheter and Device)4

 • **Access Sheath**.....4

 • **Dilator**4

 Figure 2. WATCHMAN Access System (Access Sheath/Dilator)4

 Implantation Procedure4

 TABLE 1. WATCHMAN LAA Closure Device Selection.....5

 Figure 3. WATCHMAN Access Sheath and Pigtail Catheter5

 Figure 4: WATCHMAN Access Sheath Marker Bands5

 Figure 5. WATCHMAN LAA Closure Device Position and Size6

 TABLE 2. WATCHMAN LAA Closure Device Diameter6

WARRANTY6

WATCHMAN™

1 2 F

Left Atrial Appendage Closure Device with Delivery System

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

DEVICE DESCRIPTION

The WATCHMAN LEFT ATRIAL APPENDAGE (LAA) Closure Technology consists of the Access System (Access Sheath and Dilator) and Delivery System (Delivery Catheter and LAA Closure Device). The Access System and Delivery System permit Device placement in the LAA via femoral venous access and inter-atrial septum crossing into the left atrium. The WATCHMAN Device is a self-expanding nitinol structure with a porous membrane on the proximal face. The Device is constrained within the Delivery System until deployment in the LAA. The Device is available in 5 sizes from 21 to 33 mm. Device selection is determined by LAA measurements using Fluoroscopy (Fluoro) and Transesophageal Echocardiography (TEE).

The WATCHMAN LAA Closure Device is designed to be permanently implanted at or slightly distal to the ostium (opening) of the LAA to trap potential emboli before they exit the LAA. The placement procedure can be done under local or general anesthesia in a catheterization laboratory setting.

Contents

Quantity	Description
1	WATCHMAN Left Atrial Appendage Closure Device with Delivery System

INTENDED USE/ INDICATIONS FOR USE

The WATCHMAN LAA Closure Technology is intended to prevent thrombus embolization from the left atrial appendage and reduce the risk of life-threatening bleeding events in patients with non-valvular atrial fibrillation who are eligible for anticoagulation therapy or who have a contraindication to anticoagulation therapy.

CONTRAINDICATIONS

Do not use the WATCHMAN LAA Closure Device if:

- Intracardiac thrombus is visualized by echocardiographic imaging.

Note: If thrombus is identified in the LAA, dissolve with anticoagulation therapy before attempting to implant the WATCHMAN Device.

- An atrial septal repair or closure device is present.
- The LAA anatomy will not accommodate a Device. See Table 1
- Any of the customary contraindications for other percutaneous catheterization interventions eg. patient size (i.e. too small for TEE probe, catheter size, etc.) or condition (i.e. active infection, bleeding disorder, untreated ulcer, etc.) are present.

WARNINGS

Implantation of the WATCHMAN LAA Closure Device should only be performed by physicians trained in percutaneous and transseptal procedures who have completed the WATCHMAN training program.

- **The LAA is a thin walled structure. Use caution when accessing the LAA and deploying the Device.**
- The WATCHMAN Access and Delivery Systems are sterile and intended for single use only. Do not reuse or resterilize. Reuse could result in product damage and/

or breakage that could lead to clinical complications, possibly requiring prolonged hospitalization. Resterilization could result in product contamination resulting in infection (e.g. endocarditis/sepsis/local infection), possibly requiring antibiotics or prolonged hospitalization.

- Careful consideration should be given to use of the Device in pregnant and/or breastfeeding women due to the risk of significant exposure to x-rays and the use of strong anticoagulation medication.
- The WATCHMAN device has not been studied in patients under the age of 18.
- Device selection should be based on accurate LAA measurements obtained using Fluoro and TEE in multiple angles (e.g. 0°, 45°, 90°, 135°).
- Aspirin should be started one day prior to scheduled procedure and continued daily.
- Patients should be fully heparinized throughout the procedure with an active clotting time (ACT) of 200 - 300 seconds after transseptal puncture.
- Fluoro and TEE should be used when implanting the Device.
- Do not release (unscrew) the Device unless release criteria (step 14) are satisfied.
- Potential for Device embolization exists with cardioversion < 30 days following Device implantation, verify Device position post cardioversion.
- Post- procedure warfarin therapy is required in ALL patients receiving a Device who are eligible for warfarin therapy or other equivalent oral anticoagulant per institution's protocol. Patients should remain on 81-100 mg of aspirin and warfarin for a minimum of 45 days post implant (INR 2.0-3.0). At 45 days post implant perform Device assessment with TEE. Cessation of warfarin is at physician discretion. Patients ceasing warfarin should begin clopidogrel 75mg daily and increase aspirin dosage to 300-325mg daily for 6 months post-implant and remain on aspirin 300-325mg indefinitely. For patients contraindicated to anticoagulation therapy, begin clopidogrel 75mg daily and aspirin 300-325mg daily for 6 months post-implant and remain on aspirin 300-325mg indefinitely. A visual assessment with TEE should be performed 45 days post implant to assess the device.
- Administer appropriate endocarditis prophylaxis for 6 months following device implantation. The decision to continue endocarditis prophylaxis beyond 6 months is at physician discretion.

MAGNETIC RESONANCE IMAGING

The WATCHMAN LAA Closure Device was determined to be MR Conditional per: (International, Designation: F2503-05. Standard Practice for Marking Medical Devices for Safety in the MR Environment. ASTM International, 100 Barr Harbor Drive, PO Box C700, West Conshohocken, Pennsylvania, 19428, 2005).

Non-clinical testing has demonstrated the WATCHMAN device is MR Conditional. A patient with the Device can be scanned safely under the following conditions:

- Static magnetic fields of 3 Tesla or 1.5 Tesla
- Spatial gradient field of 2500 Gauss/cm or less
- The maximum whole body averaged specific absorption rate (SAR) shall be limited to 2.0 W/kg (normal operating mode only) for 15 minutes of continuous application of RF energy during a scan
- Normal operating mode of the MRI scanner

The WATCHMAN LAA Closure Device should not migrate in this MRI environment. MR imaging within these conditions may be performed immediately following the implantation of the device. This device has not been evaluated to determine if it is MR Conditional beyond these parameters.

3.0 Tesla Temperature Information

In non-clinical testing, the WATCHMAN LAA Closure Device produced a temperature rise of < 1.1°C at a maximum MR system-reported SAR of 2.0 W/kg as measured by calorimetry for 15 minutes of continuous MR scanning in a 3 Tesla MR system (Excite, Software G3.0-052B, GE Healthcare, Milwaukee, WI).

The actual in vivo temperature rise is expected to be less than these values as the calculations did not include the cooling effects due to blood perfusion in the tissue outside the WATCHMAN device. In vivo, local SAR depends on MR field strength and may be different than the estimated whole body averaged SAR, due to body composition, device position within the imaging field, and scanner used, thereby affecting the actual temperature rise. No tests have been performed on possible nerve or other tissue stimulation activated by strong gradient magnetic fields and the resulting induced voltages.

1.5 Tesla Temperature Information

Non-clinical testing of RF-induced heating in the WATCHMAN LAA Closure Device was performed at 64 MHz in a 1.5 Tesla whole body coil MR scanner (Intera, Software Release 10.6.2.4, 2006-03-10, Philips Medical Systems, Andover, MA) and produced a temperature rise of < 1.5°C at an MR extrapolated SAR of 2.0 W/kg for 15 minutes of continuous MR scanning.

The actual in vivo temperature rise is expected to be less than these values as the calculations did not include the cooling effects due to blood perfusion in the tissue

outside the WATCHMAN™ device. In vivo, local SAR depends on MR field strength and may be different than the estimated whole body averaged SAR, due to body composition, device position within the imaging field, and scanner used, thereby affecting the actual temperature rise. No tests have been performed on possible nerve or other tissue stimulation activated by strong gradient magnetic fields and resulting induced voltages.

Image Artifact Information

MR image quality may be compromised if the area of interest is relatively close to the WATCHMAN device. Optimization of MR imaging parameters is recommended.

PRECAUTIONS

None known.

ADVERSE EVENTS

Adverse events reported during the clinical study were adjudicated by the Clinical Events Committee and are as follows:

- Pericardial Effusion
- Excessive Bleeding
- Gastrointestinal Bleeding
- Stroke - Ischemic
- Death
- Device Thrombus
- Hypotension
- Major Bleed Requiring Transfusion
- Bruising - Hematoma
- Inability to Move or Retrieve Device
- Anemia Requiring Transfusion
- Allergic Reaction to Contrast Media/Meds
- Infection/Pneumonia
- Device Embolization
- Pseudoaneurysm
- TEE Complications (throat pain, bleeding)
- Systemic Embolism
- Pleural Effusion
- AV Fistula
- Transient Ischemic Attack (TIA)
- Arrhythmias
- Thrombus at Septal Puncture
- Post Procedure Anesthesia Effects
- Thrombosis
- Groin Puncture Bleed
- Pulmonary Edema
- Vasovagal Reactions
- Cranial Bleed
- Air Embolism
- Stroke – Hemorrhagic

Some additional events that may be expected in catheterization procedures include:

- Pneumothorax
- Pulmonary Vein Obstruction
- Valvular or vascular damage

HOW SUPPLIED

- The WATCHMAN LAA Closure Device is pre-loaded in the Delivery System.
- The WATCHMAN Access System is packaged separately.
- The WATCHMAN LAA Closure products are supplied STERILE using an ethylene oxide (EO) process.
- Do not use if package is opened or damaged.
- Do not use if labeling is incomplete or illegible.

Note: Contents of inner package are STERILE.

HANDLING AND STORAGE

Store in a cool, dry, dark place.

OPERATIONAL INSTRUCTIONS

Pre-Procedural Instructions

A baseline TEE should be performed to verify that a WATCHMAN LAA Closure Device may be implanted.

1. Assess the following through multiple imaging planes (0° - 135° sweep):
 - LAA size /shape, number of lobes in LAA, and location of lobes to ostium.
 - Confirm the absence of thrombus (use Color Doppler and echo contrast as necessary).
2. Record LAA ostium and LAA length measurements (0° - 135° sweep). Measure the LAA ostium at approximately these angles.
 - at 0° measure from coronary artery marker to a point 2 cm from tip of the “limbus”
 - at 45° measure from top of the mitral valve annulus to a point 2 cm from tip of the “limbus”
 - at 90° measure from top of the mitral valve annulus to a point 2 cm from tip of the “limbus”
 - at 135° measure from top of the mitral valve annulus to a point 2 cm from tip of the “limbus”

Measured maximum LAA ostium width must be $\geq 17\text{mm}$ or $\leq 31\text{mm}$ to accommodate available device sizes.

Note: The maximum LAA ostium and LAA length measurements determine device size selection.

PROCEDURAL INSTRUCTIONS

Equipment Needed for Implantation Procedure

- Venous Introducer (optional)
- Standard transseptal access system
- 0.035” guidewire (exchange length extra support)
- 6Fr Pigtail Catheter




Figure 1. WATCHMAN Delivery System (Delivery Catheter and Device)

• Access Sheath

• Dilator

Figure 2. WATCHMAN Access System (Access Sheath/Dilator)

Implantation Procedure

WARNING: The use of echocardiographic imaging is required (TEE is recommended as an aid in placing the WATCHMAN Device).

Note: Patients should be fully heparinized throughout the procedure with a recommended minimum active clotting time (ACT) of 200-300 seconds after transseptal puncture.

1. Use standard practice to puncture vessel and insert 0.035” guidewire and vessel dilator. Use a standard transseptal access system to cross inter-atrial septum.

2. Exchange crossing sheath with exchange length extra support 0.035" guidewire. Position guidewire in left upper pulmonary vein (LUPV) or loop in left atrium.
3. Prepare WATCHMAN™ Access System.
 - A. Remove Access Sheath and Dilator under sterile conditions.
 - B. Inspect prior to use to ensure no damage.

Note: Inspect sterile package and WATCHMAN Access System prior to use. If sterile barrier has been compromised in any way, DO NOT USE.

- C. Flush Access Sheath and Dilator with saline prior to use.
- D. Insert Dilator into hemostasis valve of Access Sheath.

Note: Audible snap upon connection of Dilator and hemostasis valve of Access Sheath.

Note: Do not tighten the hemostasis valve while the Dilator is inserted in the WATCHMAN Access System. The Dilator by itself will occlude the lumen of the WATCHMAN Access System creating hemostasis. Tightening the valve onto the Dilator may damage the valve threads, which can lead to subsequent difficulty in closing the valve and an incomplete seal, once the Dilator is removed.

4. Advance WATCHMAN Access System over guidewire into left atrium (LA). As Access Sheath nears center of LA, hold Dilator and advance Access Sheath into initial position in LA or ostium of LUPV.

Note: Use caution when introducing WATCHMAN Access System to prevent damage to cardiac structures.

5. Remove Dilator and guidewire, leaving Access Sheath. Allow back bleed to minimize potential for introducing air before tightening valve. Flush with saline. If continued back bleed is observed from the valve after the Dilator is removed despite attempting to close it, loosen the valve cap (counter-clockwise rotation) until the cap spins freely. Then re-attempt closure of the valve while exerting gentle forward pressure on the valve cap during closure (clockwise rotation) to ensure proper engagement of the valve thread. While these steps are being undertaken, manual occlusion of the valve opening using a gloved finger is recommended to minimize blood loss.

Note: These steps may be repeated if necessary. However, if this does not mitigate the blood leak, the user should remove and replace the WATCHMAN Access Sheath before proceeding with the procedure.

6. Confirm LAA size and select appropriate WATCHMAN LAA Closure Device.
 - A. Using TEE, measure LAA ostium width and LAA length in 3-4 views (0°, 45°, 90°, 135°).
 - B. Choose a device based on **maximum** LAA ostium width recorded. Use Table 1 as a guide.

Note: LAA anatomy should accommodate a Device as described in Table 1.

TABLE 1. WATCHMAN LAA Closure Device Selection

Max LAA Ostium (mm)	Device Diameter (mm)
17 – 19	21
20 – 22	24
23 – 25	27
26 – 28	30
29 – 31	33

Note: Record multiple angles on cine with contrast prior to advancing Access Sheath into LAA. Use fluoro guidance while advancing pigtail catheter, or Access Sheath. Stop if resistance is felt.

- C. Carefully advance pigtail catheter through Access Sheath into distal LAA under fluoro guidance (Figure 3). Carefully advance Access Sheath over pigtail catheter until Access Sheath marker band corresponding to Device size (see Figure 4) is at or just distal to LAA ostium. Slowly remove pigtail catheter.


Figure 3. WATCHMAN Access Sheath and Pigtail Catheter


Figure 4. WATCHMAN Access Sheath Marker Bands

7. Prepare WATCHMAN Delivery System.
 - A. Remove Delivery System under sterile conditions.
 - B. Inspect prior to use to ensure no damage to handle, catheter connections and Device (through Delivery System).

Note: If sterile barrier has been compromised in any way, or Delivery System appears damaged DO NOT USE.

- C. Confirm that the distal tip of the Device is aligned with marker band on Delivery System.

CAUTION: Do not allow WATCHMAN Device to protrude to prevent damage to Delivery Catheter.

- D. Flush system with saline removing all air and maintaining fluid throughout Delivery System. Open and flush proximal valve.

Note: To avoid introducing air, apply pressurized saline bag to sideport of Access Sheath, or submerge Access Sheath hub in saline. Saline may be dripped from Delivery System during introduction into Access Sheath by injecting through flush port.

8. Loosen proximal valve of Access Sheath allowing bleed back before inserting Delivery System. Note: Hemostasis valve should spin freely (fully open).

Note: Tightening the valve onto the WATCHMAN Delivery System may damage the valve threads, which can lead to subsequent difficulty in closing the valve and an incomplete seal, once the WATCHMAN Delivery System is removed.

9. To avoid introduction of air, slowly advance Delivery System into Access Sheath under fluoro guidance.

Note: Use caution when introducing Delivery System to prevent damage to cardiac structures.

10. On fluoro, align most distal marker band on Delivery System with most distal marker band on Access Sheath. Once bands are aligned, stabilize Delivery System, retract Access Sheath and snap together as Access Sheath/Delivery System Assembly.

11. Using fluoro and TEE confirm position of Delivery System tip before deploying the Device.

Note: To inject contrast, flush catheter or measure power injector pressure while inserting Delivery System into Access Sheath. Contrast syringe or manifold must be attached to flush port of Delivery System. If using a power injector, the maximum pressure **should not** exceed 100 psi.

12. If repositioning is required, unsnap and slowly remove Delivery System from Access Sheath. If necessary reinsert pigtail catheter to reposition Access Sheath. Reinsert Delivery System as described in Steps 9 and 10.
13. Deploy WATCHMAN™ Device by loosening valve on Delivery System and holding deployment knob stationary while retracting Assembly to completely deploy Device. Leave core wire attached.
14. Device release criteria:
- A. **Position:** Plane of maximum diameter is at or just distal to and spans entire LAA Ostium (See Figure 5).


Figure 5. WATCHMAN LAA Closure Device Position and Size

- B. **Anchor:** Gently pull back then release deployment knob to visualize movement of Device and LAA together.
- C. **Size (compression):** Measure plane of maximum diameter of Device (See Figure 5). Use Table 2 as a guide.
- D. **Seal:** Ensure all lobes are distal to Device and sealed.

TABLE 2. WATCHMAN LAA Closure Device Diameter

Original Diameter (mm)	Deployed Diameter (80-92% of original) (mm)
21	16.8-19.3
24	19.2-22.1
27	21.6-24.8
30	24.0-27.6
33	26.4 -30.4

15. Partial Device Recapture

Note: Partially recapture and redeploy WATCHMAN Device if too distal to LAA ostium

- A. Advance tip of Access/Delivery System Assembly up to Device (do not unsnap). Fix deployment knob position with right hand and gently advance Access/Delivery System Assembly over shoulders of Device. Position right thumb against Delivery System hub for stability. Resistance will be felt as Device shoulders collapse. Continue to advance Assembly up to but not past fixation anchors. When resistance is felt a second time (anchor contact), stop, tighten hemostasis valve.

Note: If Device is retrieved past fixation anchors, recapture fully and replace Delivery System. Refer to Step 16. The WATCHMAN Device and Delivery System are for single use only. Do not reuse or resterilize.

- B. Reposition Access Delivery/System Assembly proximally and re-deploy by holding deployment knob and retracting Access Sheath until Device is completely deployed. Leave core wire attached.

WARNING: Do not release the WATCHMAN Device from the core wire if the Device does not meet release criteria (Step 14).

16. Full Device recapture.

Note: Fully recapture the Device if too proximal or does not meet release criteria

- A. Advance tip of Access/Delivery System Assembly up to face of Device (do not unsnap).
- B. Fix deployment knob with right hand and gently advance Access/Delivery System Assembly over shoulders of Device. Position right thumb against Delivery System for stability. Resistance will be felt as Device shoulders collapse. Continue to advance Assembly until Device is completely collapsed and recaptured (past anchors).
- C. Withdraw Device until distal tines are proximal to marker band then tighten hemostasis valve.
- D. Unsnap Delivery System from Access Sheath while maintaining position. Slowly remove Delivery System.
- E. Insert pigtail catheter to reposition Access Sheath in LAA if necessary.
- F. Repeat Steps 7-14 with new Delivery System.
17. WATCHMAN Device Release: Confirm proper position, anchor, size, and seal, and then advance Assembly to face of Device. Rotate deployment knob counter clockwise 3-5 full turns. Confirm core wire is disconnected.
18. Remove Access Sheath and Delivery System based on parameters for hemostasis.
19. Use standard of care for post procedure bleeding at access site.
20. Post Procedure Information

- A. Post-procedure warfarin therapy is required in ALL patients receiving a Device who are eligible for warfarin therapy or other equivalent oral anticoagulant per institution's protocol. Patients should remain on 81-100 mg of aspirin and warfarin for a minimum of 45 days post implant (INR 2.0-3.0). At 45 days post implant perform Device assessment with TEE. Cessation of warfarin is at physician discretion. Patients ceasing warfarin should begin clopidogrel 75mg daily and increase aspirin dosage to 300-325mg daily for 6 months post-implant and remain on aspirin 300- 325mg indefinitely. For patients contraindicated to anticoagulation therapy, begin clopidogrel 75mg daily and aspirin 300-325mg daily for 6 months post-implant and remain on aspirin 300-325mg indefinitely. A visual assessment with TEE should be performed 45 days post implant to assess the device.
- B. At 45 days assess WATCHMAN Device with TEE.
- Confirm absence of intra-cardiac thrombus.
 - Perform color Doppler assessment to include the device/LAA border at the following approximate TEE angles (0°, 45°, 90° and 135°). Measure any residual jet around the device if necessary.
- C. Prescribe appropriate endocarditis prophylaxis for 6 months following Device implantation. The decision to continue endocarditis prophylaxis beyond 6 months is at physician discretion.

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or reesterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

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© 2015 Boston Scientific Corporation or its affiliates. All rights reserved.
Boston Scientific Limited
Boston Scientific (Australia) Pty Ltd