IDC™ Interlocking Detachable Coil

Directions for Use	2
Instrucciones de uso	6
Mode d'emploi	11
Gebrauchsanweisung	16
Istruzioni per l'uso	20
Gebruiksaanwijzing	25
Instruções de Utilização	29

TABLE OF CONTENTS

WARNING	2
DEVICE DESCRIPTION	2
Figure 1. IDC Coil	2
Contents.....	2
INTENDED USE/INDICATIONS FOR USE	2
CONTRAINDICATIONS	2
WARNINGS	3
GENERAL PRECAUTIONS	3
MAGNETIC RESONANCE IMAGING (MRI) SAFETY INFORMATION	3
Medical Registration	3
ADVERSE EVENTS	3
COIL SIZE SELECTION	3
HOW SUPPLIED	3
Handling and Storage.....	3
DEVICE PREPARATION	3
Figure 2. Example of Continuous Flush Setup.....	4
DIRECTIONS FOR USE	4
Figure 3. Introducer Sheath with Twist-Lock Mechanism	4
Figure 4. Unlock via Counter-Clockwise Rotation.....	4
Figure 5. Advancing IDC Coil to Pre-Release Position (1-RO marker microcatheter)	4
Figure 6. Delivering IDC™ Coil (1-RO marker microcatheter).....	5
WARNINGS	5
PRECAUTIONS	5
IDC COIL REMOVAL PROCEDURE	5
Figure 7. Lock via clockwise rotation.....	5
WARRANTY	5

IDC™

Interlocking Detachable Coil

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

Do not use catheters, delivery wires or IDC Coils that have been damaged in any way. Damaged catheters may rupture causing vessel trauma or tip detachment during steering maneuvers. Damaged delivery wires and IDC Coils may cause vessel trauma or unpredictable distal tip response during steering maneuvers. Take care not to puncture gloves or sterile drapes when handling delivery wire.

DEVICE DESCRIPTION

The IDC Coil includes a standard or soft coil (manufactured from platinum-tungsten alloy) that is mechanically attached to a coil delivery wire. This assembly is contained within an introducer sheath. The IDC Coil is designed to be delivered under fluoroscopy with a 0.53 mm (0.021 in) inner diameter (I.D.) microcatheter (e.g. Renegade™ Microcatheter) with one radiopaque (RO) tip marker. The interlocking delivery wire design allows the coil to be advanced and retracted before final placement in the vessel, thus aiding in more controlled delivery including the ability to withdraw the coil prior to deployment.


Figure 1. IDC Coil

Contents

The IDC Coil includes an embolic coil with interlocking delivery wire and introducer sheath.

INTENDED USE/INDICATIONS FOR USE

The IDC Coil is a modified interlocking detachable coil indicated to obstruct or reduce rate of blood flow in the peripheral vasculature. This device is not intended for neurovascular use.

CONTRAINDICATIONS

None known.

WARNINGS

Recanalization has been observed with the usage of some coils. Angiographic follow-up is recommended to ensure continued occlusion.

Do not advance the delivery wire once the coil has been placed. Perforation or damage of the vessel wall could occur.

GENERAL PRECAUTIONS

Carefully read all instructions prior to use. Observe all warnings and precautions noted throughout these and other instructions relevant to the procedure. Failure to do so may result in complications.

To prevent damage to the twist-lock function of the introducer sheath, take care not to excessively tighten or loosen the twist-lock mechanism.

Do not apply excessive force while seating the introducer sheath in the microcatheter hub. Introducer sheath tip deformation and difficulty with coil delivery into the microcatheter could result.

Multiple embolization procedures may be required to achieve the desired occlusion of some vessels. To position another IDC™ Interlocking Detachable Coil return to Steps 1-9 in Directions For Use.

Replace microcatheters periodically during delivery of multiple coils or if increased resistance is noted during coil delivery.

Axial compression or tension forces may be stored in the microcatheter shaft during IDC Coil delivery, and coil release may lead to catheter tip movement. Verify repeatedly during the procedure that the distal shaft of the microcatheter is not under stress prior to IDC Coil detachment by slightly repositioning the microcatheter, delivery wire, or entire assembly simultaneously.

Carefully remove the delivery wire after coil deployment so that the delivery arm does not catch on the valve in the RHV thumbscrew.

Do not rotate the delivery wire more than one turn (360 degrees) during delivery of the IDC Coil. **Excessive rotation of the delivery wire may damage the coil or may result in premature detachment of the interlocking arms within the microcatheter.**

Do not advance the IDC Interlocking Detachable Coil if it **becomes lodged** within the microcatheter. Determine the cause of the resistance and replace the microcatheter and coil if necessary. See IDC Coil Removal Procedure for further instructions.

Advance and retract the IDC Coil smoothly, especially in tortuous anatomy. Replace the coil if excessive friction is noted within the microcatheter. If friction is noted in any successive coil, carefully examine both coil and microcatheter for possible damage. Replace both if necessary.

Do not retract the IDC Coil too quickly or against resistance. Doing so may result in a stretched coil or damage to the interlocking mechanism.

MAGNETIC RESONANCE IMAGING (MRI) SAFETY INFORMATION

Non-clinical testing has demonstrated that the IDC Coil is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 Tesla and 3.0 Tesla only
- Maximum spatial gradient magnetic field of 2500 gauss/cm (25 T/m) or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode)
- The IDC Coil should not migrate in this MRI environment. Non-clinical testing at field strengths other than 1.5 Tesla and 3 Tesla has not been performed to evaluate coil migration or heating. MR Imaging may be performed immediately following the implantation of the IDC coil(s).

Under the scan conditions defined above, the IDC Coil is expected to produce a maximum temperature rise of the following after 15 minutes of continuous scanning:

1.5 Tesla: For vessels in the body the calculated temperature rise was 5.2 °C with an uncertainty upper bound temperature of 7.1 °C for a whole body average SAR value of 2.0 W/kg and a continuous scan time of 15 minutes.

3.0 Tesla: For vessels in the body the calculated temperature rise was 2.8 °C with an uncertainty upper bound temperature of 3.8 °C for a whole body average SAR value of 2.0 W/kg and a continuous scan time of 15 minutes.

In non-clinical testing, the image artifact caused by the device extends approximately 2 mm from the IDC Coil when imaged with a Spin Echo pulse sequence and a 3 Tesla system. With a Gradient Echo pulse sequence, the image artifact extends approximately 5 mm beyond the IDC Coil with images made in a 3 Tesla system.

Medical Registration

It is recommended that patients register the conditions under which the implant can be scanned safely with the Medic Alert Foundation (www.medicalert.org) or equivalent organization.

ADVERSE EVENTS

The complications that may result from a peripheral embolization procedure include, but are not limited to:

- Complications related to catheterization (e.g., hematoma at the site of entry, vessel injury, etc.)
- Death
- Emboli
- Foreign body reactions necessitating medical intervention
- Hemorrhage
- Infection necessitating medical intervention
- Ischemia
- Pain
- Recanalization
- Temporary neurological deficit
- Tissue necrosis
- Undesirable clot formation of the vasculature
- Vasospasm

COIL SIZE SELECTION

Coil selection is a matter of physician preference and the clinical situation. The shape and diameter of the vessel to be occluded as well as proximity to branch vessels generally govern selection of the coil diameter and length. Coil diameter should approximate the vessel diameter. Selection of a coil diameter >1 mm larger than the vessel diameter may result in coil elongation and a non-compact placement with less effective reduction of blood flow. Selection of a coil diameter smaller than the vessel diameter may result in coil migration.

HOW SUPPLIED

The IDC Coil includes an embolic coil with interlocking delivery wire and introducer sheath.

Do not use if package is opened or damaged.

Do not use if labeling is incomplete or illegible.

Handling and Storage

Store in cool, dry, dark place.

DEVICE PREPARATION

Prior to use, ensure that sterile packaging is intact. Return the device if sterility appears to have been compromised.

Remove the IDC Coil from its packaging and inspect for product damage. Ensure that the interlocking arms are interlocked within the introducer sheath. Do not advance the coil outside the introducer sheath. Do not use the IDC Coil System if it is damaged. Do not use if labeling is incomplete or illegible.

In order to achieve excellent performance of the IDC™ Coil and reduce the risk of thromboembolic complications, it is critical that a continuous flow of appropriate flush solution be maintained between a) the microcatheter and guiding catheter, and b) the microcatheter and any intraluminal device. Continuous flushing:

- Reduces retrograde blood flow into the microcatheter and introducer sheath during coil delivery.
- Reduces contrast crystal formation and/or thrombosis on the delivery wire and in the guiding catheter and microcatheter lumens.
- Reduces premature coil thrombosis.

An example of a continuous flush setup is shown in Figure 2.


Figure 2. Example of Continuous Flush Setup

DIRECTIONS FOR USE

1. Place the microcatheter in the area to be embolized per standard technique. Take care to position the microcatheter tip parallel with, not perpendicular to, the vessel wall to facilitate delivery of the coil.
2. Slowly withdraw the IDC Interlocking Detachable Coil from its dispenser coil and inspect assembly. Discard if there is any evidence of damage. Ensure that the IDC Coil arms are interlocked inside the introducer sheath. Do not remove the IDC Coil assembly from the introducer sheath.
3. Release the IDC Coil inside its introducer sheath by gently pinching the sheath on both sides of the twist-lock mechanism and rotating the proximal side counter-clockwise 2-3 rotations (Figure 3 and 4).
4. Once the IDC Coil is freed, carefully advance to within 2 mm of introducer sheath tip, then re-engage twist-lock (Figure 7).

Caution: To prevent damage to the twist-lock function of the introducer sheath, take care not to excessively tighten or loosen the twist-lock mechanism.


Figure 3. Introducer Sheath with Twist-Lock Mechanism


Figure 4. Unlock via Counter-Clockwise Rotation

5. Attach a Rotating Hemostatic Valve (RHV) to the proximal luer adapter on the hub of the microcatheter. Begin continuous flow of an appropriate flush solution. In general, one drop of flush solution every 1-3 seconds from a pressure bag containing the flush solution is recommended.
6. Open the thumbscrew of the RHV and carefully insert the IDC Coil until the distal tip of the introducer sheath is firmly seated in the microcatheter hub.

Caution: Do not apply excessive force while seating the introducer sheath in the microcatheter hub. Introducer sheath tip deformation and difficulty with coil delivery into the microcatheter could result.

7. Tighten the RHV thumbscrew just enough to prevent retrograde flow but not so tight as to pinch the introducer sheath and inhibit forward movement of the delivery wire. Maintain in-line pressure of the continuous flush to prevent retrograde flow once the sheath is removed.
8. Disengage the twist-lock mechanism, then transfer the IDC Coil and delivery wire from the introducer sheath into the microcatheter by advancing the delivery wire in a smooth, continuous manner. **Ensure that the introducer sheath remains firmly seated in the microcatheter hub to prevent premature deployment.**
9. Gently withdraw and remove the introducer sheath from the microcatheter once the proximal end of the delivery wire is within 10 cm of the proximal end of the sheath. **Do not discard the sheath in case it is necessary to remove the IDC Coil prior to deployment.**

1-RO marker microcatheter: Maneuver the IDC Coil under fluoroscopy until the coil detachment zone is approximately 1 cm proximal to the microcatheter radiopaque tip marker (Figure 5). This positions the interlocking arms approximately 1 cm proximal to the microcatheter tip. **Do not advance delivery wire further until ready to release coil.**


Figure 5. Advancing IDC Coil to Pre-Release Position (1-RO marker microcatheter)


Figure 6. Delivering IDC™ Coil (1-RO marker microcatheter)

If IDC Coil repositioning is necessary, gently retract the coil under fluoroscopy. If repositioning is unsuccessful, remove and discard the IDC Coil. (See IDC Coil Removal Procedure for further instructions.)

To deploy the coil, slowly advance the delivery wire under fluoroscopy until interlocking arms pass microcatheter's tip marker (Figure 6).

WARNINGS

Recanalization has been observed with the usage of some coils. Angiographic follow-up is recommended to ensure continued occlusion.

Do not advance the delivery wire once the coil has been placed. Perforation or damage of the vessel wall could occur.

PRECAUTIONS

Multiple embolization procedures may be required to achieve the desired occlusion of some vessels. To position another IDC Interlocking Detachable Coil return to Steps 1-9 in Directions For Use.

Replace microcatheters periodically during delivery of multiple coils or if increased resistance is noted during coil delivery.

Axial compression or tension forces may be stored in the microcatheter shaft during IDC Coil delivery, and coil release may lead to catheter tip movement. Verify repeatedly during the procedure that the distal shaft of the microcatheter is not under stress prior to IDC Coil detachment by slightly repositioning the microcatheter, delivery wire, or entire assembly simultaneously.

Carefully remove the delivery wire after coil deployment so that the delivery arm does not catch on the valve in the RHV thumbscrew.

IDC COIL REMOVAL PROCEDURE

An IDC Coil must be removed if the coil is determined to be the incorrect size. If resistance is noted and repositioning is unsuccessful, remove and discard. Make sure that the introducer sheath twist-lock mechanism is disengaged to facilitate threading of the sheath over the delivery wire (Figure 4).

1. Gently begin to retract the IDC Coil under fluoroscopy. If resistance is encountered, retract the microcatheter and delivery wire simultaneously to facilitate movement.
2. Once the IDC Coil is withdrawn to approximately the mid-shaft point in the microcatheter, gently thread the distal end of the introducer sheath over the proximal end of the delivery wire.
3. Open the thumbscrew of the RHV and carefully advance the introducer sheath until it is firmly seated in the proximal luer adapter of the microcatheter.
4. Tighten the RHV thumbscrew just enough to prevent retrograde flow but not so tight as to inhibit backward movement of the delivery wire through the microcatheter.
5. Holding the introducer sheath in place, gently withdraw the IDC Coil until the interlocking arms and distal coil tip are visible inside the sheath.

6. Lock the IDC Coil into position by gently pinching the introducer sheath on both sides of the twist-lock mechanism and rotating the proximal side clockwise (Figure 7).


Figure 7. Lock via clockwise rotation

7. Withdraw the introducer sheath/delivery wire assembly from the microcatheter/RHV assembly.

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

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© 2015 Boston Scientific Corporation or its affiliates. All rights reserved.
Boston Scientific Limited
Boston Scientific (Australia) Pty Ltd