TABLE OF CONTENTS

DEVICE DESCRIPTION	2
INTENDED USE/INDICATIONS FOR USE	3
CONTRAINDICATIONS.....	3
ADVERSE EVENTS	4
HOW SUPPLIED.....	4
SYSTEM CONTROLS, DISPLAYS AND INDICATORS	5
SYSTEM INSTALLATION.....	10
OPERATION.....	17
SERVICE AND MAINTENANCE	20
PRODUCT SPECIFICATIONS	21
OPERATIONAL MESSAGES.....	32
DIAGNOSTIC MESSAGES	33
LIMITED WARRANTY AND DISCLAIMER.....	36

THESE DOCUMENTS ARE THE PROPERTY OF BOSTON SCIENTIFIC CORP. AND SHALL NOT BE REPRODUCED, DISTRIBUTED, DISCLOSED OR USED FOR MANUFACTURE OR SALE OF APPARATUS WITHOUT THE EXPRESS WRITTEN CONSENT OF BOSTON SCIENTIFIC CORP.

WARNING: Do not attempt to operate the Maestro 4000™ Cardiac Ablation System prior to reading this entire Operator's Manual. All instructions should be understood, and followed carefully. For future reference, keep this manual in a convenient, readily accessible place.

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

DEVICE DESCRIPTION

Maestro 4000 Cardiac Ablation System

The Maestro 4000 Cardiac Ablation System is a non-sterile product that supplies radiofrequency energy to catheters during ablation procedures. The System is compatible with all BSC cardiac ablation catheters including BSC open-irrigated catheters used in conjunction with the MetriQ™ Pump and MetriQ Irrigation Tubing Set.

Note: Throughout this manual, the term "System" refers to the Maestro 4000 Cardiac Ablation System.

Maestro 4000 Controller

The Maestro 4000 Controller is an RF Generator specifically designed for cardiac ablation. It produces user-selectable power-controlled or temperature-controlled RF power output in the range of 0 to 150 watts into a nominal tissue impedance of 100 ohms. It delivers RF power via a monopolar method driving current between a single active electrode at the tip of the ablation catheter and one or two dispersive pads applied on the skin. When used with BSC open-irrigated catheters, the RF Generator communicates with the MetriQ Pump to coordinate delivery of RF energy with irrigation flow to the catheter tip.

Note: Throughout this manual, the term RF Generator, or Generator, refers to the Maestro 4000 Controller.

Maestro 4000 Pod

The Pod, which is connected to the RF Generator, allows connection to all BSC catheters and provides connections for dispersive pad(s) to complete the RF circuit. The Pod also connects to EP recording systems and provides RF filtering to allow continuous electrogram recording during RF delivery. The Pod model varies by operating territory and determines the maximum power setting allowed by the RF Generator.

Maestro 4000 Remote Control

The optional Remote allows the user to control the RF Generator with up to 75 feet between the user interface and the sterile field.

Maestro Foot Switch

The optional Foot Switch provides hands-free control to start/stop RF delivery.

Dispersive Pads

A Dispersive Pad provides external patient contact to complete the RF circuit. It disperses current over a large area to minimize damage due to heating of skin and underlying tissue. It is also known as the DIP (dispersive indifferent pad), ground pad, or return pad.

User Profile

The Maestro 4000™ Cardiac Ablation System should be used only by physicians fully trained in cardiac electrophysiology. EP lab staff members may prepare the Maestro 4000 Cardiac Ablation System for use and assist with operation of the RF Generator.

Contents

The RF Generator and Pod are packaged separately as:

- (1) Maestro 4000 Controller
- (1) Maestro 4000 Pod, 100W or (1) Maestro Pod, 150W

Optional Accessories

The MetriQ™ Pump, Maestro 4000 Remote Control, Maestro 4000 Foot Switch, and communication cables required for optional equipment such as the MetriQ Pump, Maestro 4000 Remote Control, or connection to an EP recording system are packaged and sold separately.

- MetriQ Pump, Model 4100
- Maestro 4000 Remote Control, Model 4020
- Maestro 4000 Foot Switch, Model 21850
- Cable, Generator to Pump or Remote
 - Model 661, 20 ft.
 - Model 662, 25 ft.
 - Model 663, 50 ft.
 - Model 664, 75 ft.
- Pod to EGM Cable Model 653S

INTENDED USE/INDICATIONS FOR USE

The Maestro 4000 Cardiac Ablation System is intended for use with BSC cardiac ablation catheters in cardiac ablation procedures.

Note: Refer to the individual catheter Directions for Use for catheter compatibility to the Maestro 4000 Cardiac Ablation System. It is also important to carefully review the specific indications, contraindications, warnings, precautions and adverse events included with each catheter, prior to use of the catheter with the Maestro 4000 Cardiac Ablation System.

CONTRAINDICATIONS

There are no specific contraindications for use of the Maestro 4000 Cardiac Ablation System itself. However, users should read and understand the specific indications, contraindications, warnings, and precautions included with any cardiac ablation catheter used in conjunction with the System.

Note: The contraindications listed in the catheter Directions For Use also apply to the use of the Maestro 4000 Cardiac Ablation System. Carefully review the specific indications, contraindications, warnings, precautions, and adverse events included with each catheter, prior to use of the catheter with the Maestro 4000 Cardiac Ablation System.

ADVERSE EVENTS

Users should also read and understand the specific indications, contraindications, warnings, and precautions included with any catheter used in conjunction with the System.

Potential adverse events associated with the use of the Maestro 4000™ Cardiac Ablation System are, but not limited to, the following:

- Additional intervention required
- Arrhythmia
- Burns
- Cardiac Arrest
- Cardiac Tamponade
- Cerebral Vascular Accident (CVA)
- Complete Heart Block
- Conduction Pathway Injury
- Congestive Heart Failure
- Death
- Discomfort
- Edema
- Electrical Shock
- Embolism
- Esophagitis
- Exposure to Biohazardous Material
- Fistula
- Hematoma
- Infection
- Injury (Not Otherwise Specified) -
- Laceration
- Myocardial Infarction
- Myocardial Trauma
- Necrosis
- Nerve Injury
- Perforation
- Pericardial Effusion
- Pericarditis
- Pleural Effusion
- Prolonged Procedure
- Renal damage/failure
- Respiratory Distress/Insufficiency
- Swallowing Disorders
- Tissue Damage
- Transient Ischemic Attack (TIA)
- Vasospasm
- Vessel Occlusion
- Vessel Trauma

HOW SUPPLIED

The Maestro 4000 Cardiac Ablation System components are supplied as non-sterile in corrugated boxes. Do not use if package is opened or damaged. Do not use if labeling is incomplete or illegible.

Handling and Storage

Operating Environment

Ambient Temperature: 10 °C to 40 °C

Relative Humidity: 30% to 75%

Atmospheric Pressure: 70 kPa to 106 kPa

Transport Environment

Temperature: -29 °C to 60 °C

Relative Humidity: 30% to 85%

Atmospheric Pressure: Uncontrolled

Storage Environment

Ambient Temperature: 20 °C to 30 °C

Relative Humidity: Uncontrolled

Atmospheric Pressure: Uncontrolled

SYSTEM CONTROLS, DISPLAYS AND INDICATORS

The controls and displays for the Maestro 4000™ Cardiac Ablation System are located on the front panel of the RF Generator and are duplicated on the front panel of the optional Remote. These two control panels are identical except that the RF Generator's front panel also has an Isolated Patient Connection for the System's Pod. Since the front panels are the same, only the RF Generator front panel is illustrated in Figure 1.


Figure 1. RF Generator Front Panel

A. Power Setting Adjustment.....	6	L. Mains Power Switch	7
B. Power Display.....	6	M. System Fault Indicator	9
C. Temperature Setting Adjustment	6	N. Memory Buttons.....	7
D. Temperature Display	6	O. Alphanumeric Display	8
E. Impedance Setting Adjustment	6	P. Active/Inactive Panel	9
F. Impedance Display	6	R. Generator-Remote Button	8
G. Time Setting Adjustment.....	7	S. Control Mode Button	8
H. Time Display	7	T. Control Mode Indicators.....	9
I. Clear Counter Button.....	7	U. Average Button	9
J. View Settings Button.....	7	V. Maximum Button.....	9
K. RF Power Control Button	7	W. Isolated Patient Connector	9

Power Setting Adjustment

- A single press of a Power Setting Adjustment Button (▲/▼) will increase or decrease the RF power setting in 1-watt increments in both STANDBY and DELIVER states.
- Pressing and holding the Power Setting Adjustment Buttons (▲/▼) allows rapid scrolling of the RF power setting in STANDBY state only.
- The adjustment range varies by catheter/pod type (See Table 1 Adjustment Ranges for Power and Temperature).

Power Display

- In STANDBY state, the RF Power setting is displayed.
- In DELIVER state, the measured RF Power is displayed.
- Temperature Setting Adjustment: A single press of the Temperature Setting Adjustment Button (▲/▼) will increase or decrease the catheter tip temperature setting by 1 °C increments in both STANDBY and DELIVER states.
- Pressing and holding the Temperature Setting Adjustment Buttons (▲/▼) allows rapid scrolling of the catheter tip temperature setting in STANDBY state only.
- The adjustment range varies by catheter/pod type and the control mode selected (See Table 1 Adjustment Ranges for Power and Temperature).

Temperature Display

- In STANDBY state, the temperature setting is displayed.
- In DELIVER state:
 - Measured catheter tip temperature is displayed (between 15 °C and 95 °C).
 - Three dashes (---) will be displayed when the measured temperature is below 15 °C and above 95 °C.

Impedance Upper Limit Adjustment

- A single press of the Impedance Setting Adjustment Button (▲/▼) will increase or decrease the impedance upper limit setting by 1 Ω increments in both STANDBY and DELIVER states.
- Pressing and holding the Impedance Setting Adjustment Buttons (▲/▼) allows rapid scrolling of the impedance upper limit setting in STANDBY state only.
- The adjustment range for the impedance upper limit is 150 Ω to 300 Ω.

Impedance Display

- In STANDBY state, the impedance upper limit setting is displayed.
- In DELIVER state, the measured impedance is displayed.
- The impedance display will flash “HI” and RF power will shut off if the measured impedance rises above the upper limit during RF delivery.
- The impedance display will flash “LO” and RF power will shut off if the measured impedance falls below the low impedance limit for the catheter during RF delivery.
 - The low impedance limit for high-power catheters is 25 Ω.
 - The low impedance limit for non-high-power catheters is 50 Ω.

Time Setting Adjustment

- The Time Setting Adjustment Button increases or decreases the maximum duration for RF power delivery in one-second increments during STANDBY state only.
- A single press of the Time Setting Adjustment Button (▲/▼) will increase or decrease the time by 1 second.
- Pressing and holding the Time Adjust Buttons (▲/▼) allows rapid scrolling of the time setting.
- The range for the maximum duration setting varies by catheter type.
 - The range for open-irrigated catheters is 0 – 999 seconds.
 - The range for all other catheter types is 0 – 120 seconds.

Note: The maximum duration must be set to at least 4 seconds to start RF delivery.

Time Display

- In STANDBY state, the setting for maximum duration of RF power is displayed.
- In DELIVER state, the display changes to zero seconds at the onset of RF delivery, and then counts up to the maximum duration (unless delivery is manually terminated).

Clear Counter Button

- The Clear Counter Button resets the count of valid RF deliveries to zero and resets the total RF delivery time to zero.
- This button is inactive during RF delivery.

View Settings Button

- During RF delivery, pressing and holding the View Settings Button temporarily displays the settings for Power, Temperature, Impedance, and Time instead of the measured values.

RF Power Control Button

- Press to start RF delivery while in STANDBY state.
- Press to stop RF delivery while in DELIVER state.
- Press to clear a DIAGNOSTIC message and exit DIAGNOSTIC state.

Mains Power Switch

- The Mains Power Switch turns the unit ON and OFF.

Memory Buttons

- Memory Buttons - allow RF delivery settings to be stored and recalled.
 - Select the desired Control Mode, Power, Temperature, Impedance, Time, and status of the Maximum and Average Buttons (ON or OFF)
 - To store these values, press and hold a Memory button for at least two seconds until audible tone confirms that settings have been saved.
 - To recall the settings, press and release the appropriate Memory button.
- The Memory buttons are inactive during RF power delivery.

Alphanumeric Display

In STANDBY state

- Measured temperature is displayed between 15 °C – 95 °C (Outside this range, the display reads “D06-TEMP OUT OF RANGE [15C-95C]”).
- Measured impedance is displayed between 0-300 Ω (Over 300 Ω, the display reads “HI”).
- Operational messages (M06 to M10) alert the operator to special conditions (without entering DIAGNOSTIC state). See OPERATIONAL MESSAGES Section for additional information.

In DELIVER state

- Operational messages (M02 to M05 and M11) alert the operator to special conditions affecting the delivery of RF (without entering DIAGNOSTIC state). See OPERATIONAL MESSAGES Section for additional information.

In Both STANDBY and DELIVER states

- The count of valid RF deliveries and total RF delivery time are displayed.
 - A valid delivery is defined as an RF application greater than 4 seconds.
 - After 999 deliveries, the counter will roll over to zero.
 - The total RF delivery time is displayed in MM:SS format.
- The counter should be cleared between patients.
 - Upon initial power ON, the counter clears automatically
 - Reconnecting a catheter within 20 minutes will not reset the counter. If a catheter is reconnected after 20 minutes, the counter will be reset.
- Average and maximum values for power, temperature, and impedance are displayed if the feature is selected.
- Diagnostic codes (D01 – D12) alert the operator to a diagnostic condition. See DIAGNOSTIC MESSAGES Section for additional information.

Generator-Remote Button

- The Generator-Remote buttons toggle the Generator and Remote Control panels between active and inactive states.
- All controls are functional on the active panel.
- The following controls remain fully or partially functional on the “inactive” panel:
 - The RF Power Control button can be used to stop RF power.
 - The Generator-Remote and View Settings buttons are fully functional.
- The Generator-Remote Button is inactive during RF power delivery.

Control-Mode Button

- Pressing the Control-Mode Button toggles the unit between the Power-Control and Temperature-Control modes of operation.
- The button is inactive during RF power delivery.

Average Button

- When the Average Button is pressed in STANDBY state, the average values of power, temperature and impedance for the previous ablation are shown in the alphanumeric display.
- When the Average Button is pressed in DELIVER state, the average values of power, temperature and impedance for the current ablation are shown in the alphanumeric display.

Maximum Button

- When the Maximum Button is pressed in STANDBY state, the maximum values of power, temperature and impedance for the previous ablation are shown in the alphanumeric display.
- When the Maximum Button is pressed in DELIVER state, the maximum values of power, temperature and impedance for the current ablation are displayed in the alphanumeric display.

Isolated Patient Connector

Provides for connection of the Pod to the RF Generator.

Front Panel Indicators

The front panel indicators include the following:

System Fault	Red	Indicates a System Fault has occurred.
Temperature Mode	Green	Indicates System is in Temperature-Control Mode
Power Mode	Green	Indicates System is in Power-Control Mode
RF Power Control	Flashing Amber	Indicates System is in STANDBY state
RF Power Control	Solid Amber	Indicates System is in RF DELIVER state
Memory	Green	Indicates settings are loaded from memory
Panel Active	Green	Indicates the Control Panel is active
Panel Inactive	Yellow	Indicates the Control Panel is inactive
Average	Green	Indicates the Average feature is active
Maximum	Green	Indicates the Maximum feature is active

Remote (Optional Accessory)

The Remote allows an operator access to all controls, displays, and indicators while at a distance from the RF Generator. The Remote can be connected to the RF Generator with cables of various lengths. All controls, displays, and indicators on the Remote have the same appearance and function as those on the RF Generator.

Foot Switch (Optional Accessory)

A Foot Switch with a 10 foot cable can be connected to the RF Generator. The Foot Switch provides hands-free control to start/stop RF delivery. The Foot Switch functions similarly to the RF Power Control button except that the operator must continuously hold the Foot Switch down for RF power to be delivered. RF power delivery is immediately terminated when the operator's foot is lifted off the Foot Switch.

Pressing and releasing the Foot Switch can also be used to clear diagnostic messages and exit diagnostic state.

SYSTEM INSTALLATION

This section provides installation instructions for the Maestro 4000™ Cardiac Ablation System.

- **WARNING:** No modification of this equipment is allowed other than the specified fuses, as this may result in electrical shock and/or other unexpected consequences.
- **WARNING:** To avoid risk of electric shock, this equipment must only be connected to a supply mains with protective earth.
- **WARNING:** A Hospital Grade power cord (such as those supplied with the Maestro) must be used to connect the RF Generator or Remote Mains Power Inlet to an AC wall outlet designated as “Hospital Grade” or “Hospital Only.”
- **WARNING:** Do not use an extension cord with the System and do not connect the System’s power supply cord to an additional multiple portable socket.
- **WARNING:** Equipment connected to the analog and digital interfaces of the System must be certified to the respective IEC standards (i.e. IEC 950 for data processing equipment and IEC 60601-1 for medical equipment). Furthermore, all configurations shall comply with the system standard IEC 60601-1-1 (or 60601-1 3rd edition). Any user who connects additional equipment to the signal input ports or signal output ports configures a medical system, and is therefore responsible for the compliance of that system with the requirements of the system standard IEC 60601-1-1 (or 60601-1 3rd edition).
- **WARNING:** Electromagnetic interference (EMI) produced by the RF Generator during delivery of RF power may adversely affect the performance of other equipment. The user is encouraged to try to correct the interference by one or more of the following measures: Reorient or relocate the receiving device, increase the separation between the equipment, connect the equipment into an outlet on a circuit different from that to which the other device(s) are connected. Consult your local BSC field service technician for help. The use of accessories, transducers and cables other than those specified may result in increased electromagnetic emissions or decreased electromagnetic immunity of the System. EMC installation and service information is provided in the Maestro 4000 Cardiac Ablation System Safety Specifications section.
- **WARNING:** Components of the Maestro 4000 Cardiac Ablation System should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the Maestro 4000 Cardiac Ablation System should be observed to verify normal operation in the configuration in which it will be used.
- **WARNING:** When physiological monitoring equipment is used on the same patient, any monitoring electrodes should be placed as far as possible from the ablation electrodes. Needle monitoring electrodes are not recommended. Monitoring systems incorporating high frequency current-limiting devices are recommended.
- **WARNING:** The cables connecting the catheter to the System should be positioned in such a way that contact with the patient or other leads is avoided.
- **WARNING:** During power delivery, the patient should not come into contact with metal parts which are electrically connected to earth ground or which have an appreciable capacitance to earth (such as operating table supports, etc.). The use of antistatic sheeting is recommended.
- **WARNING:** Skin-to-skin contact (for example between the arms and body of the patient) should be avoided, for example by insertion of dry gauze.

- **WARNING:** There is a possibility of skin burns to the patient during ablation therapy. The use of dispersive pads, which meet or exceed IEC 60601-2-2 requirements, is required. The entire area of the dispersive pad should be reliably attached to the patient's body as close to operating field as possible. Refer to the manufacturer's directions for use for proper application. When using BSC high-power catheters, it is required that two dispersive pads be used.
- **WARNING:** The Pod cable is permanently attached to the Pod. Do not attempt to loosen or remove this cable from the Pod. Do not attempt to rotate or twist the connector itself.
- **Caution:** Pins of connectors identified with the ESD warning symbol should not be touched and connections should not be made to these connectors unless ESD precautionary procedures are used. It is recommended that all staff involved in the assembly and/or installation of the Maestro 4000™ Cardiac Ablation System receive an explanation of the ESD warning symbol and training in ESD precautionary procedures.
- **Caution:** Portable and mobile RF communications equipment can affect the Maestro 4000 Cardiac Ablation System. It is advised not to use this equipment in proximity to the Maestro 4000 Cardiac Ablation System.
- **Caution:** Adequate filtering must be used to allow continuous monitoring of the surface electrocardiogram (ECG) during RF power applications.
- **Note:** The RF Generator, Pod, and if used, the Remote Control and Foot Switch, are intended for installation within the Patient Environment.
 Note: The Patient Environment is described per IEC 60601-1 as any volume in which intentional or unintentional contact can occur between Patients and parts of the System, or between Patient and other persons touching parts of the System.

RF Generator Rear Panel

This section describes the connections and controls on the rear panel of the cardiac RF Generator shown in Figure 2.


Figure 2. RF Generator Rear Panel

Mains Power Inlet

The Mains Power Inlet connects to an AC wall outlet designated as “Hospital Grade” or “Hospital Only” using a “Hospital Grade” power cord.

Potential Equalization Terminal

- This terminal can be used to verify compliance with 60601-1.
- This terminal is used for connecting dedicated ground leads.

Foot Switch Connector

This connector is used for the optional Foot Switch.

Peristaltic Pump Connector

This connector can be used to communicate with the optional MetriQ™ Pump.

Remote Connector

This connector can be used to communicate with the optional Remote.

EP Recording System Connector

This connector allows the RF Generator to communicate the ablation parameters to the EP Recording System.

Speaker Volume Knob

This knob adjusts the volume of the RF Generator audio speaker.

Remote Rear Panel

This section describes the connections and controls on the rear panel of the Remote shown in Figure 3.


Figure 3. Remote Rear Panel

Mains Power Inlet

The Mains Power Inlet must be connected to an AC wall outlet designated as “Hospital Grade” or “Hospital Only” using a “Hospital Grade” power cord.

Potential Equalization Terminal

- This terminal can be used to verify compliance with 60601-1.
- This terminal is used for connecting dedicated ground leads.

Remote Control Connector

Provides an interface between the Remote and the RF Generator.

Speaker Volume Knob

This knob adjusts the volume of the Remote audio speaker.

Pod Connections

This section describes the connectors on the Pod shown in Figure 4.


Figure 4. Pod Connections

Standard/High-Power Catheter Connector

One nine-pin receptacle labeled “STD/XP” is provided for connection to catheters.

Dispersive Pad Connectors

- Two male two-pin receptacles are provided for connection to dispersive pads (also known as the dispersive indifferent pads, ground pads, or return pads).
- Dispersive pads complete the electrical circuit, allowing RF power to be delivered to the patient.
- Connect (2) dispersive pads to the Pod when using high-power catheters.
- A single dispersive pad is sufficient when using non-high-power catheters.
- The RF Generator will cease RF output to the patient if the current through a dispersive pad exceeds 1.1 amps.
- Read the manufacturer’s instructions before installing dispersive pads.

EP Recorder Connector

Connects the Pod to the EGM pins of the EP Recording System using a BSC Model 653S cable.

Pod-to-Generator Cable

The Pod-to-Generator cable is used to connect the Pod to the “Isolated Patient Connector” on the front panel of the RF Generator.

Compatible Catheters

- **WARNING:** Refer to the Directions for Use (DFU) provided with individual catheters to determine if the catheters are compatible with the Maestro 4000™ Cardiac Ablation System. Carefully review the specific indications, contraindications, warnings, precautions and adverse events included in a catheter's DFU prior to using the catheter with the Maestro 4000 Cardiac Ablation System.

Catheters tested for compatibility for use with the Maestro 4000 Cardiac Ablation System include the BSC Models listed below.

Note: some catheters may not be available in all geographic areas.

Standard Catheters:

- Blazer II™
- Blazer II HTD™
- Blazer Prime™ HTD

High-Power Catheters:

- Blazer II XP™
- Blazer Prime™ XP
- IntellaTip MiFi XP

Closed-Irrigation Catheters:

- Chilli II Cooled Ablation Catheter

Open-Irrigated Catheters:

- Blazer Open-Irrigated Ablation Catheter

The RF Generator recognizes the catheter/pod combination and automatically adjusts the ranges of power and temperature settings available for the selected control mode.

Table 1 Adjustment Ranges for Power and Temperature

Catheter Type	Pod Type	Temperature Range Temp-Control Mode	Temperature Range Power-Control Mode	Power Range
Standard	100W	30 – 90 °C	30 – 95 °C	0 - 50 W
Standard	150W	30 – 90 °C	30 – 95 °C	0 - 100 W
High-Power	100W	30 – 80 °C	30 – 85 °C	0 - 100 W
High-Power	150W	30 – 90 °C	30 – 95 °C	0 - 150 W
Closed-Irrigation	100W	30 – 90 °C	30 – 95 °C	0 - 50 W
Closed-Irrigation	150W	30 – 90 °C	30 – 95 °C	0 - 50 W
Open-Irrigated	100W	30 – 50 °C	30 – 50 °C	0 - 50 W
Open-Irrigated	150W	30 – 50 °C	30 – 50 °C	0 - 50 W

OPERATION

- **WARNING:** Pacemakers, implantable cardioverter/defibrillators, and leads can be adversely affected by radiofrequency energy. It is important to refer to the manufacturer's instruction for use prior to performing ablation procedures.
- **WARNING:** The impedance display of the RF Generator should be continuously monitored during RF power delivery. If a sudden rise in impedance is noted, RF power delivery should be discontinued.
- **WARNING:** Failure of the RF Generator could result in an unintended increase of output power.
- **WARNING:** The risk of igniting flammable gases or other materials is inherent in the application of RF power and precautions must be taken accordingly. Some materials, for example cotton, wool and gauze, when saturated with oxygen may be ignited by sparks produced in normal use of the HF surgical equipment.

Note: Boston Scientific relies on the physician to determine, assess and communicate to the individual patient all foreseeable risks of the cardiac ablation procedure.

This section provides basic operation instructions for the Maestro 4000™ Cardiac Ablation System.

Start-Up

Turn the RF Generator ON by toggling the Mains Power Switch on the front panel to the "I" position.

If a Remote is connected to the RF Generator:

- Turn the Remote ON by toggling its Mains Power Switch to the "I" position.
- Press the Generator Remote button to activate the Remote Control Panel.

The RF Generator and Remote (if applicable) will initiate a self-test.

- "SELF TEST" appears on the alphanumeric display, a continuous audio tone is generated for approximately two seconds, and all front panel displays/indicators are lit.
- The Pod and a Catheter must be connected to complete the self-test.
- If no System malfunction is detected, the RF Control Button will flash to indicate the RF Generator has entered STANDBY state.
- If a System malfunction is detected, the RF Generator or Remote (if applicable) will not operate. The red SYSTEM FAULT indicator will be lit. The displays for Power, Temp, Impedance and Time will all read 888.

Standby State

STANDBY state initiates automatically after System start-up and whenever Deliver state is terminated. Power, Temperature, Impedance, and Time settings may be viewed and adjusted in this state.

When STANDBY state is initiated after System start-up:

- The default settings shown in Table 2 will appear in the Power, Temperature, Impedance, and Time displays.

Table 2. Default Setting for Power, Temperature, Impedance, Time

Catheter Type	Power	Temperature	Impedance	Time
Standard, High-Power, Closed-Irrigation	0 Watts	30 °C	300Ω	0 Sec
Open-Irrigated	0 Watts	50 °C	250Ω	0 Sec

- Measured impedance and temperature values appear in the alphanumeric display.
- Ablation count and total time (also in the alphanumeric display) will be reset to zero.

Select RF Delivery Settings

Use one of the methods below to select RF Delivery settings for a procedure.

- Recall RF Delivery settings stored with the Memory Buttons.
- Manually set the RF Delivery settings as described below.

Control-Mode Selection

- **WARNING:** With non-irrigated catheters, start with a low temperature setting in temperature-control mode, and then slowly titrate the temperature setting up to achieve the desired effect. With irrigated catheters, select power-control mode with a temperature setting intended to ensure that adequate fluid is irrigating the tip of the catheter.

The mode of operation can be toggled between Power-Control Mode and Temperature-Control Mode by pressing the Control Mode Button while the unit is in STANDBY mode (The control mode button is inactive during RF power delivery).

The control-mode selection determines the effect of the Power and Temperature settings.

- In Power-Control Mode:
 - The amount of RF power delivered will equal the power setting unless the measured temperature exceeds the temperature setting.
 - The temperature setting is a limit for catheter tip temperature. If the limit is exceeded, power will be reduced or halted to lower the measured temperature below this limit.

Note: When power is reduced or halted, an operational message or diagnostic code will appear in the front panel display. See OPERATIONAL MESSAGES Section, or DIAGNOSTIC MESSAGES Section for additional information.

- In Temperature-Control Mode:
 - The power setting is a limit. The RF Generator is allowed to adjust RF power delivery up to this limit to achieve the target temperature.
 - The temperature setting is the target temperature.

Note: The target temperature may not be reached if the power setting (limit) is reached.

The default control modes for BSC Catheters are listed in Table 3. When a catheter is initially connected, the RF Generator automatically switches to the default control mode for that catheter. However, the user may switch control modes by pressing the control-mode button.

Table 3. Catheter Default Control Modes

Catheter Types	Default Control Mode
High-Power and Standard Power	Temperature Control
Open-Irrigated and Closed-Irrigated	Power Control

Power Setting

- **WARNING:** To minimize the potential for thrombus formation, inadvertent damage to cardiac tissues, and collateral damage to adjacent tissue not intended for ablation, begin by using a low power setting and gradually increase the power output if necessary, especially in areas where low blood flow and correspondingly low convective cooling may be present.

Temperature Setting

- **WARNING:** To minimize the potential for thrombus formation, inadvertent damage to cardiac tissues, and collateral damage to adjacent tissue not intended for ablation, begin by using a low temperature setting and gradually increase the power output if necessary, especially in areas where high blood flow and correspondingly high convective cooling may be present.
- **WARNING:** The displayed temperature is not the maximum tissue temperature. The measured temperature may be influenced by the degree of tissue contact and variations in blood flow. The difference between the max tissue temperature and the displayed temperature increases when using irrigated catheters, and the correlation with lesion formation is greatly reduced. The temperature displayed with irrigated catheters is not intended to guide lesion formation. The temperature displayed with irrigated catheters is intended to verify that the tip is being sufficiently irrigated. This behavior will minimize the potential for thrombus formation and/or inadvertent damage to nearby tissue.

Impedance Setting

The impedance upper limit may be adjusted between 150 Ω to 300 Ω .

Time Setting

Set the timer for the maximum duration of RF Power delivery.

Selecting Average and Maximum Display

Press the “Maximum” and/or “Average” Buttons, if those values are to be displayed.

RF Power Delivery with Open-Irrigated BSC Catheters

Use of the RF Generator with BSC open-irrigated catheters requires proper set-up and preparation of the MetriQ™ Pump and the MetriQ Irrigation Tubing Set. The MetriQ Pump Operator’s Manual provides complete instructions for set-up and operation of the MetriQ Pump and MetriQ Irrigation Tubing Set .

Starting RF Power Delivery:

- Press the RF Power Control button once, or hold down the Foot Switch pedal.
- The MetriQ Pump will switch from standby flow to ablation flow.
- Delivery of RF Power will start when the Pre-RF Delay is complete.
- The RF Power Control button will illuminate and remain lit until RF power delivery stops.
- RF Power will stop automatically when the setting for maximum duration of RF Power Delivery is reached.
- The Pump will switch to standby flow when the Post-RF Delay is complete.

Stopping RF Power Delivery:

- To stop delivery of RF Power before reaching the setting for maximum duration, press the RF Power Control Button, or release the Foot Switch pedal.
- The RF Power Control button will flash to indicate the RF Generator is in STANDBY state.
- The Pump will switch to standby flow when the Post-RF Delay is complete.

RF Power Delivery with all other Catheters

Starting RF Power Delivery:

- Press the RF Power Control button once, or hold down the Foot Switch pedal.
- Delivery of RF Power will start immediately.
- The RF Power Control button will illuminate and remain lit until RF power delivery stops.
- RF Power will stop automatically when the setting for maximum duration of RF Power Delivery is reached.

Stopping RF Power Delivery:

- To stop delivery of RF Power before reaching the setting for maximum duration, press the RF Power Control Button, or release the Foot Switch pedal.
- The RF Power Control button will flash to indicate that the RF Generator is in STANDBY state.

SERVICE AND MAINTENANCE

Preventative Inspection

- **WARNING:** Damage such as frayed cords or cables and cracks or dents on the equipment may result in electrical shock.
- **Caution:** BSC recommends that the RF Generator and Remote be powered off at the end of each procedure in order to ensure that the self-test is performed before the next procedure.

During the useful life of the equipment, maintain close watch for damage such as frayed cords or cables and cracks or dents on the equipment. If damage is identified, take the equipment out of service and contact Boston Scientific Corporation for service requirements.

Functional Self-Test

System does not require periodic preventative maintenance and calibration. Upon start up, the unit performs self-test per OPERATION Start-up Section of this manual. For further details refer to that section.

System Servicing

None of the System components are user-serviceable. Contact Boston Scientific Corporation for all service requirements.

Cleaning/Disinfecting

- **WARNING:** Do not immerse the Generator, the Remote, or accessories in any liquid.
- **WARNING:** Use of non-flammable cleaning and disinfection agents is recommended. If used, flammable agents or solvents should be allowed to evaporate before high-frequency surgery. There is a risk of pooling of flammable solutions under the PATIENT or in body depressions such as the umbilicus, and in body cavities such as the vagina. Any fluid pooled in these areas should be mopped up before HF SURGICAL EQUIPMENT is used.
- **Caution:** The Generator, Pod, Remote, Footswitch, power cords, and communication cables are not intended to be sterilized and should remain outside of the sterile field.

The outer surfaces of the Generator, the Pod, the Remote and their accessories may be wiped clean with a mild soapy solution. Isopropyl alcohol may be used to clean the outer surfaces. Avoid caustic or abrasive cleaners.

ESD Training and Precautionary Procedures

Prior to assembly, installation, or interconnection of the Maestro 4000™ Cardiac Ablation System, it is recommended that any staff (i.e. clinical/biomedical engineers and health care staff) that could touch connectors identified with the ESD warning symbol undergo ESD training. At minimum, ESD training should include an introduction to the physics of electrostatic charge, the voltage levels that can occur in normal practice, and the damage that can be done to electronic components if they are touched by an operator who is electrostatically charged. Further, an explanation should be given of methods to prevent build-up of electrostatic charge, and how and why to discharge one's body to earth or to the frame of the equipment or System, or bond oneself by means of a wrist strap to the equipment or System or to earth prior to making a connection. Finally, staff must be made aware that accessible pins of connectors identified with the ESD warning symbol should not be touched with the fingers or with a hand-held tool unless proper precautionary procedures have been followed. ESD precautionary procedures should include:

- Methods to prevent build-up of electrostatic charge (e.g. air conditioning, humidification, conductive floor coverings, non-synthetic clothing);
- Discharging one's body to the frame of the equipment or System or to earth or a large metal object;
- Bonding oneself by means of a wrist strap to the equipment or System or to earth.

End of Useful Life

When the equipment reaches the end of its useful life, dispose of the Generator, Remote and all accessories in accordance with hospital, administrative and, local government policy. Contact your BSC representative or BSC field service engineer (1.800.949.6708 in the US) prior to disposal.

The Maestro 4000 Controller components are expected to have a useful life up to 7 years.

PRODUCT SPECIFICATIONS


General Specifications

Description	Specification
Power Specifications	
Line Power	
• RF Generator	100-120/220-240V _{AC} , 50/60Hz, 300VA
• Remote	100-120/220-240V _{AC} , 50/60Hz, 25VA
Current Rating	
• RF Generator	4 A @120V _{AC}
• Remote	1 A @120V _{AC}
Fuses	
• RF Generator	T4AL250V
• Remote	T1AL250V
Power Cord	
Length	10 feet (3.0 meters)
Foot Switch (Cable)	
Length	10 feet (3.0 meters)

Description	Specification
Pod to RF Generator Cable Length Connector	15 feet (4.6 meters) 14-pin Quick Connector
Pod Dimensions (excluding cable) <ul style="list-style-type: none"> • Height • Width • Depth • Weight Pod Connectors <ul style="list-style-type: none"> • Recorder • Catheter • Indifferent Electrodes Pod Recorder Filters <ul style="list-style-type: none"> • Low Pass Filters • Low Frequency Cutoff 	2.3 in. (5.8 cm) 6.8 in. (17.3 cm) 6.2 in. (15.7 cm) 2.2 lb. (1.0 kg) Quick Connect Connector 9-pin Quick Connect Connector - Type CF-Defibrillation-Proof Standard male 2-pin for commercial pads - Type CF-Defibrillation-Proof Referenced to the Indifferent Electrode –3 dB at 5 kHz ± 1 kHz
High-Current Shutoff Mechanism	RF output ceases to the patient if more than 1.1 Amp flows in either Dispersive Pad
RF Generator Power Output @ 460 kHz	<ul style="list-style-type: none"> • Maximum Output Power: <ul style="list-style-type: none"> • 100 W for 100-W Pod • 150 W for 150-W Pod • Maximum Output Voltage: <ul style="list-style-type: none"> • 132 Vrms (187 Vpk) with 100-W pod • 162 Vrms (229 Vpk) with 150-W pod • Maximum Output Current (under normal operation): 1 A per dispersive pad • Maximum RF power setting based on catheter type used and Control Mode • Use catheters that have a rated voltage higher than the RF Generator maximum output voltage rating.
Impedance	In STANDBY Mode: The impedance is measured in the 0-300 Ω range In DELIVER Mode: High-power catheters display measured impedance in the 25-300 Ω range and display “LO” or “HI” outside that range. Non-high-power catheters display measured impedance in the 50-300 Ω range and displays “LO” or “HI” outside that range.

Description	Specification
<p>Temperature</p> <p>In Temperature-Control Mode</p> <p>In Power-Control Mode</p> <p>In either control mode:</p>	<p>The catheter tip temperature (temperature set point) can be selected within the ranges listed in Table 1.</p> <p>The catheter tip temperature (upper limit for continued RF power delivery at the power set point) can be selected (within the ranges listed in Table 1).</p> <p>The RF Generator measures and displays temperature in the 15 °C to 95 °C range and displays “LO” or “HI” outside that range.</p>
<p>RF Generator Dimensions</p> <ul style="list-style-type: none"> • Width • Height • Depth • Weight 	<p>13.0 inches (33.1 cm)</p> <p>7.3 inches (18.6 cm)</p> <p>16.5 inches deep (41.9 cm)</p> <p>22 lbs. (10 kilograms)</p>
<p>Time</p>	<p>For open-irrigated catheters, the time can be set between 0 – 999 seconds in increments of 1 second.</p> <p>For all other catheters, the time can be set between 0 – 120 seconds in increments of 1 second.</p>
<p>Counter</p>	<p>0 to 999 RF power deliveries</p>
<p>Generator Sound Output</p> <ul style="list-style-type: none"> • Power Up – Sign On Tone • Keyboard Key Click • RF Delivery • Three Tone Beep 	<p>500 Hz – 1.5 kHz – 2.5 kHz Square Wave</p> <p>1 kHz Square Wave</p> <p>500 Hz Square Wave</p> <p>300 Hz Square Wave</p>
<p>Remote Dimensions</p> <ul style="list-style-type: none"> • Height • Width • Depth • Weight 	<p>13.0 inches (33.1 cm)</p> <p>7.3 inches (18.6 cm)</p> <p>7.5 inches deep (19.1 cm)</p> <p>11 lb. (5 kg)</p>

Maestro 4000™ Cardiac Ablation System Safety Specifications

Device Description
Class I, Defibrillation proof Type CF Equipment, IPX0, not AP/APG
Mode of Operation: Continuous
EMC Emissions and Susceptibility: The Maestro 4000 Cardiac Ablation System has been tested and found to comply with the limits for medical devices to the IEC 60601-1-2:2007. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. This System generates, uses, and can radiate radiofrequency energy and, if not installed and used in accordance with the instruction given below, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation.
TUV Rheinland of North America Certified.

Electrical Isolation
Leakage current conforms to IEC 60601-1
Dielectric withstand conforms to IEC 60601-1

Guidance and manufacturer's declaration – electromagnetic emissions		
The Maestro 4000 Cardiac Ablation System is intended for use in the electromagnetic environment specified below. The customer or the user of the Maestro 4000 Cardiac Ablation System should assure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The Maestro 4000 Cardiac Ablation System must emit electromagnetic energy in order to perform its intended function. Nearby electronic equipment may be affected.
RF emissions CISPR 11	Class A	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	The Maestro 4000 Cardiac Ablation System is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Guidance and manufacturer's declaration – electromagnetic immunity			
The Maestro 4000 Cardiac Ablation System is intended for use in the electromagnetic environment specified below. The customer or the user of the Maestro 4000 Cardiac Ablation System should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.

Guidance and manufacturer's declaration – electromagnetic immunity			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Maestro 4000™ Controller requires continued operation during electrical disturbances on the power mains, it is recommended that the MetriQ™ Pump be connected to an uninterruptible power supply with surge and electrical fast transient-suppression filtering/device embedded or connected in series.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Maestro 4000 Controller requires continued operation during electrical disturbances on the power mains, it is recommended that the MetriQ Pump be connected to an uninterruptible power supply with surge and electrical fast transient-suppression filtering/device embedded or connected in series.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycles 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 seconds	<5% U_T (>95% dip in U_T) for 0.5 cycles 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Maestro 4000 RF Generator requires continued operation during electrical disturbances on the power mains, it is recommended that the MetriQ Pump be connected to an uninterruptible power supply with surge and electrical fast transient-suppression filtering/device embedded or connected in series.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of typical location in a typical commercial or hospital environment.
NOTE: U_T is the AC mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration – electromagnetic immunity

The Maestro 4000™ Cardiac Ablation System is intended for use in the electromagnetic environment specified below. The customer or the user of the Maestro 4000 Cardiac Ablation System should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 Vrms 150 kHz to 80 MHz 3 V/m 80 MHz to 2.5 GHz	3 Vrms 3 V/m	<p>Portable and mobile RF communications equipment should be used no closer to any part of the Maestro 4000 Cardiac Ablation System, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance:</p> $d = 1.17\sqrt{P} \quad 150 \text{ kHz to } 80 \text{ MHz}$ $d = 1.17\sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2.33\sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Maestro 4000 Cardiac Ablation System or any of its components are used exceeds the applicable RF compliance level above, the Maestro 4000 Cardiac Ablation System should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating components or the entire Maestro 4000 Cardiac Ablation System.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Guidance and Manufacturer's Declaration – electromagnetic immunity			
The Maestro 4000™ Cardiac Ablation System is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Maestro 4000 Cardiac Ablation System can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Maestro 4000 Cardiac Ablation System as recommended below, according to the maximum output power of the communications equipment.			
Rated Maximum Output Power of Transmitter W	Separation Distance According to Frequency of Transmitter M		
	150 kHz to 80 MHz $d = 1.17\sqrt{P}$	80 MHz to 800 MHz $d = 1.17\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.33\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.70	3.70	7.37
100	11.7	11.7	23.3

For transmitters rated at a maximum output power not listed above, the recommended separation distance *d* in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where *P* is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Power Delivery

Refer to Table 1 for maximum power settings. In the event of unusually high or low tissue impedance, the System limits the maximum power. Figure 6 through Figure 9 depict the maximum output power of the System (as a function of tissue impedance) as measured at the output of the catheter for each catheter-pod configuration. Figure 10 depicts the nominal, minimum, and maximum output power of the Generator as a function of the control setting at nominal impedance of 100 Ω.

Note: Due to the inductive and capacitive loads presented by each Catheter configuration, the power delivered at the Catheter tip may be less than the power on the System display.


Figure 6. Power Output for High-Power Catheters with 150-W Pods


Figure 7. Power Output for OI/Chilli II Catheters with 100-W and 150-W Pods and for Standard Catheters with 100-W Pods


Figure 8. Power Output for Standard Catheters with 150-W Pods.


Figure 9. Power Output for High-Power Catheters with 100-W Pods


Figure 10. Maestro 4000™ Controller (RF Generator) Power Output Versus Power Setting

Device Labels

Remote Control Output EP Recording System Foot Switch Speaker Volume

VA	POWER INPUT	300 VA
Hz	FREQUENCY	50/60 Hz
~	AC VOLTAGE	100-120 / 220-240V

2X Fuse T1AL 250V

W	RATED OUTPUT POWER	50, 100, 150
Ω	RATED LOAD	100
kHz	OPERATING FREQUENCY	460

Electrostatic Sensitive Device

Legal Manufacturer: Boston Scientific Corporation, One Boston Scientific Place, Natick, MA 01760-1537, USA. USA Customer Service 888-272-1001.

Made In USA: 560 Cottonwood Drive, Milpitas, CA 95035, USA.

CE 0344

EU Authorized Representative: Boston Scientific International S.A., 55 avenue des Champs Pierreux, TSA 51101, 92729 NANTERRE CEDEX, FRANCE.

Non-Ionizing Electromagnetic Radiation

Equipotentiality

CAUTION Attention: Consult ACCOMPANYING DOCUMENTS.

Follow Instructions For Use

871-0186-000 REV A 90818618-01A

Boston Scientific

Maestro 4000™ Controller

RF Cardiac Ablation Controller

UPN Product No. M00440000

REF Catalog No. 4000

Date of Manufacture YYYY-MM-DD

SN Serial Number 12345678

VA	POWER INPUT	25 VA
Hz	FREQUENCY	50/60 Hz
~	AC VOLTAGE	100-120 / 220-240V

Electrostatic Sensitive Device Separate Collection

Legal Manufacturer: Boston Scientific Corporation, One Boston Scientific Place, Natick, MA 01760-1537, USA. USA Customer Service 888-272-1001.

Made In USA: 560 Cottonwood Drive, Milpitas, CA 95035, USA.

CE 0344

EU Authorized Representative: Boston Scientific International S.A., 55 avenue des Champs Pierreux, TSA 51101, 92729 NANTERRE CEDEX, FRANCE.

CAUTION Attention: Consult ACCOMPANYING DOCUMENTS.

Follow Instructions For Use

871-0189-000 REV A 90818619-01A

Boston Scientific

Maestro 4000™ Remote Control

RF Cardiac Ablation Remote Control

UPN Product No. M00440200

REF Catalog No. 4020

Date of Manufacture YYYY-MM-DD

SN Serial Number 12345678

Boston Scientific

CAUTION Attention: Consult ACCOMPANYING DOCUMENTS.

Follow Instructions For Use

Maestro 4000™ Pod, 100W

RF Cardiac Ablation Pod, 100W

UPN Product No. M00440100

REF Catalog No. 4010

Date of Manufacture YYYY-MM-DD

SN Serial Number 12345678

Legal Manufacturer: Boston Scientific Corporation, One Boston Scientific Place, Natick, MA 01760-1537, USA. USA Customer Service 888-272-1001.

Made In USA: 560 Cottonwood Drive, Milpitas, CA 95035, USA.

CE 0344

EU Authorized Representative: Boston Scientific International S.A., 55 avenue des Champs Pierreux, TSA 51101, 92729 NANTERRE CEDEX, FRANCE.

CAUTION Attention: Consult ACCOMPANYING DOCUMENTS.

Follow Instructions For Use

871-0188-000 REV A 90818620-01A

Boston Scientific

CAUTION Attention: Consult ACCOMPANYING DOCUMENTS.

Follow Instructions For Use

Maestro 4000™ Pod, 150W

RF Cardiac Ablation Pod, 150W

UPN Product No. M004EPT40100

REF Catalog No. EPT4010

Date of Manufacture YYYY-MM-DD

SN Serial Number 12345678

Legal Manufacturer: Boston Scientific Corporation, One Boston Scientific Place, Natick, MA 01760-1537, USA. USA Customer Service 888-272-1001.

Made In USA: 560 Cottonwood Drive, Milpitas, CA 95035, USA.

CE 0344

EU Authorized Representative: Boston Scientific International S.A., 55 avenue des Champs Pierreux, TSA 51101, 92729 NANTERRE CEDEX, FRANCE.

CAUTION Attention: Consult ACCOMPANYING DOCUMENTS.

Follow Instructions For Use

871-0187-000 REV A 90818622-01A

OPERATIONAL MESSAGES

Table 4. Operational Messages

Code	Display	Description	Action
M02	POWER LIMITED, POWER SETTING	This message is displayed as long as the power setting is reached in temperature-control mode.	Reduce temperature setting or increase power setting, as appropriate.
M03	POWER LIMITED, LOW IMPEDANCE	This message is displayed as long as the power is limited due to reaching the current limit due to low impedance.	See power curves in RF delivery section of this document.
M04	POWER LIMITED, HIGH IMPEDANCE	This message is displayed as long as the power is limited due to reaching the voltage limit due to high impedance.	WARNING: Apparent low power output or higher than typical impedance measurements may be indicative of faulty dispersive pad application or failure of an electrical lead. Check the application of the dispersive pad and all electrical connections before continuing or selecting higher power outputs. See power curves in RF delivery section of this document.
M05	WAIT FOR PUMP	This message is displayed during pre-RF and post-RF delay.	Wait for MetriQ™ Pump to complete pre-RF delay or post-RF delay (or adjust the delay duration, and minimum temperature drop setting on MetriQ Pump, as appropriate). See the MetriQ Pump Operator's Manual for additional information on pre-RF delay.
M06	CLOSE PUMP MENU	This message is displayed when RF is requested but the MetriQ Pump's menu is open (with an OI catheter connected and Maestro 4000™ Controller is operating with MetriQ Pump in Automatic mode).	Close the MetriQ Pump's menu prior to requesting RF delivery.
M07	SET TIME	This message is displayed when RF is requested but the time setting is still 0 seconds.	Adjust time setting appropriately.
M08	SET POWER	This message is displayed when RF is requested but the power setting is still 0 Watts.	Adjust power setting appropriately.
M09	SET TEMPERATURE	This message is displayed when RF is requested but the temperature setting is still 30 degrees.	Adjust temperature setting appropriately.

Code	Display	Description	Action
M10	RECONNECT PUMP	This message is displayed when communication is lost between the Maestro 4000™ Controller and MetriQ™ Pump.	Reconnect Pump and tighten connector thumb screws if the Pump was accidentally disconnected. Inspect generator-to-pump cable and/or connectors for damage.
M11	RF REDUCED 50%, MAX TEMP REACHED	This message is displayed when a non-OI catheter is used and the Maestro 4000 Controller measured temperature exceeds the temperature setting in power-control mode (and the power setting is automatically reduced by 50% to facilitate continued lesion formation).	Verify that the power-control mode is the desired mode of operation. Decrease the power setting and/or increase the temperature setting, as appropriate.

DIAGNOSTIC MESSAGES

Table 5. Diagnostic Messages and Corresponding Corrective Actions

Code	Display	Description	Action
D01	LOW IMPEDANCE	Measured impedance below the low impedance threshold of 25 Ω when using a high power catheter and 50 Ω for all other catheters.	Inspect generator-to-pod cable and connections and replace components which may appear damaged. If error recurs, replace Pod. If error recurs, replace the Generator.
D02	HIGH IMPEDANCE	Measured impedance above the impedance setting (150 Ω – 300 Ω).	Check application of dispersive pad and abide to all recommendations in the dispersive pad DFU/IFU.
D03	LIMIT EXCEEDED	$V > 1.1 \times V_{\max}$ or $I > 1.1 \times I_{\max}$ or $P > 1.2 \times P_{\max}$	Do not disconnect catheter, catheter cable, dispersive pad(s), or pull catheter back into sheath during RF delivery. If error recurs, restart Generator. Adding a second dispersive pad may help (especially if also experiencing M03 operational message). If error recurs, replace Pod. If error recurs, replace the Generator.
D04	TEMP ABOVE SET POINT	Measured temperature is greater than set point by at least 5 °C for more than four seconds.	Check temperature setting. If error recurs, restart the Generator. If error recurs, replace the Generator.

Code	Display	Description	Action
D05	EXCESSIVE TEMPERATURE	Measured temperature greater than 95 °C for more than one second or measured temperature exceeds setting in Power Control mode.	<p>Check temperature and power settings.</p> <p>Check control-mode setting.</p> <p>Check cables and connectors in the RF pathway for possible damage or loose connection.</p> <p>Verify that temperature in standby is appropriate and doesn't change with agitating connectors/cables in the RF pathway. If temperature changes with agitation, replace suspected components one at a time.</p> <p>If error recurs, replace catheter cable.</p> <p>If error recurs, replace Pod.</p> <p>If error recurs, replace the Generator.</p>
D06	TEMP OUT OF RANGE [15C-95C]	Measured temperature <15 °C or >95 °C	<p>This condition may have been caused by switching catheters during a procedure.</p> <p>Check cables and connectors in the RF pathway for possible damage or loose connection.</p> <p>Verify that temperature in standby is appropriate and doesn't change with agitating connectors/cables in the RF pathway. If temperature changes with agitation, replace suspected components one at a time.</p> <p>If error recurs, replace catheter cable.</p> <p>If error recurs, replace catheter.</p> <p>If error recurs, replace Pod.</p> <p>If error recurs, replace the Generator.</p>
D07	TEMPERATURE TOO HIGH	An attempt to RF delivery with a measured temperature > 43 °C or greater than the temperature setting within 30 seconds of the previous RF delivery.	<p>Check temperature setting.</p> <p>Check cables and connectors in the RF pathway for possible damage or loose connection.</p> <p>Verify that temperature in standby is appropriate and doesn't change with agitating connectors/cables in the RF pathway. If temperature changes with agitation, replace suspected components one at a time.</p> <p>If error recurs, replace catheter cable.</p> <p>If error recurs, replace Pod.</p> <p>If error recurs, replace the Generator.</p>
D08	CHECK POD	Generator detects a pod failure or an issue with the generator-to-pod connection.	<p>Check generator-to-pod connection.</p> <p>If error recurs, replace pod.</p>
D09	N/A	Invalid Code	Contact BSC

Code	Display	Description	Action
D10	INVALID CATHETER ID	Invalid catheter ID	Check catheter cable for damage. If error recurs, replace catheter cable. If error recurs, replace catheter. If error recurs, replace catheter Pod.
D11	CHECK REMOTE	There is a communications problem between the Remote and the Generator or the Remote was disconnected or connected during RF.	Do not disconnect or connect the Remote during RF. Check cable and connections between the Generator and remote, tightening the thumb screws on the connectors, or replacing component if damage is suspected.
D12	CHECK PUMP	There is a communication error with the MetriQ™ Pump or there is a diagnostic condition at the Pump.	Clear diagnostic condition at Pump. Check cable and connections between the Generator and Pump, tightening the thumb screws on the connectors, or replacing component if damage is suspected.
N/A	SYSTEM FAULT	An unrecoverable error has occurred on the Generator/Remote (as indicated by which unit displayed the "SYSTEM FAULT" LED). The Generator/Remote enters SAFE mode and will not allow any operations.	Power off the Generator/Remote and restart. If error recurs, contact BSC field service for further assistance.

LIMITED WARRANTY AND DISCLAIMER

Limited Warranties

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this System. When maintained in such conditions as specified by BSC, it will be free from defects in material and workmanship at buyer's location for 12 months from the date of deliver. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose. Handling, storage, and cleaning of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond BSC's control directly affect the System and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument.

Service is limited to exchanging faulty System components. No in-field repairs will be performed. For service, contact the BSC authorized service representative. The user must pay all freight charges for all parts returned to BSC. BSC will pay freight for shipping the repaired or replaced parts back to the user. BSC extends, to the registered user, all warranties offered by third party software upon which the System depends.

Extended warranties that include both the hardware and software can be purchased any time after the current warranty period has expired. If the System falls under the initial 12 month manufacturer's warranty, then an extended warranty can be purchased anytime during the first year (12 months) or after the manufacturer's warranty expires. Contact BSC for more information.

All catheters used with the System are for single use only. Do not reuse, reprocess, or re-sterilize. Re-use, reprocessing, or re-sterilization may compromise the structural integrity of the catheter and/or lead to device failure, which, in turn, may result in patient injury, illness, or death. Reuse, reprocessing or re-sterilization may also create a risk of contamination of the catheter and/or cause patient infection or cross-infection including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient. BSC assumes no liability with respect to single use instruments that are reused, reprocessed, or re-sterilized and makes no warranties, express or implied, including, but not limited to the warranties of merchantability or fitness for intended use with respect to such instrument.

Disclaimer and Exclusion of Other Warranties

There are no warranties of any kind that extend beyond the description of the warranties above. Boston Scientific Corporation disclaims and excludes all warranties, whether express or implied, of merchantability of fitness for a particular use or purpose.

Limitation of Liability for Damages

In any claim or lawsuit for damages arising from alleged breach of warranty, breach of contract, negligence, product liability or any other legal or equitable theory, the buyer specifically agrees that Boston Scientific Corporation shall not be liable for damages for loss of profits or revenues, loss of use of the product, loss of facilities or services, any downtime costs, or for claims of buyer's customers for any such damages. Boston Scientific Corporation's sole liability for damages shall be limited to the cost to buyer of the specified goods sold by Boston Scientific Corporation to buyer which give rise to the claim for liability. The buyer's use of this product shall be deemed acceptance of the terms and conditions of these limited warranties, exclusions, disclaimers and limitations of liability for money damages.

© 2014 Boston Scientific Corporation or its affiliates. All rights reserved.

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Legal Manufacturer
 Manufactured for:
 Boston Scientific Corporation
 One Boston Scientific Place
 Natick, MA 01760-1537
 USA
 USA Customer Service 888-272-1001

Made In USA:
 560 Cottonwood Drive
 Milpitas, CA 95035
 USA

CE 0344

Caution: Attention: Consult ACCOMPANYING DOCUMENTS.

UPN Product No. M00440000
REF Catalog No. 4000
Date of Manufacture YYYY-MM-DD
SN Serial Number 12345678

maestro 4000

VA	POWER INPUT	25 VA
Hz	FREQUENCY	50/60 Hz
~	AC VOLTAGE	100-120 / 220-240V

Legal Manufacturer
 Manufactured for:
 Boston Scientific Corporation
 One Boston Scientific Place
 Natick, MA 01760-1537
 USA
 USA Customer Service 888-272-1001

Made In USA:
 560 Cottonwood Drive
 Milpitas, CA 95035
 USA

CE 0344

Caution: Attention: Consult ACCOMPANYING DOCUMENTS.

UPN Product No. M00440200
REF Catalog No. 4020
Date of Manufacture YYYY-MM-DD
SN Serial Number 12345678

maestro 4000

Caution: Attention: Consult ACCOMPANYING DOCUMENTS.

Follow Instructions For Use

Boston Scientific

Maestro 4000™ Pod, 100W
 RF Cardiac Ablation Pod, 100W

UPN Product No. M00440100
REF Catalog No. 4010
Date of Manufacture YYYY-MM-DD
SN Serial Number 12345678

Legal Manufacturer
 Manufactured for:
 Boston Scientific Corporation
 One Boston Scientific Place
 Natick, MA 01760-1537
 USA
 USA Customer Service 888-272-1001

Made In USA:
 560 Cottonwood Drive
 Milpitas, CA 95035
 USA

CE 0344

Caution: Attention: Consult ACCOMPANYING DOCUMENTS.

UPN Product No. M004EPT40100
REF Catalog No. EPT4010
Date of Manufacture YYYY-MM-DD
SN Serial Number 12345678

Legal Manufacturer
 Manufactured for:
 Boston Scientific Corporation
 One Boston Scientific Place
 Natick, MA 01760-1537
 USA
 USA Customer Service 888-272-1001

Made In USA:
 560 Cottonwood Drive
 Milpitas, CA 95035
 USA

CE 0344

Caution: Attention: Consult ACCOMPANYING DOCUMENTS.

UPN Product No. M004EPT40100
REF Catalog No. EPT4010
Date of Manufacture YYYY-MM-DD
SN Serial Number 12345678

Legal Manufacturer
 Manufactured for:
 Boston Scientific Corporation
 One Boston Scientific Place
 Natick, MA 01760-1537
 USA
 USA Customer Service 888-272-1001

Made In USA:
 560 Cottonwood Drive
 Milpitas, CA 95035
 USA

CE 0344

Caution: Attention: Consult ACCOMPANYING DOCUMENTS.

Caution: Attention: Consult ACCOMPANYING DOCUMENTS.

Follow Instructions For Use

Boston Scientific

Maestro 4000™ Pod, 150W
 RF Cardiac Ablation Pod, 150W

UPN Product No. M004EPT40100
REF Catalog No. EPT4010
Date of Manufacture YYYY-MM-DD
SN Serial Number 12345678

Legal Manufacturer
 Manufactured for:
 Boston Scientific Corporation
 One Boston Scientific Place
 Natick, MA 01760-1537
 USA
 USA Customer Service 888-272-1001

Made In USA:
 560 Cottonwood Drive
 Milpitas, CA 95035
 USA

CE 0344

Caution: Attention: Consult ACCOMPANYING DOCUMENTS.

UPN Product No. M004EPT40100
REF Catalog No. EPT4010
Date of Manufacture YYYY-MM-DD
SN Serial Number 12345678

Legal Manufacturer
 Manufactured for:
 Boston Scientific Corporation
 One Boston Scientific Place
 Natick, MA 01760-1537
 USA
 USA Customer Service 888-272-1001

Made In USA:
 560 Cottonwood Drive
 Milpitas, CA 95035
 USA

CE 0344

Caution: Attention: Consult ACCOMPANYING DOCUMENTS.

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
 Remote Control
 EP Output EP Recording System
 Foot Switch
 Speaker Volume

VA	POWER INPUT	300 VA
Hz	FREQUENCY	50/60 Hz
~	AC VOLTAGE	100-120 / 220-240V

 2X T1AL 250V
 Electrostatic Sensitive Device
 Non-ionizing Electromagnetic Radiation
 EQUIPOTENTIALITY
 CE 0344
 Separate Collection
 Rx ONLY
 CAUTION Attention: Consult ACCOMPANYING DOCUMENTS.
 Follow Instructions For Use

Legal Manufacturer
 Manufactured for:
 Boston Scientific Corporation
 One Boston Scientific Place
 Natick, MA 01760-1537
 USA
 USA Customer Service 888-272-1001
Made In USA:
 560 Cottonwood Drive
 Milpitas, CA 95035
 USA

EU Authorized Representative
 Boston Scientific International S.A.
 55 avenue des Champs Pierreux
 TSA 51101
 92729 NANTERRE CEDEX
 FRANCE

UPN Product No. M00440000
 REF Catalog No. 4000
 Date of Manufacture YYYY-MM-DD
 SN Serial Number 12345678

Boston Scientific
Maestro 4000™ Controller
RF Cardiac Ablation Controller

871-0188-000 REV A 90818618-01A

VA	POWER INPUT	25 VA
Hz	FREQUENCY	50/60 Hz
~	AC VOLTAGE	100-120 / 220-240V

 Electrostatic Sensitive Device
 Separate Collection
 Speaker Volume
 Rx ONLY
 NON-STERILE
 CAUTION Attention: Consult ACCOMPANYING DOCUMENTS.
 Follow Instructions For Use


Legal Manufacturer
 Manufactured for:
 Boston Scientific Corporation
 One Boston Scientific Place
 Natick, MA 01760-1537
 USA
 USA Customer Service 888-272-1001
Made In USA:
 560 Cottonwood Drive
 Milpitas, CA 95035
 USA

EU Authorized Representative
 Boston Scientific International S.A.
 55 avenue des Champs Pierreux
 TSA 51101
 92729 NANTERRE CEDEX
 FRANCE

UPN Product No. M00440200
 REF Catalog No. 4020
 Date of Manufacture YYYY-MM-DD
 SN Serial Number 12345678

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871-0189-000 REV A 90818619-01A

Boston Scientific
 CAUTION Attention: Consult ACCOMPANYING DOCUMENTS.
 Follow Instructions For Use

Maestro 4000™ Pod, 100W
RF Cardiac Ablation Pod, 100W

Legal Manufacturer
 Manufactured for:
 Boston Scientific Corporation
 One Boston Scientific Place
 Natick, MA 01760-1537
 USA
 USA Customer Service 888-272-1001
Made In USA:
 560 Cottonwood Drive
 Milpitas, CA 95035
 USA


EU Authorized Representative
 Boston Scientific International S.A.
 55 avenue des Champs Pierreux
 TSA 51101
 92729 NANTERRE CEDEX
 FRANCE

UPN Product No. M00440100
 REF Catalog No. 4010
 Date of Manufacture YYYY-MM-DD
 SN Serial Number 12345678

 Rx ONLY
 NON-STERILE
 AP
 Separate Collection

 CE 0344

871-0188-000 REV A 90818620-01A

Boston Scientific
 CAUTION Attention: Consult ACCOMPANYING DOCUMENTS.
 Follow Instructions For Use

Maestro 4000™ Pod, 150W
RF Cardiac Ablation Pod, 150W

Legal Manufacturer
 Manufactured for:
 Boston Scientific Corporation
 One Boston Scientific Place
 Natick, MA 01760-1537
 USA
 USA Customer Service 888-272-1001
Made In USA:
 560 Cottonwood Drive
 Milpitas, CA 95035
 USA

EU Authorized Representative
 Boston Scientific International S.A.
 55 avenue des Champs Pierreux
 TSA 51101
 92729 NANTERRE CEDEX
 FRANCE

UPN Product No. M004EPT40100
 REF Catalog No. EPT4010
 Date of Manufacture YYYY-MM-DD
 SN Serial Number 12345678

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EC REP EU Authorized Representative
 Boston Scientific International S.A.
 55 avenue des Champs Pierreux
 TSA 51101
 92729 NANTERRE CEDEX
 FRANCE

CE 0344

Electrostatic Sensitive Device
 Non-Ionizing Electromagnetic Radiation
 Equipotentiality
 Separate Collection

ONLY CAUTION Attention: Consult ACCOMPANYING DOCUMENTS.

Follow Instructions For Use

Boston Scientific
Maestro 4000™ Controller
RF Cardiac Ablation Controller

UPN Product No. M00440000
REF Catalog No. 4000
MANUFACTURE Date of Manufacture YYYY-MM-DD
SN Serial Number 12345678

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EC REP EU Authorized Representative
 Boston Scientific International S.A.
 55 avenue des Champs Pierreux
 TSA 51101
 92729 NANTERRE CEDEX
 FRANCE

CE 0344

Electrostatic Sensitive Device
 Separate Collection
 Equipotentiality

ONLY CAUTION Attention: Consult ACCOMPANYING DOCUMENTS.

Follow Instructions For Use

Boston Scientific
Maestro 4000™ Remote Control
RF Cardiac Ablation Remote Control

UPN Product No. M00440200
REF Catalog No. 4020
MANUFACTURE Date of Manufacture YYYY-MM-DD
SN Serial Number 12345678

maestro 4000

Boston Scientific

Maestro 4000™ Pod, 100W
RF Cardiac Ablation Pod, 100W

UPN Product No. M00440100
REF Catalog No. 4010
MANUFACTURE Date of Manufacture YYYY-MM-DD
SN Serial Number 12345678

ONLY CAUTION Attention: Consult ACCOMPANYING DOCUMENTS.

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EC REP EU Authorized Representative
 Boston Scientific International S.A.
 55 avenue des Champs Pierreux
 TSA 51101
 92729 NANTERRE CEDEX
 FRANCE

CE 0344

Electrostatic Sensitive Device
 Non-Ionizing Electromagnetic Radiation
 Equipotentiality
 Separate Collection

maestro 4000

871-0188-000 REV A 90818620-01A

Boston Scientific

Maestro 4000™ Pod, 150W
RF Cardiac Ablation Pod, 150W

UPN Product No. M004EPT40100
REF Catalog No. EPT4010
MANUFACTURE Date of Manufacture YYYY-MM-DD
SN Serial Number 12345678

ONLY CAUTION Attention: Consult ACCOMPANYING DOCUMENTS.

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EC REP EU Authorized Representative
 Boston Scientific International S.A.
 55 avenue des Champs Pierreux
 TSA 51101
 92729 NANTERRE CEDEX
 FRANCE

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Legal Manufacturer
Manufactured for:
Boston Scientific Corporation
One Boston Scientific Place
Natick, MA 01760-1537
USA
USA Customer Service 888-272-1001

Made In USA:
560 Cottonwood Drive
Milpitas, CA 95035
USA

CE 0344

cTUVus Mark indicates compliance to UL 60601-1 and CAN/CSA 22.2 601.1 M90 covering electrical safety requirements for the US and Canada.

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Equipotentiality

Rx ONLY

CAUTION
Attention: Consult ACCOMPANYING DOCUMENTS.

Follow Instructions For Use

Boston Scientific

Maestro 4000™ Controller
RF Cardiac Ablation Controller

UPN Product No. M00440000

REF Catalog No. 4000

Date of Manufacture YYYY-MM-DD

SN Serial Number 12345678

VA	POWER INPUT	25 VA
Hz	FREQUENCY	50/60 Hz
~	AC VOLTAGE	100-120 / 220-240V

Legal Manufacturer
Manufactured for:
Boston Scientific Corporation
One Boston Scientific Place
Natick, MA 01760-1537
USA
USA Customer Service 888-272-1001

Made In USA:
560 Cottonwood Drive
Milpitas, CA 95035
USA

CE 0344

cTUVus Mark indicates compliance to UL 60601-1 and CAN/CSA 22.2 601.1 M90 covering electrical safety requirements for the US and Canada.

871-0189-000 REV A

Rx ONLY

CAUTION
Attention: Consult ACCOMPANYING DOCUMENTS.

Follow Instructions For Use

Boston Scientific

Maestro 4000™ Remote Control
RF Cardiac Ablation Remote Control

UPN Product No. M00440200

REF Catalog No. 4020

Date of Manufacture YYYY-MM-DD

SN Serial Number 12345678

Boston Scientific

Maestro 4000™ Pod, 100W
RF Cardiac Ablation Pod, 100W

UPN Product No. M00440100

REF Catalog No. 4010

Date of Manufacture YYYY-MM-DD

SN Serial Number 12345678

Rx ONLY

CAUTION
Attention: Consult ACCOMPANYING DOCUMENTS.

Follow Instructions For Use

Legal Manufacturer
Manufactured for:
Boston Scientific Corporation
One Boston Scientific Place
Natick, MA 01760-1537
USA
USA Customer Service 888-272-1001

Made In USA:
560 Cottonwood Drive
Milpitas, CA 95035
USA

CE 0344

cTUVus Mark indicates compliance to UL 60601-1 and CAN/CSA 22.2 601.1 M90 covering electrical safety requirements for the US and Canada.

871-0188-000 REV A

Rx ONLY

CAUTION
Attention: Consult ACCOMPANYING DOCUMENTS.

Follow Instructions For Use

Boston Scientific

Maestro 4000™ Pod, 150W
RF Cardiac Ablation Pod, 150W

UPN Product No. M004EPT40100

REF Catalog No. EPT4010

Date of Manufacture YYYY-MM-DD

SN Serial Number 12345678

Rx ONLY

CAUTION
Attention: Consult ACCOMPANYING DOCUMENTS.

Follow Instructions For Use

Legal Manufacturer
Manufactured for:
Boston Scientific Corporation
One Boston Scientific Place
Natick, MA 01760-1537
USA
USA Customer Service 888-272-1001

Made In USA:
560 Cottonwood Drive
Milpitas, CA 95035
USA

CE 0344

cTUVus Mark indicates compliance to UL 60601-1 and CAN/CSA 22.2 601.1 M90 covering electrical safety requirements for the US and Canada.

871-0187-000 REV A

Boston Scientific

Maestro 4000™ Pod, 150W
RF Cardiac Ablation Pod, 150W

UPN Product No. M004EPT40100

REF Catalog No. EPT4010

Date of Manufacture YYYY-MM-DD

SN Serial Number 12345678

Rx ONLY

CAUTION
Attention: Consult ACCOMPANYING DOCUMENTS.

Follow Instructions For Use

Legal Manufacturer
Manufactured for:
Boston Scientific Corporation
One Boston Scientific Place
Natick, MA 01760-1537
USA
USA Customer Service 888-272-1001

Made In USA:
560 Cottonwood Drive
Milpitas, CA 95035
USA

CE 0344

cTUVus Mark indicates compliance to UL 60601-1 and CAN/CSA 22.2 601.1 M90 covering electrical safety requirements for the US and Canada.

871-0187-000 REV A

Rx ONLY

CAUTION
Attention: Consult ACCOMPANYING DOCUMENTS.

Follow Instructions For Use

Boston Scientific

Maestro 4000™ Pod, 150W
RF Cardiac Ablation Pod, 150W

UPN Product No. M004EPT40100

REF Catalog No. EPT4010

Date of Manufacture YYYY-MM-DD

SN Serial Number 12345678

Rx ONLY

CAUTION
Attention: Consult ACCOMPANYING DOCUMENTS.

Follow Instructions For Use

Legal Manufacturer
Manufactured for:
Boston Scientific Corporation
One Boston Scientific Place
Natick, MA 01760-1537
USA
USA Customer Service 888-272-1001

Made In USA:
560 Cottonwood Drive
Milpitas, CA 95035
USA

CE 0344

cTUVus Mark indicates compliance to UL 60601-1 and CAN/CSA 22.2 601.1 M90 covering electrical safety requirements for the US and Canada.

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