

**Boston  
Scientific**

Advancing science for life™

# Programming Wand

## Directions for Use

**Rx ONLY**

CAUTION: Federal law restricts this device to sale, distribution and use by or on the order of a physician.

91008426-02 REV A

## Guarantees

Boston Scientific Corporation reserves the right to modify, without prior notice, information relating to its products in order to improve their reliability or operating capacity.

Drawings are for illustration purposes only.

## Trademarks

All trademarks are the property of their respective holders.

## Additional Information

For Indications and related information, see the Indications DFU. For contraindications, warnings, precautions, adverse events summary, physician instructions, sterilization, component disposal, and contact information for Boston Scientific, refer to the Information for Prescribers DFU for your spinal cord stimulator system. For other device-specific information not included in this manual, labeling symbols, and warranty, refer to the appropriate DFU as listed on your Reference Guide.

**USA** For FCC rules and for clinical studies supporting the clinical use of the neurostimulation system, refer to the Information for Prescribers manual for your SCS system.


## Product Model Numbers

Model Number	Description
NM-7190	Programming Wand

## Description

The Programming Wand is used to communicate data between the Clinician Programmer and the Stimulator (Implantable Pulse Generator or Trial Stimulator).

### Programming Wand



## Directions for Use

1. Connect the Programming Wand to the Clinician Programmer using a USB cable. Plug the USB cable into the USB port on the Clinician Programmer (CP) and the USB port (#3 on diagram) on the side of the Programming Wand (Wand).

*NOTE: Clinician Programmer is IEC 60950 compliant.*

The Wand will perform a short self-test. At the end of the self-test, the Wand will beep.

- If the Power Light (#1) turns green, the Wand is ready for use.
- If the Power Light remains red, contact a Boston Scientific representative, or call Customer Service.

2. During normal operation:

- The CP Communication Indicator (#5) will flash whenever data is successfully transmitted between the CP and the Wand.
- The Stimulator Communication Indicator (#2) will flash whenever data is successfully transmitted between the Stimulator and the Wand.
- If communication fails, the Wand will beep up to three times. Move the Wand closer to the Stimulator and/or change the orientation of the Wand to improve wireless communication.

*NOTE: Common sources of interference, such as computers and their power adapters, can reduce the wireless range between the Stimulator and the Wand. To increase range, move the Wand six inches or more from these devices.*

- After communication begins, the Signal Strength bars (#4) will indicate the strength of the wireless connection between the Stimulator and the Wand.
3. After programming is complete, you may disconnect the Wand.
  4. Optional: The Wand is designed to be used with a standard mount (for example, standard camera tripod). See #6 on the diagram.

### **Customer Service**

To order replacement parts, or if you have any questions and need to contact Boston Scientific, choose your locality from the contact list in the *Information for Prescribers* manual for your SCS System.

***This page intentionally left blank.***

***This page intentionally left blank.***

***This page intentionally left blank.***

# Boston Scientific

Advancing science for life™



## Legal Manufacturer

Boston Scientific Neuromodulation  
25155 Rye Canyon Loop  
Valencia, CA 91355 USA  
(866) 789-5899 in US and Canada  
(661) 949-4000, (661) 949-4022 Fax  
(866) 789-6364 TTY  
[www.bostonscientific.com](http://www.bostonscientific.com)  
Email: [neuro.info@bsci.com](mailto:neuro.info@bsci.com)



## Australian Sponsor Address

Boston Scientific (Australia) Pty Ltd  
PO Box 332  
BOTANY  
NSW 1455  
Australia  
Free Phone 1800 676 133  
Free Fax 1800 836 666



## EU Authorized Representative

Boston Scientific Limited  
Ballybrit Business Park  
Galway, Ireland  
T: +33 (0) 1 39 30 97 00  
F: +33 (0) 1 39 30 97 99

© 2015 Boston Scientific Corporation  
or its affiliates. All rights reserved

91008426-02 REV A 2015-09