AMS™ Disposable Dilators

Instructions for Use

English	AMS™ Disposable Dilators Instructions for Use.....	3
Deutsch	AMS™ Einweg-Dilatatoren Gebrauchsanweisung.....	7
Español	Dilatadores desechables AMS™ Instrucciones de uso	11
Français	Dilateurs jetables AMS™ Mode d'emploi	15
Italiano	Dilatatori monouso AMS™ Istruzioni per l'uso.....	19
Português	Dilatadores Descartáveis da AMS™ Instruções de Utilização.....	23

Rx ONLY

AMS™

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Instructions for Use

PRODUCT DESCRIPTION

The AMS Disposable Dilators are sterile, single-use, disposable dilators.

The Dilator is constructed of a smooth, solid, biocompatible thermoplastic material (Ultem™), a material chosen to provide high strength and rigidity.


The 22.1 cm (8.7 in) total length of each dilator is composed of a 7.98 mm (0.314 in) diameter shaft that tapers up to the wider dilation ends and tapers down to the tips of the dilator.

The dilation ends differ by 1 mm in diameter (e.g., one end of the dilator is 9 mm, the opposite end is 10 mm). The Dilator diameters are marked on each shaft near the dilation ends.

Four dilators provide eight dilation dimensions: 9/10, 11/12, 13/14, and 15/16 mm sizes.

Two package configurations accommodate the surgeon's choice of prosthesis.

- For use with three-piece inflatable penile prostheses, two dilators (9/10 and 11/12) are packaged together.
- For use with two-piece inflatable penile prostheses and non-inflatable penile prostheses, four dilators (9/10, 11/12, 13/14 and 15/16) are packaged together.

INDICATIONS FOR USE

The AMS Disposable Dilators are intended for use as manual, surgical instruments to progressively dilate the corpora cavernosa prior to insertion of a penile prosthesis.

Variations in use may occur in specific procedures due to individual technique and patient anatomy. The 1 cm incremental markings along the Dilator shaft may be used as a corporal length measuring tool.

WARNINGS AND PRECAUTIONS

- The Dilators should only be used by physicians who have been trained in dilator use and application.
- Do not use this device in a manner which may cause tissue trauma.
- This device is not designed to be implanted.
- Use caution when using this device with other surgical instruments that may damage the tool.
- Do not use beyond expiration date indicated on the package.
- Inspect the package for damage. If the package integrity is compromised, do not use the Dilators as sterility may have been compromised.
- Dilators are provided sterile and intended for **single use only. DO NOT RESTERILIZE.**
- Use of the Dilators will make impossible latent natural or spontaneous erections or the use of other interventional treatment options for erectile restoration.
- Pre-existing penile or corporal scarring may make surgical use of the Dilators more complicated or impractical.
- Failure to evaluate and promptly treat tunical trauma and perforations may result in a substantial worsening of the condition and may lead to infection and/or loss of tissue.

POTENTIAL RISKS OR ADVERSE EVENTS

- Allergic Reaction
- Altered Therapeutic Response
- Infection
- Prolonged Procedure
- Tissue Damage
- Tunical trauma/perforation

- Urethral trauma/perforation

Note: Refer to the Instructions for Use for the selected penile prosthesis for specific prosthesis risks and surgical-related Adverse Events.

DILATOR USE

1. After incising the tunica albuginea of the corporal body, introduce the 9 mm end of the 9/10 mm Dilator into the distal portion of the corpora, using gentle pressure as the Dilator is advanced. Use care during insertion and dilation to avoid urethral trauma and/or perforation. Keep the Dilator lateral and observe the position of the Dilator tip as it advances against the lateral wall of the tunica albuginea.
2. Palpate the glans to detect when the Dilator reaches the end of the distal portion of the corpora.
3. Remove the Dilator and insert the 9 mm end proximally, using care to dilate the tissue to the point of corporal insertion at the ischial tuberosity.
4. Repeat the process, increasing the Dilator size as needed to progressively expand the corpora tissue.

HOW SUPPLIED AND STORAGE

WARNING: Contents supplied STERILE. Do not use if sterile barrier is damaged. If damage is found, call your AMS representative.

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

Store product in a clean, dry, dark area at room temperature.

PACKAGING

The AMS Disposable Dilators are packaged in a formed polycarbonate inner tray within a formed polycarbonate outer tray. Both inner and outer trays are sealed with a Tyvek™ lid. After sterilization, the double tray package is placed in a dustcover box and labeled.

This page intentionally left blank.
Diese Seite wurde absichtlich freigelassen.
Esta página se deja en blanco intencionadamente.
Page laissée vierge intentionnellement.
Questa pagina è stata lasciata intenzionalmente vuota.
Esta página foi intencionalmente deixada em branco.

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All rights reserved. All trademarks are property of the
Boston Scientific (Australia) Pty Ltd
American Medical Systems, Inc.
U.S.A.
Europe B.V.**
©2017 Boston Scientific Corporation or its affiliates.