## Technical Spotlight Vol.19

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- 1) Larghi A, Rimbaş M, Tringali A, et al. Endoscopic radiofrequency biliary ablation treatment: A comprehensive review. *Dig Endosc.* 2019;31:245–55.
- 2) Ortner ME, Caca K, Berr F, et al. Successful photodynamic therapy for nonresectable cholangiocarcinoma: a randomized prospective study. *Gastroenterology.* 2003;125:1355–63.
- 3) Sofi AA, Khan MA, Das A, et al. Radiofrequency ablation combined with biliary stent placement versus stent placement alone for malignant biliary strictures: a systematic review and meta-analysis. *Gastrointest Endosc.* 2018;87:944–51.
- 4) Cui W, Wang Y, Fan W, et al. Comparison of intraluminal radiofrequency ablation and stents vs. stents alone in the management of malignant biliary obstruction. *Int J Hyperthermia.* 2017;33:853–61.
- 5) Wu TT, Li WM, Li HC, et al. Percutaneous intraductal radiofrequency ablation for extrahepatic distal cholangiocarcinoma: a method for prolonging stent patency and achieving better functional status and quality of life. *Cardiovasc Intervent Radiol* 2017;40:260–9.
- 6) Wang J, Zhao L, Zhou C, et al. Percutaneous intraductal radiofrequency ablation combined with biliary stent placement for nonresectable malignant biliary obstruction improves stent patency but not survival. *Medicine.* 2016; 95: e3329.
- 7) Kadayifci A, Atar M, Forcione DG, et al. Radiofrequency ablation for the management of occluded biliary metal stents. *Endoscopy* 2016;48:1096-101.
- 8) Yoon WJ, Kim YT, Daglellar ES, et al. Evaluation of bipolar radiofrequency ablation for occluded self-expandable metal stents in the bile duct: in vivo and in vitro study. *Endoscopy* 2015;47:1167-70.
- 9) Yang J, Wang J, Zhou H, et al. Efficacy and safety of endoscopic radiofrequency ablation for unresectable extrahepatic cholangiocarcinoma: a randomized trial. *Endoscopy.* 2018;50:751–60.
- 10) Kallis Y, Phillips N, Steel A, et al. Analysis of endoscopic radiofrequency ablation of biliary malignant strictures in pancreatic cancer suggests potential survival benefit. *Dig Dis Sci.* 2015;60:3449–55.
- 11) Inoue T, Ibusuki M, Kitano R, et al. Endobiliary radiofrequency ablation combined with bilateral metal stent placement for malignant hilar biliary obstruction. *Endoscopy* 2020 (in press)
- 12) Tal AO, Vermehren J, Friedrich-Rust M, et al. Intraductal endoscopic radiofrequency ablation for the treatment of hilar non-resectable malignant bile duct obstruction. *World J Gastrointest Endosc.* 2014;6:13–9.

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製品の詳細に関しては添付文書等でご確認いただくか、弊社営業担当へご確認ください。  
© 2020 Boston Scientific Corporation or its affiliates. All rights reserved.  
All trademarks are the property of their respective owners.

**Boston  
Scientific**  
Advancing science for life™

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