

WOLVERINE™

Coronary Cutting Balloon | The ShortCUT Trial




The ShortCUT Trial demonstrated WOLVERINE Cutting Balloon is non-inferior to Shockwave IVL with regard to MSA at the site of maximum calcium in moderate-severe calcified coronary lesions within the total study population.

Presented by Dr. Suzanne Baron, TCT 2025

ShortCUT Study Design

Randomised, non-inferiority trial of N=413 patients at 21 US sites with moderate-severe calcified lesions



Primary Endpoint (non-inferiority):
MSA at site of max calcification

Key Secondary Endpoint (superiority):
Total procedural cost

Key Lesion Characteristics

| | WOLVERINE N=212 | Shockwave N=218 |
|---------------------------|--------------------|--------------------|
| Reference Vessel Diameter | 3.4 mm | 3.5 mm |
| Target Lesion Stenosis | 89% | 87% |
| Median Max Calcium Arc | 338° | 347° |
| Nodular Calcium | 25% | 25% |
| CTO | 13.5% | 12.5% |

WOLVERINE demonstrated non-inferiority to Shockwave in regards to MSA as site of max calcium. The ShortCUT Trial showed no differences in stent expansion and presence of calcium fractures in the total patient cohort.

MSA at site of max calcium in total cohort



P-value (non-inferiority) = 0.007

Secondary Endpoints

Stent Expansion at site of maximum calcium

WOLVERINE
97.7%

Shockwave
97.7%

P-value = 0.99

Calcium Fracture

WOLVERINE
78.9%

Shockwave
77.4%

P-value = 0.72

► Cost

Total procedural costs were significantly lower with WOLVERINE, showing it is a cost-effective strategy for IVUS guided PCI procedures.

► Safety

There were no significant differences in complication associated with randomised device, intraprocedural adverse events or MACCE through 30 days, validating WOLVERINE's comparable safety to Shockwave.


| Endpoint | WOLVERINE N=206 Patients | Shockwave N=207 Patients | P-Value |
|------------------------------------------------|--------------------------------|--------------------------------|---------|
| Complication associated with randomised device | 5.5% | 4.0% | 0.34 |
| In-Hospital MACCE | 1.0% | 1.0% | 1.00 |
| 30-Day MACCE | 2.9% | 2.9% | 1.00 |
| Intraprocedural adverse events | 1.0% | 1.0% | 1.00 |

- ShortCUT trial protocol called for WOLVERINE Cutting Balloon to be sized 0.5mm smaller than reference vessel diameter and inflated to high pressures (16–20 atm). Average pressure used for WOLVERINE based on protocol was 17 atm.
- HD-IVUS was strongly recommended prior to randomisation and mandatory after randomised device use (pre-stent implantation) & at procedure completion.

► Wolverine's Unique Mechanism of Action

Superior outcomes at lower focused pressures

3x the stress, focused at the atherotomes for controlled cracking



Boston Scientific
Advancing science for life™

www.bostonscientific.eu

CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings, and instructions for use can be found in the product labelling supplied with each device or at www.IFU-BSCI.com. Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. This material not intended for use in France.

© 2026 Boston Scientific Corporation or its affiliates. All rights reserved.