


# MASTER THE COMPLEX™

Optimizing revascularization through innovation, training, and education.



Small balloon in the proximal LCx



Crossing the plaque with 1.5 mm Burr



Angiography



Pre-dilatation with a 2.75 mm NC Emerge™ Balloon



Inflation of NC Emerge™ Balloon



Placement of guide wire in LAD



## CASE STUDY

### MASTER PLAQUE MODIFICATION TO TREAT DIFFUSED CALCIFIED LESIONS WITH IVUS SUPPORT

#### Prof. Sławomir Dobrzycki

The Clinic of Invasive Cardiology with Intensive Care and Hemodynamics

University Clinic Hospital in Białystok, Poland

This case illustrates the need for Opticross™ Coronary Imaging Catheter and Rotablator™ Rotational Atherectomy System to treat long diffused calcified lesions.

#### Patient History

- 59 year old patient
- One year earlier the LAD was treated with a 3.0 mm DES
- This time the patient was admitted to the neighboring hospital due to the new symptoms

#### Diagnosis

- After an unsuccessful attempt to open a calcified lesion in the proximal LCx with a small balloon, the patient was transferred to the University Clinic Hospital for a rotablation PCI procedure (**click on video 1**)

#### Procedure

- After crossing with a RotaWire, the plaque was modified with the use of a 1.5 mm Burr (**click on video 2**)
- A smooth channel was created in the proximal LCx and the result was confirmed by angio (**click on video 3**)
- The RotaWire was exchanged with a regular guidewire for better support and then the lesion was pre-dilatated with a 2.75 mm NC Emerge™ balloon (**click on video 4**)
- The NC Emerge™ 2.75 mm balloon was inflated for the second time in the proximal segment of the LCx and the pre-dilatation was successfully completed (**click on video 5**)
- After pre-dilatation a second guide wire was placed in LAD (**click on video 6**)

## KEY TAKEAWAYS

- **Exceptional deliverability of the Opticross™ Coronary Imaging Catheter is the key to assessing vessels with long diffused calcified lesions**
- **Rotablator™ Rotational Atherectomy System was used to modify calcified lesion for optimal stent delivery and implantation**
- **IVUS confirmed presence of diffused calcified lesion in the LCx and LMCA with a need for stenting a long segment of the vessel**
- **Then, IVUS helped optimize stent implantation and confirm the final result**





IVUS Pre Run in the LCx with Opticross™ Catheter

## IVUS Pre Run in the LCx

- The LCx was checked with an OptiCross™ Imaging Catheter
- IVUS imaging confirmed presence of a long diffused calcified lesion in the LCx and LM with a need for stenting a long part of the vessel
- Frames 650-680 reveal a narrowing
- Beginning with frame 850 the second wire becomes visible – the IVUS catheter enters the LM (**click on video 7**)

## IVUS Pre Run in the LAD

- Another run with the use of the Opticross™ Imaging Catheter was performed in the LM and LAD to check for any asymmetry or malapposition of the previously implanted stent
- Frames 440-520 (ostium of the LAD) reveal irregularities in shape and need to stent
- In frames 530 – 970 (LM) both wires are visible and there is a lot of diffused calcium (**click on video 8**)



IVUS Pre Run in the LM and LAD with Opticross™ Catheter



Left Main measurement

## Left Main MLA Measurement by IVUS

- The Minimum Lumen Area (MLA) was 5.03 mm<sup>2</sup>; below the cut-off limit of 6 mm<sup>2</sup> required for the LMCA
- IVUS confirmed the need for stenting the LM and the decision was made to cover the entire diseased area with stents (**click on video 9**)

## Stent

- The first 3.0 x 22 mm DES was implanted in the mid segment of LCx (**click on video 10**)  
Next, Promus Premier™ 3.5 x 16 mm stent was implanted in the Left Main in the direction of the LAD covering the ostium of the LCx (**click on video 11**)
- A specialty wire was selected to cross a cell of the Promus Premier™ stent entering the LCx and the NC Emerge™ 2.75 mm balloon was introduced there (**click on video 12**)
- The NC Emerge™ 2.75 mm balloon was inflated successfully and the ostium of the LCx opened (**click on video 13**)
- A 3.0 x 15 mm DES was introduced in the proximal LCx across the cell of the Promus™ Premier stent (**click on video 14**)
- TAP (T-Stenting and Protrusion) technique was applied to stent the proximal part of LCx (**click on video 15**)




Implantation DES in mid segment of LCx



Implantation of 3.5 x 16 mm DES in Left Main




Introduction of NC Emerge™ 2.75 mm balloon



Inflation of NC Emerge™ Balloon



Placement of a 3.0 x 15 mm DES



TAP technique



Application of two NC Emerge™ Balloon in the bifurcation



POT technique

## Proximal Optimization Technique (POT)

- The final stage was the application of two NC Emerge™ balloons in the sizes 3.0 x 15 mm (LM and LAD) and 2.75 x 15 mm (LM and LCx) in the bifurcation (**click on video 16**)
- Next, POT technique with a short 4.0 x 8 mm NC Emerge™ balloon was applied in the LM area to achieve the optimal stent expansion and strut apposition (**click on video 17**)

## RESULT

- Optimal result in the region of the LM bifurcation LAD and LCx confirmed by angio (**click on video 18**)
- This post-IVUS run in the LM confirms that the stent is well expanded and apposed to the vessel wall (**click on video 19**)



Optimal result of LM, LAD and LCx



Well-apposition of the stent to the vessel wall

**Boston Scientific**  
Advancing science for life™