

LATITUDE™ NXT Patient Management System

Important Safety Information

LATITUDE™ NXT Patient Management is a remote monitoring system that gives your health care provider access to your implanted device data. The LATITUDE Patient Management system is not intended to assist with medical emergencies. If you are not feeling well, call your health care provider or your local emergency service. The Communicator does not provide continuous monitoring.

The Communicator is designed to operate on standard telephone lines like those found in most homes. The Communicator may work on other telephone systems, such as Digital Subscriber Line (DSL) and Voice Over IP (VoIP) Internet systems, if those systems provide an analog interface for connecting the Communicator.

The Communicator is designed to work only with the implanted device of the patient for whom it was prescribed. It will not work with other patients' implanted devices and should be used only as authorized by the prescribing physician. The Communicator is not for use with any pulse generator other than a Boston Scientific device. Ask your physician if you have questions about any risks with using the Communicator or your implanted device.

It is very important that the Communicator remain plugged into the power outlet. Your communicator should remain connected to a telephone line, ethernet adaptor or cellular adaptor. Some household appliances and other sources of electromagnetic energy could interfere with the communication between your Communicator and your implanted device. You should be at least 1 m away from televisions, VCRs, DVD players, personal computers, and other electronic equipment, when you are using the Communicator.

Additional Resources

Visit www.bostonscientific.eu/mylatitude for more information about the MyLATITUDE™ Patient App.

Boston Scientific LATITUDE Customer Support is available by calling [<local hotline>](tel:local_hotline) or by email: latitude.international@bsci.com

All cited trademarks are the property of their respective owners. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device. Information for the use only in countries with applicable health authority product registrations. Material not intended for use in France.


CRM-916601-AA Printed in Germany by medicalvision.

**Boston
Scientific**
Advancing science for life™

www.bostonscientific.eu

© 2020 Boston Scientific Corporation
or its affiliates. All rights reserved.
DINCRM994EA

**Boston
Scientific**
Advancing science for life™



MyLATITUDE™
PATIENT APP



The MyLATITUDE™ Patient App is easy to download on a mobile device, such as a smartphone or tablet.

It shows basic information about the status of your LATITUDE™ NXT Communicator and your implanted Boston Scientific device, providing you with:

- Peace of mind that your LATITUDE Communicator is working properly
- Answers to common questions about your implanted cardiac device
- Setup and troubleshooting for your LATITUDE Communicator

The MyLATITUDE Patient App shows you:

- How to set up your LATITUDE Communicator
- The connection status of your LATITUDE Communicator
- What to do in the event your Communicator is not working
- Information about your scheduled transmission dates
- The status of your implanted device battery
- Answers to common questions about your device and remote monitoring

GET STARTED

- ① Open the App Store within your mobile device



- ② Search for "MyLATITUDE"



- ③ Select "MyLATITUDE Patient App"

Scan the QR code to learn more about the MyLATITUDE patient app

OR

