



# ACURATE Prime XL Clinical Combined Cohort

## 30 Day Outcomes

Boston  
Scientific


**Design:** Combined dataset from clinical studies enrolling subjects with symptomatic severe native aortic stenosis indicated for TAVR, who have a native annulus diameter between 26.5 mm and 29 mm (See Analysis Population graphic).

**Objective:** Evaluate clinical outcomes of the ACURATE Prime XL (29mm).

### Positive, 30-Day results demonstrated:

- No deaths and disabling stroke at 30 Days (low non-disabling stroke rate (2%) through 30 Days)
- Favorable hemodynamics: improved aortic valve area and single-digit gradients were sustained through 30 Days
- Zero patients had moderate or greater paravalvular regurgitation 30 Days

### Analysis Population




### Safety Outcomes [N=102]

	Discharge	30 Days
Death	0.0% (0)	0.0% (0)
Stroke	1.0% (1)	2.0% (2)
Disabling stroke	0.0% (0)	0.0% (0)
Rehospitalization <sup>†</sup>	0.0% (0)	1.0% (1)
Myocardial infarction	0.0% (0)	0.0% (0)
Bleeding (life-threatening/disabling)	0.0% (0)	0.0% (0)
Acute kidney injury (Stage 2/3)	0.0% (0)	0.0% (0)
Major vascular complication	0.0% (0)	0.0% (0)
Repeat procedure for valve-related dysfunction	0.0% (0)	0.0% (0)
New permanent pacemaker	5.9% (6)	12.7% (13)
Prosthetic valve malposition	0.0% (0)	0.0% (0)
Prosthetic valve thrombosis	0.0% (0)	0.0% (0)
Prosthetic valve endocarditis	0.0% (0)	0.0% (0)


<sup>†</sup>Hospitalization for valve-related issues or worsening CHF, per VARC-2 definition

### Valve Hemodynamics



Baseline echo data are site-reported; discharge & 30d echo data are core laboratory-adjudicated

### Paravalvular Regurgitation



Presented by Raj Makkar, MD at PCRLV 2024

CAUTION: The ACURATE Prime Aortic Valve System is CE Marked. In the USA, it is an investigational device. And is restricted under federal law to investigational use only.

Not available for sale. CE 0344

SH-2040607 AA Slide 1 of 3


# Clinical Experience with ACURATE Prime XL

Michael Reardon, MD<sup>1</sup>; Ravi Ramana, DO<sup>2</sup>; Scott Lim, MD<sup>3</sup>; Srinivasa Potluri, MD<sup>4</sup>; Robert Gooley, MBBS (Hons), PhD, MBA<sup>5</sup>; Dale Murdoch, BSc, MBBA<sup>6</sup>; Steven Yakubov, MD<sup>7</sup>; Sanjay Sami, MD<sup>8</sup>; Rodrigo Modolo, MD, PhD<sup>9</sup>; **Raj Makkar, MD<sup>10</sup>**


<sup>1</sup>Methodist DeBakey Heart and Vascular Center, Houston, Texas, USA; <sup>2</sup>Advocate Christ Medical Center, Oak Lawn, Illinois, USA; <sup>3</sup>Heart Care Centers of Illinois, Palos Park, Illinois, USA; <sup>4</sup>University of Virginia, Charlottesville, VA, USA; <sup>5</sup>Baylor Scott and White The Heart Hospital, Plano, TX, USA; <sup>6</sup>Victorian Heart Hospital Monash Health, Clayton, Victoria, Australia; <sup>7</sup>The Prince Charles Hospital, Brisbane, Queensland, Australia; <sup>8</sup>OhioHealth Riverside Methodist Hospital, Columbus, Ohio, USA; <sup>9</sup>Albany Medical College, Albany, NY, USA; <sup>10</sup>Boston Scientific Corporation, Marlborough, MA, USA; <sup>10</sup>Cedars-Sinai Medical Center, Los Angeles, California, USA.

## The ACURATE Prime XL Valve

Leverages ACURATE *neo2* design, with adaptations for compatibility with a larger annulus diameter (26.5mm–29mm)



## Analysis Population




## Baseline Characteristics

ACURATE Prime XL Combined Cohort (N=102)		
Age (years)	79.6 ± 6.8	
Gender	Male: 94.1%	Female: 5.9%
STS Score (%)	3.4 ± 3.1	
Operative Risk (per Case Review Committee)	High/Extreme: 40.2%	
	Intermediate: 42.2%	
	Low: 17.6%	
Mean AV gradient	41.3 ± 10.1 mmHg	
Aortic valve area	0.79 ± 0.2 cm <sup>2</sup>	


Baseline echocardiographic measures are site reported

## Valve Hemodynamics



Baseline echocardiographic data are site-reported; discharge & 30d echocardiographic data are adjudicated by a single core laboratory

## Paravalvular Regurgitation



## Safety Outcomes

	Discharge	30 Days
Death	0.0% (0)	0.0% (0)
Stroke	1.0% (1)	2.0% (2)
Disabling stroke	0.0% (0)	0.0% (0)
Rehospitalization <sup>†</sup>	0.0% (0)	1.0% (1)
Myocardial infarction	0.0% (0)	0.0% (0)
Bleeding (life-threatening/disabling)	0.0% (0)	1.0% (1)
Acute kidney injury (Stage 2/3)	0.0% (0)	0.0% (0)
Major vascular complication	0.0% (0)	0.0% (0)
Repeat procedure for valve-related dysfunction	0.0% (0)	0.0% (0)
New permanent pacemaker	5.9% (6)	12.7% (13)
Prosthetic valve malposition	0.0% (0)	0.0% (0)
Prosthetic valve thrombosis	0.0% (0)	0.0% (0)
Prosthetic valve endocarditis	0.0% (0)	0.0% (0)

<sup>†</sup>Hospitalization for valve-related issues or worsening CHF, per VARC-2 definition

As Presented by Raj Makkar, MD at PCRLV 2024; CAUTION: The ACURATE Prime Aortic Valve is CE-Marked. In the USA, the ACURATE Prime Aortic Valve System is an investigational device and is restricted under federal law to investigational use only. Not available for sale

**Conclusions:** Patients treated with the ACURATE Prime XL valve exhibited good early clinical outcomes at discharge and 30 days post procedure. Although the sample size is small, clinical results in this cohort were in line with those observed in other TAVR studies of patients with extra-large annuli and are encouraging for this new size iteration.

SH-2040607 AA Slide 2 of 3



## Disclaimer

Boston  
Scientific

**CAUTION:** In Europe, the ACURATE neo2™ Aortic Valve System and the ACURATE Prime™ Aortic Valve System are CE-marked. In the USA, the ACURATE neo2 Aortic Valve System and the ACURATE Prime Aortic Valve System are investigational devices and are restricted under federal law to investigational use only. Not available for sale.

**CAUTION:** The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings, and instructions for use can be found in the product labelling supplied with each device or at [www.IFU-BSCI.com](http://www.IFU-BSCI.com). Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. This material not intended for use in France.

2024 Copyright © Boston Scientific Corporation or its affiliates. All rights reserved.

**CE 0344**