


Laparoscopic Transcystic Common Bile Duct Exploration (LCBDE)

**Boston
Scientific**
Advancing science for life™

SpyGlass™ Discover Digital Catheter Technique Spotlight

Presented by:

Jordan Wood, MD and Adam Reid, MD, FACS
Southern Illinois University School of Medicine
Memorial Medical Center
Springfield, Illinois



Patient Presentation:

A 44-year-old male with history of Roux-en-Y gastric bypass presented with epigastric abdominal pain, leukocytosis, and hyperbilirubinemia. On abdominal CT and ultrasound, he was found to have evidence of acute cholecystitis. Subsequent MRCP revealed choledocholithiasis.

Procedure:

He was taken for laparoscopic cholecystectomy and intraoperative cholangiogram revealed a filling defect. (figure 1) He had been consented for laparoscopic assisted ERCP. However, his remnant stomach was difficult to access due to intraabdominal adhesions. Instead, laparoscopic transcystic common bile duct exploration (LCBDE) was performed. (figure 2) No dilation of the cystic duct was required, and the scope was passed easily into the common bile duct. A calculus was readily identified and removed using basket extraction. (figure 3) Post stone extraction we were able to use SpyGlass Discover Digital Catheter to identify the ampulla and confirm ductal clearance. (figure 4)




Defect (figure 1).



Performed (figure 2).



Extraction (figure 3).



Clearance (figure 4).

Discover Impact

Post-operatively, his diet was advanced and antibiotics were discontinued. His symptoms had improved the following morning and bilirubin was stable. He was discharged home on post-operative day one.

Discussion:

Choledocholithiasis frequently requires management with ERCP as well as surgical cholecystectomy for prevention of recurrence and often concurrent cholecystitis. This is further complicated in patients with Roux-en-Y anatomy, making a traditional endoscopic approach ERCP impossible. The Choledochoscope enabled us to remove the common bile duct stone at the same operation, using the cystic ductotomy created for the cholangiogram. This allowed us to forego the risk and time required for dissection of the remnant stomach, additional incision and port placement in the stomach, as well as the need to involve a gastroenterology team. In this patient's case, the recovery was fast and significant risk to the patient was avoided.

IMPORTANT INFORMATION: These materials are intended to describe common clinical considerations and procedural steps for the use of referenced technologies but may not be appropriate for every patient or case. Decisions surrounding patient care depend on the physician's professional judgment in consideration of all available information for the individual case. Boston Scientific (BSC) does not promote or encourage the use of its devices outside their approved labeling. Case studies are not necessarily representative of clinical outcomes in all cases as individual results may vary.

Images provided courtesy of Dr. Adam Reid.

CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labelling supplied with each device. Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. This material not intended for use in France. Rx Only.

**Boston
Scientific**
Advancing science for life™

Boston Scientific Corporation
300 Boston Scientific Way
Marlborough, MA 01752-1234
www.bostonscientific.com

©2021 Boston Scientific Corporation or
its affiliates. All rights reserved.

ENDO-926309-AA