



CRE™ WIREGUIDED BALLOON DILATOR

Technique Spotlight



Fadlallah Habr, MD

Assistant Professor of Medicine at Brown University
Rhode Island Hospital, Providence, RI




Figure 1



Figure 2



Figure 3




Figure 4

PATIENT HISTORY

A 71-year-old female with a history of chronic heartburn and dysphagia was recently diagnosed with a distal esophageal peptic stricture. She had undergone two endoscopic dilations with recurrence of her symptoms. She was referred to our center for further therapy.

PROCEDURE

After appropriate sedation, the endoscope was advanced into the esophagus which showed a tight stricture of about 2cm in length at the gastroesophageal junction (Figure 1). The endoscope was advanced with difficulty into the stomach. Due to the recurrent nature of the stricture, it was decided to place a Polyflex™ Stent after dilation. A CRE™ Wireguided Balloon Dilator was then advanced into the scope's working channel and positioned in the middle of the stricture. Incremental dilation was performed from 10-13.5mm, first using a 10-11-12mm CRE Balloon and then a 12-13.5-15mm CRE Balloon. While slowly rotating the camera lens, the inflated balloon was pulled against the tip of the endoscope. This allowed the mucosal lining of the stricture to be clearly visualized through the end of the balloon during dilation (Figure 2). After dilating the stricture (Figure 3), a 21mm Polyflex Stent was placed (Figure 4). The patient tolerated the procedure very well and was discharged home.

DISCUSSION

Benign strictures can involve any segment in the GI tract. Balloon dilation has proven to be a safe and effective technique to treat those strictures. It is usually performed under direct endoscopic visualisation, although it can also be performed under fluoroscopic guidance. I prefer to use the CRE Balloon rather than other modalities because I can visualise the area that is being dilated. Pulling the balloon snug against the tip of the endoscope after inflation allows me to see the strictured mucosa while performing the dilation. This may help prevent perforation, especially in very tight and friable strictures, or help identify if a large tear is occurring. The rounded ends of the CRE Balloon provide me good visibility of the stricture, but I also have the ability to feel the stricture as well. In order to do this, I will move the scope through the stricture with the inflated balloon pulled against its tip. This back and forth motion, while the balloon is inflated, also aids in dilating the stricture.

This case reflects the opinion of and is the responsibility of the author.

European Distribution Centre –

The Netherlands
T: +31 45 54 67 700 F: +31 45 54 67 800

European Headquarters –

Paris
T: +33 1 57 66 80 00 F: +33 1 57 66 84 99

Argentina (Freephone)
T: 0800 555 2678 F: +5411 4896 8517

Australia/New Zealand
T: +61 2 8063 8100 F: +61 2 9330 1404

Austria
T: +43 1 60 810 F: +43 1 60 810 60

Belgium (Freephone)
T: 0800 94 494 F: 0800 93 343

Brazil
T: +55 11 5502 8500 F: +55 11 5103 2212

Canada
T: +1 888 359 9691 F: +1 888 575 7396

Chile
T: +562 445 4904 F: +562 445 4915

China – Beijing
T: +86 10 8525 1588 F: +86 10 8525 1566

China – Guangzhou
T: +86 20 8767 9791 F: +86 20 8767 9789

China – Shanghai
T: +86 21 6141 5959 F: +86 21 6141 5900

Colombia

T: +57 1 629 5045 F: +57 1 612 4761

Czech Republic

T: +420 296 331 901 F: +420 296 331 935

Denmark (Freephone)
T: 80 30 80 02 F: 80 30 80 05

Eire (Freephone)
T: 1800 882 969 F: 1800 882 968

Finland
T: +358 20 762 8882 F: +358 20 762 8883

France
T: +33 1 39 30 49 00 F: +33 1 39 30 49 01

Germany (Freephone)
T: 08000 723300 F: 08000 723319

Greece
T: +30 210 9542 300 F: +30 210 9542 310

Hong Kong
T: +852 2960 7100 F: +852 2563 5276

Hungary
T: +36 1 456 30 40 F: +36 1 456 30 41

India – Bangalore
T: +91 80 2212 4928/9 F: +91 80 2207 5153

India – Chennai
T: +91 44 2220 1879 F: +91 44 2220 0716

India – Delhi
T: +91 11 4243 2222 F: +91 11 2610 0808

India – Mumbai

T: +91 22 4030 9165 F: +91 22 4040 9199

Italy

T: +39 010 60 60 1 F: +39 010 60 60 200

Korea
T: +82 2 3476 2121 F: +82 2 3476 1776

Mexico
T: +52 55 5687 63 90 F: +52 55 5687 62 28

Middle East/Gulf/North Africa
T: +961 1 805 410 F: +961 1 805 445

The Netherlands
T: +31 30 602 55 44 F: +31 30 602 55 05

Norway (Freephone)
T: 800 104 04 F: 800 101 90

Poland
T: +48 22 435 14 14 F: +48 22 435 14 10

Portugal
T: +351 1 381 25 40 F: +351 21 381 25 58

South Africa
T: +27 11 840 8600 F: +27 11 463 6077

South East Asia – Malaysia
T: +60 3 2283 3813 F: +60 3 2284 3813

South East Asia – Philippines
T: +63 2 687 6994 F: +63 2 687 3047

South East Asia – Singapore
T: +65 6418 8888 F: +65 6418 8899

South East Asia – Thailand

T: +66 2 6380 100 F: +66 2 6380 400

South East Europe

T: +30 210 9542 300 F: +30 210 9542 310

Spain – Barcelona
T: +34 93 444 72 00 F: +34 93 405 90 45

Spain – Madrid
T: +34 901 11 12 15 F: +34 91 319 50 03

Sweden
T: +46 42 25 69 00 F: +46 42 25 69 69

Switzerland
T: +41 32 626 57 00 F: +41 32 626 57 01

Taiwan
T: +886 2 2747 7278 F: +886 2 2747 7270

Turkey – Istanbul
T: +90 216 464 36 66 F: +90 216 464 36 67

UK (Freephone)
T: 0844 800 4512 F: 0844 800 4513

Uruguay
T: +59 82 900 6212 F: +59 82 900 6212

Venezuela
T: +58 212 959 6275 F: +58 212 959 5328

**Boston
Scientific**

Delivering what's next.™

www.bostonscientific.com
www.bostonscientific-international.com

RCS Nanterre B420 668 420
© 2008 Boston Scientific Corporation
or its affiliates. All rights reserved.
DINEND2209EA