

Ultraflex Single-Use Uncovered and Covered Tracheobronchial Stent System (Each)

UNCOVERED		COVERED		Expanded Stent OD (mm)	Expanded Stent Length (mm)	Cover Length (if applicable) (MM)	Tip Max OD (mm)
Uncovered Distal Release	Uncovered Proximal Release	Covered Distal Release	Covered Proximal Release				
-	M00564640	-	-	8	20	-	4.1
-	M00564650	M00564740	-	8	40	25	4.1
-	M00564660	-	-	10	20	-	4.1
M00564500	-	M00564750	-	10	30	15	4.1
-	M00564670	M00564760	-	10	40	25	4.1
-	M00564680	-	-	12	20	-	4.1
M00564510	-	M00564770	-	12	30	15	4.1
-	M00564690	M00564780	-	12	40	25	4.1
-	M00564700	-	-	14	20	-	4.1
M00564520	-	M00564790	-	14	30	15	4.1
-	M00564710	M00564800	-	14	40	25	4.1
-	M00564720	M00564810	-	14	60	45	4.1
-	-	M00564820	-	14	80	65	4.1
M00564530	-	M00564830	-	16	40	25	5.3
M00564540	-	M00564840	-	16	60	45	5.3
-	-	M00564850	-	16	80	65	5.3
M00564560	-	M00564860	-	18	40	25	5.3
M00564570	-	M00564870	-	18	60	45	5.3
-	-	M00564880	-	18	80	65	5.3
M00564590	-	M00564890	-	20	40	25	5.3
M00564600	-	M00564900	-	20	60	45	5.3
-	-	M00564910	-	20	80	65	5.3

Pulmonary Jagwire™ Guidewire (Box 2)


Order Number	O.D. (in)/(mm)	Total Length (cm)
M00515171	0.035 / 0.89	180

CRE™ Pulmonary Balloon Dilatation Catheter (Each)

Order Number	Diameter @3 ATM	Diameter (mm) @ Intermediate ATM	Diameter (mm) @Max. Inflation ATM	Balloon Length (cm)	Catheter Length (cm)
M00550300	12	13.5 @ 4.5	15 @ 8	5.5	110
M00550310	15	16.5 @ 4.5	18 @ 7	5.5	110
M00550320	18	19 @ 4.5	20 @ 6	5.5	110
M00550330	8	9 @ 5.5	10 @ 9	3.0	110
M00550340	10	11 @ 5	12 @ 8	3.0	110
M00550350	12	13.5 @ 4.5	15 @ 8	3.0	110

All cited trademarks are the property of their respective owners. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device. Information for the use only in countries with applicable health authority product registrations. Material not intended for use in France. Product available in the European Economic Area (EEA) only. Please check availability with your local sales representative or customer service.

ENDD-460902-AA JUN2017 Printed in Germany by medicalvision.



The Ultraflex Tracheobronchial Stent System is indicated for use in the treatment of tracheobronchial strictures produced by malignant neoplasms.

Accommodate Varying Airway Anatomy without Kinking

Knitted Nitinol Design

Stent geometry is designed to adapt to anatomical contours and exert constant, gentle radial pressure to maintain patency

Wide Range of Sizes

Variety of lengths and diameters in both covered and uncovered designs is intended to allow for complete bridging of stricture

Clear Secretions

Flexible Open Loop Design

Epithelisation of uncovered stent may promote mucociliary clearance

Resist Migration

Uncovered Ends

Epithelisation of ends may limit migration

Resist Tumor Ingrowth

Polyurethane Covering

On the covered version, covering helps resist tumor ingrowth.

Delivery System

Low Profile

The compressed stent and delivery system have between a 5-9 mm outer diameter. The system is designed to facilitate advancement across tumors and may be placed via flexible or rigid bronchoscopy

Flexibility


The flexible delivery catheter is designed to enhance the ease of navigation through the airway

Radiopaque Markers


Radiopaque markers on the delivery catheter are designed to target the deployed position of the stent

Distal or Proximal Release

Different release systems are designed to allow the physician greater control over stent deployment



Overview of Procedural Steps*



1 Location of Stricture and Pre-Dilation

- If necessary, airway is pre-dilated to 75% of normal size

2 Determination of Stricture Length and Diameter

- Via Bronchoscopic examination
- Via Fluoroscopic examination
- Via CT Imaging

3 Stricture Examination

- Bronchial Branches are identified
- Stricture Margins are identified (radiopaque markers may be used)
- Intraluminal tumor or granulation tissue is removed if necessary

4 Selection of Stent Size

Length


- Stricture must be completely bridged
- Normal mucosa should be overlapped by 1-2 cm

Diameter

- Should be equal to normal proximal lumen

5 Guidewire Insertion

- 0.035" Jagwire™ Guidewire is placed
- Bronchoscope is removed



6 Delivery Catheter is Inserted Over Guidewire

7 Advancement of Stent to Stricture

- Bronchoscope is re-placed into the airway
- Under bronchoscopic visualisation, catheter is advanced over the guidewire
- Using fluoroscopic guidance and radiopaque markers, stent is centered across stricture

8 Deployment of Stent

- Holding catheter steady, ring is gently pulled with finger to deploy stent

9 Adjustment of Stent Position – OPTIONAL

- Catheter is pulled toward the operator to reposition†
- Deployment is then continued

10 Complete Deployment Confirmation and System Removal

Deployment Confirmed

- Via Bronchoscopic examination
- Via Fluoroscopic examination

Delivery System Removed

- Delivery system is removed carefully by rotating catheter
- Bronchoscope is removed.

* This overview is provided for illustrative purposes only and is intended only as a brief summary of how placement procedures for the Ultraflex Tracheobronchial Stent System are generally performed. It is not a substitute for exercise of proper medical judgment by a physician in the care of specific patients, or for the Instructions for Use included with the device.

† For distal and proximal release systems, the inner catheter and stent can be pulled towards the operator and away from the operator, respectively. Repositioning may be limited by scope position or tortuous anatomy, among other factors.