



OBJECTIVE

- ▶ Post-market clinical follow-up study to evaluate the safety and effectiveness of the POLARxTM Cryoablation Balloon for pulmonary vein isolation (PVI) to treat paroxysmal atrial fibrillation (PAF) using real-world data.

METHODS

- ▶ Prospective, non-randomized, multicenter international registry ([NCT04250714](#)).
- ▶ 399 patients were enrolled across 19 centers between August 2020 and May 2021.


PROCEDURAL CHARACTERISTICS¹

- ▶ Data on 372 de novo PVI procedures (n = 2190 ablations) were collected.
- ▶ Acute PV isolation occurred in 96.8% of PVs.
- ▶ The procedure time was 68.2 ± 24.6 , left atrial dwell time was 46.6 ± 18.3 minutes and the fluoroscopy time was 15.6 ± 9.6 minutes.
- ▶ Grade 3 or 4 occlusion was achieved in 98.2% of PVs with a 71.2% rate of single-shot isolation.
- ▶ The average nadir ablation temperature was $-56.3 \pm 6.5^{\circ}\text{C}$.

SAFETY²

- ▶ The safety endpoint event rate was 4.6% (**Figure 1**). Endpoint events included serious vascular access complications (2.6%), cardiac tamponade/perforation (0.5%), thromboembolism/air embolism (0.5%), myocardial infarction (0.3%), persistent gastroparesis/vagus nerve injury (0.3%) persistent phrenic nerve injury (0.3%), and stroke/cerebrovascular accident (0.3%).


Figure 1. Freedom from Primary Safety Event



EFFICACY²

- ▶ The 12-month freedom from any arrhythmia was 83.5% and freedom from atrial fibrillation 88.1% (**Figures 2, 3**).

Figure 2. Freedom from Documented Arrhythmia Recurrence




✓ **83.5%**
Arrhythmia
Recurrence Free



EFFICACY (cont.)

Figure 3. Freedom from AF Recurrence



✓ **88.1%**
AF Recurrence Free

- ▶ Redo procedures were conducted on 19 patients. In 14/19 patients reconnection of at least one PV could be identified (RSPV n=8, RIPV n=12, LSPV n=9, LIPV n=11).
- ▶ Freedom from any arrhythmia was associated with lower nadir temperature (p=0.008) and longer time to thaw (p=0.05) during the index procedure.

CONCLUSION

- ▶ In the POLAR ICE real-world registry there was a low safety event rate (4.6%).
- ▶ The one-year arrhythmia recurrence free rates were 83.5%, with an AF recurrence free rate of 88.1%.
- ▶ Lower nadir temperatures and a longer thawing time were acute predictors of long-term clinical success.



[POLARxTM Cryoablation System
Indications, Safety, and Warnings](#)

Boston Scientific
Advancing science for lifeTM

1. Martin CA, Tilz RRR, Anic A, et al. Acute procedural efficacy and safety of a novel cryoballoon for the treatment of paroxysmal atrial fibrillation: results from the POLAR ICE study. *J Cardiovasc Electrophysiol.* 2023;34:833-840.
2. Luik, et al. Long-term success rates of a stable, low pressure cryoballoon for the treatment of paroxysmal atrial fibrillation: results of the prospective, international, multicenter POLAR-ICE study. Presented at: ESC 2023; August 2023; Amsterdam, NL.

All trademarks are the property of their respective owners.

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician. Rx only. Prior to use, please see the complete "Instructions for Use" for more information on Indications, Contraindications, Warnings, Precautions, Adverse Events, and Operator's Instructions.

Cardiology
300 Boston Scientific Way
Marlborough, MA 01752-1234
www.bostonscientific.com

Medical Professionals:
1.800.CARDIAC (227.3422)
Customer Service:
1.888.272.1001

© 2025 Boston Scientific Corporation
or its affiliates. All rights reserved.

EP-1668804-AA