




**FARAPULSE™**  
Pulsed Field Ablation System  
**Clinical Compendium**



# TABLE OF CONTENTS

<b>Publication Listing by Topic</b>	<b>6-8</b>
<b>Case Studies</b>	<b>9-10</b>
<b>Clinical Publications 2023</b>	<b>11-19</b>
<b>Comparison of Pulsed-Field Ablation versus Very High-Power Short Duration-Ablation for Pulmonary Vein Isolation</b> Wörmann J, Schipper J, Lüker J, et al.....	<b>11</b>
<b>Pulsed-Field Ablation Versus Single Catheter High-Power Short-Duration Radiofrequency Ablation for Atrial Fibrillation: Procedural Characteristics, Myocardial Injury and Midterm Outcomes</b> Badertscher P, Weidlich S, Serban T, et al.....	<b>11</b>
<b>Quantitative Assessment of Transient Autonomic Modulation after Single-Shot Pulmonary Vein Isolation with Pulsed-Field Ablation</b> Del Monte A, Cespón Fernández M, Vetta G, et al.....	<b>11-12</b>
<b>Left Atrial Posterior Wall Isolation with Pulsed Field Ablation in Persistent Atrial Fibrillation</b> Gunawardene M, Frommeyer G, Ellermann C, et al.....	<b>12</b>
<b>Pulsed-Field vs. Cryoballoon vs. Radiofrequency Ablation: A Propensity Score Matched Comparison of One-Year Outcomes after Pulmonary Vein Isolation in Patients with Paroxysmal Atrial Fibrillation</b> Maurhofer J, Kueffer T, Madaffari A, et al.....	<b>12</b>
<b>Long-Term Clinical Outcomes of Pulsed Field Ablation in the Treatment of Paroxysmal Atrial Fibrillation</b> Musikantow D, Neuzil P, Anic A, et al.....	<b>13</b>
<b>Early Recurrences Predict Late Therapy Failure after Pulsed Field Ablation of Atrial Fibrillation</b> Plank K, Bordignon S, Urbanek L, et al.....	<b>13</b>
<b>Characterization of Durability and Reconnection Patterns at Time of Repeat Ablation after Single-Shot Pulsed Field Pulmonary Vein Isolation</b> Ruwald M, Haugdal M, Worck R, et al.....	<b>13</b>
<b>Pulsed Field or Conventional Thermal Ablation for Paroxysmal Atrial Fibrillation</b> Reddy VY, Gerstenfeld EP, Natale A, et al.....	<b>14</b>
<b>Comparison of Pulsed Field Ablation and Cryoballoon Ablation for Pulmonary Vein Isolation</b> Schipper H, Steven D, Lüker J, et al.....	<b>14</b>
<b>Pulsed Field Ablation-Based Pulmonary Vein Isolation Using a Simplified Single-Access Single-Catheter Approach – The Fast and Furious PFA Study</b> Tilz R, Vogler J, Kirstei B, et al.....	<b>15</b>
<b>Pulsed-Field Ablation on Mitral Isthmus in Persistent Atrial Fibrillation - Preliminary Data on Efficacy and Safety</b> Davong B, Adeliño R, Delasnerie H, et al.....	<b>15</b>
<b>Pulmonary Vein Isolation Durability and Lesion Regression in Patients with Recurrent Arrhythmia after Pulsed Field Ablation</b> Kueffer T, Stefanova A, Madaffari A, et al.....	<b>15</b>
<b>Acute Lesion Extension Following Pulmonary Vein Isolation with Two Novel Single Shot Devices: Pulsed Field Ablation versus Multielectrode Radiofrequency Balloon</b> My I, Lemoine M, Butt M, et al.....	<b>16</b>
<b>Pulsed Field Versus Cryoballoon Pulmonary Vein Isolation for Atrial Fibrillation: Efficacy, Safety, and Long-Term Follow-Up in a 400-Patient Cohort</b> Urbanek L, Bordignon S, Schaack D, et al.....	<b>16</b>
<b>European Real-World Outcomes with Pulsed Field Ablation in Patients with Symptomatic Atrial Fibrillation - Lessons from the Multicenter EU-PORIA Registry</b> Schmidt B, Bordignon S, Neven K, et al.....	<b>16</b>
<b>Electrophysiological Findings during Re-Do Procedures after Single-Shot Pulmonary Vein Isolation for Atrial Fibrillation with Pulsed Field Ablation</b> Magni F, Scherr D, Manninger M, et al.....	<b>17</b>
<b>Lesion Formation Following Pulsed Field Ablation for Pulmonary Vein and Posterior Wall Isolation</b> Sohns C, Fink T, Braun M, et al.....	<b>17</b>
<b>Safety and Effectiveness of Pulsed Field Ablation to Treat Atrial Fibrillation: One-Year Outcomes From the MANIFEST-PF Registry</b> Turagam MK, Neuzil P, Schmidt B, et al.....	<b>17</b>
<b>Bronchial Safety After Pulsed-Field Ablation for Paroxysmal Atrial Fibrillation</b> Füting A, Reinsch N, Brokkaar L, et al.....	<b>18</b>
<b>Pulsed Field Ablation to Treat Atrial Fibrillation: Autonomic Nervous System Effects</b> Musikantow DR, Neuzil P, Petru J, et al.....	<b>18</b>

# TABLE OF CONTENTS

<b>Clinical Publications 2023</b>		<b>11-19</b>
<b>PFA for treatment of AF in patients with Congenital Anomalies of Cardiac Veins</b> Castiglione A, Küffer T, Gräni C, et al.....		<b>18</b>
<b>Effects of Pulsed Field Ablation on Autonomic Nervous System in Paroxysmal Atrial Fibrillation: A Pilot Study</b> Guo F, Wang J, Deng Q, et al.....		<b>19</b>
<b>Visualization of Fibroblast Activation Using 68Ga- FAPI PET/CT after Pulmonary Vein Isolation with Pulsed Field Compared with Cryoballoon Ablation</b> Kupusovic J, Kessler L, Bruns F, et al.....		<b>19</b>
<b>A Randomized Controlled Trial of Pulsed Field Ablation versus Standard-of-Care Ablation for Paroxysmal Atrial Fibrillation: The ADVENT Trial Rationale and Design</b> Reddy VY, Lehmann JW, Gerstenfeld EP, et al.....		<b>19</b>
<b>Pulsed Field Ablation in Real-World Atrial Fibrillation Patients: Clinical Recurrence, Operator Learning Curve and Re-Do Procedural Findings</b> Ruwald MH, Johannessen A, Lock Hansen M, et al.....		<b>19</b>
<b>Clinical Publications 2022</b>		<b>20-24</b>
<b>Pulsed-Field Ablation for the Treatment of Left Atrial Reentry Tachycardia</b> Kueffer T, Seiler J, Madaffari A, et al.....		<b>20</b>
<b>Findings from Repeat Ablation using High-Density Mapping after Pulmonary Vein Isolation with Pulsed Field Ablation</b> Tohoku S, Chun J, Bordignon S, et al.....		<b>20</b>
<b>Pulsed Field Ablation-Based Pulmonary Vein Isolation in Atrial Fibrillation Patients with Cardiac Implantable Electronic Devices: Practical Approach and Device Interrogation (PFA in CIEDs)</b> Chen S, Chun J, Bordignon S, et al.....		<b>20</b>
<b>Initial Experience with Pulsed Field Ablation for Atrial Fibrillation</b> Magni F, Mulder B, Groenveld H, et al.....		<b>21</b>
<b>Pulsed Field Ablation in Patients with Complex Consecutive Atrial Tachycardia in Conjunction with Ultra-High-Density Mapping: Proof of Concept</b> Gunawardene M, Schaeffer B, Jularic M, et al.....		<b>21</b>
<b>Pulsed-Field Ablation-Based Pulmonary Vein Isolation: Acute Safety, Efficacy and Short-Term Follow-up in a Multi-Center Real World Scenario</b> Lemoine MD, Fink T, Mencke C, et al.....		<b>21</b>
<b>Cerebral Safety After Pulsed Field Ablation for Paroxysmal Atrial Fibrillation</b> Reinsch N, Fütting A, Höwel D, et al.....		<b>22</b>
<b>Catheter Ablation Induced Phrenic Nerve Palsy by Pulsed Field Ablation—Completely Impossible? A Case Series</b> Pansera F, Bordignon S, Bologna F, et al.....		<b>22</b>
<b>Multi-National Survey on the Methods, Efficacy, and Safety on the Post-Approval Clinical Use of Pulsed Field Ablation (MANIFEST-PF)</b> Ekanem E, Reddy VY, Schmidt B, et al.....		<b>22</b>
<b>Pulsed Field Ablation for Pulmonary Vein Isolation: Real-World Experience and Characterization of the Antral Lesion Size Compared with Cryoballoon Ablation</b> Blockhaus C, Guelker J, Feyen L, et al.....		<b>23</b>
<b>Validation of a Multipolar Pulsed-Field Ablation Catheter for Endpoint Assessment in Pulmonary Vein Isolation Procedures</b> Kueffer T, Baldinger S, Servatius H, et al.....		<b>23</b>
<b>5S Study: Safe and Simple Single Shot Pulmonary Vein Isolation with Pulsed Field Ablation Using Sedation</b> Schmidt B, Bordignon S, Tohoku S, et al.....		<b>23</b>
<b>Characterization of Circumferential Antral Pulmonary Vein Isolation Areas Resulting from Pulsed-Field Catheter Ablation</b> Bohnen M, Weber R, Minners J, et al.....		<b>24</b>
<b>First Experience with Pulsed Field Ablation as Routine Treatment for Paroxysmal Atrial Fibrillation</b> Fütting A, Reinsch N, Höwel D, et al.....		<b>24</b>
<b>Troponin Release after Pulmonary Vein Isolation Using Pulsed Field Ablation Compared to Radiofrequency and Cryoballoon Ablation</b> Krisai P, Knecht S, Badertscher P, et al.....		<b>24</b>
<b>Pulsed Field Ablation Combined with Ultra-High-Density Mapping in Patients Undergoing Catheter Ablation for Atrial Fibrillation: Practical and Electrophysiological Considerations</b> Gunawardene M, Schaeffer B, Jularic M, et al.....		<b>24</b>

# TABLE OF CONTENTS

<b>Clinical Publications 2021</b>	<b>25-26</b>
<b>Does Pulsed Field Ablation Regress Over Time? A Quantitative Temporal Analysis of Pulmonary Vein Isolation</b> Kawamura I, Neuzil P, Shivamurthy P, et al.....	<b>25</b>
<b>Pulsed Field Ablation Prevents Chronic Atrial Fibrotic Changes and Restrictive Mechanics After Catheter Ablation for Atrial Fibrillation</b> Nakatani Y, Sridi-Cheniti S, Cheniti G, et al.....	<b>25</b>
<b>Pulsed Field Ablation of Paroxysmal Atrial Fibrillation: 1-year Outcomes of IMPULSE, PEFCAT, and PEFCAT II</b> Reddy VY, Dukkupati SR, Neuzil P, et al. ....	<b>25</b>
<b>How Does the Level of Pulmonary Venous Isolation Compare between Pulsed Field Ablation and Thermal Energy Ablation (Radiofrequency, Cryo, or Laser)?</b> Kawamura I, Neuzil P, Shivamurthy P, et al. ....	<b>26</b>
<b>Pulsed Field Ablation Selectively Sparing the Oesophagus During Pulmonary Vein Isolation for Atrial Fibrillation</b> Cochet H, Nakatani Y, Sridi-Cheniti S, et al.....	<b>26</b>
<b>Pulsed Field Ablation: A Promise That Came True</b> Ante A, Breskovic T, Sikiric I.....	<b>26</b>
<b>Clinical Publications 2020</b>	<b>27</b>
<b>Pulsed Field Ablation in Patients with Persistent Atrial Fibrillation</b> Reddy VY, Anic A, Koruth J, et al.....	<b>27</b>
<b>Ostial Dimensional Changes After Pulmonary Vein Isolation: Pulsed Field Ablation vs. Radiofrequency Ablation</b> Kuroki K, Whang W, Eggert C, et al.....	<b>27</b>
<b>Clinical Publications 2019</b>	<b>28</b>
<b>Pulsed Field Ablation for Pulmonary Vein Isolation in Atrial Fibrillation</b> Reddy VY, Neuzil P, Koruth JS, et al.....	<b>28</b>
<b>Clinical Publications 2018</b>	<b>28</b>
<b>Ablation of Atrial Fibrillation with Pulsed Electric Fields</b> Reddy VY, Koruth J, et al.....	<b>28</b>

# TABLE OF CONTENTS

<b>Preclinical Publications 2023</b>	<b>29</b>
<b>Electrophysiology, Pathology, and Imaging of Pulsed Field Ablation of Scarred and Healthy Ventricles in Swine</b> Kawamura I, Reddy V, Santos-Gallego C, et al. ....	29
<b>Preclinical Publications 2022</b>	<b>29</b>
<b>Effect of Epicardial Pulsed Field Ablation Directly on Coronary Arteries</b> Higuchi S, Im S, Stillson C, et al. ....	29
<b>Pulsed Field Ablation of Left Ventricular Myocardium in a Swine Infarct Model</b> Im S, Higuchi S, Lee A, et al. ....	29
<b>Preclinical Publications 2020</b>	<b>30</b>
<b>Pulsed Field Ablation vs. Radiofrequency Ablation: Esophageal Effects in a Novel Preclinical Model</b> Koruth JS, Kuroki K, Kawamura I, et al. ....	30
<b>Preclinical Publications 2019</b>	<b>30</b>
<b>Preclinical Evaluation of Pulsed Field Ablation: Electrophysiological and Histological Assessment of Thoracic Vein Isolation</b> Koruth JS, Kuroki K, Iwasawa J, et al. ....	30
<b>Endocardial Ventricular Pulsed Field Ablation: A Proof-of-Concept Preclinical Evaluation</b> Koruth JS, Kuroki K, Iwasawa J, et al. ....	30

# PUBLICATION LISTING BY TOPIC

## Clinical Trial Outcomes Data

- Reddy, et al., [Pulsed Field Ablation of Paroxysmal Atrial Fibrillation: 1-Year Outcomes of IMPULSE, PEFCAT, and PEFCAT II](#)

## 2023

- Mustkantow, et al., [Long-Term Clinical Outcomes of Pulsed Field Ablation in the Treatment of Paroxysmal Atrial Fibrillation](#)
- Reddy, et al., [Pulsed Field or Conventional Thermal Ablation for Paroxysmal Atrial Fibrillation](#)

## Safety

- *Registries:*
  - Schmidt, et al., [European Real-World Outcomes with Pulsed Field Ablation in Patients with Symptomatic Atrial Fibrillation - Lessons from the Multicenter EU-PORIA Registry](#)
  - Turagam, et al., [Safety and Effectiveness of Pulsed Field Ablation to Treat Atrial Fibrillation: One-Year Outcomes From the MANIFEST-PF Registry](#)
- *Autonomic Nervous System:*
  - Del Monte, et al., [Quantitative Assessment of Transient Autonomic Modulation after Single-Shot Pulmonary Vein Isolation with Pulsed-Field Ablation](#)
  - Guo, et al., [Effects of Pulsed Field Ablation on Autonomic Nervous System in Paroxysmal Atrial Fibrillation: A Pilot Study](#)
- *Bronchial:*
  - Fütting, et al., [Bronchial Safety After Pulsed-Field Ablation for Paroxysmal Atrial Fibrillation](#)
- *Cerebral:*
  - Guo, et al., [Effects of Pulsed Field Ablation on Autonomic Nervous System in Paroxysmal Atrial Fibrillation: A Pilot Study](#)
  - Reinsch, et al., [Cerebral Safety after Pulsed Field Ablation for Paroxysmal Atrial Fibrillation](#)
- *Esophageal:*
  - Cochet, et al., [Pulsed Field Ablation Selectively Sparing the Oesophagus During Pulmonary Vein Isolation for Atrial Fibrillation](#)
- *Phrenic Nerve:*
  - Pansera, et al., [Catheter Ablation Induced Phrenic Nerve Palsy by Pulsed Field Ablation—Completely Impossible? A Case Series](#)
- *Stenosis:*
  - Kuroki, et al., [Ostial Dimensional Changes After Pulmonary Vein Isolation: Pulsed Field Ablation vs Radiofrequency Ablation](#)
  - Reddy, et al., [Pulsed Field or Conventional Thermal Ablation for Paroxysmal Atrial Fibrillation](#)

## Real-World Clinical Data

- Chen, et al., [Pulsed Field Ablation-Based Pulmonary Vein Isolation in Atrial Fibrillation Patients with Cardiac Implantable Electronic Devices: Practical Approach and Device Interrogation \(PFA in CIEDs\)](#)
- Ekanem, et al., [Multi-National Survey on the Methods, Efficacy, and Safety on the Post-Approval Clinical Use of Pulsed Field Ablation \(MANIFEST-PF\)](#)
- Fütting, et al., [First Experience with Pulsed Field Ablation as Routine Treatment for Paroxysmal Atrial Fibrillation](#)
- Gunawardene, et al., [Pulsed Field Ablation in Patients with Complex Consecutive Atrial Tachycardia in Conjunction with Ultra-High-Density Mapping: Proof of Concept](#)
- Kueffer, et al., [Validation of a Multipolar Pulsed-Field Ablation Catheter for Endpoint Assessment in Pulmonary Vein Isolation Procedures](#)
- Kueffer, et al., [Pulsed-Field Ablation for the Treatment of Left Atrial Reentry Tachycardia](#)
- Lemoine, et al., [Pulsed-Field Ablation-Based Pulmonary Vein Isolation: Acute Safety, Efficacy and Short-Term Follow-Up in a Multi-Center Real World Scenario](#)
- Magni, et al., [Initial Experience with Pulsed Field Ablation for Atrial Fibrillation](#)
- Schmidt, et al., [SS Study: Safe and Simple Single Shot Pulmonary Vein Isolation with Pulsed Field Ablation Using Sedation](#)

# PUBLICATION LISTING BY TOPIC

## 2023

- Castiglione, et al., [PFA for Treatment of AF in Patients with Congenital Anomalies of Cardiac Veins](#)
- Plank, et al., [Early Recurrences Predict Late Therapy Failure after Pulsed Field Ablation of Atrial Fibrillation](#)
- Ruwald, et al., [Pulsed Field Ablation in Real-World Atrial Fibrillation Patients: Clinical Recurrence, Operator Learning Curve and Re-do Procedural Findings](#)
- Schmidt, et al., [European Real-World Outcomes with Pulsed Field Ablation in Patients with Symptomatic Atrial Fibrillation - Lessons from the Multicenter EU-PORIA Registry](#)
- Tiltz, et al., [Pulsed Field Ablation-Based Pulmonary Vein Isolation Using a Simplified Single-Access Single-Catheter Approach — The Fast and Furious PFA Study](#)
- Turagam, et al., [Safety and Effectiveness of Pulsed Field Ablation to Treat Atrial Fibrillation: One-Year Outcomes From the MANIFEST-PF Registry](#)

## FARAPULSE Versus Other Ablation Modalities

- Blockhaus, et al., [Pulsed Field Ablation for Pulmonary Vein Isolation: Real World Experience and Characterization of the Antral Lesion Size Compared with Cryoballoon Ablation](#)
- Kawamura, et al., [How Does the Level of Pulmonary Venous Isolation Compare Between Pulsed Field Ablation and Thermal Energy Ablation \(Radiofrequency, Cryo, or Laser\)?](#)

## 2023

- Badertscher, et al., [Pulsed-Field Ablation Versus Single Catheter High-Power Short-Duration Radiofrequency Ablation for Atrial Fibrillation: Procedural Characteristics, Myocardial Injury and Midterm Outcomes](#)
- Maurhofer, et al., [Pulsed-Field vs. Cryoballoon vs. Radiofrequency Ablation: A Propensity Score Matched Comparison of One-Year Outcomes after Pulmonary Vein Isolation in Patients with Paroxysmal Atrial Fibrillation](#)
- My, et al., [Acute Lesion Extension Following Pulmonary Vein Isolation with Two Novel Single Shot Devices: Pulsed Field Ablation versus Multielectrode Radiofrequency Balloon](#)
- Urbanek, et al., [Pulsed Field Versus Cryoballoon Pulmonary Vein Isolation for Atrial Fibrillation: Efficacy, Safety, and Long-Term Follow-Up in a 400-Patient Cohort](#)
- Schipper, et al., [Comparison of Pulsed Field Ablation and Cryoballoon Ablation for Pulmonary Vein Isolation](#)
- Wörmann, et al., [Comparison of Pulsed-Field Ablation versus Very High-Power Short Duration-Ablation for Pulmonary Vein Isolation](#)

## Lesion Characterization

- Blockhaus, et al., [Pulsed Field Ablation for Pulmonary Vein Isolation: Real World Experience and Characterization of the Antral Lesion Size Compared with Cryoballoon Ablation](#)
- Bohnen, et al., [Characterization of Circumferential Antral Pulmonary Vein Isolation Areas Resulting from Pulsed-Field Catheter Ablation](#)
- Gunawardene, et al., [Pulsed Field Ablation Combined with Ultra High-Density Mapping in Patients Undergoing Catheter Ablation for AF: Practical and Electrophysiological Considerations](#)
- Kawamura, et al., [Does Pulsed Field Ablation Regress Over Time? A Quantitative Temporal Analysis of Pulmonary Vein Isolation](#)
- Kawamura, et al., [How Does the Level of Pulmonary Venous Isolation Compare Between Pulsed Field Ablation and Thermal Energy Ablation \(Radiofrequency, Cryo, or Laser\)?](#)
- Nakatani, et al., [Pulsed Field Ablation Prevents Chronic Atrial Fibrotic Changes and Restrictive Mechanics After Catheter Ablation for Atrial Fibrillation](#)
- Reddy, et al., [Pulsed Field Ablation of Paroxysmal Atrial Fibrillation: 1-Year Outcomes of IMPULSE, PEFCAT, and PEFCAT II](#)

# PUBLICATION LISTING BY TOPIC

## 2023

- Kueffer, et al., [Pulmonary Vein Isolation Durability and Lesion Regression in Patients with Recurrent Arrhythmia after PulsedField Ablation](#)
- Magni, et al., [Electrophysiological Findings during Re-Do Procedures after Single-Shot Pulmonary Vein Isolation for Atrial Fibrillation with Pulsed Field Ablation](#)
- My, et al., [Acute Lesion Extension Following Pulmonary Vein Isolation with Two Novel Single Shot Devices: Pulsed Field Ablation versus Multielectrode Radiofrequency Balloon](#)
- Ruwald, et al., [Characterization of Durability and Reconnection Patterns at Time of Repeat Ablation after Single-Shot Pulsed Field Pulmonary Vein Isolation](#)
- Schmidt, et al., [European Real-World Outcomes with Pulsed Field Ablation in Patients with Symptomatic Atrial Fibrillation - Lessons from the Multicenter EU-PORIA Registry](#)
- Tohoku, et al., [Findings from Repeat Ablation Using High-Density Mapping after Pulmonary Vein Isolation with Pulsed Field Ablation](#)

## Posterior Wall and/or Mitral Isthmus Ablation

\*Ablation beyond pulmonary vein isolation is outside the use of labeled indication of the FARAWAVE™ PFA Catheter with the FARAPULSE PFA System

## 2023

- Davong, et al., [Pulsed-Field Ablation on Mitral Isthmus in Persistent Atrial Fibrillation - Preliminary Data on Efficacy and Safety](#)
- Gunawardene, et al., [Left Atrial Posterior Wall Isolation with Pulsed Field Ablation in Persistent Atrial Fibrillation](#)
- Ruwald, et al., [Characterization of Durability and Reconnection Patterns at Time of Repeat Ablation after Single-Shot Pulsed Field Pulmonary Vein Isolation](#)
- Sohns, et al., [Lesion Formation Following Pulsed Field Ablation for Pulmonary Vein and Posterior Wall Isolation](#)

## Biomarkers

- Krisai, et al., [Troponin Release After Pulmonary Vein Isolation Using Pulsed Field Ablation Compared to Radiofrequency and Cryoballoon Ablation](#)

## 2023

- Badertscher, et al., [Pulsed-Field Ablation Versus Single Catheter High-Power Short-Duration Radiofrequency Ablation for Atrial Fibrillation: Procedural Characteristics, Myocardial Injury and Midterm Outcomes](#)
- Guo, et al., [Effects of Pulsed Field Ablation on Autonomic Nervous System in Paroxysmal Atrial Fibrillation: A Pilot Study](#)
- Kupusovic, et al. [Visualization of Fibroblast Activation using 68Ga- FAPI PET/CT after Pulmonary Vein Isolation with Pulsed Field Compared with Cryoballoon Ablation](#)
- My, et al., [Acute Lesion Extension Following Pulmonary Vein Isolation with Two Novel Single Shot Devices: Pulsed Field Ablation versus Multielectrode Radiofrequency Balloon](#)

## Clinical Trial Design

- Reddy, et al., [A Randomized Controlled Trial of Pulsed Field Ablation versus Standard-of-Care Ablation for Paroxysmal Atrial Fibrillation: The ADVENT Trial Rationale and Design](#)



## CASE STUDIES

- Adeliño, et al., [Mitral Isthmus Ablation with Pulsed-Field Technology: The Flower Power](#)  
\*Ablation beyond pulmonary vein isolation is outside the use of labeled indication of the FARAWAVE PFA Catheter with the FARAPULSE PFA System
- Adeliño, et al., [Pulsed-Field Ablation of Recurrent Right Atrial Tachycardia: Expanding the Use of Electroporation Beyond Atrial Fibrillation](#)  
\*Ablation beyond pulmonary vein isolation is outside the use of labeled indication of the FARAWAVE PFA Catheter with the FARAPULSE PFA System
- Adragão, et al., [Pulsed-Field Ablation vs Radiofrequency Ablation for Ventricular Tachycardia: First in-Human Case of Histologic Lesion Analysis](#)  
\*Ablation beyond pulmonary vein isolation is outside the use of labeled indication of the FARAWAVE PFA Catheter with the FARAPULSE PFA System
- Ali Ellejmi, et al., [Superior Vena Cava Isolation using a Multielectrode Pulsed Field Ablation Catheter](#)  
\*Ablation beyond pulmonary vein isolation is outside the use of labeled indication of the FARAWAVE PFA Catheter with the FARAPULSE PFA System
- Ascione, et al., [A Posterior Wall Resistant to Electroporation Finally Blocked with Vein of Marshall Ethanol Infusion](#)  
\*Ablation beyond pulmonary vein isolation is outside the use of labeled indication of the FARAWAVE PFA Catheter with the FARAPULSE PFA System
- Bianchini, et al., [Pulsed-Field Ablation of Pulmonary Vein and Left Atrial Posterior Wall Combined with Left Atrial Appendage Occlusion as Single Procedure](#)  
\*Ablation beyond pulmonary vein isolation is outside the use of labeled indication of the FARAWAVE PFA Catheter with the FARAPULSE PFA System
- Chen, et al., [Pulsed Field Ablation as First-Line Treatment to Reduce Atrial Fibrillation Burden Documented by Pacemaker](#)  
\*PRECAUTION: Implantable pacemakers and implantable cardioverter/defibrillators may be adversely affected by irreversible electroporation current
- Chen, et al., [Pulsed Field Ablation as First Line “Efficient” Rhythm Control for Atrial Fibrillation Complicated with Heart Failure: Proof of Concept](#)
- Chen, et al., [Pulsed Field Ablation of Incessant Superior Vena Cava-Triggered Atrial Fibrillation: Watch Out for the Sinoatrial Node](#)  
\*Ablation beyond pulmonary vein isolation is outside the use of labeled indication of the FARAWAVE PFA Catheter with the FARAPULSE PFA System
- Della Rocca, et al., [Transient Inferior ST-Segment Elevation and Ventricular Fibrillation After Cavotricuspid Isthmus Pulsed-Field Ablation](#)  
\*Ablation beyond pulmonary vein isolation is outside the use of labeled indication of the FARAWAVE PFA Catheter with the FARAPULSE PFA System
- Gardziejczyk, et al., [Pulse-Field Ablation using Penta-Spline Catheter as a Bail-Out Strategy for Peri-Mitral Flutter Related to the Left Atrium Anterior Wall Scar](#)  
\*Ablation beyond pulmonary vein isolation is outside the use of labeled indication of the FARAWAVE PFA Catheter with the FARAPULSE PFA System
- Iacopino, et al., [Lesion Effects in Terms of Local Impedance Variations after Pulsed-Field Ablation During Pulmonary Vein Isolation: A Case Report](#)

## CASE STUDIES

- Koruth, et al., [Selective Sparing of Purkinje fibers with Pulsed-Field Myocardial Ablation](#)
- Maury, et al., [Intrapulmonary Haemorrhage during Pulsed Field Ablation](#)
- Maury, et al., [Transient Loss of Capture after Pulse Field Ablation due to Pacing Threshold Elevation](#)  
\*PRECAUTION: Implantable pacemakers and implantable cardioverter/defibrillators may be adversely affected by irreversible electroporation current
- Miraglia, et al., [Unexpected Fused Posterior Wall Lesions after Pulsed-Field Pulmonary Vein Isolation](#)
- Mittal, et al., [Pulsed Field Ablation in Common Inferior Pulmonary Trunk](#)
- Mol, et al., [A Superior Right Jugular Approach to Perform Pulmonary Vein Isolation using FARAPULSE Pulsed-Field Ablation](#)
- Ouss, et al., [First in Human Pulsed Field Ablation to Treat Scar Related Ventricular Tachycardia in Ischemic Heart Disease: A Case Report](#)  
\*Ablation beyond pulmonary vein isolation is outside the use of labeled indication of the FARAWAVE PFA Catheter with the FARAPULSE PFA System
- Ruwald, et al., [Pulsed Field Ablation of the Cavotricuspid Isthmus using a Multispline-Electrode Pulsed Field Ablation Catheter](#)  
\*Ablation beyond pulmonary vein isolation is outside the use of labeled indication of the FARAWAVE PFA Catheter with the FARAPULSE PFA System
- Schmidt, et al., [Single Shot Electroporation of Premature Ventricular Contractions from the Right Ventricular Outflow Tract](#)  
\*Ablation beyond pulmonary vein isolation is outside the use of labeled indication of the FARAWAVE PFA Catheter with the FARAPULSE PFA System
- Tokohu, et al., [Pulsed Field Ablation for Persistent Superior Vena Cava: New Solution for an Old Problem](#)  
\*Ablation beyond pulmonary vein isolation is outside the use of labeled indication of the FARAWAVE PFA Catheter with the FARAPULSE PFA System
- Urbanek, et al., [First Pulse Field Ablation of an Incessant Atrial Tachycardia from the Right Atrial Appendage](#)  
\*Ablation beyond pulmonary vein isolation is outside the use of labeled indication of the FARAWAVE PFA Catheter with the FARAPULSE PFA System

## 2023 CLINICAL PUBLICATIONS

### **Comparison of Pulsed-Field Ablation versus Very High-Power Short Duration-Ablation for Pulmonary Vein Isolation**

Wörmann J, Schipper J, Lüker J, et al.

*Journal of Cardiovascular Electrophysiology* (October, 2023), available [here](#)

- Study that compared the procedural outcome data for PVI between FARAWAVE and very high-power short duration (vHPSD) defined as 70W/7 sec lesions or 70W/5 sec for posterior wall.
- There were 57 patients in each group.
- The FARAWAVE group had significantly shorter procedure duration ( $65 \pm 17$  min) versus the vHPSD ( $95 \pm 23$  min) with longer fluoroscopy times ( $15 \pm 5$  min) vs  $12 \pm 3$  min for vHPSD.
- The freedom from arrhythmia recurrence at a median of 125 days was 80.7% in the FARAWAVE arm versus 77.2% in the vHPSD group.
- Safety event rates were low with 2 tamponades occurring in the FARAWAVE group and 2 groin bleeds in the vHPSD group. One clinically non-significant PV stenosis occurred in the vHPSD group.

---

### **Pulsed-Field Ablation Versus Single Catheter High-Power Short-Duration Radiofrequency Ablation for Atrial Fibrillation: Procedural Characteristics, Myocardial Injury and Midterm Outcomes**

Badertscher P, Weidlich S, Serban T, et al.

*Heart Rhythm* (September, 2023), available [here](#)

- Compared FARAPULSE to high-power short-duration (HPSD) RF looking at efficiency, safety, myocardial injury and midterm outcomes.
- 115 patients (56% paroxysmal) underwent ablation, 52 patients had FARAPULSE ablation and 63 had HPSD RF ablation.
- PFA procedures were significantly shorter (PFA, 58 [53-71] minutes vs HPSD, 83 [71-99] minutes with significantly longer fluoroscopy times (PFA 13 [10-16] minutes vs HPSD 2.2 [1.3-3.6]).
- The postoperative troponin levels were significantly higher in the PFA group (1540 ng/l [1010-1980]) vs HPSD (897 ng/l [725-1240]).
- The AF recurrence free rate at 6 months was 85% for the PFA group and 65% for the HPSD group.
- PFA procedures were shorter, there were higher cardiac troponin levels, and the AF-free survival during mid-term follow-up was similar.

---

### **Quantitative Assessment of Transient Autonomic Modulation after Single-Shot Pulmonary Vein Isolation with Pulsed-Field Ablation**

Del Monte A, Cespon Fenandez M, Vetta G, et al.

*Journal of Cardiovascular Electrophysiology* (September, 2023), available [here](#)

- Assessed the effects of FARAPULSE ablation on the ganglionated plexi and autonomic nervous system (ANS) by looking at the degree of acute vagal modulation induced immediately following FARAPULSE ablation.
- De novo PVI patients treated with FARAPULSE (n=40) or cryoballoon (n=36) were assessed with extracardiac vagal simulation (ECVS) to capture the effects of ablation. To capture any transient effects, the subgroup was assessed before PVI, immediately after PVI and 10 minutes after the last ablation application.

## 2023 CLINICAL PUBLICATIONS

- Baseline values were similar, but the vagal response induced by ECVS almost disappeared in the thermal group but persisted in the FARAPULSE group. Intraprocedural vagal reactions occurred more frequently with FARAPULSE than thermal. The heart rate 24-hour post ablation increased more with thermal than PFA ablation.
- In the subgroup with repeated ANS modulation assessment, PFA had a significant acute suppression of vagal response immediately after ablation which recovered almost completely within a few mins after ablation.
- FARAPULSE was found to be associated with only transitory, short vagal effects on the ANS.

---

### **Left Atrial Posterior Wall Isolation with Pulsed Field Ablation in Persistent Atrial Fibrillation**

Gunawardene M, Frommeyer G, Ellermann C, et al.

\*Ablation beyond pulmonary vein isolation is outside the use of labeled indication of the FARAWAVE PFA Catheter with the FARAPULSE PFA System

*Journal of Clinical Medicine (September, 2023)*, available [here](#)

- Persistent AF patients were treated with PVI + (n=16) or PVI ++ posterior wall isolation (n=59) with FARAWAVE with 32 patients being de novo and 43 patients were repeat ablation patients.
- In the redo cohort, 67% of all PVs were isolated.
- PVI + PWI had an average procedure time of 91 ± 30 min and two minor complications occurred.
- The 354 ± 197-day freedom from atrial arrhythmias (allowing AADs) in the PVI + PWI cohort was 79.3%.
- PWI guided by FARAPULSE had favorable outcomes with a low number of complications.

---

### **Pulsed-Field vs. Cryoballoon vs. Radiofrequency Ablation: A Propensity Score Matched Comparison of One-Year Outcomes after Pulmonary Vein Isolation in Patients with Paroxysmal Atrial Fibrillation**

Maurhofer J, Kueffer T, Madaffari A, et al.

*Journal of Interventional Cardiac Electrophysiology (September, 2023)*, available [here](#)

- CBA and RFA AF patients were propensity matched to PFA, (PFA, n=40), (CBA, n=80) and (RFA, n=80).
- Median procedure times were the shortest with CBA (75 min), followed by PFA (94 min) and RFA (182 min), with RFA having the lowest fluoroscopy dose.
- After 1-year of follow-up, freedom from any atrial arrhythmia was 85% for PFA, 66.2% for CBA, and 73.8% for RFA.
- With propensity matched patients, the results were favorable for the initial use of PFA versus CBA and RFA.

## 2023 CLINICAL PUBLICATIONS

### **Long-Term Clinical Outcomes of Pulsed Field Ablation in the Treatment of Paroxysmal Atrial Fibrillation**

Musikantow D, Neuzil P, Anic A, et al.

*JACC: Clinical Electrophysiology* (September, 2023), available [here](#)

- The first long-term safety and recurrence outcomes for the FARAPULSE PFA system in clinical trial patients.
- 121 PAF patients were treated during these feasibility studies (IMPULSE, PEFCAT, PEFCAT II), of which 49 patients were treated with the optimized waveform (“Biphasic II”). DOI: 10.1016/j.jacep.2021.02.014
- 116 patients were included in long term follow-up with a mean follow-up duration of ~4 years [49 +/- 7 months].
- No new adverse events were reported.
- All Follow-Up Results (Years 1-5) - With the optimized biphasic waveform, there was an 81% (38/47) freedom from AF/AFL recurrence.
- Late Recurrence Follow-Up Analysis (Years 2-5) - 95% freedom from AF/AFL/AT (optimized biphasic waveform).

---

### **Early Recurrences Predict Late Therapy Failure after Pulsed Field Ablation of Atrial Fibrillation**

Plank K, Bordignon S, Urbanek L, et al.

*Journal of Cardiovascular Electrophysiology* (September, 2023), available [here](#)

- 231 AF patients (55% paroxysmal) were analyzed for a medial follow-up of 367 days.
  - 46 (21%) experienced early recurrence of atrial tachyarrhythmia (ERAT) after a median of 23 days post-ablation.
  - The KM estimated freedom from AF/AT was 74.2% at 1 year, 81.8% for paroxysmal and 64.8% for persistent AF.
- Multivariate analysis found that ERAT and female sex were independent predictors of late recurrence.

---

### **Characterization of Durability and Reconnection Patterns at Time of Repeat Ablation after Single-Shot Pulsed Field Pulmonary Vein Isolation**

Ruwald M, Haugdal M, Worck R, et al.

*Journal of Interventional Cardiac Electrophysiology* (September, 2023), available [here](#)

- The pulmonary vein durability rate was 69% in repeat ablation patients (n=26) that had a FARAPULSE procedure an average of 292 ± 119 days after the de novo ablation.
- Patients who underwent posterior wall isolation had a durable PW isolation rate of 80% (4/5).
- Reconnection was observed in the LSPV (27%), LIPV (19%), RSPV (35%), RIPV (42%) with the gaps significantly clustered in the right sided anterior carina compared to other regions.

## 2023 CLINICAL PUBLICATIONS

### **Pulsed Field or Conventional Thermal Ablation for Paroxysmal Atrial Fibrillation**

Reddy VY, Gerstenfeld EP, Natale A, et al.

*New England Journal of Medicine (August, 2023)*, available [here](#), supplement available [here](#)

- The ADVENT Pivotal Trial was the first randomized clinical trial that directly compared FARAPULSE™ PFA to standard-of-care thermal ablation devices (force-sensing radiofrequency (RFA) or cryoballoon ablation (CBA)), for the treatment of paroxysmal atrial fibrillation (PAF).
- It included an experienced group of thermal ablators with limited clinical experience with the novel FARAPULSE technology.
- In this RCT, FARAPULSE demonstrated:
  - Non-inferiority for both the primary safety and effectiveness outcomes compared to thermal ablation technology (posterior probability > .999).
  - Significantly less pulmonary vein cross-sectional narrowing compared to thermal ablation (posterior probability > .999).
  - Significantly shorter procedure times, reduced LA dwell time and total ablation time versus thermal ablation. Lower standard deviations across these characteristics also indicate less variability within the PFA procedures.

---

### **Comparison of Pulsed Field Ablation and Cryoballoon Ablation for Pulmonary Vein Isolation**

Schipper H, Steven D, Lüker J, et al.

*Journal of Cardiac Electrophysiology (August, 2023)*, available [here](#)

- Retrospective analysis of de novo paroxysmal or persistent AF PVI with FARAWAVE (PFA) (n=54) and the POLARx Cryoballoon (CBA) (n=54).
- The total procedure times excluding the LA mapping were significantly shorter for the PFA group (58.0 ± 12.5 min) vs CBA (73.0 ± 24.8 min). Fluoroscopy time was significantly longer in the PFA arm. Subgroup analysis showed a significant reduction in procedure time with continued use of FARAPULSE.
- At 273 ± 129 days, the arrhythmia recurrence free rate was similar for both devices, 74% for PFA and 72% for CBA.
- HR changes between baseline and 3 month follow up did not differ between both groups (PFA: 4 ± 8 beats/min, CBA: 4 ± 11 beats/min).

## 2023 CLINICAL PUBLICATIONS

### **Pulsed Field Ablation-Based Pulmonary Vein Isolation Using a Simplified Single-Access Single-Catheter Approach — The Fast and Furious PFA Study**

Tilz R, Vogler J, Kirstei B, et al.

*Circulation Journal* (August, 2023), available [here](#)

- 50 paroxysmal (56%) and persistent AF patients underwent wide area circumferential ablation (WACA) with FARAPULSE.
  - The mean procedure time was  $27.4 \pm 6.6$  min with a mean LA dwell time of  $14.4 \pm 5.5$  min.
  - The mean time to ambulation was  $3.3 \pm 3.1$  hours with a low rate of periprocedural complications.
  - At a mean follow-up of  $6.5 \pm 2.1$  months, 82% (41/50) patients remained in sinus rhythm.
- .....

### **Pulsed-Field Ablation on Mitral Isthmus in Persistent Atrial Fibrillation - Preliminary Data on Efficacy and Safety**

Davong B, Adeliño R, Delasnerie H, et al.

\*Ablation beyond pulmonary vein isolation is outside the use of labeled indication of the FARAWAVE PFA Catheter with the FARAPULSE PFA System

*JACC: Clinical Electrophysiology* (July, 2023), available [here](#)

- PVI, posterior wall (PW) and mitral isthmus (MI) ablation were performed in 45 patients with persistent AF.
  - The acute success of PVI, PW isolation, and MI block was 100%.
  - There were 2 (4.4%) coronary artery spasms which were reversible after intravenous nitrate infusion.
  - During a mean follow-up of  $107 \pm 59.5$  days, there was a 20% rate of arrhythmia recurrence.
- .....

### **Pulmonary Vein Isolation Durability and Lesion Regression in Patients with Recurrent Arrhythmia after Pulsed Field Ablation**

Kueffer T, Stefanova A, Madaffari A, et al.

*Journal of Interventional Cardiac Electrophysiology* (July, 2023), available [here](#)

- Redo ablation was performed on 29/341 (8.5%) of patients for arrhythmia recurrence.
- At 6-months post index ablation, mapping identified 69/110 (63%) durable PV isolation. In 6 (21%) all PVs were durability isolated.
- PV reconnections were often found on the right sided veins and on the anterior aspects of the upper veins.
- Importantly, only minor regression was observed between the index and redo procedures (median of 3 mm).

## 2023 CLINICAL PUBLICATIONS

### **Acute Lesion Extension Following Pulmonary Vein Isolation with Two Novel Single Shot Devices: Pulsed Field Ablation versus Multielectrode Radiofrequency Balloon**

My I, Lemoine M, Butt M, et al.

*Journal of Cardiovascular Electrophysiology (July, 2023)*, available [here](#)

- Compared lesion formation and lesion extent (measured with mapping and biomarkers) between FARAPULSE and HELIOSTAR (multi-electrode RF balloon).
- 60 paroxysmal patients (28 PFA, 32, RF balloon) underwent PVI, high density mapping and Troponin I was quantified.
- The posterior wall ablation area was significantly larger in the PFA group.
- In a subset of 38 patients, the serum Troponin was significantly higher in the PFA group, likely due to it creating larger lesions.

---

### **Pulsed Field Versus Cryoballoon Pulmonary Vein Isolation for Atrial Fibrillation: Efficacy, Safety, and Long-Term Follow-Up in a 400-Patient Cohort**

Urbanek L, Bordignon S, Schaack D, et al.

*Circulation: Arrhythmia and Electrophysiology (July, 2023)*, available [here](#)

- 400 patients were treated with FARAPULSE (n=200) or cryoballoon ablation (CBA) (n=200).
- The mean procedure times were significantly shorter in the FARAPULSE group (34.5 [29-40] mins) vs CBA (50 [45-60] mins) with similar fluoroscopy times.
- The overall procedural complication rates were 6.5% in the CBA and 3.0% in the FARAPULSE group driven by a higher rate of phrenic nerve palsy in the CBA group.
- The 1-year freedom from arrhythmia recurrence rates in paroxysmal AF were similar with 83.1% in the CBA group and 80.3% in the FARAPULSE group.

---

### **European Real-World Outcomes with Pulsed Field Ablation in Patients with Symptomatic Atrial Fibrillation - Lessons from the Multicenter EU-PORIA Registry**

Schmidt B, Bordignon S, Neven K, et al.

*EURPOACE (July, 2023)*, available [here](#)

- Registry to study the real-world adoption, workflow, acute and long-term outcomes after pulsed field ablation (PFA) in an all-comer atrial fibrillation (AF) patient population in high-volume European centers, inclusive of learning curve.
- This registry demonstrated consistent, short procedure times with a median of 58 minutes despite a large number of operators with varied experience and workflow.
- There was a low rate of safety events (3.6%) and promising one-year efficacy rate (74%) in a large spectrum of AF patients.
- Operator experience and previous primary ablation modality did not have an effect on the one-year AF/AT recurrence rates showing a rapid adoption of the technology by new operators and prior RF and cryo users.
- A small subset of 149 patients (12%) returned for repeat ablation during follow-up. In these patients, EAM revealed a high rate of PVI with 72% of pulmonary veins being durably isolated.



## 2023 CLINICAL PUBLICATIONS

### **Electrophysiological Findings during Re-Do Procedures after Single-Shot Pulmonary Vein Isolation for Atrial Fibrillation with Pulsed Field Ablation**

Magni F, Scherr D, Manninger M, et al.

*Journal of Interventional Cardiac Electrophysiology* (May, 2023), available [here](#)

- Patients who had a de novo procedure with FARAWAVE that had recurrence and subsequent repeat ablation (14/447) procedures were analyzed. The mean time to recurrence was  $4.9 \pm 1.9$  months.
- PV reconnection was found in zero (35.7%), one (21.4%), two (14.3%) or three (28.6%) of patients.
- Durable PVI was observed in over 1/3 of redo patients. The most common arrhythmia recurrence following PVI only was AF. Concomitant (35.7%) or isolated AFL/AT (14.3%) recurrence was observed in 50% of patients.

---

### **Lesion Formation Following Pulsed Field Ablation for Pulmonary Vein and Posterior Wall Isolation**

Sohns C, Fink T, Braun M, et al.

\*Ablation beyond pulmonary vein isolation is outside the use of labeled indication of the FARAWAVE PFA Catheter with the FARAPULSE PFA System

*PACE* (May, 2023), available [here](#)

- Lesion formation was assessed with late gadolinium enhancement CMR (LGE-CMR) 3-months after FARAPULSE ablation.
- In 10 patients, PVI and posterior wall isolation (PWI) was performed with FARAWAVE. The mean procedure duration was  $62 \pm 7$  min with a mean LA dwell time of  $13 \pm 2$  min.
- The mean LA scar burden was  $8.1 \pm 2.1\%$  with a mean scar width of  $12.8 \pm 2.1$  mm. At 7 months, 9/10 (90%) of patients were recurrence free.
- LGE CMR analysis found homogenous and continuous lesion patterns with no evidence of PV stenosis or collateral damage to adjacent structures.

---

### **Safety and Effectiveness of Pulsed Field Ablation to Treat Atrial Fibrillation: One-Year Outcomes From the MANIFEST-PF Registry**

Turagam MK, Neuzil P, Schmidt B, et al.

*Circulation* (May, 2023), available [here](#)

- Multi-national retrospective survey of all patients treated with FARAPULSE from 24 EU centers (77 operators), 1,568 patients.
- Low complication rates; 1.9% major complication rate and 4.0% minor complication rate with no reported esophageal damage or PV stenosis.
- There was an 81.6% 1-year freedom from AF/AFL/AT for paroxysmal AF patients with no difference in recurrence free outcomes based on the procedural volume (PFA procedure numbers).

## 2023 CLINICAL PUBLICATIONS

### **Bronchial Safety After Pulsed-Field Ablation for Paroxysmal Atrial Fibrillation**

Füting A, Reinsch N, Brokkaar L, et al.

*Circulation: Arrhythmia and Electrophysiology (April, 2023)*, available [here](#)

- Respiratory tract CT scans were performed on 60 patients post FARAPULSE ablation to look for bronchial damage with either straight-tip (n=30) or J-tip (n=30) guidewires.
  - In 12/30 patients with the straight-tip, extra-stiff guidewire, small amounts of old blood without active bleeding were detected with no evidence of thermal lesions. There was no clinical relevance at 30 days post-procedure.
  - Use of the straight-tip guidewire may lead to asymptomatic bronchial damage which was not detected when the J-tip guidewire was used.
- .....

### **Pulsed Field Ablation to Treat Atrial Fibrillation: Autonomic Nervous System Effects**

Musikantow DR, Neuzil P, Petru J, et al.

*JACC: Clinical Electrophysiology (April, 2023)*, available [here](#)

- Heart rate was assessed pre and post PVI using FARAPULSE (n=40), Cryoablation (n=40) and radiofrequency (n=40) PVI ablation to understand the impact of pulsed field ablation on the ganglionated plexi (GP).
  - Between baseline and 3 months, heart rates increased by  $8.9 \pm 11.4$  (RF),  $11.1 \pm 9.4$  (CB), and  $-0.1 \pm 9.2$  (PFA) beats/min.
  - Unlike thermal ablation, FARAPULSE PFA had minimal effects on the GPs.
- .....

### **PFA for Treatment of AF in Patients with Congenital Anomalies of Cardiac Veins**

Castiglione A, Küffer T, Gräni C, et al.

\*Ablation beyond pulmonary vein isolation is outside the use of labeled indication of the FARAWAVE PFA Catheter with the FARAPULSE PFA System

*Journal of Cardiovascular Electrophysiology (March, 2023)*, available [here](#)

- Five patients with congenital anomalies were treated with FARAPULSE.
- PVs were isolated with no phrenic nerve palsy or other complications.
- Pre-procedural imaging and 3D mapping was found to be well suited, efficient, and versatile in AF patients with anomalous cardiac veins.

## 2023 CLINICAL PUBLICATIONS

### **Effects of Pulsed Field Ablation on Autonomic Nervous System in Paroxysmal Atrial Fibrillation: A Pilot Study**

Guo F, Wang J, Deng Q, et al.

*Heart Rhythm (March, 2023)*, available [here](#)

- Nerve injury biomarkers and DW-MRI were conducted on 18 patients in a pilot study.
- Serum nerve injury biomarkers did not differ between pre- and post--ablation. Heart rate variability did not differ and there were no acute cerebral microemboli events.
- FARAPULSE PVI did not induce nerve injury in this study.

---

### **Visualization of Fibroblast Activation Using 68Ga- FAPI PET/CT after Pulmonary Vein Isolation with Pulsed Field Compared with Cryoballoon Ablation**

Kupusovic J, Kessler L, Bruns F, et al.

*Journal of Nuclear Cardiology (March, 2023)*, available [here](#)

- Fibroblast activation was used as a surrogate for ablation damage after FARAPULSE(n=15) and CBA (n=11) ablation.
- Fibroblast activation tissue response was less pronounced in the PFA patient cohort vs CBA.

---

### **A Randomized Controlled Trial of Pulsed Field Ablation versus Standard-of-Care Ablation for Paroxysmal Atrial Fibrillation: The ADVENT Trial Rationale and Design**

Reddy VY, Lehmann JW, Gerstenfeld EP, et al.

*Heart Rhythm O2 (March, 2023)*, available [here](#)

- The ADVENT (Randomized Controlled Trial for Pulsed Field Ablation versus Standard of Care Ablation for Paroxysmal Atrial Fibrillation) trial was a multicenter, prospective, single-blind, randomized controlled trial comparing PVI using PFA vs conventional thermal (cryoballoon and contact force radiofrequency) ablation for the treatment of drug-resistant paroxysmal AF.

---

### **Pulsed Field Ablation in Real-World Atrial Fibrillation Patients: Clinical Recurrence, Operator Learning Curve and Re-Do Procedural Findings**

Ruwald MH, Johannessen A, Lock Hansen M, et al.

*Journal of Interventional Cardiac Electrophysiology (February, 2023)*, available [here](#)

- 121 patients underwent PVI with FARAPULSE. The mean procedure time was significantly reduced from the initial cases from 85 ± 34 min to 72 ± 18 min.
- There was one phrenic nerve palsy with partial remission at follow-up. The KM event-free estimate at 365 days was 80% (88% paroxysmal, 69% persistent).
- In 5/8 re-do procedures, the gaps were primarily located in the right pulmonary veins.

## 2022 CLINICAL PUBLICATIONS

### **Pulsed-Field Ablation for the Treatment of Left Atrial Reentry Tachycardia**

Kueffer T, Seiler J, Madaffari A, et al.

\*Ablation beyond pulmonary vein isolation is outside the use of labeled indication of the FARAWAVE PFA Catheter with the FARAPULSE PFA System

*Journal of Interventional Cardiac Electrophysiology (December, 2022)*, available [here](#)

- Left atrial reentry tachycardia were treated with FARAPULSE (n=22).
- Lesion used to treat the ATs included, 20 roof lines, 13 anterior lines, and 6 mitral isthmus lines with no reported complications.

---

### **Findings from Repeat Ablation using High-Density Mapping after Pulmonary Vein Isolation with Pulsed Field Ablation**

Tohoku S, Chun J, Bordignon S, et al.

*EUROPACE (November, 2022)*, available at [here](#)

- In redo patients initially treated with FARAPULSE using the 5S strategy, the incidence of pulmonary vein (PV) reconnection was assessed (inclusive of learning curve).
- Among the 360 patients, 25 patients (19 paroxysmal) underwent a redo procedure in  $6.1 \pm 4$  months.
- The PV durable isolation rate was 90.9% as assessed by high-density mapping.
- The mechanism of all but one atrial tachyarrhythmia was macro-reentry.
- The mean % of isolated posterior wall surface area was  $72.7 \pm 19.0\%$ .
- There was a low rate of PV reconnection (9.1%) in redo patients and the unique features of the FARAWAVE catheter design and optimized workflow enabled wide antral lesion creation without regression over time.

---

### **Pulsed Field Ablation-Based Pulmonary Vein Isolation in Atrial Fibrillation Patients with Cardiac Implantable Electronic Devices: Practical Approach and Device Interrogation (PFA in CIEDs)**

Chen S, Chun J, Bordignon S, et al.

\*PRECAUTION: Implantable pacemakers and implantable cardioverter/defibrillators may be adversely affected by irreversible electroporation current

*Journal of Interventional Cardiac Electrophysiology (November, 2022)*, available [here](#)

- A pilot patient cohort (n=20) underwent PFA ablation for AF (PVI) with different CIEDs.
- CIEDs included pacemaker, implantable cardioverter-defibrillators (ICD), or cardiac resynchronization therapy plus defibrillator (CRT-D).
- CIED pre- and post-PFA interrogation of the devices showed no significant alterations to the parameters or function of the CIEDs and no lead dislodgement.

## 2022 CLINICAL PUBLICATIONS

### **Initial Experience with Pulsed Field Ablation for Atrial Fibrillation**

Magni F, Mulder B, Groenveld H, et al.

*Frontiers in Cardiovascular Medicine* (November, 2022), available [here](#)

- 100 subjects (80% paroxysmal AF) underwent AF ablation with FARAWAVE.
  - The learning curves of 2 operators (junior/senior) who performed >20 procedures showed no difference in procedure time, senior ( $46.9 \pm 9.7$  min) and junior ( $45.9 \pm 9.9$  min).
  - The 2 complications that occurred were bleeding at the access site.
- 

### **Pulsed Field Ablation in Patients with Complex Consecutive Atrial Tachycardia in Conjunction with Ultra-High-Density Mapping: Proof of Concept**

Gunawardene M, Schaeffer B, Jularic M, et al.

\*Ablation beyond pulmonary vein isolation is outside the use of labeled indication of the FARAWAVE PFA Catheter with the FARA-PULSE PFA System

*Journal of Cardiovascular Electrophysiology* (Sept, 2022), available [here](#)

- Fifteen patients with atrial tachycardia (AT) underwent high density mapping to ID critical sites for AT maintenance.
  - FARAWAVE ablation was performed with 100% success, 63% terminated with the first application and 2 ATs in the right atrial requiring RF ablation.
  - No procedure-related complications occurred.
- 

### **Pulsed-Field Ablation-Based Pulmonary Vein Isolation: Acute Safety, Efficacy and Short-Term Follow-up in a Multi-Center Real World Scenario**

Lemoine MD, Fink T, Mencke C, et al.

*Clinical Research in Cardiology* (Sept, 2022), available [here](#)

- 138 patients (62% persistent AF) from 2 centers were treated with FARAWAVE.
- Mean procedure time was  $78 \pm 22$  min including pre- and post-procedure HD voltage mapping. FARAWAVE LA dwell time was  $23 \pm 9$  min with a fluoroscopy time of  $16 \pm 7$  min.
- There were 3 groin complications (2.2%), 1 pericardial tamponade (0.7%) and 1 transient ST-elevation (0.7%).
- The one-year freedom from recurrence rate was 90% in paroxysmal patients (n = 47) and 60% in persistent AF patients (n = 82).

## 2022 CLINICAL PUBLICATIONS

### **Cerebral Safety After Pulsed Field Ablation for Paroxysmal Atrial Fibrillation**

Reinsch N, Fütting A, Höwel D, et al.

*Heart Rhythm (Sept, 2022)*, available [here](#)

- In 30 patients treated with FARAWAVE, Nation Institute of Health Stroke Scale (NIHSS) scores were assessed 2- and 30-days post PVI. One day after PVI, DW-MRI and FLAIR imaging was done to document the occurrence of silent cerebral events (SCE)/silent cerebral lesions (SCL).
- NIHSS scores were 0 for all patients. Cerebral MRI scans were normal in 29/30 (97%) of patients. In one patient (3%), a single cerebral lesion was observed. 40-days post-procedure, a follow-up MRI cerebral scan showed complete lesion regression.

---

### **Catheter Ablation Induced Phrenic Nerve Palsy by Pulsed Field Ablation—Completely Impossible? A Case Series**

Pansera F, Bordignon S, Bologna F, et al.

*European Journal Case Report (Sept, 2022)*, available [here](#)

- Case series on three patients that had FARAWAVE PFA-induced phrenic nerve (PN) injury during PVI. Cases 1 and 3 had PAF without evidence of structural heart disease and case 2 had Pers AF and ischemic cardiomyopathy with preserved ejection fraction.
- Transient right hemidiaphragm palsy was seen during PFA delivery in the RSPV (Cases 1 and 2) and the RIPV (Case 3).
- The palsy lasted < 1 min and was followed by spontaneous full recovery in all cases (Case 1, 40 sec, Cases 2 and 3 lasted a few seconds).
- Transient PN palsy fully recovered rapidly suggesting PN hyperpolarization of neuronal cells or depletion of acetylcholine in the motoric endplate. Further studies are needed to understand the mechanism.

---

### **Multi-National Survey on the Methods, Efficacy, and Safety on the Post-Approval Clinical Use of Pulsed Field Ablation (MANIFEST-PF)**

Ekanem E, Reddy VY, Schmidt B, et al.

*Europace (August, 2022)*, available [here](#)

- The MANIFEST-PF registry was a retrospective survey of 24 centers with 90 operators, 1758 patients that assessed the real-world performance (use case, acute effectiveness, safety) of FARAPULSE.
- Procedure time was 65 min, fluoroscopy time was 13.7 min. There was a 99.9% mean acute PVI success rate.
- There were no esophageal complications reported, no phrenic nerve injury persisting beyond hospital discharge and no reported PV stenosis. There was a 1.6% rate of major complications, a 3.87% rate of minor complications and 0.46% rate of energy specific adverse events.
- Root cause analysis showed that most of the pericardial tamponades and stroke were attributable to catheter workflow and manipulation, independent of energy modality. Complications were plotted on a timeline, and it indicated an improvement in complication rate over time.

## 2022 CLINICAL PUBLICATIONS

### **Pulsed Field Ablation for Pulmonary Vein Isolation: Real-World Experience and Characterization of the Antral Lesion Size Compared with Cryoballoon Ablation**

Blockhaus C, Guelker J, Feyen L, et al.

*Journal of Interventional Cardiac Electrophysiology* (August, 2022), available [here](#)

- Single-center study looking at procedural characteristics and the size of acute PVI antral lesions with high-density mapping in 43 patients treated with PFA compared to 20 patients treated with cryoballoon ablation.
- All patients had 100% acute vein isolation with no early reconnections. The acute antral lesion size of PFA lesions ( $67.03 \pm 12.69\%$ ) were significantly larger compared to cryoballoon ( $57.39 \pm 10.91\%$ ).
- In the PFA group there was no acute phrenic nerve injury, and 1 (4.34%) patient stroke.

---

### **Validation of a Multipolar Pulsed-Field Ablation Catheter for Endpoint Assessment in Pulmonary Vein Isolation Procedures**

Kueffer T, Baldinger S, Servatius H, et al.

*EUROPACE* (June, 2022), available [here](#)

- In 56 patients undergoing PVI with FARAWAVE, the accuracy of FARAWAVE to detect residual PV connections was assessed with high-density mapping.
- Acute PVI was achieved in 100% of PVs.
- The accuracy of the PV assessment with FARAWAVE was 91%. In 14/213 (6.6% of veins), FARAWAVE incorrectly indicated residual PV conduction due to high-output pace-capture.
- Lowering the output to 5 V/1 ms reduced this observation to 0.9% (2/213) and increased the accuracy to 97%.
- FARAWAVE offered reliable endpoint assessment for PVI and lowering the pacing output increased the accuracy from 91% to 97%.
- At a median of 3.2 months, 3/56 (5.4%) underwent a redo procedure. The durable PV isolation rate was 10/12 (83%).

---

### **5S Study: Safe and Simple Single Shot Pulmonary Vein Isolation with Pulsed Field Ablation Using Sedation**

Schmidt B, Bordignon S, Tohoku S, et al.

*Circulation: Arrhythmia and Electrophysiology* (June, 2022), available [here](#)

- Single-center study looking at the adoption and the process of streamlining the procedure in the first 191 patients treated with FARAPULSE PFA. Electrogram validation was performed with a circular mapping catheter (CMC) in the first 25 patients, cerebral MRI was performed in 53 patients and esophageal endoscopy was performed in 52 patients.
- Electrogram information was 100% congruent between the CMC and FARAWAVE. PVI rate was 100%. No esophageal temperature rise or esophageal thermal injuries were observed. Two minor strokes occurred in the first 25 patients, likely due to air embolism during catheter exchanges.
- After the first 25 patients, the procedure times were significantly reduced from an average of  $46 \pm 14$  min to  $38 \pm 13$  min. During short term follow-up, 9% (17/191) of patients had atrial arrhythmia recurrence.

## 2022 CLINICAL PUBLICATIONS

### **Characterization of Circumferential Antral Pulmonary Vein Isolation Areas Resulting from Pulsed-Field Catheter Ablation**

Bohnen M, Weber R, Minners J, et al.

*Europace (June, 2022)*, available [here](#)

- In 40 patients, pre- and post-procedure 20-pole circular mapping catheter voltage mapping was done to evaluate PV isolation and area of isolation.
- Isolation gaps were located most frequently in the anterior antral PV segments of the left PVs.
- Additional areas of isolation beyond the antral PV segments were found on the posterior wall and roof regions.

---

### **First Experience with Pulsed Field Ablation as Routine Treatment for Paroxysmal Atrial Fibrillation**

Füting A, Reinsch N, Höwel D, et al.

*Europace (May, 2022)*, available [here](#)

- Single-center 30 patient study looking at phrenic nerve injury and high-density mapping pre-and post-ablation.
- Acute PVI rate was 100%, the median procedure time was 116 min and the FARAWAVE catheter dwell time was 29 min. There was no esophageal or phrenic nerve injury.
- 97% of patients were in sinus rhythm after 90 days.

---

### **Troponin Release After Pulmonary Vein Isolation Using Pulsed Field Ablation Compared to Radiofrequency and Cryoballoon Ablation**

Krisai P, Knecht S, Badertscher P, et al.

*Heart Rhythm (May, 2022)*, available [here](#)

- Troponin T is a measure of myocardial cell death. Troponin T was measured in 60 patients one day before and the morning after PVI ablation with FARAWAVE, radiofrequency or cryoballoon ablation. No additional lesion sets were performed.
- Post-procedure Troponin T levels with PFA were 1.6x and 1.9x higher vs. RF and Cryo, respectively with no significant difference between the RF and cryo groups.

---

### **Pulsed Field Ablation Combined with Ultra-High-Density Mapping in Patients Undergoing Catheter Ablation for Atrial Fibrillation: Practical and Electrophysiological Considerations**

Gunawardene M, Schaeffer B, Jularic M, et al.

*Journal of Cardiovascular Electrophysiology (March, 2022)*, available [here](#)

- 20 consecutive patients underwent PVI with FARAWAVE. Additional ablations were performed off-label in a sub-set of patients. PFA lesion size and decrease in voltage were assessed with high-density voltage mapping.
- High density mapping showed PV reconnection in 5 cases (6.25%). Gaps were located at the anterior-superior PV ostia and were successfully closed with additional PFA. Voltage was significantly decreased following PFA with almost no complex electrogram fractionation at the lesion border zones.
- High-density mapping for FARAWAVE PFA lesion showed wide, antral, circumferential lesion with significantly decreased atrial tissue voltage and little evidence of fraction in the lesion border zones.



## 2021 CLINICAL PUBLICATIONS

### **Does Pulsed Field Ablation Regress Over Time? A Quantitative Temporal Analysis of Pulmonary Vein Isolation**

Kawamura I, Neuzil P, Shivamurthy P, et al.

*Heart Rhythm (June, 2021)*, available [here](#)

- Patients with PAF underwent PVI with FARAWAVE. A comparison of voltage maps immediately after PFA and at a median of 84 days (interquartile range 69–90 days) later revealed that there was no significant difference in either the left and right-sided PV antral isolation areas or nonablated posterior wall area.
- The distances between low-voltage edges on the posterior wall were also not significantly different between the 2 time points.
- The level of PV antral isolation after PFA with FARAWAVE persisted without regression.

---

### **Pulsed Field Ablation Prevents Chronic Atrial Fibrotic Changes and Restrictive Mechanics After Catheter Ablation for Atrial Fibrillation**

Nakatani Y, Sridi-Cheniti S, Cheniti G, et al.

*Europace (May, 2021)*, available [here](#)

- Cardiac magnetic resonance was performed pre-ablation, acutely (< 3 h), and 3 months post-ablation in 41 patients with PAF undergoing PVI with PFA (n = 18) or thermal ablation (n = 23, 16 radiofrequency ablations, 7 cryoballoon ablations).
- Tissue changes were more homogeneous after PFA than after thermal ablation, with no sign of microvascular damage or intramural hemorrhage. In the chronic stage, the majority of acute LGE had disappeared after PFA, whereas most LGE persisted after thermal ablation.
- The maximum strain on PV antra, the LA expansion index, and LA active emptying fraction declined acutely after both PFA and thermal ablation but recovered at the chronic stage only with PFA.
- In this study, PFA induced large acute LGE lesions which mostly disappeared in the chronic stage, suggesting a reparative process involving less chronic fibrosis.

---

### **Pulsed Field Ablation of Paroxysmal Atrial Fibrillation: 1-Year Outcomes of IMPULSE, PEFCAT, and PEFCAT II**

Reddy VY, Dukkipati SR, Neuzil P, et al.

*JACC-EP (May, 2021)*, available [here](#)

- In 3 multicenter studies (IMPULSE, PEFCAT and PEFCAT II), PAF patients underwent PVI using a basket and flower PFA catheter.
- Invasive remapping was performed at 2 to 3 months, and reconnected PVs were reisolated with PFA or radiofrequency ablation. After a 90-day blanking period, arrhythmia recurrence was assessed over 1-year follow-up.
- In 121 patients, acute PVI was achieved in 100% of PVs with PFA alone.
- PV remapping, performed in 110 patients at  $93.0 \pm 30.1$  days, demonstrated durable PVI in 84.8% of PVs (64.5% of patients), and 96.0% of PVs (84.1% of patients) treated with the optimized biphasic energy PFA waveform.
- The 1-year Kaplan-Meier estimates for freedom from any atrial arrhythmia for the entire cohort and for the optimized biphasic energy PFA waveform cohort were  $78.5 \pm 3.8\%$  and  $84.5 \pm 5.4\%$ , respectively.

## 2021 CLINICAL PUBLICATIONS

### **How Does the Level of Pulmonary Venous Isolation Compare between Pulsed Field Ablation and Thermal Energy Ablation (Radiofrequency, Cryo, or Laser)?**

Kawamura I, Neuzil P, Shivamurthy P, et al.

*Europace (May, 2021)*, available [here](#)

- In a clinical trial (NCT03714178), PAF patients under-went PVI with FARAWAVE using a biphasic waveform, and after 75 days, detailed voltage maps were created.
- Comparative voltage mapping data were retrospectively collected from consecutive PAF patients who (i) underwent PVI using thermal energy, (ii) underwent re-ablation for recurrence, and (iii) had durably isolated PVs. The left and right PV antral isolation areas and non-ablated posterior wall were quantified.
- There was no significant difference between the PFA and thermal ablation cohorts in either the left- and right-sided PV isolation areas, or the non-ablated posterior wall area.

---

### **Pulsed Field Ablation Selectively Spares the Oesophagus During Pulmonary Vein Isolation for Atrial Fibrillation**

Cochet H, Nakatani Y, Sridi-Cheniti S, et al.

*Europace (February, 2021)*, available [here](#)

- Cardiac magnetic resonance (CMR) imaging was performed before, acutely (< 3 h) and 3 months post-ablation in 41 PAF patients undergoing PVI with PFA (N = 18, FARAPULSE) or thermal methods (N = 23, 16 radiofrequency, 7 cryoballoon).
- Esophageal and aortic injuries were assessed by using late gadolinium-enhanced (LGE) imaging. Phrenic nerve injuries were assessed from diaphragmatic motion on intra-procedural fluoroscopy.
- Acutely, thermal methods induced high rates of esophageal lesions (43%), all observed in patients showing direct contact between the esophagus and the ablation sites.
- Esophageal lesions were observed in no patient ablated with PFA (0%, P < 0.001 vs. thermal methods), despite similar rates of direct contact between the esophagus and the ablation sites (P = 0.41).
- Acute lesions were detected on CMR on the descending aorta in 10/23 (43%) after thermal ablation, and in 6/18 (33%) after PFA (P = 0.52). CMR at 3 months showed a complete resolution of esophageal and aortic LGE in all patients.

---

### **Pulsed Field Ablation: A Promise That Came True**

Ante A, Breskovic T, Sikiric I.

*Current Opinion in Cardiology (Jan, 2021)*, available [here](#)

- Pulsed field ablation is a nonthermal ablative modality that uses short living, strong electrical field created around catheter to create microscopic pores in cell membranes (electroporation). When adequately dosed/ configured it shows a preference for myocardial tissue necrosis.
- First in human series using pulsed field ablation for atrial fibrillation ablation have been completed and data published for several platforms. Acute safety outcomes are similar across the platforms with a low complication rate for complications typically reported for thermal ablation methods (esophageal injury, pulmonary vein stenosis, phrenic nerve palsy).
- Promising acute data on pulmonary vein isolation had been corroborated with satisfactory 1-year clinical follow-up for a single platform (i.e., FARAPULSE), whereas reports are pending for the rest. Research efforts are being expanded to a development of focal catheters, and therefore, pulsed field ablation application for ventricular arrhythmias.

## 2020 CLINICAL PUBLICATIONS

### **Pulsed Field Ablation in Patients with Persistent Atrial Fibrillation**

Reddy VY, Anic A, Koruth J, et al.

*JACC (Sept, 2020)*, available [here](#)

- PersAFOne was a single-arm study evaluating biphasic, bipolar PFA with FARAWAVE for PVI and LAPW ablation to assess the safety and lesion durability of pulsed field ablation (PFA) for both PVI and LAPW ablation in persistent AF.
  - In 25 patients, acute PVI (96 of 96 pulmonary veins) were 100% acutely successful with the FARAWAVE catheter. Using the focal PFA catheter, acute cavotricuspid isthmus block was achieved in 13 of 13 patients.
  - Post-procedure EGD and repeat cardiac computed tomography revealed no mucosal lesions or PV narrowing, respectively.
  - Invasive remapping at 2 to 3 months demonstrated durable isolation (defined by entrance block) in 82 of 85 PVs (96%) and 21 of 21 LAPWs (100%) treated with the pentaspline catheter.
- 

### **Ostial Dimensional Changes After Pulmonary Vein Isolation: Pulsed Field Ablation vs. Radiofrequency Ablation**

Kuroki K, Whang W, Eggert C, et al.

*Heart Rhythm (May, 2020)*, available [here](#)

- Data were analyzed from 4 PAF ablation trials using either PFA or RFA.
- Baseline and 3-month cardiac computed tomography scans were reconstructed into 3-dimensional images, and the long and short axes of the PV ostia were quantitatively and qualitatively assessed in a randomized blinded manner.
- PV ostial diameters decreased significantly less with PFA than with RFA (% change; long axis:  $0.9\% \pm 8.5\%$  vs.  $-11.9\% \pm 16.3\%$ ;  $P < .001$  and short axis:  $3.4\% \pm 12.7\%$  vs.  $-12.9\% \pm 18.5\%$ ;  $P < .001$ ).
- PV narrowing/stenosis was present in 0% and 0% vs. 12.0% and 32.5% of PVs and patients who underwent PFA and RFA, respectively.
- In this study, unlike after RFA, the incidence and severity of PV narrowing/stenosis after PV isolation was virtually eliminated with PFA.

## 2019 CLINICAL PUBLICATIONS

### **Pulsed Field Ablation for Pulmonary Vein Isolation in Atrial Fibrillation**

Reddy VY, Neuzil P, Koruth JS, et al.

*JACC (July, 2019)*, available [here](#)

- Two trials were conducted to determine whether PFA allows durable pulmonary vein (PV) isolation without damage to collateral structures, in patients with PAF.
  - Ablation was performed using proprietary bipolar PFA waveforms: either monophasic with general anesthesia and paralytics to minimize muscle contraction, or biphasic with sedation because there was minimal muscular stimulation. No esophageal protection strategy was used. Invasive electrophysiological mapping was repeated after 3 months to assess the durability of PV isolation.
  - 81 patients, all PVs were acutely isolated by monophasic (n = 15) or biphasic (n = 66) PFA. With successive waveform refinement, durability at 3 months improved from 18% to 100% of patients with all PVs isolated. Beyond 1 procedure-related pericardial tamponade no additional primary adverse events over the 120-day median follow-up, including: stroke, phrenic nerve injury, PV stenosis, and esophageal injury.
- 

## 2018 CLINICAL PUBLICATIONS

### **Ablation of Atrial Fibrillation with Pulsed Electric Fields**

Reddy VY, Koruth J, et al.

*JACC-EP (April, 2018)*, available [here](#)

- The first acute clinical experience of AF ablation with PFA, both epicardial box lesions during cardiac surgery, and catheter-based PVI.
- PFA was performed using a custom over-the-wire endocardial catheter for percutaneous transseptal PV isolation, and a linear catheter for encircling the PVs and posterior left atrium during concomitant cardiac surgery.
- Catheter PV ablation was successful in 15 patients (100%) 57 PVs Using 3.26 lesions/PV and surgical box lesions were successful in 6 of 7 patients (86%) 2 lesions/patient. No complications.

## 2023 PRECLINICAL PUBLICATIONS

### **Electrophysiology, Pathology, and Imaging of Pulsed Field Ablation of Scarred and Healthy Ventricles in Swine**

Kawamura I, Reddy V, Santos-Gallego C, et al.

*Circulation: Arrhythmia and Electrophysiology* (January, 2023), available [here](#)

- 6 swine were infarcted to assess penetration of scar, risk of arrhythmias and lesion imaging evaluation.
- FARAPULSE PFA successfully penetrated scar without significant differences in the lesion depth of infarcted tissue ( $5.9 \pm 1.0$  mm) vs healthy ( $5.7 \pm 1.3$  mm) myocardium.
- In ungated QRS PFA applications, sustained ventricular arrhythmias requiring defibrillation occurred in 4/187 (2.1%) applications with zero occurring during gated applications.
- Dark-blood late-gadolinium-enhanced sequences allowed for improved endocardial border detection.

---

## 2022 PRECLINICAL PUBLICATIONS

### **Effect of Epicardial Pulsed Field Ablation Directly on Coronary Arteries**

Higuchi S, Im S, Stillson C, et al.

*JACC: Clinical Electrophysiology* (Dec, 2022), available [here](#)

- 4 swine, FARAWAVE lesions were delivered directly to the left anterior descending artery, left circumflex artery or normal myocardium.
- Angiography was performed to quantify the degree of coronary artery narrowing and histology was performed at 4 and 8 weeks.
- Acute luminal narrowing immediately after PFA was 47% which gradually resolved over 30 minutes.
- Epicardial lesions had a median depth of 4.1 mm and 87.5% of the arteries had minimal to mild stenosis via neointimal hyperplasia.

---

### **Pulsed Field Ablation of Left Ventricular Myocardium in a Swine Infarct Model**

Im S, Higuchi S, Lee A, et al.

*JACC: Clinical Electrophysiology* (June, 2022), available [here](#)

- 10 swine were infarcted to evaluate how PFA and RF perform in areas of myocardial scar.
- In myocardial scar, lesion depth was not different between the FARAWAVE or the FOCAL PFA catheter.
- In myocardial scar, lesion depth was significantly greater for PFA vs. RF.
- In a pre-clinical animal model, unlike RF, FARAPULSE PFA was able to effectively ablate surviving islands of myocardium in infarct-related ventricular substrate.

## 2020 PRECLINICAL PUBLICATIONS

### **Pulsed Field Ablation vs Radiofrequency Ablation: Esophageal Effects in a Novel Preclinical Model**

Koruth JS, Kuroki K, Kawamura I, et al.

*Circulation: Arrhythmia and Electrophysiology* (January, 2020), available [here](#)

- A novel preclinical model was created to nonsurgically assess the response to esophageal injury. This was accomplished by delivering the energy source from within the inferior vena cava, against the esophagus (which was purposefully mechanically deviated towards the IVC).
  - Biphasic pulsed field ablation induced no chronic histopathologic esophageal changes, whereas radiofrequency catheter ablation demonstrated a spectrum of esophageal lesions including esophageal ulcers, abscess, and fistula.
- 

## 2019 PRECLINICAL PUBLICATIONS

### **Preclinical Evaluation of Pulsed Field Ablation: Electrophysiological and Histological Assessment of Thoracic Vein Isolation**

Koruth JS, Kuroki K, Iwasawa J, et al.

*Circulation: Arrhythmia and Electrophysiology* (December, 2019), available [here](#)

- In this study, the safety, efficacy, and durability of achieving catheter-based electrical isolation of PVI using optimized monophasic and biphasic PFA waveforms and describe procedural and histological characteristics of PFA in swine atrial tissue.
  - Both waveforms created confluent myocardial lesions that demonstrated a myocardial-specific ablative effect.
  - Biphasic PFA was more durable than monophasic PFA and radiofrequency ablation lesions.
- 

### **Endocardial Ventricular Pulsed Field Ablation: A Proof-of-Concept Preclinical Evaluation**

Koruth JS, Kuroki K, Iwasawa J, et al.

*EP Europace* (December, 2019), available [here](#)

- Assessment of safety and feasibility of FARAPULSE PFA in swine ventricles with a prototype steerable endocardial catheter.
- Gross measurements, available for 28 of 30 ablation sites, revealed average lesion dimensions to be  $6.5 \pm 1.7$  mm deep and  $22.6 \pm 4.1$  mm, with a maximum depth and width of 9.4 mm and 28.6 mm respectively. In PFA lesions, fibrous tissue homogeneously replaced myocytes. When present in the lesion zone, nerve fascicles and vasculature were preserved.



**Cardiology**

300 Boston Scientific Way  
Marlborough, MA 01752-1234  
**[www.bostonscientific.com](http://www.bostonscientific.com)**

*Medical Professionals:*  
1.800.CARDIAC (227.3422)  
*Customer Service:*  
1.888.272.1001

© 2023 Boston Scientific Corporation  
or its affiliates. All rights reserved.

EP-1426609-AC

CAUTION: Investigational device. Limited by Federal (or US) law to investigational use only. Not available for sale.  
Product not FDA approved. Not approved for sale in the US. For educational purposes only.