


ANAESTHESIA CHAPTER

# EMBLEM™ MRI S-ICD SYSTEM

## Implant Best Practices

---

EMA Version 2



## Important information

Prior to use, please review Physician Manual and Directions for Use. These materials are intended to describe common clinical considerations and procedural steps for the use of referenced technologies, but may not be appropriate for every patient or case. Decisions surrounding patient care depend on the physician's professional judgment in light of all available information for the case at hand.



Find all relevant S-ICD Best Practices videos on EDUCARE

Your personalised education and training platform:

**[educare.bostonscientific.eu](https://educare.bostonscientific.eu)**

*Boston Scientific extend their thanks to Dr Reinoud Knops, Prof. Jürgen Kuschyk and Dr Andrea Droghetti for their support, and expertise in the development of this document. Also to Dr Daniel Rodríguez Muñoz for providing additional S-ICD implant images, and Dr Michael Essandoh for review of the anaesthesia content.*

↓ To access the full S-ICD Implant Best Practices guide please [click here](#)

## 1 | S-ICD Labelling and Recommendations

### 1.1 User Manuals

## 2 | S-ICD Learning Curve

### 2.1 Improving Outcomes

## 3 | Commonly Used Accessories

### 3.1 Checklist

## 4 | Pre-implant Considerations

### 4.1 Checklist

## 5 | Screening

### 5.1 Positioning the Surface ECG Electrodes

### 5.2 Obtaining a Stable ECG Baseline

## 6 | System Placement

### 6.1 X-Ray Landmarks and System Placement

### 6.2 Fluoroscopy

### 6.3 Mark-Up

### 6.4 Female-Specific Considerations

## 7 | Anaesthesia Options when Implanting S-ICD

### 7.1 Anaesthesia Management

### 7.2 Monitored Anaesthesia Care (MAC)

### 7.3 Non-Anaesthesia Administered Sedation and Analgesia (NASA)

### 7.4 Truncal Plane Blocks

## 8 | Draping and Patient Preparation

### 8.1 Patient Prep

### 8.2 Draping

### 8.3 Arm Positioning

## 9 | Implant

### 9.1 Initial Incision and Pocket Placement

⬇ To access the full S-ICD Implant Best Practices guide please [click here](#)

- 9.2 Intermuscular Device Placement
- 9.3.1 Xiphoid Incision – Process and Position
- 9.3.2 Xiphoid Incision – Width and Depth
- 9.4 Lateral Tunnel – Two-Incision Technique
- 9.5.1 Superior Tunnel – Process
- 9.5.2 Superior Tunnel – Position
- 9.6.1 Electrode Position
- 9.6.2 Electrode Positioning – Examples
- 9.7 Securing Electrode at Xiphoid
- 9.8 Connecting the Electrode to the Can
- 9.9 Positioning Can in the Pocket
- 9.10 Closing

## 10 | VF Induction

- 10.1 General Considerations
- 10.2.1 Recommended Testing Protocol
- 10.2.2 Ensuring Proper Device Position
- 10.2.3 Ensuring Electrode Position

## 11 | S-ICD Device Replacement

- 11.1 Improving Outcomes
- 11.2 Final Programming

## 12 | Challenges

- 12.1 High BMI Patients
- 12.2 Additional Tips for High BMI Patients
- 12.3 Considerations for Low BMI Patients
- 12.4 High DFT/Failed DFT
- 12.5 Air Entrapment

## 13 | Post-implant Consideration

- 13.1 Post-Operative Care


## 14 | References

- 14.1 Bibliography

# 7

---

## Anaesthesia Options when Implanting S-ICD



## 7 | Anaesthesia Options when implanting S-ICD

### 7.1 Anaesthesia Management


An expert panel recommended several anaesthesia options and best practices that allow implant of the S-ICD without the need for general anaesthesia.<sup>14</sup> In addition to adequate sedation and pain control during the implant itself, these best practices provide superior post-operative pain control and in many instances allow patients to be discharged home the same day of surgery.<sup>14,15</sup>

S-ICD pocket creation requires substantial dissection within the highly innervated intercostal region (T2-T9); this pocket creation, together with tunnelling and defibrillation testing requires sedation and pain control practices that are different to those used for a transvenous ICD implant.

#### TIP

Partner with the anaesthesia team to determine the best anaesthesia approach to utilise during the learning curve phase (first 5-10 S-ICD implants).<sup>14</sup>

If local anaesthesia options are utilised, it is important for the implanter to communicate with the anaesthesia provider the parts of the implant where there is a need to have deeper sedation: creating the pocket, tunnelling, and defibrillation testing.



#### General Anaesthesia

If general anaesthesia is employed, consider also using a long acting local anaesthetic, as this has been shown to reduce the need for post-operative opioids.<sup>14,16-18</sup>

Using general anaesthesia may increase costs and case time, and can cause clinically significant bradycardia and hypotension in patients with heart failure, and may prolong post-op recovery.<sup>14</sup>

#### VIDEOS

For further information or to see related procedural videos please visit our EDUCARE website by clicking [here](#).

## 7 | Anaesthesia Options when implanting S-ICD

### 7.2 Monitored Anaesthesia Care (MAC)


Monitored anaesthesia care (MAC) is a service provided by an anaesthesiologist or nurse anaesthetist (not general anaesthesia).

The benefits of using MAC during the S-ICD implant is that it allows the implanter to completely focus on the technical aspects of the S-ICD implant while the anaesthesia provider controls sedation and pain, and manages haemodynamics during the most stimulating points in the procedure.

This diagram displays a list of options to control post-operative pain. Some of these medications are given preoperatively, such as preoperative acetaminophen or gabapentin. Some medications are given during the implant (e.g. the long acting local anaesthetic bupivacaine), but last longer and help to prevent pain in the postoperative period. Experts believe preventing pain is a more effective approach than treating post-implant pain.<sup>14</sup>

#### Postoperative Pain Control

*Enhanced Recovery for S-ICD Patients*



## 7 | Anaesthesia Options when implanting S-ICD

### 7.2 Monitored Anaesthesia Care (MAC) (continued)

#### BEST PRACTICES



When utilised with a long acting local anaesthetic, such as bupivacaine, and pre-operative acetaminophen, patients can often go home the same day as the S-ICD implant.<sup>15</sup>

#### TIP



In place of 1% lidocaine, a long acting local anaesthetic (0.25% bupivacaine) at a maximum dose of 2.5 mg/kg may be used. Irrespective of body weight, the maximum dose of bupivacaine administered should NOT exceed 175 mg. Intralipid 20% should be available for intravenous infusion in the event of local anaesthetic systemic toxicity.<sup>14</sup>

While MAC combined with long-acting local anaesthetic and other pre-operative and post-operative multi-modal therapies is safe and very effective for S-ICD implants, it may not be the ideal anaesthesia protocol for S-ICD implants when anaesthesia services are not available.

#### VIDEOS



**For further information or to see related procedural videos please visit our EDUCARE website by clicking [here](#).**



## 7 | Anaesthesia Options when implanting S-ICD

### 7.3 Non-Anaesthesia Administered Sedation and Analgesia (NASA)

While NASA may be more convenient to utilise, care must be taken to assure patient safety due to the risks associated with respiratory or cardiovascular complications. Especially when utilising non-reversible sedative, and hypnotic agents, such as Propofol (which has a narrow therapeutic window, and can rapidly progress from sedation to general anaesthesia).

#### BEST PRACTICES



- The nurse responsible for patient sedation and monitoring should be exempt from additional duties, to ensure rapid recognition and treatment of sedation-related complications.
- Patient selection must be carefully considered.
- The use of reversible sedatives and analgesics such as midazolam and opioids should be considered.
- As with a TV-ICD implant, consider GA or MAC versus NASA when the patient:
  - Is morbidly obese, has obstructive sleep apnea, severe heart failure, severe pulmonary hypertension, and/or has a challenging airway.

- Has psychiatric illness for which psychoactive drugs are prescribed or opiates for chronic pain, as both scenarios may lead to increased sedation and analgesia requirements.

- It is important that the NASA approach only be considered once the implanter has acquired sufficient S-ICD implant experience, as NASA requires that the sedation be administered by the implanter, or by a nurse under the direct supervision of the implanter.
- The use of a long acting local anaesthetic such as 0.25% bupivacaine in the pocket and preoperative acetaminophen and/or oxycodone is recommended when utilising the NASA approach.<sup>14</sup>

#### TIP



For implanters who have deep sedation privileges at their institution, the use of small bolus doses of etomidate 0.1 mg/kg is recommended over the use of propofol due to the potential for cardiovascular and respiratory depression with propofol.<sup>14</sup>

#### VIDEOS



**For further information or to see related procedural videos please visit our EDUCARE website by clicking [here](#).**

## 7 | Anaesthesia Options when implanting S-ICD

### 7.4 Truncal Plane Blocks

Truncal Plane Blocks, including the serratus anterior plane block (SAPB) are quick and easy to administer (5 – 10 minutes),<sup>19</sup> and should be administered up to a minimum of 40 minutes before the implant procedure.


SAPB provides long acting anaesthesia (480 minutes) to the lateral branches of intercostal nerves from the T2 to T9 level with excellent results.<sup>16,17,19-21</sup>

#### EQUIPMENT NEEDED<sup>21</sup>


- 50 – 120 mm, 21 – 22 gauge needle
- High-frequency linear ultrasound (US) probe

1. With the patient in a supine position, starting at the clavicle, count the ribs inferiorly and laterally, down to the fifth rib in the midaxillary line. The latissimus dorsi (superficial and posterior), teres major (superior) and serratus muscles (deep and inferior) are then easily identifiable by ultrasound overlying the fifth rib.
2. Infiltrate the skin at the puncture site with local anaesthesia.
3. Insert the needle into the plane under US guidance with the probe in an anteroposterior direction.


a. Ultrasound-guided SAPB<sup>21</sup>



b. Injection of anaesthesia into the superficial interfascial plane between latissimus dorsi and serratus anterior<sup>21</sup>



c. Injection of anaesthesia between serratus anterior and intercostal muscles<sup>21</sup>



#### TIP



Ensure that the needle is in the same plane as the probe so that it can be visualised.

## 7 | Anaesthesia Options when implanting S-ICD

### 7.4 Truncal Plane Blocks (continued)


#### BEST PRACTICES




Communication with the patient is important, particularly when using local anaesthesia techniques, as it puts the patient at ease and may improve the overall experience and acceptance of the device.<sup>22</sup>

1. Identify the two areas the anaesthesia is to be delivered:
  - Superficial interfascial plane between the serratus anterior and latissimus dorsi muscle.
  - Deeper block between serratus anterior and ribs plane.
2. Inject saline solution to open the fascia, followed by injection of anaesthetic solution.

Figure 1: Deep serratus anterior plane block<sup>22</sup>



L.A. = local anaesthetic  
L.D. = latissimus dorsi  
S.a.m. = serratus anterior muscle



Figures 2: Shows serratus anterior plane block demonstrating the injection site on the chest wall and the corresponding ultrasound images of the injection above and below the serratus anterior muscle.<sup>17</sup>

## 7 | Anaesthesia Options when implanting S-ICD

### 7.4 Truncal Plane Blocks (continued)

#### NOTE



Ensure that the anaesthesia is administered between the muscle layers and not directly into the muscle.

"Discuss the option of plane blocks with your anaesthesia provider and determine if this might be an appropriate option for patients who are not able to tolerate general anaesthesia or MAC, or in obese patients who may require more local anaesthetic."<sup>23</sup>

#### TIP



Ensure that sufficient volume (appropriate dose according to the patient weight) of anaesthesia is used to fully permeate the layers.

Implantation of S-ICD using SAPB may also result in less post-operative pain, and reduce post-operative opioid use.<sup>17</sup>

#### VIDEOS




**For further information or to see related procedural videos please visit our EDUCARE website by clicking [here](#).**

# 14

---

## References



# 14 | References

## 14.1 Bibliography

1. 359481-001 EN Europe 2015-11, EMBLEM™ MRI S-ICD Pulse Generator.
2. 360168-001 EN Europe 2016-09, 3501 EMBLEM™ S-ICD Electrode.
3. PN 104749-001 EN 2013/05, 3200 EMBLEM™ S-ICD Programmer.
4. 360210-001 EN Europe 2017-08, Users Manual EMBLEM™ S-ICD Electrode Delivery System.
5. Knops, R.E. *et al.* The learning curve associated with the introduction of the subcutaneous implantable defibrillator. *Europace* (2016) 18, 1010-1015.
6. Knops, R.E. & Brouwer, T.F. Implantation of the Subcutaneous Implantable Cardioverter-Defibrillator, Chapter 28, Ellenborgers *et al.*, ISBN: 978-0-323-37804-8; PII: B978-0-323-37804-8.00028-6.
7. Brouwer, T.F. *et al.* Implantation of the subcutaneous implantable cardioverter-defibrillator: an evaluation of 4 implantation techniques. *Circ Arrhythm Electrophysiol.* 2017; 10.1161/CIRCEP.116.004663.
8. Shukla, H.H. *et al.* High defibrillation thresholds in transvenous biphasic implantable defibrillators: clinical predictors and prognostic implications. *Pace*, 2003; 26.
9. Okamura, H. *et al.*, Right Parasternal Lead Placement Increases Eligibility for Subcutaneous Implantable Cardioverter Defibrillator Therapy in Adults With Congenital Heart Disease. *Circ J* 2016; 80: 1328 – 1335.
10. S-ICD White Paper The Three Ps: Planning, Preparation and Process (7 step) Seven Step Process for S-ICD Implantation, CRM-385307-AA APR2016.
11. Knops, R.E. *et al.* Two-incision technique for implantation of the subcutaneous implantable cardioverter-defibrillator. *Heart Rhythm*, 2013; 10: 1240 – 1243.
12. Migliore, F. *et al.* Intermuscular two-incision technique for subcutaneous implantable cardioverter defibrillator implantation: results from a multicenter registry. *PACE*, 2017, 00:1 – 8.
13. Amin, A.K. *et al.* Factors Associated With High-Voltage Impedance and Subcutaneous Implantable Defibrillator Ventricular Fibrillation Conversion Success. *Circ Arrhythm Electrophysiol.* 2019;12:e006665.
14. Essandoh, MK, Mark, GE, Aasbo, JD, *et al.*, Anesthesia for subcutaneous implantable cardioverter-defibrillator implantation: Perspectives from the clinical experience of a U.S. panel of physicians. *Pacing Clin Electrophysiol*, 2018. 41: p. 807-816.
15. Mark, G, *Implanter Experience and Considerations.* EP Lab Digest, 2015. 15(7 <http://www.eplabdigest.com/articles/Implanter-Experience-and-Considerations-S-ICD%E2%84%A2-System-Interview-George-Mark-MD>).
16. Miller, MA, Bhatt, HV, Weiner, M, *et al.*, Implantation of the subcutaneous implantable cardioverter-defibrillator with truncal plane blocks. *Heart Rhythm*, 2018. 15(7): p. 1108-1111.
17. Miller, MA, Garg, J, Salter, B, *et al.*, Feasibility of subcutaneous implantable cardioverter-defibrillator implantation with opioid sparing truncal plane blocks and deep sedation. *J Cardiovasc Electrophysiol.* 2019;30:141-148.
18. Kaya, E, *et al.* Conscious sedation during subcutaneous implantable cardioverter-defibrillator implantation using the intermuscular technique. *J Interv Card Electrophysiol*, 2018. Sep 14. doi: 10.1007/s10840-018-0445-9. [Epub ahead of print].
19. Ueshima, H, Hara, E, and Otake, H, Successful cases of S-ICD implantation performed under the serratus plane block. *J Clin Anesth*, 2016. 33: p. 147-8.
20. Blanco, R. *et al.* Serratus plane block: a novel ultrasound-guided thoracic wall nerve block. *Anaesthesia* 2013, 68, 1107-1113.
21. Droghetti, A, *et al.* Ultrasound-guided serratus anterior plane block combined with the two-incision technique for subcutaneous ICD implantation. *Pacing Clin Electrophysiol*, 2018. May;41(5):517-523.

# 14 | References

## 14.1 Bibliography (continued)

22. Droghetti, A. Ultrasound-guided serratus anterior plane block and parasternal block in cooperative sedation for S-ICD implantation. *PACE*, 2019, DOI: 10.1111/pace.13727.
23. Essandoh, M, K. Anesthesia and Pain Control Considerations for the EMBLEM™ S-ICD Device Implant: White Paper, 2019.
24. Droghetti, A. *et al.* Totally submuscular implantation of subcutaneous implantable cardioverter defibrillator: a safe and effective solution for obese or oversized patients. *Clinical Case Reports*, 2016.
25. Ferrari, P. *et al.* Intermuscular pocket for subcutaneous implantable cardioverter defibrillator: Single-center experience. *J. Arrhythmia* (2016).
26. Subzposh, F. *et al.* Successful Implant of the Subcutaneous Implantable Cardioverter Defibrillator in Patients with Low Body Mass Index. *JACC*. 2014; 63:12.
27. Winter, J. *et al.*, Intermuscular technique for implantation of the subcutaneous implantable cardioverter defibrillator: long-term performance and complications. *Europace*, 2016, doi:10.1093/europace/euw297.
28. Emblem S-ICD System Implant Two Incision Intermuscular technique, Dr Andrea Droghetti: <https://www.youtube.com/watch?v=14FTn2nN5zg>.
29. Di Giralamo, E. *et al.* Is shock impedance value alone to be considered a good predictor for shock efficacy in subcutaneous implantable cardioverter defibrillator? *Clinical Case Reports* 2018; 6(4): 574-577.
30. Bettin, M. *et al.* Right versus left parasternal electrode position in the entirely subcutaneous ICD. *Clinical Research in Cardiology*, 2017; <https://doi.org/10.1007/s00392-017-1194-y>.
31. Sasaki *et al.* Subcutaneous Implantable Cardioverter Defibrillator Lead Repositioning for Preventing Inappropriate Shocks Due to Myopotential Oversensing in a Post-Fulminant Myocarditis Patient. *Int Heart J Advance Publication*, 2019; doi: 10.1536/ihj.18-324.
32. Hohnloser, S.H. *et al.* Effect of Amiodarone and Sotalol on Ventricular Defibrillation Threshold The Optimal Pharmacological Therapy in Cardioverter Defibrillator Patients (OPTIC) Trial. *Circulation*. 2006;114:104-109.
33. Verma, A., *et al.*, Incidence of very high defibrillation thresholds (DFT) and efficacy of subcutaneous (SQ) array insertion during implantable cardioverter defibrillator (ICD) implantation. *Journal of Interventional Cardiac Electrophysiology*, 2010. 29(2): p. 127-133.
34. Heist, K. *et al.* Determinants of Subcutaneous Implantable Cardioverter-defibrillator Efficacy: A Computer Modeling Study. *JACC Clinical Electrophysiology*, 2017.
35. Nayak, *et al* 10 Joule Test Shock Impedance Measurement May Be Useful in Predicting Successful or Unsuccessful Defibrillation in Patients Undergoing Subcutaneous Defibrillator Implantation, presented at HRS 2018, Boston.
36. Weiss, R. *et al.* Safety and Efficacy of a Totally Subcutaneous Implantable-Cardioverter Defibrillator. *Circulation*. 2013;128:944-953.
37. Saini, A. & Ellenbogen, K.A. Oversensing and inappropriate shock in a patient with subcutaneous implantable cardioverter defibrillator: What is the mechanism? *Pacing Clin Electrophysiol*. 2017;40:897-899.
38. Zipse, M.M. *et al.* Inappropriate Shocks due to Subcutaneous Air in a Patient With a Subcutaneous Cardiac Defibrillator. *Circ Arrhythm Electrophysiol*. 2014;7:768-770.
39. Kutyla, *et al.* The Need for Pacing in patients who qualify for and ICD: Clinical Implications. *ESC abstract* 2014.
40. Boriani, G. *et al.* Battery drain in daily practice and medium-term projections on longevity of cardioverter-defibrillators: an analysis from a remote monitoring database. *Europace*, <http://dx.doi.org/10.1093/europace/euv436>.
41. Brouwer, T.F. *et al.* Surgical Management of Implantation-Related Complications of the Subcutaneous Implantable Cardioverter-Defibrillator. *J Am Coll Cardiol EP* 2016;2:89-96.

## 14 | References

### 14.1 Bibliography (continued)

42. Yap, S-C. *et al.* Air entrapment causing early inappropriate shocks in a patient with a subcutaneous cardioverter-defibrillator. *Heart Rhythm Case Reports*, 2015;1:156-158).
43. Chinitz, J.S. *et al.* Inappropriate Shocks within 24 hours after Implantation of a Subcutaneous Defibrillator with a Two-incision Technique. *The Journal of Innovations in Cardiac Rhythm Management*. 2016; 7, 2295-2298.



## Emblem™ MRI S-ICD System

### INDICATIONS FOR USE

The S-ICD System is intended to provide defibrillation therapy for the treatment of life-threatening ventricular tachyarrhythmias in patients who do not have symptomatic bradycardia, incessant ventricular tachycardia, or spontaneous, frequently recurring ventricular tachycardia that is reliably terminated with anti-tachycardia pacing.

### CONTRAINDICATIONS

Unipolar stimulation and impedance-based features are contraindicated for use with the S-ICD System.

### WARNINGS

Concomitant use of the S-ICD System and implanted electromechanical devices (for example implantable neuromodulation/neurostimulation systems, ventricular assist device (VAD), or implantable insulin pump or drug pump) can result in interactions that could compromise the function of the S-ICD, the co-implanted device, or both. The S-ICD is intended as lifesaving therapy and should be seen as priority in the decision and evaluation of concomitant system implants over non-lifesaving applications. Electromagnetic (EMI) or therapy delivery from the co-implanted device can interfere with S-ICD sensing and/or rate assessment, resulting in inappropriate therapy or failure to deliver therapy when needed. In addition, a shock from the S-ICD pulse generator could damage the co-implanted device and/or compromise its functionality. Verify sensing configuration, operation modes, surgical considerations and existing placement of all involved devices prior to any co-implant. To help prevent undesirable interactions, test the S-ICD system when used in combination with the co-implanted device, and consider the potential effect of a shock on the co-implanted device. Induction testing is recommended to ensure appropriate detection and time to therapy for the S-ICD and appropriate post-shock operation of the co-implanted device. Failure to ensure appropriate detection and time to therapy delivery of the S-ICD system could result in patient injury or death. Following completion of the interaction testing, thorough follow-up evaluation of all co-implanted devices should be performed to ensure that device functions have not been compromised. If operational settings of the co-implanted devices change or if patient conditions changes which may affect S-ICD sensing and therapy performance, re-evaluation of the co-implanted devices may be required. Do not expose a patient with an implanted S-ICD System to diathermy. EMBLEM S-ICD devices are considered MR Conditional. Unless all MRI Conditions of Use are met, MRI scanning of the patient does not meet MR Conditional requirements for the implanted system. The Programmer is MR Unsafe and must remain outside the MRI site Zone III. During MRI Protection Mode the Tachycardia therapy is suspended. MRI scanning after ERI status has been reached may lead to premature battery depletion, a shortened device replacement window, or sudden loss of therapy. The Beeper may no longer be usable following an MRI scan.

Do not expose a patient with an implanted S-ICD System to diathermy. EMBLEM S-ICD devices are considered MR Conditional. Unless all MRI Conditions of Use are met, MRI scanning of the patient does not meet MR Conditional requirements for the implanted system. The Programmer is MR Unsafe and must remain outside the MRI site Zone III. During MRI Protection Mode the Tachycardia therapy is suspended. MRI scanning after ERI status has been reached may lead to premature battery depletion, a shortened device replacement window, or sudden loss of therapy. The Beeper may no longer be usable following an MRI scan. The pulse generator may be more susceptible to low frequency electromagnetic interference at induced signals greater than 80  $\mu$ V. The S-ICD System has not been evaluated for pediatric use.

### PRECAUTIONS

For specific information on precautions, refer to the following sections of the product labeling: clinical considerations, sterilization and storage, implantation, device programming, environmental and medical therapy hazards, hospital and medical environments, home and occupational environments, follow up testing, explant and disposal, supplemental precautionary information.

### POTENTIAL ADVERSE EVENTS

Potential adverse events related to implantation of the S-ICD System may include, but are not limited to, the following: Acceleration/induction of atrial or ventricular arrhythmia, adverse reaction to induction testing, allergic/adverse reaction to system or medication, bleeding, conductor fracture, cyst formation, death, delayed therapy delivery, discomfort or prolonged healing of incision, electrode deformation and/or breakage, electrode insulation failure, erosion/extrusion, failure to deliver therapy, fever, hematoma/seroma, hemothorax, improper electrode connection to the device, inability to communicate with the device, inability to defibrillate or pace, inappropriate post shock pacing, inappropriate shock delivery, infection, injury to or pain in upper extremity, including clavicle, shoulder and arm, keloid formation, migration or dislodgement, muscle/nerve stimulation, nerve damage, pneumothorax, post-shock/post-pace discomfort, premature battery depletion, random component failures, stroke, subcutaneous emphysema, surgical revision or replacement of the system, syncope, tissue redness, irritation, numbness or necrosis.

Patients who receive an S-ICD System may develop psychological disorders that include, but are not limited to, the following: depression/anxiety, fear of device malfunction, fear of shocks, phantom shocks.

*Refer to the product labeling for specific indications, contraindications, warnings/precautions and adverse events. Rx only. (Rev. E) 046774 AI*

All cited trademarks are the property of their respective owners. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labelling supplied with each device. Information for the use only in countries with applicable health authority product registrations. Materials not intended for use in France.

CRM-627910-AA.

**Boston  
Scientific**  
Advancing science for life™

[www.bostonscientific.eu](http://www.bostonscientific.eu)

© 2021 Boston Scientific Corporation  
or its affiliates. All rights reserved.