

Improving Infection Prevention in Endoscopy and the Role of Single-Use Technologies

Endoscopy in the Era of COVID-19

Current Challenges in Endoscope Hygiene and Processing

Improving Infection Prevention in Endoscopy and the Role of Single-Use Technologies

Highlights of Expert Panel Meeting on Improving Infection Prevention in Endoscopy and the Role of Single-Use Technologies

The Use of Single-Use Duodenoscopes in Endoscopy – a Current Perspective



Advantages, Environmental Impact and Sustainability of Single-Use Endoscopes

Future Outlook: Emerging Trends and Technologies in Endoscopy



Protecting the Environment


The science is clear. Climate change is a threat to the planet and human health. **Boston Scientific became one of the first medical device manufacturers to pledge to achieve carbon neutrality by 2030** in all manufacturing and key distribution sites.



What's Next?

Boston Scientific joined the **United Nations Race to Zero and Business Ambition for 1.5°C campaign**, and is committed to set a science-based target through the Science Based Targets initiative that will put us on a path to **achieve net-zero emissions across our entire value chain by 2050**.

1. Scope 1 & 2. Inclusive of all manufacturing and key distribution sites only. Pages 39-43, 2021 Boston Scientific Performance Report. https://www.bostonscientific.com/content/dam/bostonscientific/corporate/corporate-responsibility/performance-report/2021_boston_scientific_performance_report.pdf. 2. Solid waste diverted from landfills. Inclusive of all manufacturing and key distribution sites only. Page 44, 2021 Boston Scientific Performance Report. https://www.bostonscientific.com/content/dam/bostonscientific/corporate/corporate-responsibility/performance-report/2021_boston_scientific_performance_report.pdf.



Published by Global Business Media
Global Business Media Limited
62 The Street Ashted Surrey KT21 1AT
United Kingdom

Switchboard: +44 (0)1737 850 939
Fax: +44 (0)1737 851 952
Email: info@globalbusinessmedia.org
Website: www.globalbusinessmedia.org

Publisher
Kevin Bell
Business Development Director
Marie-Anne Brooks

Editor
Dr. Michael A. James PhD
Senior Project Manager
Steve Banks
Advertising Executives
Michael McCarthy
Abigail Coombes
Production Manager
Paul Davies

For further information visit:
www.globalbusinessmedia.org

The opinions and views expressed in the editorial content in this publication are those of the authors alone and do not necessarily represent the views of any organisation with which they may be associated.

Material in advertisements and promotional features may be considered to represent the views of the advertisers and promoters. The views and opinions expressed in this publication do not necessarily express the views of the Publishers or the Editor. While every care has been taken in the preparation of this publication, neither the Publishers nor the Editor are responsible for such opinions and views or for any inaccuracies in the articles.

© 2022. The entire contents of this publication are protected by copyright. Full details are available from the Publishers. All rights reserved. No part of this publication may be reproduced, stored in a retrieval system or transmitted in any form or by any means, electronic, mechanical photocopying, recording or otherwise, without the prior permission of the copyright owner.

Contents

Foreword <i>Michael A James, PhD</i>	2
Endoscopy in the Era of COVID-19 <i>Michael A James, PhD</i>	3
Precautions for the Prevention of SARS-CoV-2 Causing Longer Wait Times Enhanced Decontamination Features The Need for New Technologies for Efficient and Safe Endoscopy	
Current Challenges in Endoscope Hygiene and Processing <i>Michael A James, PhD</i>	6
Risk of Infection During Endoscopy Insufficiency of Decontamination of Endoscopes Reprocessing Delays Conclusions	
Improving Infection Prevention in Endoscopy and the Role of Single Use Technologies <i>Enrique Vazquez-Sequeiros MD, PhD</i>	9
Potential Solutions to Reduce the Risk of Duodenoscope Transmission of Infections Questions Regarding Single-Use Duodenoscopes	
Improving Infection Prevention in Endoscopy and the Role of Single-Use Technologies – Expert Panel Meeting <i>Panel members listed at beginning of article</i>	14
Introduction Patient and Device Factors Related to Duodenoscope-Transmitted Infections Considerations for Patient Selection: Center Experience with Single-Use Duodenoscopes Patient Selection: Perspectives from Middle Eastern Centers	
The Use of Single-Use Duodenoscopes in Endoscopy – a Current Perspective <i>Prof. Dr. Mark Ellrichmann, Claudio C. Conrad</i>	21
Introduction Clinical Functionality of Single-Use Duodenoscopes Cost Assessment of Single-Use Endoscopes Patient Selection for Single-Use Duodenoscopes Conclusion	
Advantages, Environmental Impact and Sustainability of Single-Use Endoscopes <i>Michael A James, PhD</i>	25
Advantages Performance Comparison Environmental Impact Economic Sustainability Conclusion	
Future Outlook: Emerging Trends and Technologies in Endoscopy <i>Michael A James, PhD</i>	28
New Procedures and Training Advancements in Wireless Capsules Imaging Technologies ERCP/Duodenoscopy Advancements Conclusion	

Foreword

The utility of endoscopy and clinical outcomes related to endoscopic procedures continues to evolve with new circumstances and new technologies. Improved imaging, artificial intelligence and safety innovations have taken endoscopy into new interventional realms and improved diagnostic yield. Safety in regard to healthcare-acquired infection has recently been an area of concern, particularly that involving multidrug-resistant organism transmission during endoscopy.

Safety of endoscopic procedures for both patients and healthcare workers has come to the forefront in recent years because of the COVID-19 pandemic, accompanied by increased cost, stress on staff and requirements for sanitation and scope reprocessing. New technologies, including sterilizable duodenoscopes, disposable components and single-use endoscopes have been introduced as solutions to these problems. In this Report, we discuss (i) the implications of endoscopy in the COVID-19 era, including workflow, wait times and reprocessing requirements; (ii) emerging issues, insufficiencies and

challenges in successful reprocessing and high-level disinfection of reusable endoscopes in the context of recent outbreaks of multidrug-resistant organisms; (iii) the utility of single-use endoscopes to mitigate the risk of infection from endoscopic procedures, such as endoscopic retrograde cholangiopancreatography (ERCP) and endoscopy of the pancreatobiliary tree; (iv) performance, acceptance and patient selection for single-use endoscopy; (v) economic sustainability and environmental impact of single-use scopes, and (vi) emerging trends and advancements in endoscopic procedures, imaging, computer-aided diagnosis and prevention of infection. Continued innovation and technology is broadening the utility of endoscopic procedures. These innovations are bringing new challenges in assessment, implementation and training related to new technologies and procedures but are also improving our ability to ensure positive outcomes and safety for patients.

**Michael James
Editor**

Dr. Michael A. James PhD, is a medical writer, biotech entrepreneur/founder in the fields of oncology and virology, and former faculty of Surgery and Pharmacology/Toxicology at the Medical College of Wisconsin. He holds a PhD in microbiology from the University of Iowa and was trained in cancer cell biology and molecular biology at Washington University in St. Louis.

Endoscopy in the Era of COVID-19

Michael A James, PhD

THE COVID-19 pandemic has had significant impact on the safety and efficiency of endoscopy procedures. The causal virus, SARS-CoV-2, is known to be transmitted through respiratory droplets, surfaces and aerosols^[1]. In addition to consideration of potential exposure to body fluids and contamination of surfaces during endoscopy, it is appropriate to consider endoscopy an aerosol-generating procedure. With over 50% of COVID-19 patients exhibiting gastrointestinal (GI) symptoms^[2] and detection of virus in feces (although not correlative to GI symptoms^[3]), oral-fecal transmission should also be considered.

SARS-CoV-2 has been found to be stable for 3 hours as an aerosol and 72 hours on certain surfaces, including plastic and stainless steel^[4]. An early study of SARS-CoV-2 infection in Wuhan showed that healthcare workers comprised 29% of patients^{[5],[6]}. Accordingly, the World Health Organization offered interim guidance on the use of personal protective equipment (PPE) during endoscopy^[7], including the use of hairnets, face shields, goggles, N-95 or FFP-2/3 masks, full sleeve gowns, gloves and shoe covers. In high-risk cases of fever, COVID-19 symptoms, positive contact or travel to endemic areas, endoscopy was recommended to be limited to designated rooms or isolation rooms if available. Decontamination procedures included surface disinfection of procedure rooms and standard high-level disinfection of endoscopy equipment. Despite fecal detection of virus, lower endoscopy is viewed as lower risk by US and UK guiding organizations, leading to less strict PPE guidelines for these procedure^[8]. GI endoscopy in general is considered to be relatively safe when appropriate PPE and decontamination are observed^[9]. However, compliance with PPE use, decontamination procedure and safety precautions cannot be assumed. While PPE use and decontamination are important for the prevention of SARS-CoV-2 transmission to healthcare workers, documented cases of transmission to patients through endoscopy are rare. However, with increasing transmissibility of modern variants,



Endoscopy room and PPE. "File:Mani Zadeh MD Endoscopic Sinus Surgery.jpg" by BestInPlastics is licensed under CC BY-SA 3.0.

this may be changing. There is also special risk conferred by asymptomatic infection, which is on the rise^[10].

Precautions for the Prevention of SARS-CoV-2 Causing Longer Wait Times

Perhaps a more significant impact has been on efficiency of endoscopic work flow and wait times. Endoscopy procedures are taking longer because of precautions for the prevention of SARS-CoV-2 transmission, causing concern among physicians that cancer patients and those with other GI illnesses are not getting diagnoses and treatment in a timely manner. GI endoscopy wait times may have already been a problem before the pandemic as illustrated by a Canadian study in 2008 showing that 78.6% of patients were not meeting standards for wait times^[11]. In a simulation model, workflow changes brought on by the COVID-19 pandemic significantly and negatively affected operation and productivity measures for endoscopy, thereby

In addition to consideration of potential exposure to body fluids and contamination of surfaces during endoscopy, it is appropriate to consider endoscopy an aerosol-generating procedure.

With over 50% of COVID-19 patients exhibiting gastrointestinal (GI) symptoms and detection of virus in feces (although not correlative to GI symptoms), oral-fecal transmission should also be considered

A factor that has significantly affected turnover time for endoscopy procedures has been enhanced decontamination procedures that have been implemented to combat COVID-19 infection

affecting financial metrics^[12]. The British Society of Gastroenterology published guidance in 2020 recommending that non-essential endoscopic procedures be reduced and that FIT screening and Bowel Scope be paused, limiting procedures to emergency and essential procedures^[13]. Emergency and essential procedures were limited to acute bleeding, foreign bodies, obstructions, hepato-pancreatobiliary, nutrition, vacuum therapy and ongoing GI bleeds. Globally, similar guidance has left surveillance and cancer staging procedures with extended wait times and decreased availability. In a global survey of GI endoscopy units, 85% reported decreasing procedure volumes by over 50% because of the pandemic with 2.45% completely halting procedures^[14].

Components of COVID-19-related measures that have affected workflow and wait times for endoscopy have included the need for dedicated or isolated rooms (depending on COVID-19-status), COVID-19 screening and assessment before procedures, staffing issues and enhanced sanitation procedures^[12]. Site structure of endoscopy units has significantly changed as a result of COVID-19. These changes have included modification of unit layouts, incorporating risk based dedication of spaces in waiting, recovery and procedure rooms, implementation of new checkpoint and PPE donning areas and the creation of separate pathways and processes according to risk assignments^[15]. Staffing issues have included training of staff in COVID-19 precautions, reassignment of at-risk staff, illness with COVID-19 among staff, increasing staff furloughs and greater staff needs for screening and assessment^[16].

Enhanced Decontamination Features

A factor that has significantly affected turnover time for endoscopy procedures has been

enhanced decontamination procedures that have been implemented to combat COVID-19 infection^[17]. Surface disinfection of rooms and surfaces is a large part of enhanced procedures that contributes to increased turnover time. Another is disinfection of endoscopy equipment. Endoscope disinfection procedures are to include high-level disinfection (HLD), which is typically performed on an automatic endoscope reprocessor. Such reprocessors employ cycles of exposure to a chemical disinfectant solution at specific temperatures and for a recommended duration^[18]. Although guidelines have not advised significant changes to endoscope disinfection protocols, strict adherence to existing recommendations is emphasized^[19]. Proper disinfection of endoscopy equipment is highly dependent on operator training. Standard procedures for endoscope reprocessing include pre-cleaning at the bedside, testing for leaks, manual cleaning and HLD, which have been shown to be effective at eliminating SARS-CoV-2 contamination^[19]. Novel approaches to endoscope sterilization have included techniques such as low-temperature plasma-activated gas, which produces ultraviolet light and free radicals to kill a broad range of microorganisms^[19]. Additional advances have been made in the use of infection barriers and the disposable endoscopes or endoscope components.

The Need for New Technologies for Efficient and Safe Endoscopy

Effects of the COVID-19 pandemic on endoscopy in terms of both safety and efficient work flow are having a significant impact on patients, health care workers and endoscopy units. To best ensure safe procedures and provide effective and timely services to patients, new perspectives on operations in these units and new technologies for efficient and safe endoscopy will need to be implemented.

In a simulation model, workflow changes brought on by the COVID-19 pandemic significantly and negatively affected operation and productivity measures for endoscopy, thereby affecting financial metrics

References:

- [1] C. del Rio and P. N. Malani, "COVID-19—New Insights on a Rapidly Changing Epidemic," JAMA, vol. 323, no. 14, pp. 1339–1340, Apr. 2020, doi: 10.1001/jama.2020.3072.
- [2] L. Pan et al., "Clinical Characteristics of COVID-19 Patients With Digestive Symptoms in Hubei, China: A Descriptive, Cross-Sectional, Multicenter Study," Am J Gastroenterol, vol. 115, no. 5, pp. 766–773, May 2020, doi: 10.14309/ajg.0000000000000620.
- [3] J. Ong, B. E. Young, and S. Ong, "COVID-19 in gastroenterology: a clinical perspective," Gut, vol. 69, no. 6, pp. 1144–1145, Jun. 2020, doi: 10.1136/gutjnl-2020-321051.
- [4] N. van Doremalen et al., "Aerosol and Surface Stability of SARS-CoV-2 as Compared with SARS-CoV-1," N Engl J Med, vol. 382, no. 16, pp. 1564–1567, Apr. 2020, doi: 10.1056/NEJMc2004973.
- [5] R. Soetikno et al., "Considerations in performing endoscopy during the COVID-19 pandemic," Gastrointestinal Endoscopy, vol. 92, no. 1, pp. 176–183, Jul. 2020, doi: 10.1016/j.gie.2020.03.3758.
- [6] D. Y. Wang and D. B. Johnson, "When to consider alternatives to front-line immune therapies in metastatic melanoma," Melanoma Management, vol. 4, no. 2, pp. 71–74, May 2017, doi: 10.2217/mmt-2017-0005.
- [7] W. H. Organization, "Rational use of personal protective equipment for coronavirus disease (COVID-19): interim guidance, 27 February 2020," Art. no. WHO/2019-nCoV/PCPPE_use/2020.1, 2020, Accessed: May 05, 2022. [Online]. Available: <https://extranet.who.int/iris/restricted/handle/10665/331215>
- [8] J. Ong, G. B. Cross, and Y. Y. Dan, "Prevention of nosocomial SARS-CoV-2 transmission in endoscopy: international recommendations and the need for a gold standard," Gut, vol. 69, no. 6, pp. 1145–1148, Jun. 2020, doi: 10.1136/gutjnl-2020-321154.
- [9] A. Repici et al., "Low risk of COVID-19 transmission in GI endoscopy," Gut, vol. 69, no. 11, pp. 1925–1927, Nov. 2020, doi: 10.1136/gutjnl-2020-321341.
- [10] R. Subramanian, Q. He, and M. Pascual, "Quantifying asymptomatic infection and transmission of COVID-19 in New York City using observed cases, serology, and testing capacity," Proceedings of the National Academy of Sciences, vol. 118, no. 9, p. e2019716118, Mar. 2021, doi: 10.1073/pnas.2019716118.
- [11] D. Yu, W. M. Hopman, and W. G. Paterson, "Wait time for endoscopic evaluation at a Canadian tertiary care centre: Comparison with Canadian Association of Gastroenterology targets," Can J Gastroenterol, vol. 22, no. 7, pp. 621–626, Jul. 2008.
- [12] A. Das, "Impact of the COVID-19 pandemic on the workflow of an ambulatory endoscopy center: an assessment by discrete event simulation," Gastrointestinal Endoscopy, vol. 92, no. 4, pp. 914–924, Oct. 2020, doi: 10.1016/j.gie.2020.06.008.
- [13] "Endoscopy activity and COVID-19: BSG and JAG guidance," The British Society of Gastroenterology, Apr. 03, 2020. <https://www.bsg.org.uk/covid-19-advice/endoscopy-activity-and-covid-19-bsg-and-jag-guidance/> (accessed May 05, 2022).
- [14] M. Alborae et al., "The global impact of COVID-19 on gastrointestinal endoscopy units: An international survey of endoscopists," Arab Journal of Gastroenterology, vol. 21, no. 3, pp. 156–161, Sep. 2020, doi: 10.1016/j.ajg.2020.08.008.
- [15] V. Cennamo et al., "Redesign of a GI endoscopy unit during the COVID-19 emergency: A practical model," Digestive and Liver Disease, vol. 52, no. 10, pp. 1178–1187, Oct. 2020, doi: 10.1016/j.dld.2020.05.007.
- [16] A. Perisetti, H. Goyal, and N. Sharma, "Gastrointestinal Endoscopy in the Era of COVID-19," Frontiers in Medicine, vol. 7, 2020, Accessed: May 06, 2022. [Online]. Available: <https://www.frontiersin.org/article/10.3389/fmed.2020.587602>
- [17] A. Abdelqader et al., "Impact of the SARS-CoV-2 pandemic on turnover time and revenue in the endoscopy unit: single-center experience," Endosc Int Open, vol. 9, no. 11, pp. E1680–E1685, Nov. 2021, doi: 10.1055/a-1546-8302.
- [18] K.-W. Chiu, L.-S. Lu, and S.-S. Chiou, "High-level disinfection of gastrointestinal endoscope reprocessing," World J Exp Med, vol. 5, no. 1, pp. 33–39, Feb. 2015, doi: 10.5493/wjem.v5.i1.33.
- [19] T. Chua, N. Halim, and S. Reicher, "Recent Advances in Endoscope Disinfection: Where Do We Stand in the COVID era?," Tech Innov Gastrointest Endosc, vol. 23, no. 2, pp. 190–198, 2021, doi: 10.1016/j.tige.2020.10.001.

Perhaps a more significant impact has been on efficiency of endoscopic work flow and wait times. Endoscopy procedures are taking longer because of precautions for the prevention of SARS-CoV-2 transmission, causing concern among physicians that cancer patients and those with other GI illnesses are not getting diagnoses and treatment in a timely manner

Current Challenges in Endoscope Hygiene and Processing

Michael A James, PhD




Figure 1. E. coli on EMB agar. Credit: Adonofrio (Biology101.org). Licensed under CC BY 2.0.

Effective reprocessing of endoscopy equipment is challenging and may not be regularly accomplished, thus increasing the risk for cross-contamination and infection

Risk of Infection During Endoscopy

Infection related to endoscopic procedures is a long-standing issue, with nosocomial infections being a common complication following endoscopic retrograde cholangiopancreatography (ERCP)^[1]. Infection in endoscopy patients has two sources, endogenous infections which come from the patient's own microflora^[2] and exogenous, which originate from environmental sources, typically a contaminated endoscope^[2]. The most commonly isolated pathogens from contaminated duodenoscopes are *Escherichia coli*, *Staphylococcus aureus*, *Klebsiella pneumoniae*, *Enterococcus faecalis* and *Pseudomonas aeruginosa*^{[2],[3]}. The 2020 British Society of Gastroenterology guidelines for endoscope decontamination have identified four microorganisms of concern during the past 20 years^[4]. These include mycobacteria, which can be spread during bronchoscopy, particularly to immunocompromised or HIV-positive patients; bacterial spores from *Bacillus* or *Clostridium* species, which have been isolated from endoscopes but not linked to infections; multidrug-resistant gram-negative bacilli; and prions, which are highly resistant to decontamination procedures. Additionally, recent outbreaks of multidrug-resistant organism infections have been associated with duodenoscopy^[2], including those occurring after ERCP^[5]. Many of these events occurred despite proper

reprocessing of endoscopy equipment^[6]. A 2019 review of duodenoscope contamination rates published between 1998 and 2018 showed that contamination rates range between 0.69% and 60%^[6]. Because of concerns with duodenoscope contamination data and ongoing challenges with endoscope reprocessing efficacy the FDA has

recommended a transition from re-usable duodenoscopes with fixed endcaps to fully disposable duodenoscopes or reusable duodenoscopes with disposable components^[7].

Insufficiency of Decontamination of Endoscopes

Effective reprocessing of endoscopy equipment is challenging and may not be regularly accomplished, thus increasing the risk for cross-contamination and infection. A major factor leading to ineffective reprocessing and increased risk of infection is non-adherence to reprocessing guidelines coupled with the complex design of endoscopes^{[7],[8]}. In a prospective study of endoscope reprocessing effectiveness, 1.4% of manual cleaning procedures and 75.4% of automatic HLD procedures were found to adhere to documented guidelines^[9]. There is frequent residual microbial presence on reusable endoscopes following HLD, and disinfectants commonly fail tests of minimum effective concentration (MEC) prior to the expiration of their maximum use periods^[10]. The US Food and Drug Administration (FDA) monitored adherence to endoscope reprocessing protocols for multiple major manufacturers of duodenoscopes and found that protocols were not followed correctly 27–94% of the time^{[11],[12]}. In addition, unacceptably high contamination rates of endoscopic equipment even after reprocessing following guidelines were reported in this study.

Complex instructions for the use and testing of HLD and human failures to properly conduct these procedures were suggested to contribute to these issues. Non-adherence to guidelines, use of insoluble products (e.g. simethicone), contaminated water, inadequate drying and damaged endoscopes were cited as specific factors influencing the effectiveness of HLD.

Endoscope decontamination procedures recommended by 2020 British Society of Gastroenterology guidelines start with bedside pre-cleaning to remove visible organic material^[4]. This step prevents blockage of the channels due to dried soil and is followed by a leak test to rule out damage that might allow internal endoscope exposure to fluids. Manual cleaning is next, which includes brushing and flushing of the channels with a compatible detergent to remove organic material and clear blockages. Recommendations for the use of enzymatic detergents during this step have been reduced because of related hypersensitivities^[13]. Finally, the endoscope is rinsed and dried before automatic high-level disinfection in an endoscope washer/disinfector (EWD). A recent systematic review of the literature published from 2008 to 2020 sought to identify failures in endoscope processing that have led to insufficient decontamination^[14]. The most common failure was the absence of or inadequacies in drying. The next most common failure was manual cleaning without brushing all appropriate channels. Other gaps included failure to test for leaks, improper storage, improper use of disinfectant solution, insufficient time of immersion, insufficient monitoring of MEC, failure to pre-clean and improper transportation to processing facilities.

Aside from issues of adherence to reprocessing guidelines, contamination rates have been estimated as high as 60% on duodenoscopes even after proper decontamination^{[6],[15]}. Outbreaks of multidrug-resistant organisms have been attributable to ERCP in cases where endoscope reprocessing has adhered to manufacturer's guidelines^[5]. This residual risk despite strict reprocessing guideline adherence highlights the need for technologies that minimize the risk of infection.

Reprocessing Delays

Delays in endoscope reprocessing occur for a wide variety of reasons and may result in a disruption of patient care. In a survey of healthcare workers, the majority of respondents reported pressure to work faster during duodenoscope cleaning procedures^[16]. Greater than 75% of these healthcare workers experienced discomfort or fatigue when conducting cleaning procedures. When asked to identify practices to best improve reprocessing comfort and efficiency, respondents

identified mentoring and retention of experienced staff as most important. Delays in endoscopy procedures were recently found to negatively impact quality-of-care and increase the length of hospital stays for inpatients by a median of 2 days^[17]. This study also identified inpatient endoscopy delays as an independent risk factor for readmission within 30 days. Workflow limitations are a significant factor in these delays as endoscopy unit availability was one of the most common causes of delays in inpatient endoscopy in the study. Positive cultures or use on high-risk of COVID-19-positive patients requires quarantine of endoscopes, which contributes to downtime and equipment availability. Enhanced surveillance and reprocessing techniques have been developed as a result of recent outbreaks, which further increase scope downtime, requiring 3.4 times the number of scopes to maintain the same volume of procedures^[18]. These recommended procedures include microbial culture, liquid chemical sterilization, repeat HLD or ethylene oxide sterilization in addition to manufacturer-recommended reprocessing. The available clinical data indicate that these supplemental measures to enhance duodenoscope reprocessing do not effectively eliminate duodenoscope contamination^[19].

It has been proposed that a shift to sterilization or the use of sterile, disposable endoscopes is required to ensure the safety of patients^[20]. Disposable, single-use endoscopy equipment represents a technology that may minimize downtime and improve safety for endoscopic procedures. Single-use endoscopes eliminate the need for reprocessing and related costs, staff requirements and training. In an early study comparing a disposable sheath for fiber optic sigmoidoscopy to reprocessing of endoscopes found that overall downtime was reduced from 32.8 minutes to 8.1 minutes^[21]. This advantage is in addition to potentially improved staff and patient safety. Disposable technology may also minimize cost and downtime associated with scope repair. It was determined that colonoscopies require repair over 3 times each year^[22].

Conclusions

It is clear that there are ongoing challenges in effective disinfection of endoscopy equipment and the efficiency of endoscopic procedures in terms of workflow, staffing, training, downtime and cost. As such, improved endoscopy unit design, protocols and technologies are required to ensure quality patient care and the safety and efficient employment of healthcare workers. Implementation of these measures will require analysis of benefits to safety and efficiency, financial considerations and education of physicians and healthcare workers.

Improved endoscopy unit design, protocols and technologies are required to ensure quality patient care and the safety and efficient employment of healthcare workers

References:

- [1] G. G. Balan, C. V. Sfarti, S. A. Chiriac, C. Stanciu, and A. Trifan, "Duodenoscope-associated infections: a review," Eur J Clin Microbiol Infect Dis, vol. 38, no. 12, pp. 2205–2213, Dec. 2019, doi: 10.1007/s10096-019-03671-3.
- [2] M. R. Rahman, A. Perisetti, R. Coman, P. Bansal, R. Chhabra, and H. Goyal, "Duodenoscope-Associated Infections: Update on an Emerging Problem," Dig Dis Sci, vol. 64, no. 6, pp. 1409–1418, Jun. 2019, doi: 10.1007/s10620-018-5431-7.
- [3] "Antimicrobial resistance in the EU/EEA (EARS-Net) - Annual Epidemiological Report for 2019," European Centre for Disease Prevention and Control, Nov. 18, 2020. <https://www.ecdc.europa.eu/en/publications-data/surveillance-antimicrobial-resistance-europe-2019> (accessed May 17, 2022).
- [4] "2020 Guidance on Decontamination of Equipment for Gastrointestinal Endoscopy," The British Society of Gastroenterology. <https://www.bsg.org.uk/clinical-resource/guidance-on-decontamination-of-equipment-for-gastrointestinal-endoscopy/> (accessed May 10, 2022).
- [5] J. T. Higa, "Duodenoscope-related infections: Overview and epidemiology," Techniques in Gastrointestinal Endoscopy, vol. 21, no. 4, p. 150623, Oct. 2019, doi: 10.1016/j.tgie.2019.150623.
- [6] G. Thornhill and M. David, "Endoscope-associated infections: A microbiologist's perspective on current technologies," Techniques in Gastrointestinal Endoscopy, vol. 21, no. 4, p. 150625, Oct. 2019, doi: 10.1016/j.tgie.2019.150625.
- [7] C. for D. and R. Health, "Infections Associated with Reprocessed Duodenoscopes," FDA, Apr. 2022, Accessed: Jun. 23, 2022. [Online]. Available: <https://www.fda.gov/medical-devices/reprocessing-reusable-medical-devices/infections-associated-reprocessed-duodenoscopes>
- [8] T. Chua, N. Halim, and S. Reicher, "Recent Advances in Endoscope Disinfection: Where Do We Stand in the COVID era?," Tech Innov Gastrointest Endosc, vol. 23, no. 2, pp. 190–198, 2021, doi: 10.1016/j.tgie.2020.10.001.
- [9] C. L. Ofstead, H. P. Wetzler, A. K. Snyder, and R. A. Horton, "Endoscope Reprocessing Methods: A Prospective Study on the Impact of Human Factors and Automation," Gastroenterology Nursing, vol. 33, no. 4, pp. 304–311, Aug. 2010, doi: 10.1097/SGA.0b013e3181e9431a.
- [10] C. L. Ofstead, K. M. Hopkins, B. L. Buro, J. E. Eiland, and H. P. Wetzler, "Challenges in achieving effective high-level disinfection in endoscope reprocessing," American Journal of Infection Control, vol. 48, no. 3, pp. 309–315, Mar. 2020, doi: 10.1016/j.ajic.2019.09.013.
- [11] "522 Postmarket Surveillance Studies Database." https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfPMA/pss.cfm?t_id=354&c_id=3692 (accessed Jul. 12, 2022).
- [12] "522 Postmarket Surveillance Studies Database." https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfPMA/pss.cfm?t_id=355&c_id=3693 (accessed Jul. 12, 2022).
- [13] A. Adisesh, E. Murphy, C. M. Barber, and J. G. Ayres, "Occupational asthma and rhinitis due to detergent enzymes in healthcare," Occup Med (Lond), vol. 61, no. 5, pp. 364–369, Aug. 2011, doi: 10.1093/occmed/kqr107.
- [14] R. Madureira and A. Oliveira, "Endoscopic processing: what are the gaps in clinical practice?," Revista Eletrônica de Enfermagem, vol. 23, pp. 1–13, Jun. 2021, doi: 10.5216/ree.v23.66550.
- [15] C. for D. and R. Health, "Safety Communications - Supplemental Measures to Enhance Duodenoscope Reprocessing: FDA Safety Communication." <http://wayback.archive-it.org/7993/20170722150658/https://www.fda.gov/MedicalDevices/Safety/AlertsandNotices/ucm454766.htm> (accessed May 11, 2022).
- [16] A. D. Sivek et al., "Healthcare worker feedback on duodenoscope reprocessing workflow and ergonomics," American Journal of Infection Control, Jan. 2022, doi: 10.1016/j.ajic.2022.01.012.
- [17] C. C. Jacobs et al., "Factors Associated With Inpatient Endoscopy Delay and its Impact on Hospital Length-of-Stay and 30-Day Readmission," Clinical Gastroenterology and Hepatology, vol. 19, no. 12, pp. 2648–2655, Dec. 2021, doi: 10.1016/j.cgh.2021.06.009.
- [18] S. Bomman et al., "Economic burden of enhanced practices of duodenoscopes reprocessing and surveillance: balancing risk and cost containment," Endosc Int Open, vol. 09, no. 9, pp. E1404–E1412, Sep. 2021, doi: 10.1055/a-1515-2591.
- [19] L. W. Day et al., "Multisociety guideline on reprocessing flexible GI endoscopes and accessories," Gastrointestinal Endoscopy, vol. 93, no. 1, pp. 11-33, e6, Jan. 2021, doi: 10.1016/j.gie.2020.09.048.
- [20] W. A. Rutala, H. Kanamori, E. E. Sickbert-Bennett, and D. J. Weber, "What's new in reprocessing endoscopes: Are we going to ensure 'the needs of the patient come first' by shifting from disinfection to sterilization?," American Journal of Infection Control, vol. 47, pp. A62–A66, Jun. 2019, doi: 10.1016/j.ajic.2019.01.017.
- [21] R. I. Rothstein and B. Littenberg, "Disposable, sheathed, flexible sigmoidoscopy: A prospective, multicenter, randomized trial," Gastrointestinal Endoscopy, vol. 41, no. 6, pp. 566–572, Jun. 1995, doi: 10.1016/S0016-5107(95)70192-3.
- [22] S. Larsen, A. Kalloo, and S. Hutfless, "The hidden cost of colonoscopy including cost of reprocessing and infection rate: the implications for disposable colonoscopes," Gut, vol. 69, no. 2, pp. 197–200, Feb. 2020, doi: 10.1136/gutjnl-2019-319108.

Delays in endoscope reprocessing occur for a wide variety of reasons and may result in a disruption of patient care

Improving Infection Prevention in Endoscopy and the Role of Single Use Technologies

This paper describes the potential risk of infection after endoscopy of the pancreatobiliary tree and discusses about alternatives to solve this problem. This review focusses more in the new single-use duodenoscopes as a promising alternative.

Author

Enrique Vazquez-Sequeiros MD, PhD evazquezse@gmail.com

Author Affiliation

1. Head of Endoscopy Unit. Gastroenterology and Hepatology Service. University Hospital Ramon y Cajal, IRYCIS, Madrid. Spain.

Conflict of Interest Disclosure

EVS: Has received money from lectures and/or consultancy by Ella-Biomed, Boston Scientific, Olympus and Ambu.

Key Words (3-5 MeSH)

Single use duodenoscope; disposable; cost; infectious disease; safety.

Corresponding Author

Enrique Vazquez-Sequeiros MD, PhD

Email: evazquezse@gmail.com

Endoscopy Unit. Gastroenterology and Hepatology Service.

University Hospital Ramon y Cajal, IRYCIS, Madrid. Spain.

Carretera de Colmenar Km 9,200. CP 28034. Madrid.

To perform an upper or lower gastrointestinal endoscopy, an endoscope has to be introduced through a natural orifice (mouth or anus). Although a minimally invasive procedure, any endoscopy performance associates a non-neglectable risk of infection. This risk is significantly higher when a side-viewing endoscope, either a duodenoscope or an echoendoscope, more complex in design and therefore difficult to clean, is used^[1]. These types of endoscopes, with the camera and working channel exit located at the side of the distal end of the endoscope and an elevator nail for needle or catheter redirection at the same level, may certainly be more difficult to clean. A suboptimal disinfection of these side-viewing endoscopes will result in a higher risk for bacterial infection. The magnitude of this problem increases in patients with a higher risk for infection, such as oncologic patients and those

receiving chemotherapy or immune suppressive drugs, i.e. transplanted patients. If the transmitted bacteria is a so-called "superbug" or multidrug-resistant (MDR) organism (e.g. carbapenem-resistant Enterobacteriaceae), the risk for a lethal outbreak of infection related to endoscopic retrograde cholangiopancreatography (ERCP) will certainly increase, as has been shown in recent years^[2,3].

Potential Solutions to Reduce the Risk of Duodenoscope Transmission of Infections

Adequate disinfection of endoscopes has become a major problem in recent years, especially for side-viewing endoscopes, and a number of potential solutions have been proposed to solve or at least minimize this problem:

Although a minimally invasive procedure, any endoscopy performance associates a non-neglectable risk of infection

Adequate disinfection of endoscopes has become a major problem in recent years, especially for side-viewing endoscopes

1.) Sterilization: Sterilization of side-viewing duodenoscopes has been proposed by experts in the field as an effective alternative. However, this approach requires a longer time to reprocess the duodenoscope, and the elevated temperatures needed for sterilization may damage the endoscope and reduce the lifecycle of the expensive duodenoscopes^[4]. These shortcomings of sterilization would result in a need to buy a larger number of duodenoscopes to do the same number of ERCPs, and endoscopes would need to be replaced for damage more frequently. For these reasons, sterilization does not seem to be an adequate option and is only reserved for selected cases, such as contaminated duodenoscopes despite repeated high-level disinfection^[4].

2.) Modifications of design of reusable duodenoscopes: Manufacturers of reusable duodenoscopes have also been very active at trying to solve the problem of infections related to the use of their products and difficulties encountered in the cleaning of their non-disposable duodenoscopes. Some companies have modified the protocols and material used (e.g. specifically designed cleaning brushes) for manual cleaning of the distal end of the duodenoscope. Furthermore, all manufacturers of duodenoscopes (Fujifilm Corporation, Olympus Medical Systems and Pentax Medical) have redesigned the distal part of the endoscope by adding disposable/ single-use components such as disposable end-caps or single-use elevator nails in an attempt to minimize the risk for infection^[5-7]. These initiatives directed to the most critical parts of the duodenoscope disinfection process are very valuable but cannot completely prevent an infection.

3.) Fully disposable (single-use) duodenoscopes: A final alternative to avoid iatrogenic infections related to the use of a duodenoscope, and probably the more ambitious and safer one, would be the use of single-use duodenoscopes. This futuristic approach has become a reality in the very recent years and is the fruit of intensive research at leading companies in the field of endoscopy. The Boston Scientific Corporation EXALT Model D (Figure1) is the first single-use duodenoscope that reached the market and has been fully available for commercial use since that moment, while the Ambu Innovation GmbH (aScope) has also introduced a disposable duodenoscope in more recent dates. Endoscopists should be aware that single-use duodenoscopes, despite differences between manufacturers, are essentially similar in design to reusable ones (working channel length: 1240

mm; insertion tube outer diameter: 11.3 mm; angulation range: up: 120°; down: 90°; right: 110°; left: 90°; working channel inner diameter: 4.2 mm; field of view: 108-130°). The concept of single-use duodenoscopes would eliminate the transmission of microorganisms from one patient to another due to ineffective processing^[7,8] and appears to be the more reasonable approach to pursue in coming years if we want to definitively avoid risks for crossover infection related to ERCP performance.

Questions Regarding Single-Use Duodenoscopes

Although we believe that the single-use duodenoscope approach is the way to go in the future, there are still some questions that need to be further clarified:

1.) Performance characteristics of single-use vs reusable duodenoscopes: Before we decide to switch to single-use duodenoscopes, we need to be sure that they are equally effective, both technically and clinically, when compared with current state of the art "reusable duodenoscopes". At present time we can only find published data reporting the efficacy of single-use duodenoscopes manufactured by the Boston Scientific Corporation (EXALT Model D) (Figures 1, 2 and 3).^[7,8] These initial reports have been very promising, demonstrating that ERCP performance with these new endoscopes is feasible and provides reasonable clinical results^[7,8]. This was elegantly shown in a recent multicenter study conducted at 6 academic medical centers^[8]. The authors demonstrated that 96.7% of 60 patients evaluated could be successfully treated with a disposable duodenoscope, and only 2 of the 60 patients (3.3%) required a rescue ERCP



Figure 1. Single-use duodenoscope EXALT Model D with image processor (white box).




Figure 2. Endoscopy image of the major papilla.



Figure 3. A self-expandable metal stent is placed in the bile duct.

with a conventional duodenoscope^[8]. These results, along with other subsequent studies, support a similar efficacy of non-reusable duodenoscopes^[9-10].

2.) Single-use duodenoscopes are equally effective in low and in high-grade complexity procedures?: From a clinical perspective, this is an important point. One may question whether these good results of single-use duodenoscopes will remain when difficult

cases are performed. Although results in this area are still limited, we have strong evidence in the literature demonstrating in a prospective randomized controlled trial that single-use duodenoscopes have no significant differences in technical performance and safety compared with conventional reusable duodenoscopes in cases of low complexity ERCP (Grade 1-2 of the American Society for Gastrointestinal Endoscopy (ASGE) Complexity Scale)(9-10). Other non-comparative studies showed that 3.3% of cases

There are 4 specific groups of potential indications that would be wise to bear in mind when a decision between single-use or reusable duodenoscopes has to be taken

required a switch from single-use duodenoscope to a reusable duodenoscope because of ERCP failure in patients with a grade of complexity of 2 or 3 (8-9). Although some high complexity ERCP cases (Grade 4) have been reported with remarkable success, the performance characteristics of single-use duodenoscopes in such cases is yet to be proven.

3.) Clinical indications for single-use duodenoscopes: Although still in the process of being fully defined, there are 4 specific groups of potential indications that would be wise to bear in mind when a decision between single-use or reusable duodenoscopes has to be taken:

3.a.) Patients with an elevated risk of having a MDR bacteria:^[11]

- ✓ Previous colonization of the patient by a MDR.
- ✓ Patient admitted or with a recent admission to the Intensive Care Unit (ICU).
- ✓ Patient receiving or having received recently mechanical ventilation.
- ✓ Patients with exposure to broad-spectrum antimicrobial within 30 days (e.g. fluoroquinolones, 3rd generation cephalosporin and/or carbapenem).
- ✓ Recent placement of intravascular devices or urinary probes.
- ✓ Patient who has recently developed nosocomial infections.
- ✓ Patient who has undergone recent surgery or another invasive procedure.
- ✓ Patient with neutropenia (< 500 PMN/mm³).
- ✓ Presence of a local outbreak of MDR.

3.b.) Patients with a higher risk of morbidity-mortality if infected by MDR or other bacteria:

- ✓ Aged-patients.
- ✓ Patients with a long-hospitalization.
- ✓ Oncologic patients.
- ✓ Immune suppressed patients.
- ✓ Transplanted patients.
- ✓ Patients affected by an autoimmune disease.

3.c.) Patients with a higher risk of developing an infection after ERCP performance:

- ✓ ERCP placement of a biliary stent.
- ✓ Patients with a diagnosis of cholangiocarcinoma.
- ✓ Incomplete drainage of the bile duct.
- ✓ Hospitalized patients.

3.d.) Logistic indications: (Potential indications)

- ✓ ERCP has to be performed outside the Endoscopy Unit (Operation Room, Emergency Department, ICU,...).
- ✓ ERCP performed in non-laboring hours (cleaning of endoscopes not available).
- ✓ Reusable duodenoscope under repair or in quarantine.
- ✓ Low-volume ERCP centers (avoids large investment associated with buying a reusable duodenoscope).

4.) Alternative strategies of use for disposable duodenoscopes: Three potential alternative strategies for informed use of duodenoscopes are presented:

4.a.) Replacement of all reusable

duodenoscopes by single-use endoscopes:

This approach could be the best from a patient's perspective. Initial cost analysis studies (including durability of the duodenoscope, cost of the duodenoscope, washer and cleaning supplies, labor costs and cost of treatment of infections) suggest an extra cost of \$612 per procedure in high-volume institutions (a little more in low volume centers) if all reusable endoscopes are replaced by single-use duodenoscopes^[12]. The affordability of this strategy by the healthcare system appears to be difficult at present time.

4.b.) Selective use of disposable

duodenoscopes:

It seems a wise decision to pursue this more economically affordable strategy, reserving the economic resources available for those patients at an increased risk of being infected and for those with a major probability of being carriers of MDR bacteria. A marginal increase in healthcare cost (e.g. 5–10% of ERCP patients) may be better assumed by sanitary authorities.

4.c.) Use of duodenoscopes with a disposable cup:

Improvements in duodenoscope design are likely to reduce the risk of infection, and partially disposable duodenoscopes should replace current non-disposable ones. However, this alternative appears to be insufficient as supported by the 37% rate of persistent organic residue contamination after reprocessing these partially disposable duodenoscopes^[13].

5.) Environmental impact of disposable duodenoscopes:

Although not the subject of this review, it appears clear that any modification in our clinical practice, such as transitioning into single-use endoscopes, may affect our environmental impact, as suggested by recent publications^[14]. A reasonable balance between potential patient benefits, cost and environmental impact appears to be mandatory.

In summary, single-use duodenoscopes are certainly a key stone to avoid infections related to ERCP. Preliminary reports support high efficacy and safety, but efforts are still needed to better define indications and reduce the cost associated with incorporation into routine clinical practice.

References:

- [1] Ha J, Son BK. Current issues in duodenoscope-associated infections: now is the time to take action. Clin Endosc 2015;48:361-363. doi: 10.5946/ce.2015.48.5.361.
- [2] Aumeran C, Poincloux L, Souweine B, et al. Multidrug-resistant klebsiella pneumoniae outbreak after endoscopic retrograde cholangiopancreatography. Endoscopy 2010;42:895-899. doi: 10.1055/s-0030-1255647.
- [3] Epstein L, Hunter JC, Arwady MA, et al. New Delhi metallo- β -lactamase-producing carbapenem-resistant escherichia coli associated with exposure to duodenoscopes. JAMA 2014;312:1447-1455. doi: 10.1001/jama.2014.12720.
- [4] Naryzhny I, Silas D, Chi K. Impact of ethylene oxide gas sterilization of duodenoscopes after a carbapenem-resistant Enterobacteriaceae outbreak. Gastrointest Endosc 2016;84:259-262. doi: 10.1016/j.gie.2016.01.055.
- [5] Vickery K, Ngo Q-D, Zou J, Cossart YE. The effect of multiple cycles of contamination, detergent washing, and disinfection on the development of biofilm in endoscope tubing. Am J Infect Control 2009;37:470-475. doi: 10.1016/j.ajic.2008.09.016.
- [6] Verfaillie CJ, Bruno MJ, Voor in 't Holt AF, et al. Withdrawal of a novel-design duodenoscope ends outbreak of a VIM-2-producing pseudomonas aeruginosa. Endoscopy 2015;47:493-502. doi: 10.1055/s-0034-1392080.
- [7] Trindale AJ, Copland A, Bhatt A, et al. Single-use duodenoscopes and duodenoscopes with disposable end caps. Gastrointest Endosc 2021;93:997-1005. doi: 10.1016/j.gie.2020.12.033.
- [8] Muthusamy VR, Bruno MJ, Kozarek RA, et al. Clinical evaluation of a single-use duodenoscope for endoscopic retrograde cholangiopancreatography. Clin Gastroenterol Hepatol 2020;18:2108-2117.e3. doi: 10.1016/j.cgh.2019.10.052.
- [9] Complexity of ERCP: Cotton PB, Eisen G, Romagnuolo J, et al. Grading the complexity of endoscopic procedures: results of an ASGE working party. Gastrointest Endosc 2011;73:868-874. doi: 10.1016/j.gie.2010.12.036.
- [10] Bang JY, Hawes R, Varadarajulu S. Equivalent performance of single-use and reusable duodenoscopes in a randomised trial. Gut 2021;70:838-44. doi: 10.1136/gutjnl-2020-321836.
- [11] Di Franco S, Alfieri A, Pace MC, et al. Blood Stream Infections from MDR Bacteria. Life (Basel) 2021;11(6):575. doi: 10.3390/life11060575.
- [12] Bang JY, Sutton B, Hawes R, Varadarajulu S. Concept of disposable duodenoscope: at what cost?. Gut 2019;68:1915-1917. doi: 10.1136/gutjnl-2019-318227.
- [13] Ridditid W, Pakvisal P, Chatsuwat T, et al. A newly designed duodenoscope with detachable distal cap significantly reduces organic residue contamination after reprocessing. Endoscopy 2020;52:754-60. doi: 10.1055/a-1145-3562.
- [14] Namburam S, von Renteln D, Damianos J, Bradish et al. Estimating the environmental impact of disposable endoscopic equipment and endoscopes. Gut 2021;0:1–6. doi:10.1136/gutjnl-2021-324729

In summary, single-use duodenoscopes are certainly a key stone to avoid infections related to ERCP. Preliminary reports support high efficacy and safety, but efforts are still needed to better define indications and reduce the cost associated with incorporation into routine clinical practice

Improving Infection Prevention in Endoscopy and the Role of Single-Use Technologies – Expert Panel Meeting

Evidence, Experience and Expert Opinion

Highlights of the Expert Panel Meeting for the Middle East and North Africa, 15 September 2021*

**Alessandro Repici, Director of Digestive Endoscopy Unit; Professor of Gastroenterology¹

V Raman Muthusamy, Director of Endoscopy; Professor of Clinical Medicine²

Bader Ahmad S Alajlan, Consultant in Therapeutic and Interventional Endoscopy³

Ali Abdulatif Alali, Assistant Professor, Gastroenterology and Hepatology, Specialist Gastroenterology, Hepatology and Therapeutic Endoscopy⁴

Sameer AbdulRaheim Al-Awadhi, Head of Gastroenterology Unit and Senior Lecturer⁵

Abdulrahman Abdulaziz Al Fadda, Director, Advanced Therapeutic Endoscopy Program Gastroenterology Section, Department of Medicine; Director, Endoscopy Unit; Chairman, Gastroenterology and Hepatology Exam Board Committee; Lecturer, College of Medicine⁶

Ahmad Waseil Al Harbi, Consultant Gastroenterologist and Advanced Therapeutic Endoscopist⁷

Emad Suliman AlJahdali, Assistant Professor of Internal Medicine, Gastroenterology Fellowship Program Director⁸

Abed Al-Lehibi, Director, Advanced Therapeutic Endoscopy Fellowship Program, Chairman Gastroenterology / Hepatology Department⁹

Mohammad Al Beshir, Consultant Gastroenterologist & Therapeutic Endoscopist, Director, Gastroenterology Shared Training Program - SCFHS¹⁰

Mohamed Elkady, Consultant and Head of Gastroenterology, Hepatology and Endoscopy Department; Lecturer in Gastroenterology and Hepatology Department¹¹

Ahmed Abd El-Azeem El-Mikkawy, Consultant of Hepato-gastroenterology and Liver Transplantation; Associate Professor Gastroenterology & Hepatology Department¹²

Khalid Osman Elamin Elsayed, Consultant Gastroenterology and Director of Endoscopy¹³

Ayman Yosry Abdel-Rahim Hamad, Emeritus Professor of Gastroenterology and Hepatology; Chairman Faculty of Medicine¹⁴

Ibrahim Mostafa, Founder and CEO; Lecturer of Gastroenterology & Hepatology¹⁵

Data from global registries and reports have shown that patients with IBD do not present a higher risk of infection or a worse course of COVID-19 compared to the general population

*Virtual meeting sponsored by Boston Scientific

**Corresponding author

¹ Digestive Endoscopy Unit, Division of Gastroenterology, Humanitas Research Hospital; Humanitas University, Rozzano 20093, Italy

² UCLA Health System; Vatche and Tamar Manoukian Division of Digestive Diseases; David Geffen School of Medicine, University of California, Los Angeles, Los Angeles, California, USA

³ King Faisal Specialist Hospital & Research Centre, Riyadh, Saudi Arabia

⁴ Department of Medicine, Faculty of Medicine, Kuwait University, Kuwait

⁵ Rashid Hospital, Dubai Health Authority, Dubai Medical College, Dubai, United Arab Emirates

⁶ King Faisal Specialist Hospital and Research Center (KFHS&RC); Saudi Commission of Health Specialties (SCHS); Alfaisal University, Riyadh, Kingdom of Saudi Arabia

⁷ King Faisal Specialist Hospital & Research Center Jeddah, Saudi Arabia

⁸ King Abdul Aziz University Hospital, Jeddah, Kingdom of Saudi Arabia

⁹ King Fahad Medical City, Riyadh, Kingdom of Saudi Arabia

¹⁰ Division of Gastroenterology, Dept of Medicine, King Fahad Specialist Hospital, Chief, Clinical Training, Chief, Leadership Training Academy, Eastern Health Cluster, Dammam, Saudi Arabia

¹¹ Dar Alfouad Hospital, Nasr City, Cairo; Theodor Bilharz Research Institute, Giza, Egypt


¹² Dar Alfouad Hospital, Nasr City, Cairo; Theodor Bilharz Research Institute, Giza, Egypt

¹³ Burjeel Hospital, Abu Dhabi, United Arab Emirates

¹⁴ Cairo University; Cairo University Center for Hepatic Fibrosis Research (CUC-HF), Cairo, Egypt

¹⁵ ROEYA Hepatology, Gastroenterology and Endoscopy Center; Theodor Bilharz Research Institute, Giza, Egypt

Requisite events for infection transmission via duodenoscope and associated probabilities/rates



CRE, carbapenem-resistant *Enterobacteriaceae*; ERCP, endoscopic retrograde cholangiopancreatography
Figure adapted from Kim S, Russell D, Mohamadnejad M, et al. Risk factors associated with the transmission of carbapenem-resistant *Enterobacteriaceae* via contaminated duodenoscopes. *Gastrointest Endosc*. 2016;83(6):1121–9

Introduction

The burden of healthcare-acquired infections (HAIs) or nosocomial infections (NIs) is largely unrecognized and under-reported. HAIs are the sixth leading cause of death, surpassing the combined deaths from HIV/AIDS, cancer and traffic accidents^[1,2].

Data from the USA show that more than 100,000 hospitalized patients die every year due to NIs, with the direct economic burden of HAIs estimated at \$28–\$45 billion^[1]. A report by the Organization for Economic Cooperation and Development (OECD) mentions that a single resistant infection has cost of about €8,500 to €34,000 more than a non-resistant infection^[3].

While reporting from Western countries is more consistent, there are studies indicating that bacterial profiles in the Middle East are similar^[4,5], which suggests a comparable situation. In one study at King Abdulaziz Specialist hospital in Saudi Arabia, *Klebsiella pneumoniae* was identified in 30 patients^[6]. *K. pneumoniae* is an opportunistic pathogen responsible for a significant proportion of nosocomial and community-acquired infections, and its genotypic variation is a major barrier to controlling public health risks associated with pathogens^[6].

Infection outbreaks have a significant impact on patients and institutions. Those involving multidrug-resistant organisms (MDRO) have the greatest clinical consequence and are also the most easily tracked. Such outbreaks are currently rare, but the era of MDRO is only beginning and is likely to worsen in future. For an affected unit, any outbreak can create poor morale, temporarily reduce interventional capacity and surgical volume, increase the costs for additional equipment, cleaning and legal fees; and damage the reputation of the institution^[7].

Patient and Device Factors Related to Duodenoscope-Transmitted Infections

Up to 8% of patients undergoing endoscopic retrograde cholangiopancreatography (ERCP) may develop an infection due to a contaminated duodenoscope^[7].

The main patient-related risk factors for duodenoscope-transmitted infections were inpatient status and biliary stent placement in a study comparing 15 patients with active infection or colonization against 89 who had been exposed but had no infection. Additionally, biliary stent placement and male sex were associated with active infection compared to colonization^[7]. In an earlier study of 39 patients, of whom 35 had duodenoscope exposure, risk factors for transmission included stent placement and brushings^[8].


Infection outbreaks are also associated with a combination of device-related factors^[9,10]. Firstly, the duodenoscope has a complex design, with the distal tip and elevator making it difficult to achieve complete clearance of debris and contaminated fluids^[2,11,12].

Secondly, devices used during ERCP create higher levels of mechanical stress than gastroscopes or colonoscopes, with subsequent permanent damage that favors bacterial contamination and creates the conditions for biofilm formation^[10]. Foreign debris such as tissue, fibres, and metal clips are frequently retrieved from the working channel of the scope – in 48% of inspected duodenoscopes (N=50) and 59% of all flexible endoscopes (N=126) in one survey^[13].

Following large outbreaks in three major Dutch medical facilities, Rauwers et al^[14] identified a range of device abnormalities, including: brown staining on the distal tips; cracking at the distal end on the inside of the biopsy channel; inadequate connection of the replaced protective cap, leaving space between the cap and the tip frame; and white and brown oxidation stains on the elevator^[14].

Single-use technologies have become available and may provide a solution to the complex problem of device-associated infection outbreaks^[15,16]. The world's first single-use model, EXALT D (Boston Scientific Corporation, Marlborough Mass, USA) reaches equivalent performance to reusable duodenoscopes while eliminating the risk of bacterial colonization due to ineffective processing, and the need for reprocessing and repairs^[16].

Cleaning duodenoscopes is a complex process, with no single point to confirm that infection risk is truly controlled¹²



The impact of single-use devices on the environment is also a consideration, but must be balanced against the waste produced by reusable devices and the disinfection process, as well as the overall risk reduction in specific patients. Importantly, in a study comparing single-use and reusable ureteroscopes, Davis et al (2018) showed that the environmental impact of a single-use flexible endoscope is comparable to that of the reusable flexible endoscope^[17].

Considerations for Patient Selection: Center Experience with Single-Use Duodenoscopes


There are now over 400 published cases for the Exalt D single-use duodenoscope showing levels of satisfaction and adverse event rates comparable to reusable duodenoscopes. The device is under constant review and a new iteration is now available, with a particular focus on ergonomics and hand experience. However, a key question remains around patient selection for single-use technologies.

The Humanitas University Hospital, Milan, Italy

According to the protocols at Humanitas University Hospital, Exalt D is routinely used in three select indications:

- sick and frail patients with advanced cancer, or cancer suitable for surgery, who are jaundiced and require simple, complication-free decompression
- post-transplant patients, including those who have had a stem cell transplant
- patient cases involving the use of SpyGlass™ DS System.

For every case where a single-use duodenoscope might be indicated, the team discusses the advantages and disadvantages of the procedure, as well as its potential impact on the patient and budget in the context of risk for post-ERCP infection complications.



UCLA Health System, Los Angeles, USA

At UCLA Health System, the team uses four criteria for deciding on single-use technologies with the hospital analysis committee on price:


- to facilitate logistics/operations, e.g. weekend/night cases; operating room cases; backup fleet on busy days; low-volume/off-site centers
- to protect the fleet, e.g. in patients with known MDROs^[22] to avoid the risk of transmission, which has been estimated around 14%^[7]
- in critically ill patients who would not be able to withstand a duodenoscope-acquired infection; these include hospitalized patients, in particular those receiving transplantation and intensive care unit patients with pressor support
- in procedures at high risk for bacteremia/infection transmission, such as cholangioscopy, biliary radiofrequency ablation, and stent placement

The experience to date at UCLA Health System suggests that patient preference may also come to play a role in the choice of single-use duodenoscope. Emerging data on cost considerations related to an enhanced

An expert panel was convened to share experiences with the Exalt D, review important factors impacting its use, and identify those patients who will benefit the most from a single-use versus reusable duodenoscope

Exalt D Single-Use Duodenoscope: Published Data on Efficacy, Comparability and Safety

Publication	Number of patients	Key outcomes
Muthusamy et al. 2020 (18)	60	97% success with SUD; AE rate 6.6%; median overall satisfaction 9/10
Bang et al. 2021 (19)	48 (all native papilla)	96% success; no difference in cannulation success, procedure time or AE rates (4% for SUD) compared to reusable devices
Slivka et al. 2021 (20)	200 (n=19 endoscopists)	40.5% ASGE Grade 3-4; crossover rate 11.3% vs. 2.5% (N=19); success 97%; AE rate 6.5%; similar time and satisfaction between experts and "less-experts"
Napoléon et al. 2021 (21)	60	95% success; median satisfaction 9/10; AE rate 5% (not SUD related)
Ehrlich et al. 2021 (16)	36	91.7% success rate; AE rate 3.2%



duodenoscope reprocessing protocol suggest a cost efficacy for single-use duodenoscopes that may even lead to savings in certain settings^[23-25].

Patient Selection: Perspectives from Middle Eastern Centers

Defining the Issue of Duodenoscope-Associated Infection

The group agreed that the risk of duodenoscope-associated infection is increasingly apparent globally. However, many centers do not track the issue, while others are unable to do so because the patients return to other medical facilities after undergoing ERCP in the center.

Demonstrating the true extent of the problem will be challenging until active screening is implemented. In the Middle East, there is a general acceptance that MDROs are a growing issue, but the link has not yet been made to ERCP or duodenoscopes.

Experts suggested there should be more active screening of duodenoscopes and greater follow-up of ERCP patients, especially oncology patients. The group expressed enthusiasm for a multi-center prospective study to routinely culture duodenoscopes and understand the contamination rates. Protocols could be drawn up and a central laboratory used to ensure consistency and enable the collection of meaningful region-specific data.

Experience

There was some experience within the group with the Exalt D single-use duodenoscope, which was reportedly comparable with the current gold standard.

Some participants raised concerns that Exalt D might handle differently, but it was noted that Pentax™ and Olympus™ do not handle in the same way either, and the more experienced in the group advised that after 5–10 cases with the Exalt D single-use duodenoscope, they were no longer conscious of the difference.

It was observed that different centers may have different processing protocols and disinfection standards, and legal challenges were also discussed. Single-use technology may head off any legal challenges arising around reprocessing.

The group agreed that infection risk with duodenoscopes must be addressed in clinical practice. Since contamination was shown to occur in the duodenoscope channel as well as in the tip^[26], single-use tips may not provide a full solution. A 'clean' device is also challenging to prove. All these factors point to a strong argument for eliminating infection risk by investing in single-use duodenoscopes.

Indications for the Exalt D Single-Use Duodenoscope

All centers face budget limitations, which may affect how often single-use duodenoscopes are deployed. However, the group felt there is a clear indication for their use in post-transplant patients, those with compromised immunity and/or known MDRO colonization.

In some private centers in the Middle East, patients have requested a duodenoscope solely for their use, and patient preference may become a strong driver for single-use devices.

There is a need for a convincing argument for reimbursing single-use devices. Currently under a procedure tariff system, the levels of funding are insufficient. In all cases, until appropriate guidelines are developed, there needs to be an individualized discussion with payers to argue the case for single-use devices by highlighting patients with comorbidities, cancer or undergoing multiple medical therapies.

Results from clinical studies are not predictive of results in other studies. Results in other studies may vary.

All trademarks are the property of their respective owners. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labelling supplied with each device. Products shown for INFORMATION purposes only and may not be approved for sale in certain countries. This material is not intended for use in France. 2022 Copyright © Boston Scientific Corporation or its affiliates. All rights reserved. ENDO-12424022-AA

Expert Consensus

The experts agreed fully that the single-use device should be used after organ transplantation, in patients at high risk of transmission due to known MDRO colonization, and in those with the highest risk of infection, e.g. immunocompromised patients. Consensus was also strong that the evidence and experience from the Middle East region should be collated, and considered together with the identified indications for single-use duodenoscope deployment, to strengthen the case for reimbursement.

Acknowledgements

Medical writing and editorial support was provided by Carys Thomas Ampofo and Ileana Stoica and was funded by Boston Scientific.

Authors' Contributions

All authors contributed to the interpretation of findings and writing and reviewing of this manuscript. All authors meet the International Committee of Medical Journal Editors (ICMJE) criteria for authorship for this article, take responsibility for the integrity of the work as a whole, and have given their approval for this version to be published.

Funding

The funding of this consensus paper was provided by Boston Scientific, Massachusetts, USA.

It is hard to argue against a single-use duodenoscope in the context of infection, which can be a human as well as a financial disaster

It was observed that different centers may have different processing protocols and disinfection standards, and legal challenges were also discussed. Single-use technology may head off any legal challenges arising around reprocessing

Data from global registries and reports have shown that patients with IBD do not present a higher risk of infection or a worse course of COVID-19 compared to the general population

References:

[1] Haque M, Sartelli M, McKimm J, Abu Bakar M. Health care-associated infections - an overview. Infect Drug Resist. 2018;11:2321-33.

[2] Rubin ZA, Kim S, Thaker AM, Muthusamy VR. Safely reprocessing duodenoscopes: current evidence and future directions. Lancet Gastroenterol Hepatol. 2018;3(7):499-508.

[3] OECD. HEALTH AT A GLANCE: EUROPE 2018. Available from: https://www.oecd-ilibrary.org/social-issues-migration-health/health-at-a-glance-europe-2018_health_glance_eur-2018-en; Accessed November 2021.

[4] Oldenkamp R, Schultz C, Mancini E, Cappucco A. Filling the gaps in the global prevalence map of clinical antimicrobial resistance. Proc Natl Acad Sci U S A. 2021;118(1).

[5] Dandachi I, Chaddad A, Hanna J, Matta J, Daoud Z. Understanding the Epidemiology of Multi-Drug Resistant Gram-Negative Bacilli in the Middle East Using a One Health Approach. Front Microbiol. 2019;10:1941.

[6] Lagha R, Ben Abdallah F, AAH AL, Amor N, Hassan MM, Mabrouk I, et al. Molecular characterization of multidrug resistant Klebsiella pneumoniae clinical isolates recovered from King Abdulaziz Specialist Hospital at Taif City, Saudi Arabia. J Infect Public Health. 2021;14(1):143-51.

[7] Kim S, Russell D, Mohamadnejad M, Makker J, Sedarat A, Watson RR, et al. Risk factors associated with the transmission of carbapenem-resistant Enterobacteriaceae via contaminated duodenoscopes. Gastrointest Endosc. 2016;83(6):1121-9.

[8] Epstein L, Hunter JC, Arwady MA, Tsai V, Stein L, Gribogiannis M, et al. New Delhi metallo-beta-lactamase-producing carbapenem-resistant Escherichia coli associated with exposure to duodenoscopes. JAMA. 2014;312(14):1447-55.

[9] Rauwers AW, Voor In't Holt A, De Groot R. ERCP duodenoscopes in Dutch ERCP centres: high prevalence of bacterial contamination despite reprocessing. Presented at ESCMID O3842016.

[10] Beilenhoff U, BIERING H, Blum R, Brljak J, Cimbrow M, Dumonceau JM, et al. Prevention of multidrug-resistant infections from contaminated duodenoscopes: Position Statement of the European Society of Gastrointestinal Endoscopy (ESGE) and European Society of Gastroenterology Nurses and Associates (ESGENA). Endoscopy. 2017;49(11):1098-106.

[11] Beilenhoff U, BIERING H, Blum R, Brljak J, Cimbrow M, Dumonceau JM, et al. Reprocessing of flexible endoscopes and endoscopic accessories used in gastrointestinal endoscopy: Position Statement of the European Society of Gastrointestinal Endoscopy (ESGE) and European Society of Gastroenterology Nurses and Associates (ESGENA) - Update 2018. Endoscopy. 2018;50(12):1205-34.

[12] Boston Scientific. Eliminate cross contamination due to ineffective reprocessing. 2021. Available from: <https://www.bostonscientific.com/en-US/medical-specialties/gastroenterology/exalt/infection-prevention.html>; Accessed December 2021.

[13] Stephenson K. Visual inspection of flexible endoscope working channels. Vend Vantage. 2016. Available from: <http://inspektor.com/assets/visual-inspection-of-flexible-endoscope-working-channels---communique---may--june-2016.pdf>; Accessed November 2021.

[14] Rauwers AW, Troelstra A, Fluit AC, Wissink C, Loeve AJ, Vleggaar FP, et al. Independent root-cause analysis of contributing factors, including dismantling of 2 duodenoscopes, to investigate an outbreak of multidrug-resistant Klebsiella pneumoniae. Gastrointest Endosc. 2019;90(5):793-804.

[15] Kozarek RA. The Past, Present, and Future of Endoscopic Retrograde Cholangiopancreatography. Gastroenterol Hepatol (N Y). 2017;13(10):620-2.

[16] Ehrlich D, Muthusamy VR. Device profile of the EXALT Model D single-use duodenoscope for endoscopic retrograde cholangiopancreatography: overview of its safety and efficacy. Expert Rev Med Devices. 2021;18(5):421-7.

[17] Davis NF, McGrath S, Quinlan M, Jack G, Lawrentschuk N, Bolton DM. Carbon Footprint in Flexible Ureteroscopy: A Comparative Study on the Environmental Impact of Reusable and Single-Use Ureteroscopes. J Endourol. 2018;32(3):214-7.

[18] Muthusamy VR, Bruno MJ, Kozarek RA, Petersen BT, Pleskow DK, Sejpal DV, et al. Clinical Evaluation of a Single-Use Duodenoscope for Endoscopic Retrograde Cholangiopancreatography. Clin Gastroenterol Hepatol. 2020;18(9):2108-17 e3.

[19] Bang JY, Hawes R, Varadarajulu S. Equivalent performance of single-use and reusable duodenoscopes in a randomised trial. Gut. 2021;70(5):838-44.

[20] Slivka A, Ross AS, Sejpal DV, Petersen BT, Bruno MJ, Pleskow DK, et al. Single-use duodenoscope for ERCP performed by endoscopists with a range of experience in procedures of variable complexity. Gastrointest Endosc. 2021;94(6):1046-55.

[21] Napoléon B, Gonzalez J, Grandval P, et al. Evaluation of the performance of a single-use duodenoscope: a prospective multicentre French study (the Exaltes study). Endoscopy. 2021;53(S01):S16.

[22] Thornhill G, Dunkin B. The ERCP patient: risk factors for infection. 2021. Available from: <https://www.bostonscientific.com/en-EU/medical-specialties/gastroenterology/exalt/exalt-clinical1/patient-risk-factors/whitepaper.html>; Accessed October 2021.

[23] Bomman S, Kozarek RA, Thaker AM, Kodama C, Muthusamy VR, Ross AS, et al. Economic burden of enhanced practices of duodenoscopes reprocessing and surveillance: balancing risk and cost containment. Endosc Int Open. 2021;9(9):E1404-E12.

[24] Das A, Cangelosi M, Muthusamy V. Techniques and innovations in GI endoscopy. 2021. To be submitted for publication.

[25] Muthusamy V, Ross A. Economic burden of emergent practices of Duodenoscopes reprocessing and surveillance: balancing risk- and cost containment. Gastrointestinal Endoscopy. 2018;87:AB167.

[26] Boston Scientific. Duodenoscope hot spots. 2020. Available from: https://www.bostonscientific.com/content/dam/bostonscientific/endo/general/exalt/Duodenoscope_hot_spots.pdf; Accessed October 2021.

The Use of Single-Use Duodenoscopes in Endoscopy – a Current Perspective

Author

Mark Ellrichmann, Claudio C. Conrad

Corresponding Author

Prof. Dr. Mark Ellrichmann

Direct Department of Interdisciplinary Endoscopy

University Hospital Schleswig-Holstein, Campus Kiel

Arnold-Heller-Str. 3, Haus C

24159 Kiel, Germany

Mail: mark.ellrichmann@uksh.de

Tel: +49-431-500-22371

Fax: +49-431-500-22378

Key words

Single-use duodenoscopes, infection prevention, ERCP

Introduction

Clinically relevant infections caused by the use of contaminated, flexible endoscopes have become the focus of clinical and scientific interest in the literature. In this regard, endoscopy-associated infections can be differentiated as (i) endogenous and (ii) exogenous infections. Endogenous infections are caused by the spread of the patient's own local flora and mainly involve *Escherichia coli*, *Klebsiella*, *Enterobacter* and *Enterococci* species. These infections can cause issues, such as pneumonia due to aspiration of secretions in sedated patients and the spread of microorganisms from the oral cavity during endoscopic procedures. Furthermore, cholangitis and bacteremia have been reported in patients with biliary obstruction after endoscopic retrograde cholangiopancreatography (ERCP)^[1-3]. Although endogenous infections that are directly related to an endoscopic procedure can occur, they cannot be prevented by strictly controlled disinfection, sterilization or the use of disposable materials. In contrast, exogenous infections from consumables or endoscopes that are contaminated with microorganisms originating from

other patients can be significantly reduced or even prevented by the highest hygienic standards as well as the use of single use endoscopes and consumables.

There is currently no credible estimate of infection rate associated with endoscopic procedures. Despite the very high number of flexible endoscopies worldwide, the published number of iatrogenic, endoscopy-associated infections is rather low. However, the actual transmission rate is vastly underestimated because of inadequate surveillance strategies and lack of clinical symptoms^[4]. In the early 1990s, Gorse and Messner estimated an infection rate of 6% after gastrointestinal endoscopies^[5]. Furthermore, in the period between 1974 and 2004, 30 infection outbreaks from flexible endoscopies involving a total of 251 patients were identified^[6]. A recent meta-analysis by Kovaleva et al. showed that flexible bronchoscopy and ERCP confer the highest risk of exogenous infection^[4]. The decisive factor that is emerging is the complex geometry of the aforementioned endoscopes with small working channels, particularly that of the distal end, resulting in a residual biofilm that favors further microbial growth even after proper high-level disinfection (HLD).

There is currently no credible estimate of infection rate associated with endoscopic procedures

Hard costs comprise the cost of the SUD itself and the dedicated processor. Furthermore, costs for reprocessing equipment, cleaning solutions and repair have to be taken into account

A recent meta-analysis including 15 studies showed a contamination rate of duodenoscopes of 15.2% after adequate reprocessing. A subgroup analysis including only studies with HLD or sterilization demonstrated a residual contamination rate of 9.2%. However, the proportion of clinically relevant infections in these studies remained unclear^[7]. In a prospective study in 73 ERCP centers, 155 duodenoscopes were tested for residual microorganisms after HLD. Contamination was defined as the detection of >20 colony forming units (CFU)/20ml of sample solution for all microorganisms or >1 CFU/20ml of sample solution for microorganisms originating from the oral cavity or the gastrointestinal tract. Despite adequate HLD, the authors found a contamination rate of reprocessed duodenoscopes of 15%^[10]. Based on the published outbreaks of infection from Europe and the USA over the last five years, an infection rate of about 1.22% can be extrapolated, although the number of unreported cases is presumed to be high^[7-11]. Based on these worrisome data and published case series, the FDA (Federal Drug Administration of the US) recommended in a recent safety warning (24.07.2020; <https://fda.gov/medical-devices/safety-communications/fda-recommending-transition-duodenoscopes-innovative-designs-enhance-safety-fda-safety-communication>) that "Hospitals and endoscopy facilities should transition to innovative duodenoscope designs that include disposable components such as disposable endcaps, or to fully disposable duodenoscopes when they become available.". Consequently, the world's first single-use duodenoscope (SUD), Exalt Model D from Boston Scientific (Boston Scientific, Marlborough, MA, USA) was licensed in December 2019 in the US and January 2020 in Europe. Irrespective of potential hygienic advantages of SUD, three topics need to be discussed to implement SUDs in clinical routine: clinical functionality, cost effectiveness and patient selection for an SUD.

Clinical Functionality of Single-Use Duodenoscopes

When using the SUD for the first time, the endoscope feels much lighter with ergonomics of the handle being similar to the standard reusable Olympus duodenoscopes that are widely available. Of note, the distal end of the scope does not straighten by itself after releasing the locks of the wheels. Therefore, the endoscopist has to actively straighten the scope before withdrawing it from the duodenum (personal experience).

In a first clinical multicenter trial with seven ERCP experts performing 60 ERCPs using

the Exalt Model D, the ERCP was successfully completed in 96.7% of patients. The satisfaction score was rated 9/10 on average without additional adverse events being observed^[12]. A prospective, randomized study evaluated the use of an SUD (EXALT Model D) in comparison to standard reusable duodenoscopes. In 98 cases (48 SUD), no significant differences in the advancement time to the Papilla Vateri or in the successful cannulation of the target area were observed. However, the mean number of cannulation attempts was significantly lower with 6.7 attempts using a SUD compared to 12.7 with reprocessable duodenoscopes, which resulted in a significantly shorter time to successful cannulation. The authors characterized disposable duodenoscopes as offering a more favorable anatomical position in front of the papilla Vateri with the EXALT Model D^[13]. However, a final evaluation cannot yet be made because of the small number of cases. Clinical functionality was confirmed in a multicenter European study in 60 patients undergoing ERCP with an SUD. A technical success rate of 95% with a median operators' satisfaction rating of 9 (range 7–9) was observed^[14].

Cost Assessment of Single-Use Endoscopes

To assess total costs of SUDs, two groups of cost have to be evaluated: (i) hard costs and (ii) soft costs. Hard costs comprise the cost of the SUD itself and the dedicated processor. Furthermore, costs for reprocessing equipment, cleaning solutions and repair have to be taken into account. Soft costs primarily occur during the use of reusable duodenoscopes, including personnel for reprocessing, costs of endoscopy associated infections (prolonged hospital or ICU stay) and costs for surveillance programs (culturing, scope quarantine, sterilization).

In a theoretical cost calculation model by Bang et al., cost effectiveness of SUD was assessed. They showed that the overall cost per ERCP directly depended on the annual volume of ERCPs per center and the respective infection rates^[15]. Given an infection rate of 1%–1.2%, total costs of an ERCP in a high volume center (>350 ERCP/year) ranged from \$1,110–\$1,338, whereas in low volume centers (<350 ERCP/year), ERCP costs were \$1,220–\$2,685. In cases where SUD is not additionally reimbursed, the cut-off for cost-effectiveness is not currently met.

Due to the differences in international health care systems, reimbursement strategies for ERCP and SUD cannot be directly compared and need to be evaluated in further detail from a national perspective. In the US, 100% of the SUD costs are reimbursed in an outpatient setting by a transitional pass-through payment (TPT). In Germany, the

situation is completely different and is based on the DRG system (diagnosis related groups), where hospitals receive a fixed reimbursement depending on the level of clinical complexity of the procedure independently of the costs of the deployed consumables. So far, only an additional procedure code for SUD has been established, which currently does not result in an additional reimbursement.

Patient Selection for Single-Use Duodenoscopes

Based on the authors' experience and recently supported by a publication by Gromski et al., three tiers for patient selection for SUD can be defined^[17]: recommend the use of SUD, consider the use of SUD and do not recommend the use of SUD.

1. Recommend the use of SUD:

In this group, patients colonized with multi-drug resistant organisms (MDRO) or at risk of MDRO including those with (i) prior carbapenem-resistant Enterobacteriaceae (CRE) infection anywhere; (ii) prior multi-drug resistant infection with bile e.g. CRE, vancomycin-resistant Enterobacteriaceae (VRE), extended spectrum beta-lactamase-producing bacteria or (iii) high risk of recurrent cholangitis, e.g. primary sclerosing cholangitis (PSE), hilar malignancies with recurrent cholangitis or incomplete drainage of high grade strictures.

2. Consider the use of SUD

The second group consists of proven or potentially immunocompromised patients including (i) patients on immunosuppressive medication, e.g. those with liver transplant, PSC with ulcerative colitis or on biological therapy in general; (ii) patients with malignancies undergoing chemotherapy or (iii) patients with special circumstances of immune suppression, e.g. those with HIV infection, neutropenia or bone marrow transplantation.

3. Do not recommend the use of SUD

In all patients that do not belong to group 1

or 2, the use of SUD is not recommended. Furthermore, the use of SUD in weekend procedures and in low-volume centers should be discussed on an individual basis^[19].

Conclusion

SUDs are, by definition, always sterile, which eliminates the risk of contamination due to ineffective reprocessing, minimizes the risk of infection with exogenous pathogens and thus contributes to improving patient safety in flexible endoscopy. Assuming that supply chains from manufacturers to the end users reliably function, disposable endoscopes are always available. Thus, depending on the structure of the endoscopy unit, the use of disposable endoscopes frees up considerable staff resources because of elimination of reprocessing. These staff members are available for the core activity of direct patient care and assistance, which contributes to optimizing process quality.

The additional costs of consumables for disposable endoscopes are seen as a significant disadvantage, especially with the currently varying reimbursement systems. Moreover, the use of additional consumables produces a significant amount of additional contaminated hospital waste. Optimization of the recycling of this waste must be established to ensure sustainability. Focusing on the economic burden of SUD, the total costs of SUD is "optimized" by including the costs of clinically relevant infections. However, the threshold for cost neutrality has not yet been reached for the majority of hospitals. When it comes to the question of "paying it off", economic factors alone do not make up the total bill. Rather, the maximum safety and well-being of our patients should be at the center of our medical and ethical interest.

The funding of this editorial was provided by Boston Scientific.

Assuming that supply chains from manufacturers to the end users reliably function, disposable endoscopes are always available

When it comes to the question of "paying it off", economic factors alone do not make up the total bill. Rather, the maximum safety and well-being of our patients should be at the center of our medical and ethical interest

References:

- [1] Spach, D. H., Silverstein, F. E. & Stamm, W. E. Transmission of infection by gastrointestinal endoscopy and bronchoscopy. *Annals of Internal Medicine* (1993) doi:10.7326/0003-4819-118-2-199301150-00008.
- [2] Srinivasan, A. Epidemiology and prevention of infections related to endoscopy. *Current Infectious Disease Reports* (2003) doi:10.1007/s11908-003-0088-5.
- [3] Nelson, D. B. & Muscarella, L. F. Current issues in endoscope reprocessing and infection control during gastrointestinal endoscopy. *World Journal of Gastroenterology* (2006) doi:10.3748/wjg.v12.i25.3953. *Gastrointest. Endosc.* (1993) doi:10.1016/S0016-5107(93)70316-8.
- [4] Kovaleva, J., Peters, F. T. M., van der Mei Mei, H. C. & Degener, J. E. Transmission of infection by flexible gastrointestinal endoscopy and bronchoscopy. *Clinical Microbiology Reviews* (2013) doi:10.1128/CMR.00085-12.
- [5] Gorse, G. J. & Messner, R. L. Infection control practices in gastrointestinal endoscopy in the united states: A national survey. *Gastroenterol. Nurs.* (1991) doi:10.1097/00001610-199110000-00003.
- [6] Seoane-Vazquez, E., Rodriguez-Monguió, R., Visaria, J. & Carlson, A. Exogenous endoscopy-related infections, pseudo-infections, and toxic reactions: Clinical and economic burden. *Current Medical Research and Opinion* (2006) doi:10.1185/030079906X121048.
- [7] Larsen, S. et al. Rate and impact of duodenoscope contamination: A systematic review and meta-analysis. *EClinicalMedicine* (2020) doi:10.1016/j.eclinm.2020.100451.
- [8] Rauwers, A. W. et al. High prevalence rate of digestive tract bacteria in duodenoscopes: A nationwide study. *Gut* (2018) doi:10.1136/gutjnl-2017-315082.
- [9] Rauwers, A. W. et al. Independent root-cause analysis of contributing factors, including dismantling of 2 duodenoscopes, to investigate an outbreak of multidrug-resistant *Klebsiella pneumoniae*. *Gastrointest. Endosc.* (2019) doi:10.1016/j.gie.2019.05.016.
- [10] Forbes, N. et al. Characteristics and Outcomes of ERCP at a Canadian Tertiary Centre: Initial Results from a Prospective High-Fidelity Biliary Endoscopy Registry. *J. Can. Assoc. Gastroenterol.* (2021) doi:10.1093/jcag/gwaa007.
- [11] Ross, A. S., Baliga, C., Verma, P., Duchin, J. & Gluck, M. A quarantine process for the resolution of duodenoscope-associated transmission of multidrug-resistant *Escherichia coli*. *Gastrointest. Endosc.* (2015) doi:10.1016/j.gie.2015.04.036.
- [12] Muthusamy, V. R. et al. Clinical Evaluation of a Single-Use Duodenoscope for Endoscopic Retrograde Cholangiopancreatography. *Clin. Gastroenterol. Hepatol.* (2020) doi:10.1016/j.cgh.2019.10.052.
- [13] Bang, J. Y., Hawes, R. & Varadarajulu, S. Equivalent performance of single-use and reusable duodenoscopes in a randomised trial. *Gut* (2021) doi:10.1136/gutjnl-2020-321836.
- [14] Napoléon, B. et al. Evaluation of the performances of a single-use duodenoscope: Prospective multi-center national study. *Dig. Endosc. Off. J. Japan Gastroenterol. Endosc. Soc.* 34, 215–221 (2022).
- [15] Bang, J. Y., Sutton, B., Hawes, R. & Varadarajulu, S. Concept of disposable duodenoscope: At what cost? *Gut* (2019) doi:10.1136/gutjnl-2019-318227.
- [16] Barakat, M. T., Ghosh, S. & Banerjee, S. Cost utility analysis of strategies for minimizing risk of duodenoscope related infections. *Gastrointest. Endosc.* (2022) doi:10.1016/j.gie.2022.01.002.
- [17] Gromski, M. A. & Sherman, S. Technological review: developments in innovative duodenoscopes. *Gastrointest. Endosc.* 95, 42–50 (2022).

Hospitals and endoscopy facilities should transition to innovative duodenoscope designs that include disposable components such as disposable endcaps, or to fully disposable duodenoscopes when they become available

Advantages, Environmental Impact and Sustainability of Single-Use Endoscopes

Michael A James, PhD

Interest in the use of disposable, single-use endoscopes has grown because of evident risk of nosocomial infection from reusable endoscopes. While reusable endoscopes are most commonly used, single-use products have been gaining traction since the 1970s. Single-use duodenoscopes (Figure 1) were FDA-cleared for clinical use in 2019. In addition to reduction in risk of infection, elimination of downtime due to reprocessing of reusable endoscopes, decreased requirements for staffing, training and consumables for reprocessing, immediate availability and quality of function are attractive aspects of single-use endoscopes. Even given these attractors to the use of disposable endoscopes, significant thought is required regarding their environmental and economic sustainability.

Advantages

Risk of infection from endoscopic devices continues to be prevalent, particularly from duodenoscopes during endoscopic retrograde cholangiopancreatography (ERCP). This is evidenced by a large number of published studies highlighting outbreaks of multidrug-resistant *Klebsiella pneumoniae*, *Escherichia coli* and *Pseudomonas aeruginosa* infections from duodenoscopy since 2014^[1]. A recent study of ERCP in the Netherlands involving 73 endoscopy centers showed that 22% of patient-ready duodenoscopes were contaminated with bacteria, 15% with microorganisms of gastrointestinal origin^[2]. It is estimated that there are 670,000 infections with drug-resistant bacteria resulting in 33,000 deaths annually in Europe^[3]. Since disposable, single-use duodenoscopes are sterile, it can be assumed that risk of infection due to contamination is eliminated with their use instead of reusable endoscopes. In addition to mitigating risk of nosocomial spread of multidrug-resistant organisms, single-use endoscopes may reduce risk of disease transmission and complication in high-risk



Figure 1. Example of a single-use, disposable duodenoscope. Source: Boston Scientific

patients, such as the immunocompromised and COVID-19-positive patients.

Since turnover time, efficiency, availability and waiting periods for endoscopy procedures are a significant problem for patients, and enhanced sterilization technique recommendations post-COVID-19 have increased reprocessing burden, eliminating reprocessing downtime is a clear advantage of single-use endoscopes. A systematic review of 22 studies across gastrointestinal (GI), urologic and respiratory procedures determined that disposable sheath endoscopes resulted in faster turnover times than reusable components since high-level disinfection (HLD) or sterilization was not required^[4].

Performance Comparison

Both disposable distal endcaps and disposable duodenoscopes have comparable or superior specifications to reusable endoscopes including field-of-view, diameter, bending angles and working length^[5]. Complementary metal oxide semiconductor (CMOS) imaging has become

Multiple studies have demonstrated ≥ 95% success rates with single-use duodenoscopy for ERCP, and one reported equivalent time requirements and satisfaction ratings regardless of level of expertise of the operator

The impact of endoscopy on the environment and climate change is significant. The environmental impact of single-use endoscopes is equivocal in some cases, but depends on the balance of increased manufacturing impact versus decreased impact of reprocessing procedures

relatively inexpensive and improved greatly in terms of resolution in recent decades. Combined with LED lighting, these sensors, employed in single-use endoscopes, provide clear imaging at high resolution.

A single-use duodenoscope for ERCP was found to achieve comparable completion times and comparable performance ratings for overall performance, tip control and image quality versus reusable duodenoscopes^[6]. Navigation/pushability was rated slightly lower for the single-use endoscope. Another clinical study found similar overall safety and technical performance and decreased attempts required for successful cannulation with single-use duodenoscopes compared to reusable endoscopes^[7]. However, ease of passage into the stomach, image stability and quality and air-water button functionality were inferior. Multiple studies have demonstrated $\geq 95\%$ success rates with single-use duodenoscopy for ERCP^[7-10], and one reported equivalent time requirements and satisfaction ratings regardless of level of expertise of the operator^[9]. Continued innovation in imaging and functionality has and will continue to improve visualization using single-use endoscopes. This is true for bronchoscopy as well, which requires clear imaging of lungs and small airways. Several late-generation bronchoscopes feature improved degree of angulation and image quality to overcome these limitations^[11].

Environmental Impact

Reducing or eliminating healthcare-acquired infections (HAIs) is a major benefit of single-use endoscopy; however, potential differences in the impact of these procedure on the environment compared to reusable endoscopy should be considered. The World Gastroenterology Organization's Climate Change Working Group has expressed trepidations about the impact of endoscopy on the environment and set a goal of further assessing the impact of reusable endoscopes^[12]. While disposable endoscopes themselves create consumable waste, this burden is partially mitigated by a decrease in waste from disinfecting consumables needed for the reprocessing of reusable endoscopes. Some studies suggest that full transition to single-use endoscopes would increase net waste^[13], while others indicate that it would result in a lower carbon footprint^[14]. The overall carbon footprint of current endoscopy practice in the US has been estimated to be 86,000 tons of CO₂ equivalents^[15]. A study comparing total lifetime carbon emissions reported 4.43 kg of CO₂ equivalent per case with single-use ureteroscopy and 4.47 kg per case with reusable uteroscopy^[14]. However, another study indicated a significantly higher carbon footprint of a single-use duodenoscope

compared to reusable duodenoscopes, mostly attributable to manufacturing^[16]. Manufacturers are taking steps to improve that impact though transition to renewable energy recycling solutions and goals for carbon neutral manufacturing and net-zero emissions^[17].

Economic Sustainability

Enhanced reprocessing and surveillance practices for duodenoscopes significantly increases reprocessing costs. Doubling HLD procedures increased costs by 47% as reported by a recent US study in 2 high-volume institutions^[18]. In this study, ethylene oxide sterilization increased cost by 270%, and the surveillance practice of culture and quarantine increased cost by 160%. These increases translated into an additional annual cost of \$406,000 for high-volume centers a necessitated a 3.4-fold increase in the number of endoscopes.

There is also the economic impact of HAIs to consider. US data estimates the HAI-related death of 90,000 hospitalized patients yearly with an estimated cost of HAIs of \$28–\$45 billion^[19]. An analysis of ERCP cases in US institutions estimated a cost per procedure of \$612–\$1362 with reusable duodenoscopes including the cost of procedure-related infection at a rate of 0.4–1.0%^[20]. While the upfront cost of disposable endoscopes was higher, depending on the volume of procedure performed at the institution, the price per single-use procedure was reduced to \$112 without the cost of procedure-related HAIs.

Since high-volume endoscopy centers can distribute the cost of investment in equipment, maintenance and treatment of infection, they have a lower break-even price per procedure as opposed to small-volume centers, whose capital investment per case is higher^[20]. For this reason, the use of disposable, single-use endoscopes may be most economically viable in centers with smaller procedure volumes. The insurer Medicare has created new ICD-10-PCS procedure codes and approved applications for New Technology Add-on Payments and Transitional Pass-through Payments for single-use duodenoscopes making 40% of ECRPs in the US eligible for device reimbursement^[21].

Conclusion

The major advantages of single-use endoscopes are reduction or elimination of endoscopy-related infections and elimination of downtime and resources required for reprocessing reusable endoscopes. The impact of endoscopy on the environment and climate change is significant. The environmental impact of single-use endoscopes is equivocal in some cases, but depends on the balance of increased

manufacturing impact versus decreased impact of reprocessing procedures. Improvements in manufacturing impact are ongoing. Accessing the economic sustainability of single-use

endoscopes is complex; however, considering the cost of HAIs, single-use endoscopes may present an economically viable technology, particularly for low-volume or high-risk-case use.

References:

- [1] W. A. Rutala, H. Kanamori, E. E. Sickbert-Bennett, and D. J. Weber, "What's new in reprocessing endoscopes: Are we going to ensure 'the needs of the patient come first' by shifting from disinfection to sterilization?," *American Journal of Infection Control*, vol. 47, pp. A62–A66, Jun. 2019, doi: 10.1016/j.ajic.2019.01.017.
- [2] A. W. Rauwers et al., "High prevalence rate of digestive tract bacteria in duodenoscopes: a nationwide study," *Gut*, vol. 67, no. 9, pp. 1637–1645, Sep. 2018, doi: 10.1136/gutjnl-2017-315082.
- [3] A. Cassini et al., "Attributable deaths and disability-adjusted life-years caused by infections with antibiotic-resistant bacteria in the EU and the European Economic Area in 2015: a population-level modelling analysis," *The Lancet Infectious Diseases*, vol. 19, no. 1, pp. 56–66, Jan. 2019, doi: 10.1016/S1473-3099(18)30605-4.
- [4] C. L. Ofstead, K. M. Hopkins, M. R. Quick, K. B. Brooks, J. E. Eiland, and H. P. Wetzler, "A Systematic Review of Disposable Sheath Use During Flexible Endoscopy," *AORN Journal*, vol. 109, no. 6, pp. 757–771, 2019, doi: 10.1002/aorn.12699.
- [5] M. T. Elghannam* et al., "Single-use endoscopes: A narrative review," *International Journal of Gastrointestinal Intervention*, vol. 11, no. 1, pp. 1–4, Jan. 2022, doi: 10.18528/ijgi210055.
- [6] A. S. Ross et al., "Novel single-use duodenoscope compared with 3 models of reusable duodenoscopes for ERCP: a randomized bench-model comparison," *Gastrointest Endosc*, vol. 91, no. 2, pp. 396–403, Feb. 2020, doi: 10.1016/j.gie.2019.08.032.
- [7] J. Y. Bang, R. Hawes, and S. Varadarajulu, "Equivalent performance of single-use and reusable duodenoscopes in a randomised trial," *Gut*, vol. 70, no. 5, pp. 838–844, May 2021, doi: 10.1136/gutjnl-2020-321836.
- [8] V. R. Muthusamy et al., "Clinical Evaluation of a Single-Use Duodenoscope for Endoscopic Retrograde Cholangiopancreatography," *Clin Gastroenterol Hepatol*, vol. 18, no. 9, pp. 2108–2117.e3, Aug. 2020, doi: 10.1016/j.cgh.2019.10.052.
- [9] A. Slivka et al., "Single-use duodenoscope for ERCP performed by endoscopists with a range of experience in procedures of variable complexity," *Gastrointestinal Endoscopy*, vol. 94, no. 6, pp. 1046–1055, Dec. 2021, doi: 10.1016/j.gie.2021.06.017.
- [10] B. Napoléon et al., "Evaluation of the performances of a single-use duodenoscope: Prospective multi-center national study," *Dig Endosc*, vol. 34, no. 1, pp. 215–221, Jan. 2022, doi: 10.1111/den.13965.
- [11] "Single-Use (Disposable) Flexible Bronchoscopes: The Future of Bronchoscopy? | SpringerLink." <https://link.springer.com/article/10.1007/s12325-020-01495-8> (accessed May 17, 2022).
- [12] "Single-Use Duodenoscopes for ERCP: Rationale, Feasibility, Cost, and Environmental Impact – Gastroenterology & Hepatology." <https://www.gastroenterologyandhepatology.net/archives/may-2022/single-use-duodenoscopes-for-ercp-rationale-feasibility-cost-and-environmental-impact/> (accessed May 20, 2022).
- [13] S. Namburur et al., "Estimating the environmental impact of disposable endoscopic equipment and endoscopes," *Gut*, Nov. 2021, doi: 10.1136/gutjnl-2021-324729.
- [14] N. F. Davis, S. McGrath, M. Quinlan, G. Jack, N. Lawrentschuk, and D. M. Bolton, "Carbon Footprint in Flexible Ureteroscopy: A Comparative Study on the Environmental Impact of Reusable and Single-Use Ureteroscopes," *J Endourol*, vol. 32, no. 3, pp. 214–217, Mar. 2018, doi: 10.1089/end.2018.0001.
- [15] R. Baddeley, L. Aabakken, A. Veitch, and B. Hayee, "Green Endoscopy: Counting the Carbon Cost of Our Practice," *Gastroenterology*, vol. 162, no. 6, pp. 1556–1560, May 2022, doi: 10.1053/j.gastro.2022.01.057.
- [16] L. V. Hernandez, N. N. T. Le, C. Patnode, O. Siddiqui, and O. Jolliet, "ID: 3526786 COMPARING THE IMPACT OF REUSABLE AND SINGLE-USE DUODENOSCOPES USING LIFE CYCLE ASSESSMENT," *Gastrointestinal Endoscopy*, vol. 93, no. 6, p. AB29, Jun. 2021, doi: 10.1016/j.gie.2021.03.123.
- [17] "EXALT Model D Sustainability Solutions," www.bostonscientific.com. <https://www.bostonscientific.com/en-US/medical-specialties/gastroenterology/exalt/sustainability.html> (accessed May 20, 2022).
- [18] S. Bomman et al., "Economic burden of enhanced practices of duodenoscopes reprocessing and surveillance: balancing risk and cost containment," *Endosc Int Open*, vol. 9, no. 9, pp. E1404–E1412, Aug. 2021, doi: 10.1055/a-1515-2591.
- [19] P. W. Stone, "Economic burden of healthcare-associated infections: an American perspective," *Expert Rev Pharmacoecon Outcomes Res*, vol. 9, no. 5, pp. 417–422, Oct. 2009, doi: 10.1586/erp.09.53.
- [20] J. Y. Bang, B. Sutton, R. Hawes, and S. Varadarajulu, "Concept of disposable duodenoscope: at what cost?," *Gut*, vol. 68, no. 11, pp. 1915–1917, Nov. 2019, doi: 10.1136/gutjnl-2019-318227.
- [21] "Economics - EXALT Model D Single-Use Duodenoscope," www.bostonscientific.com. <https://www.bostonscientific.com/en-US/medical-specialties/gastroenterology/exalt/economics.html> (accessed May 20, 2022).

Future Outlook: Emerging Trends and Technologies in Endoscopy

Michael A James, PhD

New Procedures and Training

Continually evolving endoscopy technologies are improving procedural safety and broadening the diagnostic and interventional repertoire of the field, but also bringing increasing complexity and training demands. While training programs previously focused on endoscopic retrograde cholangiopancreatography (ERCP) and endoscopic ultrasound (EUS), education and training now needs to include stenting, ablation procedures, third space endoscopy and endoscopic bariatric therapy (EBT). Regulation of training through accredited programs and assessment of competence are yet to catch up to increased training complexity and demand. This increased complexity and demand have led to a focus on competency-based assessment rather than that which is based on the number of training procedures accomplished^[1,2].

Several advanced procedures that have recently been more broadly implemented contribute to increasing complexity and training demands in endoscopy. Submucosal endoscopy, or third-space endoscopy, involves access to deeper gastrointestinal tract layers in the submucosal space. These include endoscopic submucosal dissection (ESD) and peroral endoscopic myotomy (POEM). ESD has the advantage over mucosal dissection in that en block lesion resection is possible, making large lesions resectable and decreasing the risk of recurrence locally^[3]. With extensive training required and risk of serious complications, ESD expertise and training have not been sufficient in Western countries as opposed to Japan, which has enjoyed specialized training and success in terms of safety and efficacy^[2]. The European Society of Gastrointestinal Endoscopy (ESGE) has accordingly introduced recommendations for training in ESD requiring expert supervision and structured fellowships^[4]. Training requires a priori proficiency in endoscopic classification and lesion characterization and skill in hemostasis, adverse event management, resection and injection. POEM is a technique that provides

a minimally invasive procedure for managing motility disorders of the esophagus with superior outcomes and safety compared to surgery^[5]. These disorders include esophageal achalasia and esophageal motility disorders. A paucity of clinical data and small sample sizes have limited the establishment of recommendations for these procedures^[6]. POEM is another complex procedure requiring considerable skill. Finally, EBT is emerging as an adjunctive therapy for obesity, such as sleeve gastrectomy, and bariatric surgery-related adverse events. With a variety of EBT procedures, complexity and training requirements also vary. Those with greater complexity will require supervision and training programs that are structured for rigor.

Advancements in Wireless Capsules

Wireless capsules have been used for small bowel endoscopy since 2001 and since expanded to colon and esophagus endoscopy. Capsules have been instrumental in managing gastrointestinal bleeding, Crohn's disease, polyposis and celiac disease^[7]. Limitations of these newer applications in the esophagus, stomach and large bowel stem from the increased surface area and volume of these spaces^[8]. Capsule endoscopy (CE) has limited utility for identifying gastric lesions for similar reasons. However, advancements have been made in attempt to overcome these limitations in the form of locomotion systems, air insufflation techniques, improved imaging technologies and AI-assisted image interpretation. Reduced cost and improved technological advancements present the potential for broader application of CE to the gastrointestinal tract.

Despite recent improvements in imaging resolution and field-of-view, identification of disease using CE remains challenging. Imaging algorithms, adoptive frame rate and suspected blood indicators have somewhat improved interpretation; however, these improvements have been incremental with modest or insignificant




Figure 1. Output from CADx identification of pre-malignant components of colorectal polyps. Benign-green, pre-malignant-red. Fonollà, R.; E. W. van der Zander, Q.; Schreuder, R.M.; Masclee, A.A.M.; Schoon, E.J.; van der Sommen, F.; de With, P.H.N. A CNN CADx System for Multimodal Classification of Colorectal Polyps Combining WL, BLI, and LCI Modalities. Appl. Sci. 2020, 10, 5040. <https://doi.org/10.3390/app10155040> <https://creativecommons.org/licenses/by/4.0/>

improvement in diagnostic interpretation and misidentification^[7]. Recently, the use of artificial intelligence for the interpretation of CE imaging has been a subject of interest and technologic development to overcome limitations in human detection and diagnosis. Terms used for AI in CE have included computer-aided detection (CADe) and computer-aided diagnosis (CADx). AI has been investigated for use in polyp identification, optical biopsy and inspection of the colonic mucosa^[9]. Since colon polyps can be difficult to identify, CADe has been applied to their detection. Initial studies of this application demonstrated a detection rate of 97.7% of polyps, although false positives were common (36 per colonoscopy)^[10]. Optical biopsy allows the diagnosis of adenomas among polyps without pathologic review, which is useful given that small polyps are unlikely to exhibit invasive or advanced histology^[9]. Dependence on the expertise of operators has led to the application of CADx to optical biopsy, employing magnifying and non-magnifying narrow band imaging, laser-induced fluorescence spectroscopy and endocytoscopy^[9] (Figure 1). The accuracy of specific computer-aided magnifying narrow band imaging, endocytoscopy and laser-induced fluorescence spectroscopy technologies have been quantified at 93%, 90% and 85%, respectively^[11-13]. These AI technologies are burgeoning and promise to improve the diagnostic yield of CE.

Imaging Technologies

In addition to AI technologies, which can be applied to both CE and traditional endoscopy, several advancements in imaging have been made to improve detection and diagnosis. Hypoxia imaging using oxygen-sensitive fluorescent stains have been applied to the detection of neoplasias in the pharynx, esophagus and colon and of inflammatory bowel disease^[14]. In cases of the above malignancies, hypoxia imaging was able to detect significant changes in oxygen saturation mapping of neoplastic versus non neoplastic tissues^[15].

In inflammatory bowel

disease, hypoxia of the intestinal epithelium is increased, making its detection using hypoxia imaging plausible^[14]. This application requires further research. Three-dimensional imaging, which combines optical axes to create depth is an older technology that has expanded in application to endoscopy. A study of detection of non-polypoid lesion indicated a 25% increase in trainee detection in simulations over detection using two-dimensional images^[16].

The most prominent technological imaging advancement in recent years has been AI-assisted visual sensing, as described in the context of CE above. This technology also applies to other endoscopic techniques, including duodenoscopy for gallbladder stones^[17]. In a recent clinical study, AI-assisted duodenoscopic imaging of gallbladder stones using a convolutional neural network resulted in 95% precision and less procedural, recovery and hospital time^[17]. Further, bleeding and residual stones were greatly decreased.

ERCP/Duodenoscopy Advancements

The long established role for ERCP in diagnostic procedures has expanded to interventions over time. A common example is biliary obstruction management. ERCP can be used to place stents for strictures and radiofrequency ablation for hilar structures. Improvements in duodenoscopy,

Recently, the use of artificial intelligence for the interpretation of CE imaging has been a subject of interest and technologic development to overcome limitations in human detection and diagnosis

With safety as a priority, duodenoscope technologies are also now addressing the significant concern of healthcare acquired infection. Such technological advancements are certain to result in more well-equipped clinicians and increased wellness for patients

in addition to the applicable imaging and AI integrations, as discussed above, have included robotic manipulation. For example, a flexible mechanical robotic system for coaxial catheter intervention using cholangiopancreatography was recently described with demonstrated feasibility in models^[18]. Given better outcomes in modern duodenoscopy, there has also been a trend in outpatient ERCP, previously predominantly done in an inpatient setting^[19].

The main focus of advancements in duodenoscopy and ERCP in recent years has been in infection prevention. The resurgence of multidrug-resistant organism infections related to duodenoscopy has largely driven this focus. While such infections are not considered common, they are of considerable public health concern, and rates contamination with MDROs on endoscopes are surprisingly high^[20]. Solutions to this problem that have been put forth include more rigorous reprocessing training and assessment, further research into ERCP-associated infection risk and epidemiology, mandatory inspection and servicing of duodenoscopes and novel technologies, including sterilizable duodenoscopes, disposable scope components and single-use disposable duodenoscopes^[20]. Single-use scopes have been demonstrated to be tantamount to reusable scopes without the risk of transmission of infectious disease between patients^[21].

The clinical utility of endoscopy has greatly broadened in recent years while simultaneously becoming safer. Better technology, techniques, imaging and computer-assisted interpretation have resulted in shorter, safer procedures that can more commonly be done in outpatient setting

Conclusion

The clinical utility of endoscopy has greatly broadened in recent years while simultaneously becoming safer. Better technology, techniques, imaging and computer-assisted interpretation have resulted in shorter, safer procedures that can more commonly be done in outpatient setting. These advancements have also improved both disease-specific and adverse-event related outcomes. With the use of algorithms to more reliably identify disease and robotics for more precise manipulation, there is promise for continued improvement in outcomes. With increasing procedure complexity, the risk of infection may increase, making the prevention of infection particularly timely with the development of advanced techniques. With safety as a priority, duodenoscope technologies are also now addressing the significant concern of healthcare acquired infection. Such technological advancements are certain to result in more well-equipped clinicians and increased wellness for patients.

References:

- [1] T. J. Nasca, I. Philibert, T. Brigham, and T. C. Flynn, "The Next GME Accreditation System — Rationale and Benefits," *New England Journal of Medicine*, vol. 366, no. 11, pp. 1051–1056, Mar. 2012, doi: 10.1056/NEJMs1200117.
- [2] D. Yang, M. S. Wagh, and P. V. Draganov, "The status of training in new technologies in advanced endoscopy: from defining competence to credentialing and privileging," *Gastrointestinal Endoscopy*, vol. 92, no. 5, pp. 1016–1025, Nov. 2020, doi: 10.1016/j.gie.2020.05.047.
- [3] D. Yang, M. Othman, and P. V. Draganov, "Endoscopic Mucosal Resection vs Endoscopic Submucosal Dissection For Barrett's Esophagus and Colorectal Neoplasia," *Clinical Gastroenterology and Hepatology*, vol. 17, no. 6, pp. 1019–1028, May 2019, doi: 10.1016/j.cgh.2018.09.030.
- [4] P. Pimentel-Nunes et al., "Curriculum for endoscopic submucosal dissection training in Europe: European Society of Gastrointestinal Endoscopy (ESGE) Position Statement," *Endoscopy*, vol. 51, no. 10, pp. 980–992, Oct. 2019, doi: 10.1055/a-0996-0912.
- [5] E. Akintoye, N. Kumar, I. Obaitan, Q. A. Alayo, and C. C. Thompson, "Peroral endoscopic myotomy: a meta-analysis," *Endoscopy*, vol. 48, no. 12, pp. 1059–1068, Dec. 2016, doi: 10.1055/s-0042-114426.
- [6] N. Di Lorenzo et al., "Clinical practice guidelines of the European Association for Endoscopic Surgery (EAES) on bariatric surgery: update 2020 endorsed by IFSO-EC, EASO and ESPCOP," *Surg Endosc*, vol. 34, no. 6, pp. 2332–2358, Jun. 2020, doi: 10.1007/s00464-020-07555-y.
- [7] M. F. Byrne and F. Donnellan, "Artificial intelligence and capsule endoscopy: Is the truly 'smart' capsule nearly here?," *Gastrointestinal Endoscopy*, vol. 89, no. 1, pp. 195–197, Jan. 2019, doi: 10.1016/j.gie.2018.08.017.
- [8] S.-J. Nam, H. S. Lee, and Y. J. Lim, "Evaluation of Gastric Disease with Capsule Endoscopy," *Clin Endosc*, vol. 51, no. 4, pp. 323–328, Jul. 2018, doi: 10.5946/ce.2018.092.
- [9] M. F. Byrne, N. Shahidi, and D. K. Rex, "Will Computer-Aided Detection and Diagnosis Revolutionize Colonoscopy?," *Gastroenterology*, vol. 153, no. 6, pp. 1460–1464.e1, Dec. 2017, doi: 10.1053/j.gastro.2017.10.026.
- [10] Y. Wang, W. Tavanapong, J. Wong, J. H. Oh, and P. C. de Groen, "Polyp-Alert: near real-time feedback during colonoscopy," *Comput Methods Programs Biomed*, vol. 120, no. 3, pp. 164–179, Jul. 2015, doi: 10.1016/j.cmpb.2015.04.002.
- [11] Y. Kominami et al., "Computer-aided diagnosis of colorectal polyp histology by using a real-time image recognition system and narrow-band imaging magnifying colonoscopy," *Gastrointestinal Endoscopy*, vol. 83, no. 3, pp. 643–649, Mar. 2016, doi: 10.1016/j.gie.2015.08.004.
- [12] M. Misawa et al., "Characterization of Colorectal Lesions Using a Computer-Aided Diagnostic System for Narrow-Band Imaging Endocytoscopy," *Gastroenterology*, vol. 150, no. 7, pp. 1531–1532.e3, Jun. 2016, doi: 10.1053/j.gastro.2016.04.004.
- [13] T. Rath, G. E. Tontini, M. Vieth, A. Nagel, M. F. Neurath, and H. Neumann, "In vivo real-time assessment of colorectal polyp histology using an optical biopsy forceps system based on laser-induced fluorescence spectroscopy," *Endoscopy*, vol. 48, no. 6, pp. 557–562, Jun. 2016, doi: 10.1055/s-0042-102251.
- [14] S. Gulati, M. Patel, A. Emmanuel, A. Haji, B. Hayee, and H. Neumann, "The future of endoscopy: Advances in endoscopic image innovations," *Digestive Endoscopy*, vol. 32, no. 4, pp. 512–522, 2020, doi: 10.1111/den.13481.
- [15] K. Kaneko et al., "Hypoxia imaging endoscopy equipped with laser light source from preclinical live animal study to first-in-human subject research," *PLoS One*, vol. 9, no. 6, p. e99055, 2014, doi: 10.1371/journal.pone.0099055.
- [16] S. Sakata, P. M. Grove, A. R. L. Stevenson, and D. G. Hewett, "The impact of three-dimensional imaging on polyp detection during colonoscopy: a proof of concept study," *Gut*, vol. 65, no. 5, pp. 730–731, May 2016, doi: 10.1136/gutjnl-2016-311507.
- [17] D. Li, B. Du, Y. Shen, and L. Ge, "Artificial Intelligence-Assisted Visual Sensing Technology under Duodenoscopy of Gallbladder Stones," *Journal of Sensors*, vol. 2021, p. e5158577, Oct. 2021, doi: 10.1155/2021/5158577.
- [18] W. Wu et al., "Design of coaxial intervention catheter robot for Endoscopic Retrograde Cholangio-Pancreatography," in 2020 Chinese Automation Congress (CAC), Nov. 2020, pp. 6712–6716. doi: 10.1109/CAC51589.2020.9327900.
- [19] A. Almuhaideb, D. Olson, and A. A. Adam, "Advancements in Endoscopic Biliary Interventions by Gastroenterology," *Semin Intervent Radiol*, vol. 38, no. 3, pp. 280–290, Aug. 2021, doi: 10.1055/s-0041-1731266.
- [20] C. L. Ofstead, B. L. Buro, K. M. Hopkins, J. E. Eiland, H. P. Wetzler, and D. R. Lichtenstein, "Duodenoscope-associated infection prevention: A call for evidence-based decision making," *Endosc Int Open*, vol. 08, no. 12, pp. E1769–E1781, Dec. 2020, doi: 10.1055/a-1264-7173.
- [21] M. T. Elghannam* et al., "Single-use endoscopes: A narrative review," *International Journal of Gastrointestinal Intervention*, vol. 11, no. 1, pp. 1–4, Jan. 2022, doi: 10.18528/ijgi210055.

Notes:

Notes:

Sponsored by

**Boston
Scientific**

Advancing science for life™