


**FARAPULSE™**  
Pulsed Field Ablation System



**The global leader in  
PFA clinical research**





# FARAPULSE™

## Pulsed Field Ablation System

## Clinical leadership

FARAPULSE™ Pulsed Field Ablation System has been extensively researched and is supported by more published clinical evidence than any other PFA system.



The  
**MOST ROBUST**  
clinical trial  
portfolio

## SELECTED STUDIES

### COMPLETED

#### IMPULSE/PEFCAT I - II

Safety and Feasibility Study

#### PersAFOne I

Safety and Feasibility Study in Persistent AF

#### ADVENT

Randomized Trial vs. RFA and CBA

#### FARA-FREEDOM

Single-Arm Post-Market Trial (EU)

#### EU-PORIA

Multicenter Registry (n = 1,233)

#### MANIFEST-PF

Multicenter Registry (n = 1,568)

#### MANIFEST-17K

Multicenter Registry (n = 17,642)

### INITIATED/PLANNED

#### ADVANTAGE I - II

PVI and PWA in Persistent AF with FARAPOINT™ for CTI

#### PersAFOne II - III

Includes FARAPOINT focal PFA for CTI

#### FARADISE

Prospective Registry (n = 1,000+)

#### AVANT GUARD

First-line PFA vs. AAD in Persistent AF with ILR

#### DISRUPT-AF

Registry and Randomized Workflow Study

#### NAVIGATE PF

FARAVIEW™ Nav-enabled Mapping Integration

#### OPTION-A

Concomitant WATCHMAN™ and FARAPULSE

#### REMATCH Trial

Redo AF with FARAWAVE™ and FARAPOINT



**120+ clinical publications and counting.**

Scan the code to view all the latest FARAPULSE research and findings.

# Demonstrated safety

FARAPULSE enables effective ablation while eliminating risks associated with thermal procedures.

## No reported thermal complications<sup>1</sup>


In MANIFEST-17K, the largest PFA registry study conducted to date (106 centers, 413 operators, 17,642 pts), patients treated with FARAPULSE reported zero thermal complications across all procedures.



No esophageal  
fistula



No pulmonary  
vein stenosis



No persistent  
phrenic nerve injury

# 0.98%

major adverse event rate


## Major adverse event rate <1%<sup>1</sup>

In the same registry, FARAPULSE demonstrated a low adverse event rate—less than one percent of all patients (173 of 17,642 pts).

Della Rocca et al., 2023

## Lower minor adverse event rate vs. thermal<sup>2</sup>

In a propensity score-matched study between PFA, cryoballoon ablation (CBA) and radiofrequency ablation (RFA), FARAPULSE had a significantly lower minor complication rate.



Authors attributed the performance in large part to the ability of PFA to

**SELECTIVELY TARGET MYOCARDIAL CELLS**

while sparing surrounding tissue


# Proven PVI durability

FARAPULSE™ Pulsed Field Ablation System offers a proven dosing strategy that delivers lasting pulmonary vein (PV) isolation.

Della Rocca et al., 2023

## Higher PV durability vs. thermal<sup>2</sup>


Propensity-matched analysis demonstrated that the rate of PV reconnection at repeat procedures was significantly lower among patients receiving PFA at index procedure.



## Multiple studies

## Freedom from recurrence vs. thermal

Across multiple studies, FARAPULSE patients have trended towards a lower rate of arrhythmia recurrence at 12 months compared to thermal ablation.





# Significantly reduced atrial arrhythmia (AA) burden vs. thermal<sup>9,10</sup>




### Quality of Life

Significantly greater QoL improvement in patients with AA Burden <0.1% vs. ≥10%



### Healthcare Utilization

Significantly lower risk for redo ablation, cardioversion and hospitalization with AA Burden <0.1% vs. ≥0.1%



### Ablation Modality


FARAPULSE patients significantly more likely to have AA Burden <0.1% vs. RFA or CBA

## Multiple studies


# High real-world freedom from recurrence (1 year)

FARAPULSE has demonstrated a high one-year freedom from recurrence rate, consistent across numerous clinical settings and trial sizes.

## PAROXYSMAL AF



## PERSISTENT AF



# Predictable procedures

FARAPULSE™ Pulsed Field Ablation System delivers a consistent, predictable operator experience—for any workflow.

Reddy et al., 2023

## Optimized for consistency

FARAPULSE PFA System excels at achieving optimal placement, regardless of pulmonary vein anatomy. This maneuverability translates to limited LA dwell time, with FARAPULSE requiring 42% less time to perform PVI compared with thermal.<sup>9</sup>


|                | PFA             | Thermal         |
|----------------|-----------------|-----------------|
| Ablation time  | 29.2<br>± 14.3  | 50.0<br>± 24.6  |
| LA dwell time  | 59.4<br>± 18.3  | 83.7<br>± 30.3  |
| Procedure time | 105.8<br>± 29.4 | 123.1<br>± 42.1 |
| Fluoro. time   | 21.1<br>± 11.0  | 13.9<br>± 12.8  |

Procedure time and LA dwell time include the mandated 20-minute wait period, per study protocol

Multiple studies

## Real-world procedure times (large registries)


Across multiple large-registries (>1,000 pts), FARAPULSE's procedure times remained consistent—allowing for more predictable, efficient workflows.





1. Ekanem, E., Neuzil, P., Reichlin, T. et al. Safety of pulsed field ablation in more than 17,000 patients with atrial fibrillation in the MANIFEST-17K study. *Nat Med* (2024). doi:10.1038/s41591-024-03114-3.
2. Della Rocca DG, Marcon L, Magnocavallo M, et al. Pulsed electric field, cryoballoon, and radiofrequency for paroxysmal atrial fibrillation ablation: A propensity score-matched comparison. *Europace*. 2023;26(1):euac016.
3. Badertscher P, Weidlich S, Knecht S, et al. Efficacy and safety of pulmonary vein isolation with pulsed field ablation vs. novel cryoballoon ablation system for atrial fibrillation. *Europace* 25.12 (2023): euad329.
4. Maurhofer, Jens, et al. Pulsed-field vs. cryoballoon vs. radiofrequency ablation: A propensity score matched comparison of one-year outcomes after pulmonary vein isolation in patients with paroxysmal atrial fibrillation. *J Interv Card Electrophysiol*. 67.2 (2024): 389-397.
5. Ekanem E, Reddy VY, Schmidt B, et al., Multi-national survey on the methods, efficacy, and safety on the post-approval clinical use of pulsed field ablation (MANIFEST-PF) [published correction appears in *Europace*. 2023 Feb 16;25(2):449. doi:10.1093/europace/euac250]. *Europace*. 2022;24(8):1256-1266. doi:10.1093/europace/euac050.
6. Urbanek L, Bordignon S, Schaack D, et al. Pulsed field versus cryoballoon pulmonary vein isolation for atrial fibrillation: Efficacy, safety, and long-term follow-up in a 400-patient cohort. *Circulation: Arrhythmia and Electrophysiology* 16.7 (2023): 389-398.
7. Chaumont C, McDonnell E, Boveda S, et al. Prospective 1-year results of atrial fibrillation ablation using the pentaspline pulsed field ablation catheter: The initial French experience. *Arch Cardiovasc Dis*. Mar 5 2024;doi:10.1016/j.acvd.2024.01.005.
8. Rattka M, Mavrikis E, Vlachopoulou D, et al. Pulsed field ablation and cryoballoon ablation for pulmonary vein isolation: Insights on efficacy, safety and cardiac function. *J Interv Card Electrophysiol*. (2024): 1-8.
9. Reddy VY, Gerstenfeld EP, Natale A, et al. Pulsed field or conventional thermal ablation for paroxysmal atrial fibrillation. *New England Journal of Medicine*. 2023; 389(18):1660-1671. doi:10.1056/NEJMoa2307291.
10. Reddy VY, Mansour M, Calkins H, et al. Pulsed Field vs Conventional Thermal Ablation for Paroxysmal Atrial Fibrillation: Recurrent Atrial Arrhythmia Burden. *J Am Coll Cardiol*. Published online May 6, 2024. doi:10.1016/j.jacc.2024.05.001.
11. Schmidt B, Bordignon S, Nevenet K, et al. European real-world outcomes with pulsed field ablation in patients with symptomatic atrial fibrillation: Lessons from the multi-centre EU-PORIA Registry. *Europace*. 2023;25(7):euad185.
12. Boresma, et al., Real World Data Collection in Subjects Treated with the FARAPULSE Pulsed Field Ablation System (FARADISE). HRS 2024.

CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings, and instructions for use can be found in the product labelling supplied with each device or at [www.IFU-BSCI.com](http://www.IFU-BSCI.com). Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. This material not intended for use in France. ©2024 Boston Scientific Corporation or its affiliates. All rights reserved



**Boston  
Scientific**  
Advancing science for life™

[www.bostonscientific.eu](http://www.bostonscientific.eu)

© 2024 Boston Scientific Corporation or its affiliates. All rights reserved.

EP-1813001-AA

CE 0459