

IMAGEREADY™ MRI CHECKLIST FOR VERCISE™ DIRECTIONAL* DBS SYSTEMS

Refer to ImageReady MRI Guidelines for Boston Scientific Deep Brain Stimulation (DBS) Systems on www.bostonscientific.com/ manuals for labeling and safety conditions.

Patient Name: _____ Date: _____

Physician's Name: _____

Office Address: _____

Phone: _____

A. TYPE OF MR CONDITIONAL DBS SYSTEM			
<input type="checkbox"/> Lead-Only System (Leads, Lead Boots and/or Burr Hole Covers)			
<input type="checkbox"/> Full System** (Leads, Lead Extensions, Stimulator and/or Burr Hole Covers)			
B. MR CONDITIONAL SYSTEM COMPONENTS	MODEL NUMBER	MRI ELIGIBLE	NOT MRI ELIGIBLE
VERCISE GENUS™ DBS SYSTEMS – LEADS			
Vercise Cartesia™ 30 cm 8-Contact DBS Directional Lead	DB-2202-30	<input type="checkbox"/>	<input type="checkbox"/>
Vercise Cartesia 45 cm 8-Contact DBS Directional Lead	DB-2202-45	<input type="checkbox"/>	<input type="checkbox"/>
Vercise 30 cm 8-Contact DBS Lead	DB-2201-30-DC	<input type="checkbox"/>	<input type="checkbox"/>
Vercise 45 cm 8-Contact DBS Lead	DB-2201-45-DC	<input type="checkbox"/>	<input type="checkbox"/>
VERCISE GENUS DBS SYSTEM – EXTENSIONS			
55 cm 1x8 Contact Lead Extension	NM-3138-55	<input type="checkbox"/>	<input type="checkbox"/>
VERCISE GENUS DBS SYSTEM – STIMULATORS			
Vercise Genus R16 Implantable Pulse Generator	DB-1216	<input type="checkbox"/>	<input type="checkbox"/>
Vercise Genus P16 Implantable Pulse Generator	DB-1416	<input type="checkbox"/>	<input type="checkbox"/>
Vercise Genus P8 Implantable Pulse Generator	DB-1408	<input type="checkbox"/>	<input type="checkbox"/>
VERCISE GENUS DBS SYSTEM – FIXATION AND ACCESSORIES			
SureTek™ Burr Hole Cover	Provided in SureTek Burr Hole Cover Kit (DB-4600-C) and SureTek Burr Hole Cover Spares Kit (DB-4605-C)	<input type="checkbox"/>	<input type="checkbox"/>
Lead Boot	Provided in the Vercise Physician Spares Kit DB-2500-C and with leads (see above)	<input type="checkbox"/>	<input type="checkbox"/>
Silicone Suture Sleeves	Provided in the Vercise Physician Spares Kit DB-2500-C and with leads (see above)	<input type="checkbox"/>	<input type="checkbox"/>
Port Plug	Provided in the Port Plug Spares Kit, SC-4401, IPG Kits, and Lead Extension Kits, DB-3128	<input type="checkbox"/>	<input type="checkbox"/>

*A System that includes the Vercise™ PC, Vercise Gevia™, or Vercise Genus™ IPG and Vercise Cartesia™ Directional Lead(s) forms the Vercise Directional System.

**MRI Conditional under specific conditions.

VERCISE GEVIA™ DBS SYSTEM – LEADS			
Vercise Cartesia™ 30 cm 8-Contact DBS Directional Lead	DB-2202-30	<input type="checkbox"/>	<input type="checkbox"/>
Vercise Cartesia 45 cm 8-Contact DBS Directional Lead	DB-2202-45	<input type="checkbox"/>	<input type="checkbox"/>
30 cm 8-Contact DBS Lead	DB-2201-30-DC	<input type="checkbox"/>	<input type="checkbox"/>
45 cm 8-Contact DBS Lead	DB-2201-45-DC	<input type="checkbox"/>	<input type="checkbox"/>
VERCISE GEVIA DBS SYSTEM – EXTENSION			
55 cm 1x8 Contact Lead Extension	NM-3138-55	<input type="checkbox"/>	<input type="checkbox"/>
VERCISE GEVIA DBS SYSTEM – STIMULATOR			
Vercise Gevia Implantable Pulse Generator	DB-1200	<input type="checkbox"/>	<input type="checkbox"/>
VERCISE GEVIA DBS SYSTEM – FIXATION AND ACCESSORIES			
SureTek™ Burr Hole Cover	Provided in SureTek Burr Hole Cover Kit (DB-4600-C) and SureTek Burr Hole Cover Spares Kit (DB-4605-C)	<input type="checkbox"/>	<input type="checkbox"/>
Lead Boot	Provided in the Vercise Physician Spares Kit DB-2500-C and with leads (see above)	<input type="checkbox"/>	<input type="checkbox"/>
Silicone Suture Sleeves	Provided in the Vercise Physician Spares Kit DB-2500-C and with leads (see above)	<input type="checkbox"/>	<input type="checkbox"/>
Port Plug	Provided in the Port Plug Spares Kit, SC-4401, IPG Kit, DB-1200	<input type="checkbox"/>	<input type="checkbox"/>
C. DBS IMPLANT CONFIGURATION AND SYSTEM INTEGRITY (Check all that apply for Lead-Only System or Full System)			
LEAD-ONLY SYSTEM		EXTERNALIZED LEAD-ONLY SYSTEM	
Fully Implanted Lead-Only System		Externalized Lead-Only System	
MRI ELIGIBLE		MRI ELIGIBLE	
<input type="checkbox"/> Stimulator NOT implanted.		<input type="checkbox"/> Stimulator NOT implanted.	
<input type="checkbox"/> Lead extensions NOT implanted.		<input type="checkbox"/> Lead extensions NOT implanted.	
<input type="checkbox"/> Leads capped with lead boot.		<input type="checkbox"/> Partially implanted lead(s) extending out of the patient are straight with no loops and stylet has been removed from externalized lead.	
<input type="checkbox"/> Lead(s) fully implanted under the scalp of the skull.		<input type="checkbox"/> The external portion of the partially implanted lead(s) is NOT in contact with either the patient or any part of the scanner.	
<input type="checkbox"/> Patient has up to two leads implanted.		<input type="checkbox"/> Patient has up to two leads implanted.	
<input type="checkbox"/> NO evidence of fractured leads.		<input type="checkbox"/> NO evidence of fractured leads.	

Note: All of the conditions above must be met for the patient to be MRI eligible.

FULL SYSTEM**	
MRI ELIGIBLE	
<input type="checkbox"/> The stimulator must be implanted under the skin in a location near the clavicle (pectoral region).	<input type="checkbox"/> Unused stimulator ports have a port plug inserted.
<input type="checkbox"/> The stimulator must be implanted on the same side of the body as the implanted lead extension(s).	<input type="checkbox"/> For a bilateral implant where two leads and lead extensions are connected to a single stimulator, both lead extensions are routed on the same side of the body as the stimulator.
<input type="checkbox"/> Patient has up to two leads implanted.	<input type="checkbox"/> No evidence of fractured leads or compromised stimulator-lead system integrity. as of: _____/_____/_____
<input type="checkbox"/> Extensions directly connected to stimulator. No adapters present.	<input type="checkbox"/> Single stimulator is implanted.
<input type="checkbox"/> For Genus Systems: Any system with a single lead and unused lead extension port(s) have a port plug inserted into the unused lead extension port(s). <input type="checkbox"/> For Gevia Systems: The system does NOT have any unused lead extension ports (all lead extensions present have a lead inserted).	

Note: All of the conditions above must be met for the patient to be MRI eligible.

D. INSTRUCTIONS FOR THE PATIENT OR MRI CENTER PRIOR TO THE MRI SCAN (Full System only)

1. Rechargeable stimulators must be fully charged (stimulator battery level on the remote control must be at three bars) before the MRI scan. Patient must bring their charger to the MRI center in case additional charging is necessary.
2. MRI mode must be enabled on the stimulator using the patient remote control before performing an MRI scan. Patient must bring their remote control to the MRI center.



This information is not meant as a replacement for the ImageReady™ MRI Guidelines for Boston Scientific Deep Brain Stimulation Systems.



The Vercise Genus™ DBS System, Vercise Gevia™ DBS System, and Vercise™ DBS Lead-only system (before Stimulator is implanted) provide safe access to full-body MRI scans when used with specific components and the patient is exposed to the MRI environment under specific conditions defined in the supplemental manual ImageReady™ MRI Guidelines for Boston Scientific DBS Systems.

Indications for Use: The Boston Scientific Deep Brain Stimulation Systems are indicated for use in:

- Bilateral stimulation of the subthalamic nucleus (STN) as an adjunctive therapy in reducing some of the symptoms of moderate to advanced levodopa responsive Parkinson's disease (PD) that are not adequately controlled with medication.
- Bilateral stimulation of the internal globus pallidus (GPI) as an adjunctive therapy in reducing some of the symptoms of advanced levodopa responsive Parkinson's disease (PD) that are not adequately controlled with medication.

Contraindications, warnings, precautions, side effects: The Deep Brain Stimulation Systems or any of its components, is contraindicated for: Diathermy as either a treatment for a medical condition or as part of a surgical procedure, Electroconvulsive Therapy (ECT) and Transcranial Magnetic Stimulation (TMS) as the safety of these therapies in patients implanted with the Vercise™ DBS System has not been established, patients who are unable to operate the system, patients who are poor surgical candidates or who experience unsuccessful test stimulation. Patients implanted with Boston Scientific Deep Brain Stimulation Systems without ImageReady™ MRI Technology should not be exposed to Magnetic Resonance Imaging (MRI). Patients implanted with Vercise Gevia™ or Vercise Genus™ or Vercise DBS Lead-only system (before Stimulator is implanted) with ImageReady MRI Technology are Full Body MR Conditional only when exposed to the MRI environment under the specific conditions defined in ImageReady MRI Guidelines for Boston Scientific Deep Brain Stimulation Systems. Assess patients for the risks of depression and suicide. This assessment should consider both the risk of depression and suicide as well as the potential clinical benefits of DBS therapy. Monitor patients for new or worsening symptoms of depression, suicidal thoughts or behaviors, or changes in mood or impulse control and manage appropriately. Refer to the Instructions for Use provided with the Vercise DBS System or BostonScientific.com for potential adverse effects, warnings, and precautions prior to using this product.

Caution: U.S. Federal law restricts this device to sale by or on the order of a physician.

**Boston
Scientific**
Advancing science for life™

25155 Rye Canyon Loop
Valencia, CA 91355 USA

Copyright ©2021 by
Boston Scientific Corporation
or its affiliates. All rights reserved.

NM-566609-AD