



TheraSphere™ Y-90 Glass Microspheres | HEPATOLOGY

Liver Transplantation Following Y-90 for HCC 15-year, 207 patient experience evaluated TheraSphere™ as treatment for Bridging or Downstaging HCC to Liver Transplant therapy with a median overall survival* of 12.5 years

Liver Transplantation Following Yttrium-90 Radioembolization: 15-year Experience in 207-Patient Cohort - Ahmed Gabr, Laura Kulik, Samdeep Mouli, Ahsun Riaz, Rehan Ali, Kush Desai, Ronald A Mora, Daniel Ganger, Haripriya Maddur, Steven Flamm, Justin Boike, Christopher Moore, Bartley Thornburg, Ali Alasadi, Talia Baker, Daniel Borja-Cacho, Nitin Katariya, Daniela P Ladner, Juan Carlos Caicedo, Robert J Lewandowski, Riad Salem


METHODS

A multidisciplinary team comprised of hepatology, oncology, transplant surgery, and interventional radiology retrospectively reviewed data from 207 patients with unresectable HCC who underwent liver transplant after being treated with Y-90 as part of a bridging or downstaging care pathway.


RESULTS

Of the 207 patients included in these analyses, 38 (19%) patients were downstaged to within Milan transplant criteria and 169 (82%) bridged to transplant with TheraSphere either using lobar (18%, median dose of 124 Gy [IQR: 132-146]) or radiation segmentectomy (82%, median dose 260 Gy [IQR: 235-350]) administration.

169 Patients were Bridged and 38 Patients were Downstaged to T2 for Liver Transplant
Median Time to LT was 7.5 Months



Median Overall Survival after Liver Transplant was 12.5 years



* Post transplant


HEPATOLOGY: Liver Transplantation Following Y-90

RESULTS & CONCLUSIONS

Three, five and ten-year Overall Survival (OS) rates were 84%, 77%, and 60% respectively. 94 (45%), 60 (29%) and 53 (26%) of patients showed complete (no viable tumor), extensive (50-99% necrosis) and partial tumor necrosis (<50% necrosis) on histopathology of which 2%, 7% and 34% had recurrence in each group, respectively. There was a trend towards better OS for patients achieving complete/extensive tumor necrosis (p=0.056). Median recurrence free survival (RFS) post transplant was 120 Months (95%CI: 69-150). There was no differences in OS or RFS for bridged versus downstaged patients.

74% of Patients Showed Complete or Extensive Tumor Necrosis on Histopathology with an average 10 Year Recurrence-Free Survival following Liver Transplant

Histopathology and Recurrence Rate Post TheraSphere



SURVIVAL AND RECURRENCE OUTCOMES (N = 207)

	Median	3-year	5-year	10-year
Overall Survival from Y-90	157 mo. (13.1 years) [CI: 120-157]	87%	80%	62%
Overall Survival from LT	150 mo. (12.5 years) [CI: 120-150]	84%	77%	60%
Recurrence-Free Survival from LT	120 mo. (10.0 years) [CI: 69-150]	77%	65%	43%
Disease-Specific Mortality Rate	Not Reached	6%	11%	16%
Time-to-Recurrence (Recurrence-Free Probability)	Not Reached	88%	79%	76%
Overall Survival <65 years	Not Reached at 150 mo.	88%	85%	71%
Overall Survival ≥65 years P=0.003	12.5 years	73%	58%	43%

Abbreviations:

HCC: hepatocellular carcinoma

LT: Liver transplantation

Y90: Yttrium-90 radioembolization

MELD: Model of endstage liver disease

BCLC: Barcelona Clinic Liver Cancer

cTACE: conventional chemoembolization

LRT: locoregional therapy

OS: Overall survival

RFS: Recurrence-free survival

TTP: time-to-progression

TTR: Time-to-recurrence

DSM: Disease-specific-mortality

MRI: gadolinium-enhanced magnetic resonance imaging

CT: triphasic contrast-enhanced computerized tomography

CP: Child-Pugh

IQR: Interquartile range

KM: Kaplan-Meier analysis

CI: 95% Confidence Interval

ECOG: Eastern Cooperative Oncology Group

UNOS: United Network for Organ Sharing

AFP: Alpha fetoprotein

ETOH: Alcoholic cirrhosis

NASH: Non-alcoholic steatohepatitis

PBC: Primary biliary cirrhosis

PSC: Primary sclerosing cholangitis

HCV: Hepatitis C virus infection

HBV: Hepatitis B virus Infection

TNM: Tumor, Node, Metastasis

HEPATOLOGY: Liver Transplantation Following Y-90

Baseline Characteristics at Y-90

	Median [IQR]	N (%)
Age (years)	60 [56-65]	
	0	145 (70%)
ECOG	1	61 (29.5%)
	2	1 (0.5%)
Child-Pugh	A	99 (48%)
	B	91 (44%)
	C	17 (8%)
BCLC	A	106 (51%)
	B	20 (10%)
	C	64 (31%)
	D	17 (8%)
UNOS TNM	T1	9 (4%)
	T2	160 (77%)
	T3	22 (11%)
	T4a	12 (6%)
	T4b	4 (2%)
AFP (ng/dL)	<13 (normal)	114 (55%)
	13-100	48 (23%)
	>100	45 (22%)
	Range	0.8-15735
Prior Liver therapy	Surgical Resection	8 (3.5%)
	Prior HCC LRT	35 (17%)
	Treatment Naive	164 (79.5%)
Listing	Prior to Y-90	117 (57%)
	After Y-90	90 (43%)
Y-90 Administration	Lobar	37 (18%)
	Segmental	170 (82%)
Y-90 Dose (Gy)	Lobar	124 [132-146]
	Segmental	260 [235-350]

Baseline Characteristics at Liver Transplant

	Median [IQR]	N (%)
Age (years)	62 [57-66]	
MELD-Na Score	13 [10-17]	
Wait-list time (months)	7 [4-10]	
Time from Y-90 (months)	7.5 [4.4-10.3]	
Etiology of HC	Autoimmune hepatitis	3 (1.5%)
	Alpha 1 antitrypsin	1 (0.5%)
	Biliary Atresia	1 (0.5%)
	Cryptogenic	13 (6%)
	ETOH	30 (14%)
	HCV + ETOH	11 (5%)
	HV	102 (49%)
	HBV	22 (10%)
	NASH	13 (6%)
	PBC	7 (3%)
	Wilson's	1 (0.5%)
	PSC	1 (0.5%)
	Hemochromatosis	2 (1%)
AFP (ng/dL)	<13 (normal)	132 (64%)
	13-100	62 (30%)
	>100	13 (6%)
	Range	0.8-13774
Liver Parenchyma	Cirrhosis	202 (97.5%)
	Bridging Fibrosis	5 (2.5%)
	Grade 1	37 (18%)
	Grade 2	69 (33%)
	Grade 3	6 (3%)
Tumor Grade	Fibromellar	1 (0.5%)
	Mixed HCC - cholangiocarcinoma	4 (2%)
	Unable to identify due to extensive necrosis	90 (43.5%)
	Complete (100%)	94 (45%)
Tumor Necrosis	Extensive (51-99%)	60 (29%)
	Partial (<50%)	53 (26%)

Gabr, A., Kulik, L., Mouli, S., Riaz, A., Ali, R., Desai, K., Mora, R.A., Ganger, D., Maddur, H., Flamm, S., Boike, J., Moore, C., Thornburg, B., Alasadi, A., Baker, T., Borja-Cacho, D., Kariyana, N., Ladner, D.P., Caicedo, J.C., Lewandowski, J.P., and Salem, R. (2020). Liver Transplantation Following Yttrium-90 Radioembolization: 15-year Experience in 207-Patient Cohort. *Hepatology*. Accepted Author Manuscript. doi:10.1002/hep.31318

Boston Scientific is not responsible for the collection, analysis or reporting of the investigator-sponsored research output which is the sole responsibility of the investigators. Boston Scientific's involvement in investigator-sponsored research is limited to providing financial support for research that advances medical and scientific knowledge about our products. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device.

TheraSphere™ Yttrium-90 Glass Microspheres

INDICATION FOR USE: TheraSphere is indicated for use as selective internal radiation therapy (SIRT) for local tumor control of solitary tumors (1-8 cm in diameter) in patients with unresectable hepatocellular carcinoma (HCC), Child-Pugh Score A cirrhosis, well-compensated liver function, no macrovascular invasion, and good performance status. **CONTRAINDICATIONS:** TheraSphere is contraindicated in patients: whose Tc-99m macroaggregated albumin (MAA) hepatic arterial perfusion scintigraphy shows any deposition to the gastrointestinal tract that may not be corrected by angiographic techniques • who show shunting of blood to the lungs that could result in delivery of greater than 16.5 mCi (0.61 GBq) of Y-90 to the lungs. Radiation pneumonitis has been seen rarely in patients receiving doses to the lungs greater than 30 Gy in a single treatment • in whom hepatic artery catheterization is contraindicated, such as patients with vascular abnormalities or bleeding diathesis • who have pulmonary insufficiency (conventionally defined by an arterial oxygen pressure (PaO₂) of < 60 mmHg or oxygen saturation (SaO₂) of < 90%) or severe liver dysfunction, including hepatic encephalopathy, clinically evident ascites or treatment with diuretics for ascites • with portal vein thrombosis (PVT) Type 4 involvement and lack of Tc-99m MAA deposition on the PVT seen on the Tc-99m MAA imaging with >70% tumor replacement in the liver • with comorbidities or poor overall health (e.g., ECOG performance status rating ≥ 2) which may make the patient a poor candidate for locoregional radiation treatment • who are pregnant. **WARNINGS:** The following pre-treatment, high-risk factors (disease characteristics) have been associated with serious adverse events deemed possibly related to use of the device: infiltrative tumor type • tumor nodules too numerous to count • AST or ALT > 5 times ULN • bilirubin > 2 mg/dL • tumor volume > 50% combined with albumin < 3 g/dL. Keep the TheraSphere dose vial upright and stored in its lead pot before and during patient treatment, except as required for radiation measurement. Do not open the dose vial acrylic shield prior to patient treatment. Post-treatment, waste materials require caution to prevent contamination and beta shielding due to residual glass microspheres. **PRECAUTIONS: GENERAL PRECAUTIONS:** As in any intra-arterial procedure, aseptic technique should be practiced, and care should be taken to ensure minimum patient anesthesia exposure extraneous to therapeutic objective. • Consideration of patient comorbidities should be used when determining the type and volume of fluid to infuse via catheter to avoid electrolyte imbalance, fluid shift, and hyperglycemia. • It is important to avoid any aggressive arterial procedure that may lead to arterial spasm that impairs TheraSphere distribution into the perfused liver target volume which may lead to underdosing or non-target deposition of TheraSphere. **PRECAUTION IN PATIENTS WITH IMPAIRED LIVER FUNCTION:** No efficacy or safety data from the LEGACY study are available to support the use of the device in patients with Child-Pugh score B or C cirrhosis. **PRECAUTION IN VULNERABLE PATIENTS:** No effectiveness or safety data are available to support the use of the device in children or breast-feeding women. **ENDOCRINE DISRUPTION, CARINOGENICITY, MUTAGENICITY, TOXICITY TO REPRODUCTION:** Ideally the use of this radioactive device in women of childbearing capability should be performed during the first few (approximately 10) days following the onset of menses. **RADIATION SAFETY:** Radioactive products should be used only by healthcare professionals who are qualified by specific training in the safe use and handling of radionuclides and whose experience and training have been approved by the appropriate government agency authorized to license the use of radionuclides. • As in the use of any radioactive material, ensure minimum radiation exposure to the patient extraneous to the therapeutic objective, and to minimize radiation exposure to workers and others in contact with the patient. **RELEASE AND POST-TREATMENT PRECAUTIONS:** Post-treatment patient care: use universal precautions for body fluid contact. Take Y-90 may be detectable in blood and urine; handle with gloves and dispose as normal body fluids. The radiation field is expected to be less than 1 mrem/h (10 µSv/h) at 3 ft (1 m) from the patient's abdomen. Supplemental shielding and segregation of the patient are not required to maintain exposure to others below regulated limits. • Release instructions: The patient should follow good hygiene (e.g., proper hand washing). Caregivers, family, and others do not require restrictions on patient contact; however, they can minimize their radiation exposure by avoiding prolonged time (>12 hours per day) within 1 ft (0.3 m) of the patient's abdomen for the first week post therapy. Patients should be advised that radiation emitted from the patient may be detectable at security screening (e.g., international travel). • Special precautions post-administration: If the patient requires hospitalization, surgery, medical assessment or treatment regarding any part of their thorax or abdomen within first 2 weeks of treatment, the patient should advise the hospital and treating physician of the Y-90 TheraSphere implant. The physician should consult their radiation safety staff for handling and disposal of liver tissue. • Special liver tissue handling: Special liver tissue handling may be required for post-treatment surgery, explant, or transplant since the glass microspheres remain permanently implanted in the liver tissue. Disclosure of the treatment will be required if cremation is considered. **POTENTIAL ADVERSE EVENTS:** The use of this product leads to irradiation of both tumorous and normal liver tissue. As a result, patients with compromised liver function may be at greater risk of liver function impairment and hence could experience complications. Clinical side effects usually occur within the first 4 to 6 weeks after treatment. Based on clinical trial data, literature reviews and post-market surveillance, adverse events potentially associated with treatment using Y-90 microspheres, including TheraSphere, may include the following: Allergic reaction • Altered liver function, acute or chronic • Anorexia • Anxiety • Ascites • Bile Duct Injury • Bleeding/hemorrhage • Chills • Rigors • Cholecystitis (inflammatory or infectious) • Colitis • Death • Dehydration • Diarrhea • Dizziness • Dyspnea • Edema (any location) • Electrolyte abnormalities • Elevated BUN/creatinine • Fall • Fatigue • Fever • Gastrointestinal bleeding / Hemorrhage • Gastrointestinal ulcer or ulceration • Hepatic encephalopathy • Hepatorenal failure • Hiccups • Hypertension • Hypotension • Infection (any location) • Liver failure, acute or chronic • Lymphopenia • Malaise • Mood alteration • Muscle weakness • Nausea • Neutropenia • Pain (any location) • Pancreatitis • Platelet count abnormalities • Pleural effusion • Portal hypertension • Pre-existing chronic liver disease decompensation • Pulmonary edema • Pulmonary fibrosis • Radiation hepatitis • Radiation induced disease, acute • Radio Embolization Induced Liver Disease (REILD) • Sepsis • Supraventricular arrhythmias • Thrombosis (arterial or venous) • Tumor inflammation (including tumor edema) • Tumor-lysis syndrome • Vomiting • Weight loss. Complications related to the administration procedure itself may include: Allergic reaction; Arterial injury including vessel dissection • Aspiration pneumonia • Bruising/bleeding/hematoma at site • Constipation/abdominal distension • Fatigue • Flushing • Infection • Nausea • Nerve damage. **CAUTION:** Federal (USA) law restricts this device to sale by or on order of a physician. PI-992004-AA. **Note:** Dose to the liver does not exceed 150 Gy. TheraSphere is a registered trademark of Theragnostics Corporation used under license by Biocompatibles UK Ltd. All other trademarks are property of their respective owners.



Peripheral Interventions
300 Boston Scientific Way
Marlborough, MA 01752-1234
www.bostonscientific.com

To order product or for more information contact customer service at 1.888.272.1001.

© 2022 Boston Scientific Corporation or its affiliates. All rights reserved.

PI-823306-AA