



Simplicit^{90Y}™ Personalized Dosimetry Software | KEY FEATURES

Personalized Treatment Simplifying SIRT Workflow

Dosimetry planning is standardized and improves consistency, while post-treatment verification allows assessment of the absorbed dose delivered for each patient.



KEY FEATURES

Offers a variety of tools and functions to help improve Selective Internal Radiation Therapy (SIRT) dosimetry workflow:

- Multimodal image fusion and registration
- Automated tools for liver segmentation (CT, MRI & CBCT)
- Intuitive Lung Shunt Fraction (LSF) calculation tools
- Dosimetric estimates with single-compartment, multi-compartment, and voxel-wise techniques
- Pre- and post-treatment dosimetry
- Vial selector feature allows you to populate the TheraSphere™ order form directly from the software
- Ability to both import and export contours as DICOM RTSS objects
- Saves data and generates comprehensive reports

NEW FEATURES FOR Simplicit^{90Y} 2.6.2


- Updated lung dose tracking for advanced Lung Shunt Calculation
- Ability to load multiple CBCT's into one session
- Calculate normal tissue volume in case of multiple tumors
- Vertical monitor display option
- TheraSphere PDF order form incorporated to allow you to populate, save and use with your standard ordering process

Personalized Dosimetry

Simplicit^{90Y} is a comprehensive software solution for dosimetry planning, allowing you to:

Enhance Consistency And Efficiency

Simplicit^{90Y} allows you to easily incorporate multimodal images with a range of high-performance registration tools, all with one common interface. Coupled with automated and semi-automated tools for segmentation, this reduces inter-user variability and improves dosimetry consistency and planning confidence.



High performance rigid and deformable image registration

Analyze with Ease

Manipulate your data and expedite analysis with tools for the following:


- Image-based (2D, 3D) dosimetry assessment
- Rapid LSF calculation
- Advanced image registration quality control tools

Confirm Y-90 Treatment Quality


Simplicit^{90Y} provides the capability to visualize prospective dose distribution and assess the absorbed dose delivered to structures of interest. By allowing for pre- and post-treatment dosimetry, this software can help determine the effectiveness of a patient's Y-90 SIRT with confidence.

Personalize Dosimetry

Simplicit^{90Y} can be used to interactively tailor the absorbed dose per perfused volume by adjusting the injected activity. The software tools can be customized to a patient's specific tumor presentation and anatomy.



View absorbed dose distribution with isodose contour line display



Show heterogeneity of absorbed dose distribution in critical structures, and assess tumor coverage

Intended Use (US Only)

Simplicit^{90Y} is intended to be used by trained medical professionals for TheraSphere[™] pre-treatment dosimetry planning and post-treatment dosimetry evaluation following Y90 treatment. Simplicity^{90Y} is a medical image and information management system that is intended to receive, transmit, store, retrieve, display and process digital medical images, as well as create, display and print reports from those images. The medical modalities of these medical imaging systems include, but are not limited to, CT, MRI, SPECT and PET. Simplicity^{90Y} provides the user with the means to display, register and fuse medical images from multiple modalities. Simplicity^{90Y} provides tools to create, transform, and modify contours for the user to define objects in medical image volumes for use in TheraSphere pre-treatment dosimetry planning and for post-treatment dosimetry. The objects include, but are not limited to, tumors and normal tissues. For post-Yttrium-90 (Y-90) treatment, Simplicity^{90Y} should only be used for the retrospective determination of dose and should not be used to prospectively calculate dose or for the case where there is a need for retreatment using Y90 microspheres.

Indication for Use (US Only)

Simplicit^{90Y} is a standalone software device that is used by trained medical professionals as a tool to aid in evaluation and information management of digital medical images. Simplicity^{90Y} supports the reading, rendering and display of a range of DICOM compliant imaging and related formats including but not limited to CT, PT, NM, SPECT, MR, SC, RTSS. Simplicity^{90Y} enables the saving of sessions in a proprietary format as well as the export of formats including CSV and PDF files. Simplicity^{90Y} is indicated, as an accessory to TheraSphere, to provide pre-treatment dosimetry planning support including Lung Shunt Fraction estimation (based on planar scintigraphy) and liver single-compartment MIRD schema dosimetry in accordance with TheraSphere labelling. Simplicity^{90Y} provides tools to create, transform, and modify contours/Regions of Interest for calculation of Lung Shunt Fraction and Perfused Volume. Simplicity^{90Y} includes features to aid in TheraSphere dose vital selection, dose vital ordering and creation of customizable reports. Simplicity^{90Y} is indicated for pre-treatment dosimetry and evaluation following Yttrium-90 (Y-90) microsphere treatment. Simplicity^{90Y} provides tools to create, transform, and modify contours/Regions of Interest for the user to define objects in medical image volumes to support TheraSphere post-Y90 treatment calculation and evaluation. The objects include, but are not limited to, tumors and normal tissues, and liver volumes. Simplicity^{90Y} is indicated for registration, fusion display and review of medical images allowing medical professionals to incorporate images, such as CT, MRI, PET, CBCT and SPECT in TheraSphere Yttrium-90 (Y-90) microspheres pre-treatment planning and post-Y90 treatment evaluation. For post-Yttrium-90 (Y-90) treatment, Simplicity^{90Y} should only be used for the retrospective determination of dose and should not be used to prospectively calculate dose or for the case where there is a need for retreatment using Y-90 microspheres. PI-994110-AA

TheraSphere[™] Yttrium-90 Glass Microspheres

INDICATION FOR USE: TheraSphere is indicated for use as selective internal radiation therapy (SIRT) for local tumor control of solitary tumors (1-8 cm in diameter), in patients with unresectable hepatocellular carcinoma (HCC), Child-Pugh Score A cirrhosis, well-compensated liver function, no macrovascular invasion, and good performance status. **CONTRAINDICATIONS:** TheraSphere is contraindicated in patients: whose Tc-99m macroaggregated albumin (MAA) hepatic arterial perfusion scintigraphy shows any deposition to the gastrointestinal tract that may not be corrected by angiographic techniques • who show shunting of blood to the lungs that could result in delivery of greater than 16.5 mCi (0.61 GBq) of Y-90 to the lungs. Radiation pneumonitis has been seen early in patients receiving doses to the lungs greater than 30 Gy in a single treatment • in whom hepatic artery catheterization is contraindicated, such as patients with vascular abnormalities or bleeding diathesis • who have pulmonary insufficiency (constrictive pericarditis, pulmonary hypertension, or severe pulmonary hypertension) • who have portal vein thrombosis (PVT) Type 4 involvement and lack of Tc-99m MAA deposition on the PVT seen on the Tc-99m MAA imaging with >70% tumor replacement in the liver (with comorbidities or poor overall health (e.g., ECOG performance status rating >2) which may make the patient a poor candidate for locoregional radiation treatment • who are pregnant. **WARNINGS:** The following pre-treatment, high-risk factors (disease characteristics) have been associated with serious adverse events deemed possibly related to use of the device: infiltrative tumor type • tumor nodules too numerous to count • AST or ALT > 5 times ULN • bilirubin > 2.0 mg/dL • tumor volume > 50% combined with albumin < 3 g/dL. Keep the TheraSphere dose vial upright and stored in its lead pot before and during patient treatment, except as required for radiation measurement. Do not open the dose vial acrylic shield prior to patient treatment. Post-treatment, waste materials require caution to prevent contamination and beta shielding due to residual glass microspheres. **PRECAUTIONS: GENERAL PRECAUTIONS:** As in any intra-arterial procedure, aseptic technique should be practiced, and care should be taken to ensure minimum patient anesthesia exposure extraneous to therapeutic objective. • Consideration of patient comorbidities should be used when determining the type and volume of fluid to infuse via catheter to avoid electrolyte imbalance, fluid shift, and hyperglycemia. • It is important to avoid any aggressive arterial procedure that may lead to arterial spasm that impairs TheraSphere distribution into the perfused liver target volume which may lead to underdosing or non-target deposition of TheraSphere. **PRECAUTION IN PATIENTS WITH IMPAIRED LIVER FUNCTION:** No efficacy or safety data from the LEGACY study are available to support the use of the device in patients with Child-Pugh score B or C cirrhosis. **PRECAUTION IN VULNERABLE PATIENTS:** No effectiveness or safety data are available to support the use of the device in children or breast-feeding women. **ENDOCRINE DISRUPTION, CARCINOGENICITY, MUTAGENICITY, TOXICITY TO REPRODUCTION:** Ideally the use of this radioactive device in women of childbearing capability should be performed during the first few (approximately 10) days following the onset of menses. **RADIATION SAFETY:** Radioactive products should be used only by healthcare professionals who are qualified by specific training in the safe use and handling of radionuclides and whose experience and training have been approved by the appropriate government agency authorized to license the use of radionuclides. • As in the use of any radioactive material, ensure minimum radiation exposure to the patient, extant to the therapeutic objective, and to minimize radiation exposure to workers and others in contact with the patient. **RELEASE AND POST-TREATMENT PRECAUTIONS:** Post-treatment patient care: use universal precautions for body fluid contact. Tare Y-90 may be detectable in blood and urine; handle with gloves and dispose as normal body fluids. The radiation field is expected to be less than 1 mrem/h (0.01 μSv/h) at 3 ft (1 m) from the patient's abdomen. Supplemental shielding and segregation of the patient are not required to maintain exposure to others below regulated limits. • Release instructions: The patient should follow good hygiene (e.g., proper hand washing). Caregivers, family, and others do not require restrictions on patient contact; however, they can minimize their radiation exposure by avoiding prolonged time (>12 hours per day) within 1 ft (0.3 m) of the patient's abdomen for the first week post therapy. Patients should be advised that radiation emitted from the patient may be detectable at security screening (e.g., international travel). • Special precautions post-administration: The Y-90 TheraSphere implant. The physician should consult their radiation safety staff for handling and disposal of liver tissue. • Special liver tissue handling: Special liver tissue handling may be required for post-treatment surgery, explant, or transplant since the glass microspheres remain permanently implanted in the liver tissue. Disclosure of the treatment will be required if reoperation is considered. **POTENTIAL ADVERSE EVENTS:** The use of this product leads to irradiation of both tumorous and normal liver tissue. As a result, patients with compromised liver function may be at greater risk of liver function impairment and hence could experience complications. Clinical side effects usually occur within the first 4 to 6 weeks after treatment. Based on clinical trial data, literature reviews and post-market surveillance, adverse events potentially associated with treatment using Y-90 microspheres, including TheraSphere, may include the following: Allergic reaction • Altered liver function, acute or chronic • Anorexia • Anxiety • Ascites • Bleeding/hemorrhage • Chills • Rigors • Cholecystitis (inflammatory or infectious) • Colitis • Death • Dehydration • Diarrhea • Dizziness • Dyspnea • Edema (retention) • Electrolyte abnormalities • Elevated BUN/creatinine • Fall • Fatigue • Fever • Gastrointestinal bleeding / hemorrhage • Gastrointestinal iliac or ulceration • Hepatic encephalopathy • Hepatorenal failure • Hypertension • Hypertension • Infection (any location) • Liver failure, acute or chronic • Lymphopenia • Malaise • Mood alteration • Muscle weakness • Nausea • Neurotoxicity • Pain (any location) • Pancreatitis • Platelet count abnormalities • Pleural effusion • Portal hypertension • Pre-existing chronic liver disease decompensation • Pulmonary edema • Pulmonary fibrosis • Radiation hepatitis • Radiation induced disease, acute • Radio Embolization Induced Liver Disease (REILD) • Sepsis • Supraventricular arrhythmia • Thrombosis (arterial or venous) • Tumor irradiation (including tumor edema) • Tumor lysis syndrome • Vomiting • Weight loss. Complications related to the administration procedure itself may include: Allergic reaction: Arterial injury including vessel dissection • Aspiration pneumonia • Bruising/bleeding/hematoma at site • Constipation/abdominal distension • Fatigue • Flushing • Infection • Nausea • Nerve damage. **CAUTION:** Federal (USA) law restricts this device to sale by or on order of a physician. PI-992004-AA. **Note:** Dose to the liver does not exceed 150 Gy. TheraSphere is a registered trademark of Theragnostics Corporation used under license by Biocompatibles UK Ltd. All other trademarks are property of their respective owners.



Peripheral Interventions
300 Boston Scientific Way
Marlborough, MA 01752-1234
www.bostonscientific.com

To order product or for more information
contact customer service at 1.888.272.1001.

© 2022 Boston Scientific Corporation
or its affiliates. All rights reserved.

PI-988103-AB