

# IMPERIAL RCT SUMMARY


The world's first head-to-head DES SFA Trial, evaluating Boston Scientific Corporation's Eluvia™ Drug-Eluting Vascular Stent System and Cook Medical's Zilver™ PTX™ Stent




## Sustaining strong results through five years

The results of the IMPERIAL RCT show that Eluvia Drug-Eluting Stent (DES) is clinically effective and safe in treating patients with symptomatic SFA disease both in the short-term during the height of restenosis risk, and long-term out to five years.

Eluvia DES demonstrated superiority over Zilver PTX<sup>1</sup> with a statistically significant primary patency through 1-Year



Eluvia DES showed lower revascularisation rates than Zilver PTX through 5 years<sup>2</sup> with statistical significance<sup>3</sup> at 2-Years



## IMPERIAL TRIAL 2-YEAR CLINICAL RESULTS

Excellent Patient Follow-up at 24-Months (~90%)

	IMPERIAL RCT <sup>4</sup> (n = 309)	IMPERIAL Long Lesions <sup>5,6</sup> (n = 56)	Diabetic Subgroup <sup>7,8</sup> (n = 116)	Severe/ Moderate Calcium Subgroup <sup>8</sup> (n = 193)	CTO Subgroup <sup>8</sup> (n = 96)
Study Design	RCT, global multicenter	Single arm, global multicenter	RCT, global multicenter	RCT, global multicenter	RCT, global multicenter
<b>24-month primary patency rate*</b>	<b>83.0%</b>	<b>77.2%</b>	<b>85.7%</b>	<b>85.0%</b>	<b>76.4%</b>
Lesion length (mm)	86.5	162.8	87.0	89.9	94.4
Severe calcification	40%	28%	46%	n/a	n/a
Total occlusions	31%	32%	25%	n/a	100%

**Durable, consistent outcomes in long and complex lesions**

Highest primary patency ever reported at 2 years\*\*

Highly durable outcomes in ~16cm lesions at 2 years

TLR (12%) in line with overall cohort and low stent thrombosis rate (0.9%)

Remarkable primary patency and <10% TLR in heavy calcium

Highly durable outcomes in CTOs at 2 years

**Eluvia DES patients on average avoided reintervention 6 months longer than Zilver PTX patients at 3-Years<sup>2†</sup>**

Zilver PTX

13 months

p=0.0058

Eluvia DES

19 months



**IMPERIAL TRIAL OBJECTIVE**

Evaluate the safety and effectiveness of the Boston Scientific Corporation Eluvia™ Drug-Eluting Vascular Stent System for treating Superficial Femoral Artery (SFA) and/or Proximal Popliteal Artery (PPA) lesions up to 140mm in length.

**IMPERIAL TRIAL DESIGN**

Global multi-center, 2:1 randomisation against Cook Medical's Zilver™ PTX™ Stent, controlled, single-blind, non-inferiority trial; core lab adjudicated.

- 465 (RCT) patients across 64 sites
- 5-year follow-up
- Degree of stenosis ≥ 70% (visual angiographic assessment)
- Vessel diameter ≥ 4mm and ≤ 6mm
- Total lesion length ≥ 30mm and ≤ 140mm

**BASELINE CHARACTERISTICS**

Patient Demographics	Eluvia (n=309)	Zilver PTX (n=156)
Age (Years)	68.5±9.5	67.8±9.4
Male Gender	66.0%	66.7%
Diabetes Mellitus	41.7%	43.6%
History of Smoking	86.1%	84.0%

Lesion Characteristics	Eluvia (n=309)	Zilver PTX (n=156)
Target Lesion Length (mm)	86.5±36.9	81.8±37.3
Severely Calcified	40.1%	32.3%
Total Occlusions	31.2%	30.3%
Extending into Distal SFA	66.3%	65.4%

\* Intention to treat. Kaplan-Meier estimate utilising time-to-event of clinically-driven TLR up to 730 days and Duplex Ultrasound data at 24 months. Primary patency defined as duplex ultrasound PSVR

≤2.4, in the absence of clinically-driven target lesion revascularisation or bypass of the target lesion, as assessed by the DUS core lab.

\*\*Highest-two-year primary patency based on 24-month Kaplan-Meier estimates reported for IMPERIAL, IN.PACT SFA, ILLUMENATE, LEVANT II and Primary Randomisation for Zilver PTX RCT.

† Among patients who underwent a CD-TLR within 3 years of the index procedure

1. IMPERIAL Trial: A global randomized controlled multi-center trial with 2:1 randomization of the Eluvia™ Drug-Eluting Stent against Cook Medical's Zilver™ PTX™ Stent, single-blind, non-inferiority design; independent core lab adjudication. Superiority determined in a post hoc analysis that was specified prior to unblinding. 12-Month Primary Patency rate of 86.8% in the Eluvia arm vs. 77.5% in the Zilver PTX arm (p-value = 0.0144). Gray WA, Lancet. 2018 Sep 24, pii: S0140-6736(18)32262-1.

2. Gray W. 5-year Results from the IMPERIAL Randomized Study of Eluvia and Zilver PTX Drug-eluting Stents and Long Lesion Substudy for Femoropopliteal Artery Disease; CRT 2023, Washington DC Feb 27, 2023.

3. Intention to treat. Iida O, VIVA 2019. RCT, randomized controlled trial; TLR, target lesion revascularization.

4. In IMPERIAL RCT, Eluvia K-M Primary Patency was 83% vs. 77.1% for Zilver PTX at 24 months, p=0.1008.

5. Golzaar, J. et al, Journal of Endovascular Therapy, Jan 2020. <https://doi.org/10.1177/1526602820901723>.

6. Vermassen, F. VIVA Late-Breaking Clinical Trials June 2020.

7. In IMPERIAL Diabetic Subgroup, Eluvia K-M Primary Patency was 95.2% vs. 81.5% for Zilver PTX at 12 months. Diabetic = Medically Treated Diabetes.

8. Dr. Gray LINC presentation -2. Gray, W. 2 year Outcomes from the IMPERIAL Randomized Head to Head Study of Eluvia DES and ZilverPTX. LINC 2020.

CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings, and instructions for use can be found in the product labelling supplied with each device or at [www.IFU-BSCI.com](http://www.IFU-BSCI.com). Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. This material not intended for use in France.



**bostonscientific.eu**

© 2023 Boston Scientific Corporation or its affiliates. All rights reserved.

PI-1542203-AA

CE 0344