



# Salvage Focal Prostate Cryotherapy of 1.2cc Anterior Lesion


**Taimur T Shah**

Charing Cross Hospital, Imperial College Healthcare NHS Trust and Imperial College London (ICL)




## PRESENTATION

- 72-year-old man presented with a possible recurrence having had brachytherapy for Gleason 3+3 prostate cancer ten years previously
  - PSA nadir after radiotherapy was 0.3ng/ml and had risen to 7ng/ml
  - Otherwise fit and well (had undergone an open appendectomy during his thirties)
- The patient underwent staging investigations consisting of a prostate mpMRI, choline PET-CT and a NM bone scan. Scans showed no evidence of metastatic disease and so the patient underwent a template mapping and targeted biopsy
- 3T mpMRI showed a 1.2cc, PI-RADS2, 5/5 right anterior lesion **A**
- Transperineal prostate mapping and targeted biopsy confirmed right-sided Gleason 4+3, 8mm maximum cancer core length (MCCL) prostate cancer **B**
- Diagnosis of recurrent prostate cancer T2aNOM0; the patient was offered all salvage treatment options and chose focal cryotherapy



**Presentation:** Multiparametric MRI reveals 1.2cc right anterior lesion. Top: T2 axial view (left) and ADC map (right). Bottom: b2000 diffusion-weighted image (left) and dynamic contrast-enhanced image (right)



**Pre-cryotherapy:** MRI diagram of right anterior lesion


Diagrams adapted from Villers A et al. Diagnostic and Interventional Imaging (2012) 93, 262–267

//

*Focal salvage cryotherapy has managed to provide oncological control and maintain this patient's quality of life in his later years. //*

# Salvage Focal Prostate Cryotherapy of 1.2cc Anterior Lesion

Doctor Taimur T Shah Charing Cross Hospital, Imperial College Healthcare NHS Trust and Imperial College London (ICL)



## TREATMENT

- Right anterior focal cryotherapy performed using five IceSeed™ 1.5 needles
- Two thermocouples used for temperature monitoring within the lesion and at Denonvilliers' fascia
- Urethral warmer inserted over a wire after a flexible cystoscopy confirmed no needles present within the urethra or bladder
- Two 10-minute freeze cycles




\* Each 10-minute freeze cycle must include  $\geq 3$  minutes at  $< -40^{\circ}\text{C}$  at the most lateral border of the tumour margin




## OUTCOME

- The PSA post treatment dropped significantly to 0.05ng/ml and has remained stable in the two years since, indicating a very good response
- Early and late post-operative MRI images showed very good tumour coverage and no residual cancer **C D**
- Repeated annual mpMRI has not shown any evidence of residual or recurrent cancer



Post cryotherapy: Six-week follow-up MRI dynamic contrast-enhanced T1 images



Post cryotherapy: One-year follow-up MRI. Left to right: T2 axial view; ADC map; b2000 diffusion-weighted image; dynamic contrast-enhanced image



## CONCLUSION

- Salvage cases are generally regarded as being difficult, particularly in patients who have had previous brachytherapy. Outcomes from salvage surgery – whether it be a prostatectomy or a minimally invasive therapy such as cryotherapy – are generally worse than in patients with primary disease. Salvage prostatectomy, however, carries significant morbidity for the patient with very high rates of incontinence and erectile dysfunction and a risk of rectal injury, along with a long post-operative recovery period<sup>1-4</sup>
- This patient had a very well-localised anterior recurrence after previous brachytherapy. He underwent salvage focal cryotherapy as a day-case procedure and had a catheter *in situ* for one week. He had no major adverse events and although he suffered from minor incontinence, this did not require him to wear a pad. His PSA dropped significantly to 0.05ng/ml and has remained stable
- Focal salvage cryotherapy has managed to provide oncological control and also maintain this patient's quality of life in his later years

### References:

1. Kimura M *et al.* BJU International 2009; 105:191-201.
2. Shah T *et al.* Anticancer Ther 2014; 14(11):1337-47. 3. Chade DC *et al.* Eur Urol 2012; 61(5): 961-71. doi: 10.1016/j.eururo.2012.01.022.
4. Pearce SM *et al.* Urol Oncol 2015; 33(4):163.e1-6. doi: 10.1016/j.urolonc.2015.01.016.

ADC: Apparent diffusion coefficient  
MCCL: Maximum cancer core length  
3T mpMRI: 3 Tesla multiparametric magnetic resonance imaging  
MRI: Magnetic resonance imaging  
NM: Nuclear medicine

PET-CT: Positron emission tomography-computed tomography  
PI-RADS: Prostate Imaging Reporting and Data System  
PSA: Prostate specific antigen

Results from case studies are not predictive of results in other cases. Results in other cases may vary. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labelings supplied with each device. Information for the use only in countries with applicable health authority product registrations. Material not intended for use in France.

PI-880008-AA

**Boston Scientific**  
Advancing science for life™

[www.bostonscientific.eu](http://www.bostonscientific.eu)

© 2020 Boston Scientific Corporation or its affiliates. All rights reserved.