

SUMMARY

Boston Scientific pacemakers, defibrillators, and S-ICDs include an identifier in the header of the pulse generator that is visible via X-ray or fluoroscopy. The identifier is useful in determining the device manufacturer and model, allowing a clinician to determine the appropriate programmer to communicate with the device. This article includes a reference table of Boston Scientific device identifiers, and includes images of the X-ray identifier location.

Products Referenced

All referenced Boston Scientific pacemakers, ICDs, CRT-Ds, CRT-Ps, and S-ICDs

Products referenced are unregistered or registered trademarks of Boston Scientific Corporation or its affiliates. All other trademarks are the property of their respective owners.

For comprehensive information on device operation, reference the full instructions for use found at: www.bostonscientific-elabeling.com.

CAUTION: Law restricts this device to sale by or on the order of a physician. Indications, contraindications, precautions and warnings can be found with product labeling supplied with each device. Products referenced herein may not be approved in all geographies. Information is for the use in countries with applicable Health Authority product registrations.

All graphics produced by Boston Scientific Corporation, unless otherwise noted.

CRT-D: Cardiac Resynchronization Therapy Defibrillator
CRT-P: Cardiac Resynchronization Therapy Pacemaker
ICD: Implantable Cardioverter Defibrillator
S-ICD: Subcutaneous Implantable Defibrillator

Contact Information

www.bostonscientific.com

Americas

Technical Services

LATITUDE™ Customer Support

1.800.CARDIAC (227.3422)

+1.651.582.4000

Patient Services

1.866.484.3268

Europe, Middle East, Africa

Technical Services

+32 2 416 7222

inttechservice@bsci.com

LATITUDE Customer Support

latitude.europe@bsci.com

Japan

Technical Services

japantechservice@bsci.com

LATITUDE Customer Support

japan.latitude@bsci.com

Asia Pacific

Technical Services

+61 2 8063 8299

aptechservice@bsci.com

LATITUDE Customer Support

latitudeasiapacific@bsci.com

© 2018 by Boston Scientific Corporation or its affiliates. All Rights Reserved.

X-ray Identification in Boston Scientific Pulse Generators

X-ray Identifier

Boston Scientific pulse generators, including pacemakers, defibrillators, and S-ICDs, have an identifier that is visible on X-ray film or under fluoroscopy. This identifier provides noninvasive confirmation of the manufacturer and consists of letters (manufacturer identification) and a number (cross-references to a programmer software application) that is required in order to communicate with the device.

The X-ray identifier is embedded in the header of the device. In general, the X-ray identifier will be visible by X-ray or fluorography at the approximate location identified in Figure 1. Note that the location may vary slightly based on the device family or the header design.




Figure 1. X-ray identifier

The tables below identify current generations of Boston Scientific pacemakers, defibrillators and S-ICDs, the programmer software application used by the device, the X-ray identifier (manufacturer letters + software identification number), and an actual X-ray image.

Identification via a Boston Scientific Programmer

The pulse generator can also be identified by interrogating it with a Boston Scientific Programmer. The model number is stored in device memory and is shown on the programmer upon interrogation. If the pulse generator is not manufactured by Boston Scientific or the model is not supported by the programmer being used, a message will be displayed when an attempt is made to interrogate the device. For example, the message displayed on the Model 3120 ZOOM™ LATITUDE™ Programmer is shown in Figure 2.





Figure 2. PG Not Identified Message

Table 1. X-ray identifier - Boston Scientific Pacemakers and CRT-Ps

Note that images below include a representative sample from each device family - the location of the X-ray identifier may vary slightly based on the device family or the design of the header.



**ACCOLADE™, PROPONENT™, ESSENTIO™, ALTRUA™ 2 pacemakers
VISIONIST™, VISIONIST™ X4, VALITUDE™, VALITUDE™ X4 CRT-Ps**


Uses Programmer Software Application: Model 2869/3869

X-ray identifier: **BSC012**

ACCOLADE MRI, PROPONENT MRI, ESSENTIO MRI pacemakers

Uses Programmer Software Application: Model 2869/3869

X-ray identifier: **▲BSC012**



**VITALIO™, FORMIO™, INGENIO™, ADVANTIO™ pacemakers
INLIVEN™, INTUA™, INVIVE™ CRT-Ps**


Uses Programmer Software Application: Model 2869/3869

X-ray identifier: **BSC011**

VITALIO MRI, FORMIO MRI, INGENIO MRI, ADVANTIO MRI pacemakers

Uses Programmer Software Application: Model 2869/3869

X-ray identifier: **▲BSC011**




ALTRUA™ pacemakers

Uses Programmer Software Application: Model 2892


X-ray identifier: **BOS003**

Table 2. X-ray identifier - Boston Scientific ICDs and CRT-Ds

Note that images below include a representative sample from each device family - the location of the X-ray identifier may vary slightly based on the device family or the design of the header.



RESONATE™, PERCIVA™, CHARISMA™, VIGILANT™, MOMENTUM™ ICDs
RESONATE, CHARISMA, VIGILANT, MOMENTUM CRT-Ds
AUTOGEN™, DYNAGEN™, INOGEN™, ORIGEN™ ICDs
AUTOGEN, DYNAGEN, INOGEN, ORIGEN CRT-Ds
 Uses Programmer Software Application: Model 2868/3868
 X-ray identifier: **BSC140**




INCEPTA™, PUNCTUA™, ENERGEN™ ICDs
INCEPTA, PUNCTUA, ENERGEN CRT-Ds
 Uses Programmer Software Application: Model 2868/3868
 X-ray identifier: **BSC120**




TELIGEN™ ICD and COGNIS™ CRT-D
 Uses Programmer Software Application: Model 2868/3868
 X-ray Identifier: **BOS112**

Table 3. X-ray identifier - Boston Scientific S-ICDs

Note that images below include a representative sample from each device family - the location of the X-ray identifier may vary slightly based on the device family or the design of the header.



EMBLEM and EMBLEM MRI S-ICD
Uses Programmer Software Application: Model 2877
X-ray identifier: **BSC507**



SQ-RX S-ICD
Uses Programmer Software Application: Model 2877
X-ray identifier: **CH1010**