
Swiss LithoClast® Select System


DIRECTIONS FOR USE
MODE D'EMPLOI

© Copyright EMS SA
FB-411/US ed. 2015/06
(part of BSC ref. 840-302)

Rx ONLY Caution! Federal (USA) law restricts this device
to sale by or on the order of a physician

**Boston
Scientific**


COMPONENTS

- 1 Control unit
- 2 Pinch valve
- 3 Holder for Stone Catcher
- 4 Rotary knob for ultrasound lithotripsy power setting
- 5 Ultrasound standby/ON switch with control lamp
- 6 Ultrasound handpiece cord connector
- 7 Pulse mode selection [single pulse/ multiple pulse]
- 8 Toggle mode switch for increase/decrease buttons
- 9 Decrease button (in selected mode)
- 10 Increase button (in selected mode)
- 11 Pneumatic handpiece cord connector
- 12 Rotary knob for pneumatic lithotripsy power setting
- 13 Pneumatic standby/ON switch with control lamp
- 14 Ultrasound mode display screen
 - a Ultrasound power level [%]
 - b Ultrasound duty cycle setting [%]
 - c Ultrasound use since last stand by [sec/hr]
- 15 Pneumatic mode display screen
 - d Pneumatic power level [%]
 - e Number of impulses since last stand by [sec/hr]
 - f Pulse frequency [Hz]

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- 16 Power switch ON/OFF
- 17 Power socket
- 18 Fuse box
- 19 Connector for the foot pedal
- 20 Equipotential connector
- 21 Connector for compressed air supply
- 22 Manufacturer identification plate
- 23 Dual foot pedal (combined ultrasound/pneumatic)
- 24 Stone Catcher single use with roller clamp
- 25 Spare fuses
- 26 Power cord
- 27 Compressed air hose

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DEAR CUSTOMERS,

The Swiss LithoClast® Select System is a combination ultrasonic and pneumatic urological lithotripter used for the fragmentation of urinary tract calculi in the kidney, ureter, and bladder.

These operational instructions are applicable only for control units equipped with software version "rel. 3.x", which includes the Swiss LithoClast® Select System and the Swiss LithoClast® Ultra Upgraded System. The software version is displayed when the control unit is turned on. Prior to use, please verify that "rel. 3.x" is shown briefly on the display panel. Please note that units with software "rel. 3.x" are hereafter referred to as the "Swiss LithoClast Select System".

Contraindications:

Use of the product is contraindicated in patients with the following:

- Active bleeding disorders
- Solitary functioning kidney
- Creatinine greater than or equal to 3 µg %
- During pregnancy
- Under the age of 18
- Stricture and obstruction problems
- With an electrical stimulator implanted e.g. pacemaker.

Potential Complications:

Potential complications associated with fragmentation of urinary tract calculi by pneumatic and/or ultrasound energy are:

- Perforation
- Hemorrhage
- Lesion
- Stone migration
- Pain/colic
- Macroscopic hematuria
- Infection
- Ureteral obstruction

Instructions concerning the specific components (handpieces, probes) are provided separately. Read the instruction manual delivered with each of the components before using the Swiss LithoClast® Select System.

Do not hesitate to contact your BSC representative for further information.

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PLEASE READ THIS BEFORE BEGINNING WORK!

About this manual

Note that the English version of this manual is the master from which translations derive. In case of any discrepancy, the binding version is the English text.

These operating instructions are to ensure the correct installation and use of this product. Always keep this instruction close at hand.

 Please read these operating instructions carefully as they explain all the most important details and procedures. Please pay special attention to the safety precautions.

 This instruction manual describes combined ultrasonic and pneumatic lithotripsy, as well as ultrasonic or pneumatic lithotripsy application modes.

To prevent possible patient or user injury or damage to the Swiss LithoClast® Select System, please note the corresponding directives:

 **Danger**
Risk of patient or user injury

 **Caution**
Risk of damage to the device or environmental harm

 **Please note**
Useful additional information and hints

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Autorisé

CONTENTS


Installation	3	Maintenance.....	28
Connections	3	Fuses replacement	30
Electrical Safety precautions	3	Safety precautions	32
Disposable Stone Catcher device	4	Storing of the product	34
Operating the unit	6	Accessories.....	34
Suction control	6	Product warranty.....	34
Ultrasound lithotripsy	8	Warranty	34
Preparation.....	8	Customer Service	36
Ultrasound power.....	10	Electromagnetic compatibility.....	36
Duty cycle.....	10	Symbols	38
Activation.....	10	Technical data.....	44
Pneumatic lithotripsy.....	12	Troubleshooting	46
Preparation.....	12		
Pulse modes.....	14		
Pulse frequencies	14		
Pneumatic power.....	14		
Activation.....	14		
Ultrasound and pneumatic combined lithotripsy mode	16		
Preparation.....	16		
Endoscopy	18		
Treatment recommendations.....	20		
General.....	20		
Recommended settings	20		
Selection between individual and combined modes	22		
Warnings	24		
Post-treatment procedure.....	24		
Cleaning, disinfecting and sterilizing	26		
Preparation and cleaning	26		
Handpieces and accessories.....	26		

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INSTALLATION

Clean, disinfect and sterilize all the components before first use and after all subsequent uses as described in the "Reprocessing Instructions" manual. Assemble the application parts under sterile conditions.

Perform function testing according to the directions for use supplied with the respective accessory.

Do not cover the unit. Leave a free space of at least 5 cm around the device for cooling.

CONNECTIONS

Connect the power supply cord to the unit. To prevent damage to the unit, check that its rated voltage meets the local line voltage. Connect only to a FI protected main supply (FI = Residual current protection).

Connect the quick connector of the compressed air tube to the unit, by pushing the connector until it clicks on.

Use only medical-quality, dry compressed air, filtered at a minimum of 5 µm. This supply of compressed air can be either:

- the central compressed air supply system (wall socket)
- a separate compressor
- a cylinder (ISO 7396)

Connect the foot pedal cord to the unit. Pay attention to the connector orientation and fasten the securing nut.

Be sure that the connectors are perfectly dry before making the connection. To disconnect, pull only the connector sleeve; never pull on the cord itself.

When applicable connect the equipotential plug to other medical devices used.

Electrical Safety precautions

To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective earth.

No modification shall be made on this product.

The main switch of the unit must be accessible at any time.

When the main switch is on the "0" position, the unit is disconnected from the supply network.

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

Disposable Stone Catcher device

The disposable Stone Catcher device is provided sterile. Read the instruction manual delivered with it before performing installation.

Connect the short silicone tubing end with the conical connector to the external suction source.

 The pinch valve default state is closed and opens when the ultrasound (left) foot pedal is depressed halfway.

 Install a roller clamp on the suction tubing to enable control of the suction flow independently of the lower pressure.

 The sterile disposable Stone Catcher device cannot be reprocessed.

Stone Catcher device tubing

 The use of other tubing than the one supplied with the Stone Catcher device can either create an incomplete occlusion of the tubing in the pinch valve, or the tubing can collapse during use. Use only the Stone Catcher device tubing provided by BSC.

The suction tube (standard suction tubing with 3 m length and 5 mm internal diameter) must be sterile and single use. Use of a transparent tubing allows for observation of fragment transport to visualize a possible obstruction.

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OPERATING THE UNIT

Suction control

 Recommendations on suction settings for lithotripsy using BSC handpieces and probes are provided separately.

 Suction flow can be adjusted by using the roller clamp on the suction tubing.

 Excessive suction can impair the endoscopic view, collapse the organ, or damage the mucosa. Always check that the suction tubing is inserted into the pinch valve and that the pinch valve is functioning properly.

 Make sure that suction is always operational. Continuous suction during ultrasound lithotripsy is required to evacuate fragments and liquids as well as to ensure the cooling of the probe and the handpiece.

 Interruption of the suction must be avoided. It is operator's responsibility to continuously monitor the suction flow to check that it is operational and that the level of lower pressure is correctly set.

 Fragments blocked in the lumen of the probe and the handpiece may lead to loss of suction and heating of the probe. If blockage occurs, stop lithotripsy to remove fragments from the probe and from the handpiece lumen before continuing.

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
ULTRASOUND LITHOTRIPSY

Preparation

 Please refer to the Ultrasound Handpiece Directions for use for handpiece preparation and assembly details. Never connect a partially assembled handpiece to the control unit.

 Be sure that the control unit is off while connecting the handpiece.

The Standby/On button light will turn orange 5 seconds after turning on the unit. Press the Standby/On button to start the corresponding mode.


 The cumulative ultrasound usage time (seconds/hours) is displayed briefly on the left display when the unit is turned on.

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Ultrasound power

The left rotary knob is used to adjust the ultrasound power from 10 % to 100 % (in 10 % increments).

Duty cycle

The duty cycle is the ratio between ultrasound ON and ultrasound OFF over an interval of time of 80 milliseconds. Set the duty cycle using the arrow keys.

The last setting of the duty cycle is kept in memory at power off.

The handpiece and probe error detection is possible only with a duty cycle set at 100 % (see Troubleshooting chapter).

Activation

- Press the ultrasound (left) foot pedal halfway to activate the suction (pinch valve opens).
- Press the ultrasound (left) foot pedal completely to activate the ultrasound handpiece and probe with suction.

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
PNEUMATIC LITHOTRIPSY

Preparation

 Refer to the Pneumatic Handpiece Directions for use for handpiece preparation and assembly details. Never connect a partially assembled handpiece to the control unit.

 Be sure that the control unit is off while connecting the handpiece.

The Standby/On button light will turn orange 5 seconds after turning on the unit. Press the Standby/On button to start the corresponding mode.


 The cumulative number of pulses is displayed briefly on the right display when the unit is turned on.

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Pulse modes

Select either single or multiple pulse mode with the "Litho mode" button.

Pulse frequencies

When the pulse mode is set to multiple, the frequency can be selected with the arrow keys. Frequency is adjustable on a 1Hz to 12Hz range. The last setting of the pulse mode and frequency is kept in memory at power off.

Pneumatic power

The right rotary knob is used to adjust the pneumatic power in percent from 0 % to 100 % (in 10 % increments).

Activation

Press the pneumatic (right) foot pedal to activate the pneumatic pulses.

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
ULTRASOUND AND PNEUMATIC COMBINED LITHOTRIPSY MODE

Preparation

 Refer to the handpieces manuals for handpiece preparation and assembly details. Never connect a partially assembled handpiece to the control unit.

 Be sure that the control unit is off while connecting the handpieces.


The Standby/On button light will turn orange 5 seconds after turning on the unit. Press the Standby/On button to start the corresponding mode.

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 For preparation and settings, refer to the chapters on ultrasound and pneumatic lithotripsy.

Activate the ultrasound mode as well as the pneumatic mode on the control unit.

The "Select" button toggles between the adjustment of the duty cycle and the setting of the pneumatic pulse frequency. The function in use is indicated by a dot either in the duty cycle display or in the frequency display. The arrow keys are used to adjust the values.

 Press both foot pedals to activate pneumatic and ultrasound lithotripsy mode with suction simultaneously (place the foot over both pedals).

 The individual lithotripsy modes can be selectively controlled by pressing only the related foot pedal (left for ultrasound / right for pneumatic).

ENDOSCOPY

 It may happen that the probe tip is hidden in folds of the edematous tissue. Endoscopes with a 5° to 15° viewing angle offer a better visualization of the probe tip in the near range of the endoscope.

 Ensure that the tip of the pneumatic or ultrasound probe does not extend beyond the tip of the endoscope by more than 10 to 20 mm when it is fully inserted in the endoscope working channel.

 Risk of mucosal perforation: When using a 0.8 mm or 1.0 mm pneumatic probe with an ureteroscope, it is recommended to adjust (if needed) the length of the probe to the endoscope length.

 Adjusting the length of the probe will help limit the risk of perforation and reduce the energy loss, particularly with thin probes.

 To avoid damage to the endoscope, the probes should not be activated with their tip still inside the working channel of the endoscope.

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TREATMENT RECOMMENDATIONS

General

Intracorporeal lithotripsy by use of ultrasound and pneumatic energies, used simultaneously or individually, requires experience and skill in endoscopic technique. A good balance of irrigation through the endoscope versus suction is important for a proper endoscopic view. Excessive and badly timed flushing should be avoided in order not to flush away the stone.

The stone can only be fragmented if the ultrasound probe and/or the pneumatic probe is/are in direct contact with it. Before releasing a pulse or activating the ultrasound, touch the stone with the probe tip and move it slightly to verify contact.

Due to the push back effect, mobile ureteric stones should not be treated from below with the probe positioned in the center of the stone. Best practice is to position the probe tangentially at the stone so that the lithotripsy pulses will not push it back.

The stone should be hit at the corner or at the edges depending on its structure.

It is useful to treat the stone by pushing it slightly against the ureter wall and working with a low pulse rate. Use of a stone migration prevention device may be indicated.

It is possible to treat impacted stones directly "in situ" using only pneumatic lithotripsy, because pneumatic lithotripsy does not generate heat, therefore no risk of thermal damage to surrounding tissues.

Recommended settings

 For each lithotripsy method, whether used individually or in combination, it is recommended to start fragmentation at medium power settings of 50 % to 60 %.

The frequency of the pulses should not be higher than 7Hz in individual mode with a duty cycle set at 100 % at the beginning whether used alone or in combination with pneumatic lithotripsy, and can be increased or decreased gradually.

For combined application, it is recommended to set the pneumatic pulses frequency at 2~3Hz to help controlling the fragmentation.

 Working with a duty cycle of less than 100 % can be helpful to reduce warm-up of the ultrasound instrument if e.g. the endoscopic procedure requires to reduce suction through the ultrasound instrument.

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Selection between individual and combined modes

The Swiss LithoClast® Select System offers two modes of lithotripsy; ultrasound lithotripsy and pneumatic direct impact lithotripsy, known as the Swiss LithoClast® principle.

These two lithotripsy modes can operate in combination or individually. The selection of the treatment mode depends on the stone location and on the stone properties.

Stone fragmentation testing has been conducted to verify that simultaneous use of ultrasound and shock wave pneumatic lithotripsy can reduce the time needed for stone fragmentation and evacuation. In situations where it is time consuming and hard to fragment the stones by using ultrasound energy alone, the simultaneous application of the pneumatic probe guided through the ultrasound probe can considerably improve the process. This is particularly true when treating hard stones, such as mono hydrate stones.

The combined use is not indicated when the fragmentation process works fine with the ultrasound probe alone. The flow characteristic of the ultrasound probe is limited when the pneumatic probe is inside. Fragment evacuation is more efficient when the suction channel is kept free.

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WARNINGS

 **Danger of urothelial and mucosal perforation!** Due to the tendency of the ultrasound probe to perforate the stone rather than fragment it, there is a risk of injuring the mucosa of the kidney and bladder behind the stone. During the entire treatment, the user must keep the probe tips under endoscopic view (the probe tip is extended 10-20 mm beyond the endoscope tip).

 Do not apply pressure with the distal part of the probes against tissues as damage to tissue or perforation could occur. Use particular caution with thin probes.

 Do not touch the ultrasound probe when ultrasound is activated. The ultrasound vibrations of the probes may cause heat and instant burns.

 The user's hand may have a tendency to push forward under the influence of the repeated pneumatic pulses. This tendency must be corrected by the user who will continuously control the relative position of the probe to reduce risks of perforations.

 When using the unit in multiple pulse mode, be sure to check the progress of fragmentation and the probe position on a frequent basis (every 1 to 2 seconds).

Sterile backup ultrasound and pneumatic lithotripsy probes must be available in the event of a probe failure during treatment.

If a probe breaks distally, use sterile grasping forceps to remove probe pieces from the urinary tract.

POST-TREATMENT PROCEDURE

Emptying of suction tubes and Stone Catcher device

- Remove the silicone tubing of the Stone Catcher device from the pinch valve,
- Suck air in to empty the suction tubes,
- Tilt the Stone Catcher device in order to empty the receptacle by suction.

Elimination or conservation of Stone Catcher device content

- If fragments are not to be kept for analysis, dispose of the stone fragments per hospital protocol.
- If the fragments are to be kept, close the receptacle with the yellow transport closing cap.

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CLEANING, DISINFECTING AND STERILIZING

 Please conform to the recommendations of the "Reprocessing Instructions" manual delivered with your product regarding cleaning, disinfecting, packaging and sterilizing of the components.

 Please follow present day rules about sterilization applicable in your country.

Preparation and cleaning

Unit

 The housing of the unit is not watertight.

Please proceed as follows:

- Turn off the device
- Disconnect the power supply connector
- Disconnect the air supply connector
- Disconnect the foot pedal connector
- Disconnect the handpiece(s) and attach the caps onto the connectors

 Clean the unit only with an alcohol-based, commercially available (ethanol, isopropanol), colorless disinfectant. The use of scouring powder or an abrasive sponge will damage its surface.

Handpieces and accessories

 Clean the foot pedal only with an alcohol-based, commercially available (ethanol, isopropanol), colorless disinfectant. The use of scouring powder or an abrasive sponge will damage its surface.

 Dispose of the disposable Stone Catcher device.

 For the other components, please refer to related instruction manual.

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MAINTENANCE

Should legal provisions in your country demand maintenance intervals, these must be observed.

Display of the software version

When the device is turned on, the installed software version is displayed for few seconds.

Display of cumulated operation time

During power up, the total time of ultrasound usage as well as the total number of pneumatic impulses are displayed in the related displays on front panel.

Ultrasound time

The cumulative ultrasound usage time (seconds/hours) is shown briefly on the left display at start-up.

Remark: If the number is superior to 9999 (seconds), the hours mode is shown, for example: 012h = 12 hours.

If the number is superior to 999h (hours), the exponential mode is used: i.e. 45E3 = 45×10^3 = 45'000 hours

Pneumatic pulse counter

The cumulative number of pulses is shown briefly on the right display at start-up.

Remark: If the number is superior to 9999, the exponential mode is used
For example: 56E4 = 56×10^4 = 560'000

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
Fuses replacement

Disconnect the power supply cord at the back of the unit.

Replace the fuses with those specified on the identification plate at the back of the unit.

Test the control unit. If the fuse "burns" again, please contact your BSC sales representative.

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SAFETY PRECAUTIONS

Boston Scientific (distributor) or EMS (manufacturer) assumes no liability for direct or consequential injury or damage resulting from improper use, arising in particular through non-observance of the operating instructions, or improper preparation and maintenance.

 Read this instruction manual before using this device. This product may only be applied for its intended use by qualified personnel and for the applications described in this manual. If the unit is used in combination with other instruments, please refer to their instruction.

 Use for the intended purpose only. Before using the product, make sure the operating instructions have been carefully studied. This also applies to any equipment used with this product. Failure to observe the operating instructions may result in the patient or user suffering serious injury or the product being damaged.

 Risk of explosion: Do not use this product in the presence of inflammable anaesthetics or gases.

 This product must be used only by trained and qualified personnel.

 Always examine the product for damage before commencing treatment. Damaged accessories or a damaged unit must not be used and must be replaced. Use original Swiss LithoClast® Select System spare parts and accessories only.

 Proper care and maintenance are required to ensure safe operation of medical devices. Therefore it is recommended that devices and instruments be checked regularly for proper functioning and completeness, upon receipt and prior to each application, to prevent injury to patients and users.

 Product repairs are to be executed exclusively through Boston Scientific. After any repair or modification, the technical safety of the device and accessories must meet the corresponding standard applicable for such equipment. Each repair, modification or test must be recorded in a report.

 Before each use, always clean, disinfect and sterilize the different pieces and accessories of this product. Please refer to the information provided in the "Reprocessing Instructions" manual delivered with your product. Non-sterile pieces and accessories may cause bacterial or viral infections.

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STORING OF THE PRODUCT

 Keep the original packaging until the product is to be disposed of permanently. You can use it for shipping or storing at any time.

Should you wish to store the product for an extended period of time:

- Proceed as described in the chapter "Cleaning, disinfecting and sterilizing",
- Pack the product and all accessories in the original packaging.

Please refer to the storage and transport conditions in the "Technical data".

ACCESSORIES

Accessories are available from Boston Scientific Customer Service.

 For safety reasons this product may only be used with genuine Swiss LithoClast® Select System accessories. The user assumes all risks of using this product with non-genuine accessories or otherwise in non-compliance with the instructions manual.

PRODUCT WARRANTY

The warranty is valid for one year from the date of purchase.

Damages due to non-adherence to the operating instructions or wear out of parts are excluded from warranty.

 The warranty of the product will be cancelled if you try to open or disassemble it.

Warranty

The manufacturer warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond the manufacturer's and/or BSC's control directly affect the instrument and the results obtained from its use. The manufacturer's obligation under this warranty is limited to the repair or replacement of this instrument if deemed by the manufacturer to be defective in material or workmanship, and neither the manufacturer nor BSC shall be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this instrument. The user assumes responsibility for any consequences which may arise pertaining to the use of the products. BSC and the manufacturer neither assume, nor authorize any other person to assume for it, any other or additional liability or responsibility in connection with this instrument.

Swiss LithoClast® is a registered trademark.

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CUSTOMER SERVICE

Please contact your BSC representative for any servicing or repairs.

Boston Scientific and EMS decline responsibility for the safety of the device and declare the warranty null and void if service or repair is carried out by an unauthorized third party or if non-genuine spare parts are used.

It is best to ship your product in the original packaging. It protects your product against damage during shipment.

 Before dispatching your product including all accessories, please clean, disinfect and sterilize as described in the operating instructions.

 Risk of transmitting micro-organisms! To protect the personnel of the repair center and for safety reasons during transport and shipment, all products and accessories returned to the factory for repair or revision must be cleaned, disinfected and sterilized in accordance with the instructions manual. Repair can be refused for product or accessories received in a contaminated condition.

ELECTROMAGNETIC COMPATIBILITY

 This product needs special precautions regarding electromagnetic compatibility and needs to be installed and put into service according to the note "Electromagnetic compatibility".

 Certain types of mobile telecommunication equipment could potentially interfere with this product. The separation distances recommended in the note "Electromagnetic compatibility" must be taken into account.

 This product should not be used adjacent to or stacked with an other unit. If adjacent or stacked use is necessary, functioning should be verified prior to operation in the configuration in which it will be used.

 The use of accessories and cords other than those specified or sold by BSC as replacement parts, may result in increased emissions or decreased immunity of this product.

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SYMBOLS

EMS ⁺	Manufacturer's logo
	Caution!
CE ₀₁₂₄	CE symbol refers to directive 93/42/EEC, including EN 60601-1 and EN 60601-1-2
	Power switch "ON"
○	Power switch "OFF"
	Equipotential plug
	Foot pedal connection
IP..	Degree of protection against water permeability
	Disposal of Old Electrical & Electronic Equipment (Applicable in the European Union and other European countries with separate collection systems)

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

SYMBOLS

	Standby / ON
Litho Mode	Switch for toggling between single and multiple pulse modes for the Swiss LithoClast pneumatic module
	Single pulse mode
	Multiple pulse mode
	Pneumatic handpiece connection
	Applied part, type BF
	Fuse
	Input
	Date of manufacture

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

SYMBOLS

	CSA marking with "C" identification for products in conformance with Canadian standards and "US" for products in conformance with US standards
	Ultrasound handpiece connection
Select	Switch for toggling between ultrasound and pneumatic modules in order to change settings
	Setting increase
	Setting decrease
	Rotation direction indicator to increase the ultrasound and pneumatic power
	Grounding
	Read the operating instructions
	GOST R marking for products in conformance with Russian standards
	DEKRA INMETRO identification for products in conformance with Brazilian electrical standards

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TECHNICAL DATA

DESCRIPTION

Distributor	Boston Scientific Corporation Marlborough, MA USA 01752
Manufacturer	EMS SA, CH – 1260 Nyon, Switzerland
Model	Swiss LithoClast® Select System
Power supply (voltage)	100 – 240 VAC
Power supply (frequency)	50 – 60 Hz
Power consumption	500 VA
Ultrasound frequency	23 – 26.4 kHz
Ultrasound power consumption	78 W
Energy output	Appro. 85 mJ
Classification	EN-60601-1: Class I Applied part BF IP 20 (Foot pedal: IP X8) MDD 93/42 EEC: Class IIb
Operating mode	Continuous modes
Primary fuse	Ø 5 X 20 mm T4A 250 V, time lag, low breaking capacity
Compressed air supply	3.5 to 6.5 bar (3500 to 6500 hPa) (51 psi to 94 psi)
Weight	11 kg (24.3 lb)
Housing dimensions	Height: 135 mm (5.3 in) Width: 360 mm (14.2 in) Depth: 420 mm (16.5 in)
Operating conditions	Temp. +10°C to +30°C (50°F to 86°F) Rel. humidity 30% to 75% Atm. pressure 700 hPa to 1060 hPa (10.1 psi to 15.3 psi), Altitude max. 3000 m (9800 ft)

Storage and transport conditions

-10°C to +40°C (14°F to 104°F)
10% to 95% relative humidity
500 hPa to 1060 hPa air pressure
(7.2 psi to 15.3 psi)

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TROUBLESHOOTING

Error codes

The Swiss LithoClast® Select System provides an error code system in case of malfunction or incorrect handling. The code is made of two blocks of characters like **Er 12**. On the display, the code is on the left for ultrasound, on the right for pneumatics, and is split left-right for general. Codes of the A type disappear with the cause. The B type code must be confirmed following solution by pressing and releasing the foot pedal. The ultrasound duty cycle must be at 100% for the error message to be detected.

No.	type	linked to	cause/description	solution
Er 10	A	push button or pedal	Defective foot pedal or foot pedal accidentally pressed while switching the unit on.	Pay attention not to press a foot pedal when switching on the unit. Check that nothing pushes on or impair the normal functioning of the pedal. If failure persists call BSC customer service.
Er 11			Defective push button or push button accidentally pressed while switching the unit on.	Pay attention not to press a push button when switching on the unit. If failure persists call BSC customer service.
Er 20	B	Ultrasound Handpiece	The internal ultrasound power value cannot be obtained. Defective ultrasonic probe or handpiece.	Please proceed as follows: 1 Verify that the ultrasound probe is tightened, 2 Press the foot pedal to try again, 3 If the error appears again, replace the probe, 4 Press the foot pedal to try again, 5 If the error appears again, reduce the ultrasound power to 80 %, 6 Press the foot pedal to try again, 7 If the error persists, Replace the ultrasound handpiece, 8 Press the foot pedal to try again, 9 If the error is still present, call BSC customer service.
Er 21			Defective ultrasound handpiece or cord (i.e. a damaged insulation).	
Er 23			Loss of ultrasound effectiveness, the control unit cannot adjust the resonance frequency. The handpiece is loose or not connected, the ultrasound probe is damaged or untightened.	
Er 24			Ultrasound efficiency is not stable. The resonance frequency is lost (similar to code Er 23 but this code does not come up because the ultrasound functions intermittently).	
Er 30	B	Ultrasound mode	Internal power supply failure.	Turn the device off, then on again. If failure persists call BSC customer service.
Er 40	A	Pneumatic mode	Internal working pressure is not in the specified range. Input pressure can be too high or too low. Activation of the pneumatic pulses is stopped.	1 Verify and adjust the air pressure according to the value specified on the identification plate on the unit back panel. 2 Reduce working pressure using the pneumatic rotary knob. Press the foot pedal to try again. 3 Do not increase pressure beyond 100 % again. If the error is still present, call BSC customer service.

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