
Swiss LithoClast® Select System

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
DIRECTIONS FOR USE
MODE D'EMPLOI

© Copyright EMS SA
FB-389/US ed. 2015/06
(part of BSC ref. 840-302)

Rx ONLY Caution! Federal (USA) law restricts this device
to sale by or on the order of a physician

**Boston
Scientific**

COMPONENTS

- 1 Ultrasound Handpiece kit
- 2 Ultrasound Handpiece
- 3 5 mm wrench
8 mm wrench
- 4 Angled suction connector for combination use ultrasound /
pneumatic
- 5 Straight suction connector for ultrasound only use
- 6 Electrical connector
- 7 Protective sterilization cap for the electrical connector
- 8 Protective cap for the Ultrasound Handpiece
- 9 Ultrasound Probe

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Technical modifications

BSC reserves the right to modify the technique, accessories, operating instructions or contents of the set due to technical or scientific improvements.

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DEAR CUSTOMERS,

This Vario Ultrasound Handpiece is used with the Swiss LithoClast® Select System for the fragmentation of urinary tract calculi in the kidney, ureter, and bladder.

The Vario Ultrasound Handpiece is compatible with Swiss LithoClast® Select System control units equipped with software version "rel. 3.x". Control units having a previous version of software installed require an upgrade to enable controlling the Vario Ultrasound Handpiece.

Contraindications

Use of the product is contra-indicated in patients with the following:

- Active bleeding disorders
- Solitary functioning kidney
- Creatinine greater than or equal to 3 µg %
- During pregnancy
- Under the age of 18
- Stricture and obstruction problems
- With an electrical stimulator implanted e.g. pacemaker

Instructions concerning the specific components are provided separately in the relative Direction for use. Read the Directions for use delivered with each of the components before using the Swiss LithoClast® Select System. Do not hesitate to contact your BSC representative for further information.

PLEASE READ THIS BEFORE BEGINNING WORK!

About this manual

Note that the English version of this manual is the master from which translations derive. In case of any discrepancy, the binding version is the English text.

These operating instructions are to ensure the correct installation and use of this product.

Please read these operating instructions carefully as they explain all of the most important details and procedures. Please pay special attention to the safety precautions.

Always keep this instruction close at hand.

To prevent possible patient or user injury or damage to the Swiss LithoClast® Select System, please note the corresponding directives.

Danger

Risk of patient or user injury

Caution

Risk of damage to the device or environmental harm

Please note

Useful additional information and hints

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Autorisé

CONTENTS

Preparing for treatment.....	2
Assembly of the Vario Ultrasound Handpiece.....	2
Connecting the Vario Ultrasound Handpiece	4
Ultrasound Lithotripsy Mode.....	6
Combined ultrasound and pneumatic lithotripsy	8
Connecting the combined handpieces to the unit	12
Preparation.....	14
Application	16
Selection of lithotripsy mode.....	18
Change from combined ultrasound and pneumatic lithotripsy to ultrasound mode only	20
Operating the handpiece	22
Cleaning, disinfecting and sterilizing	24
Preparation of the Ultrasound Handpiece for reprocessing	24
Preparation of the combination handpiece for reprocessing	26
Maintenance.....	32
O-ring	32
Safety precautions	34
Storing the product	36
Accessories.....	36
Product warranty.....	36
Warranty	36
Symbols	38
Customer Service	40
Troubleshooting	40
Technical data.....	40

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
PREPARING FOR TREATMENT

Assembly of the Vario Ultrasound Handpiece

 The equipment has been provided non-sterile and must be cleaned and sterilized before use, please refer to the "Reprocessing instructions" manual.

 The Ultrasound Probe must be tightened onto the Ultrasound Handpiece using the supplied wrench. Keep the wrench sterile. If the probe loosens during treatment, retighten the probe with the wrench. Even a slightly loose probe connection will cause a total loss of fragmentation performance.


 To thread the suction connector in properly, align it carefully with the handpiece and apply a light axial pressure while screwing it in place.

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
Connecting the Vario Ultrasound Handpiece

The single use Stone Catcher device is provided sterile. Read the instruction manual delivered with it before performing installation.


 Connect the long tube of the sterile single use Stone Catcher device to the suction connector of the Ultrasound Handpiece.

 The electrical connector must be perfectly dry before connecting it to the unit.

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Ultrasound Lithotripsy Mode

Preparation

All components (including wrenches) are provided non sterile, therefore they must be sterilized and dried per the recommendations prior to assembly. In a sterile manner, attach the Ultrasound Probe to the front end of the handpiece with clockwise turns until finger tight, then securely fasten with the 5mm wrench.

In a sterile manner, attach the straight suction connector to the back end of the handpiece with slight axial pressure and a ¼ clockwise turn. Once complete, connect the long end of the tubing from the disposable Stone Catcher device (packaged sterile) to the straight suction connector.

Remove the protective sterilization cap and connect the plug of the Vario Ultrasound Handpiece to the control unit by aligning the red dots. The connector must be completely dry before attachment.

Hand the Stone Catcher device off the sterile field and place the short portion of tubing from the Stone Catcher device into the pinch valve on the control unit. Also connect the roller clamp to the short tubing to help regulate the flow of suction. Once completed, place the canister in the Stone Catcher device Holder.

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
Combined ultrasound and pneumatic lithotripsy

Assembly of the Vario Ultrasound Handpiece for combined use with the Pneumatic Handpiece pn3

 The equipment has been provided non-sterile and must be cleaned and sterilized before use, please refer to the "Reprocessing instructions" manual.

 To thread the suction connector in properly, align it carefully with the handpiece and apply a light axial pressure while screwing it in place.


 The Ultrasound Probe must be tightened onto the Ultrasound Handpiece by using the supplied wrench. Store the wrench sterile. If the probe loosens during treatment, retighten the probe by using the wrench. Even a slightly loose probe connection will cause a total loss of fragmentation performance.

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

Adjust the length of the Pneumatic Probe by turning the wheel of the adjustment interface so that the tip of the Pneumatic Probe is flush with the tip of the Ultrasound Probe.

 If the Pneumatic Probe cannot be adjusted flush with the Ultrasound Probe, this indicates the wrong ultrasound / Pneumatic Probe combination is being used. Combination probes are marked by an identical ring pattern on the proximal shaft of both probes.

 Both lithotripters work only in direct contact with calculi. The length adjustment of the Pneumatic Probe relative to the Ultrasound Probe is important, when using both lithotripters simultaneously.

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
Connecting the combined handpieces to the unit

The single use Stone Catcher device is provided sterile. Read the instruction manual delivered with it before performing installation.



 Connect the long tube of the sterile single use Stone Catcher device to the suction connector of the Ultrasound Handpiece.

 Both connectors must be perfectly dry before connecting them to the unit.

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PREPARATION

The Vario Handpiece, Pneumatic Handpiece, Ultrasound Probe wrenches, the suction connectors, and Ultrasound Probes are provided non sterile, therefore they must be sterilized per the recommendations prior to assembly. The Pneumatic Probes and Stone Catcher devices are single use and provided sterile.

In a sterile manner, attach the Ultrasound Probe to the front end of the handpiece with clockwise turns until finger tight, then securely fasten with the 5 mm wrench. Wrenches should be sterilized with handpieces per sterility guidelines. Attach the angled suction connector to the back end of the handpiece with slight axial pressure and a ¼ clockwise turn.

Once completed, insert the Pneumatic Probe down the center lumen of the angled suction adapter, through the Vario Ultrasound Handpiece and Ultrasound Probe, then screw the Pneumatic Probe with clockwise turns into the angled suction adapter until finger tight.

Attach the pneumatic handpiece to the Pneumatic Probe with slight axial pressure and a ¼ clockwise turn until the probe 'clicks' into place. Once attached, adjust the white dial on the Pneumatic Probe to ensure the Pneumatic Probe is flush or even with the Ultrasound Probe. Clockwise turns move the Pneumatic Probe out and counter clockwise turn pull the Pneumatic Probe in.

Once the probe has been adjusted attach the long end of the tubing from the Stone Catcher device to the Vario Handpiece and then hand the Stone Catcher device off the sterile field. Place the short portion of tubing from the Stone Catcher device into the pinch valve on the control unit. Also connect the roller clamp to the short tubing to regulate the flow of suction. Once completed, place the canister in the Stone Catcher device Holder.

Remove the protective sterilization caps. Connect the plug of the Vario Ultrasound Handpiece to the control unit by aligning the red dots and connect the pneumatic handpiece by just pushing the connector into the control unit. Both connectors must be completely dry before attachment.

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APPLICATION

Coding system for the ultrasound and pneumatic probe sets

	Set 330 mm serie	Set 403 mm serie
Ultrasound Probes	 Ø 3.8 mm  Ø 3.3 mm	 Ø 3.8 mm  Ø 3.3 mm
Pneumatic Probes	 Ø 1.0 mm	 Ø 1.0 mm

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Selection of lithotripsy mode

The following table shows which Ultrasound Probe size and Pneumatic combination Probe size is used for various stone types and locations:

Lithotripsy Mode

Stone Type	Ultrasound only	Pneumatic & Ultrasound combination
Hard kidney stones through percutaneous access (monohydrate, cystine stones)	Ultrasound probe: 3.3 mm or 3.8 mm Use ultrasound lithotripsy only for hard stones after initial fragmentation with the pneumatic / ultrasound combination (see directions for intra-operative change between combination hand-piece setup and ultrasound only)	Ultrasound probe: 3.3 mm or 3.8 mm Pneumatic probe: 1.0 mm Pneumatic probe 1.3 mm (optional) for use with 3.8 mm Ultrasound Probe only
Soft kidney stones (struvite, matrix stones)	Ultrasound probe: 3.3 mm or 3.8 mm	
Ureteral stones	Ultrasound probe: 1.5 mm	
Hard bladder stones (monohydrate stones)	Ultrasound probe: 3.3 mm or 3.8 mm This method can be used for fragment removal following the pneumatic fragmentation	Ultrasound probe: 3.3 mm or 3.8 mm Pneumatic probe: 1.0 mm Pneumatic probe 1.3 mm (optional) for use with 3.8 mm Ultrasound Probe only
Soft bladder stones	Ultrasound probe: 3.3 mm or 3.8 mm	
Small kidney stones	Ultrasound probe 1.9 mm for use with mini nephroscopes with 6 Fr. work channel size	

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
Change from combined ultrasound and pneumatic lithotripsy to ultrasound mode only

Disconnect the Pneumatic Handpiece from the control unit and attach the sterilization cap to the end of the hose.

Disconnect the Pneumatic Probe from the Pneumatic Handpiece by applying slight axial pressure and turning the handpiece $\frac{1}{4}$ turn counter clockwise.

Remove the angled suction connector and Pneumatic Handpiece from the Ultrasound Handpiece.

Attach the straight connector to the Ultrasound Handpiece and attach the Stone Catcher device tubing to the straight connector.


 It is necessary to change the angled suction connector with the straight suction connector.

 Always use the straight suction connector when working in ultrasound mode only.

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OPERATING THE HANDPIECE

 Do not touch the Ultrasound Probe when the ultrasound is activated. The ultrasound vibrations in the probes may cause instant burns by frictional heating.

 Danger of urothelial and mucosal perforation! There is a risk of injuring the mucosa of the kidney and bladder, which lies behind the stone. Always work under endoscopic visualization and keep the Pneumatic Probe tip flush with the Ultrasound Probe tip when using the combined modes.

 Fragments can block suction through the Ultrasound Probe. This will lead to a loss of the ultrasound suction capabilities and to heat-up of Ultrasound Probe and handpiece. Check at regular intervals that the suction is functional.

 The Ultrasound Handpiece and probe error detection is only possible when working with 100% "Duty cycle".

 Check that the suction connector is tight on the handpiece to maintain a good seal on the suction system.

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
CLEANING, DISINFECTING AND STERILIZING

Preparation of the Ultrasound Handpiece for reprocessing

 Please conform to the recommendations of the Reprocessing Instructions manual delivered with your product regarding procedure of cleaning, disinfecting, packaging and sterilizing of the components.

 Please follow present day rules about sterilization applicable in your country.

 The handpiece must be fully disassembled for each sterilization.

 Ensure sterilization cap is connected to each hose before processing.



 If the Ultrasound Probe is hard to untighten from the handpiece, use the 5 mm wrench and the 8 mm wrench as shown to unlock.

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
Preparation of the combination handpiece for reprocessing


 Please conform to the recommendations of the "Reprocessing Instructions" manual delivered with your product regarding procedure of cleaning, disinfecting, packaging and sterilizing of the components.

 Please follow present day rules about sterilization applicable in your country.

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



The Ultrasound Handpiece must be fully disassembled for each sterilization.

If the probe has become too tight during the course of the procedure then disconnect the distal protective cap by screwing counter clockwise. Then use the 8 mm wrench to provide counter tension to allow the probe to be loosened and removed.

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Please proceed as follows for proper disassembly of the handpieces.

Ultrasound Lithotripsy Mode

At the end of the procedure disconnect the Vario Ultrasound Handpiece from the control unit and attach the sterilization cap to the end of the hose. Also remove the suction tubing from the pinch valve, disconnect the suction tubing from the Stone Catcher device, and disconnect the Stone Catcher device tubing from the Vario handpiece.

Remove the Stone Catcher device cap with the tubing and replace with the Yellow protective cap so stones can be sent for analysis. Dispose of the remaining Stone Catcher device parts per hospital protocol.

Disconnect the Ultrasound Probe from the Vario Handpiece using the 5mm wrench. If the probe has become too tight during the course of the procedure then disconnect the distal protective cap by screwing counter clockwise. Then use the 8 mm wrench to provide counter tension to allow the probe to be loosened and removed. Disconnect the straight suction connector from the back end of the Vario Handpiece.

Always ensure the protective sterilization cap is on the end of the hose before reprocessing.

Pneumatic Lithotripsy Mode

At the end of the procedure disconnect the Pneumatic Handpiece from the control unit and attach the sterilization cap to the end of the hose. Disconnect the Pneumatic Probe from the Pneumatic Handpiece by applying slight axial pressure and turning the handpiece ¼ turn counter clockwise. Dispose of the probes after use.

 Never disconnect the air supply tubing from the Pneumatic Handpiece for reprocessing. Always ensure the protective sterilization cap is on the end of the hose before reprocessing.

Combination Ultrasound and Pneumatic Lithotripsy Mode

At the end of the procedure disconnect both the Pneumatic Handpiece and the Vario Handpiece from the control unit and attach the sterilization cap to the end of each hose.

Disconnect the Stone Catcher device tubing from the Vario handpiece. Remove the Stone Catcher device cap with the tubing and replace with the Yellow protective cap so stones can be sent for analysis. Disconnect the Pneumatic Probe from the handpiece by applying slight axial pressure and turning the handpiece ¼ turn counter clockwise. Dispose of the remaining Stone Catcher device parts per hospital protocol.

Unscrew the Pneumatic Probe from the Vario Handpiece and dispose of. Do not disconnect the air supply tubing from the Pneumatic Handpiece for reprocessing.

Disconnect the Ultrasound Probe from the Vario Handpiece using the 5mm wrench. If the probe has become too tight during the course of the procedure then disconnect the distal protective cap by screwing counter clockwise. Then use the 8 mm wrench to provide counter tension to allow the probe to be loosened and removed. Disconnect the angled suction connector from the back end of the Vario Handpiece.

 Never disconnect the air supply tubing from the Pneumatic Handpiece for reprocessing. Always ensure the protective sterilization cap is on the end of the hose before reprocessing. Handpiece damage will occur if moisture enters the handpiece interior.

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
MAINTENANCE

O-ring

 If the O-ring is damaged or worn out, it must be replaced immediately.


Change the O-ring or replace the suction connector.

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SAFETY PRECAUTIONS

Boston Scientific (distributor) or EMS (manufacturer) assumes no liability for direct or consequential injury or damage resulting from improper use, arising in particular through non-observance of the operating instructions, or improper preparation and maintenance.

 Read this instruction manual before using this product. This product may only be applied for its intended use by qualified personnel and for the applications described in this manual. If the unit is used in combination with other instruments, please refer to their instruction.

 Use for the intended purpose only. Before using the product, make sure the operating instructions have been carefully studied. This also applies to any equipment used with this product. Failure to observe the operating instructions may result in the patient or user suffering serious injury or the product being damaged.

 This product must be used only by trained and qualified personnel.

 Always examine the product for damage before commencing treatment. Damaged accessories or a damaged unit must not be used and must be replaced. Use original Swiss LithoClast® Select System spare parts and accessories only.

 Proper care and maintenance are required to ensure safe operation of medical devices. Therefore it is recommended that devices and instruments be checked regularly for proper functioning and completeness, upon receipt and prior to each application, to prevent injury to patients and users.

 Product repairs are to be executed exclusively through Boston Scientific. After any repair or modification, the technical safety of the device and accessories must meet the corresponding standard applicable for such equipment. Each repair, modification or test must be recorded in a report.

 Before each use, always clean, disinfect and sterilize the different pieces and accessories of this product. Please refer to the information provided in the "Reprocessing Instructions" manual delivered with the product. Non-sterile pieces and accessories may cause bacterial or viral infections.

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STORING THE PRODUCT

 Keep the original packaging until the product is to be disposed of permanently. You can use it for shipping or storing at any time.

Should you wish to store the product for an extended period of time:

- Proceed as described in the chapter "Cleaning, disinfecting and sterilizing",
- Pack the product and all accessories in the original packaging.

Please refer to the storage and transport conditions in the "Technical data".

ACCESSORIES

Accessories are available from Boston Scientific Customer Service.

 For safety reasons this product may only be used with genuine Swiss LithoClast® Select System accessories. The user assumes all risks of using this product with non-genuine accessories or otherwise in non-compliance with the instructions manual.

PRODUCT WARRANTY

The warranty is valid for one year from the date of purchase.

Damages due to non-adherence to the operating instructions or wear out of parts are excluded from warranty.

 The warranty of the product will be cancelled if you try to open or disassemble it.

Warranty

The manufacturer warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond the manufacturer's and/or BSC's control directly affect the instrument and the results obtained from its use. The manufacturer's obligation under this warranty is limited to the repair or replacement of this instrument if deemed by the manufacturer to be defective in material or workmanship, and neither the manufacturer nor BSC shall be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this instrument. The user assumes responsibility for any consequences which may arise pertaining to the use of the products. BSC and the manufacturer neither assume, nor authorize any other person to assume for it, any other or additional liability or responsibility in connection with this instrument.

Swiss LithoClast® is a registered trademark.

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SYMBOLS

EMS⁺	Manufacturer's logo
	Caution! Read the operating instructions
CE 0124	EC marking: Refers to directive 93/42 EEC, including EN 60601-1 and EN 60601-1-2
	Sterilizable at up to 135°C in the autoclave
	For single use only. Do not reuse.
LOT	Lot number
 	Disposal of Old Electrical & Electronic Equipment (Applicable in the European Union and other European countries with separate collection systems)

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CUSTOMER SERVICE

Please contact your BSC representative for any servicing or repairs.

Boston Scientific and EMS decline responsibility for the safety of the device and declare the warranty null and void if service or repair is carried out by an unauthorized third party or if non-genuine spare parts are used.

It is best to ship your product in the original packaging. It protects your product against damage during shipment.

 Before dispatching your product including all accessories, please clean, disinfect and sterilize as described in the operating instructions.

 Risk of transmitting micro-organisms! To protect the personnel of the repair center and for safety reasons during transport and shipment, all products and accessories returned to the factory for repair or revision must be cleaned, disinfected and sterilized in accordance with the instructions manual. Repair can be refused for product or accessories received in a contaminated condition.

TROUBLESHOOTING

 For all codes displayed on the unit may, please read the instruction manual of the unit.

TECHNICAL DATA

Description

Distributor	Boston Scientific Corporation Marlborough, MA USA 01752
Manufacturer	EMS SA, CH-1260 Nyon, Switzerland
Model	Vario Ultrasound Handpiece for use with the Swiss LithoClas [®] Select System
Ultrasound frequency	23.2 – 26.4 kHz
Classification EN-60601-1	Class I Applied part BF IP67 Mode: Continuous operation
Classification MDD 93/42 EEC	Class IIb
Weight	270 g (9.5 oz)
Dimensions	Length : 132-138 mm (5.2 - 5.4 in) Width : 31.5 mm (1.2 in)
Operating conditions	+10°C to +40°C (50°F to 104°F) 30% to 75% relative humidity
Validated life time	100 treatments
Storage and transport conditions	-10°C to +40°C (14°F to 104°F) 10% to 95% relative humidity 500 hPa to 1060 hPa air pressure (7.2 psi to 15.3 psi)

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