

Capio™ SLIM Suture Capturing Device

Order Number	Description	Unit
M0068318250	Capio™ SLIM Suture Capturing Device	bx 1
M0068318261	Capio™ SLIM Suture Capturing Device	bx 5

Capio Sutures

Order Number	*New Order No.	Description	Size	Unit
M0068331131	M0068332231	Non-absorbable, coated braided polyester, double armed with TC tapercut needle (dart), 48"	0	bx 12
M0068331231	M0068332241	Non-absorbable, polypropylene mono ilament, double armed with TC tapercut needle (dart), 48"	0	bx 12
M0068332131	M0068332341	Absorbable, coated braided PGA, double armed with TC tapercut needle (dart), 48"	0	bx 12
M0068331141	M0068332351	Non-absorbable, coated braided polyester, double armed with TC tapercut needle (dart) and a T 26mm 1/2 circle taper needle, 36"	0	bx 12
M0068331241	M0068332451	Non-absorbable, polypropylene mono ilament, double armed with TC tapercut needle (dart) and a T 26mm 1/2 circle taper needle, 36"	0	bx 12
M0068331371	Unchanged	Monodek™ Absorbable, mono ilament PDO, double armed with TC tapercut needle (dart) and a T 26mm 1/2 circle taper needle, 48"	0	bx 12

*Existing Order Numbers will be replaced by new ones during 2021

Your pelvic loor reconstructive tool to meet a variety of your algorithm needs:

- Native tissue repair
- Biologic graft augmentation
- And more...

Capio is a registered trademark of Boston Scientific Corporation or its affiliates. All other trademarks are the property of their respective owners. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labelling supplied with Rx Only. Information for the use only in countries with applicable health authority product registrations. Data on file. Bench test results may not necessarily be indicative of clinical performance. Results from case studies are not predictive of results in other cases. Results in other cases may vary. Not intended for use or distribution in France.

WH-653602-AA Printed in Germany by medicalvision.

**Boston
Scientific**
Advancing science for life™

www.bostonscientific.eu

© 2019 Boston Scientific Corporation or its affiliates. All rights reserved. DINURO2412EA

**Boston
Scientific**
Advancing science for life™




Capio™ SLIM
Suture Capturing Device

Capio™ SLIM

Suture Capturing Device

The Capio™ SLIM Suture Capturing device features a reduced profile to minimize the space required within the surgical field, easy suture loading and a funnel shaped catch, designed for consistent ease of use.



Comparison of Suture Capturing Device Heads

Building on the history and reducing the profile, the Capio™ SLIM Suture Capturing Device features one of the smallest device profiles and suture dart carrier diameter while keeping the same bite depth as the current Capio Open Access Suturing Device.

Capio™ SLIM
Suture Capturing Device



Digitex™
Suture Delivery System



FiXT™
Suturing Device



Figure 1: Comparison of Device Heads illustrates the difference in carrier diameter amongst these fixation devices. The Capio™ SLIM Suture Capturing Device carrier diameter is 45% less than the FiXT Suturing Device, 43% less than the Digitex Suture Delivery System.

Fixation Device Comparison

	Shaft Diameter	Head Width	Device Weight	Carrier Diameter
Capio OPEN ACCESS Suture Capturing Device**	10.1 mm	9.8 mm	45.4 g	1.2 mm
Capio™ SLIM Suture Capturing Device*	3.0 mm	6.3 mm	38.2 g	1.2 mm
Digitex Suture Delivery System*	9.5 mm	9.6 mm	133.3 g	2.1 mm
FiXT Suturing Device*	4.8 mm	6.9 mm	121.1 g	2.2 mm

Figure 2: Fixation Device Comparison illustrates the device profile differences amongst these fixation devices. The Capio™ SLIM Suture Capturing Device has the smallest device profile based on the following measurements: shaft diameter, head width and device weight.

* Measurements were recorded using one (1) of each device. Data on file at Boston Scientific.

** Head width and device weight values are an average of five (5) device measurement recordings. Shaft diameter values are an average of fifteen (15) device measurement recordings. Data on file at Boston Scientific.

The evolution continues with the Capio™ SLIM Suture Capturing Device

Since 1995, the Capio Device has provided consistent suture placement in difficult to access pelvic floor locations.



1995

The Capio Suture Capturing Device debuts, revolutionizing the field of urogynecology.



1999

Capio Device may be used with other biological graft materials offering another treatment option



2011

The Capio Device remains “the foundation” of Pelvic Floor Reconstruction.



2012

Capio™ SLIM: the evolution continues.