


AGENT™

Drug-Coated Balloon



1. DCB Bench testing performed by Boston Scientific Corporation. Data on file at Boston Scientific. Tip profile data on file for 3.00 mm DCBs. Flexibility comparison data on file for 3.50 mm DCBs excluding Selution, data on file for 3.00 mm DCB.
2. Pre-clinical pharmacokinetic studies, data on file at Boston Scientific for AGENT, SeQuent Please and Selution. Magic Touch data derived from product brochure (TCT 2020).
3. Blagosklonny et al. Paclitaxel induces primary and postmitotic G1 arrest in human arterial smooth muscle cells. Cell Cycle. 2004 Aug;3(8):1050-6. Epub 2004 Aug 24.
4. Shahrzori ZMF, et al. Drug-Coated Balloons: Recent Evidence and Upcoming Novelties. J Cardiovasc Dev Dis. 2025 May 20;12(5):194. doi: 10.3390/jcdd12050194. PMID: 40422965; PMCID: PMC12112549. Other coronary PTX DCB drug concentrations; SeQuent Please 3 µg/mm², Prevail 3.5 µg/mm², Pantera Lux 3 µg/mm². Also see respective product IFU's.
5. Market data on file at BSC as of July 2023.
6. Hamm, C, et al, AGENT ISR study, Eurointervention, 2020, 3-year comparative clinical study results on file at Boston Scientific Corporation: Data available upon request.
7. Nakamura, M, et al, Drug-coated balloon for the treatment of small vessel coronary artery disease, Circulation Japan 2022.
8. Yeh RW, Shlofmitz R, Moses J, et al. Paclitaxel-Coated Balloon vs Uncoated Balloon for Coronary In-Stent Restenosis: The AGENT IDE Randomized Clinical Trial. JAMA. 2024;331(12):1015–1024. doi:10.1001/jama.2024.1361


Please review the full IFU for indications for use and operating instructions.

CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings, and instructions for use can be found in the product labelling supplied with each device or at www.IFU-BSCI.com. Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. This material not intended for use in France.



Leave the Right AGENT Behind

► Perform complex interventions
without compromise



Effortless Deliverability
Navigate and cross complex lesions with ease




Efficient Drug Transfer
Deliver the optimal dose when it matters most,
with less balloon drug load



**Exceptional Long-Term
Outcomes**
Drive lower event rates consistently as demonstrated
by AGENT clinical trials


➤ Effortless Deliverability

Navigate and cross complex lesions with ease



➤ **Laser Bonded Tip**
Improves crossability and reduces tip catch


➤ **Bi-segment Inner Shaft**
Maximizes flexibility without compromising push for reduced resistance in tortuous anatomy



➤ **Z-Glide™**
Improves tracking with proprietary lubricious coating


➤ The AGENT Difference

Ultra-Low Tip Profile¹



LOWER IS BETTER

Uncompromised Flexibility¹



Distal End Flexibility


Overcoming the rigidity typically introduced by drugcoating matrices (drug + excipient/carrier) is key to optimal distal end flexibility. AGENT DCB demonstrated best-in-class distal end flexibility, minimizing resistance and enhancing deliverability in tortuous anatomy.

LOWER IS BETTER

[†] Measure of how much work is required to deliver the distal end of the device through a curved simulated artery – how easily does the distal assembly flex.

Efficient Drug Transfer

Deliver the optimal dose when it matters most, with less balloon drug load. Due to the efficiency of AGENT™ Drug-Coated Balloon’s proprietary TransPax™ Coating Technology, AGENT is able to use the lowest paclitaxel dosage compared to other paclitaxel-coated balloons, while still resulting in the highest concentration of drug in the tissue.



TransPax™ coating optimizes the three phases of drug delivery

AGENT DCB is designed with a novel excipient, sharp-edge structure, and uniform crystalline formulation. As a result, AGENT transfers more drug to the tissue and less of it downstream.

Design Requirement

1 Targeted Transfer
Balloon to vessel wall

Paclitaxel (PTX) Drug Choice

Hydrophobic
Paclitaxel is durable under hydration

AGENT Coating Design

Optimized Excipient
Acetyl tributyl citrate (ATBC) maximizes balloon-to-vessel wall transfer

CCCC(=O)OC(C)(OC(=O)CCCC)OC(=O)CCCC

Design Requirement

2 Rapid Absorption
Vessel wall to tissue

Paclitaxel (PTX) Drug Choice

Lipophilic
Paclitaxel has a high affinity for fatty tissue

AGENT Coating Design

Sharp-Edge Structure
Needle-like coating improves tissue penetration

Design Requirement

3 Sustained Retention
Prolonged tissue concentration

Paclitaxel (PTX) Drug Choice

Chemically Stable
Paclitaxel has a long half-life

AGENT Coating Design

Crystalline Formulation
Crystallized coating maintains therapeutic effect


²Independent studies performed using the same FDA recommended and accepted porcine coronary de novo vessel test method.

➤ Exceptional Long-Term Outcomes

Drive lower event rates consistently as demonstrated by AGENT clinical trials. **Over 4,400 patients** have been evaluated or are undergoing evaluation with AGENT DCB.⁵ This includes studies in ISR, de novo lesions, small vessels, large vessels, bifurcations, complex lesions, and patient-specific cases.



AGENT ISR Study⁶

AGENT demonstrated comparable event rates at 3-years compared to SeQuent Please



AGENT Japan SV Study⁷

AGENT demonstrated numerically lower levels of late lumen loss compared to SeQuent Please




➤ Proven clinical performance

AGENT US IDE Study⁸

This prospective, multicenter, 2:1 randomized controlled trial showed statistically lower event rates for AGENT DCB at one year.

41%
relative risk
reduction
for TLF

NNT = 9 Patients*



➤ At one-year, AGENT DCB also demonstrated statistically lower event rates:

50%
risk reduction
for TLR

(13.0% vs. 24.7%, P=0.0005)

49%
risk reduction
for TV-MI

(5.8% vs. 11.1%, P=0.023)

0%
ZERO def/prob ST

(0.0% vs. 3.2%, P=0.0004)

*Number Needed to Treat to prevent 1 TLF

➤ Modern PCI Portfolio


With the most comprehensive Modern PCI portfolio on the market, Boston Scientific has the tools needed to See, Prep, and Treat each lesion.




AVVIGO⁺
Multi-Modality Guidance System



OPTICROSS^{HD}
60 MHz Coronary Imaging Catheter



COMET^{II}
Pressure Guidewire



WOLVERINE[™]
Cutting Balloon[™] Dilatation Device



ROTAPRO[™]
Rotational Atherectomy System




GUIDEZILLA[™]
Guide Extension Catheter




AGENT[™]
Drug-Coated Balloon




SYNERGY^{XD}
Everolimus-Eluting Platinum Chromium
Coronary Stent System



SYNERGY^{MEGATRON}
Everolimus-Eluting Platinum Chromium
Coronary Stent System



AGENT[™]
Drug-Coated Balloon



One platform. Expanded potential.
AGENT DCB's 40 mm balloon is here.



Drug-Coated Balloon Length (mm)

Balloon Length Balloon Diameter	12mm	15mm	20mm	30mm	40mm
2.00	H74939222201210	H74939222201510	H74939222202010	H74939222203010	H74939222204010
2.25	H74939222221210	H74939222221510	H74939222222010	H74939222223010	-
2.50	H74939222251210	H74939222251510	H74939222252010	H74939222253010	H74939222254010
2.75	H74939222271210	H74939222271510	H74939222272010	H74939222273010	-
3.00	H74939222301210	H74939222301510	H74939222302010	H74939222303010	H74939222304010
3.50	H74939222351210	H74939222351510	H74939222352010	H74939222353010	H74939222354010
4.00	H74939222401210	H74939222401510	H74939222402010	H74939222403010	H74939222404010