




# CASE STUDY: RENAL CRYOABLATION

## Cryoablation of 7cm Renal Tumor after Previous Contralateral Nephrectomy

Alex King, MD | University Hospital Southampton | Southampton, UK

### PRESENTATION


- 78-year-old male
- Right kidney multifocal renal cell carcinoma (RCC) with vein invasion and concurrent 7 cm left anterior renal tumor **A**
  - Left tumor shown on biopsy to be conventional clear cell RCC




Presentation: Multiphase CT scans reveal multifocal RCC and a 7 cm left anterior renal tumor in coronal (left), axial (center) and sagittal (right) views

### TREATMENT

- Right nephrectomy
- Left kidney cryoablation (three months after nephrectomy)
  - Hydrodissection to protect adjacent structures **B**
  - Seven IceRod™ 1.5 PLUS needles used to sculpt the ice to match the shape of the tumor **C D E**




Cryoablation: CT (axial view) shows hydrodissection (white arrows) and hypodense iceball formation (blue arrow)




Cryoablation: Needle placement and hypodense iceball formation in sagittal view (blue arrow)

# Cryoablation of 7 cm Renal Tumor after Previous Contralateral Nephrectomy

Alex King, MD | University Hospital Southampton | Southampton, UK



**Cryoablation:** Coronal CT demonstrates the distribution of needles in the slightly 'bean-shaped' tumor, illustrating how multi-needle cryoablation allows the ice to be sculpted to match the tumor morphology



**Cryoablation:** Visibility of leading edge of iceball on axial CT during treatment (arrows) allows confirmation of tumor coverage and appropriate parenchymal 'safety' margin

## OUTCOME


- Cryotherapy required a single night in hospital with simple analgesia
- No reported adverse events
- Early follow-up imaging of left kidney confirmed adequate ablation **F**

## CONCLUSION

- This case highlights the benefit of cryoablation to avoid bilateral nephrectomy and subsequent dialysis or the morbidity of open partial nephrectomy

**|| Careful hydrodissection and precise probe placement permitted a safe and adequate ablation without procedural complication. ||**

**|| In the setting of a previous contralateral nephrectomy, percutaneous cryoablation of this 7 cm tumor offered a nephron-sparing option without the morbidity of open surgery. ||**



**Post cryoablation:** Three-month follow-up CT scans. Top: conventional venous phase imaging. Bottom: iodine map from dual energy scan showing absence of residual tumoral enhancement

CT: Computed Tomography



### Cryoablation Needles Indications, Safety, and Warnings

<https://www.bostonscientific.com/cryoneedles-indications>



### ICEfx™ Cryoablation System Indications, Safety, and Warnings

<https://www.bostonscientific.com/icefx-indications>



### Visual ICE™ Cryoablation System Indications, Safety, and Warnings

<https://www.bostonscientific.com/visualicecryo-indications>

**Boston  
Scientific**  
Advancing science for life™

**Interventional Oncology  
and Embolization**  
300 Boston Scientific Way  
Marlborough, MA 01752-1234  
[www.bostonscientific.com](http://www.bostonscientific.com)

To order product or for more information  
contact customer service at 1.888.272.1001.

© 2026 Boston Scientific Corporation  
or its affiliates. All rights reserved.

PI-997708-AB