




CASE STUDY: LUNG CYROABLATION

Cryoablation of mCRC in the Lung Located Adjacent to the Aorta

Professor Thierry de Baere | Gustave Roussy Cancer Institute | Villejuif, France

PRESENTATION

- 57-year-old female with history of lung metastases from colorectal cancer
 - Video-assisted thoracic surgery 20 months previously for wedge resection of single 1.2 cm lung metastasis in right upper lobe
- New 1 cm lung metastasis identified in the left lower lobe **A**



Presentation: CT in the prone position under general anaesthesia shows a 1 cm lung nodule close to the aorta



TREATMENT PLAN

- Tumor board decided on thermal ablation
 - Cryoablation was selected from the thermal ablation techniques available due to the proximity of the nodule to the aorta, which rendered it difficult to target with an expandable needle
 - For successful cryoablation, the metastasis would need to be moved away from the aorta to avoid both the possible damage to the aorta from the ice and the heat-sink effect that would reduce the chances of complete ablation


TREATMENT

- One IceSphere™ 1.5 needle was placed in the target tumor under CT guidance **B**
- 50 ml of CO₂ was injected into the pleural space using a spring-loaded needle **C**
- The IceSphere™ 1.5 needle was activated for 30 seconds at 20% maximum power to stick the needle to the tumor ("stick-mode") to allow the needle and tumor to be gently pulled back by 3 cm **D**
- An additional 250 ml CO₂ was then injected and the needle pulled back a few centimeters further **E**
- With the tumor isolated in the middle of CO₂-filled pleural cavity, away from vulnerable neighboring structures, ablation could be completed safely
- Freeze and thaw cycles were completed per lung protocol
- The iceball was visible on CT as a 'ground glass' opacity covering the tumor with appropriate 'safety' margin beyond the tumor edge



Cryoablation of mCRC in the Lung Located Adjacent to the Aorta

Professor Thierry de Baere | Gustave Roussy Cancer Institute | Villejuif, France



Cryoablation: CT (axial view) showing single puncture for placement of IceSphere™ 1.5 needle in the target tumor



Cryoablation: CT (axial view) showing 50 ml of CO₂ injected into pleural space. Some CO₂ visible in the pleural space (arrow)




Cryoablation: After 30 seconds of activation at 20% maximum power, the IceSphere™ needle is gently pulled back by 3 cm




Cryoablation: Additional 250 ml of CO₂ is injected and the needle pulled back by a few more centimeters to isolate the tumor

OUTCOME

- One-month follow-up CT showed condensation at the location of the cryoablation **F**
- One-year follow-up CT showed near-complete regression of the ablation zone but unfortunately also confirmed progression of metastatic disease with appearance of multiple distant nodules **G**



Post cryoablation: One-month follow-up CT (axial view) shows condensation in the location of the cryoablation (arrow)



Post cryoablation: One-year follow-up CT shows nearly complete regression of the ablation zone (arrow) and multiple distant metastases.

CONCLUSION

- The unique “stick-mode” feature offered by cryoablation allows the tumor to be moved away from vulnerable adjacent organs, and consequently for technically challenging lung tumors to be treated safely and effectively



Cryoablation Needles Indications, Safety, and Warnings

<https://www.bostonscientific.com/cryoneedles-indications>



ICEfx™ Cryoablation System Indications, Safety, and Warnings

<https://www.bostonscientific.com/icefx-indications>



Visual ICE™ Cryoablation System Indications, Safety, and Warnings

<https://www.bostonscientific.com/visualicecryo-indications>

**Boston
Scientific**
Advancing science for life™

Peripheral Interventions
300 Boston Scientific Way
Marlborough, MA 01752-1234
www.bostonscientific.com

To order product or for more information
contact customer service at 1.888.272.1001.

© 2025 Boston Scientific Corporation
or its affiliates. All rights reserved.

PI-1439406-AB